CK-MB Calibrator

Order Information

Cat. No.	Package size
105-001132-00	3×1 mL

Intended Use

CK-MB Calibrator is used for calibration of quantitative determination of Creatine Kinase-MB (CK-MB) analyte on Mindray BS series chemistry analyzers.

Summary

CK-MB Calibrator contains Creatine Kinase-MB (CK-MB) analyte. By calibrating the selected analytes, Mindray BS-series chemistry analyzers can generate valid calibration factors, and then calculate the activities of each sample automatically after analysis. Mindray BS measurement system is composed of Mindray BS-series chemistry analyzers, Mindray reagent kits, calibrators and control materials.

Components

CK-MB Calibrator is a lyophilized calibrator based on serum. The calibrator values are lot-specific with the matched models listed in the value sheet.

Storage and stability

The calibrator is stable up to expiration date indicated on the label when stored in unopened vial at 2-8°C and protected from light.

Once reconstituted, please store calibrator tightly capped when not in use.

Once thawed to room temperature, the calibrator should not be re-frozen.

The stability information is as below.

Temperatures of storage	Stability
15-25°C	1 day
2-8°C	2 days
(-25)-(-15)°C	28 days

Preparation

1. It is recommended to take the vial out of the freezer and allow reaching room temperature.
2. Tap the vertically positioned vial gently and ensure that the lyophilized material is at the bottom of the vial.
3. Remove the screw cap and rubber stopper carefully, avoiding any loss of lyophilizate.
4. Reconstitute by adding 1.0 mL distilled/deionized water to the side of the vial slowly and exactly.

- Carefully replace the rubber stopper and leave to stand for 30 minutes at room temperature.
- During the standing time, mix contents by inverting several times and swirling gently to ensure that all components of lyophilizate are dissolved. Avoid the formation of foam.
- According to the calibration procedure described in the operation manual for the analyzers, use the specific value listed in the value sheet to set calibration parameter and execute calibration.
- Calibration is recommended after reagent lot change, execution of specific maintenance, out of control or other troubleshooting procedure.
- Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

Materials required but not provided

- General laboratory materials: NaCl solution 9 g/L (saline), distilled/deionized water.
- Mindray Creatine Kinase-MB Kit and Mindray ClinChem Multi Control. The catalogue number of reagent and control can be obtained from the corresponding Mindray reagent kit and control instruction for use.
- Mindray BS series chemistry analyzers, and General laboratory equipment.

Performance Characteristics

■ Calibration

- The calibrator values are assigned by Mindray standard transfer procedure and routine method. It is lot-specific with the matched models listed in the value sheet.
- Calibration is stable for approximately 30 days in BS-800 chemistry analyzers. The calibration stability may vary in different instruments, each laboratory should set a calibration frequency in the instrument parameters appropriate to their usage pattern.
Recalibration may be necessary when the following occur:
 - As changed reagent lot.
 - As required following quality control procedures or out of control.
 - As executes specific maintenance or troubleshooting procedure of chemistry analyzers.
- Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

■ Traceability of Calibrator

The calibrator values assigned by Mindray standard transfer procedure and routine method are listed in the value sheet. The traceability process is based on ISO 17511¹ and ISO 18153², the analyte (measurand) in this calibrator is traceable to Manufacturer's selected measurement procedure.

■ Homogeneity

Homogeneity	Units
Intra-vial homogeneity	SD≤3.00 U/L or CV≤2.00%
Inter-vial homogeneity	SD≤3.35 U/L or CV≤2.24%

Quality control

Use Mindray control as sample to monitor the status of calibration. The control result should be within the defined range as shown in the control value sheet.

Result interpretation

If the control result falls outside the range, the measurement system should be checked. For instance, position of reagent or sample placed on the analyzer; expiration date or storage condition of the calibrator, reagent and control; the parameter setting in software; or performance of the analyzer.

Warnings and precaution

1. For in vitro diagnostic use only. For laboratory professional use.
2. Please take the necessary precautions for handling all laboratory reagents.
3. Please confirm the integrity of the package before use. Do not use the calibrators with damaged packages. The calibrators avoid direct exposure to sunlight. The results can't be assured when stored at inappropriate condition.
4. If unintentionally opened before used, store the calibrators tightly capped at 2-8°C and protected from light, and the stability is equally to in-use stability.
5. Do not mix calibrator with different lots and bottles.
Do not use the calibrator beyond the expiration date and the in-use date.
Do not mix new opened calibrator with in-use calibrator.
Avoid the formation of foam.
6. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration or controls do not meet the insert and/or the Mindray System criteria.
7. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.

8. Do not swallow. Avoid contact with skin and mucous membranes.
9. When the calibrator accidentally enters the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
10. Safety data sheet is available for professional user on request.
11. Disposal of all waste material should be in accordance with local guidelines.
12. The calibrator was tested with CE-marked methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen³.
13. All human material should be considered potentially infectious.
14. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable.
15. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established

References

1. ISO 17511:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and control materials.
2. ISO 18153:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values for catalytic concentration of enzymes assigned to calibrators and control materials.
3. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington:2007.

Graphical symbols

In Vitro Diagnostic
medical device

Unique device
identifier

European
Conformity

Consult Instructions
For use

Use-by
date

Authorized representative in
the European Community

Batch Code

Temperature
limit

Manufacturer

Catalogue
number

Biological Risks

Keep away from sunlight

Indicates a medical device that needs protection from light sources

© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998; **Fax:** +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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Referências

1. ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
2. ISO 18153:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values for catalytic concentration of enzymes assigned to calibrators and control materials.
3. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington:2007.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Morada: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, R. P. China

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Bibliografia

1. ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
2. ISO 18153:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values for catalytic concentration of enzymes assigned to calibrators and control materials.
3. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office;

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Indirizzo: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, Repubblica popolare cinese

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Referanslar

1. ISO 17511:2003. In vitro diagnostic medical devices- Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and Kontrol materials.
2. ISO 18153:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values for catalytic concentration of enzymes assigned to calibrators and control materials.
3. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington:2007.

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Références

1. ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
2. ISO 18153:2003. Dispositifs médicaux de diagnostic in vitro—Mesurage

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3. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.

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© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)