

DECLARATION OF CONFORMITY

1) Manufacturer (Name, department): **Monobind Inc.**

Address: **100 North Pointe, LAKE FOREST, CA 92630. UNITED STATES**

and

2) European authorized representative: **CEpartner4U BV,**

Address: **ESDOORNLAAN 13, 3951DB MAARN, THE NETHERLANDS;**

(on product labels printed as:

CEpartner4U , ESDOORNLAAN 13, 3951DB MAARN, THE NETHERLANDS. www.cepartner4u.com)

3) Product(s) (name, type or model/batch number, etc.):

Immunoassay products;

AccuBind® ELISA,

AccuLite® CLIA,

QSure® Control,

Instruments

see appendix

4) The product(s) described above is in conformity with:

<u>Document No.</u>	<u>Title</u>
98/79/EC	<i>In vitro</i> Diagnostic Medical Devices Directive

5) Additional information (Conformity procedure, Notified Body, CE certificate, Registration nr., etc.):

Conformity assessment procedure for CE marking: *In vitro* Diagnostic Medical Device Directive, Annex III

Registration nr. : **NL- CA002-22758 and NL- CA002-22762**

Lake Forest, USA; 2021-09-20



(Place & date of issue (yyyy-mm-dd))

Tony Shatola; QA Director, Monobind Inc.
(name, function and signature of manufacturer)

Appendix

Date: 2021-09-20

List of devices.

<i>Device types</i>	<i>Item# AccuBind® ELISA Microwells</i>	<i>Item# AccuLite® CLIA Microwells</i>	<i>Item# QSure® Control</i>	<i>Item# Instru- ment</i>	<i>EDMS code</i>	<i>Risk Class</i>	<i>First date of CE-marking</i>
Allergy & Anemia							
Ferritin Test System	2825-300A 2825-300B	2875-300A 2875-300B			12.07.01.02.00	Low	2005-11-11
Folate Test System	7525-300A 7525-300B	7575-300A 7575-300B			12.07.01.03.00	Low	2010-06-29
Immunoglobulin E (IgE) Test System	2525-300A 2525-300B	2575-300A 2575-300B			12.02.01.02.00	Low	2005-11-11
Transferrin Soluble Receptor (sTfR) Test System	8625-300A 8625-300B	8675-300A 8675-300B			12.07.01.06.00	Low	2010-06-29
Vitamin B-12 (Vit B12) Test System	7625-300A 7625-300B	7675-300A 7675-300B			12.07.02.04.00	Low	2011-09-26
Folate, Vitamin B-12 (Anemia Panel VAST) Test System	7825-300A 7825-300B	7875-300A 7875-300B			12.07.01.00.00	Low	2013-09-16
Autoimmune							
Anti-Cyclic Citrullinated Peptide IgG (Anti-CCP IgG) Test System	12725-300A 12725-300B	12775-300A 12775-300B			12.11.01.90.00	Low	2019-04-03
Anti-Thyroglobulin (Anti-Tg) Test System	1025-300A 1025-300B	1075-300A 1075-300B			12.10.03.04.00	Low	2005-11-11
Anti-Thyroperoxidase (Anti-TPO) Test System	1125-300A 1125-300B	1175-300A 1175-300B			12.10.03.01.00	Low	2005-11-11
Bone Metabolism & Growth							
Calcitonin Test System	9325-300A 9325-300B	9375-300A 9375-300B			12.06.03.02.00	Low	2019-04-03
Growth Hormone (hGH) Test System	1725-300A 1725-300B	1775-300A 1775-300B			12.06.04.02.00	Low	2005-11-11
Parathyroid Hormone (PTH) Test System	9025-300A 9025-300B	9075-300A 9075-300B			12.06.03.13.00	Low	2011-09-26
Parathyroid Hormone (PTH) 3rd & 2nd Gen (VAST) Test System	10025-300A 10025-300B	10075-300A 10075-300B			12.06.03.13.00	Low	2019-04-03
25(OH) Vitamin D Total Direct (Vit D-Direct) Test System	7725-300A 7725-300B	7775-300A 7775-300B			12.06.03.10.00	Low	2017-07-05
Cancer Markers							
Alpha-Fetoprotein (AFP) Test System	1925-300A 1925-300B	1975-300A 1975-300B			12.03.90.01.00	Low	2005-11-11
CA-125 Test System	3025-300A 3025-300B	3075-300A 3075-300B			12.03.01.06.00	Low	2005-11-11
CA 15-3 Test System	5625-300A 5625-300B	5675-300A 5675-300B			12.03.01.02.00	Low	2010-06-29
CA 19-9 Test System	3925-300A 3925-300B	3975-300A 3975-300B			12.03.01.03.00	Low	2005-11-11
Carcinoembryonic Antigen (CEA) Test System	1825-300A 1825-300B	1875-300A 1875-300B			12.03.01.31.00	Low	2005-11-11
Next Generation Carcinoembryonic Antigen	4625-300A	4675-300A			12.03.01.31.00	Low	2010-06-29

<i>Device types</i>	<i>Item# AccuBind® ELISA Microwells</i>	<i>Item# AccuLite® CLIA Microwells</i>	<i>Item# QSure® Control</i>	<i>Item# Instru- ment</i>	<i>EDMS code</i>	<i>Risk Class</i>	<i>First date of CE-marking</i>
(CEA-Next Gen) Test System	4625-300B	4675-300B					
Free β-Subunit Human Chorionic Gonadotropin (Free Beta hCG) Test System	2025-300A 2025-300B	2075-300A 2075-300B			12.03.01.90.00	Low	2005-11-11
Cardiac Markers							
CK-MB Test System	2925-300A 2925-300B	2975-300A 2975-300B			12.13.01.02.00	Low	2005-11-11
Digoxin (DIG) Test System	925-300A 925-300B	975-300A 975-300B			12.08.01.01.00	Low	2005-11-11
High Sensitivity CRP (hs-CRP) Test System	3125-300A 3125-300B	3175-300A 3175-300B			12.13.01.90.00	Low	2005-11-11
Myoglobin Test System	3225-300A 3225-300B	3275-300A 3275-300B			12.13.01.05.00	Low	2005-11-11
Troponin I (cTnI) Test System	3825-300A 3825-300B	3875-300A 3875-300B			12.13.01.07.00	Low	2005-11-11
Diabetes							
C-Peptide Test System	2725-300A 2725-300B	2775-300A 2775-300B			12.06.01.01.00	Low	2005-11-11
Insulin Test System	2425-300A 2425-300B	2475-300A 2475-300B			12.06.01.03.00	Low	2005-11-11
Rapid Insulin Test System	5825-300A 5825-300B				12.06.01.03.00	Low	2010-06-29
Insulin - C-Peptide (Diabetes Panel VAST)	7325-300A 7325-300B	7375-300A 7375-300B			12.06.01.03.00	Low	2005-11-11
Endocrine							
ACTH Test System	10625-300	10675-300			12.06.04.01.00	Low	2019-04-03
Aldosterone Test System	10125-300	10175-300			12.06.02.01.00	Low	2019-04-03
Leptin Test System	10925-300	10975-300			12.06.90.17.00	Low	2019-04-03
Fertility & Prenatal							
Anti-Müllerian Hormone (AMH) Test System	9725-300A 9725-300B	9775-300A 9775-300B			12.05.02.16.00	Low	2019-04-03
Folicle Stimulating Hormone (FSH) Test System	425-300A 425-300B	475-300A 475-300B			12.05.01.04.00	Low	2005-11-11
B-Human Chorionic Gonadotropin (hCG) Test System	825-300A 825-300B	875-300A 875-300B			12.05.02.05.00	Low	2005-11-11
B-Human Chorionic Gonadotropin Extended Range (hCG-XR) Test System	8825-300A 8825-300B	8875-300A 8875-300B			12.05.02.05.00	Low	2013-09-16
Rapid B-Human Chorionic Gonadotropin (Rapid hCG) Test System	3325-300A 3325-300B				12.05.02.05.00	Low	2005-11-11
Inhibin A Test System	9525-300A 9525-300B	9575-300A 9575-300B			12.05.01.90.00	Low	2019-04-03
Inhibin B Test System	9625-300A 9625-300B	9675-300A 9675-300B			12.05.01.90.00	Low	2019-04-03
Luteinizing Hormone (LH) Test System	625-300A 625-300B	675-300A 675-300B			12.05.01.05.00	Low	2005-11-11
Pregnancy Associated Plasma Protein – A Mass Units (PAPP-A Mass Units) Test System	12625-300A 12625-300B	12675-300A 12675-300B			12.05.02.10.00	Low	2017-07-05
Prolactin Hormone (PRL) Test System	725-300A 725-300B	775-300A 775-300B			12.05.01.08.00	Low	2005-11-11

<i>Device types</i>	<i>Item# AccuBind® ELISA Microwells</i>	<i>Item# AccuLite® CLIA Microwells</i>	<i>Item# QSure® Control</i>	<i>Item# Instru- ment</i>	<i>EDMS code</i>	<i>Risk Class</i>	<i>First date of CE-marking</i>
Prolactin Hormone Sequential (PRLs) Test System	4425-300A 4425-300B	4475-300A 4475-300B			12.05.01.08.00	Low	2005-11-11
Human Chorionic Gonadotropin (hCG) , Human Prolactin (hPRL), Human Luteinizing Hormone (hLH), Follicle Stimulating Hormone (FSH) (Fertility Panel VAST) Test System	8325-300B 8325-300D 8325-300E	8375-300B 8375-300D 8375-300E			12.05.01.90.00	Low	2006-08-24
Alpha-Fetoprotein (AFP), Human Chorionic Gonadotropin (hCG), Unconjugated Estiol (u-E3) Triple Screen (Triple Screen Panel VAST) Test System	8525-300A 8525-300B	8575-300A 8575-300B			12.05.01.90.00	Low	2010-06-29
Infectious Diseases							
Anti-H. Pylori IgG (H. Pylori Ab IgG) Test System	1425-300A 1425-300B	1475-300A 1475-300B			15.01.04.03.00	Low	2005-11-11
Anti-H. Pylori IgM (H. Pylori Ab IgM) Test System	1525-300A 1525-300B	1575-300A 1575-300B			15.01.04.03.00	Low	2005-11-11
Anti-H. Pylori IgA (H. Pylori Ab IgA) Test System	1625-300A 1625-300B	1675-300A 1675-300B			15.01.04.03.00	Low	2005-11-11
Anti-SARS-CoV-2 (COVID-19) IgG Test System	11925-300A 11925-300B	11975-300A 11975-300B			15.04.80.90.00	Low	2020-08-25
Anti-SARS-CoV-2 (COVID-19) IgM Test System	11725-300A 11725-300B	11775-300A 11775-300B			15.04.80.90.00	Low	2020-08-25
Anti-SARS-CoV-2 (COVID-19) IgA Test System	11825-300A 11825-300B	11875-300A 11875-300B			15.04.80.90.00	Low	2020-08-25
Anti-SARS-CoV-2 (COVID-19) S1-RBD IgG Test System	12025-300A 12025-300B	12075-300A 12075-300B			15.04.80.90.00	Low	2021-09-20
D-Dimer Test System	9225-300A 9225-300B	9275-300A 9275-300B			13.02.05.03.00	Low	2020-08-25
Procalcitonin (PCT) Test System	1425-300A 1425-300B	1475-300A 1475-300B			12.06.90.16.00	Low	2017-07-05
Neonatal							
Neonatal 17OHP (N-17OHP) Test System	5525-300A 5525-300B				12.05.01.07.00	Low	2008-02-01
Neonatal (N-T4) Thyroxine Test System	2625-300A 2625-300B				12.04.01.12.00	Low	2005-11-11
Neonatal TBG (N-TBG) Test System	8925-300A 8925-300B				12.04.01.09.00	Low	2013-09-16
Neonatal TSH (N-TSH) Test System	3425-300A 3425-300B 3425-300D 3425-300E				12.04.01.90.00	Low	2005-11-11
Steroid							
Androstenedione (ANST) Test System	12425-300A 12425-300B	12475-300A 12475-300B			12.05.01.01.00	Low	2021-09-20
Cortisol Test System	3625-300A 3625-300B	3675-300A 3675-300B			12.06.02.04.00	Low	2005-11-11
Dehydroepiandrosterone (DHEA) Test System	7425-300A 7425-300B	7475-300A 7475-300B			12.05.01.02.00	Low	2011-09-26
Dehydroepiandrosterone Sulfate (DHEA-S) Test System	5125-300A 5125-300B	5175-300A 5175-300B			12.05.01.02.00	Low	2010-06-29
Estrone (E1) Test System	10325-300A 10325-300B	10375-300A 10375-300B			12.05.02.04.00	Low	2019-04-03

<i>Device types</i>	<i>Item# AccuBind® ELISA Microwells</i>	<i>Item# AccuLite® CLIA Microwells</i>	<i>Item# QSure® Control</i>	<i>Item# Instru- ment</i>	<i>EDMS code</i>	<i>Risk Class</i>	<i>First date of CE-marking</i>
Estradiol (E2) Test System	4925-300A 4925-300B	4975-300A 4975-300B			12.05.01.03.00	Low	2010-06-29
Unconjugated Estiol (u-E3) Test System	5025-300A 5025-300B	5075-300A 5075-300B			12.05.02.02.00	Low	2010-06-29
Progesterone Test System	4825-300A 4825-300B	4875-300A 4875-300B			12.05.01.06.00	Low	2010-06-29
17-OH Progesterone (17-OHP) Test System	5225-300A 5225-300B	5275-300A 5275-300B			12.05.01.07.00	Low	2010-06-29
17-OH Progesterone SI (17-OHP-SI) Test System	9925-300A 9925-300B	9975-300A 9975-300B			12.05.01.07.00	Low	2010-10-18
Sex Hormone Binding Globulin (SHBG) Test System	9125-300A 9125-300B	9175-300A 9175-300B			12.05.01.09.00	Low	2013-09-16
Testosterone Test System	3725-300A 3725-300B	3775-300A 3775-300B			12.05.01.10.00	Low	2007-11-01
Free Testosterone Test System	5325-300A 5325-300B	5375-300A 5375-300B			12.05.01.10.00	Low	2010-06-29
Thyroid							
Total Triiodothyronine (tT3) Test System	125-300A 125-300B 125-300D 125-300E	175-300A 175-300B 175-300D 175-300E			12.04.01.05.00	Low	2005-11-11
Free Triiodothyronine (fT3) Test System	1325-300A 1325-300B 1325-300A 1325-300B	1375-300A 1375-300B 1375-300D 1375-300E			12.04.01.01.00	Low	2005-11-11
Total Triiodothyronine (tT3 SBS) Test System	8125-300A 8125-300B	8175-300A 8175-300B			12.04.01.01.00	Low	2010-06-29
Rapid Total Triiodothyronine (Rapid -tT3) Test System	11225-300A 11225-300B				12.04.01.01.00	Low	2017-07-05
T3-Uptake (T3U) Test System	525-300A 525-300B	575-300A 575-300B			12.04.01.06.00	Low	2005-11-11
Thyroxine (tT4) Test System	225-300A 225-300B 225-300D 225-300E	275-300A 275-300B 275-300D 275-300E			12.04.01.07.00	Low	2005-11-11
Free Thyroxine (fT4) Test System	1225-300A 1225-300B 1225-300D 1225-300E	1275-300A 1275-300B 1275-300D 1275-300E			12.04.01.02.00	Low	2005-11-11
Total Thyroxine (tT4 SBS) Test System	8225-300A 8225-300B	8275-300A 8275-300B			12.04.01.01.00	Low	2010-06-29
Rapid Total Thyroxine (Rapid -tT4) Test System	11125-300A 11125-300B				12.04.01.01.00	Low	2017-07-05
Thyrotropin (TSH) Test System	325-300A 325-300B 325-300D 325-300E	375-300A 375-300B 375-300D 375-300E			12.04.01.11.00	Low	2005-11-11
Rapid TSH Test System	6025-300A 6025-300B	6075-300A 6075-300B			12.04.01.11.00	Low	2010-06-29
Thyroxine-Binding Globulin (TBG) Test System	3525-300A 3525-300B	3575-300A 3575-300B			12.04.01.09.00	Low	2005-11-11
Thyroglobulin (Tg) Test System	2225-300A	2275-300A			12.04.01.08.00	Low	2005-11-11

Device types	Item# AccuBind® ELISA Microwells	Item# AccuLite® CLIA Microwells	Item# QSure® Control	Item# Instru- ment	EDMS code	Risk Class	First date of CE-marking
	2225-300B	2275-300B					
Total Thyroxine (tT4), Total Triiodothyronine (tT3) & Thyroid Stimulating Hormone (TSH) (Thyroid Panel VAST) Test System	8025-300B 8025-300D 8025-300E	8075-300B 8075-300D 8075-300E			12.04.01.01.00	Low	2005-11-11
Free Thyroxine (fT4), Free Triiodothyronine (fT3) & Thyroid Stimulating Hormone (TSH) (Free Thyroid Panel VAST) Test System	7025-300B 7025-300D 7025-300E	7075-300B 7075-300D 7075-300E			12.04.01.01.00	Low	2010-06-29

Miscellaneous Controls							
Anti-H. Pylori Control (IgA, IgG, IgM) – Positive & Negative			HPC-300		12.50.01.16.00	Low	2013-09-16
Anti-Tg & Anti-TPO Control – Positive & Negative			AIT-101		12.50.01.16.00	Low	2010-06-29
Maternal Control – (AFP, uE3, hCG, Free beta hCG) Tri Level			MC-300		12.50.01.16.00	Low	2010-06-29
TBG Control – Tri-Level			TBG-300		12.50.01.16.00	Low	2013-09-16
Tg Control – Tri-Level			TG-300		12.50.01.16.00	Low	2010-06-29
Tumor Marker Control – (CA 125, CA 15-3, CA 19-9) Tri-Level			TMC-300		12.50.01.16.00	Low	2013-09-16

Miscellaneous Instruments							
Autoplex® ELISA & CLIA Analyzer				IN006	21.02.10.01	Low	2010-06-29
Autoplex® G2 ELISA & CLIA Analyzer				IN006-2	21.02.10.01	Low	2013-09-16
Autoplex® G3 ELISA & CLIA Analyzer				IN006-3	21.02.10.01	Low	2017-07-05
NeoEldex® ELISA Analyzer				IN009	21.02.10.01	Low	2011-09-26
Impulse® 3 CLIA Analyzer				IN007	21.02.10.01	Low	2010-06-29
NeoLumax® CLIA Analyzer				IN010	21.02.10.01	Low	2011-09-26
LuMatic® CLIA Analyzer				IN008	21.02.10.01	Low	2011-09-26
PrisMatic® ELISA Analyzer				IN013	21.02.10.01	Low	2013-09-16
PlateWash - Immunoassay Washer				IN002	21.02.10.01	Low	2010-06-29
TITIN® ELISA & CLIA Analyzer				IN015-EC	21.02.10.01	Low	2017-07-05
TITIN® ELISA Analyzer				IN015-E	21.02.10.01	Low	2017-07-05
TITIN-s® ELISA & CLIA Analyzer				IN016-EC	21.02.10.01	Low	2017-07-05
TITIN-s® ELISA Analyzer				IN016-E	21.02.10.01	Low	2017-07-05



Certificate of Registration of Quality Management System to ISO 13485:2016

Brazil - RDC ANVISA n. 16/2013 RDC ANVISA n. 23/2012 RDC ANVISA n. 67/2009

Canada - Medical Devices Regulations – Part 1- SOR 98/282

United States- 21 CFR 803, 21 CFR 806, 21 CFR 807 – Subparts A to D,

☒ 21 CFR 820 - Quality System Regulation,

The National Standards Authority of Ireland is an MDSAP Recognized Auditing Organization and certifies that:

Monobind Inc.

100 North Pointe Drive

Lake Forest, CA 92630

USA

Facility ID: F002818

has been assessed and deemed to comply with the requirements of the above standard and regulations in respect of the scope of operations given below:

The Design, Manufacture, and Distribution of In-Vitro Diagnostic Medical Device Immunoassays and Related Reagents and Controls. The Distribution of Related Washers and Analyzers.

Additional sites covered under this multi-site certification are listed on the Annex (File No. MP19.4585)

Approved by:
Kevin Mullaney
Director of Certification

Certificate Number: MP19.4585 / Rev 2
Certification Granted: 2019/09/25
Effective Date: 2022/09/25
Expiry Date: 2025/09/24



National Standards Authority of Ireland, 1 Swift Square, Northwood, Santry, Dublin 9, Ireland T +353 1 807 3800

National Standards Authority of Ireland, 20 Trafalgar Square, Nashua, New Hampshire, NH 03063, USA T +1 603 882 4412

All valid certifications are listed on NSAI's website – www.nsaiinc.com The continued validity of this certificate may be verified under "Approved Client Listing"



NSAI

Annex to Certificate Number: MP19.4585 / Rev 2

Scope of Registration:

The Design, Manufacture, and Distribution of In-Vitro Diagnostic Medical Device Immunoassays and Related Reagents and Controls. The Distribution of Related Washers and Analyzers.

Activity

Location

Headquarters, Design,
Manufacture

Monobind Inc.
100 North Pointe Drive
Lake Forest, CA 92630
USA
File No.: MP19.4585
Facility ID: F002818

Manufacture, Distribution

Monobind Inc.
103 North Pointe Drive
Lake Forest, CA 92630
USA
File No.: MP19.4585/A
Facility ID: F002818

Verified by:
Director of Certification



β-Human Chorionic Gonadotropin (hCG) Test System Product Code: 825-300

1.0 INTRODUCTION

Intended Use: The Quantitative Determination of Chorionic Gonadotropin (hCG) Concentration in Human Serum by a Microplate Enzyme Immunoassay, Colorimetric

2.0 SUMMARY AND EXPLANATION OF THE TEST

Human chorionic gonadotropin (hCG) concentration increases dramatically in blood and urine during normal pregnancy. hCG is secreted by placental tissue, beginning with the primitive trophoblast, almost from the time of implantation, and serves to support the corpus luteum during the early weeks of pregnancy. hCG or hCG similar glycoproteins can also be produced by a wide variety of trophoblastic and nontrophoblastic tumors. The measurement of hCG, by assay systems with suitable sensitivity and specificity has proven great value in the detection of pregnancy and the diagnosis of early pregnancy disorders.

According to the literature, hCG is detectable as early as 10 days after ovulation, reaching 100 mIU/ml by the first missed period. At the time for the next ovulation, the hCG level is 200 mIU/ml (approximately 28 days after conception).¹ A peak of 50,000 or even 100,000 mIU/ml is attained by the third month, then a gradual decline is observed.^{2,3}

In this method, hCG calibrator, patient specimen or control is first added to a streptavidin coated well. Biotinylated monoclonal and enzyme labeled antibodies (directed against distinct and different epitopes of hCG) are added and the reactants mixed. Reaction between the various hCG antibodies and native hCG forms a sandwich complex that binds with the streptavidin coated to the well.

After the completion of the required incubation period, the enzyme-chorionic gonadotropin antibody bound conjugate is separated from the unbound enzyme-chorionic gonadotropin conjugate by aspiration or decantation. The activity of the enzyme present on the surface of the well is quantitated by reaction with a suitable substrate to produce color.


The employment of several serum references of known chorionic gonadotropin levels permits construction of a dose response curve of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with chorionic gonadotropin concentration.

3.0 PRINCIPLE


Immunoassay (TYPE 3):

The essential reagents required for an immunoassay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitope recognition, **in excess**, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well

through the interaction of streptavidin coated on the well and exogenously added biotinylated monoclonal anti-hCG antibody. Upon mixing monoclonal biotinylated antibody, the enzyme-labeled antibody and a serum containing the native antigen, reaction results between the native antigen and the antibodies without competition or steric hindrance to form a soluble sandwich complex. The interaction is illustrated by the following equation:



Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:



After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

4.0 REAGENTS

Materials Provided:

A. hCG Calibrators – 1 ml/vial - Icons A-F

Six (6) vials of references for hCG Antigen at levels of 0(A), 5(B), 25(C), 50(D), 100(E) and 250(F) mIU/ml. Store at 2-8°C. A preservative has been added.

Note: The calibrators, human serum based, were calibrated using a reference preparation, which was assayed against the WHO 3rd IS (75/537).

B. hCG Enzyme Reagent – 13 ml/vial - Icon

One (1) vial containing enzyme labeled affinity purified antibody, biotinylated monoclonal mouse IgG in buffer, dye, and preservative. Store at 2-8°C.

C. Streptavidin Coated Plate – 96 wells – Icon

One 96-well microplate coated with streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

D. Wash Solution Concentrate – 20 ml/vial - Icon

One (1) vial containing a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.

E. Substrate A – 7 ml/vial - Icon S^A

One (1) vial containing tetramethylbenzidine (TMB) in buffer. Store at 2-8°C.

F. Substrate B – 7 ml/vial - Icon S^B

One (1) vial containing hydrogen peroxide (H₂O₂) in buffer. Store at 2-8°C.

G. Stop Solution – 8 ml/vial - Icon

One (1) vial containing a strong acid (1N HCl). Store at 2-8°C.

H. Product Instructions.

Note 1: Do not use reagents beyond the kit expiration date.

Note 2: Avoid extended exposure to heat and light. **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on the label.**

Note 3: Above reagents are for a single 96-well microplate

4.1 Required But Not Provided:

- Pipette(s) capable of delivering 0.025 and 0.050ml (25 & 50µl) volumes with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.100 and 0.350ml (100 & 350µl) volumes with a precision of better than 1.5%.
- Microplate washers or a squeeze bottle (optional).
- Microplate Reader with 450nm and 620nm wavelength absorbance capability.
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplate cover for incubation steps.
- Vacuum aspirator (optional) for wash steps.

- Timer.
- Quality control materials

5.0 PRECAUTIONS

For In Vitro Diagnostic Use
Not for Internal or External Use in Humans or Animals

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA licensed reagents. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-8395.

Safe Disposal of kit components must be according to local regulatory and statutory requirement.

6.0 SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood, serum in type and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain redtop venipuncture tube without additives or anti-coagulants. Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells.

In patients receiving therapy with high biotin doses (i.e. >5mg/day), no sample should be taken until at least 8 hours after the last biotin administration, preferably overnight to ensure fasting sample.

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.05 ml (50µl) of the specimen is required.

7.0 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, normal and elevated range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8.0 REAGENT PREPARATION

- Wash Buffer**
Dilute contents of wash concentrate to 1000ml with distilled or deionized water in a suitable storage container. Store diluted buffer at 2-30°C for up to 60 days.
- Working Substrate Solution** – Stable for one year
Pour the contents of the amber vial labeled Solution 'A' into the clear vial labeled Solution 'B'. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2 - 8°C.

Note1: Do not use the working substrate if it looks blue.

Note 2: Do not use reagents that are contaminated or have bacteria growth.

9.0 TEST PROCEDURE

Before proceeding with the assay, bring all reagents, serum reference calibrators and controls to room temperature (20-27°C).

****Test Procedure should be performed by a skilled individual or trained professional****

- Format the microplate wells for each serum reference calibrator, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C**

- Pipette 0.025 ml (25µl) of the appropriate serum reference calibrator, control or specimen into the assigned well.
- Add 0.100 ml (100µl) of hCG-Enzyme Reagent to all wells.
- Swirl the microplate gently for 20-30 seconds to mix and cover.
- Incubate 60 minutes at room temperature.
- Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
- Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100µl) of working substrate solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells**
- DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION**
- Incubate at room temperature for fifteen (15) minutes.
- Add 0.050ml (50µl) of stop solution to each well and gently mix for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells**
- Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. **The results should be read within thirty (30) minutes of adding the stop solution.**

10.0 CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of Human chorionic gonadotropin (hCG) in unknown specimens.


- Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
- Plot the absorbance for each duplicate serum reference versus the corresponding hCG concentration in mIU/ml on linear graph paper (do not average the duplicates of the serum references before plotting).
- Draw the best-fit curve through the plotted points.
- To determine the concentration of hCG for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in mIU/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (1.745) intersects the dose response curve at (157 mIU/ml) hCG concentration (See Figure 1).

Note: Computer data reduction software designed for ELISA assays may also be used for the data reduction. **If such software is utilized, the validation of the software should be ascertained.**

EXAMPLE 1

Sample I.D.	Well Number	Abs (A)	Mean Abs (B)	Value (mIU/ml)
Cal A	A1	0.002	0.004	0
	B1	0.005		
Cal B	C1	0.073	0.071	5
	D1	0.069		
Cal C	E1	0.340	0.350	25
	F1	0.360		
Cal D	G1	0.637	0.650	50
	H1	0.663		
Cal E	A2	1.223	1.212	100
	B2	1.199		
Cal F	C2	2.518	2.502	250
	D2	2.486		
Ctrl 1	E2	0.075	0.076	5.8
	F2	0.077		
Ctrl 2	G2	0.280	0.290	21.9
	H2	0.301		
Patient	A3	1.736	1.745	157
	B3	1.754		

*The data presented in Example 1 and Figure 1 are for illustration only and **should not** be used in lieu of a dose response curve prepared with each assay.



11.0 Q.C. PARAMETERS

In order for the assay results to be considered valid the following criteria should be met:

- The absorbance (OD) of calibrator 'F' should be ≥ 1.3 .
- Four out of six quality control pools should be within the established ranges.

12.0 RISK ANALYSIS

The MSDS and Risk Analysis Form for this product are available on request from Monobind Inc.

12.1 Assay Performance

- It is important that the time of reaction in each well is held constant to achieve reproducible results.
- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Patient specimens with hCG concentrations above 250 mIU/ml may be diluted with normal male serum (hCG < 1 mIU/ml) and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind IFU may yield inaccurate results.
- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
- Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from Monobind@monobind.com.

12.2 Interpretation

- Measurements and interpretation of results must be performed by a skilled individual or trained professional.**
- Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- The reagents for the test system have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays (Boscato LM, Stuart MC. 'Heterophilic antibodies: a problem for all immunoassays' Clin. Chem.

- 1988:3427-33). For diagnostic purposes, the results from this assay should be in combination with clinical examination, patient history and all other clinical findings. For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability**.
 - If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.
 - False positive results may occur in the presence of a wide variety of trophoblastic and nontrophoblastic tumors that secrete hCG. Therefore, the possibility of an hCG secreting neoplasia should be eliminated prior to diagnosing pregnancy.
 - Also, false positive results may be seen when assaying specimens from individuals taking the drugs Pergonal* and Clomid**. Additionally Pergonal will often be followed with an injection of hCG.
 - Spontaneous microabortions and ectopic pregnancies will tend to have values which are lower than expected during a normal pregnancy while somewhat higher values are often seen in multiple pregnancies.^{5,6,7}
 - Following therapeutic abortion, detectable hCG may persist for as long as three to four weeks. The disappearance rate of hCG, after spontaneous abortion, will vary depending upon the quantity of viable residual trophoblast.^{4,5,6,7}
 - A hCG value alone is not of diagnostic value** and should only be used in conjunction with other clinical manifestations (observations) and diagnostic procedures.

*Pegonal is a registered trademark of Serono Laboratories, Inc.
 **Clomid is a registered trademark of Merriell-National Laboratories

13.0 EXPECTED RANGES OF VALUES

A study of an apparent normal adult population was undertaken to determine expected values for the HCG AccuBind® ELISA Test System. The mean (X) values, standard deviations (σ) and expected ranges ($\pm 2\sigma$) are presented in Table 1.

TABLE 1
Expected Values for the hCG ELISA Test System
(In mIU/ml - 3rd IS 75/537)

Number	25
Mean	2.9
Standard Deviation	1.4
Expected Ranges ($\pm 2\sigma$)	0.1 - 5.7

Expected levels for hCG during normal pregnancy (3) are listed in Table 2.

TABLE 2
Expected Values for hCG levels (3rd IS 75/537)
during normal pregnancy (in mIU/ml)

	10 - 30
1 st week	30 - 100
2 nd week	100 - 1000
3 rd week	1,000 - 10,000
4 th week	30,000 - 100,000
2 nd & 3 rd month	10,000 - 30,000
2 nd trimester	5,000 - 15,000
3 rd trimester	

Values for hCG for a normal, healthy population and pregnant women, during gestation cycle, are given in Table 3. The values depicted below represent limited in house studies in concordance with published literature.^{8,9,10}

TABLE 3
Median Values during Gestation.

Gestation (Week)	hCG (IU/ml)
15	40.88
16	33.87
17	28.71
18	26.74
19	18.76
20	19.24
21	23.46

It is important to keep in mind that establishment of a range of values which can be expected to be found by a given method for a population of "normal"-persons is dependent upon a multiplicity of factors: the specificity of the method, the population tested and

the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

14.0 PERFORMANCE CHARACTERISTICS

14.1 Precision

The within and between assay precisions of the hCG AccuBind® ELISA were determined by analyses on three different levels of control sera. The number (N), mean value (X), standard deviation (σ) and coefficient of variation (C.V.) for each of these control sera are presented in Table 4 and Table 5.

TABLE 4
Within Assay Precision (Values in mIU/ml)

Sample	N	X	σ	C.V.
Level 1	20	4.4	0.22	4.9%
Level 2	20	18.7	0.75	4.0%
Level 3	20	214.8	14.59	6.8%

TABLE 5
Between Assay Precision* (Values in mIU/ml)

Sample	N	X	σ	C.V.
Level 1	20	5.4	0.52	9.6%
Level 2	20	22.4	1.97	8.8%
Level 3	20	213.1	15.16	7.1%

*As measured in ten experiments in duplicate.

14.2 Sensitivity

The hCG AccuBind® ELISA test system has a sensitivity of 0.003 mIU/well. This is equivalent to a sample containing 0.102 mIU/ml hCG concentration. The analytical sensitivity (detection limit) was ascertained by determining the variability of the '0 mIU/ml' calibrator and using the 2σ (95% certainty) statistic to calculate the minimum dose.

14.3 Accuracy

This hCG AccuBind® ELISA test system was compared with a reference radioimmunoassay. Biological specimens from normal and pregnant populations were assayed. The total number of such specimens was 110. The least square regression equation and the correlation coefficient were computed for the hCG ELISA in comparison with the reference method. The data obtained is displayed below.

TABLE 6

Method	Mean (x)	Least Square Regression Analysis	Correlation Coefficient
Monobind	14.8	$y = 0.081 + 0.93(x)$	0.989
Reference	15.1		

Only slight amounts of bias between the hCG ELISA method and the reference method are indicated by the closeness of the mean values. The least square regression equation and correlation coefficient indicates excellent method agreement.

14.4 Specificity

The cross-reactivity of the hCG AccuBind® ELISA to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of chorionic gonadotropin needed to produce the same absorbance.

Substance	Cross Reactivity	Concentration
Chorionic Gonadotropin (hCG)	1.0000	----
β -hCG subunit	< 0.0001	1000ng/ml
Follitropin (FSH)	< 0.0001	1000ng/ml
Lutropin Hormone (LH)	< 0.0001	1000ng/ml
hrotropin (TSH)	< 0.0001	1000ng/ml

14.5 Hook Effect

The test shows no hook effect up to concentrations of > 150,000 mIU/ml.

15.0 REFERENCES

- Kosasa TS, "Measurement of Human Chorionic Gonadotropin", *Journal of Reproductive Medicine*, 26, 201-6 (1981).
- Danzer H, Braunstein GD, et al, "Maternal Serum Human Chorionic Gonadotropin Concentrations and Fetal Sex Predictions", *Fertility and Sterility*, 34, 336-40 (1980).
- Braunstein G.D., et al., "Serum Human Chorionic Gonadotropin Levels through Normal Pregnancy", *American Journal of Obstetrics and Gynecology* 126, pg. 678-81 (1976).
- Goldstein DP, and Kosasa TS, "The Subunit Radioimmunoassay for HCG Clinical Application", *Gynecology*, 6, 145-84 (1975).
- Batzer F, "Hormonal Evaluation of Early Pregnancy", *Fertility and Sterility*, 34, 1-12 (1980).
- Braunstein GD, et al, "First-Trimester Chorionic Gonadotropin Measurements as an Aid to the Diagnosis of Early Pregnancy Disorders", *American Journal of Obstetrics and Gynecology*, 131, 25-32 (1978).
- Lenton E, Neal L and Sulaiman R, "Plasma Concentrations of Human Gonadotropin from the time of Implantation until the Second Week of Pregnancy", *Fertility and Sterility*, 37, 773-78 (1982).
- Canick JA, Rish S. 'The accuracy of assigned risks in maternal serum screening', *Prenatal Diagnosis*; 18:413-415 (1998).
- NIH State-of-the-Science Conference Statement on Management of Menopause-Related Symptoms. NIH Consensus State Sci Statements. Mar 21-23; 22(1), 1-38 (2005).
- Tietz NW, Ed: *Clinical Guide to Laboratory Tests* 3rd Ed, Philadelphia, WA Saunders Co (1995).

Revision: 5 Date: 2021-Sep-23 DCO: 1509
 MP825 Product Code: 825-300

Size	96(A)	192(B)
Reagent (fill)	A) 1ml set	1ml set
	B) 1 (13ml)	2 (13ml)
	C) 1 plate	2 plates
	D) 1 (20ml)	1 (20ml)
	E) 1 (7ml)	2 (7ml)
	F) 1 (7ml)	2 (7ml)
	G) 1 (8ml)	2 (8ml)

For Orders and Inquiries, please contact


Monobind Inc.
 100 North Pointe Drive
 Lake Forest, CA 92530 USA

Tel: +1 949.951.2665 Mail: info@monobind.com
 Fax: +1 949.951.3539 Fax: www.monobind.com



Please visit our website to learn more about our products and services.

Glossary of Symbols
 (EN 960/ISO 15223)





Dehydroepiandrosterone Sulfate (DHEA-S) Test System Product Code: 5125-300

1.0 INTRODUCTION

Intended Use: The Quantitative Determination of Dehydroepiandrosterone Sulfate Concentration in Human Serum or Plasma by a Microplate Enzyme Immunoassay, Colorimetric

2.0 SUMMARY AND EXPLANATION OF THE TEST

Dehydroepiandrosterone sulfate (DHEA-S) is the major C19 steroid secreted by the adrenal cortex, and is a precursor in testosterone and estrogen biosynthesis. DHEA-S, the sulfate ester of DHEA, is derived from sulfated precursors and by enzymatic conversion of DHEA in adrenal and extradrenal tissues. Due to the presence of a 17-oxo [rather than hydroxyl] group, DHEA-S possesses relatively weak androgenic activity, which for unsulfated DHEA has been estimated at ~10% that of testosterone.¹ However, the bioactivity of DHEA-S may be increased by its relatively high serum concentrations, approximately 100 to 1000-fold higher than DHEA or testosterone, and its weak affinity for sex-hormone binding globulin.²

The physiologic role of DHEA-S is not well-defined. Serum levels are relatively high in the fetus and neonate, low during childhood, and increase during puberty.^{3,4} Increased levels of DHEA-S during adrenarche may contribute to the development of secondary sexual hair. DHEA-S levels show a progressive decline after the third decade of life.⁵ Unlike DHEA, DHEA-S levels do not show significant diurnal variation and little day-to-day variation. DHEA-S levels are not responsive to acute corticotropin administration,⁴ and do not vary significantly during the normal menstrual cycle.² This may be due to the slower metabolic clearance rate of DHEA-S as compared to DHEA.⁶

Measurement of serum DHEA-S is a useful marker of adrenal androgen synthesis. Abnormally low levels have been reported in hypoadrenalism,³ while elevated levels occur in several conditions; including virilizing adrenal adenoma and carcinoma,⁷ 21-hydroxylase and 3β-hydroxysteroid dehydrogenase deficiencies⁸ and some cases of female hirsutism.² Since very little DHEA-S is produced by the gonads,^{8,9} measurement of DHEA-S may aid in the localization of the androgen source in virilizing conditions. Methods for measurement of DHEA-S include gas-liquid chromatography, double-isotope derivative techniques, competitive protein-binding assays, and radioimmunoassay. Although significant cross-reactivity occurs with DHEA, androstenedione and androsterone, the relative concentrations of these competing substances in most normal and pathologic samples predicts a minimal effect on assay performance.


The Monobind DHEA-S ELISA Kit uses a specific anti-DHEA-S antibody, and does not require prior sample extraction of serum or plasma. Cross-reactivity to other naturally occurring and structurally related steroids is low. The employment of several serum references of known DHEA-S concentration permits

construction of a graph of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with DHEA-S concentration.

3.0 PRINCIPLE

Competitive Enzyme Immunoassay (TYPE 7):

The essential reagents required for an enzyme immunoassay include antibody, enzyme-antigen conjugate and native antigen. Upon mixing biotinylated antibody, enzyme-antigen conjugate and a serum containing the native antigen, a competition reaction results between the native antigen and the enzyme-antigen conjugate for a limited number of antibody binding sites. The interaction is illustrated by the following equation:



Ab_{Bt} = Biotinylated x-DHEA-S IgG Antibody (Constant Quantity)

Ag = Native Antigen (Variable Quantity)

$\text{Enz} \text{ Ag}$ = Enzyme-antigen Conjugate (Constant Quantity)

AgAb_{Bt} = Antigen-Antibody Complex

$\text{Enz} \text{ AgAb}_{\text{Bt}}$ = Enzyme-antigen Conjugate -Antibody Complex

k_a = Rate Constant of Association

k_{-a} = Rate Constant of Disassociation

$K = k_a / k_{-a}$ = Equilibrium Constant

A simultaneous reaction between the biotin attached to the antibody and the streptavidin immobilized on the microwell occurs. This effects the separation of the antibody bound fraction after decantation or aspiration.

$\text{AgAb}_{\text{Bt}} + \text{Enz} \text{ AgAb}_{\text{Bt}} + \text{Streptavidin}_{\text{CW}} \rightarrow \text{immobilized complex}$

$\text{Streptavidin}_{\text{CW}}$ = Streptavidin immobilized on well

Immobilized complex = sandwich complex bound to the solid surface

The enzyme activity in the antibody bound fraction is inversely proportional to the native antigen concentration. By utilizing several different serum references of known antigen concentration, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

4.0 REAGENTS

Materials Provided:

A. DHEA-S Calibrators – 1ml/vial - Icons A-F

Six (6) vials of serum reference for DHEA-S at concentrations of 0 (A), 0.2 (B), 1.0 (C), 2.0 (D), 4.0 (E) and 8.0 (F) in µg/ml. Store at 2-8°C. A preservative has been added. The calibrators can be expressed in molar concentrations (nM/L) by using 2.71 as a conversion factor.

For example: 1 µg/ml x 2.71 = 2.71 µM/L

B. DHEA-S Enzyme Reagent – 6.0 ml/vial ☒

One (1) vial of DHEA-S (Analog)-horseradish peroxidase (HRP) conjugate in a protein-stabilizing matrix with red dye. Store at 2-8°C.

C. DHEA-S Biotin Reagent – 6.0 ml - Icon ▽

One (1) bottle of reagent contains anti-DHEA-S biotinylated purified rabbit IgG conjugate in buffer, blue dye and preservative. Store at 2-8°C.

D. Streptavidin Coated Plate – 96 wells -Icon ▴

One 96-well microplate coated with 1.0 µg/ml streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

E. Wash Solution Concentrate – 20ml/vial - Icon ♀

One (1) vial contains a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.

F. Substrate A – 7ml/vial - Icon S^A

One (1) vial contains tetramethylbenzidine (TMB) in buffer. Store at 2-8°C.

G. Substrate B – 7ml/vial - Icon S^B

One (1) vial contains hydrogen peroxide (H₂O₂) in buffer. Store at 2-8°C.

H. Stop Solution – 8ml/vial - Icon ⏹

One (1) vial contains a strong acid (1N HCl). Store at 2-8°C.

I. Product Instructions

Note 1: Do not use reagents beyond the kit expiration date.

Note 2: Avoid extended exposure to heat and light. **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on the label.**

Note 3: Above reagents are for a single 96-well microplate.

4.1 Required But Not Provided:

1. Pipette capable of delivering 0.010ml (10µl) and 0.050ml (50µl) with a precision of better than 1.5%.
2. Dispenser(s) for repetitive deliveries of 0.100ml (100µl) and 0.350ml (350µl) volumes with a precision of better than 1.5%.
3. Adjustable volume (200-1000µl) dispenser(s) for conjugate.
4. Microplate washer or a squeeze bottle (optional).
5. Microplate Reader with 450nm and 620nm wavelength absorbance capability.
6. Absorbent Paper for blotting the microplate wells.
7. Plastic wrap or microplate cover for incubation steps.
8. Vacuum aspirator (optional) for wash steps.
9. Timer.
10. Quality control materials.

5.0 PRECAUTIONS

**For In Vitro Diagnostic Use
Not for Internal or External Use in Humans or Animals**

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA required tests. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-8395.

Safe Disposal of kit components must be according to local regulatory and statutory requirement.

6.0 SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood serum or heparanized plasma in type, and taken with the usual precautions in the collection of venipuncture samples. For accurate comparison to establish normal values, a fasting morning serum sample should be obtained. The blood should be collected in a redtop veni-puncture tube with or without additives or anti-coagulants (for serum) or evacuated tube(s) containing EDTA or heparin (for plasma). Allow the blood to clot for serum samples. Centrifuge the specimen to separate the serum or plasma from the cells.

In patients receiving therapy with high biotin doses (i.e. >5mg/day), no sample should be taken until at least 8 hours after the last biotin administration, preferably overnight to ensure fasting sample.

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.020ml (20µl) of the specimen is required.

7.0 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, normal and high range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8.0 REAGENT PREPARATION

1. Wash Buffer

Dilute contents of wash solution to 1000ml with distilled or deionized water in a suitable storage container. Diluted buffer can be stored at 2-30°C for up to 60 days.

2. Working Substrate Solution - Stable for 1 year

Pour the contents of the amber vial labeled Solution 'A' into the clear vial labeled Solution 'B'. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2 - 8°C.

Note1: Do not use the working substrate if it looks blue.

Note 2: Do not use reagents that are contaminated or have bacteria growth.

9.0 TEST PROCEDURE

Before proceeding with the assay, bring all reagents, serum reference calibrators and controls to room temperature (20-27°C).

*****Test Procedure should be performed by a skilled individual or trained professional*****

1. Format the microplates' wells for each serum reference calibrator, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
2. Pipette 0.010 ml (10 µl) of the appropriate serum reference calibrator, control or specimen into the assigned well.
3. Add 0.050 ml (50µl) of the DHEA-S Enzyme Reagent to all wells.
4. Swirl the microplate gently for 20-30 seconds to mix.
5. Add 0.050 ml (50µl) of Anti- DHEA-S Biotin Reagent to all wells.
6. Swirl the microplate gently for 20-30 seconds to mix.
7. Cover and incubate for 30 minutes at room temperature.
8. Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
9. Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
10. Add 0.100 ml (100µl) of working substrate solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells.**
11. **DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION**
12. Incubate at room temperature for fifteen (15) minutes.
13. Add 0.050ml (50µl) of stop solution to each well and gently mix for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
14. Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm). **The results should be read within thirty (30) minutes of adding the stop solution.**

Note: Dilute the samples suspected of concentrations higher than 8.0 ug/ml 1:5 and 1:10 with DHEA-S 0' µg/ml calibrator or patient serum pools with a known low value for DHEA-S.

10.0 CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of DHEA-S in unknown specimens.


1. Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
2. Plot the absorbance for each duplicate serum reference versus the corresponding DHEA-S concentration in ug/ml on linear graph paper (do not average the duplicates of the serum references before plotting).
3. Connect the points with a best-fit curve.
4. To determine the concentration of DHEA-S for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ug/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance in the patient sample (1.078) intersects the dose response curve at (1.21 µg/ml) DHEA-S concentration (See Figure 1).

Note: Computer data reduction software designed for ELISA assays may also be used for the data reduction. **If such**

software is utilized, the validation of the software should be ascertained.

EXAMPLE 1

Sample I.D.	Well Number	Abs (A)	Mean Abs (B)	Value (µg/ml)
Cal A	A1	2.562	2.572	0.0
	B1	2.582		
Cal B	C1	1.865	1.847	0.2
	D1	1.829		
Cal C	E1	1.186	1.163	1.0
	F1	1.140		
Cal D	G1	0.855	0.850	2.0
	H1	0.845		
Cal E	A2	0.555	0.556	4.0
	B2	0.557		
Cal F	C2	0.355	0.349	8.0
	D2	0.344		
Cont 1	G2	1.394	1.387	0.62
	H2	1.380		
Pat# 1	A3	1.065	1.078	1.21
	B3	1.091		



*The represented in Example 1 and Figure 1 is for illustration only and **should NOT** be used in lieu of a dose response curve prepared with each assay.

11.0 Q.C. PARAMETERS

In order for the assay results to be considered valid the following criteria should be met:

- The absorbance (OD) of calibrator 0 ug/ml should be ≥ 1.3
- Four out of six quality control pools should be within the established ranges.

12.0 RISK ANALYSIS

The MSDS and Risk Analysis Form for this product are available on request from Monobind Inc.

12.1 Assay Performance

- It is important that the time of reaction in each well is held constant to achieve reproducible results.
- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Patient specimens with DHEA-S concentrations above 8.0 µg/mL may be diluted (1/5, 1/10 or higher) with DHEA-S '0' calibrator and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential.

Any deviation from Monobind' IFU may yield inaccurate results.

- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
- Risk Analysis: as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from Monobind@monobind.com.

12.2 Interpretation

1. Measurements and interpretation of results must be performed by a skilled individual or trained professional.

- Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- The reagents for the test system have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays (Boscato LM, Stuart MC. 'Heterophilic antibodies: a problem for all immunoassays' Clin. Chem. 1988;34:27-33). For diagnostic purposes, the results from this assay should be in combination with clinical examination, patient history and all other clinical findings.
- For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability**.
- If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.
- Clinically, a **DHEA-S value alone is not of diagnostic value** and should only be used in conjunction with other clinical manifestations (observations) and diagnostic procedures.

13.0 EXPECTED RANGES OF VALUES

In agreement with established reference intervals for a "normal" adult population, the expected ranges for the DHEA-S AccuBind® ELISA Test System are detailed in Table 1.

TABLE 1 Expected Values for the DHEA-S Test System	
POPULATION	RANGE (µg/ml)
Male	0.06 – 4.58
Female	0.03 – 5.88

It is important to keep in mind that establishment of a range of values which can be expected to be found by a given method for a population of "normal" persons is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

14.0 PERFORMANCE CHARACTERISTICS

14.1 Precision

The within and between assay precision of the DHEA-S AccuBind® ELISA Test System were determined by analyses on three different levels of pool control sera. The number, mean values, standard deviation and coefficient of variation for each of these control sera are presented in Table 2 and Table 3.

TABLE 2 Within Assay Precision (Values in µg/ml)				
Sample	N	X	σ	C.V.
Low	16	0.66	0.06	9.8%
Normal	16	1.14	0.05	4.9%
High	16	4.84	0.21	4.3%

TABLE 3

Between Assay Precision (Values in µg/ml)				
Sample	N	X	σ	C.V.
Low	10	0.61	0.06	9.5%
Normal	10	1.36	0.04	3.1%
High	10	4.73	0.16	3.4%

*As measured in tens experiments in duplicate over a ten day period.

14.2 Sensitivity

The DHEA-S AccuBind® ELISA Test System has a sensitivity of 0.042 ug/ml. The sensitivity was ascertained by determining the variability of the 0 ug/ml serum calibrator and using the 2σ (95% certainty) statistic to calculate the minimum dose.

14.3 Accuracy

The DHEA-S AccuBind® ELISA Test System was compared with a chemiluminescence immunoassay method. Biological specimens from low, normal and relatively high DHEA-S level populations were used (The values ranged from 0.2 ug/ml – 7.7 ug/ml). The total number of such specimens was 77. The least square regression equation and the correlation coefficient were computed for this DHEA-S EIA in comparison with the reference method. The data obtained is displayed in Table 4.

TABLE 4			
Method	Mean (x)	Least Square Regression Analysis	Correlation Coefficient
Monobind (y)	1.12	y = 0.1448 + 0.986x	0.983
Reference (X)	1.18		

Only slight amounts of bias between this method and the reference method are indicated by the closeness of the mean values. The least square regression equation and correlation coefficient indicates excellent method agreement.

14.4 Specificity

The % cross reactivity of the DHEA-S antibody to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of DHEA-S needed to displace the same amount of labeled analog.

Substance	Cross Reactivity
DHEA-S	1.0000
DHEA	0.0004
Androstenedione	0.0003
Dihydrotestosterone	0.0008
Cortisone	<0.0001
Corticosterone	<0.0001
Cortisol	0.0004
Spirolactone	<0.0001
Estrilol	<0.0001
Estradiol	<0.0001
Estrone	<0.0001
Testosterone	<0.0001

15.0 REFERENCES

- Dorfman RI, Shipley RA: *Androgens*, John Wiley and Sons, New York, 1956, pp. 116-128.
- Pang S, Riddick L: Hirsutism. IN Lifshitz F (Ed): *Pediatric Endocrinology*, A Clinical Guide, second edition. Marcel Dekker, Inc., New York, 1990, pp. 259-291.
- De Peretti E, Forest MG: Pattern of plasma dehydroepiandrosterone sulfate levels in humans from birth to adulthood: evidence for testicular production. *J Clin Endocrinol Metab*, 47:572-577 (1978).
- Lashansky G, Saenger P, Fishman K, Gautier T, Mayes D, Berg G, Di Martino-Nardi J, Reiter E, "Normative data for adrenal steroidogenesis in a healthy pediatric population: age and sex-related changes after adrenocorticotropin stimulation", *J Clin Endocrinol Metab*, 73:674-686 (1991).
- Zumoff B, Roenfeldt RS, Stain GW, Levin J, Fukushima DK, "Sex differences in twenty-four hour mean plasma concentration of dehydroepiandrosterone (DHEA) and dehydroepiandrosterone sulfate (DHEAS) and the DHEA to DHEAS ratio in normal adults", *J Clin Endocrinol Metab* 51:330-333 (1980).

- Pang S, Lerner A, Stoner E, Oberfield S, Engle I, New M, "Late-onset adrenal steroid 3'-hydroxysteroid dehydrogenase deficiency: A cause of hirsutism in pubertal and postpubertal women", *J Clin Endocrinol Metab*, 60:428-439 (1985).
- Lee PDK, Winter RJ, and Green OC, "Virilizing adrenocortical tumors in childhood: eight cases and a review of the literature", *Pediatrics*, 76:437-444 (1985).
- Rittmaster RS, "Differential suppression of testosterone and estradiol in hirsute women with the superactive gonadotropin-releasing hormone agonist leuprolide" *J Clin Endocrinol Metab*, 67:651-655,(1988).
- Beer NA, Jakubowicz DJ, Beer RM, Arocha IR, and Nestler JE "Effects of nifedipine on glucose tolerance and serum insulin dehydroepiandrosterone sulfate levels in insulin-resistant obese and hypertensive men", *J Clin Endocrinol Metab*, 76:178-183 (1993).
- Tietz NW, *Textbook of clinical chemistry*, 2nd ed. Philadelphia: W.B. Saunders, 1994.

Revision: 5 Date: 2019-Jul-16 DCO: 1353
MP5125 Product Code: 5125-300

Size		96(A)	192(B)
Reagent (fill)	A)	1ml set	1ml set
	B)	1 (6ml)	2 (6ml)
	C)	1 (6ml)	2 (6ml)
	D)	1 plate	2 plates
	E)	1 (20ml)	1 (20ml)
	F)	1 (7ml)	2 (7ml)
	G)	1 (7ml)	2 (7ml)
	H)	1 (8ml)	2 (8ml)

For Orders and Inquires, please contact

Monobind Inc.
100 North Pointe Drive
Lake Forest, CA 92630 USA

Tel: +1 949.951.2665 Mail: info@monobind.com
Fax: +1 949.951.3539 Fax: www.monobind.com



Please visit our website to learn more about our products and services.

Glossary of Symbols (EN 980/ISO 15223)

