

FT13 CALIBRATION DETERGENT CALIBRATION

- Connect the 78401-444 hookup (self-disinfection hookup) into the basin
- Place the detergent container on the electronic scale, reset it and start the SET DETERGENT cycle (available on the technical menu only).

INSTRUMENT				CYCLE TYPE	
CATEGORY	MODEL	SERIAL NUMBER	ID	CYCLE NAME	
Sigmoidoscope	ES-3840	1	1	double clean	
Gastroscope	QIF-Q165	A012345	1	TEST	
Colonoscope	CF-100L	123		COMPLETE DISINFECTION	
Duodenoscope	JF-130	A123		FAST DISINFECTION	
TEST				COMPLETE STERILIZATION	
				FAST STERILIZATION	
				CLEANING SAMPLE	
				CALIBRATION	
				cal	
				CALIBRATION WATER	
				CALIBRATION SOLUTION A	
				CALIBRATION SOLUTION B	
				CALIBRATION CLEANER	
				calib_detergent	
DOCTOR				SAMPLE	
PATIENT					
SELF-DISINFECTION		MENU		START	

- Wait the end of the cycle (green screen with OK button), and push OK button.
- Check the scale and FT13 values and if necessary adjust FT13 value in the CALIBRATION screen (CALIBRATED columns)
- The allowed scales range is min 16ml (ISACLEAN) or 34 ml (INTERCEPT PLUS)
- The allowed FT13 flowmeter is 16ml or 34 ml
- Repeat the procedure from point at least three times with correct values

SENSOR	ZERO POINT	MIN	MAX	CONSTANT	MIN	MAX	READING	CALIBRATED
MAX_V11	0	0	0	10	1	20		
MAX_V12	0	0	0	10	1	20		
SGR	0	0	0	1	1	1	0	0
FT1	0	0	0	18,865	11	14	0	0
FT2	0	0	0	17,085	14,5	17,5	0	0
FT5	0	0	0	1	1	1	0	0
FT6	0	0	0	1	1	1	0	0
FT7	0	0	0	1	1	1	0	0
FT8	0	0	0	1	1	1	0	0
FT9	0	0	0	1	1	1	0	0
FT10	0	0	0	1	1	1	0	0
FT13	0	0	0	0,35	0,25	0,4	0	0
FT14	0	0	0	17,089	10	19	0	0
FT19	0	0	0	1	0,2	0,35	0	0
T1	0	0	0	1,3	1	1	18,25	21,125
T2	0	0	0	1,1	1	1	19,25	21,725
TP1	650	625	675	3	3	5	1580	2790
TP2	650	625	675	0,66	0,6	0,85	653	1,98
TP3	650	625	675	3,9	3	4,5	682	124,8
TP4	650	625	675	3,9	3	4,5	1084	1602,6
TP5	650	625	675	3,9	3,5	5	1927	4980,3
TP7	650	625	675	3,9	3	4,5	716	257,4
PRINT		TECHNICAL MENU		SAVE				

- | INSTRUMENT | | | | CYCLE TYPE |
|---------------|---------|---------------|----|-------------------------------|
| CATEGORY | MODEL | SERIAL NUMBER | ID | |
| Sigmoidoscope | ES-3040 | 1 | 1 | CYCLE NAME |
| Gastroscope | GR-Q165 | A012345 | 1 | FAST DISINFECTION |
| Colonoscope | CF-100L | 123 | | COMPLETE STERILIZATION |
| Duodenoscopy | JF-130 | A123 | | FAST STERILIZATION |
| | | | | CLEANING SAMPLE |
| TEST | | | | SETTING |
| | | | | SET SOLUTION A |
| | | | | SET WATER |
| | | | | SET CLEANER |
| | | | | calc. detergent |
| | | | | FAST REPROCESSING |
| | | | | CALIBRATION SOLUTION A |
| | | | | CALIBRATION SOLUTION B |
| | | | | SET SOLUTION B |
| | | | | |

DOCTOR
SAMPLE

PATIENT
MENU
START

- | SENSOR | ZERO POINT | MIN | MAX | CONSTANT | MIN | MAX | READING | CALIBRATED |
|---------|------------|-----|-----|----------|------|------|---------|------------|
| MAX V11 | 0 | 0 | 0 | 10 | 1 | 20 | | |
| MAX V12 | 0 | 0 | 0 | 10 | 1 | 20 | | |
| SGR | 0 | 0 | 0 | 1 | 1 | 1 | 0 | 0 |
| FT1 | 0 | 0 | 0 | 18,865 | 11 | 14 | 0 | 0 |
| FT2 | 0 | 0 | 0 | 17,085 | 14,5 | 17,5 | 0 | 0 |
| FT5 | 0 | 0 | 0 | 1 | 1 | 1 | 0 | 0 |
| FT6 | 0 | 0 | 0 | 1 | 1 | 1 | 0 | 0 |
| FT7 | 0 | 0 | 0 | 1 | 1 | 1 | 0 | 0 |
| FT8 | 0 | 0 | 0 | 1 | 1 | 1 | 0 | 0 |
| FT9 | 0 | 0 | 0 | 1 | 1 | 1 | 0 | 0 |
| FT10 | 0 | 0 | 0 | 1 | 1 | 1 | 0 | 0 |
| FT13 | 0 | 0 | 0 | 0,35 | 0,25 | 0,4 | 0 | 0 |
| FT14 | 0 | 0 | 0 | 17,600 | 16 | 19 | 0 | 0 |
| FT19 | 0 | 0 | 0 | 1 | 0,2 | 0,35 | 0 | 0 |
| T1 | 0 | 0 | 0 | 1,3 | 1 | 1 | 16,25 | 21,125 |
| T2 | 0 | 0 | 0 | 1,1 | 1 | 1 | 19,75 | 21,725 |
| TP1 | 650 | 625 | 675 | 3 | 3 | 5 | 1580 | 2790 |
| TP2 | 650 | 625 | 675 | 0,66 | 0,6 | 0,85 | 653 | 1,98 |
| TP3 | 650 | 625 | 675 | 3,9 | 3 | 4,5 | 682 | 124,8 |
| TP4 | 650 | 625 | 675 | 3,9 | 3 | 4,5 | 1084 | 1602,6 |
| TP5 | 650 | 625 | 675 | 3,9 | 3,5 | 5 | 1927 | 4960,3 |
| TP7 | 650 | 625 | 675 | 3,9 | 3 | 4,5 | 716 | 257,4 |