Easypump® II ST/LT

en	Instruction for use	10
de	Gebrauchsanweisung	12
bg	Инструкция за употреба	14
cs	Návod k použití	16
da	Brugsanvisning	18
el	Οδηγίες χρήσης	20
es	Instrucciones de uso	22
et	Kasutusjuhised	24
fi	Käyttöohje	26
fr	Mode d'emploi	28
hr	Upute za uporabu	30
hu	Használati utasítás	32
id	Petunjuk penggunaan	34
it	Istruzioni per l'uso	36
kk	Қолдану нұсқаулығы	38
lt	Naudojimo instrukcija	40
lv	Lietošanas instrukcija	42
ms	Arahan Penggunaan	44
nl	Gebruikersinformatie	46
no	Bruksanvisning	48
pl	Instrukcja użytkowania	50
pt	Instruções de utilização	52
ro	Instrucțiuni de utilizare	54
ru	Инструкция по применению	56
sk	Návod na použitie	58
sl	Navodila za uporabo	60
sv	Bruksanvisning	62
tr	Kullanım Kılavuzu	64
uk	Інструкції для застосування	66
vi	Hướng dẫn sử dụng	68
zh	使用说明	70

B | BRAUN

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OVER- AND UNDERFILLING OF EASYPUMP® II ST-LT

PRODUCT SPECIFICATIONS

LT-VARIANTS															
LT-Variants		60 - 12	500 - 125	80 - 16	125 - 25	270 - 27	60 - 30	120 - 30	400 - 40	100 - 50	270 - 54	400 - 80	270 - 68	400 - 100	270 - 135
Code No. (REF)*		4540002	4540003	4540004	4540006	4540008	4540010	4540012	4540014	4540016	4540018	4540022	4540026	4540028	4540032
Nominal Flow Rate (ml/h)		5	40	5	5	10	2	4	10	2	5	5	4	4	2
Nominal Fill Volume (ml)		60	500	80	125	270	60	120	400	100	270	400	270	400	270
Minimal Fill Volume (ml)		30	240	50	60	120	30	60	240	50	135	240	135	240	135
Maximal Fill Volume (ml)		65	560	125	125	335	65	125	560	125	335	560	335	560	335
Maximum Retained Volume (ml)		≤2	≤10	≤3	≤3	≤8	≤2	≤3	≤10	≤3	≤8	≤10	≤8	≤10	≤8

*Suffix of article code may vary according to sales region

STORAGE

TIME TO REACH ROOM TEMPERATURE (APPROX. +23°C ±2°C)														
Hours From Refrigerator (approx. +2°C to +8°C)	6	12	6	12	12	6	12	12	6	12	12	12	12	12
Hours From Freezer (approx. -18°C)	12	18	12	18	18	12	18	18	12	18	18	18	18	18

OVER- AND UNDERFILLING DATA

APPROXIMATE DELIVERY TIME		Fill Volume (ml)														
Hours (h:mm)		Days														
6 h	-	-	240	-	-	-	-	-	-	-	-	-	-	-	-	-
8 h	-	-	40	320	-	-	-	-	-	-	-	-	-	-	-	-
10 h	-	-	50	400	50	-	-	-	-	-	-	-	-	-	-	-
12 h	-	-	60	480	60	60	120	-	-	-	-	-	-	-	-	-
18 h	-	-	-	90	90	180	36	72	-	-	-	-	-	-	-	-
24 h	1 d	-	-	120	120	240	48	96	240	-	-	-	-	-	-	-
30 h	-	-	-	-	-	-	60	120	300	60	150	-	-	-	-	-
48 h	2 d	-	-	-	-	-	-	-	480	96	240	240	-	-	-	-
60 h	2.5 d	-	-	-	-	-	-	-	-	120	-	300	240	240	-	-
72 h	3 d	-	-	-	-	-	-	-	-	-	-	360	288	288	144	-
96 h	4 d	-	-	-	-	-	-	-	-	-	-	480	-	384	192	-
120 h	5 d	-	-	-	-	-	-	-	-	-	-	-	-	480	240	-
144 h	6 d	-	-	-	-	-	-	-	-	-	-	-	-	-	288	-
168 h	7 d	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
192 h	8 d	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
216 h	9 d	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
240 h	10 d	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
264 h	11 d	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
288 h	12 d	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

PRODUCT SPECIFICATIONS

ST-VARIANTS										
ST-Variants	100 - 0.5	250 - 0.5	50 - 1	100 - 1	250 - 1	250 - 1.5	400 - 2	500 - 2	100 - 2	400 - 4
Code No. (REF)*	4540040	4540042	4540044	4540046	4540048	4540050	4540052	4540054	4540056	4540058
Nominal Flow Rate (ml/h)	200	500	50	100	250	175	200	250	50	100
Nominal Fill Volume (ml)	100	250	50	100	250	250	400	500	100	400
Minimal Fill Volume (ml)	50	125	30	50	135	135	240	240	50	240
Maximal Fill Volume (ml)	125	295	65	125	335	335	560	560	125	560
Maximum Retained Volume (ml)	≤3	≤8	≤2	≤3	≤8	≤8	≤10	≤10	≤3	≤10

STORAGE

*suffix of article code may vary according to sales region

TIME TO REACH ROOM TEMPERATURE (APPROX. +23°C ±2°C)										
Hours From Refrigerator (approx. +2°C to +8°C)	6	12	6	6	12	12	12	12	6	12
Hours From Freezer (approx. -18°C)	12	18	12	12	18	18	18	18	12	18

OVER- AND UNDERFILLING DATA

APPROXIMATE DELIVERY TIME										
Hours (h:mm)	Fill Volume (ml)									
00:15 h	50	125	-	-	-	-	-	-	-	-
00:30 h	100	250	-	-	-	-	-	-	-	-
00:45 h	-	-	38	75	188	-	-	-	-	-
01:00 h	-	-	50	100	250	175	-	250	50	-
01:15 h	-	-	63	125	-	219	-	313	63	-
01:30 h	-	-	-	-	-	263	300	375	75	-
01:45 h	-	-	-	-	-	-	350	438	88	-
02:00 h	-	-	-	-	-	-	400	500	100	-
02:15 h	-	-	-	-	-	-	450	-	113	-
02:30 h	-	-	-	-	-	-	500	-	125	250
02:45 h	-	-	-	-	-	-	550	-	-	275
03:00 h	-	-	-	-	-	-	-	-	-	300
03:30 h	-	-	-	-	-	-	-	-	-	350
04:00 h	-	-	-	-	-	-	-	-	-	400
04:30 h	-	-	-	-	-	-	-	-	-	450
05:00 h	-	-	-	-	-	-	-	-	-	500
05:30 h	-	-	-	-	-	-	-	-	-	550

B. Braun Melsungen AG | Hospital Care | 34209 Melsungen | Germany
Tel. +49 5661 71-0 | www.bbraun.com

Note: In order to deliver the Nominal Fill Volume, please add the Maximum Retained Volume:

Example LT 60-12: 60 ml + 2 ml = 62 ml (Nominal Fill Volume + Maximum Retained Volume = Total Filling Volume)

en Instructions for Use

Description

Easypump® II ST/LT disposable elastomeric infusion pump system is a non-electrically driven portable infusion device, enabling patients to be treated in an ambulatory manner. Medication is delivered to the patient by positive pressure applied by the elastomeric membrane. The flow rate is determined by the combination of the flow regulation device (flow restrictor) and the positive pressure of the elastomeric membrane. This pressure delivers the solution against the back-pressure of catheters and blood pressure in the veins.

Materials of the elastomeric pumps that constitute the fluid path are: MVQ, PVC (DEHP-free), PMMA, PES, ABS, glass, e-PTFE.

Easypump® II ST/LT offers two Luer Lock compatible interfaces for other (medical) devices and infrastructure:

- Fill port (for filling the pump i.e. by syringes, automatic filling devices, etc.)
- Patient connector (for patient access devices such as vascular access devices, epidural catheters via filter or catheter connector, catheter/cannula placed subcutaneously, etc.)

Sterility

Do not use if package is damaged. Visually inspect the product packaging to ensure that the sterile barrier system is intact.

Intended purpose

Single use only.

Indication

Easypump® II ST/LT is indicated for delivering a pre-determined amount of medication to the patient via intravenous, subcutaneous or epidural routes (according to pump model and SPCs of drugs) in a continuous and accurate manner. Chemotherapy, antibiotic therapy, anaesthesia and IV pain management are the most common therapies in which Easypump® II ST/LT can be used.

Medications and/or fluids must be administered as per the instructions provided by the drug manufacturer (i.e. as per the Summary of product characteristics (SPCs) of the medication/fluid).

Patient population

Easypump® II ST/LT is intended for patients who receive a treatment in an ambulatory manner. The use of Easypump® II ST/LT depends on the intended infusion therapy as well as the physical and mental conditions of the patient.

Easypump® II ST/LT can be used for both adult and paediatric patients (from 3.5 kg upwards).

Restrictions for use: known contraindications for intended infusion therapy according to the drug/solution SPC and devices.

Intended User

Easypump® II ST/LT should only be used by healthcare worker who have been adequately trained within the educational background in this technique.

Upon medical assessment and adequate training, patients and caregivers may also be allowed to take over defined handling steps according to national guidelines. Healthcare workers are responsible for ensuring that patients and/or caregivers are instructed in the correct use of Easypump® II ST/LT and that the instruction of caregiver and/or patient is documented by the healthcare worker according to local documentation guidelines.

Contraindications

Easypump® II ST/LT is contraindicated for the delivery of blood, blood products, insulin, total parenteral nutrition, lipid emulsions and life-supporting or -sustaining medications as well as for intra-articular use. Easypump® II ST/LT should not be used in patients with known hypersensitivity to any of the materials of the device.

Residual Risks/Side effects

Common complications associated with the use of elastomeric pumps for continuous infusion are:


- Catheter-related complications (such as catheter migration, dislodgement, obstruction, insertion site infection, penetration of the vessel, nerve injury, needle trauma, perineuraxial injection). Note: After initial correct positioning of an intravenous catheter, for example, it may get dislodged. If this is the case, the infusion may be delivered paravenously, i.e. into the tissue. In case of paravenous infusion, severe tissue reactions and necrosis might occur. In such situations, please stop the infusion immediately, leave the catheter in place, contact your physician and check the SPC of the drug/solution. Intravenous catheters may also lead to infusion site infections.
- Tubing-related complications (such as kinking)
- Infusion and infusion pump-related complications (such as inaccurate flow, leakage, obstruction). Note: In case of over-infusion, severe side effects may occur, depending on the infused medication. In case of under-delivery or stopped delivery, the treatment may be compromised. In case of leakage, being exposed to a medication not intended for such use might lead to side effects.
- Drug toxicity. Any drug may lead to side effects and toxicity. Please refer to the specific SPC of the drug manufacturer.

Precautions

- For detailed information about the range of drugs typically administered via elastomeric pumps, please refer to the Drug Stability List.
- Temperature dependency: Easypump® II ST/LT is designed to work at room temperature $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($73^{\circ}\text{F} \pm 3.6^{\circ}\text{F}$). For every 1°C above or below this temperature, the flow rate will increase or decrease by approximately 3%. An increase in temperature results in an increase in flow rate and, in contrast, a decrease in temperature results in a decrease in flow rate. The flow restrictor is calibrated to work at 31°C (88°F). To maintain a stable flow rate the flow restrictor should be in close contact with the patient's skin at all times (31°C).
- If Easypump® II ST/LT needs to be stored in the refrigerator or freezer, allow the unit to warm up to room temperature before usage. **Storage should not exceed a maximum duration of 30 days**, deducting the pump-specific nominal infusion time (Example: Maximum 30 days – nominal infusion time 5 days = 25 days maximum storage).
- If storage in a freezer is intended, the line must not be primed since the filter can be damaged.
- Easypump® II ST/LT is MRI compatible.
- Underfilling and overfilling: Filling the pump to less than the nominal volume generally results in a shorter delivery time. Filling the pump to more than the nominal volume results in a longer delivery time.
- Delivery time can increase significantly as a result of extended storage time.
- Diluent dependency: Easypump® II ST/LT flow rates are calculated on the basis of using 0.9% NaCl. Using dextrose (D5W) as diluent or the addition of any drug of a higher viscosity than normal saline will increase delivery time (i.e. by 10% in case of dextrose [D5W]).
- If the pump is intended to be stored, and dextrose is used as a diluent, the system should only be primed shortly prior to the start of the infusion, as the flow rate may otherwise decrease significantly.
- An uneven inflation of the pump during the filling process has no effect on the flow rate.
- Air bubbles in the silicone sleeve would not affect the functionality of the pump. There is an air vent in the filter which is capable of filtering out the trapped air bubbles. The silicone sleeve has gas permeability properties, therefore the trapped air in the silicone sleeve can be diffused out of the pump.
- Ambient pressure dependency: Easypump® II ST/LT should be used within an ambient pressure between 86 kPa and 106 kPa.
- Back-pressure of catheters and blood pressure affects the flow rate.
- During epidural application the flow rate may increase.

- Easypump® II ST/LT should be positioned at approx. the same level as the venous access device. Placing the pump at a higher or lower level may result in an altered flow rate.

When filled to nominal volume and under the environmental conditions described here, the flow accuracy is $\pm 15\%$ of the labelled flow rate for Easypump® II ST/LT.

Warning

- Do not use if package has been opened or is damaged.
- Single use only. Do not re-sterilize or re-use.
- Re-use of single-use devices creates a potential risk for the patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Do not use alcohol to disinfect the fill port as it may impair product performance.
- Blockage or slow flow can occur due to crystallization/precipitation of certain drugs.
- There is no alarm or alert when flow changes occur to advise the user of the potential for serious injury or death due to under-/over delivery or even stop of infusion.
- External pressure such as squeezing or lying on the pump increases the flow rate.
- Monitor all connections for tightness over the whole infusion procedure.

Operating Instructions

Directions for filling

Elastomeric pumps are to be filled with reasonable care and only by qualified personnel.

Prior to usage, please inspect for any damage. Do not use if damaged.

Elastomeric pumps can be filled with a Luer Lock syringe or other similar filling devices using an aseptic technique. The tubing should be primed with diluent before adding medication.

Priming (use aseptic technique):

1. Unfold the outer layer properly prior to starting the filling process.
2. Close the clamp.
3. Flip open the comfort cap.
4. Twist off the closing cone from the filling port and place it on a sterile surface.
5. Attach the filling device to the filling port.
6. Inject diluent. Repeat as often as necessary.
7. Remove the filling device from the filling port, reattach the closing cone on the filling port and close all caps.
8. Open the closing cone of the patient connector.
9. Open the clamp and prime the system.
10. Close the clamp and reattach the closing cone of the patient connector. Do not re-open the clamp until the pump is connected to the patient.

Specific priming instructions for drugs prone to crystallization/precipitation (i.e. 5-Fluorouracil):

- Follow steps 1–5 mentioned under "Priming".

- Initially fill Easypump® II with approximately 10 ml of diluent.
- Follow steps 7–10 listed under "Priming".

Note: If not primed correctly, blockage or slow flow can occur due to crystallization/precipitation of certain drugs in the area highlighted below.

Filling – Adding medication (use aseptic technique):

Always refer to drug manufacturer's package insert for drug reconstitution/dilution, storage procedures and handling.

Repeat steps 2–6 until required volume is achieved.

After adding the medication, flush the filling port with 0.9% sodium chloride solution to avoid high concentration of the drug inside the filling area.

Note: Please make sure that the closing cones of the filling port and patient connector are properly closed after filling and prior to use!

Upon completion, the tubing will be filled exclusively with diluent, to prevent crystallization/precipitation, while the pump reservoir will contain the medication until the infusion is started.

Patient Connection:

1. Open the closing cone of the patient connector.
2. Connect the patient connector of the pump to the patient's access device.
3. Make sure that the flow restrictor is completely taped on the patient's skin.

Note: Length of flow restrictor tubing may vary, because each flow restrictor is individually calibrated.

4. Make sure, that the filter is not covered by any dressing.
5. Open the clamp to start the infusion.

End of Therapy:

1. Make sure to close the clamp at the end of therapy and disconnect the patient connector.

Note: Due to higher pump pressure, blood backflow cannot happen, however blood can diffuse with the solutions in the tubing.

Duration of use

See over- and underfilling table in the Appendix.

Storage and handling conditions

Keep away from sunlight. Keep dry.

Disposal

The product Easypump® II ST/LT is made from raw materials that, both in combination and as individual components, do not need to be disposed of separately. If used correctly according to the instructions for use, no thermal or chemical influences can alter the used materials in a manner that would require separate disposal. Disposal of used devices must be ensured by the healthcare provider to be in accordance with clinical waste management regulations in each respective country where the medical devices are used.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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[CH REP] B. Braun Medical AG