

ASPARTATE AMINOTRANSFERASE (AST/GOT)

COD 12531 5 x 40 mL + 5 x 10 mL

Only for in vitro use in the clinical laboratory



ASPARTATE AMINOTRANSFERASE (AST/GOT)

IFCC

INTENDED USE

Reagent for the measurement of aspartate aminotransferase (AST or GOT) concentration in human serum or plasma. The obtained values are useful as an aid in the diagnosis and monitoring of disorders of liver diseases, especially acute ones.

This reagent is for use in the BioSystems A25 and A15 analyzers or in other analyzer with similar performance characteristics.

CLINICAL SIGNIFICANCE

The aminotransferases catalyze the formation of glutamic acid from 2-oxoglutarate by transfer of amino groups. AST is found in highest concentration in the liver and heart muscle but it is also abundant in skeletal muscle, kidney and pancreas.


The serum concentration of AST is elevated in hepatitis and other forms of hepatic disease associated with necrosis: infectious mononucleosis, cholestasis, cirrhosis, metastatic carcinoma of the liver, delirium tremens, and after administration of various drugs^{1,2}.

Serum AST concentration is also elevated after myocardial infarction, in skeletal muscle disease (as progressive muscular dystrophy), in acute pancreatitis or hemolytic disease and other^{1,2}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Aspartate aminotransferase (AST or GOT) catalyzes the transfer of the amino group from aspartate to 2-oxoglutarate, forming oxalacetate and glutamate. The catalytic concentration is determined from the rate of decrease of NADH, measured at 340 nm, by means of the malate dehydrogenase (MDH) coupled reaction³⁻⁶.



CONTENTS AND COMPOSITION

A. Reagent: 5 x 40 mL. Tris 121 mmol/L, L-aspartate 362 mmol/L, malate dehydrogenase > 460 U/L, lactate dehydrogenase > 660 U/L, pH 7.8.

WARNING: H315: Causes skin irritation. H319: Causes serious eye irritation. P280: Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P332+P313: If skin irritation occurs: Get medical advice/attention.

B. Reagent: 5 x 10 mL. NADH 1.9 mmol/L, 2-oxoglutarate 75 mmol/L, sodium hydroxide 148 mmol/L, sodium azide 9.5 g/L.

WARNING: H302: Harmful if swallowed. EUH031: Contact with acids liberates toxic gas. P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. P330: Rinse mouth.

For further warnings and precautions, see the product safety data sheet (SDS).

STORAGE AND STABILITY

Store at 2-8°C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 15 days.

Indications of deterioration: Absorbance of the blank below the limit indicated in "Test Parameters".

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

– Pyridoxal phosphate (BioSystems cod. 11666):

C. Reagent: Pyridoxal phosphate AST 10 mmol/L. 5 mL.

– Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044).

REAGENT PREPARATION

Reagents are provided ready to use.

Working Reagent: Pour the contents of the Reagent B into the Reagent A bottle. Mix gently. Other volumes can be prepared in the proportion: 4 mL Reagent A + 1 mL Reagent B. Stable for 1 month at 2-8°C.

Working Reagent with Pyridoxal Phosphate (Note 1): Mix as follows: 10 mL of Working Reagent + 0.1 mL of Reagent C (cod 11666). Stable for 6 days at 2-8°C and for 6 days in the refrigerated compartment of the analyzer.

SAMPLES

Serum and plasma collected by standard procedures.

Aspartate aminotransferase in serum and plasma is stable for 7 days at 2-8°C. Use heparin as anticoagulant⁸.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 15 days (6 days for Working Reagent with Pyridoxal Phosphate), after reagent lot change or as required by quality control procedures.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, 18009 and 18042) and II (cod. 18007, 18010 and 18043) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

REFERENCE VALUES

| Reaction temperature | 37°C | 30°C |
|-----------------------------------|---------------------------|---------------------------|
| Without pyr-P, up to ¹ | 40 U/L = 0.67 μ kat/L | 25 U/L = 0.42 μ kat/L |
| With pyr-P, up to ^{3,4} | 50 U/L = 0.83 μ kat/L | 30 U/L = 0.50 μ kat/L |

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using an A25 analyzer. Results are similar with A15.

- Detection limit: 3.39 U/L = 0.056 μ kat/L.
- Linearity limit: 500 U/L = 8.33 μ kat/L.
- Precision:

| Mean concentration | Repeatability (CV) | Within-laboratory (CV) |
|------------------------------|--------------------|------------------------|
| 33.0 U/L = 0.549 μ kat/L | 3.0 % | 5.1 % |
| 170 U/L = 2.83 μ kat/L | 1.7 % | 4.7 % |

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

LIMITATIONS OF THE PROCEDURE

- Interferences: Lipemia (triglycerides 2 g/L) interfere. Bilirubin (20 mg/dL) and Hemolysis (hemoglobin 50 mg/dL) do not interfere. Other drugs and substances may interfere⁷.

NOTES

- The IFCC recommended method specifies the addition of pyridoxal phosphate.

BIBLIOGRAPHY

- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2005.
- Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
- Sociedad Española de Química Clínica, Comité Científico, Comisión de Enzimas. Método recomendado para la determinación en rutina de la concentración catalítica de la aspartato aminotransferasa en suero sanguíneo humano. Quim Clin 1987; 6: 235-239.
- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C, Part 5. Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002; 40:725-733.
- IFCC reference procedures for measurement of catalytic concentrations of enzymes: corrigendum, notes and useful advice. Clin Chem Lab Med 2010; 48: 615-621.
- Gella FJ, Olivella T, Cruz Pastor M, Arenas J, Moreno R, Durban R and Gómez JA. A simple procedure for routine determination of aspartate aminotransferase and alanine aminotransferase with pyridoxal phosphate. Clin Chim Acta 1985; 153: 241-247.
- Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
- World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev.2, 2002.

TEST PARAMETERS

These reagents may be used in several automatic analyzers. Specific instructions for application in many of them are available on request.

R1: Use Reagent A, or Reagent A with Pyridoxal Phosphate, R2: Use Reagent B.

| | A25 | A15 |
|------------------------------|--------------|--------------|
| GENERAL | | |
| Name | AST | AST |
| Sample type | SER | SER |
| Analysis mode | kinetic mon. | kinetic mon. |
| Units | U/L | U/L |
| Turbidimetry test | no | No |
| Decimals | 0 | 0 |
| Type of reaction | decreasing | decreasing |
| PROCEDURE | | |
| Reading mode | monoch. | monoch. |
| Main filter | 340 | 340 |
| Reference filter | - | - |
| Sample | 25 | 25 |
| Vol. R1 | 300 | 300 |
| Vol. R2 | - | - |
| Washing | 1.2 | 1.2 |
| Reading 1 (cycle) | 7 (12)* | 5 (8)* |
| Reading 2 (cycle) | 18 | 12 |
| Reagent 2 (cycle) | - | - |
| Predilution factor | - | - |
| CALIBRATION AND BLANK | | |
| Calibration type | multiple | multiple |
| Number of calibrators | - | - |
| Calibration curve | - | - |
| OPTIONS | | |
| Blank absorbance limit | 1.400 | 1.400 |
| Kinetic blank limit | - | - |
| Linearity limit | 500 | 500 |
| Substrate depletion | 0.5 | 0.5 |

*Application for Working Reagent with Pyridoxal Phosphate.

