AW Family

AW VolumeShare 7

The AW Family is GE Healthcare's ecosystem of processing and advanced visualization solutions across platforms and applications.

Advantage Workstation (AW) with VolumeShare 7 enables multi-modality image review, comparison, and processing with simplicity and power. Featuring, 64-bit technology, this workstation is a cornerstone solution to Radiology departments around the world.

This document is intended to overview the features, specifications, use cases, and other key information of AW VolumeShare 7.

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SECRET DE AFACERI

Contents

Overview	3
Key Features	3
Industry Standards.....	3
DICOM Conformance Standards:.....	3
Filming Protocols	3
Indications for Use.....	4
Regulatory Compliance	4
Product Details	5
Modalities	5
License Management.....	5
User Interface	5
Standard Applications & Capabilities.....	6
Workflow Management.....	9
Security Capabilities.....	11
System Components	11
Workstation Configuration	11
About GE HealthCare	13

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SECRET DE AFACERI

Overview

Key Features

- Simple drag and drop action for networking, media interchange, and filming.
- Multi-tasking capabilities.
- Active directory integration allows enterprise level user authentication. Directories supported include Microsoft® Active Directory® (MSAD) and other LDAP authentication services.
- Patient list to manage images from local workstation or remote systems.
- Postfetch feature, designed to optimize your reading workflow by automatically gathering a patient's prior exams according to your pre-set criteria.
- Support for Key Images.
- End Review automates routine filming and networking tasks with just one click.
- Search Advantage for fast and easy search of a patient's exam history on PACS or any other DICOM®-compliant device.
- Enhanced Quick Filters of the Patient List filters studies by Modality, Date, End Review status or Exam Description.
- DICOM CD/DVD/USB Creation Tool.
- 2D Viewer for image display, manipulation, annotation, review.
- Integrated Filmer with enhanced flexibility to perform all filming and data exporting tasks.
- Access to a wide variety of applications for greater diagnostic flexibility.

Industry Standards

AW VolumeShare 7 complies with a wide variety of industry standards to facilitate adoption of features and performance improvements as the computing and medical imaging industry evolves.

DICOM Conformance Standards:

- DICOM 3.0 Storage Service Class for RT, CT, MR, CR, X-ray (Angio and R&F), Digital X-ray (DX), MG, NM, PET, U/S, Secondary Capture, Secondary Capture Color DICOM Image Objects. (Service Class User (SCU) for image send and Service Class Provider (SCP) for image receive).
- DICOM 3.0 Query/Retrieve Service Class (SCU and SCP).
- DICOM 3.0 Storage Commitment Service Class (SCU).
- DICOM Print (Color and B&W).
- DICOM Media Interchange (CD-R, DVD+R(W)).

Filming Protocols

DICOM Print (Color and B&W) and Adobe® Postscript (Color and B&W) for supported printers.

CONFIDENTIAL

SECRET DE
AFACERI

Indications for Use

AW VolumeShare 7 is a review workstation, which allows easy selection, review, processing and filming of multi-modality DICOM images from a variety of diagnostic imaging systems. When interpreted by a trained physician, filmed or displayed images on the AW monitor may be used as a basis for diagnosis, except in the case of mammography images.

Regulatory Compliance

 This product complies with the Regulation (EU) 2017/745 of the European Parliament and the Council on the medical device (MDR).

This This product or its features may not be available in some countries or regions. Please contact your sales associate.

Rx Only

 **CONFIDENTIAL**

**SECRET DE
AFACERI**

Product Details

The Advantage Workstation VolumeShare 7 software technology makes a multi-modality advanced visualization workflow solution that helps to enhance diagnostic confidence and productivity. This solution features software that is optimized for 64-bit technology and multicore processor hardware to provide leading edge performance. It includes AW's premier 3D image analysis package, Volume Viewer, with a suite of volumetric visualization, and analysis tools for CT, MR, 3D X-ray and PET.

In addition, this package includes the multi-modality 2D Viewer with dedicated functions for review of CT, MR, X-Ray, Angio, DX, U/S and PET images. It also includes the Filmer, a multimedia export tool for creating electronic films, filming, and exporting to an internal web server, CD, or DVD.

AW VolumeShare 7 is also available for purchase without Volume Viewer for users who require only 2D image applications. Volume Viewer dependent features will not be supported with this configuration.

Modalities

Supported DICOM modalities include:

- Computed tomography (CT),
- Magnetic resonance imaging (MR),
- Radiofluoroscopy (RF),
- Xray angiography (XA),
- Computed radiography (CR),
- Digital radiography (DX),
- Mammography (MG),

- Nuclear medicine (NM),
- Positron emission tomography (PET),
- Ultrasound (US)
- Structured Reporting (SR),
- Key Objects (KO)

License Management

- Application Usage Monitor tracks how often a particular license was used, when, and by whom. You can compile and view usage reports for a specific period of time. This helps you understand how your systems contribute to departmental productivity and to plan for future license requirements based on actual usage data.
- For floating licenses, a secondary license server can be assigned for redundancy. If the primary license server becomes unavailable, the secondary can be contacted so that work can continue.

User Interface

- Displays on one or two 1280x1024 monitors with scroll or optical mouse.
- Simple drag-and-drop mechanism for networking, media interchange, filming and data export options.
- Automatic help message display briefly describes function currently under the mouse pointer.
- Shortcut keys and programmable window/level function keys for accelerated control.
- Multi-tasking capabilities that allow Patient List, 2D Viewer, Filmer and Volume Viewer to run simultaneously with Fast Switch capabilities between applications.

CONFIDENTIAL

SECRET DE
AFACERI

Standard Applications & Capabilities

Separate application licenses must be purchased to access advanced applications functionality.

2D Viewer

The 2D Viewer is an application used to display, manipulate, annotate and review 2D images by a trained physician for diagnostic interpretations.

Display Customization

Display customization allows the user to manage layouts to display data. Key benefits include:

- Allows modification of number of exams or series displayed to facilitate either single or multi-exam reviews.
- Standard layouts provide flexibility to tailor the image display from 1x1 to 8x8. The current layout is retained if the images are from the same modality and orientation.
- Annotation levels allow selection of the image information fields to display.
- Toolbar customization enables control of which buttons are displayed based on user preference.

Study Navigation

- The navigator lets users assign a series to a view on the fly.
- Cine mode also provides temporal, spatial, or manual playback loops.
- With two exams or series loaded, Cine mode supports a side-by-side display format with synchronized playback loops for more efficient comparisons.

Image Review

The initial image window and level setting is based on the DICOM header. Once displayed, several methods of adjusting image window and level are provided.

Provides routine image manipulation features:

- Flip/Rotate, Zoom, Pan, Magnifying Glass,
- Inverse Video: inverts grayscale color map,
- Display normal: lets you return the image to its default viewing parameters.

You can access the following features using a single mouse click directly on the image. Having these direct manipulation tools easily at hand gives you fast access with less distraction from your review task by eliminating the need to return to the graphical user interface controls: 2D distance, Angle, Report cursor, Box ROI, elliptical ROI and free-hand ROI.

The program continuously updates statistics on the fly.

- The Annotation feature lets you highlight areas of interest by adding text and line/arrow pointers to any image feature.
- The Copy/Paste/Eraser feature lets you copy/paste/remove any text or graphic placed on an image.
- Cross-reference indicates the position of the current slice over the scout or localizer or any other non-parallel series.
- Save lets you store a copy of the image as it appears on the screen for future review.
- The Key Images feature allows you to flag images and create key objects.

Advanced X-Ray Analysis

The user can apply shutters to X-ray images to focus on specific areas within the image. Shutters are black, opaque overlays with elliptical or

CONFIDENTIAL

SECRET DE
AFACERI

rectangular cutouts that you can size as desired. Once a shutter is applied, you can move the underlying image to shift the viewable area. Shutters can also be applied to CT and MR images with the Image Matte feature.

The user can apply several levels of edge enhancement filters to an image.

- **Image Subtraction:** A graphical user interface lets you select a mask and subtraction of all images in an associated sequence. Subtracted images may be saved as a sequence with one mouse click.
- **Landscape:** With this tool, you can introduce a percentage of the mask image into the subtracted image for anatomical reference. A graphical user interface lets you specify the percentage of mask to be reintroduced.
- **Pixel Shift:** Auto pixel shift optimizes a mask to image registration. You can also move a mask manually to optimize subtraction quality in a specific region of interest.
- **Split Pixel Shift:** You can split the screen horizontally or vertically for pixel shift.
- **Maximum/Minimum Opacifications:** You can integrate selected images to provide a resulting Max. Op./Min. Op. image.

Batch Filming

The Print Series feature lets users automatically batch film an entire series with a single keystroke.

Volume Viewer

Volume Viewer is a multi-modality AW Family advanced visualization application. It provides excellent 3D visualization and processing capabilities for reading and comparing CT, MR, 3D X-Ray, PET, PET/MR, and PET/CT datasets. It

features a breadth of high-performance analysis tools such as Multiplanar Reformat, Curvilinear Reformating, Volume Rendering and Lumen Navigation.

Volume Viewer is the pre-requisite application for what are known as Volume Viewer Foundation (VVF) applications within the AW Family Advanced Visualization Application portfolio. VVF applications share a similar user interface and have some interoperability. For detailed description of Volume Viewer features, please refer to the Volume Viewer Product Data Sheet.

Filmer

The integrated Filmer enhances the efficiency of the review station and gives users greater filming and exporting flexibility. This feature supports two modes: the Mini Filmer mode and Full-Screen mode, which provides the ability for customization and film layout template creation. Three key mechanisms give users flexibility:

- Free format filming
- DICOM structured reporting (SR)
- Data export (HTML/PDF and JPEG, PNG, MPEG, AVI, or QTVR)

With the Filmer, users can easily extract significant images from any AW Family AV application (2D Viewer, Volume Viewer, Ready View, etc.).

- You can film images individually by dragging and dropping to the on-screen Filmer, or by a single F1 keystroke.
- Multiple image formatting allows filming multiple images in a single page frame with the F2 keystroke.
- With Film MID users can send multiple images to a single Filmer frame.

CONFIDENTIAL

SECRET DE
AFACERI

- Batch filming is supported by applications that provide that capability (Print Series in the 2D Viewer, Batch Film Protocols in Volume Viewer).

Communication between applications and the Filmer is accomplished in the Mini Filmer mode, which provides a minimized footprint. The Mini Filmer mode provides the following:

- Compatibility with Batch Filming from Volume Viewer (may be optional in your configuration)
- Store and position images transferred from an application.
- Type of export (film, media, database).
- Rapid switching between Full-Screen Filmer and application.

The resultant electronic films become a quick summary of the patient study and radiological interpretation which can then be reviewed by clinicians and physicians. They can also be saved as independent files (DICOM SR and Secondary Captures) for teaching purposes. Electronic film can contain one or several pages with specific layouts for each page or all pages.

A flexible Edit Mode provides the ability to easily add, manipulate, format, or delete images from the film. Images can contain text and graphics from measurements and your annotations, and may be window/leveled, magnified, flipped, rotated, or cine. Additional annotations can be added to the image in edit mode.

Preview Mode displays the film as it will be printed or exported, considering the layout applied to each page and the compression level specified for non-DICOM exporting (JPEG/PNG and MPEG/AVI/QTVR).

Printing

AW VolumeShare 7 includes Network DICOM Print (B&W and color). For supported printers, AW VolumeShare 7 includes Network PostScript capability (B&W and color). For PostScript and DICOM printers, any printing format created in the Filmer is supported (e.g., non-square matrix formats for rectangular images such as CT run-offs).

- AW VolumeShare 7 postscript printing has been tested on the following devices: Codonics 1660M, 1660MD or Horizon, Lexmark Optra 1650N, 1855N, SC1275N, C710N, C720N, T612 or T614, Seiko 1720D, Kodak DMI3600, Quantum GL2101HD, declared with film/thick paper or with plain paper, Tally T8106, HP LaserJet and Xerox Phaser
- Digital cameras and analog cameras using a 3M-952 protocol (including DASM interface) are not supported on AW VolumeShare 7.

Exporting

Data export is integrated in the Filmer, thus providing all image processing tools needed for multimedia image export. The Filmer exports any electronic film to DICOM SR, PDF/HTML, JPEG, PNG, MPEG, AVI or QTVR format. AW VolumeShare 7 enhanced security of the export functionality with the use of SFTP connectivity. Non-DICOM data can also be saved on a multi-session CD/DVD. The Data Export capability is intended only for publishing and communication, not for diagnostic purposes. Its simplicity is reflected in the different export mechanisms available:

- CD/DVD removable media
- USB flash drive
- Network HTTP and FTP protocols

CONFIDENTIAL

SECRET DE
AFACERI

Cardiac Review and Export

Processing and review of CT, MR and PET cardiac exams with manual oblique reformatted protocols can be exported as a multi-phase Cine movie that allows referring physicians to review exams in a dynamic mode.

Workflow Management

AW VolumeShare 7 provides features designed to optimize your reading workflow:

- The Patient List provides tools that enable you to sort and filter imaging studies stored on the local workstation or on remote systems.
- Postfetch allows you to retrieve a patient's prior related DICOM exams from a remote DICOM host. Retrieval can be triggered by the arrival of a new patient study on the workstation from the network or supported media.
- Support for IHE Key Image Notes profile, allowing you to mark certain images in the 2D or 3D viewer as key images. Key images are displayed as a separate series in the exam list and accessed or sent to an IHE KIN-supporting PACS. Specific tags and Key Image Notes may be created using the 2D Viewer application.
- Support for external USB disks to serve as DICOM storage media. The USB disk can be accessed for reading and writing via the media button, the same way you would access a CD or DVD. The maximum number of images that can be stored on a USB device is restricted to 300,000 across all exams. USB media must be formatted using the FAT32 file system. See the AW VolumeShare 7 Basic Display/Viewer/Filmer User Manual for more details on USB storage.

- Preference sharing among users. Sharable preferences include: Filmer layouts, Volume Viewer custom protocols, and viewer preferences. A standard USB flash drive may also be used to share preferences between workstations.
- One-Touch protocols enable you to define an application or preset protocol to launch automatically based on DICOM elements.
- DICOM Query/Retrieve Storage Class User (SCU) and Storage Class Provider (SCP) provide seamless network integration.
- DICOM Storage Commitment SCU lets you know when exams have been archived successfully on DICOM devices, such as PACS, that support Storage Commitment SCP.

End Review

End Review automates the routine tasks required at the end of reviewing each exam. The "End Review" flag in the patient list allows you to mark exams as "Done" after post-processing has been completed. End Review allows you to automatically perform one or several of these actions:

- Print pages prepared in the Filmer to the default printer (DICOM or Postscript filmers supported), with the option to automatically clear the Filmer after printing.
- Save electronic films in the workstation's attached DICOM database.
- Automatically push the entire exam or only the series created on the AW to one or more remote hosts.

Remote Network Host Management

AW VolumeShare 7 Patient List provides a control panel to select a remote host or destination for

CONFIDENTIAL

SECRET DE
AFACERI

network transfer. Icons signify different device types (acquisition, post processing, PACS, etc.) that are accessible to the workstation via the network.

You can initiate a network transfer simply by dragging an exam, patient, series, or image(s) to the remote host icon. You can activate a remote browser to display detailed exam information from the remote host, including filtering information if supported.

Search Advantage

This advanced search engine lets you search a patient's exam history on PACS or other DICOM device with a few mouse clicks.

Worklists and Filters

You can filter the patient list by modality, date, end review status, or exam description. Most filters are available on the remote host patient list as well.

You can filter the patient list further by choosing one or any of these parameters:

- Modality.
- Patient name or patient ID.
- Exam location (hospital name) and exam description.
- Series description
- Date and time of day: today, or specified date or date ranges with a specified time or time ranges.
- Radiologist's or referring physician's name.
- Accession number.

Quick Sort (ascending or descending) and Quick Access (entry field) are available for any of the fields displayed at exam level (patient name, patient ID, exam location and description,

radiologist's and referring physician's name, date and time).

Queues

AW VolumeShare 7 manages three queues: network transfer; media interchange, filming. A menu on the patient list gives you easy access to queue status, and an animated icon gives you continuous network queue status information.

DICOM CD/DVD/USB Creation Tool

- The DICOM media creation tool offers you the flexibility to save and retrieve from supported CD, DVD, or USB storage devices.
- Add selected exam, series, or images by dragging and dropping to the pre-mastering window.
- CD and DVD usage percent is shown before you launch the writing process.
- You can optionally down sample certain 1024²x512² X-ray angiographic images during CD creation.
- CD/DVD composer lets you manage large amounts of data. If data size is larger than medium capacity, you are prompted to supply additional media storage space.
- You can specify the number of copies created during a CD/DVD save session.
- You can embed a DICOM viewer on DICOM CD/DVD so data can be reviewed on a PC running a Windows operating system.
- You get lossless JPEG reading of any DICOM CD/DVD media. Lossless writing is available for X-ray exams.
- CD/DVD drives operate at 16x write speeds and read speeds up to 48x depending on choice of media.

CONFIDENTIAL

SECRET DE
AFACERI

- 1,600,000 512*512*16 bits uncompressed images OR
- 400,000 1024*1024*16 bits uncompressed images OR
- 100,000 2048*2048*16 bits uncompressed images
- Internal DVD Writer drive for read/write of DICOM CD/DVD media, read/write of Data Export CD/DVD data and service use (DVD Install)
- Integrated dual Ethernet 10/100/1000 Mbit/s ports.
- Additional single Ethernet 10/100/1000 Mbit/s port
- **1 USB QWERTY (or regional) Keyboard and Mouse**

Footprint

- Height 38.60 cm (15.2 in.)
- Width 16.89 cm (6.65 in.)
- Depth 44.47 cm (17.5 in.)
- Approximate Weight 17 kg (38 lbs.)

Operating Environment

- Temperature: +5°C to +35°C (40° to 95° F)
- Humidity: 10% to 85% (relative non-condensing)
- Altitude: 0 to 5000m (16,404 ft.)
- Acoustics: LWAd less than 4.3 Bels
- Shock: 40 G peak, half-sine, 2-3 ms

Non-Operating Environment

- Temperature: -40°C to +60°C (-40° to 140° F)
- Humidity: 10% to 90% (relative, non-condensing)
- Altitude: 0 to 12192m (37,000 ft.)

Monitors

- **(2) 19" color Flat Panel LCD monitors**

- DICOM Part 14 factory calibrated
- Native resolution 1280 x 1024 (5:4 aspect ratio)
- Approximate Weight: 5.3 Kg (11.68 lbs.).
- AC 100 - 240 V : 50 / 60 Hz

Image Networking

- Standard 10/100/1000 Base-T Ethernet for DICOM
- 1000 Base-T dedicated network for optimal Direct Connect performance
- Protocols supported:
 - DICOM 3.0 Storage SCU/SCP and Query/Retrieve SCU/SCP
 - RSVP
 - TCP/IP network layer
- SdCNet is no longer supported
- AW VolumeShare 7 does not support the AdvantageNET network protocol.
- AW VolumeShare 7 does not support DICOM images from GE Healthcare Signa™ version 5.x 1.5T MR systems
- AW VolumeShare 7 software is supported on the previous generation AW VolumeShare 5 HP Z820 and Z440 workstations.

CONFIDENTIAL

SECRET DE AFACERI

GE HealthCare

About GE HealthCare

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 100 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from prevention and screening, to diagnosis, treatment, therapy, and monitoring. We are an \$18 billion business with 51,000 employees working to create a world where healthcare has no limits.

Follow us on [Facebook](#), [LinkedIn](#), [Twitter](#), [Instagram](#) and [Insights](#) for the latest news, or visit our website [gehealthcare.com](#) for more information.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

GE and the GE monogram are trademarks of General Electric Company. Used under trademark license.

Adobe is either a registered trademark or trademark of Adobe Systems Incorporated in the United States and/or other countries.

ClamAV is a registered trademark of Cisco Systems.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

Intel and Xeon, are either registered trademarks or trademarks of Intel Corporation in the United States and/or other countries.

Microsoft, Active Directory, and Windows, are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Novell and eDirectory are registered trademarks of Novell, Inc. in the United States and other countries.

Oracle and Java are registered trademarks of Oracle and/or its affiliates.

All other product names and logos are trademarks or registered trademarks of their respective companies.

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SECRET DE AFACERI

Volume Viewer

GE Healthcare's premium 3D Advanced Visualization and Image Processing platform brings you powerful capabilities to help you visualize and interpret your multimodality imaging data with confidence and ease.

AutoLaunch

Preprocessing enabler

Imaging modalities provide a wealth of diagnostic information, but also present a challenge to your Radiology workflow as image volumes continue to increase while your time for reading and reporting decreases.

Volume Viewer provides you with a rich 3D image processing toolset aimed at creating and displaying the views you need with little user input and streamlining interpretation and reporting by providing the visualization tools you need with minimum clicks.

With an intuitive, modern user interface, Volume Viewer helps a wide diversity of users to learn and master the expanding portfolio of sophisticated tools and applications that it provides. This customizable user interface maximizes the real estate allocated for displaying images and provides access to the users' favorite tools directly from the image viewport. In addition, a customizable toolbar provides one click access to commonly used tools to facilitate a productive review.

Overview

Volume Viewer provides excellent 3D visualization and processing capabilities for reading and comparing CT, MR, 3D X-ray, PET, PET/MR and PET/CT datasets. Volume Viewer also features a broad portfolio of high-performance analysis tools, automating routine tasks and helping to make 3D image processing a stress-free component of your routine workflow.

Volume Viewer is available on VolumeShare 7, a multi-modality advanced visualization workflow solution that helps to enhance diagnostic precision and productivity.

Highlights

- Customizable protocol-driven workflows supporting CT, MR, PET and Innova 3D XA modalities, designed to help enhance your clinical review productivity
- Compare multiple exams from different modalities, at loading² and on the fly²
- Adaptive layouts that take advantage of the dual monitor capabilities in landscape and portrait orientations
- Easy point-and-click segmentation tools for segmenting structures of interest
- CT Dual Energy and GSI⁶ support with dedicated workflows and tools
- Advanced 4D PET support with Q.Freeze⁶ tool and 4D PET Dynamic workflow
- Interactive Summary Table collecting measurements as they are deposited on the images.
- Rich set of export tools including export for 3D printing³.

General Features

- **Unified user interface**, with one click access to your favorite tools and more space for displaying images.
- **Customizable protocol-driven workflows** designed to help enhance your clinical review productivity for a wide range of care areas, with 3D and 4D real time interaction.
- **"Smart layout"** protocol, which adapts the screen layout to the type of loaded exams.
- **High resolution, real-time rendering modes:** MIP/Min IP/Average; standard volume rendering and enhanced volume rendering with Volume Illumination¹ on the fly; Navigator view for endoluminal fly-through; lumen and curved views.
- **Compare mode** for reviewing multiple exams and modalities, simplified with **Dynamic Load to load and register²** new series into the current review session.
- Rich set of **2D/3D ROI tools**, including **AutoContour** for 3D semi-automatic contouring of structures of interest (CT, MR, PET), and **ROI color coded** for tissue classification based on voxel values.
- **Advanced tools**, to take advantage of full 3D capabilities: **AutoSelect**, for easy point-and-click segmentation; **One or two clicks Quick Vessel Trace** to analyze all vessels, in curved reformat, lumen, or MPR view; **Advanced 3D processing tools**.
- **Summary Table** extends reporting capabilities, while collecting measurements as they are deposited on the images.

- **Multiple Export capabilities**, in different DICOM and non-DICOM formats, including **Key Image Notes**.
- **3D Suite³:** STL, OBJ, VRML and 3MF file formats to export single or multi mesh objects, for 3D printing. 3D Suite is not intended for clinical use.
- **Save state**, which allows to save and restore to pause and resume the review at a later time or to facilitate collaboration workflows.

Modality Features

- **Dual energy and GSI⁴ GE Healthcare CT images** supported with dedicated review protocols.
- **Bone VCAR⁵** combines spine labelling on the fly for all CT review and a dedicated layout to improve spine review.
- **MR review support:** Recognition of PSD name, and specific MR parameters at image loading, which simplifies protocol customization.
- **Dedicated review workflows for PET/CT and PET/MR images**, including **PET SUV measurements**.
- **Q.Freeze⁶** combines the quantitative benefits of 4D PET gated imaging into an image that provides both frozen patient motion and reduced image noise.
- **Dedicated protocol-driven workflows for Innova 3D XA images⁹.**

Platform features

- **"Smart Compression"** technology for fast and responsive client performance on AW Server.
- **Fast automated Pre-processing⁷**

- **Multi-sessions switch on AW workstation with AutoLaunch⁸.**
- **Support of Chinese, Japanese, Korean and Russian languages.**

System Requirements

- **AW Server 3.2 and above**, with recommended monitor resolution up to dual 2MP or single 3MP.
- **AW 4.7 Workstation and above** (some performances limitations can occur with Z800)
- **CentricityTM Universal Viewer**

Intended use / Indications for use

Volume Viewer is a medical diagnostic software that allows the processing, review, analysis and communication of 3D reconstructed images and their relationship to originally acquired images from CT, MR, X-Ray Angiography and PET Scanning devices. The combination of acquired images, reconstructed images, annotations and measurements performed by the clinician are intended to provide to the referring physician clinically relevant information for diagnosis, surgery and treatment planning.

Regulatory compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC. This product or its feature may not be available in some other countries or regions. Please contact your sale associate.

Rx Only

Required system configuration:

System Requirements may not be available in all countries or regions.

Please contact your sales associate.

Minimum hardware requirements:

Minimum 32 GB RAM (32GB)

Minimum 2000 MB free disk space

Minimum 4GB free GPU memory

Minimum 2 free cores

Operating System: Windows 7 or later with the following configuration: WindowsTM 7 (64-bit), 64-bit, 32GB RAM, 2000 MB free disk space, 4GB free GPU memory, 2 free cores, 2000 MB free disk space. Configuration is not for use in third party environments.

Product Availability: subject to 3D reconstruction.

Minimum system requirements may vary.

CONFIDENTIAL

SECRET DE AFACERI

© 2019 General Electric Company.
All rights reserved. Data subject to change.
GE, the GE Monogram, Imagination at work, Centricity and AutoBone are trademarks of General Electric company
DOC1460639
GE Medical Systems SCS

Volume Viewer


General Features

Unified user interface

- The user interface provides one click access to favorite tools and more space for clinical images.
- The page layout is organized to facilitate navigation through the review steps. Guided protocols walk the unfamiliar user through a clear workflow of creating and manipulating volumetric data, helping to reduce the learning curve by providing a consistent workflow.
- The customizable toolbar can be docked at the top, bottom, left, or right, by a simple drag & drop. This preferred location is saved for future use for each user.
- An extensive set of advanced tools are grouped into categories on the toolbar to facilitate quick retrieval when needed.
- Each category and tool are easily configurable by a simple drag & drop in the Toolbar so that they are only one click away for your future reviews
- To reduce long mouse trips, up to 6 favorite tools can be accessible from each image via a right mouse menu.
- Customization of the toolbar and tool behaviors with a user's preferences is provided in a single intuitive Preferences menu.
- Keyboard shortcuts and quick access tools for on-viewport controls

The following mouse modes are available to help you manipulate the images:

- Interactive VR adjustment lets you adjust VR opacity interactively with the mouse as Window Level/Window Width
- Direct Paging (continuous or contiguous)
- Free hand rotation on 3D and oblique views
- Percentage fusion, to easily change the transparency of objects on fused images. There are several fusion modes available for your use.

- Mag glass displays zoomed-up imagery around the cursor.

Protocol management and loading

- Customizable protocol-driven workflows for a wide range of clinical applications, including support for multiple review steps. You are guided through the review steps by a Navigation bar at the top of the User Interface.
 - Review protocols created by user with Volume Viewer can be saved as favorite protocols and accessed easily in routine.
 - Progressive Load enhances image loading performance. This feature has two modes:
 - Images are loaded in sequential mode (e.g., Reformat)
 - Images are loaded in interlaced mode (e.g., 3D/VR)
- In either case, image review can begin as soon as the first image is displayed.
- Navigation through series and exams can be performed without exiting the patient list.

Smart layout

- Volume Viewer General Review protocol presents an adaptive layout for single and multi-modality PACS-like reading of data. This protocol distributes and displays all loaded series of a single exam in an intelligent layout to maximize use of available monitor space. With multiple exams, the protocol automatically detects and distributes "like" series for optimal comparison.

GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

Advanced visualization capabilities

Volume Viewer offers various capabilities to display advanced rendering modes:

- 3D Volume Rendering / Volume Illumination¹ images
- MPR views with different slice thickness. The following rendering modes are available for all Thick Slab: MIP, Min IP, Volume Rendering, Average.
- Navigator views which display interactive endoluminal views exportable as movies.
- Curvilinear reformatting allowing display of curved, lumen and cross section views for various structures (vessels, spine, etc).
- Merge multiple volume rendered models into a single view.

Volume Rendering Volume Illumination

Predefined cut planes

Volume Viewer lets the user define cut planes to isolate specific structures in the VR model, in 16 configurations, including Left / Right / Inferior / Superior / Anterior / Posterior / Front. You can then display the voxels values on each plane.

Multiple Volume Rendering and Volume Illumination¹ objects

- Volume Viewer allows you to merge up to 8 Volume Rendered models from the same series into a single 3D view.
- You can adjust independently the threshold, colors, and transparency of each Volume Rendering or Volume Illumination

Auto-center fly through with smart cursor

- In fly-through studies (airways, colon, angiography) navigation along the centerline of the structure of interest can be recorded step by step and can be exported as a movie.
- In addition to the regular 180° or less Navigator view, Volume Viewer enables a Fish-Eye View with any wide camera angle value from 180° to 360°. This provides a view of structures both in front of and behind the users' virtual location on the same image.

Lumen view

- Lumen View provides an unfolded 3D view around a user defined centerline. The lumen view can be interactively adjusted (rotation around the centerline, width, field of view).

3D and 4D capabilities

- Volume Viewer allows 3D cursor synchronization in any orientation on the fly.

GE imagination at work

CONFIDENTIAL

SECRET DE AFACERI

- Real time reformat in oblique planes is available for your exams, as well as simple and double oblique interactive modes.
- Quick access to cross reference display gives you the ability to display the current slice location or all slices, or other planes.
- Real time interaction in 4D mode, by using Cine tool on all planes and 3D views for CT, MR, and PET multiphase data.

Compare mode


- With Volume Viewer several volumes can be loaded either from the same exam or from different exams. This is useful for exam comparison (Compare Mode) or for Multiphase examinations (of the liver for example). The Compare Mode is also available for PET/CT examinations and takes advantage of a dual screen configuration.
- Multi-phase CT and MR data are listed as separate phases in the Series Selection screen, allowing selection of all phases or a subset of phase data for review. The Series Selection Panel displays the multi-phase data split into arterial and portal phases, for example.
- Dynamic Load, compatible with both Volume Viewer and Fusion protocols, lets you drag & drop 3D volumes from CT, MR, PET, and 3D XA modalities into a desired viewport. Together with the Integrated Registration² option, Dynamic Load allows new volumes to be registered and loaded on the fly. Save State series can be restored as a separate session. For dual monitor configurations, a Save State series can be displayed on the right monitor together with a current session on the left monitor. For single monitor configurations, a user may simply switch between the two sessions.
- The zoom and pan functions are propagated to all images displayed in the same orientation.

2D/3D ROI for quantitative measurements


- Volume Viewer contains a set of basic 2D/3D tools: distance, angle, report cursor, arrow tool, annotation, free hand ROI tool, that can be placed on fused images as well.
- Easy deposition and labelling of measurements.
- All 3D ROIs can be customized, and color coded to display statistics computed on different intensity range and help assessing tissue classification based on voxel values.
- All measurements are considered bookmarks and can be accessed by navigating through the list of measurements.

AutoContour tool

- AutoContour provides consistent tools for 3D contouring ROI, with the same workflow used for CT, MR, and PET images.
- For CT and MR images, the contouring tool will suggest contour boundaries after defining 1 point or a diameter on a structure of interest, which can be edited prior to validation.

- For PET images, the contouring tool has been adapted for standard SUV calculations.

Advanced 3D tools

Segmentation tools

- AutoSelect segmentation tool helps adding or removing any continuous structures of interest by simple point and click.
- You can also access advanced segmentation tools such as Scalpel, Threshold, Dilate, Erode, Subtraction, Intersection, Addition, Filter Floaters.

Quick Vessel Trace (CT and MR images)

- You can perform vessel analysis with one or two points clicks from any CT or MR Angiography images.
- Quick Vessel Trace automatically extracts the vessel after user clicks and launches vessel tracking for fast review in curved reformat, cross-sectional, best L-section, lumen, and MRP view.

- Tracking can be extended proximally or distally for a full view of the vessel.

Summary table

- The Summary Table collects measurements and organizes them according to finding and exam date.
- The Summary Table allows interactive navigation through findings in the image viewports and can be docked or expanded.
- The Summary Table can be copied into the client clipboard and then pasted in email, reports, office software etc. (on AWS only)

Export capabilities

Volume Viewer contains multiple standard options for exporting the results of a review session:

- Save images to the database, as new DICOM series.
- Save the significant images as Key Image Note objects in the database. Key Image Notes and End Review allow you to flag images of interest as Key Image Notes (IHE profile) and push them to the archiving system when you exit the application.
- Save State objects are used to save the work in a new series containing all the post processing data (3D Model, displays, measurements, annotations, etc), for future review.
- Save STL, OBJ, VRML or 3MF file from 3D Volume Rendering / Volume Illumination¹ viewport for 3D printing purposes through 3D Suite³

Additional export tools are available to create specific sequences of images to be exported:

- The Batch tool allows creating a sequence of reformatted images or a sequence of rotating 3D views.
- The Movie tool creates a comprehensive movie including different rotations, zooms, and pan of the image, which can be exported as DICOM series or mpeg file.
- Quick Export: Exports in a single click a full batch of contiguous images at the displayed thickness for 2D images, or a batch of rotations of a 3D View.
- Cardiac Review and Export: Processing and reviewing cardiac exams for CT, MR and PET with manual oblique

reformatted protocols can be exported into a multi-phase Cine movie that allows the referring physician to review the exam in a dynamic mode.

- Capture and send images from the viewport to your computer clipboard on the fly (Ctrl+C/Ctrl+V on AW Server only).


Modality Features

CT image quality

- CT filters to enhance or smooth CT images to help improve image quality on the fly.
- The lung filter helps enhance contours of images reconstructed in standard mode for excellent visualization of lung structures

- A smoothing filter reduces noise while maintaining high image contrast.

CT Dual energy protocols

- Three protocols are available to enable review of dual energy images acquired on compatible GE Healthcare scanners and enhancing pixels with specific ratio: DE Calcium ($HU_{80}/HU_{140} \geq 1.25$), DE Uric Acid ($HU_{80}/HU_{140} \leq 1.25$) and DE Custom (custom threshold on HU_{Low}/HU_{High}).

GE imagination at work

- In addition, all Volume Viewer standard protocols are compatible with GSI[®] GE Healthcare CT images.

Bone VCAR⁵


- Bone VCAR is available in Volume Viewer toolbar and can be launched on the fly
- Vertebrae are automatically recognized and labelled (deep learning algorithm)
- A dedicated protocol provides a specific layout of the spine including curved views

Review of MR studies

Volume Viewer offers support for review of MR studies, through dedicated protocols.

- The MR-specific parameters (PSD name, image weighting, scan plane, etc) are recognized at loading and are used to customize the layout display.
- Anatomy-based protocols for standard review (e.g. spine, brain stroke...), with dedicated and guided review steps, to help review standard MR exams.
- Dedicated MR Cardiac Viewer for single and multiple Cine Review. Enables comparison between cardiac series such as Time Course and Myocardial Delayed Enhancement (MDE) images.
- Support of MR Multi Echo, Multiphase, Diffusion series.
- Support of parametric series (e.g., ADC, MTT, and other parametric maps) created in READY View, allowing direct functional measurements.


- 2D and 3D ROI propagation to other phases / series / exams allowing you to easily correlate information from multiple sources.
- Improved MR annotation consistency.
- In combination with Integrated Registration², enables direct access to MR image registration (inter/intra exam), image fusion and Whole Body MR Review protocols
- Ability to automatically bind several MR axial series corresponding to sequential axial locations into a single series. This is useful for consolidating multi-stage acquisitions for displaying Whole Body MR scanning.

SUV measurements on PET images

- Volume Viewer supports Standardized Uptake Values (SUV) for image display and measurements. Several SUV scales are available like SUVbw, SUV lbm, SUVbsa, as well as SUV Peak.
- Window/level presets may be user-defined for PET images based on SUV values.
- The SUV values are available in all the basic 2D/3D ROI tools of Volume Viewer and in the AutoContour tool.

PET/CT and PET/MR dedicated protocols

- Volume Viewer contains a list of predefined protocols allowing the review of PET/CT and PET/MR data. The screen layouts contain fused views between the morphological and functional images, as well as 3DMIP rendering of PET images.
- Factory protocols are customizable for your own review.
- PET 4D protocol allows users to load and display gated and dynamic PET series, as well as summing or reframing these series.

GE imagination at work

CONFIDENTIAL

SECRET DE AFACERI

PET IQ improvements

- Enhanced 3D visualization algorithm to fully support Q.Clear PET images.
- Q.Freeze⁶ option provides a respiratory motion-corrected Gated PET volume, statistically similar to conventional static PET with significantly reduced or eliminated blurring effects due to patient respiration. It provides also a corrected Gated PET series where counts are summed back to each respiratory phase.

Dedicated features for 3D XA images⁹

The optional Innova 3D XR application produces XA images (512x512 images or 256x256 images depending on your settings). This format opens the Volume Viewer world to Innova 3D XA images:

- Single exam: Innova Navigator, Curved Reformat, Multi Oblique, MPVR 3mm Average, Segment Structure
- Compare exam: Compare Axial, Compare Coronal, Compare Sagittal, Reformat
- Layouts and protocols customization
- The Innova options provide access to Workflow booster tools, easy multi-modality review "on-the-fly" dedicated Innova review protocols with dedicated Review Steps.

Smart compression

- The Smart Compression technology automatically displays images with full fidelity when the image is still, then uses the selected compression level for increased interaction speed during user interaction. This allows for diagnostic reads on full fidelity static images with responsive dynamic display even at low bandwidth. On-image visual indicators notify the user when compression is in effect.

Pre-processing⁷ features on AW Server

- Preprocessing is a Volume Viewer feature that automatically performs routine processing tasks and saves the results so they are ready for you. When new exams are transferred to the database, Volume Viewer recognizes user-defined keywords in the Series Description and launches the appropriate pre-processing. The results are stored in a Save State object with the original study so you can load it when you are ready for review.

AutoLaunch⁸ features on AW workstation

- AutoLaunch is a Volume Viewer feature that automatically preloads exams, so they are ready for your review.

- When you are processing a study and a new exam is transferred to the AW Workstation, this feature automatically launches the new exam with an appropriate Volume Viewer protocol in the background. When you are ready, one click in the AutoLaunch window instantly brings up the exam in the Volume Viewer application, eliminating waiting time and extra steps to load the exam into computer memory for processing.
- AutoLaunch is compatible with CT, MR and PET single volume protocols of Volume Viewer.
- When combined with optional applications⁶, this feature gives access to data already preprocessed, that you can review by just clicking on AutoLaunch window.

Summary of Operation

Volumetric models are loaded by selecting the exam or series. The user can select a protocol category from an anatomical selector or go directly to a Review Layout. In either case, images are loaded progressively in the background; this gives control to the user in just a few seconds after selecting the images. Selecting a Review Layout launches a volumetric display protocol with predefined layout preferences. Review Layouts may be customized and combined to suit the user's workflow. Selecting a protocol category unlocks a variety of visual protocols that include the layout, threshold, rendering mode and filming formats. Some of these protocols direct the user through the process providing capabilities to interactively view and manipulate the model, increasing productivity and consistency for all modalities.

References

¹Requires Volume Illumination license

Volume Illumination may not be available in some other countries or regions.

Please contact your sales associate

²Requires Integrated Registration license

³Requires 3D Suite license

⁴Requires Volume Viewer GSI license

⁵Requires Bone VCAR license

⁶Requires Q.Freeze license

⁷Preprocessing is compatible with the following purchasable options: AutoBone™ Xpress, CardIQ Xpress Reveal, CardIQ Xpress Function, CardIQ Xpress Process, Advanced CTC Pro 3D EC, Colon VCAR, and CT Perfusion 4D and FastStroke. ColonVCAR is not for sale in the US.

⁸Requires AutoLaunch license for AW workstation

⁹Requires Innova Volume Viewer license

GE Medical Systems SCS
283 rue de la Miniere
78530 Buc France

DOC1460639 rev9

General Electric reserves the right to make changes in specifications and features, or discontinue the product or service described at any time, without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. Illustrations are provided for informational purposes, and your configuration may differ.

This information does not constitute legal, financial, coding, or regulatory advice in connection with your use of the product or service. Please consult your professional advisors for any such advice. Operation of GE Healthcare products should neither circumvent nor take precedence over required patient care, including human intervention of healthcare providers. GE Healthcare products and services do not code medical procedures. Accurate coding is the responsibility of the provider or billing professional.

GE, the GE Monogram, imagination at work, Centricity and AutoBone are trademarks of General Electric company.

General Electric Company, by and through its GE Healthcare division.

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
General Electric Company, prin intermediul diviziei sale GE Healthcare.

©2019 General Electric Company

[Scrieți aici]

GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

CT Perfusion 4D

Comprehensive perfusion analysis available for neurology, oncology and cardiology. Easy-to-use, protocol-driven workflow for all organs.

Perfusion 4D has been expanded beyond stroke and oncology to now include dynamic myocardial perfusion. For any organ where blood perfusion is of interest Perfusion 4D is a complete package to quantitatively provide the necessary information to guide treatment decisions.

The protocol-driven design leads the user step-by-step through the process, reducing keystrokes and improving repeatability. Get the information you want quickly and reliably.

CONFIDENTIAL

SECRET DE AFACERI


Overview

GE's Perfusion 4D is a fast, easy-to-use automated software for analyzing CT Perfusion images related to stroke, tumor angiogenesis and dynamic myocardial perfusion. Its simple user interface and automated perfusion post-processing make it easy for you to diagnose quickly and accurately – and make treatment decisions more confidently.

GE leverages its deconvolution expertise, with its innovative delay correction algorithm. Perfusion 4D takes perfusion assessment to the next level, with the addition of CT Dynamic Myocardial Perfusion. Providing the classic color overlays with quantitative measurements to evaluate myocardial blood perfusion.

Highlights

- Improved IQ of functional maps in presence of noise
- Improved processing times
- Improved ROI mirroring for brain protocols
- Deep Learning algorithm for ventricle removal
- Automated midline detection for brain symmetry
- Absolute and relative values with customizable thresholds for tissue classification
- Permits injection rates of 4cc/sec
- Dynamic registration
- Dynamic myocardial perfusion maps

Visit us:

<https://www.gehealthcare.com/en/products/advanced-visualization/all-arterial-time-ct-perfusion-4d.html>

Visit us:

<https://www.gehealthcare.com/en/products/advanced-visualization/all-arterial-time-ct-perfusion-4d.html>

Features

- Visualize all information in true volumetric form with the ability to employ all the volumetric-based image analysis tools in the AW Volume Viewer.
- Deep Learning algorithm for brain ventricles removal¹
- Automated midline detection for brain symmetry
- Absolute and relative values with customizable thresholds for tissue classification¹
- Background processing will send processed images and functional maps in email² format to stroke team²
- Whole heart perfusion with semi-automatic guided workflow to aid in assessment of myocardial ischemia
- Smart Maps using the Intelligent 4D noise suppression algorithm improves functional map image quality in the presence of noise.
- Dynamic non-rigid registration protocols for dynamic acquisitions within the body.
- Productivity enhancements include:
 - Ability to select automated or manual vessel detection for Brain Perfusion
 - Faster processing times for functional map computations
 - Volumetric visualization of functional maps
 - Access from anywhere using AW Server.

- Time stamps on functional maps and Regions of interest (ROIs).
- Interactive large vessels exclusion.
- Many basic features, including the ability to:
 - Process double phase acquisition data.
 - Create 4D Regions of Interest (ROIs).
 - Display ROIs statistical information.
 - Display averaged time intensity information of the ROIs.
 - Simultaneously review any number of functional maps.
 - Customize protocols
 - Save parameters
- Includes CT perfusion protocols for stroke, brain & body tumor, myocardium as well as dynamic registration
- Functional maps:
 - Regional Blood Volume (BV)
 - Regional Blood Flow (BF)
 - Regional Mean Transit Time (MTT)
 - Contrast arrival delay (IRF T0)
 - Transit time to peak of impulse residue function (TMax)
 - Mean Slope of Increase
 - Capillary permeability surface area for lengthened acquisition protocols (PS)

- Average image
- Base image

Tissue Classification

System Requirements³

- AW Workstations AW 4.7 Ext. 14 or higher
- AW Server 3.2 Ext. 3.2 or higher

Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR).

This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

Indications for Use

CT Perfusion 4D is an image analysis software package that allows the user to produce dynamic image data and to generate information with regard to changes in image intensity over time. **It supports the analysis of CT Perfusion Images (in the head and body)** after the intravenous injection of contrast, in calculation of the various perfusion-related parameters (i.e. regional blood flow, regional blood volume, mean transit time and capillary permeability). The results are displayed in a user-friendly graphic format as parametric images.

This software will aid in the assessment of the extent and type of perfusion, blood volume and capillary permeability changes, which may be related to stroke or tumor angiogenesis and the treatment thereof.

CONFIDENTIAL

SECRET DE AFACERI

1. Requires CT Perfusion 4D Neuro Enhanced license. CT Perfusion 4D Neuro Enhanced may not be available in all countries or regions. Please contact your sales associate.
2. Requires Send by Email license and AW Server version with email service, not compatible with stand-alone AW systems. Send by Email may not be available in all countries or regions. Please contact your sales associate.
3. Send by email functionality requires an AW Server version with email service available.

GE imagination at work

CONFIDENTIAL

SECRET DE AFACERI

© 2021 General Electric Company.
All rights reserved. Data subject to change.
E Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
DOC1896217

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SECRET DE AFACERI

Visit us:

http://www3.gehealthcare.com/products/categories/advanced_visualization/applications/perfusion_4d_multi-organs

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SECRET DE AFACERI

GE imagination at work

VessellIQ™ Xpress with AutoBone™ Xpress

Fast and Easy to Use Vessel Analysis with Automated Real-Time Tracking, Auto-Labeling and Simplified User Interface.

Peripheral vascular disease (PVD) is a condition similar to coronary artery disease and carotid artery disease. In PVD, fatty deposits build up in the inner linings of the artery walls. These blockages restrict blood circulation, mainly in arteries leading to the kidneys, stomach, arms, legs and feet.

There are different clinical techniques used to diagnose PVD. Computed Tomography Angiography (CTA) can potentially provide sufficient information for therapeutic decision-making and reduced imaging costs for the management of PVD. One of the advantages of computed tomography (CT) is that it shows both vascular & nonvascular anatomy and also provides 3D information for best view. Also, CT may be more useful for patients with limited or no peripheral access or with stents, who may not be candidates for magnetic resonance angiography (MRA). Complementary to the ability to perform CTA is to perform clinical assessment by a dedicated post-processing tool.

Overview

VessellIQ™ Xpress with Autobone Xpress is a post-processing software package for the Advantage Workstation (AW) platform, AW Server, CT scanners and PACS reading stations. It is a tool to analyze CT angiographic data including stenosis analysis, thrombus, pre and post stent planning procedures, and directional vessel tortuosity visualization. Autobone Xpress provides zero-click bone segmentation for head, neck and other anatomy. The automation and ease of use streamlines workflow.

Highlights

- Automated detection of aorta and iliacs with auto-labeling of vessels
- Automated real-time fast tracking for all vasculature
- New contextual wheel menu including easy access to advanced options
- New centerline edit on capabilities
- Improved predefined measurements and edit capabilities
- New reporting interface which includes measurements in summary table
- New user interface

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/autobone-vessellq-xpress>

Features

- Anatomy based protocols for guided workflow in the analysis for your peripheral vascular CTA acquisitions.
- Automatic bone removal for efficient review of vasculature.
- Automated abdominal aorta CTA protocol with vessel tracking and vessel labeling.
- Fast tracking detects entire length of vessel and displays the centerline based on a single start point.
- Advanced editing tools bringing new flexibility for faster editing based on targeted anatomy.
- Thrombus detection with measurements and comparison over time.
- Dynamic AVA provides access to vessel tracking at any time, from any protocol, on any vessel.
- Advanced measurement tools for stenosis detection and color identification of HU value for quantification of different densities.
- Create a queue of cases that are pre-loaded and pre-processed in the background while you read another exam.
- Measurements are organized in the summary table allowing for easy viewing and exporting.
- Fully customizable workflows with adjustable layouts, personalized parameter and settings, custom review steps.

- Save State lets you save and restore the state of the processed images at any stage

System Requirements

- AW Workstation AW 4.7 ext 14 or higher
- AW Server 3.2 Ext. 3.2 or higher

Indications for Use

VessellQ Xpress is an optional, non-invasive, optimized, post-processing application intended to provide images and tools to analyze vascular anatomy and pathology, aiding physicians in diagnosis and determination of treatment paths, from a set of Computed Tomography (CT) Angiographic images. VessellQ Xpress is an option for the Advantage Workstation (AW) platform, CT Scanner, and/or PACS, which can be used in the analysis of 2D and 3D CT Angiography images/data for the purpose of cardiovascular and vascular disease assessment. This software-only device is designed to support physician assessment for a wide variety of clinical uses such as stenosis analysis, pre/post stent planning, pre/post valve replacement planning, and directional vessel tortuosity visualization. VessellQ Xpress' automatic visualization tools provide users the capability to segment bony structures for accurate identification of the vessels. Additional tools enable analysis of the vascular anatomy including the aorta, valves, and branching vessels for: anatomical sizing;

density and volume analysis of segmented vasculature and calcified / non-calcified plaque; and measurements of abnormalities. The TAVI Analysis option for VessellQ Xpress is a planning tool used for Trans Aortic Valve Implantation (TAVI) procedures. It automatically segments the aorta and displays the aortic valve in multiple views for measurements of anatomic structures commonly needed for aortic annulus replacement planning. TAVI Analysis provides guided workflow and semi-automated tools to aid in evaluation of appropriate access pathways for interventional procedure planning.

AutoBone Xpress option is a software package that is intended to facilitate segmentation of bony structures and calcifications from abdominal and extremity CT Angiography data.

Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, Imagination at work, AutoBone Xpress, VessellQ Xpress are
trademarks of General Electric Company.
All other third-party trademarks are the property of their respective owners
Doc1471990

Features Detail

Anatomy-based protocols

- The program offers you predefined protocols for fast, easy review of all vascular structures:
 - Intra-cerebral
 - Carotids
 - Thoracic Aorta
 - Abdominal Aorta
 - Runoff
 - Spinal
 - Extremity
- Protocols may be customized to meet your requirements and saved for quick access.

AutoBone™ Xpress

- Anatomical driven protocols with automatic bone removal.
- Improve vessel visualization by removing obstructive bony detail

- Automatic calcium removal for vessel lumen visualization

Auto Aorta

- With a single click, the program automatically removes bones and performs centerline tracking of the aorta and iliacs and labels vessels.

- Simply confirm the automated tracking or modify it by adding additional branch points, deleting a single point, or deleting multiple points from any view.

Fast Tracking

- Semi-automatic detection of complete vessel from a single starting point

- Advanced editing with variable interpolation tool for quick centerline edits adjusted for length of vessel.

- Bulls eye adjustment capability on x-section views to perfect tracking of the smallest most tortuous vessels.
- Volume measurements of the thrombus allow you to track change over time. Additional measurements include :
 - Minimum diameter
 - Maximum diameter
 - Mean diameter
 - Cross-sectional area

Vessel Label Data Base

- Each protocol is associated with a database of vessel names
- Customize vessel label database to meet your requirements for every anatomical region.
- A right mouse click provides access to the list of vessel names per anatomy.

Dynamic AVA

- Real-time tracking of the vessel center line and instant visualization of an unfolded view for quick vessel inspection.
- Immediate access to vessel analysis tools.
- Edit contours with a selection of drawing tools with variable interpolation for fast editing.

CONFIDENTIAL

SECRET DE
AFACERI

Advanced Vessel Analysis Tools

- Semi-automated vessel analysis lets you track, extract, visualize, and size vessels.
- Track a vessel and its branches for quick review in curved, oblique longitudinal or cross-sectional views.
- Multi-point deposition capabilities for improved center line tracking of challenging cases.
- When you edit a vessel contour, the program automatically interpolates the edit to neighboring slices.
- Generate and save rotating movies
 - Curved planar reformation and cross-sectional (lumen) views
 - Best L-section reformatted images
- Select vessel segments for quantitative or qualitative analysis.

- Available measurements include:
 - Distance
 - Cross-sectional area

- Single or dual reference point comparison
- Relative percent stenosis / aneurysm
- Color Identification
- 3D Angle
- Thrombus Extraction
- Volume
- Mean, min, max diameter
-
- Move freely between review steps and protocols without leaving the exam.

Stenosis Detection & ColorID

- With the customizable one- or two-point deposit technique, you can automatically quantify degree of stenosis and segment length.
- By assigning customizable color mapping to CT number ranges within an image, you can identify plaque on an axial image, curved reformat, or MPVR image. Up to ten levels of color coding are available. Colors transition smoothly from one density to another.
- Make volume and area measurements to track plaque size over time.

- Custom sizeable 3D ROI lets you apply color over all anatomy.
- Generate multiple ROIs on one vessel or multiple vessels.
- Automatically include each defined plaque, its volume and area into the report.

Quick AVA

- Dynamic AVA, provides access to vessel analysis from any protocol.
- Launch the Quick AVA feature from the main tool bar.
- Perform vessel analysis from any 3D or reformatted image with one- or two-point clicks.
- If you detect a lesion, you can analyze it by dropping a point above and below the section. The vessel analysis tool launches automatically for fast review in curved reformat, cross-sectional, lumen, and MRP view.
- Extend tracking proximally or distally for extended views at any time.

Summary Table & Report

- Measurements are automatically captured in summary table for easy review.

- Simply click on any measurement in the summary table and you will be taken to that measurement location.

CONFIDENTIAL

SECRET DE AFACERI

- Prepare report with a single click or generate an advanced report with customized options
- Send the report to database, filer or PDF.

Save State & Protocol

- Save the current state of tracking and measurements for easy reloading with a single mouse click.
Save state captures:
 - 3D images
 - Vessel tracking
 - Measurements
 - Vessel contours.
- Transfer the Save State to the PACS or other AW workstations for future review.
- Save custom protocols including anatomical labels and measurements.

Summary

VesselIQ Xpress provides user-friendly tools an intuitive interface to help streamline your vessel analysis workflow and help you generate more comprehensive reports for referring physicians.

Multiple Study Comparison

- With VesselIQ Xpress, compare a previous study with the current study. You also have the flexibility to launch comparison within Quick AVA.

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
SECRET DE AFACERI

Lung VCAR

For lung nodule analysis - automatic visualization, measurement, reporting and follow-up

Clinical and Technical Background

CT scanning is utilized as the exam of choice to visualize and analyze complex lung pathology. The detection of pulmonary nodules and assessment of their evolution with CT are of major importance in chest imaging.

Overview

Lung VCAR brings efficient CT pulmonary nodule assessment and diagnosis. The innovative Digital Contrast Agent (DCA) feature automatically visualizes lung nodules to help you confirm the presence or absence of suspicious lesions from 2 to 12 mm in size. Lung VCAR allows automated follow-up for lesion matching by the registration of two or more datasets, automatic lesion classification.

Highlights

- Synchronized 2D, Digital Contrast Agent (DCA) and segmentation analysis.
- Automatic nodule visualization.
- Automatic nodule analysis (volume, doubling time, % growth).
- Automatic follow-up.
- Reporting workflow.

CONFIDENTIAL

SECRET DE AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/aiw>

Features

Review

- Ability to synchronize multiple images for nodule comparison.
- Ability to review single or multiple exams and compare axial, sagittal, oblique, coronal, and volume-rendered images.
- Automatically propagates previous exam bookmarks to current exam.
- Automatically segments both right and left lung to reduce visual distraction.
- Digital Contrast Agent (DCA) automatically highlights spherical shapes to enhance visualization of suspicious nodules.

Analysis

- Performs automatic segmentation of all nodule types.
- Provides automatic nodule analysis, including % growth, doubling time
- Automatic Nodule Contour to verify pixels within the volume

Indications for Use

Lung VCAR/AdvantageALA is intended to provide an optimized non-invasive application to measure abnormalities in the lung (for example, nodules, lesions, etc.) from a set of computed tomography (CT) images.

The software is designed to support the physician in confirming the presence or absence of physician identified lung lesions (e.g. nodules). The software allows measurement of volume over time using a consistent standardized measurement protocol, thus providing an estimation of the volume doubling time. The Lung VCAR/AdvantageALA software allows analysis and displays statistics for nodule characterization all the different nodule types. Lung VCAR optional Digital Contrast Agent (DCA) module is an automated highlight feature for the visual identification of possible lesions. Digital Contrast Agent (DCA) is a 3D filter that produces images that highlight spherical (S) or cylindrical (C) anatomical regions, such as nodules, cysts, scars, and vessels. Images are made available to the physician to aid in characterization of suspicious nodules and thus, the patient management care decision process. Lung VCAR/AdvantageALA provides the physician with additional information, meant to complement diagnosis based on classical techniques.

System Requirements

Lung VCAR/AdvantageALA can be installed on AW Workstations AW 4.7 Ext. 14 or later and AW Server 3.2 Ext 3.2 or higher.

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
**SECRET DE
AFACERI**

GE imagination at work

© 2021 General Electric Company
All rights reserved. Data subject to change
GE and GE Monogram are trademarks of General Electric Company
* Trademark of General Electric Company
DOC0751472

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Colon VCAR EC

A comprehensive reading workflow solution for colonic lesion detection.

Colorectal cancer is a major health concern for both men and women. Recent studies conclude that early detection of colonic polyps can reduce the mortality rate of colon cancer. You need a protocol that lets you quickly and accurately detect and report results of CT colon scans; a program that optimizes your workflow and improves your reading productivity.

Overview

Colon VCAR EC gives you a comprehensive reading workflow solution for detecting colonic lesions with flexibility and efficient performance. With it you can visualize anatomy that would otherwise be obscured by tagged fluid. Plus digital contrast agent filter highlights potentially suspicious regions.

What's new

- The program allows you to read and problem solve in 2D, 3D, or 360 degree dissection views.
- Prone/Supine polyp linking provides quick navigation between datasets.
- Digital Contrast Agent (DCA) provides automatic visualization of shapes characteristic of polyps.
- Electronic Cleansing (EC) subtracts tagged stool and fluid, making it easier to identify lesions.

CONFIDENTIAL

SECRET DE AFACERI


Visit us:

www.gehealthcare.com/aw/applications/

Features

- Prone and supine synchronized reviewing for problem solving.
- You synchronize prone and supine navigation based on polyp bookmark location and distance from center line.
- Provides you with quick, easy to use polyp and center line measurement tools for complete reporting.
- You can leverage any combination of synchronized 2D, 3D, and 360 views.
- Dissection views make primary reading and problem solving more efficient.

- You can vary navigation speed.
- Lets you track entire colon from rectum to cecum for thorough, precise exams.
- You are provided with unrestricted viewing angles of all colon surfaces.
- You can apply color to bookmarked lesions for better visualization.
- The program allows you to assess luminal, mural and extracolonic anatomy and pelvic structures.

- A selection of application-specific tools ensures ease of use and thorough exams.
- Your workflow benefits from fast, interactive electronic patient reporting.

System Requirements

- AW Workstation
- AW Server

Indications for Use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. ColonVCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability.

It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.

Regulatory Compliance

This product complies with the European CE marking regulation following Medical Devices Directive: Directive 93/42/EEC.

This product has not been approved for sale in the United States by the United States Food and Drug Administration (FDA).

CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.

Features Detail

Electronic Cleansing


Electronic Cleansing, accessible in both 2D and 3D review, subtracts tagged stool and fluid from the generated data making it easier to identify lesions that may be visible in only one view.

Tagged material is identified based on Hounsfield units. Processing tasks such as Auto Segmentation, 3D flythrough, and Dissection also take advantage of tagged material subtraction.

Auto Segmentation shows the tagged material while 3D flythrough enables quick anatomy visualization without the tagged material. With this method, the tagged material is removed yet shown in a different color at an adjustable opacity level.

Similarly, the tagged material can be removed from the dissection view as seen in the images below.

User Interface

The Colon VCAR EC user interface has been optimized to facilitate fast, accurate examination of the colon's full extent. It provides seamless integration between automatic and manual interactions for a level of flexibility designed to reduce reading times.

Simultaneous, Indexed Viewing

The program can display both prone and supine

views, which can be manually synchronized for simultaneous and indexed viewing. During review, you can "Lock to Target" with a single button click to quickly problem solve during the primary read. And you are able to bookmark locations of suspected polyps.

360-degree Dissection

This powerful tool provides a 360-degree view of the entire internal large bowel simultaneously displayed in supine and prone positions. This ability speeds viewing and analysis.

This feature lets you leverage all available views: 2D, 3D, and 360-degree dissection during the primary read. For more traditional reading methods, a 2D, 3D, or Dissection view can be laid out and reserved for primary reading with reference to other views for problem solving.

Joystick Navigation

Colon VCAR EC lets you adjust the speed and direction of a flythrough simply by using the mouse. Flythroughs can also be set on autopilot for hands-off reviewing. Additionally, a single mouse click lets you lock onto a region of interest and fly around the target anomaly.

Computer Aided Reading

Colon VCAR EC's Contrast Agent (DCA), a 3D filter, highlights spherical anatomy such as polyps in color. The resulting images facilitate the characterization of potential polyps.


Virtual Biopsy

Problem solving while 3D reading is facilitated by the Virtual biopsy tool, which allows you to view areas in the navigator in color-coded or black/white scales. Different densities, determined by Hounsfield units, can be displayed using the color scale.

CONFIDENTIAL

SECRET DE AFACERI

Auto-Segmentation

Colon VCAR EC quickly performs the time-intensive task of colon segmentation for both prone and supine data sets. All you do is select the protocol.

The program automatically defines a luminal centerline in both prone and supine positions. This can be applied to a high percentage of cases, depending on the level of insufflation and collapse.

Semi-automatic and manual segmentation modes are also available. Semi-automatic segmentation provides complementary, and often necessary, assistance in cases when the colon is collapsed and operator intervention is required. This mode is typically used to connect colon segments around a collapse to form a single centerline. Manual segmentation is typically used to extend the centerline to the end of the cecum for completeness. Manual editing is also available to refine automatic tracking if needed.

Dissection (centerline) Contouring Verification

Using a color mapping technique, Dissection Contouring Verification outlines colon anatomy included in the centerline trace. You can review verification in any orthogonal plane.

Automatic High Resolution Flythrough

Real time 512-matrix resolution automatic flythrough enhances navigational image quality, and provides operational flexibility, functionality, and simplicity. Anatomy remains sharp with minimal pixilation. Steering remains centered as you navigate through the colon. You can specify your navigation speed, and 3D navigation is correlated with 2D and 360-degree Dissection views for rapid synchronized referencing in primary or secondary reading.

512 matrix resolution

Save State

This feature gives you the ability to save any centerline edits, bookmarks, and colored polyps for later reloading. It is also convenient for quick recall during consultation or exam review.

Prone/Supine Polyp Linking

After bookmarking a polyp in prone and supine datasets, the program provides you with the capability to match and label them for easy tracking. This gives you the ability to read registered prone/supine datasets for a one pass read.

Synchronized Prone/Supine 3D Dissection View Navigation

Once prone/supine polyp linkage is made, the program co-registers and links the datasets and polyps, letting you visualize opposing anatomy for a synchronized review.

Small Bowel Extraction

This feature quickly segments the small bowel, giving you an unobstructed view of the large intestine.

Analysis Mode

The program's Analysis Mode provides a variety of tools to help streamline your workflow:

- Lock Mode lets you target and possible lesion and stay locked on it at all times.
- Synchronized display can be customized to show two oblique images, axial, sagittal, coronal, volume rendered and/or 3D views within the colon in a Lock mode. Leveraging multiple synchronized planes and views in a problem-solving mode may enhance your reading productivity.

CONFIDENTIAL

SECRET DE AFACERI


GE imagination at work

Colon VCAR EC 3

Color Centerline Tracking

With Colon VCAR EC you can display a real-time 3D model of the colon with a defined centerline to quickly orient the location of the area you are viewing to the overall colon. This synchronization is always linked and available for quick reverence.

Patient Report*

Customizable, interactive reporting capability gives greater structure and flexibility in generating your reports. You can choose how many images to display per sheet, place your facility's logo, apply physician signatures, list patient demographics, and detail your comments and impressions. Export your report to CD, HTTP, or print it to paper.

The Quick Report panel gives you instant access to report input fields and provides access to automatic bookmark annotation on the colon map.

Measurement Tools


With Colon VCAR EC, you can measure polyp distance from distal rectum to polyp location in millimeters, perform 2D measurements, and obtain statistics of ROI in Hounsfield units.

Archive

You can define the number of images per polyp you want to archive to either a networked location or to magnetic optical disk.


Customizable Review and Analysis Layouts

Create layouts that match your reading style and save them as protocols for future recall. Viewports can be shaped into a variety of square or rectangular configurations.

Five of the most popular protocol layouts are pre-built for quick loading. Choose from several dissection configurations:

- 360-dissection gives you convenience and productivity benefits.
- 120-dissection lets you read around distortions if necessary.

- Use the 2D read configuration if the exam has poor insufflation and multiple collapses. 2D can also be used for extracolonic review.

- Use the 3-series layout to load a prone, supine, and decubitus series in parallel to obtain a correlated view.

- You can also customize new layouts to suit your own reading style.
- If you work on a one-monitor Advantage Workstation, a single monitor display is also available.

Movie Generator

You can generate movie clips of the area of interest simply by defining the start and end locations with a mouse click. A cross-sectional movie will be generated. You can also generate lock-to-target movies. Movies can be saved as mpeg, avi, HTTP and saved to CD.

GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

Colon VCAR EC 2

Entire Colon Movie

Generate a fly-through movie of the entire colon with one button click. Export the result as an mpeg or as screen captures to a PACS workstation for review or consultation.

Additional Features

- Automatic center line tracking
- SmartCursor™ for easier navigation
- Synchronized reformatted views
- Multiple correlations of anatomy with CT colonography
- Exploration of suspicious areas in primary 2D reformatted sections
- Primary 3D reformatted sections. Obliques, Navigator, Quick Virtual Dissection, and axial views are all correlated to the same anatomy of interest.
- Black/white or color display

Summary

Colon VCAR EC helps make colon evaluation a less traumatic procedure for patients. Preparation for the exam is less severe; the procedure is less time intensive and less expensive than optical colonoscopy. Moreover, the tools and features this program puts into your hands will streamline your workflow and potentially enable greater diagnostic confidence.

Indications for Use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. ColonVCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability. It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.

*Report feature is only available on AW Workstation

CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
*Trademark of General Electric Company.

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SmartScore 4.0

Advanced imaging software that detects, quantifies and scores cardiac calcium plaque burden. Instantly.

Cardiovascular disease remains one of the most common health issues in the world today. As with many conditions, early detection and patient risk assessment are vital to preventing or minimizing long-term negative effects. But many cases of this disease aren't diagnosed until the patient presents with symptoms. Conventional procedures to assess risk can be expensive, time consuming, and uncomfortable. An accessible, patient-friendly risk assessment method could help clinicians devise a regimen for their patients that might lessen the chance for serious cardiac events.

Overview

SmartScore 4.0 is designed to identify the presence of regional and global coronary artery calcification from a CT scan, then measure and score the results. Scores can be calculated using a standard Agatston/Janowitz (AJ) method. When correlated with a patient's personal information, the score can yield an estimation of a patient's risk for coronary artery disease.

What's new

- Non-invasive alternative to conventional assessment procedures.
- Score can be correlated to age group cohort to determine patient's risk per population
- Provides information on coronary artery wall calcium plaque buildup.
- Automatically detects calcium and highlights it in green.
- Streamlines workflow by networking patient demographics from CT scanner directly to SmartScore program.
- Free Hand Trace lets you outline specific ROIs.

CONFIDENTIAL

SECRET DE AFACERI

Visit us.

www.gehealthcare.com/aw/applications/smartscore-4/

Features

- AJ 130 scoring method uses conventional Agatston/Janowitz technique with a threshold of 130 HU which is adjusted to the appropriate image slice thickness.
- Volume scoring calculates volumes (mm^3) of calcified plaque above the 130 HU threshold.
- Mass scoring calculates mass (mg) of calcified plaque above the 130 HU threshold.
- Individual and aggregate scores are computed for each artery type.
- Provides two methods of calcium scoring.
- Customizable parameters can be predefined based on your site's preferences.
- User interface streamlines your workflow and reduces reading time.
- Report options allow you to custom tailor your reports and distribute them in a variety of formats.

System Requirements

- EKG monitor with recording device and x-ray translucent lead cable.
- Gantry hardware upgrade kit for those scanner systems already in operation.

Minimum platform release:

- Software for Advantage Workstation 4.2P or higher.
- AW Server 2.

Recommended Options

Postscript Printers:

- Codonics: NP-1660M
- Kodak: 3600 DMI
- Codonics 1660M, 1660MD or Horizon
- Lexmark Optra 1650N, 1855N, SC1275N, C710N, C72N, T612, or T614
- Seiko 1720D
- Kodak Dmi3600
- Quantium GL2101HD with film/thick paper
- Quantium GL2101HD with plain paper (see PI-008)
- Tally T8106
- HP LaserJet
- Xerox Phaser

Intended Use

SmartScore is a non-invasive software option that can be used to evaluate calcified plaques in the coronary arteries, which may be a risk factor for coronary artery disease. SmartScore may be used to monitor the progression/regression of calcium in coronary arteries over time, which may aid in the prognosis of cardiac disease.

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SECRET DE
AFACERI

GE imagination at work

© 2012 General Electric Company
All rights reserved. Data subject to change.

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CardIQ Xpress 2.0

Bringing integration and automation to your CT Cardiac post processing needs.

Cardiac disease is one of the leading health concerns worldwide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physician's approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives insight into cardiac anatomy and tools for thorough evaluation into the extent of disease.

Overview

CardIQ Xpress 2.0 is an integrated post processing imaging analysis application dedicated to cardiovascular imaging on GE Advantage Workstation (AW) and AW Server Enterprise System (AWS). The CardIQ Xpress 2.0 software option can be used to display, reformat, and analyze 2D or 3D cardiac CT images for qualitative or quantitative assessment of heart anatomy and coronary artery vessels from a single or multiple cardiac phase image data sets. Cardiac motion is a very real challenge that can occur at any heart rate. CardIQ Xpress 2.0 is designed to work with SnapShot Freeze* images to automatically process and display images generated with reduced motion blur artifact.

Highlights

- Automatically segment coronary tree across phases.
- Automatically tracks and labels coronary arteries.
- Improved centerline editing tools for faster edits.
- Right mouse wheel menu for quick access to renaming, deleting, and editing centerlines.
- Plaque ID tool assists in visualizing and quantifying plaque burden.
- Relative perfusion highlights and quantifies hypo-dense areas of myocardium.

CONFIDENTIAL

SECRET DE AFACERI

Visit us: <https://www.gehealthcare.com/products/advanced-visualization/all-applications/heart/qxpress-2.0.html>

Features

- Pre-processing automatically recognizes cardiac datasets and performs all segmentations as data arrives on the system saving processing time
- Auto Coronary analysis automatically segments the coronary tree, tracks and labels the coronary arteries
- Three pre-defined orientations provide quick access to planes for best visualization of the coronaries in 2D
- Quick AVA allows access to vessel tracking at any time
- Relative perfusion color codes and quantifies percent of hypo-dense areas of myocardium with four selectable color maps and fusion overlays with the coronary tree
- Plaque ID provides volume measurements for four distinct Hounsfield ranges to aid with identification and visualization of coronary plaque
- Automatically display SnapShot Freeze¹ processed images for reduced motion blur
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes
- Stenosis measurements, IVUS views and preset volume rendering models assist in communication of specific findings
- Measure ES and ED for ejection fraction & volume with automatic extraction of the left ventricle

- Create multiphase beating hearts
- Select oblique reformat views in the standard cath angles for easy analysis of the coronary vessels
- Display 4D valve views with a single click

System Requirements

- AW Workstation AW 4.7 Ext 14 or higher
- AW Server 3.2 Ext. 2.0 or higher

Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions.

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Rx Only

Indications for Use

CardIQ Xpress 2.0 is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images. CardIQ Xpress 2.0 is a CT, non-invasive, image analysis software package, which aids in diagnosing of cardiovascular disease to include, coronary artery disease, functional parameters of the heart,

heart structures and follow-up for stent placement, bypasses and plaque imaging. CardIQ Xpress 2.0 offers unique tools such as automatic tracking, which will pre-process the CT data into multiple viewing ports to allow for an expedited read time improving workflow. With CardIQ Xpress 2.0, the user can color code the myocardial tissue to show hypo/hyper-dense areas in the myocardial tissue of the heart. With the IVUS-like view the user can color code the HU units of the plaque to better visualize the difference between calcified and non-calcified plaque in the wall of the vessel and the lumen to determine the amount of atherosclerosis. The user can see the different valve planes along with a variety of new layouts to align the heart. The IVUS-like view is created by applying GE's Volume Rendering on a cross-section perpendicular to the detected centerline. This view merely displays a cross section as in IVUS imaging and color codes like IVUS images. No new or additional diagnostic information is added. CardIQ Xpress 2.0 is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
DOC1096477

Features Detail

Pre-Processing

- Automated processing of single or multiphase data
- Customizable menu to choose which protocols are pre-processed
- Pre-processed protocols include: Auto Coronary Analysis, Labeling of vessels, VR heart, Angiographic View, and Auto EF

Auto Launch

- Multiple exams can be loaded into the auto launch to ready them for review (AW only)
- Color-code display of ready-to-read exam listing (AW only)
- The ability to switch between exams without having to quit out of the application and reload

Automated Coronary Vessel Analysis (CVA)

- CVA allows the users to track, extract, visualize and measure coronary arteries from either single or multiple cardiac phase images
- Auto launch case selection and display of automatically tracked and labeled coronary vessels
- Simple, corrections to vessel branch tracking
- Curved, oblique, longitudinal, and cross-sectional reformatted views are automatically generated in the software for quick review

- Current-state tracking points within vessel analysis may be saved for future review and/or manipulation
- Color-coded plaque analysis with volume measurements for evaluating the change in plaque size over time
- Generate and save rotational movies from curved planar reformation, best section and cross-sectional (lumen) views

- Quantitative or qualitative coronary vessel assessment on user-selected vessel segments.
- Vessel measurements including:
 - Distance and volume
 - Cross-sectional area and mean diameter
 - Single or Dual reference point comparison
 - Relative percent stenosis
- Single-screen filming capabilities with multi-views within the screen to show the entire picture of the vessel
- Automatically display of the coronary vessel tree
- Using pre-set protocols, 2D or 3D coronary vessel tree models can be generated and displayed in an automatic fashion to give users a qualitative overview of coronary vessel structures

- New right mouse wheel menu for quick access to renaming, deleting and editing centerlines
- Improved centerline editing tools for faster edits
- The vessel tree models can be in the form of either 3D volume rendered or Maximum Intensity Projection (MIP)

2D Review

- Three pre-defined orientation protocols for easy review of the coronary vessels and chambers
- Dual reformat review allows automatic review of the coronaries from an axial image while linking to oblique views in longitudinal and cross-sectional planes
- Reformat review predefines workflow steps to automatically present thick, multi-planer views of cardiac anatomy
- One-click access to quick vessel analysis

GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

Quick/Dynamic AVA

- Single or dual click vessel analysis from any protocol; 3D, reformat or vessel tree
- Real-time tracking of the center line with instant visualization of a vessel for quick inspection
- Ability to extend tracking proximally or distally for full view of the vessel
- Ability to select a point above and below potential lesions for automated vessel analysis

Vessel Data Base

- A vessel name data base for quick access to vessel labels
- Anatomically relevant listings
- Ability to add new vessel names into the data base
- Automated processing of single or multiphase data

Relative Perfusion

- Hypo-dense areas of myocardium can be high-lighted with color to aid in the visualization of ischemic heart disease
- Pre-defined layouts for density enhancement
- A hybrid display view to show the relationship of the vessel tree with the perfusion defect
- Quantification of a hypo-perfused area as it relates to myocardial defects
- Four selectable color maps to help display the hypo-dense areas
- Exportable statistics via one touch

IVUS Views

- Interactive volume rendered images to better define the calcium, lumen and non-calcified plaque related to lumen narrowing or wall abnormalities
- Can be applied to any 2D MIP image to include; cross-section, MPR, best L-section images

- Display cross-sections perpendicular to the vessel centerline to create an IVUS-like image. No new diagnostic information is added

PlaqID

- Customizable color mapping to HU ranges for easy identification of plaque in axial, curved reformat or MPVR images
- Four distinct colors to aid in differentiating vessel lumen, non-calcified plaque and calcified plaque
- Easy modification of the color and transparency of plaque
- Volume, area to track size of plaque over time
- Smooth transition of color from one density to another
- Plaque volume and area automatically captured in summary table.

VR Heart

- One-touch automatic segmentation algorithm designed to extract the cardiac anatomy from within the chest image
- Optimize algorithms for the easy review of heart and bypass grafts
- Multiple optimized VR curve settings to enhance different structures within the heart.
- Automatic segmentation of both single or multiphase cardiac image data sets
- 4D beating heart images with the ability to page through the phases, rotate the image with real-time functionality

Multi-Phase Image Review

- Multiphase images can be reviewed in any protocol
- 10 phase images show the heart throughout the complete cycle
- Quickly edit phase images to keep only the phase or phases needed for analysis of the coronary vessels once multi-phase review is complete

GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

Transparency View

- 3D vessel tree overlaid on the heart chambers with the ability to adjust opacity of chambers of heart, myocardium vessels tree and bones

3D Angiographic View

- One-touch automatic visualization of the coronary arteries in an x-ray angiographic view
- Single click cardiac cath views to orientate the angiographic view into RAO/LAO, or CRANIAL/CAUDAL orientation
- Added or removed vessels are easily archived with the Auto Select feature
- Toggle between gray scale or inverse video views

One Touch Cath Views

- Default cardiac cath views are provided with CardIQ Xpress 2.0 software to provide 3D VR heart, 3D vessel tree or reformatted image in the orientation of a standard cardiac catheterization
- Single click on the angulation area to change between one-touch views
- Interactive RAO/LAO & CRANIAL/CAUDAL orientations allow for manually positioning the view
- Customizable orientations can be easily generated and saved for future use

3D Ejection Fraction

- Automatic extraction of the left ventricle across all phases and the automatic detection of systole and diastole delivers calculated ejection fraction and stroke volumes
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes

Layouts

- Configurable layouts
- Dual monitor support
- Configurable default layout for use at AVA launch
- Easily capture the screen layout for identical visualization on PACs or film

4D Movie Mode

- Movie mode allows the user to view and save multiple views of the heart and display as a movie sequence
- Load multiple phases from systole to diastole to create a beating heart movie
- Save movies as a DICOM image set or export in JPEG/MPEG format

SnapShot Freeze (SSF)

- Automatically display motion blur-reduced SnapShot Freeze processed images
- Reprocess SnapShot Freeze images for artifacts caused by motion blurring after manual vessel editing
- Motion blurring correction requires the use images produced by a CT using the optional SnapShot Freeze feature

Protocols

- CardIQ Xpress 2.0 is supplied with a set of pre-defined protocols that are easily adaptable to customized protocols
- The protocols have comprehensive instructions and tools which minimize the need to memorize the procedures or to refer continuously to user documentation.

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GE imagination at work

CardIQ Function Xpress

Real time cardiac review the instant you're ready to read.

Cardiac disease is one of the leading health concerns world wide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physicians approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives you insight into cardiac anatomy and function supported by quantifiable data.

Overview

Cardiac IQ Function Xpress post-processing software helps you evaluate cardiac function and diagnose cardiovascular disease with a high degree of confidence. Providing accurate and reproducible quantification of left and right ventricular volumes, ejection fractions, and myocardial mass, Cardiac IQ Function Xpress is optimized to assess cardiac function using multiphase, multi-slice cardiac CT images. The CardIQ Function Xpress option provides an easy-to-use and time-effective way for you to perform cardiovascular functional analysis.

CONFIDENTIAL

SECRET DE AFACERI

What's new

- Automatically pre-processes and loads exams.
- Automatically detects all chambers in all phases.
- Automatically performs ejection fraction and chamber volume analysis with 91% reliability on LV and RV automatic ejection fraction.
- Provides single click myocardial analysis with bull's eye images.
- Performing simultaneous case review in real time.
- Lets you load multiphase data with no load time.
- Lets you review multiple cases in real time.
- Lets you generate reports with imported analytics.

Visit us: www.gehealthcare.com/ow/applications/cardiq-function-xpress/

Features

- Apply the CardIQ Function Xpress application to standard axial CT images acquired on a GE scanner using the cardiac CT SnapShot imaging acquisition option.
- Automatically selects each heart chamber for individual volume analysis.
- Performs behind-the-scene function data process for real time ejection fraction review, volume analysis, and myocardial analysis.
- Enables you to extract, render, and display 3D volumetric models of the endocardium for ejection fraction calculation.
- Automatic detection of the endocardial and epicardial walls for wall motion, wall thickness, wall thickening and myocardial mass analysis. Lets you perform volume analysis of all chambers.
- You can automatically calculate left atrium volume while excluding the pulmonary vein.
- With a single click you can activate visual wall motion with short axis images in basal, mid, and distal orientations along with a two-chamber long axis view.
- You can perform myocardial analysis with wall motion, wall thickness, and mass calculations.
- The application's flexible reporting tool lets you include graphical representations.
- Display on screen a table of key functional parameters for instant visualization.

Intended Use

CardIQ Function Xpress is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

Image Requirements

CardIQ Function Xpress accepts standard CT image sets acquired on qualified GE CT scanners using the appropriate cardiac imaging software. These images must meet the same image requirements as those for the basic Volume Share 2 application.

System Requirements

- AW Volume Share Workstation with Volume Share 2 running 4 8400 hardware or higher (running 64 bit).
- Auto launch and preprocessing are available only on 8400 workstations with 16 GB RAM.
- 2 monitor configuration.
- Color landscape monitor.

Recommended Printers

- Codonics NP-1660M
- Kodak 3600 DMI
- Codonics 1660M, 1660MD, or Horizon
- Lexmark Optra 1650N, 1855N, SC1275N, C710N, C72N, T612, or T614
- Seiko 1720D
- Quantum GL2101HD with film/thick paper
- Quantum GL2101HD with plain paper (see PI-008)
- Tally T8106
- HP LaserJet
- Xerox Phaser

Regulatory Compliance

This product complies with the European CE Marking regulation for Medical Devices Directive: Directive 93/42/EEC, dated 14 June 1993.

CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2012 General Electric Company
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.


Overview

CardIQ Function Xpress allows you to accurately and reproducibly quantify myocardial mass, and left and right ventricular volumes and ejection fractions. The software automatically detects left ventricular endocardial and epicardial contours so that you can more easily assess left ventricular (LV) and right ventricular (RV) functional parameters. The program is optimized to perform assessment of cardiac function using GE multiphase, multi-slice cardiac CT images.

The CardIQ Function Xpress menu page consists of multiple image processing tools for single or multiphase cardiac studies. These protocols can be accessed either from the AW Volume Viewer application or directly from the application menu page on the browser.

CardIQ Function Xpress benefits to cardiac patients include:

- Minimally invasive alternative to help diagnose heart function.
- Data is acquired during a coronary CTA exam, so no additional dose or contrast is required.
- Offers a diagnostic option to patient who may not be able to undergo other non-invasive exams.

Feature Details

CardIQ Function Xpress uses single or multiphase image data sets to measure, graphically plot, and tabulate a number of ventricular function and wall motion parameters.

Pre-processing

- The program automatically processes multiphase data for functional analysis.
- You can select configurations to process only the chambers you require.
- For more streamlined reading, ejection fraction, volume analysis, and myocardial analysis are completed before you view the case.

Auto Launch

- Load multiple cases into auto launch so they are ready to review when you are.
- Switch between exams without having to quit the application and reload.
- Exams ready for review display in green.

Automatic Chamber Detection

- The program automatically loads single or multiphase data sets and selects the heart's chambers for volume analysis.
- Analyze individual chambers by automatically contouring any of the heart's four chambers.
- The system automatically detects the left and right ventricles and the left atrium and contours the endocardium for ejection fraction calculation.
- Choose single or multi-phase chamber detection.

Automatic Ejection Fraction

- The system features a protocol that, when chosen, automatically detects the end diastolic and end systolic phases and calculates an accurate ejection fraction of the left ventricle, right ventricle, and left atrium.

- Ejection fraction measurements are available for

can include logos, diagrams, signatures, patient history, and messages.

- Generate reports in PDF, HTML, or plain text.
- You can combine CardIQ Function Xpress reports with CardIQ Xpress 2.0 Elite reports.

Simplified Workflow

CardIQ Function Xpress simplifies your workflow with:

- Automatic ejection fraction protocol
- Streamlined multiphase beating images for easy visualization of cardiac wall motion.
- Single click myocardial analysis with bull's eye images.
- Simplified reporting.

Summary

CardIQ Function Xpress provides you with the tool you need to accurately analyze, evaluate, and report on cardiac anatomy and function. Its ease of use will streamline your workflow and its accuracy will give you greater diagnostic confidence.

Indication for Use

CardIQ Function Xpress is used in conjunction with CT cardiac images to automatically calculate and display various left ventricular and right ventricular functional parameters such as ejection fraction, end systolic and end diastolic volumes, stroke volumes, wall motion, wall thickening, cardiac output, myocardial mass, and systemic and pulmonary vascular resistance. Volume measurement of each chamber of the heart is also available. With CardIQ Function Xpress atrium volumes may be used to determine volume assessment of atrial disease to including, but not limited to, atrial fibrillation. CardIQ Function Xpress is a CT, non-invasive image analysis software package, which aids in the assessment of cardiac function and in determination of cardiovascular disease diagnosis and management. CardIQ Function Xpress is available on the AW Workstation, CT Scanner, PAC, or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

CardIQ Function Xpress is an extension of Volume Viewer 4 and is dedicated to cardiovascular imaging. It provides an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of CT Angiographic images.

CONFIDENTIAL

SECRET DE
AFACERI


GE imagination at work

© 2012 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

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A. Fodor

CardEP

Streamlined post-processing for enhanced electrophysiology procedures.

CT cardiac exams provide critical information to practitioners of cardiac therapy. Precise information pertaining to the left atrium's complex anatomy, the pulmonary veins, the coronary sinus, or the cardiac veins has a major impact on the efficacy of subsequent cardiac therapy: It can speed procedures and facilitate treatment. But analyzing and reporting the wealth of information CT cardiac studies provides can be time consuming. A program that automates many routine functions and gives you tools to easily quantify and qualify many aspects of cardiac function would streamline your workflow and give you greater diagnostic confidence.

Overview

CardEP is an integrated post processing image analysis software for the application of cardiovascular and electrophysiology imaging on the AW Workstation (AW) and AW Server platforms. With it, you can process display, reformat, and analyze 2D or 3D cardiac images for qualitative or quantitative assessment of heart anatomy and pulmonary veins from single or multi-phase cardiac image data sets.

Highlights

- Visualize the anatomical layout of the left atrium and pulmonary vessels.
- Visualize the origins of the pulmonary vessels.
- Automatic heart segmentation.
- Helps you visualize and quantify cardiac venous pathways.
- Offers you a variety of 2D, 3D, or reformatted protocols with which to perform image analysis.

CONFIDENTIAL

SECRET DE AFACERI

Features

CardEP offers a variety of features that give you image analysis flexibility:

- Volume render 3D heart and chamber views
- Obtain segmented 3D views of the left atrium and pulmonary veins
- View hollow 3D VR heart models
- Measure the pulmonary veins and left atrial appendage with predefined protocols
- Render and display 2D/3D views of the left chambers in single or multiple phases of the cardiac cycle with one-touch automatic extraction
- Reformat standard axial CT images of single or multiple cardiac phases into short or long axis and save the image series for further analysis.
- Automatically track, extract, and display pulmonary veins with the pulmonary vessel analysis tool
- Generate pulmonary vein measurements
- Automatically track, extract, and display coronary sinus and cardiac veins with the coronary sinus analysis tool
- Perform one-touch standard 2D reformation and 3D volume rendering, including hollow view, of cardiac anatomy with optimized Volume Viewer and 3D rendering tools

- Use predefined navigator views of the pulmonary veins to get a fly-through perspective of the vessels.


Image Requirements

CardEP accepts standard cardiac gated CT image sets acquired on CT scanners equipped with DICOM 3.0 standards. Images must meet the same image requirements as those for the basic Volume Viewer application.

System Requirements

The CardEP option can be installed on Advantage Workstation (AW) or AW Server with VolumeShare 7 or later software.

Compatible hardware:

- Z800, Z820, Z440 (and later)
- AW Server 3.2 and above

(recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048))

Indications for Use

CardEP is a post-processing software option for the Advantage Workstation (AW) Platform. This product can be used for the analysis of CT angiographic images for the assessment of the heart to include the atria, pulmonary veins and coronary sinus. It provides quantitative analysis tools which include a number of display, measurement and model export capabilities. This product can be used to aid trained physicians in the visualization and assessment of cardiac anatomy.

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SECRET DE AFACERI

GE imagination at work

© 2016 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company
DOC1117327

Feature Details

Operation

- On the menu page, multiple image processing tools for both single and multiple cardiac phase studies are provided. Access each pre-set protocol from the AW Volume Viewer application or from the application pull-down menu page on the browser.

- CardEP images present you with a detailed display of patient heart rate, cardiac scanning parameters, and each scan's temporal resolution.

3D VR Atrium

- Pre-set protocols enable you to semi-automatically generate and display 2D or 3D left atrial/ventricular views for a qualitative overview of the heart's chambers.

- The atrial/ventricular model can be either 3D volume rendered or Maximum Intensity Projection (MIP).

- With the endocardial view feature, you can view inside the atrium to analyze the relationship of the pulmonary vein ostia from the chamber's endocardial surface.
- With a one-touch protocol you can quickly change from a 3D model to an endoscopic view of the atrium and pulmonary veins.
- Anatomical landmark selection lets you deposit points on the 3D model, embed them, and then export the model for future registration use with a mapping system.
- Save predefined landmarks and models and export them to CD.
- Save current model state to the AW browser and reload for future analysis.
- Optimized VR and curve setting provides clearer contrast enhanced cardiac display.
- Choose from single or multi cardiac phase sets.
- Add/remove structure (Autoselect), scalping and brush threshold tools are all available for quick optimization of a 3D model.

3D VR Heart

- Automatically extracts cardiac anatomy from within the chest cavity.

- Optimized VR and curve setting

provides clearer contrast enhanced cardiac display.

- Automatic segmentation for single or multi-phase cardiac data.
- The program supports multiphase thresholding, filtering, and scalping techniques.
- Create movies of rotating and beating hearts with easy export capabilities.
- Export enhanced 3D models, including anatomical markers.

Image Reformation

- With batch processing tools, you can reformat either single or multi-phase cardiac images in various non-axial orientations, e.g. short axis, long axis, chamber views. Save reformatted images in separate image series for further analysis (e.g. from PACs).

- Offers predefined protocols to generate automatically short and long axis views.
- Multi-phase image capability lets you quickly and efficiently assess phases or function.
- View the cross-sectional diameter of the pulmonary vein and left atrial appendage with specifically designed protocols.
- Measurement tools integrated

CONFIDENTIAL

SECRET DE AFACERI

with protocols enable you to quickly ascertain length and cross-sectional diameter measurements.

Coronary Sinus and Pulmonary Vein Analysis

- Load single or multi-phase cardiac data sets in the vessel analysis tool.

- With the semi-automated pulmonary vein/coronary sinus analysis tool (CSA), based on the Advanced Vessel Analysis (AVA) technique, you can track, extract, visualize, and measure the veins and sinus (physical dimension and stenosis sizing) from either single or multiple cardiac phase image data sets.
- Generate curved, oblique, longitudinal, and cross-sectional reformatted views of individual vessels.
- Generate 3D models of the atrium and pulmonary veins during your vessel analysis procedure.
- Save current state tracking points within your vessel analysis for future image review or manipulation.
- Generate and save rotational movies from curved planar reformation and cross-sectional (lumen) views.
- Generate and save movies from Best L Section reformatted images.
- Perform quantitative and qualitative coronary sinus assessments on segment you

choose:

- Distance measurement
 - Cross-sectional area
 - Single or dual reference point comparison
 - Relative percentage stenosis
 - Volume measurement
 - Mean diameter
- Generate beating, rotating movies for visualization and communicating findings.

Filming

CardEP offers you filming flexibility. Batch filming capability is provided, and any image can be queued for filming.

One Touch EP Views

- Default views provide you with visualization of the atria and pulmonary veins in correct anatomical orientation. You are provided the capability to create your own presets to save and apply to future analyses.
- Use the interactive RAO/LAO and CRA/CAU annotation feature to manually orient the view.

Pulmonary Vein (PV) and Left Atrial Appendage Reformation

- Predefined protocols allow you to generate reformation of the pulmonary veins and left atrial appendage in the long axis.
- Generate cross-sectional views of the veins or appendage to determine the ostium's diameter.

EP Navigator Views

- With one-touch, the navigator displays endoscopic views of the left atrium.
- View the pulmonary vein and atrial appendage ostia from either the left or right side using

a predefined orientation layout preset.

- Optimize interior atrium views with automated threshold settings.
- Global virtual endoscopic view provides you ostial displays of the right and left pulmonary veins and the left atrial appendage.

Batch Movie Mode

- Batch movie mode lets you define first and last images for movie creation.
- Create rotating and beating movies
- Save automatically generated movies to the viewer as DICOM image sets or send them to data export in JPEG/MPEG/AVI formats.
- Use the movie mode with MIP, 3D VR, and reformatted images.
- Generate and save movies from Best L Section Reformatted images, curved planar reformatted, cross-sectional, and lumen views.

Protocols

CardEP contains predefined protocols that you can edit or modify to adapt them to your specific customized protocols and save under a new name. The CardEP protocols comprise comprehensive instructions and the tools you need to perform analyses. Tool icons are displayed in the protocol panels for fast and easy access.

GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

Summary

CardEP gives you the tools and flexibility to quickly analyze cardiac studies. The information and statistics you can generate with this program may help referring physicians with their assessment of therapy options for their patients in the visualization and assessment of cardiac anatomy.

Indications for Use

CardEP is a post-processing software option for the Advantage Workstation (AW) Platform. This product can be used for the analysis of CT angiographic images for the assessment of the heart to include the atria, pulmonary veins and coronary sinus. It provides quantitative analysis tools which include a number of display, measurement and model export capabilities. This product can be used to aid trained physicians in the visualization and assessment of cardiac anatomy.

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SECRET DE AFACERI

GE imagination at work

© 2016 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
0001117377

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