ARCHITECT**SYSTEM****ro**

HIV Ag/Ab Combo

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PROCEDURĂ

Materiale furnizate:

- 4J27 ARCHITECT HIV Ag/Ab Combo Reagent Kit

Materiale necesare dar nefurnizate:

- ARCHITECT i System
- Fișierul de testare ARCHITECT HIV Ag/Ab Combo, care poate fi obținut de pe:
 - ARCHITECT i System e-Assay CD-ROM care poate fi găsit pe www.abbottdiagnostics.com
 - ARCHITECT i System Assay CD-ROM
- 4J27-03 ARCHITECT HIV Ag/Ab Combo Calibrator
- 4J27-12 ARCHITECT HIV Ag/Ab Combo Controls
- ARCHITECT i **PRE-TRIGGER SOLUTION**
- ARCHITECT i **TRIGGER SOLUTION**
- ARCHITECT i **WASH BUFFER**
- ARCHITECT i **REACTION VESSELS**
- ARCHITECT i **SAMPLE CUPS**
- ARCHITECT i **SEPTUM**
- ARCHITECT i **REPLACEMENT CAPS**
- Pipete sau vârfuri de pipetă (opțional) pentru a dispensa volumele specificate în fereastra „pacient” sau „control” de pe ecran.

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BIBLIOGRAFIE

1. Barré-Sinoussi F, Chermann JC, Rey F, *et al.* Isolation of a T-lymphotropic retrovirus from a patient at risk for acquired immune deficiency syndrome (AIDS). *Science* 1983;220:868-71.
2. Popovic M, Sarngadharan MG, Read E, *et al.* Detection, isolation, and continuous production of cytopathic retroviruses (HTLV-III) from patients with AIDS and pre-AIDS. *Science* 1984;224:497-500.
3. Gallo RC, Salahuddin SZ, Popovic M, *et al.* Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* 1984;224:500-3.
4. Leitner T. Genetic subtypes of HIV-1. In: Myers G, Korber BT, Foley BT *et al.*, editors. *Human Retroviruses and AIDS 1996*. Los Alamos, NM: Los Alamos National Laboratory, 1996:III-28-III-40. (Available on-line at <http://hiv-web.lanl.gov>)
5. Robertson DL, Anderson JP, Bradac JA, *et al.* HIV-1 nomenclature proposal: a reference guide to HIV-1 classification. In: Kuiken CL, Foley B, Hahn B, *et al.*, editors. *Human Retroviruses and AIDS 1999*. Los Alamos, NM: Los Alamos National Laboratory; 1999:492-505. (Available on-line at <http://hiv-web.lanl.gov>)
6. Clavel F, Guétard D, Brun-Vézinet F *et al.* Isolation of a new human retrovirus from West African patients with AIDS. *Science* 1986;233:343-6.
7. Clavel F. HIV-2, the West African AIDS virus. *AIDS* 1987;1:135-40.
8. Schochetman G, George JR, editors. *AIDS Testing: A Comprehensive Guide to Technical, Medical, Social, Legal, and Management Issues*. 2nd ed. New York: NY Springer-Verlag; 1994.
9. Sarngadharan MG, Popovic M, Bruch L, *et al.* Antibodies reactive with human T-lymphotropic retroviruses (HTLV-III) in the serum of patients with AIDS. *Science* 1984;224:506-8.
10. Jackson JB, Kwok SY, Sninsky JJ, *et al.* Human immunodeficiency virus type 1 detected in all seropositive symptomatic and asymptomatic individuals. *J Clin Microbiol* 1990;28:16-9.
11. Simon F, Maucière P, Roques P, *et al.* Identification of a new human immunodeficiency virus type 1 distinct from group M and group O. *Nature Medicine* 1998;4:1032-7.
12. Ayoub A, Souquière S, Njinku B, *et al.* HIV-1 group N among HIV-1-seropositive individuals in Cameroon. *AIDS* 2000;14:2623-5.
13. Gürtler LG, Hauser PH, Eberle J, *et al.* A new subtype of human immunodeficiency virus type 1 (MVP-5180) from Cameroon. *J Virology* 1994;68(3):1581-5.
14. Peeters M, Gueye A, Mboup S, *et al.* Geographical distribution of HIV-1 group O viruses in Africa. *AIDS* 1997;11:493-8.
15. Hunt JC, Golden AM, Lund JK, *et al.* Envelope sequence variability and serologic characterization of HIV type 1 group O isolates from Equatorial Guinea. *AIDS Res Hum Retroviruses* 1997;13(12):995-1005.
16. Hackett J Jr, Zekeng L, Brennan CA, *et al.* Genetic analysis of HIV type 1 group O p24⁹ sequences from Cameroon and Equatorial Guinea. *AIDS Res Hum Retroviruses* 1997;13(13):1155-8.
17. Vanden Haesevelde M, Decourt J, De Leys RJ, *et al.* Genomic cloning and complete sequence analysis of a highly divergent African human immunodeficiency virus isolate. *J Virology* 1994;68(3):1586-96.
18. Hampl H, Sawitzky D, Stöffler-Meilicke M, *et al.* First case of HIV-1 subtype O infection in Germany. *Infection* 1995;23:369-70.
19. Loussert-Ajaka I, Chaix M-L, Korber B, *et al.* Variability of human immunodeficiency virus type 1 group O strains isolated from Cameroonian patients living in France. *J Virology* 1995;69:5640-9.
20. Soriano V, Gutiérrez M, García-Lerma G *et al.* First case of HIV-1 group

O infection in Spain. *Vox Sang* 1996;71:66.

21. Rayfield MA, Sullivan P, Bondea CI, *et al.* HIV-1 group O virus identified for the first time in the United States. *Emerging Infectious Diseases* 1996;2:209-12.
22. Sullivan PS, Do AN, Robbins K, *et al.* Surveillance for variant strains of HIV: subtype G and group O HIV-1. *JAMA* 1997;278:292.
23. Peeters M. Recombinant HIV sequences: their role in the global epidemic. In: Kuiken C, Foley B, Hahn B, *et al.*, editors. *HIV Sequence Compendium 2000*. Los Alamos, NM: Los Alamos National Laboratory, 2000:54-73. (Available on-line at <http://hiv-web.lanl.gov>)
24. Yamaguchi J, Vallari A, McArthur C, *et al.* Complete genome sequence of CG-0018a-01 establishes HIV-1 subtype L. *J Acquir Immune Defic Syndr* 2020;83(3):319-322.
25. Montavon C, Toure-Kane C, Liegeois F, *et al.* Most *env* and *gag* subtype A HIV-1 viruses circulating in West and West Central Africa are similar to the prototype AG recombinant virus IBNG. *J Acquir Immune Defic Syndr* 2000;23:363-74.
26. Brennan CA, Lund JK, Golden A, *et al.* Serologic and phylogenetic characterization of HIV-1 subtypes in Uganda. *AIDS* 1997;11:1823-32.
27. Vidal N, Peeters M, Mulanga-Kabeya C, *et al.* Unprecedented degree of human immunodeficiency virus type 1 (HIV-1) group M genetic diversity in the Democratic Republic of Congo suggests that the HIV-1 pandemic originated in Central Africa. *J Virology* 2000;74:10498-507.
28. Hussein M, Abebe A, Pollakis G, *et al.* HIV-1 subtype C in commercial sex workers in Addis Ababa, Ethiopia. *J Acquir Immune Defic Syndr* 2000;23:120-7.
29. Van Harmelen JH, Van Der Ryst E, Loubser AS, *et al.* A predominately HIV type 1 subtype C-restricted epidemic in South African urban populations. *AIDS Res Hum Retroviruses* 1999;15:395-8.
30. Couturier E, Damond F, Roques P, *et al.* HIV-1 diversity in France, 1996-1998. *AIDS* 2000;14:289-96.
31. Franssen K, Buvé A, Nkengasong JN, *et al.* Longstanding presence in Belgians of multiple non-B HIV-1 subtypes. *Lancet* 1996;347:1403.
32. Weniger BG, Takebe Y, Ou C, *et al.* The molecular epidemiology of HIV in Asia. *AIDS* 1994;8(S2):S13-28.
33. Russell KL, Carcamo C, Watts DM, *et al.* Emerging genetic diversity of HIV-1 in South America. *AIDS* 2000;14:1785-91.
34. Tanuri A, Swanson P, Devare S, *et al.* HIV-1 subtypes among blood donors from Rio de Janeiro, Brazil. *J Acquir Immune Defic Syndr* 1999;20:60-6.
35. Marlink RG, Ricard D, M'Boup S, *et al.* Clinical, hematologic, and immunologic cross-sectional evaluation of individuals exposed to human immunodeficiency virus type-2 (HIV-2). *AIDS Res Hum Retroviruses* 1988;4:137-48.
36. Simon F, Matheron S, Tamalet C, *et al.* Cellular and plasma viral load in patients infected with HIV-2. *AIDS* 1993;7:1411-7.
37. Kanki PJ, Travers KU, M'Boup S, *et al.* Slower heterosexual spread of HIV-2 than HIV-1. *Lancet* 1994;343:943-6.
38. DeCock KM, Adjuorolo G, Ekplini E, *et al.* Epidemiology and transmission of HIV-2: why there is no HIV-2 pandemic. *JAMA* 1993;270:2083-6.
39. Gao F, Yue L, Robertson DL, *et al.* Genetic diversity of human immunodeficiency virus type 2: evidence for distinct sequence subtypes with differences in virus biology. *J Virology* 1994;68:7433-47.
40. Yamaguchi J, Devare SG, Brennan CA. Identification of a new HIV-2 subtype based on phylogenetic analysis of full-length genomic sequence. *AIDS Res Hum Retroviruses* 2000;16:925-30.
41. Dawson GJ, Heller JS, Wood CA, *et al.* Reliable detection of individuals seropositive for the human immunodeficiency virus (HIV) by Competitive Immunoassays using *Escherichia coli*-expressed HIV structural proteins. *J Infect Dis* 1988;157:149-55.
42. Montagnier L, Clavel F, Krust B, *et al.* Identification and antigenicity of the major envelope glycoprotein of lymphadenopathy-associated virus. *Virology* 1985;144:283-9.
43. Barin F, McLane MF, Allan JS, *et al.* Virus envelope protein of HTLV-III represents major target antigen for antibodies in AIDS patients. *Science* 1985; 228:1094-6.
44. Schulz TF, Aschauer JM, Hengster P, *et al.* Envelope gene-derived recombinant peptide in the serodiagnosis of human immunodeficiency virus infection. *Lancet* 1986;2:111-2.
45. Allan JS, Coligan JE, Barin F, *et al.* Major glycoprotein antigens that induce antibodies in AIDS patients are encoded by HTLV-III. *Science* 1985;228:1091-4.
46. HIV-1/HIV-2/SIV Complete genomes. In: Kuiken CL, Foley B, Freed E, *et al.*, editors. *HIV Sequence Compendium* 2002. Los Alamos, NM: Los Alamos National Library; 2002:377-480.
47. Yamaguchi J, Bodelle P, Kaptué L, *et al.* Near full-length genomes of 15 HIV type 1 group O isolates. *AIDS Res Hum Retroviruses* 2003;19:979-88.
48. Hunt JC, Johnson-Paepke J, Boardway K, *et al.* Discrimination between HIV-1 and HIV-2 seropositive individuals using mouse monoclonal antibodies directed to HIV transmembrane proteins. *AIDS Res Hum Retroviruses* 1990;6:883-98.
49. Barin F, M'Boup S, Denis F, *et al.* Serological evidence for virus related to simian T-lymphotropic retrovirus III in residents of West Africa. *Lancet* 1985;2:1387-9.
50. Kanki PJ, Barin F, M'Boup S, *et al.* New human T-lymphotropic retrovirus related to simian T-lymphotropic virus type III (STLV-III_{AGM}). *Science* 1986;232:238-43.
51. Kanki PJ, M'Boup S, Ricard D, *et al.* Human T-lymphotropic virus type 4 and the human immunodeficiency virus in West Africa. *Science* 1987;236:827-31.
52. Clavel F, Guyader M, Guétard D, *et al.* Molecular cloning and polymorphism of the human immune deficiency virus type 2. *Nature* 1986;324:691-5.
53. Cot M, Poulain M, Delagneau JF *et al.* Dual HIV-1 and HIV-2 infection in West Africa supported by synthetic peptide analysis. *AIDS Res Hum Retroviruses* 1988;4:239-41.
54. Simon F, Cot M, Lesager C, *et al.* Differentiation between HIV-1 and HIV-2 infection by radioimmunoprecipitation and synthetic peptides in double reactive sera. *AIDS* 1989;3:401-4.
55. Norrby E, Biberfeld G, Chiodi F, *et al.* Discrimination between antibodies to HIV and to related retroviruses using site-directed serology. *Nature* 1987;329:248-50.
56. Gürtler LG, Zekeng L, Simon F, *et al.* Reactivity of five anti-HIV-1 subtype O specimens with six different anti-HIV screening ELISAs and three immunoblots. *J Virol Methods* 1995;51:177-84.
57. Loussert-Ajaka I, Ly TD, Chaix ML, *et al.* HIV-1/HIV-2 seronegativity in HIV-1 subtype O infected patients. *Lancet* 1994;343:1393-4.
58. Devare SG, Desai SM, Dawson GJ, *et al.* Diagnosis and monitoring of HIV-1 and HIV-2 infection. In: Khan NC, Melnick JL, editors. *Human Immunodeficiency Virus: Innovative Techniques for Isolation and Identification (Monographs in Virology, vol.18)*. Basel: S Karger; 1990; 18:105-21.
59. Kessler HA, Blaauw B, Spear J, *et al.* Diagnosis of human immunodeficiency virus infection in seronegative homosexuals presenting with an acute viral syndrome. *JAMA* 1987;258:1196-9.
60. Phair JP. Human immunodeficiency virus antigenemia. *JAMA* 1987;258:1218.
61. Allain J-P, Laurian Y, Paul DA, *et al.* Serological markers in early stages of human immunodeficiency virus infection in haemophiliacs. *Lancet* 1986;2(8518):1233-6.
62. Kenny C, Parkin J, Underhill G, *et al.* HIV antigen testing. *Lancet* 1987;1(8532):565-6.
63. Wall RA, Denning DW, Amos A. HIV antigenaemia in acute HIV infection. *Lancet* 1987;1(8532):566.
64. Stute R. HIV antigen detection in routine blood donor screening. *Lancet* 1987;1(8532):566.
65. Goudsmit J, De Wolf F, Paul DA, *et al.* Expression of human immunodeficiency virus antigen (HIV-Ag) in serum and cerebrospinal fluid during acute and chronic infection. *Lancet* 1986;2(8500):177-80.
66. von Sydow M, Gaines H, Sonnerberg A, *et al.* Antigen detection in primary HIV infection. *Brit Med J* 1988;296:238-40.
67. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
68. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
69. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
70. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
71. U.S. Department of Health and Human Services, Food and Drug Administration, Center for Biologics Evaluation and Research. Guidance for Industry Recommendations for Obtaining a Labeling Claim for Communicable Disease Donor Screening Tests Using Cadaveric Blood Specimens from Donors of Human Cells, Tissues, and Cellular

and Tissue-Based Products (HCT/PS), November 2004.

<http://www.fda.gov/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Tissue/ucm073972.htm>

Accessed July 01, 2013.

72. Primus FJ, Kelley EA, Hansen HJ, *et al.* "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibody therapy. *Clin Chem* 1988;34:261–4.
73. Schroff RW, Foon KA, Beatty SM, *et al.* Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879–85.
74. Boscato LM and Stuart MC. Heterophilic antibodies; a problem for all immunoassays. *Clin Chem* 1988;34(1):27–33.
75. Box GEP, Hunter WG, Hunter JS. *Statistics for experimenters; an introduction to design, data analysis, and model building*. New York, NY: John Wiley & Sons, Inc, 1978:510-39, 571–83.
76. SAS Institute, Inc. SAS Technical Report P-229, *SAS/STAT Software: Changes and enhancements, Release 6.07*. Cary, NC: SAS Institute, 1992:289–366.

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BIBLIOGRAFIE

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

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BIBLIOGRAFIE

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

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PRODUCT OF GERMANY

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EU Declaration of Conformity

Basic UDI-DI: 038074ARJ0427MG
Basic UDI-DI Name: ARCHITECT HIV Ag/Ab Combo
Risk Class: Class D

List Number and Size Code	Product and Trade Name	GMDN Code	EMDN Code
4J27-27	ARCHITECT HIV Ag/Ab Combo Reagent Kit	48446	W0105030404
4J27-32	ARCHITECT HIV Ag/Ab Combo Reagent Kit	48446	W0105030404
4J27-37	ARCHITECT HIV Ag/Ab Combo Reagent Kit	48446	W0105030404
4J27-03	ARCHITECT HIV Ag/Ab Combo Calibrator	48448	W0105080903
4J27-12	ARCHITECT HIV Ag/Ab Combo Controls	48449	W0105080803

Manufacturer (Name and Address)	Abbott GmbH, Max-Planck-Ring 2, 65205 Wiesbaden, Germany		
Manufacturer SRN	DE-MF-000009455		
Authorized Representative (Name and Address)	N/A		
Authorized Representative SRN	N/A		
Produced by (Site of manufacture) (Name and Address)	Abbott GmbH, Max-Planck-Ring 2, 65205 Wiesbaden, Germany		
Notified Body (Name and Identification Number)	TÜV SÜD Product Service GmbH Ridlerstraße 65, 80339 Munich, Germany Notified Body Number 0123		
Conformity Assessment Procedure	Quality Management System Annex IX Chapters I and III	EU Certificate No. V10 010051 0150	
	Technical Documentation Assessment Annex IX Chapter II	EU Certificate No. V70 010051 0155	
Common Specifications (CS)	COMMISSION IMPLEMENTING REGULATION (EU) 2022/1107 of 4 July 2022 laying down common specifications for certain class D in vitro diagnostic medical devices in accordance with Regulation (EU) 2017/746 of the European Parliament and of the Council.		

We, the undersigned, hereby declare that the in vitro diagnostic medical device(s) described above conform with the applicable provisions of the Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on In Vitro Diagnostic Medical Devices. This declaration is made in accordance with Annex IV of the IVD Regulation and is issued under the sole responsibility of the manufacturer.

Full Name: Claudia Becker	Full Name: Bridget Norton
Function: Director Quality Assurance	Function: Assoc. Director Regulatory Affairs
Signature:	Signature:
Date of Approval: 15 Jul 2025	Date of Approval: 15 July 2025
Signed for, and on behalf of: Abbott GmbH, Wiesbaden, Germany	
Date Issued: 15 Jul 2025	Place Issued: 65205 Wiesbaden, Germany
Supersedes: 28-May-2024	Effective (Date or Lot Number): 15 Jul 2025

EU Technical Documentation Assessment Certificate (IVDR)

Pursuant to Regulation (EU) 2017/746 on in Vitro Diagnostic Medical Devices,
Annex IX, Chapter II (Class D Devices)

No. V70 010051 0155 Rev. 00

Manufacturer:

Abbott GmbH

Max-Planck-Ring 2
65205 Wiesbaden
GERMANY

SRN Manufacturer - DE-MF-000009455

The Certification Body of TÜV SÜD Product Service GmbH certifies that the manufacturer has drawn up and presented a Technical Documentation according to Annex II and III of the Regulation (EU) 2017/746 on in Vitro Diagnostic Medical Devices. Details on devices covered by the Technical Documentation are described on the following page(s).

The Report referenced below summarizes the result of the assessment and includes reference to relevant CS, harmonized standards and test reports. The conformity assessment has been carried out according to Annex IX, Chapter II, of this regulation with a positive result. In order to maintain this certificate, the manufacturer shall submit Periodic Safety Update Reports at least annually to the notified body TÜV SÜD Product Service GmbH. Verification of manufactured class D devices according to Annex IX Sections 4.12 and 4.13 is applicable. In order to place the devices on the market with CE-marking, an EU Quality Management System Certificate pursuant to Annex IX Chapters I and III is necessary in addition to this EU Technical Documentation Assessment Certificate. All applicable requirements of the Testing, Certification, Validation and Verification Regulations TÜV SÜD Group have to be complied with.

For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:V70 010051 0155 Rev. 00

Report No.: 713302431

Valid from: 2024-04-09

Valid until: 2029-04-08

Marta Carnielli
Head of Certification IVD

Issue date: 2024-04-09

EU Technical Documentation Assessment Certificate (IVDR)

Pursuant to Regulation (EU) 2017/746 on in Vitro Diagnostic Medical Devices,
Annex IX, Chapter II (Class D Devices)

No. V70 010051 0155 Rev. 00

Classification: Class D
Device Group: W010503 - RETROVIRUSES (INFECT. IMMUNOLOGY/NAT)
Basic UDI-DI: 038074ARJ0427MG

Intended Purpose: Device for the automated simultaneous qualitative detection of HIV p24 antigen and antibodies to human immunodeficiency virus type 1 and/or type 2 (HIV-1/HIV-2) in human serum and plasma of individuals suspected to have HIV-1/HIV-2 infection and potential blood/cell/tissue/organ donors (incl. specimens collected post-mortem (non-heart-beating)). Aid in the diagnosis of HIV-1/HIV-2 infection and as a screening test to prevent transmission of HIV-1/HIV-2 to recipients of blood, blood components, cells, tissue and organs. HIV is the etiologic agent of AIDS. HIV is transmitted by sexual contact, exposure to blood or blood products, and prenatal infection of a fetus or perinatal infection of a newborn.
For Laboratory Professional Use Only.

Device(s): ARCHITECT HIV Ag/Ab Combo Reagent Kit - LN 4J27-27
ARCHITECT HIV Ag/Ab Combo Reagent Kit - LN 4J27-32
ARCHITECT HIV Ag/Ab Combo Reagent Kit - LN 4J27-37
ARCHITECT HIV Ag/Ab Combo Calibrator - LN 4J27-03
ARCHITECT HIV Ag/Ab Combo Controls - LN 4J27-12

Classification: Class D
Device Group: W010503 - RETROVIRUSES (INFECT. IMMUNOLOGY/NAT)
Basic UDI-DI: 038074AIP0807LB

Intended Purpose: Device for the automated simultaneous qualitative detection of HIV p24 antigen and antibodies to human immunodeficiency virus type 1 and/or type 2 (HIV-1/HIV-2) in human serum and plasma of individuals suspected to have HIV-1/HIV-2 infection and potential blood/cell/tissue/organ donors (incl. specimens collected post-mortem (non-heart-beating)). Aid in the diagnosis of HIV-1/HIV-2 infection and as a screening test to prevent transmission of HIV-1/HIV-2 to recipients of blood, blood components, cells, tissue and organs. HIV is the etiologic agent of AIDS. HIV is transmitted by sexual contact, exposure to blood or blood products, and prenatal infection of a fetus or perinatal infection of a newborn.
For Laboratory Professional Use Only.

Device(s): Alinity i HIV Ag/Ab Combo Reagent Kit - LN 08P0722
Alinity i HIV Ag/Ab Combo Reagent Kit - LN 08P0732
Alinity i HIV Ag/Ab Combo Calibrator - LN 08P0701
Alinity i HIV Ag/Ab Combo Controls - LN 08P0710

EU Technical Documentation Assessment Certificate (IVDR)

Pursuant to Regulation (EU) 2017/746 on in Vitro Diagnostic Medical Devices,
Annex IX, Chapter II (Class D Devices)

No. V70 010051 0155 Rev. 00

Classification: Class D
Device Group: W010503 - RETROVIRUSES (INFECT. IMMUNOLOGY/NAT)
Basic UDI-DI: 038074WIE0009T4

Intended Purpose: Device for the automated simultaneous qualitative detection of HIV p24 antigen and antibodies to human immunodeficiency virus type 1 and/or type 2 (HIV-1/HIV-2) in human serum and plasma of individuals suspected to have HIV-1/HIV-2 infection and potential blood/cell/tissue/organ donors (incl. specimens collected post-mortem (non-heart-beating)). Aid in the diagnosis of HIV-1/HIV-2 infection and as a screening test to prevent transmission of HIV-1/HIV-2 to recipients of blood, blood components, cells, tissue and organs. HIV is the etiologic agent of AIDS. HIV is transmitted by sexual contact, exposure to blood or blood products, and prenatal infection of a fetus or perinatal infection of a newborn.
For Laboratory Professional Use Only.

Device(s):
HIV Ag/Ab Combo Reagent Kit - LN 4J27-84
HIV Ag/Ab Combo Reagent Kit - LN 4J27-87
HIV Ag/Ab Combo Reagent Kit - LN 4J27-89
HIV Ag/Ab Combo Calibrator - LN 4J27-09
HIV Ag/Ab Combo Controls - LN 4J27-17

Classification: Class D
Device Group: W010503 - RETROVIRUSES (INFECT. IMMUNOLOGY/NAT)
Basic UDI-DI: 038074WIE0010SM

Intended Purpose: Device for the automated simultaneous qualitative detection of HIV p24 antigen and antibodies to human immunodeficiency virus type 1 and/or type 2 (HIV-1/HIV-2) in human serum and plasma of individuals suspected to have HIV-1/HIV-2 infection and potential blood/cell/tissue/organ donors (incl. specimens collected post-mortem (non-heart-beating)). Aid in the diagnosis of HIV-1/HIV-2 infection and as a screening test to prevent transmission of HIV-1/HIV-2 to recipients of blood, blood components, cells, tissue and organs. HIV is the etiologic agent of AIDS. HIV is transmitted by sexual contact, exposure to blood or blood products, and prenatal infection of a fetus or perinatal infection of a newborn.
For Laboratory Professional Use Only.

Device(s):
HIV Ag/Ab Combo Reagent Kit - LN 08P0784
HIV Ag/Ab Combo Reagent Kit - LN 08P0787
HIV Ag/Ab Combo Calibrator - LN 08P0709
HIV Ag/Ab Combo Controls - LN 08P0717

EU Technical Documentation Assessment Certificate (IVDR)

Pursuant to Regulation (EU) 2017/746 on in Vitro Diagnostic Medical Devices,
Annex IX, Chapter II (Class D Devices)

No. V70 010051 0155 Rev. 00

The validity of this certificate -none-
depends on conditions and/or
is limited to the following:

Revision History:

Rev.	Dated	Report	Description
00	2024-04-09	713302431	Initial issuance

Management Service

CERTIFICATE

Certificate Registration No.: **12 100 64551/02 TMS** / Order No.: **707151799**

**The Certification Body
of TÜV SÜD Management Service GmbH**

certifies that the organization

**Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden
Germany**

at the site

Abbott

**Abbott Diagnostics GmbH
Max-Planck-Ring 2
65205 Wiesbaden
Germany**

for the scope

**The provision of Warehousing and Distribution services of
In-Vitro Diagnostic Medical Devices and Medical Devices**

has established and applies a Quality Management System.

An audit was performed and has furnished proof
that the requirements according to

DIN EN ISO 9001:2015

are fulfilled.

The certificate is valid in conjunction
with the main certificate from **2024-10-01** until **2027-09-30**.

Fred Wenke
Head of Certification Body
Munich, 2024-07-16

Certificate

No. Q5 010051 0139 Rev. 01

Holder of Certificate: **Abbott GmbH**
Max-Planck-Ring 2
65205 Wiesbaden
GERMANY

Certification Mark:

Scope of Certificate: **Design and Development, Manufacture and Warehousing of In-vitro Diagnostic Reagents for Clinical Chemistry, Immunochemistry, Hematology and Infectious Immunology. The provision of Warehousing and Distribution services of In-vitro Diagnostic medical devices and medical devices.**

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). All applicable requirements of the Testing, Certification, Validation and Verification Regulations TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:Q5 010051 0139 Rev. 01

Report No.: 713332354_IVDR

Valid from: 2024-10-01
Valid until: 2027-09-30

Date, 2024-07-11

Christoph Dicks
Head of Certification/Notified Body

Certificate

No. Q5 010051 0139 Rev. 01

Applied Standard(s): ISO 13485:2016
(EN ISO 13485:2016/AC:2018, EN ISO 13485:2016/A11:2021)
Medical devices - Quality management systems -
Requirements for regulatory purposes

Facility(ies): **Abbott GmbH**
Max-Planck-Ring 2, 65205 Wiesbaden, GERMANY

Design and Development, Manufacture and Warehousing of
In-vitro Diagnostic Reagents for Clinical Chemistry,
Immunochemistry, Hematology and Infectious Immunology.

Abbott Diagnostics GmbH
Max-Planck-Ring 2, 65205 Wiesbaden, GERMANY

The provision of Warehousing and Distribution services of
In-vitro Diagnostic medical devices and medical devices.

./.

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© 2004, 2022 Abbott Laboratories