

Triglyceride2

FOR USE WITH
ARCHITECT

Revised May 2023.

REF 04T1020

REF 04T1030

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

For laboratory professional use only.

NAME

Triglyceride2 (also referred to as Trig2)

INTENDED USE

The Triglyceride2 assay is used for the quantitation of triglyceride in human serum or plasma on the ARCHITECT c Systems.

The Triglyceride2 assay is to be used as an aid in the diagnosis and treatment of diabetes mellitus, nephrosis, liver obstruction, other diseases involving lipid metabolism, or various endocrine disorders.

SUMMARY AND EXPLANATION OF THE TEST

Triglycerides are a family of lipids absorbed from the diet and produced endogenously from carbohydrates and fatty acids. Measurement of triglyceride is important in the diagnosis and management of hyperlipidemia. These diseases can be genetic or secondary to other disorders including nephrosis, diabetes mellitus, and endocrine disturbances. The National Cholesterol Education Program (NCEP) cites evidence that triglycerides are an independent risk factor for atherosclerosis.¹ Individuals with hypertension, obesity, and/or diabetes are at greater risk than are those without these conditions.^{2, 3}

The Adult Treatment Panel of the NCEP recommends that all adults 20 years of age and over should have a fasting lipoprotein profile (total cholesterol, LDL cholesterol, HDL cholesterol, and triglyceride) once every five years to screen for coronary heart disease risk.¹

PRINCIPLES OF THE PROCEDURE

The Triglyceride2 assay is an automated clinical chemistry assay. Triglycerides are enzymatically hydrolyzed by lipase to free fatty acids and glycerol. The glycerol is phosphorylated by adenosine triphosphate (ATP) with glycerol kinase (GK) to produce glycerol-3-phosphate and adenosine diphosphate (ADP). Glycerol-3-phosphate is oxidized to dihydroxyacetone phosphate (DAP) by glycerol phosphate oxidase (GPO) producing hydrogen peroxide (H₂O₂). In a color reaction catalyzed by peroxidase, the H₂O₂ reacts with 4-aminoantipyrine (4-AAP) and n,n-bis(4-sulfobutyl)-3-methylaniline (TODB) to produce a red colored quinoneimine dye. The absorbance of this dye at 604 nm is proportional to the concentration of triglyceride present in the sample.

Methodology: Glycerol phosphate oxidase

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Kit Contents

Triglyceride2 Reagent Kit 04T10

NOTE: Some kit sizes may not be available. Please contact your local distributor.

Volumes (mL) listed in the following table indicate the volume per cartridge.

REF	04T1020	04T1030
Tests per cartridge	200	475
Number of cartridges per kit	4	4
Tests per kit	800	1900
R1	40.0 mL	90.5 mL

R1 Active ingredients: adenosine-5'-triphosphate (ATP) 1.210 g/L, 4-aminoantipyrine (4-AAP) 0.061 g/L, glycerol kinase (GK) 0.400 KU/L, lipoprotein lipase 3.000 KU/L, L-glycerol-3-phosphate oxidase 4.000 KU/L, n,n-bis(4-sulfobutyl)-3-methylaniline (TODB) 0.212 g/L and peroxidase (POD) 4.000 KU/L. Preservative: sodium azide.

Warnings and Precautions

- IVD
- For *In Vitro* Diagnostic Use
- Rx ONLY

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials and all consumables contaminated with potentially infectious materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate regional, national, and institutional biosafety practices should be used for materials that contain, are suspected of containing, or are contaminated with infectious agents.⁴⁻⁷

The following warnings and precautions apply to: R1	
WARNING	Contains PIPES sodium salt* and sodium azide.
H316*	Causes mild skin irritation.
EUH032	Contact with acids liberates very toxic gas.
Response	
P332+P313*	If skin irritation occurs: Get medical advice / attention.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

* Not applicable where regulation EC 1272/2008 (CLP) or OSHA Hazard Communication 29 CFR 1910.1200 (HCS) 2012 have been implemented.

Follow local chemical disposal regulations based on your location along with recommendations and content in the Safety Data Sheet to determine the safe disposal of this product.

For the most current hazard information, see the product Safety Data Sheet.

Safety Data Sheets are available at www.corelaboratory.abbott or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not pool reagents within a kit or between kits.
- Do not reuse containers, caps or plugs due to the risk of contamination and the potential to compromise reagent performance.
- Upon receipt, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position.
Onboard	System Temperature	30 days	
Opened	2 to 8°C	Until expiration date	Store in upright position.

Reagents may be stored on or off the ARCHITECT c System. If reagents are removed from the system, store at 2 to 8°C (with replacement caps) in their original boxes.

For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Triglyceride2 assay file must be installed on the ARCHITECT c System prior to performing the assay.

Installation of all the required SmartWash updates on the ARCHITECT c Systems Assay Disk Version 17.00 (or higher) must be completed prior to performing the assay. See below for impacted assays:

Assay Name	Short Name	REF	Assay Number	VERSION	
				Conventional Units / Alternate Units	SI Units / Alternate Units
Albumin BCG	AlbG	7D53	1015	12	10
Lipase	Lip	7D80	1029	12	12
Magnesium	MAG	3P68	1070	6	4
Magnesium Urine	MAGU	3P68	1099	8	4
Ultra HDL	UHDL	3K33	1093	7 (c8000) 4 (c4000, c16000)	5 (c8000) 4 (c4000, c16000)

The ARCHITECT System software version 9.25 or higher must be installed on the ARCHITECT c System prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

Alternate Result Units

Conversion formula:

$$(\text{Concentration in Default result unit}) \times (\text{Conversion factor}) = (\text{Concentration in Alternate result unit})$$

Default Result Unit	Conversion Factor	Alternate Result Unit
mg/dL	0.0113	mmol/L

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen types listed below were verified for use with this assay.

Other specimen types and collection tube types have not been verified with this assay.

Specimen Types	Collection Tubes
Serum	Serum Serum separator
Plasma	Lithium heparin Lithium heparin separator Sodium heparin

- Liquid anticoagulants may have a dilution effect resulting in lower concentration values for individual specimens.

The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - pooled specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
 - specimens with fungal growth
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low-speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low-speed vortex or by inverting 10 times.

- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Recentrifugation of Specimens

- Transfer specimens to a centrifuge tube and centrifuge.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

Specimen Storage

Specimen Type	Temperature	Maximum Storage Time
Serum/Plasma	Room temperature (20 to 25°C)	2 days ⁸
	2 to 8°C	7 days ⁸
	-20°C	3 months ⁹

Avoid multiple freeze/thaw cycles.⁹

It is the responsibility of the individual laboratory to determine specific specimen stability criteria for their laboratory per their laboratory workflow.

For additional information on sample handling and processing, refer to CLSI GP44-A4.¹⁰ The storage information provided here is based on references.

Each laboratory may establish a range around -20°C from either the freezer manufacturer's specifications or your laboratory standard operating procedure(s) for specimen storage.

Stored specimens must be inspected for particulates. If present, mix with a low-speed vortex or by inversion and centrifuge the specimen to remove particulates prior to testing.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

Do not exceed the storage limitations listed above.

PROCEDURE

Materials Provided

04T10 Triglyceride2 Reagent Kit

Materials Required but not Provided

- Triglyceride2 assay file found on www.corelaboratory.abbott
- 04V1501 Consolidated Chemistry Calibrator
- Controls containing triglyceride
- Saline (0.85% to 0.90% NaCl) for specimen dilution

For information on materials required for operation of the instrument, refer to the ARCHITECT System Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the ARCHITECT System Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the ARCHITECT System Operations Manual, Section 5 to ensure sufficient specimen is present.
- Minimum sample cup volume is calculated by the system and printed on the Order List report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Minimum sample volume requirements:
 - Sample volume for single test: 1.6 µL.

NOTE: This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the ARCHITECT System Operations Manual, Section 5.

- Refer to the Consolidated Chemistry Calibrator package insert [REF] 04V1501 and/or commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Samples with a triglyceride value exceeding

1505 mg/dL (17.01 mmol/L) are flagged with the code "> 1505 mg/dL" (> 17.01 mmol/L) and may be diluted with either the Automated Dilution Protocol or the Manual Dilution Procedure. In addition, samples with a Reaction check failure are flagged with the error code "1054", and may be diluted with either the Automated Dilution Protocol or the Manual Dilution Procedure. Refer to the ARCHITECT System Operations Manual, Section 10.

Automated Dilution Protocol

The system performs a dilution of the sample, relative to the standard dilution, and automatically calculates the concentration by multiplying the result by the dilution factor.

Dilution Name	Dilution Factor
Standard	1:1.98
1:4	1:7.92

For details on configuring automated dilutions, refer to the ARCHITECT System Operations Manual, Section 2.

Manual Dilution Procedure

Dilute the sample with saline (0.85% to 0.90% NaCl).

The operator must enter the manual dilution factor in the Patient or Control order screen. The system will use this dilution factor to automatically calculate the concentration of the sample and report the result.

If the operator does not enter the manual dilution factor, the result must be manually multiplied by the appropriate manual dilution factor before reporting the result. If a diluted sample result is less than 4 mg/dL (0.06 mmol/L), do not report the result. Rerun using an appropriate dilution.

NOTE: The default Low Linearity value of the assay file corresponds to the lower limit of the reportable interval of 2 mg/dL (0.02 mmol/L). To flag values using the lower limit of the analytical measuring interval of 4 mg/dL (0.06 mmol/L), the operator must edit the Low Linearity value, adjusted by the standard dilution factor.

For detailed information on editing the result settings of assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

Calibration

For instructions on performing a calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Calibration is stable for approximately 30 days (720 hours), with a blank adjustment after 15 days (360 hours). Calibration is required with each change in reagent lot. Verify calibration with at least 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Quality Control Procedures

As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- At least 2 levels of controls (low and high) are to be run every 24 hours.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

Quality Control Guidance

Refer to "Basic QC Practices" by James O. Westgard, Ph.D. for guidance on laboratory quality control practices.¹¹

RESULTS

Calculation

The Triglyceride2 assay utilizes the Linear data reduction method to generate a calibration and results.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

Reportable Interval

Based on representative data for the limit of quantitation (LoQ) and the limit of detection (LoD), the ranges over which results can be reported are provided below according to the definitions from CLSI EP34, 1st ed.¹²

	mg/dL	mmol/L
Analytical Measuring Interval (AMI) ^a	4–1505	0.06 ^d –17.01
Extended Measuring Interval (EMI) ^b	1505–6020	17.01–68.05
Reportable Interval ^c	2–6020	0.02–68.05

^a AMI: The AMI extends from the LoQ to the upper limit of quantitation (ULoQ). This is determined by the range of values in mg/dL (mmol/L) that demonstrated acceptable performance for linearity, imprecision, and bias.

^b EMI: The EMI extends from the ULoQ to the ULoQ × sample dilution.

^c The reportable interval extends from the LoD to the upper limit of the EMI.

^d Value determined based on the dilution factor and instrument rounding.

NOTE: The default Low Linearity value of the assay file corresponds to the lower limit of the reportable interval.

LIMITATIONS OF THE PROCEDURE

- Results should be used in conjunction with other data; e.g., symptoms, results of other tests, and clinical impressions.
- Specimens with conjugated bilirubin levels greater than 2 mg/dL or unconjugated bilirubin greater than 5 mg/dL may cause falsely depressed results with the Triglyceride2 assay. Refer to the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert for further information.
- Substances that demonstrated interference with the Triglyceride2 assay are listed in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- Potential interference has not been evaluated for substances other than those described in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- SmartWashes for assays impacted by Triglyceride2 must be configured to avoid interference due to reagent carryover. See the INSTRUMENT PROCEDURE section of this package insert for the required assay file updates.

EXPECTED VALUES

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

Reference Range¹

		Range (mg/dL)	Range (mmol/L) ^a
Adult	Normal	< 150	< 1.70
	Borderline High	150–199	1.70–2.25
	High	200–499	2.26–5.64
	Very High	≥ 500	≥ 5.65

^a Alternate result units were calculated by Abbott and are not included in the citation provided.

Reference Range¹³

Pediatric Triglyceride Percentile ^a	Boys, By Age Group (Years)			Girls, By Age Group (Years)		
	5–9	10–14	15–19	5–9	10–14	15–19
75th: Acceptable	58 (0.66)	74 (0.84)	88 (0.99)	74 (0.84)	85 (0.96)	85 (0.96)
90th: Borderline	70 (0.79)	94 (1.06)	125 (1.41)	103 (1.16)	104 (1.18)	112 (1.27)
95th: High	85 (0.96)	111 (1.25)	143 (1.62)	120 (1.36)	120 (1.36)	126 (1.42)

^a Values are presented in mg/dL (mmol/L). Alternate results units (mmol/L) were calculated by Abbott and are not included in the citation provided.

Abbott has not evaluated reference ranges in the pediatric population.

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A3.¹⁴ Testing was conducted using 3 lots of the Triglyceride2 reagent, 3 lots of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 3 instruments. Two controls and 3 human serum panels were tested in a minimum of 2 replicates, twice per day on 20 days on 3 reagent lot/calibrator lot/instrument combinations, where a unique reagent lot and a unique calibrator lot are paired with 1 instrument. The performance from a representative combination is shown in the following table.

Sample	n	Mean (mg/dL)	Within-Run (Repeatability)		Within-Laboratory ^a	
			SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	80	184	1.6	0.9	3.1 (3.1–3.7)	1.7 (1.7–2.1)
Control Level 2	80	83	0.6	0.8	1.6 (1.6–1.9)	1.9 (1.9–2.3)
Panel A	80	10	0.3	2.6	1.0 (1.0–1.2)	9.8 (9.8–11.3)
Panel B	80	549	4.8	0.9	4.9 (4.9–8.0)	0.9 (0.9–1.5)
Panel C	80	1265	9.5	0.8	10.1 (10.1–20.6)	0.8 (0.8–1.7)

^a Includes within-run, between-run, and between-day variability.

^b Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot/instrument combinations.

Sample	n	Mean (mmol/L)	Within-Run (Repeatability)		Within-Laboratory ^a	
			SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	80	2.08	0.019	0.9	0.035 (0.035–0.042)	1.7 (1.7–2.0)
Control Level 2	80	0.94	0.007	0.8	0.018 (0.017–0.021)	1.9 (1.9–2.3)
Panel A	80	0.12	0.004	3.2	0.012 (0.012–0.013)	10.4 (10.1–11.3)
Panel B	80	6.21	0.054	0.9	0.056 (0.056–0.090)	0.9 (0.9–1.5)
Panel C	80	14.29	0.108	0.8	0.115 (0.115–0.231)	0.8 (0.8–1.7)

^a Includes within-run, between-run, and between-day variability.

^b Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot/instrument combinations.

Reproducibility

A study was performed based on guidance from CLSI EP05-A3.¹⁴ Testing was conducted using 1 lot of the Triglyceride2 reagent, 1 lot of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 3 instruments. Each instrument was operated by a different technician, and each technician prepared an individual sample set. Two controls and 2 human serum panels were tested in a minimum of 3 replicates at 2 separate times per day on 5 different days.

Sample	n	Mean (mg/dL)	Repeatability		Within-Laboratory ^a		Reproducibility ^b	
			SD	%CV	SD	%CV	SD	%CV
Control Level 1	90	178	1.4	0.8	1.8	1.0	3.0	1.7
Control Level 2	90	80	0.9	1.1	1.0	1.3	1.6	2.0
Panel B	90	542	7.1	1.3	7.3	1.3	12.3	2.3
Panel C	90	1265	16.4	1.3	16.5	1.3	24.8	2.0

^a Includes repeatability (within-run), between-run, and between-day variability.

^b Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

Sample	n	Mean (mmol/L)	Repeatability		Within-Laboratory ^a		Reproducibility ^b	
			SD	%CV	SD	%CV	SD	%CV
Control Level 1	90	2.02	0.016	0.8	0.020	1.0	0.034	1.7
Control Level 2	90	0.91	0.009	1.0	0.010	1.2	0.018	2.0
Panel B	90	6.13	0.080	1.3	0.082	1.3	0.139	2.3
Panel C	90	14.29	0.185	1.3	0.186	1.3	0.280	2.0

^a Includes repeatability (within-run), between-run, and between-day variability.

^b Includes repeatability (within-run), between-run, between-day, and between-instrument variability.

Accuracy

A study was performed to estimate the bias of the Triglyceride2 assay relative to standard reference material (ACS Grade Glycerol). Testing was conducted using 2 concentrations of the standard across 3 lots of the Triglyceride2 reagent, 3 lots of the Consolidated Chemistry Calibrator, and 1 instrument. The bias ranged from 0.4% to 1.7%.

Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2.¹⁵ Testing was conducted using 3 lots of the Triglyceride2 reagent on each of 2 instruments over a minimum of 3 days. The maximum observed limit of blank (LoB), limit of detection (LoD), and limit of quantitation (LoQ) values are summarized below.

	mg/dL	mmol/L
LoB ^a	1	0.01
LoD ^b	2	0.02
LoQ ^c	4	0.06

^a The LoB represents the 95th percentile from $n \geq 60$ replicates of zero-analyte samples.

^b The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on $n \geq 60$ replicates of low-analyte level samples.

^c The LoQ is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was determined from $n \geq 60$ replicates of low-analyte level samples.

Linearity

A study was performed based on guidance from CLSI EP06-A.¹⁶ This assay is linear across the analytical measuring interval of 4 to 1505 mg/dL (0.06 to 17.01 mmol/L).

Analytical Specificity

Interference

Potentially Interfering Endogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.¹⁷ Each substance was tested at 2 levels of the analyte (approximately 100 mg/dL and 250 mg/dL).

No significant interference (interference within $\pm 10\%$) was observed at the following concentrations.

No Significant Interference (Interference within $\pm 10\%$)		
Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
Bilirubin (conjugated)	2 mg/dL	24 μ mol/L
Bilirubin (unconjugated)	5 mg/dL	86 μ mol/L
Hemoglobin	1000 mg/dL	10 g/L
Total protein	15 g/dL	150 g/L

Interference beyond $\pm 10\%$ (based on 95% Confidence Interval [CI]) was observed at the concentrations shown below for the following substances.

Interference beyond $\pm 10\%$ (based on 95% Confidence Interval [CI])					
Potentially Interfering Substance	Interferent Level		Analyte Level		% Interference (95% CI)
	Default Units	Alternate Units	Default Units	Alternate Units	
Bilirubin (conjugated)	3 mg/dL	36 μ mol/L	100 mg/dL	1.13 mmol/L	-10% (-11%, -9%)
Bilirubin (unconjugated)	7 mg/dL	120 μ mol/L	100 mg/dL	1.13 mmol/L	-11% (-12%, -10%)
Total protein	17.5 g/dL	175 g/L	250 mg/dL	2.83 mmol/L	13% (12%, 13%)
Total protein	17.5 g/dL	175 g/L	100 mg/dL	1.13 mmol/L	11% (11%, 12%)

Potentially Interfering Exogenous Substances

A study was performed based on guidance from CLSI EP07, 3rd ed.¹⁷ Each substance was tested at 2 levels of the analyte (approximately 100 mg/dL and 250 mg/dL).

No significant interference (interference within $\pm 10\%$) was observed at the following concentrations.

No Significant Interference (Interference within $\pm 10\%$)		
Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
Acetaminophen	160 mg/L	1059 μ mol/L
Acetylcysteine	140 mg/L	858 μ mol/L
Acetylsalicylic acid	30 mg/L	167 μ mol/L
Alcohol (ethanol)	600 mg/dL	130 mmol/L
Ampicillin-Na	80 mg/L	215 μ mol/L
Ascorbic acid	20 mg/L	114 μ mol/L
Biotin	4250 ng/mL	17 μ mol/L
Ca-dobesilate	10 mg/L	24 μ mol/L
Cefoxitin	6600 mg/L	15 mmol/L
Cyclosporine	2 mg/L	1.7 μ mol/L
Dicycnone (ethamsylate)	15 mg/L	57 μ mol/L
Dipyrrone	100 mg/L	300 μ mol/L
Dobutamine	0.2 mg/dL	6.6 μ mol/L
Doxycycline	20 mg/L	45 μ mol/L
Hydroxyurea (hydroxycarbamide)	4 mg/dL	524 μ mol/L
Ibuprofen	220 mg/L	1067 μ mol/L
Levodopa	3 mg/L	15 μ mol/L
Methylidopa	3 mg/L	14 μ mol/L
Metronidazole	130 mg/L	759 μ mol/L
N-acetyl-4-benzoquinone imine (NAPQI)	20 mg/L	134 μ mol/L
Phenylbutazone	330 mg/L	1069 μ mol/L
Rifampicin	50 mg/L	61 μ mol/L
Sodium heparin	4 U/mL	N/A
Theophylline (1,3-dimethylxanthine)	60 mg/L	333 μ mol/L
4-Acetamidoantipyrine (4-acetylaminantipyrine)	40 mg/L	163 μ mol/L
4-Aminoantipyrine	40 mg/L	197 μ mol/L
4-Formylaminoantipyrine	40 mg/L	173 μ mol/L
4-Methylaminoantipyrine	40 mg/L	184 μ mol/L

N/A = Not Applicable

Interference beyond $\pm 10\%$ (based on 95% Confidence Interval [CI]) was observed at the concentrations shown below for the following substances.

Interference beyond $\pm 10\%$ (based on 95% Confidence Interval [CI])					
Potentially Interfering Substance	Interferent Level		Analyte Level		% Interference (95% CI)
	Default Units	Alternate Units	Default Units	Alternate Units	
Acetylcysteine	150 mg/L	920 μ mol/L	100 mg/dL	1.13 mmol/L	-10% (-11%, -9%)
Ascorbic acid	30 mg/L	170 μ mol/L	100 mg/dL	1.13 mmol/L	-21% (-22%, -21%)
Ca-dobesilate	15 mg/L	36 μ mol/L	100 mg/dL	1.13 mmol/L	-11% (-12%, -11%)
Levodopa	5 mg/L	25 μ mol/L	100 mg/dL	1.13 mmol/L	-11% (-12%, -11%)
Methylidopa	5 mg/L	24 μ mol/L	100 mg/dL	1.13 mmol/L	-10% (-11%, -10%)
Dicycnone (ethamsylate)	25 mg/L	95 μ mol/L	100 mg/dL	1.13 mmol/L	-15% (-16%, -14%)

Interferences from medication or endogenous substances may affect results.¹⁸

Method Comparison

A study was performed based on guidance from CLSI EP09-A3, Third Edition¹⁹ using the Passing-Bablok regression method.






Triglyceride2 vs. Triglyceride on the ARCHITECT c System					
	n	Units	Correlation Coefficient	Intercept	Slope
Serum	137	mg/dL (mmol/L)	1.00	1.41 (0.02)	0.99
					Concentration Range 7–1298 (0.08–14.67)

BIBLIOGRAPHY

- Cleeman JI. Executive summary of the third report of the National Cholesterol Education Program (NCEP) Expert Panel on detection, evaluation, and treatment of high blood cholesterol in adults (Adult Treatment Panel III). *JAMA* 2001;285(19):2486-2497.
- Rubins HB. Triglycerides and coronary heart disease: implications of recent clinical trials. *J Cardiovasc Risk* 2000;7(5):339-345.
- Forrester JS. Triglycerides: risk factor or fellow traveler? *Curr Opin Cardiol* 2001;16:261-264.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Cuhadar S, Atay A, Koseoglu M, et al. Stability studies of common biochemical analytes in serum separator tubes with or without gel barrier subjected to various storage conditions. *Biochem Med* 2012;22(2):202-214.
- Cuhadar S, Koseoglu M, Atay A, et al. The effect of storage time and freeze-thaw cycles on the stability of serum samples. *Biochem Med* 2013;23(1):70-77.
- Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: CLSI; 2010.
- Westgard JO. *Basic QC Practices; Training in Statistical Quality Control for Medical Laboratories*. 4th ed. Westgard QC, Inc.; 2016.
- Clinical and Laboratory Standards Institute (CLSI). *Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking*. 1st ed. CLSI Guideline EP34. Wayne, PA: CLSI; 2018.
- Miller M, Stone NJ, Ballantyne C, et al. Triglycerides and cardiovascular disease: a scientific statement from the American Heart Association. *J Am Heart Assoc* 2011;123:2292–2333.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures: Approved Guideline—Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI Guideline EP07. Wayne, PA: CLSI; 2018.
- Young DS. Laboratory test listings. In: *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. AACC Press; 2000:chap 3.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

■ Key to Symbols

ISO 15223 Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
REF	List Number
SN	Serial number

Other Symbols

CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
FOR USE WITH	Identifies products to be used together
PRODUCT OF IRELAND	Product of Ireland
R1	Reagent 1
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

The ARCHITECT c System family of instruments consists of c4000, c8000, and c16000 instruments.

ARCHITECT and related brand marks are trademarks of Abbott.

Other trademarks are the property of their respective owners.



Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000



Customer Service: Contact your local representative or find country-specific contact information on www.corelaboratory.abbott

For customers in the European Union: if, in the course of using this device, you have reason to believe that a serious incident has occurred, report it to the manufacturer and to your national authority.

Revised May 2023.

©2021, 2023 Abbott Laboratories