

Promus ELITE™ Stent System

In-Service Presentation




Promus ELITE™
Everolimus-Eluting Platinum Chromium
Coronary Stent System

Promus ELITE™ Stent System Overview

What Have We Achieved Together?

Continued innovation and advancing DES technology



Promus ELITE™ Stent


Overview

- Builds on the **proven performance** of Promus PP-DES family and PtCr platform
- **Customized stent architecture** provides:
 - Radial strength
 - Flexibility
 - Reduced recoil
 - Controlled overexpansion
- Market-leading **Everolimus and PVDF-HFP polymer**
- **Outstanding permanent polymer DES safety** as demonstrated in randomized and real-world clinical trials
 - PLATINUM Family of trials
 - Lowest relative risk of Def/Prob Stent Thrombosis among PP-DES in Kang Meta-Analysis²
 - Numerically lowest PP-DES ST rates in real-world SCAAR Registry³

More than
8M
Promus DES
implanted WW¹



Market-leading Everolimus
Drug & PVDF-HFP Polymer



Outstanding deliverability and acute performance¹


Promus PREMIER™ Stent System

Customized Architecture and Outcomes



Proprietary Laser-Cut Hypotube Technology

Improved pushability and trackability



Promus ELITE™

Deliverability and Acute Performance

What Makes the Promus ELITE™ Stent Unique?

Boston
Scientific

KEY COMPONENTS

ADVANTAGES



Enhanced
Delivery System

Improved Push and Deliverability¹



Promus PREMIER
Stent Architecture

Ideal Balance of Strength and
Flexibility



Market-leading Everolimus
Drug & PVDF-HFP Polymer

Outstanding Safety and Efficacy²

Promus ELITE™ Stent System

Outstanding Acute Performance. Proven Long-term Outcomes

Boston
Scientific




Outstanding Deliverability



Trusted Stent Platform



Proven Clinical Outcomes



Promus ELITE™ Stent System

Outstanding Acute Performance. Proven Long-term Outcomes

Boston
Scientific



Outstanding Deliverability



Trusted Stent Platform



Proven Clinical Outcomes

What makes a stent system deliverable?

Outstanding Deliverability

Does it fit?
Lesion Entry Profile
Average Stent Profile

Can it make the turn?
Trackability
System flex
Tip flexibility



Is it responsive?
Pushability
Responsiveness


All these attributes combined = deliverability

Promus ELITE™ Stent System

New Enhanced Delivery System

Boston
Scientific

Outstanding Deliverability



KEY ADVANTAGES¹

- Improved pushability
- Exceptional trackability
- Outstanding flexibility

Promus ELITE™ Stent System

New Enhanced Delivery System

Boston
Scientific


Outstanding Deliverability

Bi-Segment inner shaft provides balance of proximal push and distal flex

Laser-cut Hypotube improves pushability.
PTFE Coating to reduce friction¹

Red tip provides tip flexibility¹
and visibility when loading a guidewire

Dual-layer balloon for optimal flexibility,
compliance and controlled balloon growth¹



Promus ELITE™ Stent System


New Proprietary Laser-Cut Hypotube

Boston
Scientific

Outstanding Deliverability

NEW Laser-Cut Hypotube

- Variable cut pattern specifically designed for coronary anatomy
- ~300 cuts over 100 mm length
- Extends from midshaft to exit port to improve pushability
- Additional length maintains midshaft flexibility
- ↑ Pushability and Flexibility



End Tab

Laser Cut Section

Promus ELITE™ Stent System

Improved Pushability with Hypotube Technology

Boston
Scientific

Outstanding Deliverability

Pushability¹

+15%
Improvement

Promus PREMIER™
Stent System

Promus ELITE
Stent System

g/cm


Push Force

(How far a catheter shaft axially compresses at a fixed force)


Promus ELITE™ Stent System

Outstanding Deliverability


Pushability (g/cm)




Lesion Entry Profile (mm)



Trackability (gf.cm)



System Flex (N*mm)




Promus ELITE™ Stent System

Low Stent and Tip Profile¹

Boston
Scientific

Outstanding Deliverability



Promus ELITE™ Stent System

Outstanding Acute Performance. Proven Long-term Outcomes

Boston
Scientific



Outstanding Deliverability



Trusted Stent Platform



Proven Clinical Outcomes

Customized Stent Architecture

Ideal Balance of Strength and Flexibility

Trusted Stent Platform




KEY ADVANTAGES¹

- Exceptional Conformability
- Unmatched Radial Strength
- Most Visible Alloy
- Minimum Recoil
- Market-Leading Drug and Polymer combination

Platinum Chromium Alloy


Unique Properties and Benefits

Trusted Stent Platform




- Platinum has over twice the density of Iron or Cobalt
- Platinum provides increased strength when alloyed with stainless steel
- Lowest Nickel content (9%)
- Specifically developed for coronary stents
- Enhanced Visibility


Lesion coverage and Bailout Stenting




Platform flexibility and conformability



Incomplete stent apposition




Resistance to fracture




Platinum Chromium Alloy

Unique Properties and Benefits

In Vivo Radiopacity¹



Best-in-Class Visibility²



Promus ELITE™ Stent

Customized Stent Architecture


Trusted Stent Platform

Strength and Flexibility Where It Matters

Additional connectors on
proximal two segments
(4 or 5)



Proximal end more robust to provide
increased axial strength¹







Overall design maintains ***flexibility,
conformability and fracture resistance¹***

Customized Stent Architecture

Four Unique, Customized Stent Models

Trusted Stent Platform

Customized for Strength and Flexibility Where It Matters

Stent Model	# of Peaks	# of Connectors	Labeled Post-Dilatation Limits	Proximal End ▼	Distal End ▼
Small Vessel (2.25 mm)	8	2 throughout	2.75		
Small Workhorse (2.50-2.75 mm)	8	4 on proximal end; 2 throughout stent body	3.50		
Workhorse (3.00-3.50 mm)	8	4 on proximal end; 2 throughout stent body	4.25		
Large Vessel (4.00mm)	10	5 on proximal end; 2 throughout stent body	5.75		


Promus ELITE™ Stent

Customized Design for Outstanding Conformability

In bench testing Promus ELITE stent is more than:

- **2x** more conformable than Resolute Onyx
- **8x** more conformable than Xience Sierra

Conformability Test¹



(Amount of torque required to bend the stent; measures the ability of the stent to naturally conform to the vessel)


Outstanding Conformability

Customized Design for Outstanding Conformability

Trusted Stent Platform


Hands-on model demonstrates conformability and the way stents conform to the natural vessel shape

Promus ELITE Stent¹




Xience™ Stent¹

Promus ELITE Stent¹



Resolute Onyx Stent¹

Outstanding conformability²




Promus ELITE™ Stent

Exceptionally Low Recoil

Boston
Scientific

Trusted Stent Platform




Promus ELITE™ Stent


Excellent Radial Strength

In bench testing Promus ELITE stent shows:

- **8%** more radial strength than Resolute Onyx
- **24%** more radial strength than Xience Sierra



Radial Compression Test




Stent


Promus ELITE™ Stent System

Trusted Stent Platform


Conformability (N.mm)



Radial Strength (N/mm)



In-Vitro Stent Recoil %



Stent Platform Specifications

Trusted Stent Platform

**Promus ELITE™
Stent**



**Resolute Onyx™
Stent**



**Xience Sierra™
Stent**



Stent Material	Platinum Chromium (PtCr)	Cobalt Nickel (CoNi) w/ Platinum Core wire	Cobalt Chromium (CoCr)
Polymer	PVDF	BioLinx™	PVDF
Drug	Everolimus	Zotarolimus	Everolimus
Drug Elution Time	4 Months	6 Months	4 Months
Strut Thickness	0.0032" (81 µm)	0.0032" (81 µm)	0.0032" (81 µm)

Promus ELITE™ Stent

Low Balloon Overhang

Trusted Stent Platform

Inflated to 1114.575 kPa (11 ATM)¹

Promus ELITE™
Stent System



0.45 mm
Total: 1.01 mm

Inflated to 1823.25 kPa (18 ATM)¹

Promus ELITE™
Stent System



0.53 mm
Total: 1.14 mm

Minimal
balloon overhang to
help minimize vessel
trauma or damage
outside the stent¹

Promus ELITE™ Stent

Drug and Polymer Properties

Trusted Stent Platform

Market-Leading Everolimus Drug and PVDF-HFP Polymer Maintains Excellent Drug Distribution and Uniformity

	PtCr Everolimus Stents	
	Promus ELITE Stent	Promus PREMIER Stent
Drug	Everolimus	Everolimus
Drug Release ¹	100% by 120 days	100% by 120 days
Minimum Drug Load ¹	38 µg (2.25x8 mm)	38 µg (2.25x8 mm)
Maximum Drug Load ¹	243 µg (4.00x38 mm)	243 µg (4.00x38 mm)
Drug Dose Density ¹	100 µg/cm ²	100 µg/cm ²
Drug to Polymer Ratio ¹	1:4.9	1:4.9
Polymer	Fluorinated Co-polymer (PVDF – HFP)	Fluorinated Co-polymer (PVDF – HFP)
Total Coating Thickness ¹	7µm	7µm




Promus ELITE™ Stent

Nominal and Rated Burst Pressures*

Trusted Stent Platform

Promus ELITE Stent has exceptional compliance and low balloon growth




Promus ELITE™ Stent

Post Dilatation Limits*

Trusted Stent Platform

Promus ELITE Stent has a labeled over expansion limit up to 5.75 mm



Promus ELITE™ Stent System

Outstanding Acute Performance. Proven Long-term Outcomes

Boston
Scientific




Outstanding Deliverability



Trusted Stent Platform



Proven Clinical Outcomes



Proven Clinical Outcomes



Promus PtCr EES¹



QCA study

100 patients

30-day clinical and 9-month
angiographic & IVUS data



Performance
Endpoints Met



Workhorse Trial

1,530 patients

Non-inferiority vs. Xience V™
(PROMUS™) Stent
5-yr data



Primary
Endpoint Met



Small Vessel Trial

94 patients

Performance goal based on
TAXUS™ Express™ Stent
4-yr data



Primary
Endpoint Met



Long Lesion Trial

102 patients

Performance goal based on
TAXUS Express Stent
4-yr data



Primary
Endpoint Met

Promus ELITE™ Stent

PLATINUM Clinical Trials Key Takeaways¹

Proven Clinical Outcomes

Promus PtCr EES BEATS Xience CoCr EES²

Numerically lower event rates at 5 years²

Impressive results with 0% Incomplete Stent Apposition.³

Exceptional safety and efficacy in Small Vessels with just 2.5% TLR and 0% ARC ST as well as in Long Lesions with 0% ARC ST or MI.⁴


1. PLATINUM Clinical Trial Program studies the PROMUS Element Stent (Promus PtCr EES) and Xience V Stent (Xience CoCr EES). The principal safety and effectiveness information for the Promus PREMIER ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus PREMIER ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus PREMIER ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus PREMIER ELITE Stent System.
2. PLATINUM Workhorse Trial; presented by Stone G, MD. ACC 2015.
3. PLATINUM QCA Trial; Meredith, et al. Eurointervention 2011;7:84.
4. PLATINUM Small Vessel Trial; presented by Dominic Alocco, MD, PCR 2012. Ian T. Meredith, AM, MBBS, PhD is the PI. PLATINUM Long Lesion Trial; presented by Paul Teirstein, TCT 2012.

Promus ELITE™ Stent

PLATINUM Workhorse Trial 5-Year Results

Proven Clinical Outcomes

Numerically Lower Event Rates¹



1. Stone, GW. PLATINUM Workhorse Trial. ACC 2015. PLATINUM Clinical Trial Program studies the PROMUS Element Stent (Promus PtCr EES) and Xience V Stent (Xience CoCr EES). The principal safety and effectiveness information for the Promus PREMIER ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus PREMIER ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus PREMIER ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus PREMIER ELITE Stent System.


‡ Deaths due to unknown causes were adjudicated as cardiac death.

Promus ELITE™ Stent

PLATINUM Workhorse Trial 5-Year Results

Proven Clinical Outcomes

Numerically Lower Ischemia-Driven TLR through 5 Years¹




1. Stone, GW. PLATINUM Workhorse Trial. ACC 2015. PLATINUM Clinical Trial Program studies the PROMUS Element Stent (Promus PtCr EES) and Xience V Stent (Xience CoCr EES). The principal safety and effectiveness information for the Promus PREMIER ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus PREMIER ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus PREMIER ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus PREMIER ELITE Stent System.

Promus ELITE™ Stent

PLATINUM Workhorse Trial 5-Year Results

Proven Clinical Outcomes



1. PLATINUM Clinical Trial Program studies the PROMUS Element Stent (Promus PtCr EES) and Xience V Stent (Xience CoCr EES). The principal safety and effectiveness information for the Promus PREMIER ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus PREMIER ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus PREMIER ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus PREMIER ELITE Stent System. PLATINUM Trial: Stone, et al. JACC 2011; 57:1700.


Promus ELITE™ Stent

Outstanding Conformability


Boston
Scientific

Trusted Stent Platform

Outstanding conformability²



Significantly less vessel straightening with
Promus PtCr EES in the
PLATINUM Workhorse Trial³



1. Image courtesy of John Ormiston, MD., Mercy Hospital Auckland, New Zealand. Results from case studies are not predictive of results in other cases. Results in other cases may vary.


2. In severely angulated lesions only. Popma J, MD. Stent Design Impacts Geometric Vessel Distortion following Coronary Artery Stenting in Severely Angulated Lesions: Angiographic Analysis of the PLATINUM Workhorse Trial. ACC 2013.

Promus ELITE™ Stent

2-Year Results in Perspective


Resolute US Trial²

Resolute CoNi ZES (n = 968)
(2.50-3.50 mm cohort)
RVD = 2.63 mm³
Lesion Length = 13.06 mm³




SPIRIT III Trial⁴

Xience CoCr EES (n = 642)
RVD = 2.77 mm
Lesion Length = 14.7 mm



PLATINUM Workhorse Trial¹

Promus PtCr EES (n = 758)
RVD = 2.67 mm
Lesion Length = 13.0 mm



1. PLATINUM Clinical Trial Program studies the PROMUS Element Stent (Promus PtCr EES) and Xience V Stent (Xience CoCr EES). The principal safety and effectiveness information for the Promus PREMIER ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus PREMIER ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus PREMIER ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus PREMIER ELITE Stent System. Presented by Gregg W. Stone, MD, ACC 2012.
2. RESOLUTE US Trial studied the Resolute™ Stent (Resolute CoNi ZES). Presented by Laura Mauri MD, MSc; ACC 2012.
3. Resolute Integrity Stent System DFU. Results from different studies are not directly comparable. Information provided for educational purposes only.
4. SPIRIT III Trial studied the Xience V™ Stent. Stone, et al. Circulation 2009; 119: 680-686. The SPIRIT Clinical Trials are sponsored by Abbott.

Promus ELITE™ Stent

Small Vessel Trial 2-Year Results in Perspective

Boston
Scientific

No ARC def/prob ST
No Myocardial Infarction

1.4% TV-MI
1.4% ARC def/prob ST


PLATINUM Small Vessel Trial^{1,3}

Promus PtCr EES

(n = 89)

RVD = 2.04 mm

Lesion Length = 14.15 mm




RESOLUTE US Trial²

Resolute CoNi ZES

(n = 146)

RVD = 2.15 mm

Lesion Length = 12.40 mm



1. PLATINUM Clinical Trial Program studies the PROMUS Element Stent (Promus PtCr EES) and Xience V Stent (Xience CoCr EES). The principal safety and effectiveness information for the Promus PREMIER ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus PREMIER ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus PREMIER ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus PREMIER ELITE Stent System. PLATINUM Small Vessel Trial, Presented by Dominic Allocco MD, PCR 2012. Ian T. Meredith, AM, MBBS, PhD is the PI. There were no MIs in the PLATINUM Small Vessel Trial.

2. Presented by Laura Mauri MD, MSc, ACC 2012. RESOLUTE US Trial studied the Resolute Stent. Results from different studies are not directly comparable. Information provided for educational purposes only.

3. PROMUS Element Stent has a dedicated small vessel 2.25 mm stent model.

Promus ELITE™ Stent

Long Lesion Trial 2-Year Results in Perspective


Boston
Scientific

No ARC def/prob ST
No Myocardial Infarction

4.8% TV-MI


PLATINUM Long Lesion Trial¹

Promus PtCr EES
(n = 100)
RVD = 2.56 mm
Lesion Length = 24.38 mm



SPIRIT PRIME Long Lesion Registry²

Xience CoCr EES
(n = 104)
RVD = 2.8 mm
Lesion Length = 26.1 mm




1. PLATINUM Clinical Trial Program studies the PROMUS Element Stent (Promus PtCr EES) and Xience V Stent (XienceCoCr EES). The principal safety and effectiveness information for the Promus PREMIER ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus PREMIER ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus PREMIER ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus PREMIER ELITE Stent System. PLATINUM Long Lesion Trial, Presented by Teirstein, P, MD. TCT 2012.
2. SPIRIT Prime Trial, Presented by Costa, M et al. TCT 2012. SPIRIT PRIME Trial studied the Xience Prime Stent. Results from different studies are not directly comparable. Information provided for educational purposes only.
3. Deaths due to unknown causes were adjudicated as cardiac death.


Promus ELITE™ Stent

Incomplete Stent Apposition Results in Perspective

Post-Procedure Incomplete Stent Apposition



Late Incomplete Stent Apposition *In perspective*




1. PLATINUM Clinical Trial Program studies the PROMUS Element Stent (Promus PtCr EES) and Xience V Stent (Xience CoCr EES). The principal safety and effectiveness information for the Promus PREMIER ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus PREMIER ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus PREMIER ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus PREMIER ELITE Stent System. PLATINUM QCA Trial; Meredith, et al. Eurointervention 2011;7:84. Performance goal based on data from SPIRIT III.

2. SPIRIT III: Stone, et al. JAMA. 2008;299:1903. The SPIRIT III Trial studied the Xience V™ Stent (Xience CoCr EES). Results from different studies are not directly comparable. Information provided for educational purposes only.

Promus ELITE™ Stent

Incomplete Stent Apposition Results in Perspective



1. PLATINUM QCA Trial; Meredith, et al. Eurointervention 2011;7:84. PLATINUM Clinical Trial Program studies the PROMUS Element Stent (Promus PtCr EES) and Xience V Stent (XienceCoCr EES). The principal safety and effectiveness information for the Promus PREMIER ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus PREMIER ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus PREMIER ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus PREMIER ELITE Stent System. Results from different studies are not directly comparable. Information provided for educational purposes only.

2. ENDEAVOR II: Endeavor DFU.

3. SPIRIT III: Stone, et al. JAMA. 2008;299:1903. SPIRIT Clinical Trials are sponsored by Abbott.


4. Waseda et al. Circ Journal 2010 74(10): 2097-2102

Promus ELITE™ Stent

PtCr Everolimus-Eluting Clinical Experience

Promus PtCr EES Demonstrated Significantly Lower MACE and Death Compared to other Permanent Polymer DES

CLINICAL EVENTS AT 12 months – REWARDS Premier TLX Registry



1. Michael A. Gaglia, Jr., MD. Presented at CRT 2016. The REWARDS Premier TLX Registry is multicenter, retrospective registry to collect baseline, clinical, procedural, in hospital and 9-12 month follow-up data to compare major adverse events cardiac events in patients receiving Promus Premier drug eluting stent to data already collected from the REWARDS-TLX Registry in which patients received either the Taxus Liberté or Xience V drug eluting stent.


Promus ELITE™ Stent

PtCr Everolimus-Eluting Clinical Experience

Proven Clinical Outcomes

Incidence of major adverse cardiovascular events of Promus PtCr EES compared to Xience V™ (CoCr) Stent – REWARDS Premier TLX Registry

Statistically Significant Lower TVR-MACE† at 1 Year¹



Promus PtCR EES demonstrated a significant reduction in major adverse cardiac events¹

Adapted from CRT 2015 Presentation by Michael A. Gaglia, Jr., MD, MSc, FSCAI

1 – Results from the REWARDS Premier TLX Registry. Presented by Michael A. Gaglia, Jr., MD, MSc, FSCAI at CRT 2015. The primary objective was to compare Promus Premier DES (PtCr) to Xience V (CoCr) and Taxus Liberte in regards to the incidence of major adverse cardiovascular events at 1 year after percutaneous coronary intervention (PCI).

† – TVR-MACE is a composite endpoint of all-cause death, Q-wave MI, and TVR


Promus ELITE™ Stent

Clinical Safety

Kang Network Meta-Analysis¹


Promus PtCr EES ranked #2 for the lowest relative risk of Def/Prob Stent Thrombosis

Promus PtCr EES vs. comparators



SCAAR Registry²

Promus PtCr EES reported numerically lowest Permanent Polymer ST rates in real-world SCAAR registry



1. Caution should be taken in interpreting study results, as some stents had limited numbers of comparisons, and some of the studies had a potential risk of bias. *All Promus PtCr-EES stents, also includes PROMUS Element, PROMUS Element Plus and Promus PREMIER. Def/prob ST was available in 110 studies with 111,088 patients.

Kang, S. et al. J Am Coll Cardiol Interv. doi:10.1016/j.jcin.2016.03.038

2. Adapted from Presentation by Stefan James, MD at Euro PCR 2017 (Swedish SCAAR database, real-world outcomes of latest DES generations). Promus PtCr EES PP Stent Family includes the Promus PREMIER PP Stent.

Promus ELITE™ Stent

Stent Thrombosis Network Meta-Analysis Published in *The Lancet*



Significantly lower risk of ARC ST (Def/Prob) with Everolimus-eluting Stents compared to Zotarolimus-eluting or Bare-Metal stents in a **LARGE Network Meta-Analysis reported¹**

1-Year ARC ST (Def/Prob)

PROMUS®/Xience V® Stent vs BMS

PROMUS/Xience V Stent vs Endeavor™ Stent

2-Year ARC ST (Def/Prob)

PROMUS/Xience V Stent vs BMS

PROMUS/Xience V Stent vs Endeavor Stent


Early ARC ST (Def/Prob)

PROMUS/Xience V Stent vs BMS

PROMUS Element™ Stent vs BMS

PROMUS/Xience V Stent vs Endeavor Stent

PROMUS Element Stent vs Endeavor Stent




Late ARC ST (Def/Prob)

PROMUS/Xience V Stent vs BMS

PROMUS/Xience V Stent vs Endeavor™ Stent

PROMUS/Xience V Stent vs Resolute Integrity™ Stent



Promus ELITE™ Stent

Indications and Labeling

Promus ELITE™ Stent

Clinical Indications

CE Mark

- ✓ Concomitant Diabetes mellitus (DM)
- ✓ Ostial Lesions
- ✓ Unprotected Left Main Coronary Artery Lesions
- ✓ Chronic Total Occlusion¹
- ✓ In-Stent Restenosis²
- ✓ Acute Myocardial Infarction (AMI)
- ✓ Coronary Bifurcation³



¹ For treatment of occluded vessels, contrast visualization of the distal vessel to confirm position of guidewire within the lumen is recommended.

² For in stent restenosis, where details of the original stent are known, the expanded inner diameter of the new stent should not exceed the dilation limits of the original stent. Where details of the original stent are not known, the expanded inner diameter of the new stent should not exceed the reference vessel diameter.

³ When treating Bifurcations, care must be exercised to access the secondary vessel via the repeating geometry in the body of the stent within the primary vessel.


Promus ELITE™ Stent

Left Main Stenting

- **Left Main Vessels are Large (>5.5mm on average¹)**
 - Promus ELITE Stent is indicated for treatment of patients presenting with unprotected left main coronary artery lesions
 - Promus ELITE has a labeled overexpansion of up to 5.75mm²
- **Left Main Vessels require stents with significant radial strength**
 - Promus ELITE offers excellent radial strength³



IVUS LM in male patients



Graph 2. IVUS data on 865 LAD->LM pullbacks, male patients (Rotterdam). Average vessel Diameter (VD): 5.60 mm , Mean VD: 5.47 mm. Only in 13% <5 mm.



1. Shand J, et al. Prospective Intravascular Ultrasound Investigation of the Necessity for and Efficacy of Postdilatation Beyond Nominal Diameter of 3 Current Generation DES Platforms for the Percutaneous Treatment of the Left Main Coronary Artery. Cathet Cardiovasc Interv; 2014;84:351-358.

2. Labeled Post-Dilatation Limits. Promus PREMIER Stent, Xience Xpedition Stent, Resolute Integrity Stent and Resolute Onyx DFU

3. Bench testing performed by Boston Scientific Corporation. Data on file at Boston Scientific. All stents 2.5mm. Promus PREMIER Stent n=5, all other stents are n=3. Bench test results not necessarily indicative of clinical performance.



The expanded indications (bifurcation, ISR, LM, OL and CTO) are for CE Mark countries only

Promus ELITE™ Stent

CTO Stenting



- **CTO lesions tend to be long and in large vessels**
 - Promus ELITE has a labeled **overexpansion** of **up to 5.75mm**, allowing the physician to customize the stent to the appropriate vessel size¹
- **Stent mechanical properties are important considerations for CTO stenting**
 - The Promus ELITE Stent features a customized architecture offering exceptional **strength** and **conformability**²



For treatment of occluded vessels with the Promus ELITE stent system, contrast visualization of the distal vessel to confirm position of guidewire within the lumen is recommended.

1. Labeled Post-Dilatation Limits. Promus PREMIER Stent, Xience Xpedition Stent, Resolute Integrity Stent and Resolute Onyx DFU

2. Bench testing performed by Boston Scientific Corporation. Data on file at Boston Scientific. All stents 2.5 mm. Promus PREMIER Stent and PROMUS Element Stent n = 15, all other stents n = 3.

Bench test results not necessarily indicative of clinical performance. The expanded indications (bifurcation, ISR, LM, OL and CTO) are for CE Mark countries only

Promus ELITE™ Stent

Stenting Bifurcation

Boston
Scientific



- **Bifurcation lesions are often tapered¹**
 - Promus ELITE has a labeled overexpansion of up to 5.75mm, allowing the physician to customize the stent to the appropriate vessel size and ensure great apposition²
- **With bifurcation we want to maintain the natural vessel shape¹**
 - The Promus ELITE Stent's customized architecture offers exceptional strength and conformability³
- **Bifurcation PCI benefits from an appropriate Cell Diameter & Expansion¹**
 - The Promus ELITE Stent has large cell diameters in the body of the stent to accommodate side branch access⁴

When treating Bifurcations, care must be exercised to access the secondary vessel via the repeating geometry in the body of the stent within the primary vessel.





¹ Attributes collected from: Curtiss T. Stinis, M.D., F.A.C.C, F.S.C.A.I. -Scripps CI 2013 and ESC/EACTS GUIDELINES 2014. EHJ doi:10.1093/eurheartj/ehu278. ² Labeled Post-Dilatation Limits. Promus PREMIER Stent, Xience Xpedition Stent, Resolute Integrity Stent and Resolute Onyx DFU ³ Bench testing performed by Boston Scientific Corporation. Data on file at Boston Scientific. All stents 2.5mm. Promus PREMIER Stent n=5, all other stents are n=3. Bench test results not necessarily indicative of clinical performance.

⁴ Circular Cell Diameter: n = 5. Data on file at Boston Scientific. Size Matrix on the Promus PREMIER Stent DFUs.

The expanded indications (bifurcation, ISR, LM, OL and CTO) are for CE Mark countries only

Promus ELITE™ Stent

Side Branch Access

	2.25mm	2.50- 2.75mm	3.00- 3.50mm	4.00
Maximum Expanded Cell Diameter (MECD) in Stent Body (mm)				
Circular Cell Diameter (CCD) in Stent Body (mm)	0.63	0.75	0.91	1.06
Cell Perimeter in Stent Body (mm)	13.15	14.75	18.01	23.48
Ratio of Proximal Cell Perimeter to Body Cell Perimeter	1	0.5	0.5	0.4



Individualisation of Patient Treatment


Antiplatelet drugs should be used in combination with (Promus ELITE) drug-eluting stents. Physicians should use the information from the large body of clinical evidence for everolimus-eluting stents, coupled with current literature on drug-eluting stents, current European Society of Cardiology recommendation (or other applicable country guidelines) and the specific needs of the individual patient to determine the specific antiplatelet / anticoagulation regimen to be used for their patients in general practice.

It is very important that the patient be compliant with post-procedural antiplatelet recommendations given by their physician. In selected higher risk patients where the physician determines that the risks outweigh the benefits of continued DAPT, it may be reasonable to interrupt or discontinue therapy after 1 month of DAPT based on the low stent thrombosis rates and no observed increased risk of stent thrombosis demonstrated in the current literature. Patients who require premature discontinuation of antiplatelet therapy should be monitored closely and have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physician.

Misc Product Information

Promus ELITE™ Stent System

Key Stent Delivery System Specifications



	Monorail® Stent Delivery System
Lesion Entry Profile	0.045 mm
Average Stent Profile	1.07 mm
Guide Catheter Compatibility	5F ≥1.42 mm (0.056")
Y-Adapter Ports	Single access port to inflation lumen Guidewire exit port located 26 cm from tip Designed for guidewire ≤ 0.014" (0.36 mm)
Balloon Fold	5 wings
Marker Band Material and Length	Platinum Iridium; 1 mm
Balloon Material	PEBAX
Marker Band Placement	Nominally placed 0.4 mm (0.016") beyond stent on each end
Distal Shaft Coating	Bioslide™ hydrophilic coating
Proximal Shaft Hypotube Coating	PTFE
Shelf Life	24 months

Promus ELITE Stent System

Compliance Chart (Inner Diameter*)

Pressure Atm - kPa	Stent I.D. (mm)					
	2.25	2.50	2.75	3.00	3.50	4.00
8 - 814	---	2.29	2.50	2.72	3.24	3.72
9 - 910	2.13	2.37	2.58	2.81	3.34	3.81
10 - 1014	2.19	2.43	2.65	2.88	3.43	3.89
11 - 1117	2.24	2.50	2.72	2.95	3.51	3.96
12 - 1213	2.29	2.55	2.78	3.01	3.58	4.02
13 - 1317	2.34	2.60	2.84	3.06	3.63	4.08
14 - 1420	2.38	2.65	2.89	3.10	3.68	4.13
15 - 1517	2.42	2.68	2.93	3.14	3.73	4.17
16 - 1620	2.45	2.72	2.96	3.17**	3.77**	4.21**
17 - 1724	2.47	2.75	2.99	3.20	3.81	4.25
18 - 1827	2.50**	2.77**	3.03**	3.24	3.85	4.30
19 - 1924	2.52	2.80	3.06	3.28	3.91	4.36
20 - 2027	2.55	2.83	3.09	3.32	3.97	4.43
21 - 2130	2.57	2.87	3.13	---	---	---
22 - 2227	2.60	2.90	3.17	---	---	---

NOMINAL

RBP

Diameter (mm)	Length (mm)							
2.25	8	12	16	20	24	28	32	n/a
2.50	8	12	16	20	24	28	32	38
2.75	8	12	16	20	24	28	32	38
3.00	8	12	16	20	24	28	32	38
3.50	8	12	16	20	24	28	32	38
4.00	8	12	16	20	24	28	32	38

Promus ELITE™ Stent System

Outstanding Acute Performance. Proven Long-term Outcomes

Boston
Scientific




Outstanding Deliverability



Trusted Stent Platform



Proven Clinical Outcomes



Thank You



Interventional Cardiology
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

*To order product or for more information
contact customer service at 1.888.272.1001.*

© 2018 Boston Scientific Corporation
or its affiliates. All rights reserved.

Caution: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information contained herein is for use or distribution outside the US, France and Japan only.