CERTIFICATE OF REGISTRATION

Lorne Laboratories Ltd

Unit 1 Cutbush Park Industrial Estate
Danehill
Lower Earley
Berkshire RG6 4UT UNITED KINGDOM

UL LLC® (UL Solutions) issues this certificate to the Firm named above, after assessing the Firm's quality system and finding it in compliance with:

ISO 13485:2016
EN ISO 13485:2016

The design and manufacture of in vitro diagnostic blood grouping reagents. The purchase for resale of in vitro diagnostic serology test kits.

Authorized by

Paul Hilgeman
Senior Business Manager - Medical
CMIT – Medical Regulatory

Check Certificate Status:
[here](#)

File Number	A12241	Cycle Start	May 23, 2023
Certificate Number	1458.230523	Effective Date	May 23, 2023
Initial Issue Date	June 26, 2018	Expiry Date	May 22, 2026

This quality system registration is included in UL's Directory of Registered Firms and applies to the provision of goods and/or services as specified in the scope of registration from the address(es) shown above. By issuance of this certificate the firm represents that it will maintain its registration in accordance with the applicable requirements. This certificate is not transferable and remains the property of UL LLC® (UL Solutions).

UL Solutions
333 Pfingsten Road
Northbrook, IL 60062-2096 USA

CERTIFICATE

EC Certificate No. 1434-IVDD-074/2022

**EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices**

Polish Centre for Testing and Certification certifies
that manufactured by:

**Lorne Laboratories Ltd
Unit 1 Cutbush Park Industrial Estate
Danehill
Lower Earley
Berkshire RG6 4UT, UNITED KINGDOM**

i.e. *in vitro* diagnostic medical devices
List A

The list of medical devices covered by this certificate is provided in the Annex 1

in terms of design documentation, comply with requirements
of Annex IV (Section 4) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 28.04.2022 to 27.05.2025

The date of issue of the Certificate: 28.04.2022

The date of the first issue of the Certificate: 10.04.2019

Issued under the Contract No. MD-004/2022
Application No: 504/2022
Certificate bears the qualified signature.
Warsaw, 28/04/2022
Module H6/V1

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Digitally signed by
Aleksandra
Kostrzewa

President

ANNEX 1 TO THE CERTIFICATE

VALID ONLY WITH CERTIFICATE

No 1434-IVDD-074/2022

List of medical devices covered by the certificate:

Anti-A Monoclonal 600010
Anti-B Monoclonal 610010
Anti-A,B Monoclonal 620010
Anti-D Clone 1 Monoclonal 730010
Anti-D Clone 2 Monoclonal 710010
Anti-D Duoclone Monoclonal 740010
Anti-C Monoclonal 690005
Anti-E Monoclonal 691005
Anti-c Monoclonal 692005
Anti-e Monoclonal 693005
Anti-C+D+E Monoclonal 700010
Anti-K Monoclonal 760010

Issued under the Contract No. MD-004/2022
Application No: 504/2022
Certificate bears the qualified signature.
Warsaw, 28/04/2022

Aleksandra
Kostrzewa
President

Digitally signed by
Aleksandra
Kostrzewa

CERTIFICATE

EC Certificate No. 1434-IVDD-075/2022
Full Quality Assurance System
Directive 98/79/EC concerning
***in vitro* diagnostic medical devices**

Polish Centre for Testing and Certification certifies
that the quality assurance system in the organization:

Lorne Laboratories Ltd
Unit 1 Cutbush Park Industrial Estate
Danehill
Lower Earley
Berkshire RG6 4UT, UNITED KINGDOM

for the design, manufacture and final inspection of *in vitro* diagnostic medical device
List A

*The list of medical devices covered by this certificate is provided
in the Annex 1 to EC Design-examination Certificate No. 1434-IVDD-074/2022*

complies with requirements
of Annex IV (excluding Section 4, 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 28.04.2022 to 27.05.2025

The date of issue of the Certificate: 28.04.2022

The date of the first issue of the Certificate: 10.04.2019

Issued under the Contract No. MD-004/2022
Application No: 505/2022
Certificate bears the qualified signature.
Warsaw, 28/04/2022
Module H7

Aleksandra Kostrzewa
Digitally signed by
Aleksandra
Kostrzewa
President

CERTIFICATE

EC Certificate No. 1434-IVDD-027/2022

**Full Quality Assurance System
Directive 98/79/EC concerning
in vitro diagnostic medical devices**

Polish Centre for Testing and Certification certifies
that the quality assurance system in the organization:

Lorne Laboratories Ltd

**Unit 1 Cutbush Park Industrial Estate, Danehill,
Lower Earley, Berkshire RG6 4UT, UNITED KINGDOM**

for the design, manufacture and final inspection of *in vitro* diagnostic medical device
List B

The list of medical devices covered by this certificate is provided in the Annex 1

complies with requirements
of Annex IV (excluding Section 4, 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 23.03.2022 to 27.05.2025

The date of issue of the Certificate: 03.03.2022

The date of the first issue of the Certificate: 10.04.2019

Issued under the Contract No. MD-173/2021
Application No: 577/2022
Certificate bears the qualified signature.
Warsaw, 03/03/2022
Module H7

Aleksandra Kostrzewa Digitally signed
by Aleksandra
Kostrzewa

President

ANNEX 1 TO THE CERTIFICATE

VALID ONLY WITH CERTIFICATE

No 1434-IVDD-027/2022

List of medical devices covered by the certificate:

Anti-Jka Polyclonal 323002

Anti-Jkb Polyclonal 324002

Anti-Fyb Polyclonal 317002

AHG Elite Clear 415010

AHG Elite Green 435010

Anti-Fya Monoclonal 774002

Anti-Human IgG Clear 401010

Anti-Human IgG Green 402010

Anti-Jka Monoclonal 775002

Anti-Jkb Monoclonal 776002

CE 1434

Issued under the Contract No. MD-173/2021
Application No: 577/2021
Certificate bears the qualified signature.
Warsaw, 03/03/2022

Aleksandra Kostrzewa Digitally signed
by Aleksandra Kostrzewa

President

EC DECLARATION OF CONFORMITY

Lorne Laboratories Ltd declares that the following in vitro diagnostic reagent:

Product Name	Catalogue Number
Anti-S Monoclonal	770002

has been classified as non List A, non List B (Directive 98/79/EC, Annex II) and complies with the essential requirements and provisions of Directive 98/79/EC of the European Parliament and of the Council (also SI 2002 No.618 which transposes the requirements of Directive 98/79/EC).

and is in conformity with the national standards transposing harmonised standards:

- BS EN ISO 13485:2012
- BS EN 13612:2002
- BS EN 13641:2002
- BS EN ISO 14971:2012
- BS EN ISO 15223-1:2016
- BS EN ISO 18113-2:2011
- BS EN ISO 23640:2015

The conformity assessment procedure performed was in accordance with Annex III of Directive 98/79/EC.

This declaration of conformity is issued under the sole responsibility of Lorne Laboratories Ltd and is valid from 06 April 2017.

Eddy Velthuis
Technical Director

DIRECTIONS FOR USE

Anti-A, Anti-B and Anti-A,B Monoclonal:

For Tube, Bio-Rad-ID, Ortho BioVue, Microplate and Slide Techniques.

SUMMARY

In 1900, Landsteiner discovered the serum of some people would agglutinate the red cells of others. Four common phenotypes are now recognised: O, A, B and AB. Subgroups of A and B have since been identified.

Forward Group			Reverse Group			ABO Phenotype	Caucasians % ¹	
A	B	A,B	A ₁	A ₂	B	O		
+	0	+	0	0	+	0	A	43
0	+	+	+	+	0	0	B	9
0	0	0	+	+	+	0	O	44
+	+	+	0	0	0	0	AB	4

INTENDED PURPOSE

The ABO reagents are blood grouping reagents intended to be used to qualitatively determine the presence or absence of the A and/or B antigens on the red cells of blood donors or patients requiring a blood transfusion when tested in accordance with the recommended techniques stated in this IFU.

PRINCIPLE

The reagents contain antibodies against the appropriate A and/or B antigen on human red cells and will cause direct agglutination (clumping) of red cells that carry the corresponding ABO antigen. No agglutination generally indicates the absence of the corresponding ABO antigen on human red cells (see **Limitations**).

REAGENT

Lorne Monoclonal IgM ABO blood grouping reagents contain mouse monoclonal antibodies diluted in a phosphate buffer containing sodium chloride, EDTA and bovine albumin. The reagents do not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. Each reagent is supplied at optimal dilution for use with all the recommended techniques stated below without the need for further dilution or addition. For lot reference number and expiry date see **Vial Label**.

Product	Cell Line/Clone	Colour	Dye Used
Anti-A	9113D10	Blue	Patent Blue
Anti-B	9621A8	Yellow	Iartrazine
Anti-A,B	152D12 + 9113D10 + ES15	Colourless	None

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Blood samples can be collected into EDTA, citrate, CPDA anticoagulants or as a clotted sample. The samples should be tested as soon as possible following collection. If a delay in testing should occur, store the samples at 2-8°C. Samples displaying gross haemolysis or microbial contamination should not be used for testing. Blood samples showing evidence of lysis may give unreliable results. It is preferable (but not essential) to wash all blood samples with PBS or isotonic saline before being tested.

PRECAUTIONS

- The reagents are intended for in vitro diagnostic use only.
- If a reagent vial is cracked or leaking, discard the contents immediately.
- Do not use the reagents past the expiration date (see **Vial Label**).
- Do not use the reagents if a precipitate is present.
- Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
- The reagents have been filtered through a 0.2 µm capsule to reduce the bio-burden, but are not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date as long as there is no marked turbidity, which can indicate reagent deterioration or contamination.
- The reagents contain < 0.1% sodium azide. Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive metal azides. On disposal flush away with large volumes of water.
- No known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

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For information on disposal of the reagent and decontamination of a spillage site see Material Safety Data Sheets, available on request.

CONTROLS AND ADVICE

- It is recommended a positive control and a negative control be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.

- Since these reagents do not contain macromolecular potentiators, it is very unlikely that false positive reactions are caused with IgG coated cells.
- Blood specimens of weak A or B subgroups (e.g Ax) may give rise to false negative or weak reactions when tested using slides, microtitre plates or gel cards. It is advisable to re-test weak subgroups using tube technique.
- Individuals older than six months should have their ABO blood-grouping results confirmed by testing their serum or plasma against known group A, and B cells before their ABO blood group can be confirmed.
- Before use, let the reagent warm up to room temperature. As soon as the reagent has been used, put the reagent back in storage at 2-8°C.
- In the Recommended Techniques one volume is approximately 50µl when using the vial dropper provided.
- The use of the reagents and the interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of the country where the reagents are in use.
- The user must determine the suitability of the reagents for use in other techniques.

REAGENTS AND MATERIALS REQUIRED

- Applicator sticks.
- Automatic plate reader.
- Bio-Rad ID-Cards (NaCl, enzyme test and cold agglutinins).
- Bio-Rad ID-Centrifuge.
- Bio-Rad ID-CellStab or ID-Diluent 2.
- Glass microscope slides or white card tiles.
- Glass test tubes (10 x 75 mm or 12 x 75 mm).
- Microplate centrifuge.
- Ortho BioVue System Cassettes (Neutral).
- Ortho BioVue System Centrifuge.
- Ortho 0.8% Red Cell Diluent.
- Plate shaker.
- PBS solution (pH 6.8–7.2) or Isotonic saline solution (pH 6.5–7.5).
- Positive and negative control red cells:
- Anti-A: group A (positive control) and group O (negative control).
- Anti-B: group B (positive control) and group O (negative control).
- Anti-A,B: group A and group B (positive controls) and group O (negative control).
- Test tube centrifuge.
- Validated "U" well microplates.
- Volumetric pipettes.

RECOMMENDED TECHNIQUES

A. Tube Technique

- Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
- Place in a labelled test tube: 1 volume of Lorne Anti-ABO reagent and 1 volume of red cell suspension.
- Mix thoroughly and incubate at room temperature for 1 minute.
- Centrifuge all tubes for 10 seconds at 1000 rcf or for a suitable alternative time and force.
- Gently resuspend red cell button and read macroscopically for agglutination
- Any tubes, which show a negative or questionable result, should be incubated for 15 minutes at room temperature.
- Following incubation, repeat steps 4 and 5.

B. Bio-Rad-ID Technique (NaCl, enzyme test and cold agglutinins cards)

- Prepare a 0.8% suspension of red cells in ID-CellStab or ID-Diluent 2.
- Remove aluminium foil from as many microtubes as needed.
- Place in appropriate microtube: 50µl of red cell suspension and 25µl of Lorne Anti-ABO reagent.
- Centrifuge ID-Card(s) in the Bio-Rad gel card centrifuge.
- Read macroscopically for agglutination.

C. Ortho BioVue Technique (Neutral cassettes)

- Prepare a 0.8% suspension of red cells in 0.8% Ortho Red Cell Diluent.
- Remove aluminium foil from as many reaction chambers as needed.
- Place in appropriate reaction chamber: 50µl of red cell suspension and 40µl of Lorne Anti-ABO reagent.
- Centrifuge cassette(s) in an Ortho BioVue System Centrifuge.
- Read macroscopically for agglutination.

D. Microplate Technique, using "U" wells

- Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
- Place in the appropriate well: 1 volume Lorne Anti-ABO reagent and 1 volume red cell suspension.
- Mix thoroughly, preferably using a microplate shaker, taking care to avoid cross-well contamination.
- Incubate at room temperature for 15 minutes (time dependant on user).
- Centrifuge the microplate for 1 minute at 140 rcf or for a suitable alternative time and force.
- Resuspend the cell buttons using carefully controlled agitation on a microplate shaker

7. Read macroscopically or with a validated automatic reader.
8. Any weak reactions should be repeated by the tube technique.

E. Slide Technique

1. Prepare a 35-45% suspension of red cells in serum, plasma or PBS or Isotonic saline or use anti-coagulated whole blood (in its own plasma).
2. Place on a labelled glass slide or card tile: 1 volume of Lorne Anti-ABO reagent and 1 volume of red cell suspension.
3. Using a clean applicator stick, mix reagent and cells over an area of about 20 x 40 mm.
4. Slowly tilt the slide back and forth for 30 seconds, with occasional further mixing during the 1-minute period, maintaining slide at room temperature.
5. Read macroscopically after 1 minute over a diffuse light and do not mistake fibrin strands as agglutination.
6. Any weak reactions should be repeated by the tube technique.

INTERPRETATION OF TEST RESULTS

1. Positive: Agglutination of the red cells constitutes a positive test result and within accepted limitations of test procedure, indicates the presence of the appropriate ABO antigen on the red cells.
2. Negative: No agglutination of the red cells constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of the appropriate ABO antigen on the red cells.
3. Discrepancies: If the results obtained with reverse group don't correlate with forward group, further investigation is required.

STABILITY OF THE REACTIONS

1. Read all tube and microplate tests immediately after centrifugation.
2. Slide tests should be interpreted after a maximum of one minute to ensure specificity and to avoid the possibility a negative result may be incorrectly interpreted as positive due to drying of the reagent.
3. Caution should be exercised in the interpretation of results of tests performed at temperatures other than those recommended.

LIMITATIONS

1. ABO antigens are not fully developed at birth and so weaker reactions may therefore occur with cord or neonatal specimens.
2. When using Monoclonal Anti-A,B, blood specimens of weak A or B subgroups (e.g Ax) may give rise to false negative or weak reactions when tested using slides, microtitre plates or gel cards. It is advisable to re-test weak subgroups using the tube technique.
3. Lorne monoclonal Anti-A and monoclonal Anti-B are not validated to detect Ax and A3 or Bx and B3 antigens resp and we therefore do not claim reactivity of the monoclonal Anti-A or Anti-B reagent against these weak A and B sub-groups.
4. Stored blood may give weaker reactions than fresh blood.
5. False positive or false negative results may also occur due to:
 - Contamination of test materials
 - Improper storage, cell concentration, incubation time or temperature
 - Improper or excessive centrifugation
 - Deviation from the recommended techniques
 - Cord samples contaminated with Wharton's jelly

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1. Prior to release, each lot of Lorne ABO monoclonal reagent was tested using the recommended test methods listed in this IFU. The tests complied with the test requirements as stated in the current version/issue of the 'Guidelines for the Blood Transfusion Services in the United Kingdom' and the 'Common Technical Specifications'.
2. Specificity of source monoclonal antibodies is demonstrated using a panel of antigen-negative cells.
3. The potency of the reagents has been tested against the following minimum potency reference standards obtained from National Institute of Biological Standards and Controls (NIBSC): Anti-A reference standard 03/188 And / Or Anti-B reference standard 03/164
4. Lorne Anti-B does not react with "Acquired-B" red cells.
5. Lorne Monoclonal ABO reagents do not detect crypt antigens such as T, Tn or Cad.
6. The Quality Control of the reagents was performed using red cells with phenotypes that were verified by a UK blood transfusion centre and had been washed with PBS or Isotonic saline prior to use.

DISCLAIMER

1. The user is responsible for the performance of the reagents by any method other than those mentioned in the Recommended Techniques.
2. Any deviations from the Recommended Techniques should be validated prior to use⁸.

BIBLIOGRAPHY

1. Marion E. Reid and Christine Lomas-Francis, Blood Group Antigens and Antibodies, SBB Books, New York 2007; Page 181.
2. Issitt PD. Applied Blood Group Serology, 3rd Edition. Montgomery Scientific, Miami 1985; Chapter 6
3. Guidelines for the Blood Transfusion Service in the United Kingdom 6th Edition 2002. The Stationery Office.

4. AABB Technical Manual, 16th Edition, AABB 2008.
5. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

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Advena Ltd, Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013, Malta

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BIBLIOGRAFIE

1. Marion E.Reid & Christine Lomas-Francis, Blood Group Antigens & Antibodies, SBB Books, New York 2007; pagina 181.
2. Issitt PD. Applied Blood Group Serology, 3rd Edition. Montgomery Scientific, Miami 1985; capitolul 6.
3. Guidelines for the Blood Transfusion Service in the United Kingdom, 6th Edition 2002. The Stationary Office.
4. AABB Technical Manual, 16th edition, AABB 2008.
5. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

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

TABEL SIMBOLURI

	Batch Number		<i>In-vitro</i> Diagnostic
	Catalogue Reference		Store At
	Expiry Date		Manufacturer
	Read Pack Insert		

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Lorne Laboratories Limited

Unit 1 Cutbush Park Industrial Estate

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PRINCIPLE

The reagents contain antibodies against the D antigen on human red cells and will cause direct agglutination (clumping) of human red cells that carry the D antigen. No agglutination (no clumping) generally indicates the absence of the D antigen on human red cells (see **Limitations**).

REAGENT

Lorne Monoclonal IgM Anti-D Clone 1 and Clone 2 blood grouping reagents are low protein reagents containing a human monoclonal IgM antibody diluted with sodium chloride (0.9 g%), bovine albumin (2.0 g%) and macromolecular potentiators (1.5 g%). When typing patient samples, each reagent will directly agglutinate Rh D positive cells, including majority of variants (but not D^u) and a high proportion of weak D (D^w) phenotypes when using the recommended techniques. The reagents do not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. Each reagent is supplied at optimal dilution for use on patient samples with all recommended techniques stated below without need for further dilution or addition. For lot reference number and expiry date see Vial Label.

Product	Cell Line / Clone
Anti-D Clone 1	RUM-1
Anti-D Clone 2	MS-201

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The collective term D^w is widely used to describe red cells which have a weaker expression of the D antigen than normal. The term weak D denotes individuals with a reduced number of complete D antigen sites per red cell. The term partial D denotes individuals with missing D antigen epitopes. D^w cells is a partial D category which misses most D epitopes. Both Clone 1 and Clone 2 reagents will detect most examples of partial and weak D red cells by direct agglutination, but will not detect D^w cells.

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PRECAUTIONS

1. The reagents are intended for in vitro diagnostic use only.
2. If a reagent vial is cracked or leaking, discard the contents immediately.
3. Do not use the reagents past the expiration date (see Vial Label).
4. Do not use the reagents if a precipitate is present.
5. Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
6. The reagents have been filtered through a 0.2 µm capsule to reduce the bio-burden, but are not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date as long as there is no marked turbidity, which can indicate reagent deterioration or contamination.
7. The reagents contain < 0.1% sodium azide. Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive metal azides. On disposal flush away with large volumes of water.
8. Materials used to produce the products were tested at source and found to be negative for HIV 1+2 and HCV antibodies and HBsAg using approved microbiological tests.

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CONTROLS AND ADVICE

1. It is recommended a positive control (ideally R,r cells), and a negative control (ideally rr cells) be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
2. When typing red cells from a patient who is diagnosed with a disease that causes the red cells to become coated with antibody or other proteins (such as HDN, AIHA), it is important to test the patient's red cells using Lorne's Monoclonal D Negative Control (catalogue # 650010). Tests must be considered invalid if red cells are agglutinated using Lorne's Monoclonal D Negative Control (catalogue # 650010).
3. Weak and partial D antigen variants are poorly detected by the gel card, microtitre plate and slide technique. It is recommended that weak and partial D variants are tested using the tube test technique.
4. Before use, let the reagent warm up to room temperature. As soon as the reagent has been used, put the reagent back in storage at 2-8°C.
5. In the Recommended Techniques one volume is approximately 50µl when using the vial dropper provided.
6. The use of the reagents and the interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of the country where the reagents are in use.
7. The user must determine suitability of the reagents for use in other techniques.

REAGENTS AND MATERIALS REQUIRED

- Applicator sticks.
- Automatic plate reader.
- Bio-Rad ID-Cards (NaCl, enzyme test and cold agglutinins).
- Bio-Rad ID-Centrifuge.
- Bio-Rad ID-CellStab or ID-Diluent 2.
- Glass microscope slides or white card tiles.
- Glass test tubes (10 x 75 mm or 12 x 75 mm).
- Microplate centrifuge.
- Ortho BioVue System Cassettes (Neutral).
- Ortho BioVue System Centrifuge.
- Ortho 0.8% Red Cell Diluent.
- Plate shaker.
- PBS solution (pH 6.8-7.2) or isotonic saline solution (pH 6.5-7.5).
- Positive (ideally R,r) and negative (rr) control red cells.
- Test tube centrifuge.
- Validated "U" well microplates.
- Volumetric pipettes.

RECOMMENDED TECHNIQUES

A. Tube Technique

1. Prepare a 2-3% suspension of red cells in PBS or isotonic saline.
2. Place in a labelled test tube: 1 volume of Lorne Anti-D reagent and 1 volume of red cell suspension.
3. Mix thoroughly and centrifuge all tubes for 20 seconds at 1000 rcf or for a suitable alternative time and force.
4. Gently resuspend red cell button and read macroscopically for agglutination
5. Any tubes, which show negative or questionable result (as can happen with weak D samples), should be incubated for 15 minutes at room temperature.
6. Following incubation, repeat steps 3 and 4.

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1. Prepare a 0.8% suspension of red cells in ID-CellStab or ID-Diluent 2.
2. Remove aluminium foil from as many microtubes as needed.
3. Place in appropriate microtube: 50µl of red cell suspension and 25µl of Lorne Anti-D reagent.
4. Centrifuge ID-Card(s) in a Bio-Rad gel card centrifuge.
5. Read macroscopically for agglutination.

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1. Prepare a 0.8% suspension of red cells in 0.8% Ortho Red Cell Diluent.
2. Remove aluminium foil from as many reaction chambers as needed.
3. Place in appropriate reaction chamber: 50µl of red cell suspension and 40µl of Lorne Anti-D reagent.
4. Centrifuge cassette(s) in an Ortho BioVue System Centrifuge.
5. Read macroscopically for agglutination.

D. Microplate Technique, using "U" wells

1. Prepare a 2-3% suspension of red cells in PBS or isotonic saline.
2. Place in the appropriate well: 1 volume Lorne Anti-D reagent and 1 volume

- red cell suspension.
- Mix thoroughly, preferably using a microplate shaker, taking care to avoid cross-well contamination.
 - Incubate at room temperature for 15 minutes (time dependant on user).
 - Centrifuge the microplate for 1 minute at 140 rcf or for a suitable alternative time and force.
 - Resuspend the cell buttons using carefully controlled agitation on a microplate shaker
 - Read macroscopically or with a validated automatic reader.
 - Any weak reactions should be repeated by the tube technique.

E. Slide Technique

- Prepare a 35-45% suspension of red cells in serum, plasma or PBS or Isotonic saline or use anti-coagulated whole blood (in its own plasma).
- Place on a labelled glass slide or card tile: 1 volume of Lorne Anti-D reagent and 1 volume of red cell suspension.
- Using a clean applicator stick, mix reagent and cells over an area of about 20 x 40 mm.
- Slowly tilt the slide back and forth for 30 seconds, with occasional further mixing during the 1-minute period, maintaining slide at room temperature.
- Read macroscopically after 1 minute over a diffuse light and do not mistake fibrin strands as agglutination.
- Any weak reactions should be repeated by the tube technique.

INTERPRETATION OF TEST RESULTS

- Positive: Agglutination of the red cells constitutes a positive test result and within accepted limitations of test procedure, indicates the presence of the D antigen on the red cells.
- Negative: No agglutination of the red cells constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of the D antigen on the red cells.
- Test results of cells that are agglutinated using the reagent negative control shall be excluded, as the agglutination is most probably caused by the effect of the macromolecular potentiators in the reagent on sensitised cells.

STABILITY OF THE REACTIONS

- Read all tube and microplate tests immediately after centrifugation.
- Slide tests should be interpreted after a maximum of one minute to ensure specificity and to avoid the possibility a negative result may be incorrectly interpreted as positive due to drying of the reagent.
- Caution should be exercised in the interpretation of results of tests performed at temperatures other than those recommended.

LIMITATIONS

- Lorne Anti-D is not suitable for use with enzyme treated cells, cells suspended in LISS or for use in indirect antiglobulin (IAT) techniques.
- Stored blood may give weaker reactions than fresh blood.
- False positive agglutination may be seen due to the presence of macromolecular potentiators in the reagent when testing IgG sensitised cells, e.g. ALHA, HDN.
- False positive or false negative results may also occur due to:
 - Contamination of test materials
 - Improper storage, cell concentration, incubation time or temperature
 - Improper or excessive centrifugation
 - Deviation from the recommended techniques

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- Prior to release, each lot of Lorne Anti-D monoclonal reagent was tested using the recommended test methods listed in this IFU. The tests complied with the test requirements as stated in the current version/issue of the 'Guidelines for the Blood Transfusion Services in the United Kingdom' and the 'Common Technical Specifications'.
- Anti-D grouping reagents for D grouping of patients should not react with D⁰ cells using the method(s) recommended for use.
- Specificity of source monoclonal antibodies is demonstrated using a panel of antigen-negative cells.
- The potency of the reagents has been tested against the following minimum potency reference standard obtained from National Institute of Biological Standards and Controls (NIBSC):
 - Anti-D reference 99/836.
- The Quality Control of the reagents was performed using red cells with phenotypes that were verified by a UK blood transfusion centre and had been washed with PBS or Isotonic saline prior to use.

DISCLAIMER

- The user is responsible for the performance of the reagents by any method other than those mentioned in the Recommended Techniques.
- Any deviations from the Recommended Techniques should be validated prior to use⁹.

BIBLIOGRAPHY

- Issitt PD. Applied Blood Group Serology, 3rd Edition, Montgomery Scientific, Miami, 1985, Chapter 10.
- AABB Technical Manual, 16th Edition, AABB 2008.

- Marion E. Reid and Christine Lomas-Francis, Blood Group Antigens and Antibodies, SBB Books, New York 2007; Page 192.
- Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine 1995. 5, 171-184
- Guidelines for the Blood Transfusion Service in the United Kingdom, 6th Edition 2002. The Stationery Office.
- British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

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EC REP Advena Ltd, Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013, Malta

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BIBLIOGRAFIE

1. Issitt PD. Applied Blood Group Serology, 3rd Edition, Montgomery Scientific, Miami, 1985, capitolul 10.
2. AABB Technical Manual, 16th edition, AABB 2008.
3. Marion E.Reid & Christine Lomas-Francis, Blood Group Antigens & Antibodies, SBB Books, New York 2007; pagina 192.
4. Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine

- typing of patients and donors. Transfusion Medicine 1991; 1995. **5**, 171-184
5. Guidelines for the Blood Transfusion Service in the United Kingdom, 6th Edition 2002. The Stationary Office.
6. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, **5**, 145-150.

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EC	REP	Advena Ltd. Tower Business Centre, 2 nd Flr., Tower Street, Swatar, BKR 4013, Malta
-----------	------------	---

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SUMMARY

The Rh blood group system was discovered in 1940. The D antigen is the most clinically significant non-ABO red blood cell antigen and has been implicated in causing Haemolytic Transfusion Reactions and Haemolytic Disease of the Newborn.

Anti-D	Phenotype	Caucasians %*	Afro-Americans %*
+	Rh D +ve	83	92
0	Rh D -ve	17	8

INTENDED PURPOSE

The Anti-D reagents are blood grouping reagents intended to be used to qualitatively determine the presence or absence of the Rh D antigen on the red cells of blood donors or patients requiring a blood transfusion when tested in accordance with the recommended techniques stated in this IFU.

PRINCIPLE

The reagents contain antibodies against the D antigen on human red cells and will cause direct agglutination (clumping) of human red cells that carry the D antigen and indirect agglutination of human red cells that are Category D^o in the antiglobulin phase of testing. No agglutination (no clumping) generally indicates the absence of the D antigen on human red cells (see **Limitations**).

REAGENT

Lorne Monoclonal Anti-D Duoclon blood grouping reagent is a low protein, blended reagent containing a human monoclonal IgM and IgG anti-D, diluted in a phosphate buffer containing sodium chloride (0.9 g%), bovine albumin (2.0 g%) and macromolecular potentiators (1.5 g%). When typing patient samples, this reagent will directly agglutinate Rh D positive cells, including majority of variants (but not D^o) and a high proportion of weak D (D^w) phenotypes when using the recommended techniques. The reagents do not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. The reagent is supplied at optimal dilution for use on patient samples with all recommended techniques stated below without need for further dilution or addition. For lot reference number and expiry date see Vial Label.

IgM / IgG	Cell Line / Clone
IgM	RUM-1
IgG	MS-26

WEAKENED EXPRESSION OF THE RhD ANTIGEN

The collective term D^o is widely used to describe red cells which have a weaker expression of the D antigen than normal. The term weak D denotes individuals with a reduced number of complete D antigen sites per red cell. The term partial D denotes individuals with missing D antigen epitopes. D^o is a partial D category which misses most D epitopes. Duoclon reagent will detect most examples of partial and weak D red cells by direct agglutination, but will not detect D^o cells. This reagent will detect D^o and partial D cells in the IAT phase.

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Blood samples can be collected into EDTA, citrate, CPDA anticoagulant or as a clotted sample. The samples should be tested as soon as possible following collection. If a delay in testing should occur, store the samples at 2-8°C. Samples displaying gross haemolysis or microbial contamination should not be used for testing. Blood samples showing evidence of lysis may give unreliable results. It is preferable (but not essential) to wash all blood samples with PBS or isotonic saline before being tested.

PRECAUTIONS

1. The reagent is intended for in vitro diagnostic use only.
2. If a reagent vial is cracked or leaking, discard the contents immediately.
3. Do not use the reagent past the expiration date (see Vial Label).
4. Do not use the reagent if a precipitate is present.
5. Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
6. The reagent has been filtered through a 0.2 µm capsule to reduce the bio-burden, but is not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date as long as there is no marked turbidity, which can indicate reagent deterioration or contamination.
7. The reagent contains <0.1% sodium azide. Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive metal azides. On disposal flush away with large volumes of water.
8. Materials used to produce the reagent were tested at source and found to be negative for HIV 1+2 and HCV antibodies and HBsAg using approved microbiological tests.

9. No known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

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CONTROLS AND ADVICE

1. It is recommended that a positive control (ideally R1r cells) and a negative control (ideally rr cells) be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
2. When typing red cells from a patient who is diagnosed with a disease that causes the red cells to become coated with antibody or other proteins (such as HDN, AIHA), it is important to test the patient's red cells using Lorne's reagent negative control (Monoclonal D Negative Control, catalogue # 650010). Tests must be considered invalid if red cells are agglutinated using Lorne's Monoclonal D Negative Control (catalogue # 650010).
3. Test samples for category D^o determination by the Indirect Antiglobulin Test, Coombs Bio-Rad-ID and Coombs Ortho BioVue Techniques only.
4. Weak and variant D antigens are poorly detected by gel card, microtitre plate and slide techniques. It is recommended that weak and partial variants are tested using the tube test technique.
5. The antiglobulin tube technique can only be considered valid if all negative tests react positively with IgG sensitised red cells.
6. Before use, let the reagent warm up to room temperature. As soon as the reagent has been used, put the reagent back in storage at 2-8°C.
7. In the Recommended Techniques one volume is approximately 50µl when using the vial dropper provided.
8. The use of the reagent and the interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of the country where the reagents are in use.
9. The user must determine suitability of reagents for use in other techniques.

REAGENTS AND MATERIALS REQUIRED

- Anti-human globulin e.g. Lorne AHG Elite (Cat # 435010) or Anti-Human IgG e.g. Lorne Anti-Human IgG (Cat # 402010).
- Applicator sticks.
- Automatic plate reader.
- Coombs cell washer.
- Bio-Rad ID-Cards (LISS/Coombs) and (NaCl, enzyme test and cold agglutinins).
- Bio-Rad ID-Centrifuge.
- Bio-Rad ID-CellStab or ID-Diluent 2.
- Bio-Rad ID-Incubator equilibrated to 37°C ± 2°C.
- Glass microscope slides or white card tiles.
- Glass test tubes (10 x 75 mm or 12 x 75 mm).
- IgG sensitised red cells e.g. Lorne Coombs Control Cells (Cat # 970010).
- Microplate centrifuge.
- Ortho BioVue System Cassettes (AHG/Coombs) and (Neutral).
- Ortho BioVue System Centrifuge.
- Ortho BioVue System Heat Block equilibrated to 37°C ± 2°C.
- Ortho 0.8% Red Cell Diluent.
- Plate shaker.
- PBS solution (pH 6.8-7.2) or Isotonic saline solution (pH 6.5-7.5).
- Positive (ideally R,r) and negative (rr) control red cells.
- Test tube centrifuge.
- Validated "U" well microplates.
- Volumetric pipettes.
- Water bath or dry heat incubator equilibrated to 37°C ± 2°C.

RECOMMENDED TECHNIQUES (NOT CATEGORY D^o)

A. Tube Technique

1. Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
2. Place in a labelled test tube: 1 volume of Lorne Duoclon reagent and 1 volume of red cell suspension.
3. Mix thoroughly and centrifuge all tubes for 20 seconds at 1000 rcf or for a suitable alternative time and force.
4. Gently resuspend red cell button and read macroscopically for agglutination
5. Any tubes, which show a negative or questionable result (which can happen with D^o or weak D samples), should be incubated for 15 minutes at room temperature.
6. Following incubation, repeat steps 3 and 4.

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1. Prepare a 0.8% suspension of red cells in ID-CellStab or ID-Diluent 2.
2. Remove aluminium foil from as many microtubes as needed.
3. Place in appropriate microtube: 50µl test red cell suspension and 25µl Lorne Duoclon reagent.
4. Centrifuge the ID-Card(s) in a Bio-Rad gel card centrifuge.
5. Read macroscopically for agglutination.

C. Ortho BioVue Technique (Neutral cards)

1. Prepare a 0.8% suspension of red cells in 0.8% Ortho Red Cell Diluent
2. Remove aluminium foil from as many reaction chambers as needed.
3. Place in appropriate reaction chamber: 50µl of red cell suspension and 40µl of Lorne Duoclone reagent.
4. Centrifuge cassette(s) in an Ortho BioVue System Centrifuge.
5. Read macroscopically for agglutination.

D. Microplate Technique, using "U" wells

1. Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
2. Place in the appropriate well: 1 volume of Lorne Duoclone reagent and 1 volume of red cell suspension.
3. Mix thoroughly, preferably using a microplate shaker, taking care to avoid cross-well contamination.
4. Incubate at room temperature for 15 minutes (time dependant on user).
5. Centrifuge the microplate for 1 minute at 140 rcf or for a suitable alternative time and force.
6. Resuspend the cell buttons using carefully controlled agitation on a microplate shaker
7. Read macroscopically or with a validated automatic reader.
8. Any weak reactions should be repeated by the tube technique.

E. Slide Technique

1. Prepare a 35-45% suspension of red cells in serum, plasma or PBS or Isotonic saline or use anti-coagulated whole blood (in it's own plasma).
2. Place on a labelled glass slide or card slide: 1 volume of Lorne Duoclone reagent and 1 volume of red cell suspension.
3. Using a clean applicator stick, mix reagent and cells over an area of about 20 x 40 mm.
4. Slowly tilt the slide back and forth for 30 seconds, with occasional further mixing during the 1 minute period, maintaining slide at room temperature.
5. Read macroscopically after 1 minute over a diffuse light and do not mistake fibrin strands as agglutination.
6. Any weak reactions should be repeated by the tube technique.

RECOMMENDED TECHNIQUES (TO DETECT CATEGORY D^W)

A. Indirect Antiglobulin Technique (IAT)

1. Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
2. Place in a labelled test tube: 1 volume of Lorne Duoclone and 1 volume of red cell suspension.
3. Mix thoroughly and incubate at 37°C for 15 minutes.
4. Wash red cells at least once with PBS or Isotonic saline, taking care to decant saline between washes and resuspend each cell button after each wash. Completely decant saline after last wash.
5. Add 2 drops of AHG or anti-IgG to each dry cell button.
6. Mix thoroughly and centrifuge all tubes for 20 seconds at 1000 rcf for a suitable alternative time and force.
7. Resuspend each cell button and read macroscopically.
8. Confirm validity of all negative reactions with IgG sensitised red cells.

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1. Prepare 0.8% suspension of red cells in ID-CellStab or ID-Diluent 2.
2. Remove aluminium foil from as many microtubes as needed.
3. Place in appropriate microtube: 50µl of red cell suspension and 25µl of Lorne Duoclone.
4. Incubate the ID-Card(s) for 15 minutes at 37°C.
5. Centrifuge the ID-Card(s) in a Bio-Rad gel card centrifuge.
6. Read macroscopically for agglutination.

C. Ortho BioVue Technique (AHG/Coombs cards)

1. Prepare a 0.8% suspension of red cells in 0.8% Ortho Red Cell Diluent.
2. Remove aluminium foil from as many reaction chambers as needed.
3. Place in appropriate reaction chamber: 50µl of red cell suspension and 40µl of Lorne Duoclone.
4. Incubate the cassette(s) for 15 minutes at 37°C.
5. Centrifuge cassette(s) in an Ortho BioVue System Centrifuge.
6. Read macroscopically for agglutination.

INTERPRETATION OF TEST RESULTS

1. Positive: Agglutination of the red cells constitutes a positive test result and within accepted limitations of test procedure, indicates the presence of the D antigen on the test red cells.
2. Negative: No agglutination of the red cells constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of the D antigen on the test red cells.
3. Test results of cells that are agglutinated using the reagent negative control shall be excluded, as the agglutination is most probably caused by the effect of the macromolecular potentiators in the reagent on sensitised cells.

STABILITY OF THE REACTIONS

1. Read all tube and microplate tests immediately after centrifugation.
2. Complete washing steps without interruption and centrifuge and read tests immediately after addition of anti-human globulin because delays may result in dissociation of antigen-antibody complexes, leading to false negative or

weak positive reactions.

3. Slide tests should be interpreted after a maximum of 1 minute to ensure specificity and to avoid the possibility a negative result may be incorrectly interpreted as positive due to drying of the reagent.

4. Caution should be exercised in the interpretation of results of tests performed at temperatures other than those recommended.

LIMITATIONS

1. Lorne Anti-D is not suitable for use with enzyme treated cells or cells suspended in LISS.
2. The use of solutions for making red cell suspensions other than those described in the "Recommended Techniques" sections in the document must be validated prior to use. Some solutions may give rise to false positive or false negative reactions.
3. Stored blood may give weaker reactions than fresh blood.
4. False positive agglutination may be seen when testing IgG sensitised cells.
5. False positive or false negative results may also occur due to:
 - Contamination of test materials
 - Improper storage, cell concentration, incubation time or temperature
 - Improper or excessive centrifugation
 - Deviation from the recommended techniques

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1. Prior to release, each lot of Lorne Anti-D Duoclone monoclonal reagent was tested using the recommended test methods listed in this IFU. The tests complied with the test requirements as stated in the current version/issue of the 'Guidelines for the Blood Transfusion Services in the United Kingdom' and the 'Common Technical Specifications'.
2. Specificity of source monoclonal antibodies is demonstrated using a panel of antigen-negative cells.
3. The potency of the reagent has been tested against the following minimum potency reference standard obtained from National Institute of Biological Standards and Controls (NIBSC): Anti-D reference 99/836.
4. The Quality Control of the reagent was performed using red cells with phenotypes that were verified by a UK blood transfusion centre and had been washed with PBS or Isotonic saline prior to use.

DISCLAIMER

1. The user is responsible for the performance of the reagent by any method other than those mentioned in the Recommended Techniques.
2. Any deviations from the Recommended Techniques should be validated prior to use⁶.

BIBLIOGRAPHY

1. Issitt PD. Applied Blood Group Serology, 3rd Edition, Montgomery Scientific, Miami, 1985, Chapter 10.
2. AABB Technical Manual, 16th Edition, AABB 2008.
3. Marion E. Reid and Christine Lomas-Francis, Blood Group Antigens and Antibodies, SBB Books, New York 2007; Page 192.
4. Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine 1995, 5, 171-184
5. Guidelines for the Blood Transfusion Service in the United Kingdom, 6th Edition 2002. The Stationery Office.
6. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

AVAILABLE REAGENT SIZES

Vial Size	Catalogue Number	Tests Per Vial
10 ml	740010	200
1000 ml	740000*	20,000
5000 ml	740000X5*	100,000

*This size is for Further Manufacturing Use (FFMU) only and is therefore not CE marked.

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BIBLIOGRAFIE

- Kholer G, Milstein C. Continuous culture of fused cells secreting antibody of predefined specificity. Nature 1975, 256, 495-497.
- Race RR, Sanger R. Blood Groups in Man 6th Edition, Oxford, Blackwell Scientific Publishers 1975, Chapter 2.
- Issitt PD. Applied Blood Group Serology, 3rd Edition, Montgomery Scientific, Miami, 1985, Chapter 10.
- Mollison PL. Blood Transfusion in Clinical Medicine, 8th Edition, Oxford, Blackwell Scientific Publications, 1987, Chapter 7.
- Tippett P. Sub-divisions of the Rh (D) antigen. Medical. Laboratory Science 1988; 45, 88-93
- Thompson KM, Hughes-Jones NC. Production and characteristics of monoclonal anti-Rh. Bailliere's Clinical Haematology 1990; April
- Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine 1995. 5, 171-184
- Guidelines for the Blood Transfusion Service in the United Kingdom. H.M.S.O. Current Edition.
- British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

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



**BLOOD GROUPING REAGENTS
DIRECTIONS FOR USE**

Anti-C, Anti-E, Anti-c and Anti-e Monoclonal: For Tube, Bio-Rad-ID, Ortho BioVue, Microplate + Slide Techniques.

SUMMARY

Levine and Stetson discovered the Rh blood group system in 1940. Apart from D the other major Rh antigens are C, E, c and e. The D antigen is highly immunogenic; the C and e antigens are less immunogenic than E and c. The corresponding antibodies are all clinically significant since they may cause both Transfusion Reactions and Haemolytic Disease of the Newborn.

INTENDED PURPOSE

The Rh reagents are blood grouping reagents intended to be used to qualitatively determine the presence or absence of Rh antigens on the red cells of blood donors or patients requiring a blood transfusion when tested in accordance with the recommended techniques stated in this IFU.

PRINCIPLE

The reagents contain antibodies to the appropriate Rhesus antigen on human red cells and will cause direct agglutination (clumping) of red cells that carry the corresponding Rh antigen. No agglutination (no clumping) generally indicates the absence of the corresponding Rh antigen (see **Limitations**).

REAGENTS

Lorne Monoclonal IgM Anti-Rh blood grouping reagents are low protein reagents containing human monoclonal antibodies diluted with sodium chloride, bovine albumin and macromolecular potentiators (4.0 g%). The reagents do not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. Each reagent is supplied at optimal dilution for use with all recommended techniques stated below without need for further dilution or addition. For lot reference number and expiry date see **Vial Label**.

Reagent	Cell Line / Clone
Anti-C	MS-24
Anti-E	MS-258
Anti-c	MS-33
Anti-e	MS-16 + MS-63

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PRECAUTIONS

1. The reagents are intended for *in vitro* diagnostic use only.
2. If a reagent vial is cracked or leaking, discard the contents immediately.
3. Do not use the reagents past the expiration date (see **Vial Label**).
4. Do not use the reagents if a precipitate is present.
5. Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
6. The reagents have been filtered through a 0.2 µm capsule to reduce the bio-burden, but is not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date as long as there is no marked turbidity, which can indicate reagent deterioration or contamination.
7. The reagents contain < 0.1% sodium azide. Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive metal azides. On disposal flush away with large volumes of water.
8. Materials used to produce the products were tested at source and found to be negative for HIV 1+2 and HCV antibodies and HBsAg using approved microbiological tests.
9. No known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

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For information on disposal of the reagent and decontamination of a spillage site see **Material Safety Data Sheets**, available on request.

CONTROLS AND ADVICE

1. It is recommended a positive control (ideally heterozygous) and a negative control be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.

2. When typing red cells from patients known or suspected to have auto-antibodies, protein abnormalities or a positive Direct Antiglobulin Test (DAT), it is important that a reagent negative control is tested in parallel. For the reagent negative control, only Lorne Monoclonal Rh Control, catalogue number 640010, must be used.
3. Weak Rh antigens may be poorly detected by the gel card, microtitre plate and slide technique. It is recommended that weak Rh antigens are tested using the tube test technique.
4. Before use, let the reagent warm up to room temperature. As soon as the reagent has been used, put the reagent back in storage at 2-8°C.
5. In the **Recommended Techniques** one volume is approximately 50µl when using the vial dropper provided.
6. The use of reagents and interpretation of results must be carried out by properly trained and qualified personnel in accordance with requirements of the country where the reagents are in use. The user must determine suitability of the reagents for use in other techniques.

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Tube Technique

- Glass test tubes (10 x 75 mm or 12 x 75 mm).
- Centrifuge capable of spinning at 1000 g for 20 seconds.
- PBS solution (pH 6.8-7.2) or Isotonic saline solution (pH 6.5-7.5).
- Positive and negative control red cells:
Monoclonal Anti-C: R₁r (positive control) and rr (negative control).
Monoclonal Anti-E: R₂r (positive control) and rr (negative control).
Monoclonal Anti-c: R₁r (positive control) and R₁R₁ (negative control).
Monoclonal Anti-e: R₂r (positive control) and R₂R₂ (negative control).

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Slide Technique

- Glass microscope slides or white card tiles.
- Applicator sticks.
- Timer or stopwatch

All Techniques

- Volumetric pipettes.

RECOMMENDED TECHNIQUES

A. Tube Technique

1. Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
2. Place in a labelled test tube: 1 volume Lorne Anti-Rh reagent and 1 volume red cell suspension.
3. Mix thoroughly and centrifuge all tubes for 20 seconds at 1000 rcf or for a suitable alternative time and force.
4. Gently resuspend red cell button and read macroscopically for agglutination
5. Any tubes, which show a negative or questionable result, should be incubated for 15 minutes at room temperature.
6. Following incubation, repeat steps 3 and 4.

B. Bio-Rad ID Technique (NaCl, Enzyme and Cold agglutinins cards)

1. Prepare a 0.8% suspension of red cells in ID-CellStab or ID-Diluent 2.
2. Remove aluminium foil from as many microtubes on the NaCl/Enzyme/Cold agglutinins ID-Card as needed.
3. Place in appropriate microtube: 50µl of red cell suspension and 25µl of Lorne Anti-Rh reagent.
4. Centrifuge the ID-Card(s) in a Bio-Rad ID centrifuge.
5. Read macroscopically for agglutination.

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1. Prepare a 0.8% suspension of red cells in 0.8% Ortho Red Cell Diluent.
2. Remove aluminium foil from as many reaction chambers on the Neutral cassette as needed.

- Place in appropriate reaction chamber: 50µl of red cell suspension and 40µl of Lorne Anti-Rh reagent.
- Centrifuge cassette(s) for 5 minutes in an Ortho BioVue System Centrifuge.
- Read macroscopically for agglutination.

D. Microplate Technique, using "U" wells

- Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
- Place in the appropriate well: 1 volume Lorne Anti-Rh reagent and 1 volume red cell suspension.
- Mix thoroughly, preferably using a microplate shaker, taking care to avoid cross-well contamination.
- Incubate at room temperature for 15 minutes (time dependant on user).
- Centrifuge the microplate for 1 minute at 140 rcf or for a suitable alternative time and force.
- Resuspend the cell buttons using carefully controlled agitation on a microplate shaker
- Read macroscopically or with a validated automatic reader.
- Any weak reactions should be repeated by the tube technique.

E. Slide Technique

- Prepare a 35-45% suspension of red cells in serum, plasma or PBS or Isotonic saline. If this is not possible, whole anti-coagulated blood may also be used as the sample.
- Place on a labelled glass slide or card tile: 1 volume of Lorne Anti-Rh reagent and 1 volume of red cell suspension.
- Using a clean applicator stick, mix reagent and cells over an area of about 20 x 40 mm.
- Slowly tilt the slide back and forth for 1-minute, maintaining slide at room temperature.
- Read macroscopically after 1 minute over a diffuse light and do not mistake fibrin strands as agglutination.
- Any weak reactions should be repeated by the tube technique.

INTERPRETATION OF TEST RESULTS

- Positive:** Agglutination of the red cells constitutes a positive test result and within accepted limitations of test procedure, indicates the presence of the appropriate Rh antigen on the red cells.
- Negative:** No agglutination of the red cells constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of the appropriate Rh antigen on the red cells.
- Test results of cells that are agglutinated using the reagent negative control shall be excluded, as the agglutination is most probably caused by the effect of the macromolecular potentiators in the reagent on sensitised cells.

STABILITY OF THE REACTIONS

- Read all tube and microplate tests straight after centrifugation.
- Slide tests should be interpreted within one minute to ensure specificity and to avoid the possibility a negative result may be incorrectly interpreted as positive due to drying of the reagent.
- Caution should be exercised in the interpretation of results of tests performed at temperatures other than those recommended.

LIMITATIONS

- Lorne Anti-Rh reagents are not suitable for use with enzyme treated cells or for use in indirect antiglobulin techniques.
- Many Monoclonal human IgM anti-Rh antibodies have been shown to possess anti-i/I cold agglutinin activity, particularly with cord cells or enzyme treated cells. This may become apparent if tests are incubated below the recommended temperature.
- Some red cells express variant Rh antigens and may give weaker reactions than seen with randomly selected positive control cells. Anti-C may give weaker reactions with C antigen of R₂R₂ individuals. Similarly, Anti-e may give slightly weaker reactions in absence of C antigen, e.g. R₂r, r⁺r and rr.
- Suppressed or diminished expression of certain blood group antigens may conversely give rise to false negative reactions. For these reasons, caution should always be exercised when assigning genotypes on the basis of test results.
- False positive or false negative results may also occur due to:
 - Contamination of test materials
 - Improper storage, cell concentration, incubation time or temperature
 - Improper or excessive centrifugation
 - Deviation from the recommended techniques

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- Prior to release, each lot of Rh reagent was tested using the recommended test methods listed in this IFU. The tests complied with the test requirements as stated in the current version/issue of the "Guidelines for the Blood Transfusion Services in the United Kingdom" and the "Common Technical Specifications".
- Specificity of source Monoclonal antibodies is demonstrated using a panel of antigen-negative cells.
- The Quality Control of the reagents was performed using red cells with phenotypes that were verified by a UK blood transfusion centre and had been washed with PBS or Isotonic saline prior to use.

DISCLAIMER

- The user is responsible for the performance of the reagents by any method other than those mentioned in the **Recommended Techniques**.
- Any deviations from the **Recommended Techniques** should be validated prior to use⁵.

BIBLIOGRAPHY

- Issitt PD. Applied Blood Group Serology, 3rd Edition, Montgomery Scientific, Miami, 1985, Chapter 10.
- AABB Technical Manual, 16th edition, AABB 2008.
- Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine 1995, 5, 171-184
- Guidelines for the Blood Transfusion Service in the United Kingdom, 6th Edition 2002. The Stationary Office.
- British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

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BIBLIOGRAFIE

1. Widman FK. Technical Manual, 9th Edition. American Association of Blood Banks, Arlington, VA, 1985; capitolul 8
2. Race RR, Sanger R. Blood Groups in Man, 6th Edition. Blackwell Scientific, Oxford 1975; capitolul 2
3. Mollison PL. Blood Transfusion in Clinical Medicine, 8th Edition. Blackwell Scientific, Oxford 1987; capitolul 7
4. Issitt PD. Applied Blood Group Serology, 3rd Edition. Montgomery Scientific, Miami 1985; capitolul 6
5. Guidelines for the Blood Transfusion Service in the United Kingdom. H.M.S.O. Ediția curentă.
6. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

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
Advena Ltd. Tower Business Centre, 2nd Flr.,
 Tower Street, Swatar, BKR 4013, Malta

MONOCLONAL BLOOD GROUPING REAGENTS
DIRECTIONS FOR USE

Anti-Fy^a Monoclonal: For Indirect Antiglobulin Techniques.

SUMMARY

The Fy^a and Fy^b antigens were reported in 1950 and 1951 respectively. Anti-Fy^a and anti-Fy^b have both been implicated in immediate and delayed Haemolytic Transfusion Reactions and Haemolytic Disease of the Newborn.

Anti-Fy ^a	Anti-Fy ^b	Phenotype	Caucasians %	Afro-Americans %
+	0	Fy(a+b-)	17	9
0	+	Fy(a-b+)	34	22
+	+	Fy(a+b+)	49	1
0	0	Fy(a-b-)	Very Rare	68

INTENDED PURPOSE

The reagent is a blood grouping reagent intended to be used to qualitatively determine the presence or absence of the Fy(a) antigen (FY001) on the red cells of blood donors or patients requiring a blood transfusion when tested in accordance with the recommended techniques stated in this IFU.

PRINCIPLE

The reagent will cause indirect agglutination (clumping) of red cells, that carry the Fy(a) antigen, in the antiglobulin phase of testing. No agglutination generally indicates the absence of the Fy(a) antigen (see **Limitations**).

REAGENTS

This Monoclonal IgG blood grouping reagent contains human monoclonal antibodies diluted in a phosphate buffer containing sodium chloride and bovine albumin. The reagent does not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. The reagent is supplied at optimal dilution for use with all the recommended techniques stated below without the need for further dilution or addition. For lot reference number and expiry date see **Vial Label**.

Product	Cell Line/Clone
Anti-Fya	DG-FYA-02

STORAGE

Do not freeze. Reagent vials should be stored at 2 - 8°C on receipt. Prolonged storage at temperatures outside this range may result in accelerated loss of reagent reactivity. This reagent has undergone transportation stability studies at 37°C and -25°C as described in document BS EN ISO 23640:2015.

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Blood samples can be collected into EDTA, citrate, CPDA anticoagulants or as a clotted sample. The samples should be tested as soon as possible following collection. If a delay in testing should occur, store the samples at 2-8°C. Samples displaying gross haemolysis or microbial contamination should not be used for testing. Blood samples showing evidence of lysis may give unreliable results. It is preferable (but not essential) to wash all blood samples with PBS or Isotonic saline before being tested.

PRECAUTIONS

- The reagents are intended for *in vitro* diagnostic use only.
- If a reagent vial is cracked or leaking, discard the contents immediately.
- Do not use the reagents past the expiration date (see **Vial Label**).
- Do not use the reagents if a precipitate is present.
- Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
- The reagents have been filtered through a 0.2 µm capsule to reduce the bio-burden, but is not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date as long as there is no marked turbidity, which can indicate reagent deterioration or contamination.
- The reagents contain < 0.1% sodium azide. Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive metal azides. On disposal flush away with large volumes of water.
- Materials used to produce the reagents were tested at source and found to be negative for HIV 1+2 and HCV antibodies and HBsAg using approved microbiological tests.
- No known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

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For information on disposal of the reagents and decontamination of a spillage site see **Material Safety Data Sheets**, available on request.

CONTROLS AND ADVICE

- It is recommended a positive control (ideally heterozygous cells) and a negative control be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
- The antiglobulin techniques can only be considered valid if all negative tests react positively with IgG sensitised red cells.
- Before use, let the reagent warm up to room temperature. As soon as the reagent has been used, put the reagent back in storage at 2-8°C.
- In the **Tube Technique** one volume is approximately 50µl when using the vial dropper provided.
- The use of the reagents and the interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of the country where the reagents are in use.
- User must determine suitability of the reagents for use in other techniques.

REAGENTS AND MATERIALS REQUIRED

- Anti-human globulin i.e. Lorne AHG Elite (Cat # 435010 or 415010) or Anti-Human IgG i.e. Lorne Anti-Human IgG (Cat # 402010 or 401010).
- Coombs cell washer.
- Bio-Rad ID-Cards (LISS/Coombs or Coombs Anti-IgG).
- Bio-Rad ID-Centrifuge.
- Bio-Rad ID-CellStab or ID-Diluent 2.
- Bio-Rad ID-Incubator equilibrated to 37°C ± 2°C.
- Glass test tubes (10 x 75 mm or 12 x 75 mm).
- IgG sensitised red cells i.e. Lorne Coombs Control Cells (Cat # 970010).
- Ortho BioVue System Cassettes (AHG Polyspecific or AHG Anti-IgG).
- Ortho BioVue System Centrifuge.
- Ortho BioVue System Heat Block equilibrated to 37°C ± 2°C.
- Ortho 0.8% Red Cell Diluent.
- PBS solution (pH 6.8-7.2) or Isotonic saline solution (pH 6.5-7.5).
- Positive (ideally heterozygous) and negative control red cells.
- Volumetric pipettes.
- Water bath or dry heat incubator equilibrated to 37°C ± 2°C.

RECOMMENDED TECHNIQUES

A. Indirect Antiglobulin Technique (IAT)

- Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
- Place in a labelled test tube: 1 volume of Lorne reagent and 1 volume of red cell suspension.
- Mix thoroughly and incubate at 37°C for 15 minutes.
- Wash red cells 1 time with PBS or Isotonic saline, taking care to decant saline after the wash.
- Add 2 volumes of anti-human globulin or anti-IgG to each dry cell button.
- Mix thoroughly and centrifuge all tubes for 20 seconds at 1000 rcf or for a suitable alternative time and force.
- Gently resuspend red cell button and read macroscopically for agglutination
- Confirm validity of all negative reactions with IgG sensitised red cells.

B. Bio-Rad ID Technique (LISS/Coombs cards)

- Prepare a 0.8% suspension of red cells in ID-CellStab or ID-Diluent 2.
- Remove aluminium foil from as many microtubes on either LISS/Coombs or Coombs Anti-IgG ID-Card(s) as needed.
- Place in appropriate microtube: 50µl of red cell suspension and 25µl of Lorne reagent.
- Incubate the LISS/Coombs ID-Card(s) for 15 minutes at 37°C.
- Centrifuge the LISS/Coombs ID-Card(s) in the Bio-Rad ID-Card centrifuge.
- Read macroscopically for agglutination.

C. Ortho BioVue Technique (AHG cassettes)

- Prepare a 0.8% suspension of red cells in 0.8% Ortho Red Cell Diluent.
- Remove aluminium foil from as many reaction chambers on either AHG Polyspecific or AHG Anti-IgG cassettes as needed.
- Place in appropriate reaction chamber: 50µl of red cell suspension and 40µl of Lorne reagent.
- Incubate the cassette(s) for 15 minutes at 37°C.
- Centrifuge cassette(s) in an Ortho BioVue System Centrifuge.
- Read macroscopically for agglutination.

INTERPRETATION OF TEST RESULTS

- Positive:** Agglutination of the red cells constitutes a positive test result and within accepted limitations of test procedure, indicates the presence of the appropriate antigen on the red cells.
- Negative:** No agglutination of the red cells constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of the appropriate antigen on the red cells.

STABILITY OF THE REACTIONS

1. Washing steps should be completed without interruption and tests centrifuged and read immediately after addition of the reagent. Delays may result in dissociation of antigen-antibody complexes, causing false negative or weak positive results.
2. Caution should be exercised in the interpretation of results of tests performed at temperatures other than those **recommended**.

LIMITATIONS

1. Red cells that have a positive DAT due to a coating of IgG cannot be typed by the **Indirect Antiglobulin Technique**.
2. Suppressed or diminished expression of certain blood group antigens may conversely give rise to false negative reactions and so caution should always be exercised when assigning genotypes on the basis of test results.
3. False positive or false negative results may also occur due to:
 - Contamination of test materials
 - Improper storage, cell concentration, incubation time or temperature
 - Improper or excessive centrifugation
 - Deviation from the recommended techniques

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1. Prior to release, each lot of reagent was tested using the recommended test methods listed in this IFU. The tests complied with the test requirements as stated in the current version/issue of the "Guidelines for the Blood Transfusion Services in the United Kingdom".
2. Specificity of source monoclonal antibodies is demonstrated using a panel of antigen-negative cells
3. The Quality Control of the reagents was performed using red cells with phenotypes that were verified by a UK blood transfusion centre and had been washed with PBS or Isotonic saline prior to use.

DISCLAIMER

1. The user is responsible for the performance of the reagents by any method other than those mentioned in the **Recommended Techniques**.
2. Any deviations from the **Recommended Techniques** should be validated prior to use⁶.

BIBLIOGRAPHY

1. Widman FK. Technical Manual, 9th Edition. American Association of Blood Banks, Arlington, VA, 1985; Chapter 8
2. Race RR, Sanger R. Blood Groups in Man, 6th Edition. Blackwell Scientific, Oxford 1975; Chapter 2
3. Mollison PL. Blood Transfusion in Clinical Medicine, 8th Edition. Blackwell Scientific, Oxford 1987; Chapter 7
4. Issitt PD. Applied Blood Group Serology, 3rd Edition. Montgomery Scientific, Miami 1985; Chapter 6
5. Guidelines for the Blood Transfusion Service in the United Kingdom. H.M.S.O. Current Edition.
6. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

AVAILABLE REAGENT SIZES

	Vial Size	Catalogue Number	Tests per vial
Anti-Fya	2 ml	774002	40
Monoclonal	1000 ml	774000*	20,000

*This size is For Further Manufacturing Use (FFMU) only and is therefore not CE marked.

Lorne Laboratories Limited
 Unit 1 Cutbush Park Industrial Estate
 Danehill
 Lower Earley
 Berkshire, RG6 4UT
 United Kingdom
 Tel: +44 (0) 118 921 2264
 Fax: +44 (0) 118 986 4518
 E-mail: info@lornelabs.com

EC	REP	Advena Ltd. Tower Business Centre, 2 nd Fl., Tower Street, Swatar, BKR 4013, Malta
----	-----	--

MONOCLONAL BLOOD GROUPING REAGENTS
DIRECTIONS FOR USE

Anti-S and Anti-s Monoclonal: For Indirect Antiglobulin Techniques.

SUMMARY

The S and s antigens were reported in 1947 and 1951 respectively and form part of the MNS system. Anti-S and anti-s have both been implicated in Haemolytic Transfusion Reactions and Haemolytic Disease of the Newborn.

Anti-S	Anti-s	Phenotype	Caucasians ¹	Afro-Americans ¹
+	0	S+s-	11%	6%
+	+	S+s+	44%	25%
0	+	S-s+	45%	68%
0	0	S-s-	0%	1.5%

INTENDED USE

The reagents are blood grouping reagents intended to be used to qualitatively determine the presence or absence of the S or s antigen on the red cells of blood donors or patients requiring a blood transfusion when tested in accordance with the recommended techniques stated in this IFU.

PRINCIPLE

The reagents contain antibodies to the S or s antigen on human red cells and will cause indirect agglutination (clumping) of human red cells, that carry the corresponding specific antigen, in the antiglobulin phase of testing. No agglutination (clumping) generally indicates the absence of the corresponding specific antigen (see **Limitations**).

REAGENTS

These Monoclonal IgG blood grouping reagents contain human monoclonal antibodies diluted in a phosphate buffer containing sodium chloride and bovine albumin. The reagents do not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. Each reagent is supplied at optimal dilution for use with all the recommended techniques stated below without the need for further dilution or addition. For lot reference number and expiry date see **Vial Label**.

Product	Cell Line/Clone
Anti-S	P3S13JS123
Anti-s	P3YAN3

STORAGE

Do not freeze. Reagent vials should be stored at 2 - 8°C on receipt. Prolonged storage at temperatures outside this range may result in accelerated loss of reagent reactivity. This reagent has undergone transportation stability studies at 37°C and -25°C as described in document BS EN ISO 23640:2015.

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Blood samples can be collected into EDTA, citrate, CPDA anticoagulants or as a clotted sample. The samples should be tested as soon as possible following collection. If a delay in testing should occur, store the samples at 2-8°C. Samples displaying gross haemolysis or microbial contamination should not be used for testing. Blood samples showing evidence of lysis may give unreliable results. It is preferable (but not essential) to wash all blood samples with PBS or Isotonic saline before being tested.

PRECAUTIONS

- The reagents are intended for *in vitro* diagnostic use only.
- If a reagent vial is cracked or leaking, discard the contents immediately.
- Do not use the reagents past the expiration date (see **Vial Label**).
- Do not use the reagents if a precipitate is present.
- Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
- The reagents have been filtered through a 0.2 µm capsule to reduce the bio-burden, but is not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date as long as there is no marked turbidity, which can indicate reagent deterioration or contamination.
- The reagents contain 0.1% sodium azide. Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive metal azides. On disposal flush away with large volumes of water.
- Materials used to produce the reagents were tested at source and found to be negative for HIV 1+2 and HCV antibodies and HBsAg using approved microbiological tests.
- No known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

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For information on disposal of the reagents and decontamination of a spillage site see **Material Safety Data Sheets**, available on request.

CONTROLS AND ADVICE

- It is recommended a positive control (ideally heterozygous cells) and a negative control be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
- The antiglobulin techniques can only be considered valid if all negative tests react positively with IgG sensitised red cells.
- Before use, let the reagent warm up to room temperature. As soon as the reagent has been used, put the reagent back in storage at 2-8°C.
- In the **Tube Technique** one volume is approximately 50µl when using the vial dropper provided.
- The use of the reagents and the interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of the country where the reagents are in use.
- The user must determine the suitability of the reagents for use in other techniques.

REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

Tube Technique

- Anti-human globulin i.e. Lorne AHG Elite (Cat # 435010 or 415010) or Anti-Human IgG i.e. Lorne Anti-Human IgG (Cat # 402010 or 401010).
- Coombs cell washer.
- Glass test tubes (10 x 75 mm or 12 x 75 mm).
- PBS solution (pH 6.8-7.2) or Isotonic saline solution (pH 6.5-7.5).
- IgG sensitised red cells i.e. Lorne Coombs Control Cells (Cat # 970010).
- Positive (ideally heterozygous) and negative control red cells.
- Water bath or dry heat incubator equilibrated to 37°C ± 2°C.

Bio-Rad-ID Micro Typing Technique

- Bio-Rad ID-Cards (LISS/Coombs or Coombs Anti-IgG).
- Bio-Rad ID-Centrifuge.
- Bio-Rad ID-CellStab or ID-Diluent 2.
- Bio-Rad ID-Incubator equilibrated to 37°C ± 2°C.

Ortho BioVue Typing Technique

- Ortho BioVue System Cassettes (AHG Polyspecific or AHG Anti-IgG).
- Ortho BioVue System Centrifuge.
- Ortho BioVue System Heat Block equilibrated to 37°C ± 2°C.
- Ortho 0.8% Red Cell Diluent.

All Techniques

- Volumetric pipettes.

RECOMMENDED TECHNIQUES

A. Indirect Antiglobulin Technique (IAT)

- Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
- Place in a labelled test tube: 1 volume of Lorne reagent and 1 volume of red cell suspension.
- Mix thoroughly and incubate at 37°C for 15 minutes.
- Wash red cells 1 time with PBS or Isotonic saline, taking care to completely decant saline after the wash.
- Add 2 volumes of anti-human globulin or anti-IgG to each dry cell button.
- Mix thoroughly and centrifuge all tubes for 20 seconds at 1000 rcf or for a suitable alternative time and force.
- Gently resuspend red cell button and read macroscopically for agglutination
- Confirm validity of all negative reactions with IgG sensitised red cells.

B. Bio-Rad-ID Technique (LISS/Coombs cards)

- Prepare a 0.8% suspension of red cells in the ID-Diluent or ID-Diluent 2.
- Remove aluminium foil from as many microtubes on either LISS/Coombs or Coombs Anti-IgG ID-Card(s) as needed.
- Place in appropriate microtube: 50µl of red cell suspension and 25µl of Lorne reagent.
- Incubate the LISS/Coombs ID-Card(s) for 15 minutes at 37°C.
- Centrifuge the LISS/Coombs ID-Card(s) in the Bio-Rad ID-Card centrifuge.
- Read macroscopically for agglutination.

C. Ortho BioVue Technique (AHG cassettes)

- Prepare a 0.8% suspension of red cells in 0.8% Ortho Red Cell Diluent.
- Remove aluminium foil from as many reaction chambers on either AHG Polyspecific or AHG Anti-IgG cassettes as needed.
- Place in appropriate reaction chamber: 50µl of red cell suspension and 40µl of Lorne reagent.

4. Incubate the cassette(s) for 15 minutes at 37°C.
5. Centrifuge cassette(s) in an Ortho BioVue System Centrifuge.
6. Read macroscopically for agglutination.

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INTERPRETATION OF TEST RESULTS

1. **Positive:** Agglutination of the red cells constitutes a positive test result and within accepted limitations of test procedure, indicates the presence of the appropriate antigen on the red cells.
2. **Negative:** No agglutination of the red cells constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of the appropriate antigen on the red cells.

STABILITY OF THE REACTIONS

1. Washing steps should be completed without interruption and tests centrifuged and read immediately after addition of the reagent. Delays may result in dissociation of antigen-antibody complexes, causing false negative or weak positive results.
2. Caution should be exercised in the interpretation of results of tests performed at temperatures other than those **recommended**.

LIMITATIONS

1. Red cells that have a positive DAT due to a coating of IgG cannot be typed by the **Indirect Antiglobulin Technique**.
2. Suppressed or diminished expression of certain blood group antigens may conversely give rise to false negative reactions and so caution should always be exercised when assigning genotypes on the basis of test results.
3. False positive or false negative results may also occur due to:
 - Contamination of test materials
 - Improper storage, cell concentration, incubation time or temperature
 - Improper or excessive centrifugation
 - Deviation from the recommended techniques

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1. Prior to release, each lot of reagent was tested using the recommended test methods listed in this IFU. The tests complied with the test requirements as stated in the current version/issue of the "Guidelines for the Blood Transfusion Services in the United Kingdom".
2. Specificity of source monoclonal antibodies is demonstrated using a panel of antigen-negative cells
3. The Quality Control of the reagents was performed using red cells with phenotypes that were verified by a UK blood transfusion centre and had been washed with PBS or Isotonic saline prior to use.

DISCLAIMER

1. The user is responsible for the performance of the reagents by any method other than those mentioned in the **Recommended Techniques**.
2. Any deviations from the **Recommended Techniques** should be validated prior to use⁵.

BIBLIOGRAPHY

1. Marion E.Reid & Christine Lomas-Francis, Blood Group Antigens & Antibodies, SBB Books, New York 2007; Page 190.
2. Issitt PD. Applied Blood Group Serology, 3rd Edition. Montgomery Scientific, Miami 1985; Chapter 14.
3. AABB Technical Manual, 16th edition, AABB 2008.
4. Guidelines for the Blood Transfusion Service in the United Kingdom, 6th Edition 2002. The Stationary Office.
5. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, **5**, 145-150.

AVAILABLE REAGENT SIZES

	Vial Size	Catalogue Number	Tests per vial
Anti-S Monoclonal	2 ml	770002	40
	1000 ml	770000*	20,000
Anti-s Monoclonal	2 ml	771002	40
	1000 ml	771000*	20,000

*This size is For Further Manufacturing Use (FFMU) only is therefore not CE marked.

Lorne Laboratories Limited
Unit 1 Cutbush Park Industrial Estate
Danehill
Lower Earley
Berkshire, RG6 4UT
United Kingdom
Tel: +44 (0) 118 921 2264
Fax: +44 (0) 118 986 4518
E-mail: info@lornelabs.com

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BIBLIOGRAFIE

- Voak D, Downie DM, Moore BPL, and Engelfreit CP. Anti-Human Globulin reagent specification. The European and ISBT/ICSH View. Biotest Bulletin 1: 7-22 (1986).
- The Department of Health and Social Security. Health Services Management Antiglobulin Test. False negative results, HN (Hazard) (83) 625 Nov 1983.
- Bruce M, Watt AH, Hare W, Blue A, Mitchell R. A serious source of error in antiglobulin testing. Transfusion 1986; **26**: 177-181.
- Voak D, Downie DM, Moore BPL, Ford DS, Engelfreit CP, Case J. Replicate tests for the detection and correction of errors in AHG (AHG) tests: optimum conditions and quality control. Haematologia 1988; **21**(1): 3-16.
- Guidelines for the Blood Transfusion Service in the United Kingdom, 6th Edition 2002. The Stationary Office.
- British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, **5**, 145-150.

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
Advena Ltd. Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013, Malta

DIRECTIONS FOR USE

AHG Elite (Clear or Green): For Antiglobulin Techniques.**SUMMARY**

In 1945, Coombs, Mourant and Race described the use of anti-human globulin serum for detecting red cell-bound non-agglutinating antibodies. In 1957, Dacie et al showed that the antibodies present in antiglobulin sera were directed against certain components of complement. Anti-human globulin reagents detect non-agglutinating antibody molecules as well as molecules of complement attached to red cells following *in vivo* or *in vitro* antigen-antibody reactions.

INTENDED PURPOSE

These reagents are polyspecific blood grouping reagents intended to be used to qualitatively detect the presence or absence of sensitising IgG antibodies (all 4 subclasses) and complement factors C3d and C3b on human red cells when tested in accordance with the recommended techniques stated in this IFU.

PRINCIPLE

The reagents contain antibodies against human IgG antibodies and C3 complement factors (C3d and C3b) on human red cells and will cause direct agglutination (clumping) of red cells that are sensitised with human IgG antibodies and/or C3 complement factors (C3d and C3b). No agglutination generally indicates the absence of sensitising human IgG antibodies and C3 complement factors (C3d and C3b) on human red cells (See **Limitations**).

REAGENT

Lorne AHG Elite Clear and AHG Elite Green reagents contain anti-IgG derived from rabbits with non-specific activity removed by adsorption and mouse monoclonal IgM anti-C3d, Clone BRIC-8. The antibodies are diluted in a buffered solution containing bovine albumin. The reagents do not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. Each reagent is supplied at optimal dilution, for use with all the recommended techniques stated below without need for further dilution or addition. For lot reference number and expiry date see Vial Label.

Reagent	Cell Line/Clone	Colour	Dye Used
AHG Elite Clear	Rabbit Anti-Human IgG BRIC-8 (Anti-C3d)	Colourless	None
AHG Elite Green	Rabbit Anti-Human IgG BRIC-8 (Anti-C3d)	Green	Patent Blue and Tartrazine

STORAGE

Reagent vials should be stored at 2 - 8°C on receipt. Prolonged storage at temperatures outside this range may result in accelerated loss of reagent reactivity. This reagent has undergone transportation stability studies at 37°C and -25°C as described in document BS EN ISO 23640:2015.

SAMPLE COLLECTION AND PREPARATION

Samples should be drawn aseptically into EDTA to prevent *in vitro* complement binding and tested as soon as possible. If EDTA is unavailable, samples drawn into ACD, CPD or CPDA-1 are preferable to clotted ones. If only clotted samples are available, do not refrigerate them before testing.

PRECAUTIONS

- The reagents are intended for *in vitro* diagnostic use only.
- If a reagent vial is cracked or leaking, discard the contents immediately.
- Do not use the reagents past the expiration date (see Vial Label).
- Do not use the reagents if a precipitate is present.
- Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
- The reagents have been filtered through a 0.2 µm capsule to reduce the bio-burden, but are not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date as long as there is no marked turbidity, which can indicate reagent deterioration or contamination.
- The reagents contain < 0.1% sodium azide. Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive metal azides. On disposal flush away with large volumes of water.
- Materials used to produce the products were tested at source and found to be negative for HIV 1+2 and HCV antibodies and HBsAg using approved microbiological tests.
- No known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

DISPOSAL OF REAGENT AND DEALING WITH SPILLAGES

For information on disposal of the reagents and decontamination of a spillage site see Material Safety Data Sheets, available on request.

CONTROLS AND ADVICE

- It is recommended a positive control (weak Anti-D <0.1 IU/ml) and a negative control (an inert serum) be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
- The antiglobulin techniques can only be considered valid if all negative tests react positively with IgG sensitised red cells.

- Before use, let the reagent warm up to room temperature. As soon as the reagent has been used, put the reagent back in storage at 2-8°C.
- In the Recommended Techniques one volume is approximately 50µl when using the vial dropper provided.
- Use of the reagents and the interpretation of results must be carried out by properly trained and qualified personnel in accordance with requirements of the country where the reagents are in use. User must determine the suitability of the reagents for use in other techniques.

REAGENTS AND MATERIALS REQUIRED

- Coombs cell washer.
- Glass test tubes (10 x 75 mm or 12 x 75 mm).
- IgG sensitised red cells e.g. Lorne Coombs Control Cells (Cat # 970010).
- Inert antibody e.g. Lorne Inert AB Serum (Cat # 110010).
- Low Ionic Strength Solution (LISS): Containing 0.03M NaCl, 0.003M Na₂HPO₄, NaH₂PO₄ buffer pH 6.7 at 22°C ± 1°C and 0.24M glycine.
- PBS solution (pH 6.8-7.2) or Isotonic saline solution (pH 6.5-7.5).
- Volumetric pipettes.
- Water bath or dry heat incubator equilibrated to 37°C ± 2°C.
- Weak anti-D e.g. Lorne Precise Weak Anti-D (Cat # 209005).

RECOMMENDED TECHNIQUES**A. Direct Antiglobulin Technique (DAT)**

- Wash 1 volume of red cells (2-3% suspension in PBS or Isotonic saline) 4 times with PBS or Isotonic saline, taking care to decant saline between washes and resuspend each cell button after each wash. Completely decant saline after last wash.
- Add 2 volumes of Lorne AHG Elite to each dry cell button.
- Mix thoroughly and centrifuge all tubes for 20 seconds at 1000 rcf or for a suitable alternative time and force.
- Gently resuspend red cell button and read macroscopically for agglutination

B. Indirect Antiglobulin Technique (NISS IAT)

- Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
- Place in a labelled test tube: 2 volumes of test serum and 1 volume of red cell suspension.
- Mix thoroughly and incubate at 37°C for 15 minutes.
- Wash red cells 4 times with PBS or Isotonic saline, taking care to decant saline between washes and resuspend each red cell button after each wash. Completely decant saline after last wash.
- Add 2 volumes of Lorne AHG Elite to each dry cell button.
- Mix thoroughly and centrifuge all tubes for 20 seconds at 1000 rcf or for a suitable alternative time and force.
- Gently resuspend red cell button and read macroscopically for agglutination

C. LISS Indirect Antiglobulin Technique (LISS IAT)

- Prepare a 1.5-2% suspension of red cells in LISS.
- Place in a labelled test tube: 2 volumes of test serum and 2 volumes of red cell suspension.
- Mix thoroughly and incubate at 37°C for 15 minutes.
- Follow steps 4 to 7 of NISS IAT above.

INTERPRETATION OF TEST RESULTS

- Positive: Agglutination of test red cells constitutes a positive test result and within the accepted limitations of the test procedure, indicates the presence of IgG and/or complement (C3d/C3b) on the red cells.
- Negative: No agglutination of the test red cells constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of IgG and complement (C3d/C3b) on the red cells.

STABILITY OF THE REACTIONS

- Washing steps should be completed without interruption and tests centrifuged and read immediately after addition of the reagent. Delays may result in dissociation of antigen-antibody complexes, causing false negative or weak positive results.
- Caution should be exercised in the interpretation of results of tests performed at temperatures other than those recommended.

LIMITATIONS

- Red cells that have a positive DAT due to a coating of IgG cannot be typed by the Indirect Antiglobulin Techniques.
- A positive DAT due to complement sensitisation may not reflect *in vivo* complement fixation if test cells are from a refrigerated clotted specimen.
- Inadequate washing of red cells in the indirect antiglobulin techniques may neutralise the AHG reagent.
- Following completion of the wash phase excess residual saline may dilute the AHG Elite, reducing its potency.
- A negative direct antiglobulin test result does not necessarily preclude clinical diagnosis of ABO Haemolytic Disease of the Newborn or Auto Immune Haemolytic Anaemia. It also does not necessarily rule out HDN, especially if ABO incompatibility is suspected.
- False positive or false negative results may also occur due to:

- Contamination of test materials
- Improper storage, cell concentration, incubation time or temperature
- Improper or excessive centrifugation

SPECIFIC PERFORMANCE CHARACTERISTICS

1. Prior to release, each lot of the reagents were tested using the recommended test methods listed in this IFU against red cells coated with Anti-D, Anti-K and Anti-FyA to check suitable reactivity. The tests complied with the test requirements as stated in the current version/issue of the 'Guidelines for the Blood Transfusion Services in the United Kingdom'.
2. The anti-IgG and anti-C3d potencies have been tested against the following minimum potency reference standard obtained from the National Institute of Biological Standards and Controls (NIBSC): Anti-AHG reference standard 96/666
3. Anti-C3d potency is demonstrated in tests employing cells coated with C3d and C3b.
4. The presence of contaminating heterospecific agglutinins or antibodies to C4d has been excluded in tests employing red cells of all ABO groups and cells coated with C4d.
5. The reactivity of any Anti-IgM, Anti-IgA or Anti-light chain components that might be present has not been established.
6. The Quality Control of the reagents was performed using red cells with phenotypes that were verified by a UK blood transfusion centre and had been washed with PBS or Isotonic saline prior to use.

DISCLAIMER

1. The user is responsible for the performance of the reagents by any method other than those mentioned in the Recommended Techniques.
2. Any deviations from the Recommended Techniques should be validated prior to use¹².

BIBLIOGRAPHY

1. Voak D, Downie DM, Moore BPL, and Engelfreit CP. Anti-Human Globulin reagent specification. The European and ISBT/ICSH View. Biotest Bulletin 1; 7-22 (1986).
2. The Department of Health and Social Security. Health Services Management Antiglobulin Test. False negative results, HN (Hazard) (83) 625 Nov 1983.
3. Bruce M, Watt AH, Hare W, Blue A, Mitchell R. A serious source of error in antiglobulin testing. Transfusion 1986; 26: 177-181.
4. Voak D, Downie DM, Moore BPL, Ford DS, Engelfreit CP. Case J. Replicate tests for the detection and correction of errors in anti-human globulin (AHG) tests: optimum conditions and quality control. Haematologia 1988; 21(1): 3-16.
5. Guidelines for the Blood Transfusion Service in the United Kingdom. 6th Edition 2002. The Stationery Office.
6. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

AVAILABLE REAGENT SIZES

	Vial Size	Catalogue Number	Test Per Vial
Lorne AHG Elite	10 ml	415010	100
(Clear)	1000 ml	415000*	10,000
Lorne AHG Elite	10 ml	435010	100
(Clear)	1000 ml	435000*	10,000

*These sizes are For Further Manufacturing Use (FFMU) only and are therefore not CE marked.

Advena Ltd, Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

LORNE
LABORATORIES

Lorne Laboratories Limited

Unit 1 Cutbush Park Industrial Estate, Danehill, Lower Earley, Berkshire RG6 4UT United Kingdom

Tel: +44 (0) 118 921 2264 Fax: +44 (0) 118 986 4518 Email: info@lornelabs.com www.lornelabs.com

Blood Reagents and Diagnostic Kits

Quality blood reagents and diagnostic kits delivered worldwide

Lorne Laboratories provides high quality blood-grouping reagents and diagnostic kits. We do so at affordable prices and back this up with excellent customer service.

For over forty four years, Lorne has seen that philosophy transform the company from a small UK operation to one that spans the globe. The Lorne name is known and respected in over 110 countries around the world and it has come to mean 'quality' to the many blood transfusion professionals who use the products, both in the UK and across the five continents in which we operate.

Lorne has achieved ISO 13485 and MDSAP accreditation and our reagents and diagnostic kits are all CE marked. A significant number of our blood-grouping reagents are also registered with Health Canada.

Contents

BLOOD GROUPING REAGENTS Page

ABO system – Monoclonal Reagents.....	4
Rhesus System – Monoclonal Reagents.....	4
Grouping Reagents for M, N and S Blood Group Systems	5
Grouping Reagents for Kell Blood Group System.....	5
Grouping Reagents for Rare Blood Groups.....	6
Control Reagents.....	7
Anti-Human IgG Reagent.....	7
Anti-Human Globulin Reagent	7
Monoclonal Anti-C3d Reagent.....	7

ENZYMES AND POTENTIATORS Page

Buffered Saline Tablets	8
Papenzyme-plus.....	8
Bromelite	8
LISS Concentrate.....	8
LISS Ready for use.....	8
LISS-ADD	8
PEG-ADD.....	8
Serological Albumin	8
Accessories	9

RED CELLS Page

Red Cells – Reverse Grouping Cells	10
Red Cells – Antibody Screening Cells.....	10
Red Cells – Identicells.....	10
Red Cells – Coombs Control Cells.....	10
Red Cells Preservative – Preservacell.....	11
Red Cells Preservative – ABO Preservacell.....	11
Alsevers Solution.....	11
Red Cell Elute.....	11

DIAGNOSTIC KITS Page

Syphilis Kits	12
Latex Kits.....	13
Febrile Antigens	14
Rose Waaler.....	14

NON CE MARKED BLOOD GROUPING REAGENTS Page

ABO System - Monoclonal Reagents.....	15
Rhesus System - Monoclonal Reagents.....	15

ABO SYSTEM – MONOCLONAL REAGENTS

Lorne Monoclonal IgM ABO blood grouping reagents contain mouse monoclonal antibodies diluted in a phosphate buffer containing sodium chloride, EDTA and bovine albumin. Each reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

Item	Product Code	Size	Maximum Shelf Life
Anti-A Monoclonal	600010	10ml	36 Months
Anti-B Monoclonal	610010	10ml	36 Months
Anti-A,B Monoclonal	620010	10ml	36 Months

Lorne Anti-A₁ Lectin blood grouping reagent is prepared from an extract of *Dolichos biflorus* seeds, diluted with a sodium chloride solution containing bovine albumin. The reagent is supplied at optimal dilution for use by tube technique.
Lorne Anti-H Lectin blood grouping reagent is prepared from an extract of *Ulex europaeus* seeds, diluted with a sodium chloride solution containing bovine albumin. The reagent is supplied at optimal dilution for use by tube technique.

Anti-A ₁ Lectin	116005	5ml	24 Months
Anti-H Lectin	115002	2ml	24 Months

RHESUS SYSTEM – MONOCLONAL REAGENTS

Lorne Monoclonal IgM Anti-D Clone 1 and Clone 2 blood grouping reagents are low protein reagents containing a human monoclonal IgM antibody diluted with sodium chloride, bovine albumin and macromolecular potentiators. When typing patient samples, each reagent will directly agglutinate Rh D positive cells, including the majority of variants (but not D^v) and a high proportion of weak D (D^u) phenotypes when using the slide, tube, gel card and microplate techniques.

Item	Product Code	Size	Maximum Shelf Life
Anti-D Clone 1 Monoclonal	730010	10ml	30 Months
Anti-D Clone 2 Monoclonal	710010	10ml	30 Months

Lorne Monoclonal Anti-D Duoclone blood grouping reagent is a low protein, blended reagent containing human monoclonal IgM and IgG Anti-D, diluted in a phosphate buffer containing sodium chloride, bovine albumin and macromolecular potentiators. When typing patient samples, this reagent will directly agglutinate Rh D positive cells, including the majority of variants (but not D^v) and a high proportion of weak D (D^u) phenotypes when using the slide, tube, gel card and microplate techniques. It will agglutinate D^v cells in the IAT phase of testing.

Anti-D Duoclone Monoclonal	740010	10ml	30 Months
----------------------------	--------	------	-----------

Lorne Monoclonal IgM Anti-Rh blood grouping reagents are low protein reagents containing human monoclonal antibodies diluted with sodium chloride, bovine albumin and macromolecular potentiators. Each reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

Anti-C Monoclonal	690005	5ml	24 Months
Anti-C ^w Monoclonal	750002	2ml	24 Months
Anti-E Monoclonal	691005	5ml	24 Months
Anti-c Monoclonal	692005	5ml	24 Months
Anti-e Monoclonal	693005	5ml	24 Months
Anti-C+D+E Monoclonal	700010	10ml	24 Months

All pack inserts are available on www.lornelabs.com

GROUPING REAGENTS FOR M, N AND S BLOOD GROUP SYSTEMS

Lorne Human Anti-M blood grouping reagent is prepared from human serum diluted in a sodium chloride solution containing bovine albumin. The reagent is supplied at optimal dilution for use by tube and gel card techniques.

Item	Product Code	Size	Maximum Shelf Life
Anti-M Polyclonal	311002	2ml	24 Months

Lorne Anti-N Lectin blood grouping reagent is prepared from an extract of *Vicia*, diluted with a sodium chloride solution containing bovine albumin. The reagent is supplied at optimal dilution for use by tube and gel card techniques.

Anti-N Lectin	312002	2ml	24 Months
---------------	--------	-----	-----------

Lorne Monoclonal Anti-S and Anti-s blood grouping reagents contain human monoclonal antibodies diluted in a phosphate buffer containing sodium chloride and bovine albumin. They are supplied at optimal dilution for use by indirect Coombs tube test and gel card techniques.

Anti-S Monoclonal	770002	2ml	24 Months
Anti-s Monoclonal	771002	2ml	24 Months

GROUPING REAGENTS FOR KELL BLOOD GROUP SYSTEM

Lorne Monoclonal Anti-K blood grouping reagent is a low protein reagent containing the monoclonal IgM antibody, diluted in a phosphate buffer containing sodium chloride, bovine albumin and macromolecular potentiators. The reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

Item	Product Code	Size	Maximum Shelf Life
Anti-K Monoclonal	760010	10ml	24 Months

Lorne Monoclonal Anti-k (Cellano) blood grouping reagent contains monoclonal IgG antibody diluted in sodium chloride containing macromolecular potentiators and bovine albumin. Each reagent is supplied at optimal dilution for use by indirect Coombs tube and gel card techniques without the need for further dilution or addition.

Anti-k (Cellano) Monoclonal	325002	2ml	24 Months
-----------------------------	--------	-----	-----------

Lorne Human Anti-Kp^a and Anti-Kp^b blood grouping reagents are prepared from human serum diluted in sodium chloride containing macromolecular potentiators and bovine albumin. Each reagent is supplied at optimal dilution for use by indirect Coombs tube and gel card techniques without the need for further dilution or addition.

Anti-Kp ^a Polyclonal	321002	2ml	24 Months
Anti-Kp ^b Polyclonal	322002	2ml	24 Months

All products available in bulk quantities

GROUPING REAGENTS FOR RARE BLOOD GROUPS

Lorne Monoclonal Fy^a blood grouping reagent contains human monoclonal antibodies diluted in a phosphate buffer containing sodium chloride and bovine albumin.

Lorne Human Anti-Fy^b blood grouping reagent is prepared from human serum diluted in a sodium chloride solution containing macromolecular potentiators and bovine albumin. Each reagent is supplied at optimal dilution for use by indirect Coombs tube test and gel card techniques.

Item	Product Code	Size	Maximum Shelf Life
Anti-Fy ^a Monoclonal	774002	2ml	24 Months
Anti-Fy ^b Polyclonal	317002	2ml	24 Months

Lorne Anti-Jk^a and Anti-Jk^b Monoclonal reagents contain human monoclonal IgM antibodies, sodium chloride, phosphate buffer and bovine albumin. When typing red cells in accordance with the instructions for use, the reagents will directly agglutinate either Jk^a positive cells or Jk^b positive cells. Both reagents are supplied at optimal dilution for use by the test tube technique.

Anti-Jk ^a Monoclonal	775002	2ml	24 Months
Anti-Jk ^b Monoclonal	776002	2ml	24 Months

Lorne Human Anti-Lu^a and Anti-Lu^b blood grouping reagents are prepared from human serum diluted in a sodium chloride solution containing macromolecular potentiators and bovine albumin. Each reagent is supplied at optimal dilution for use by indirect Coombs tube test and gel card techniques.

Anti-Lu ^a Polyclonal	330002	2ml	24 Months
Anti-Lu ^b Polyclonal	331002	2ml	24 Months

Lorne Monoclonal IgM Anti-P₁ blood grouping reagent contains mouse monoclonal IgM antibodies, diluted in a solution containing sodium chloride and bovine albumin. The reagent is supplied at optimal dilution for use by tube and gel card techniques.

Anti-P ₁ Monoclonal	315002	2ml	24 Months
--------------------------------	--------	-----	-----------

Lorne Monoclonal Anti-Le^a contains human monoclonal IgM and Anti-Le^b contains mouse monoclonal IgM antibodies, diluted in a phosphate buffer containing sodium chloride, EDTA, bovine albumin and macromolecular potentiators. Each reagent is supplied at optimal dilution for use by the tube technique. The Anti-Le^a reagent can also be used for gel card techniques.

Anti-Le ^a Monoclonal	632002	2ml	24 Months
Anti-Le ^b Monoclonal	631002	2ml	24 Months

Lorne Human Anti-Di^a blood grouping reagent is prepared from human serum diluted in a sodium chloride solution containing macromolecular potentiators and bovine albumin. The reagent is supplied at optimal dilution for use by indirect Coombs tube test and gel card techniques.

Anti-Di ^a Polyclonal	328002	2ml	24 Months
---------------------------------	--------	-----	-----------

All pack inserts are available on www.lornelabs.com

CONTROL REAGENTS

Lorne Monoclonal Rh Control is formulated with the same levels of potentiators and protein as Lorne Monoclonal Rh Grouping Reagents with the blood group antibodies omitted. The reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

Item	Product Code	Size	Maximum Shelf Life
Monoclonal Rh Control	640010	10ml	24 Months

Lorne Monoclonal D Negative Control is for the control of Monoclonal Anti-D reagents and is formulated with the same concentrations of phosphate buffer, sodium chloride, bovine albumin and macromolecular potentiators as Lorne Monoclonal Anti-D reagents with the antibodies omitted. The reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

Monoclonal D Negative Control	650010	10ml	30 Months
-------------------------------	--------	------	-----------

Lorne Inert AB Serum is prepared from pooled human serum. No potentiators or any other chemicals have been added to the reagent. This reagent is supplied at optimal dilution for use by the techniques recommended for the reagent to be controlled.

Inert AB Serum	110010	10ml	24 Months
----------------	--------	------	-----------

Lorne Precise Weak Anti-D Control Reagent is prepared from pools of human serum containing low activity Anti-D. The pool is diluted in inert serum to give a final concentration of 0.09 IU/ml Anti-D. ABO antibodies are not absorbed. This polyclonal reagent is supplied at the optimal dilution, for use by the techniques recommended for the reagent to be controlled.

Precise Weak Anti-D	209005	5ml	24 Months
---------------------	--------	-----	-----------

ANTI-HUMAN IgG REAGENT

Lorne Monospecific Anti-Human IgG Clear and Anti-Human IgG Green reagents contain anti-human IgG derived from rabbits. All non-specific activity is removed by absorption. These monospecific reagents are supplied at optimal dilution, for use by spin tube technique.

Item	Product Code	Size	Maximum Shelf Life
Anti-Human IgG (Clear)	401010	10ml	24 Months
Anti-Human IgG (Green)	402010	10ml	24 Months

ANTI-HUMAN GLOBULIN REAGENT

Lorne Polyspecific Anti-Human Globulin Elite Clear and Anti-Human Globulin Elite Green reagents contain anti-human IgG derived from rabbits with non-specific activity removed by absorption and mouse monoclonal IgM Anti-C3d, Clone BRIC-8. The antibodies are diluted in a buffered solution containing bovine albumin. These polyspecific reagents are supplied at optimal dilution, for use by spin tube technique.

Item	Product Code	Size	Maximum Shelf Life
A.H.G Elite (Clear)	415010	10ml	24 Months
A.H.G Elite (Green)	435010	10ml	24 Months

MONOCLONAL ANTI-C3d REAGENT

Lorne Monoclonal IgM Anti-C3d blood grouping reagent contains mouse monoclonal Anti-C3d, Clone BRIC-8. The reagent is supplied at optimal dilution, for use by direct tube technique.

Item	Product Code	Size	Maximum Shelf Life
Anti-C3d Monoclonal	427002	2ml	24 Months

All products available in bulk quantities

ENZYMES AND POTENTIATORS

Lorne Phosphate Buffered Saline Tablets provide a safe, standardised saline solution for transfusion serology. Each tablet makes 1 litre of solution.

Item	Product Code	Size	Maximum Shelf Life
Buffered Saline Tablets	490025	25	24 Months
	490250	250	24 Months

Lorne Papanzyme-plus reagent is a ready to use liquid preparation of stabilised papain. The reagent is standardised by serological methods for use in blood group antibody investigations. The reagent is supplied at optimal dilution for use by tube technique.

Papanzyme-plus	441010	10ml	12 Months
----------------	--------	------	-----------

Lorne Bromelite reagent is a ready to use liquid preparation of stabilised bromelin. The reagent is standardised by serological methods for use in blood group antibody investigations. The reagent is supplied at optimal dilution for use by tube technique.

Bromelite	443010	10ml	12 Months
-----------	--------	------	-----------

Lorne LISS Concentrate is a solution of glycine, phosphate buffer and 0.3 M sodium chloride. The solution is supplied at a stronger concentration than needed for serological use. It must be diluted 10 times in de-ionised water before being used by all recommended techniques in the pack insert.

LISS Concentrate	460500	500ml	24 Months
	460025	2500ml	24 Months

Lorne LISS ready for use is a low ionic strength solution containing glycine, sodium chloride and phosphate buffer. The reagent is supplied at the optimal dilution ready for use by all recommended techniques in the pack insert.

LISS Ready for use	470020	20 x 250ml	12 Months
	470250	4 x 250ml	12 Months
	470025	2500ml	12 Months

Lorne LISS-ADD is a low ionic strength solution containing glycine, sodium chloride, phosphate buffer and bovine albumin. The reagent is supplied at the optimal dilution, for use by all the recommended techniques in the pack insert.

LISS-ADD	480010	10ml	24 Months
----------	--------	------	-----------

Lorne PEG-ADD is a low ionic strength solution contains glycine, a phosphate buffer and polyethylene glycol. The reagent is supplied at optimal dilution for use by all the recommended techniques stated in the pack insert without the need for further dilution or addition.

PEG-ADD	485010	10ml	24 Months
---------	--------	------	-----------

Lorne 22% and 30% Serological Albumin is prepared from a mixture of bovine serum albumin and buffered saline. No artificial avidity enhancers or high molecular weight agglutination potentiators are added to any BSA preparation. None of the BSA reagents contain sodium caprylate. Each BSA reagent is supplied at optimal dilution for use by all recommended techniques in the pack insert.

Serological Albumin 22%	451010	10ml	24 Months
Serological Albumin 30%	452010	10ml	24 Months

All products available in bulk quantities

ACCESSORIES

Item	Product Code	Size
 <p>The Essex Blockfile is a sturdy box for the storage of up to 200 cassette type histology blocks in each box. Delivered flat in packs of 25, they are easily assembled when needed and provide convenient and economical long-term storage for tissue samples. Dimensions – (L) 390mm x (D) 200mm x (H) 45mm</p>	882025	25 Boxes
 <p>The Brighton Slidestak is a sturdy box for the storage of over 1000 microscope slides upright in each box. Delivered flat in packs of 25, they are easily assembled and provide convenient and economical long-term storage. Dimensions – (L) 310mm x (D) 170mm x (H) 80mm</p>	881025	25 Boxes
 <p>The Cardtiles are suitable for a wide range of haemagglutination tests. Each card has twenty 30mm square white test sites.</p>	880100	100 Cards
 <p>The Mini-Cardtiles are specially laminated cards suitable for VDRL and similar agglutination tests requiring a white background. Each card has ten 20mm diameter white test sites that are slightly compressed to form shallow wells.</p>	880120	25 Cards
 <p>The Latex Cardtiles are suitable for all latex agglutination tests. Each card has six 30mm diameter black reaction sites.</p>	880130	25 Cards
 <p>The mini pipettes are 11.5cm in length with one sealed end. The diameter of the tube is 4mm.</p>	044000	500
 <p>We use these products in our own production. The 5ml vials are glass and the 10ml vials are plastic. Each vial is supplied with its own dropper. The 10ml vials are supplied in a sterile format in multiples of 1224 while all the other products are supplied individually and require autoclaving.</p>	LAB00002	5ml
Vials and Droppers	LAB00036	10ml
 <p>We use these products in our own production. They are not branded. We supply flatpacks that hold 5 or 10 vials</p>	LAB00027	5 Vial
Flatpacks	LAB00007	10 Vial

All pack inserts are available on www.lornelabs.com

RED CELLS - REVERSE GROUPING CELLS

These Grouping Cells are made from red cells that have been washed to remove all traces of blood group antibodies and then resuspended in a preservative solution of buffered saline containing adenine, inosine, chloramphenicol and neomycin sulphate. The preservative solution does not interfere with complement-mediated haemolysis. Each group A₁, A₂ and B vial contains a 2.5-3.0% red cell suspension derived from the blood of a single donor, whereas the group O vial contains a 2.5-3.0% suspension of pooled group O red cells from two donors in equal proportions. Antigens for which the O cells have been typed are noted on the antigen profile accompanying each lot. NB: Some antigens are present on only 50% of the cells in each lot.

Item	Product Code	Size	Minimum Shelf Life
Reverse Grouping A ₁ + B Cells	910020	2 x 10ml	30 Days
Reverse Grouping A ₁ + A ₂ + B Cells	910030	3 x 10ml	30 Days
Reverse Grouping A ₁ + A ₂ + B + O Cells	910040	4 X 10ml	30 Days
Reverse Grouping A ₂ Cells	920002	1 x 10ml	30 Days

RED CELLS - ANTIBODY SCREENING CELLS

Maxi-Screen 3 red cells are for antibody screening. Each reagent vial contains a 2.5-3.0% suspension of red cells derived from the blood of a single group O donor. The donor red cells have been washed to remove blood group antibodies and then resuspended in a preservative solution containing adenine and inosine to help preserve carbohydrate metabolism and chloramphenicol and neomycin sulphate as preservatives.

Item	Product Code	Size	Minimum Shelf Life
Maxi-Screen 3	950030	3 x 10ml	30 Days

RED CELLS - IDENTICELLS

Identicells are made up of 10 vials which each contain a 2.5-3.0% suspension of red cells derived from the blood of a single group O donor. The donor cells have been washed to remove blood group antibodies and then resuspended in a preservative solution containing adenine, inosine, chloramphenicol and neomycin sulphate as preservatives.

Item	Product Code	Size	Minimum Shelf Life
Identicells	960050	10 x 5ml	30 Days

RED CELLS - COOMBS CONTROL CELLS

Coombs Control Cells are made up of a 3.8-4.2% suspension of single donor group O red cells washed to remove all blood group antibodies and then resuspended in a preservative solution. The preservative solution contains neomycin sulphate and chloramphenicol as preservatives. The cells are then sensitised with IgG.

Item	Product Code	Size	Minimum Shelf Life
Coombs Control Cells	970010	10ml	30 Days

Deliveries take place every 28 days.

As well as the standard Red Cell products above, we can also supply other products for special standing orders only. These include 0.8% panels and screening cells for use with column agglutination systems and Papanised red cell panels. Contact us with full details of your requirements for a quotation.

All the above Red Cell products are CE marked and comply with the Red Book.

All pack inserts are available on www.lornelabs.com

They are not Lorne labelled products.

RED CELLS PRESERVATIVE – PRESERVACELL

Lorne Preservacell is a phosphate buffered solution containing glucose, calcium chloride and purine bases, with chloramphenicol, gentamycin sulphate and neomycin sulphate as antibiotics. The reagent is supplied at optimal dilution.

Item	Product Code	Size	Maximum Shelf Life
Preservacell	980500	500ml	12 Months

RED CELLS PRESERVATIVE – ABO PRESERVACELL

Lorne ABO Preservacell is a phosphate buffered solution containing glucose, calcium chloride and purine bases, with chloramphenicol, gentamycin sulphate and neomycin sulphate as antibiotics. EDTA is added to stop complement binding so that potent ABO haemolysins act as simple agglutinins. The reagent is supplied at optimal dilution.

Item	Product Code	Size	Maximum Shelf Life
ABO Preservacell	981500	500ml	12 Months

ALSEVERS SOLUTION

Alsevers Solution is an isotonic, balanced salt solution that is routinely used as an anticoagulant/blood preservative, which permits the storage of whole blood at refrigerator temperatures for approximately 10 weeks. The solution contains antibiotics and other chemicals that ensure maximum viability of stored red cells.

Item	Product Code	Size	Maximum Shelf Life
Alsevers Solution	983000	1000ml	12 Months

RED CELL ELUTE

Lorne Red Cell Elute is an acid elution kit. The kit consists of Concentrated Wash Solution, which is used to minimise antibody dissociation during washing, Acid Eluting Solution, which is a low pH glycine buffer containing a colouring agent and a Base Buffering Solution, Tris solution containing bovine albumin. The Concentrated Wash Solution requires dilution and all the other solutions are supplied at optimal dilution for use by all recommended techniques in the pack insert.

Item	Product Code	Size	Maximum Shelf Life
Red Cell Elute	930110	Kit	30 Months

All products available in bulk quantities

SYPHILIS KITS

Lorne TPHA Kit detects antibodies to *T. pallidum*. Test Cells are preserved avian erythrocytes coated with antigenic components of pathogenic *T. pallidum* (Nichol's strain). Any non-specific reactions are detected using the Control Cells; avian erythrocytes not coated with *T. pallidum* antigens. Non-specific reactions can also be absorbed out using Control Cells. Antibodies to non-pathogenic treponemes are absorbed by an extract of Reiter's treponemes in the cell suspension. Reagents are supplied at optimal dilution for use by the recommended techniques in the pack insert.

Item	Product Code	Size	Maximum Shelf Life
TPHA Microtitre Plate Kit	043100A	100 Tests	18 Months

The VDRL test is a non-treponemal slide agglutination test for the qualitative and semi-quantitative detection of plasma reagins. The antigen suspension, a lipid complex, is agglutinated when mixed with samples containing reagins of patients affected by syphilis. The reagent is supplied ready to use.

VDRL Stabilised Reagent Kit	046511A	250 Tests	30 Months
-----------------------------	---------	-----------	-----------

RPR Carbon Antigen contains micro particulate carbon, which aids in the microscopic reading of results. Lorne provides reagents, controls and kits. All the reagents are supplied at optimum dilution for use by all recommended techniques without the need for further dilution or addition.

RPR Carbon Antigen	045005A	100 Tests	30 Months
RPR Carbon Positive Control	047001A	1ml	30 Months
RPR Carbon Kit	044150A	150 Tests	30 Months
	044500A	500 Tests	30 Months

RPR Carbon Kit

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Strep Test Kit

All pack inserts are available on www.lornelabs.com

LATEX KITS

Lorne ASO Latex Kit is a serologic test for the detection of ASO antibodies. All the reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

Item	Product Code	Size	Maximum Shelf Life
ASO Latex Kit	031100A	100 Tests	30 Months

Lorne RF Latex Kit is for the detection of Rheumatoid Factor. The latex reagent is a suspension of latex particles coated with human gamma globulins, which agglutinate in the presence of Rheumatoid Factor. All latex reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

RF Latex Kit	830100A	100 Tests	30 Months
--------------	---------	-----------	-----------

Lorne CRP Latex Test Kit is for the detection of CRP. The test reagent consists of latex particles coated with either goat or rabbit Anti-CRP (IgG). All the latex reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

CRP Latex Kit	850100A	100 Tests	30 Months
---------------	---------	-----------	-----------

Lorne IM Latex Test Kit is for detection of the heterophile antibody associated with Infectious Mononucleosis. The test reagent consists of latex particles coated with partially purified glycoprotein from bovine red cells. All reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

IM Latex Kit	041050A	50 Tests	30 Months
--------------	---------	----------	-----------

Lorne LE Latex Kit is for the detection of nuclear proteins in Systemic Lupus Erythematosus (SLE or LE). The test reagent consists of DNP coated latex particles. All the reagents are supplied at optimal dilution for use by recommended techniques in the pack insert.

LE Latex Test Kit	840050	50 Tests	24 Months
-------------------	--------	----------	-----------

Lorne Strep Kit is a Latex agglutination grouping kit for the identification of Streptococci of Lancefield groups A, B, C, D, F and G by agglutination of specific antibody coated latex particles in the presence of enzymatically-extracted antigen. All the reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

Strep Test Kit	860050	6 x 50 Tests	18 Months
----------------	--------	--------------	-----------

Lorne Staph Kit is a Latex agglutination kit for the identification of *Staph. aureus*. Includes latex reagent, control reagent and agglutination slides. All the reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

Staph Test Kit	870050	50 Tests	18 Months
	870100	100 Tests	18 Months

All products available in bulk quantities

FEBRILE ANTIGENS

Lorne Stained Febrile Antigens are for the detection of certain Salmonellae, Rickettsiae and Brucellae pathogens. The antigens are suspensions of killed bacteria, stained to enhance the reading of agglutination tests. The blue stained antigens are specific to the somatic “O” antigens and the red stained antigens are specific to the flagellar “H” antigens. Suspensions of Proteus OX2, OX19 and OXK are used to detect rickettsial antibodies.

Rose Bengal is a slide agglutination test for the qualitative and semi-quantitative detection of anti-Brucella antibodies in human and animal serum. The stained bacterial suspension agglutinates when mixed with samples containing specific IgG or IgM antibodies present in the patient sample.

Item	Product Code	Size	Maximum Shelf Life
Salmonella Typhi H	502005A	100 Tests	30 Months
Salmonella Paratyphi AH	504005A	100 Tests	30 Months
Salmonella Paratyphi BH	506005A	100 Tests	30 Months
Salmonella Paratyphi CH	508005A	100 Tests	30 Months
Salmonella Typhi O	510005A	100 Tests	30 Months
Salmonella Paratyphi AO	512005A	100 Tests	30 Months
Salmonella Paratyphi BO	514005A	100 Tests	30 Months
Salmonella Paratyphi CO	516005A	100 Tests	30 Months
Brucella Abortus	518005A	100 Tests	30 Months
Brucella Melitensis	520005A	100 Tests	30 Months
Proteus OX2	522005A	100 Tests	30 Months
Proteus OX19	524005A	100 Tests	30 Months
Proteus OXK	526005A	100 Tests	30 Months
Febrile Antigen Kit + Controls	532042A	8x100 Tests and 2x1ml	30 Months
Febrile Positive Control	536001A	1ml	30 Months
Febrile Negative Control	537001A	1ml	30 Months
Rose Bengal	155050A	50 Tests	30 Months

ROSE WAALER

Rose Waaler is a technique of passive haemagglutination for qualitative and semi-quantitative detection of human serum Rheumatoid Factor. Sheep red blood cells are coated with a concentration of antiserum to sheep red blood cells that is too low to cause agglutination. The addition of serum from a patient with Rheumatoid Factor will cause agglutination.

Item	Product Code	Size	Maximum Shelf Life
Rose Waaler	156050A	50 Tests	30 Months

Febrile Antigen Kit + Controls

All pack inserts are available on www.lornelabs.com

ABO SYSTEM – MONOCLONAL REAGENTS

Lorne Monoclonal IgM ABO Standard Grade blood grouping reagents contain mouse monoclonal antibodies diluted in a phosphate buffer containing sodium chloride, EDTA and bovine albumin. Each reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

Item	Product Code	Size	Maximum Shelf Life
Anti-A Monoclonal - Standard Grade	600010E	10ml	36 Months
Anti-B Monoclonal - Standard Grade	610010E	10ml	36 Months
Anti-A,B Monoclonal - Standard Grade	620010E	10ml	36 Months

RHESUS SYSTEM – MONOCLONAL REAGENTS

Lorne Monoclonal IgM Anti-D Clone 1 Standard Grade blood grouping reagent is a low protein reagent containing a human monoclonal IgM antibody diluted with sodium chloride, bovine albumin and macromolecular potentiators. When typing patient samples, the reagent will directly agglutinate Rh D positive cells, including the majority of variants (but not D^{vi}) and a high proportion of weak D (D^u) phenotypes when using the slide, tube, gel card and microplate techniques.

Item	Product Code	Size	Maximum Shelf Life
Anti-D IgM Monoclonal - Standard Grade	730010E	10ml	30 Months

Lorne Monoclonal Anti-D Duoclone Standard Grade blood grouping reagent is a low protein, blended reagent containing human monoclonal IgM and IgG Anti-D, diluted in a phosphate buffer containing sodium chloride, bovine albumin and macromolecular potentiators. When typing patient samples, this reagent will directly agglutinate Rh D positive cells, including the majority of variants (but not D^{vi}) and a high proportion of weak D (D^u) phenotypes when using the slide, tube, gel card and microplate techniques. It will agglutinate D^{vi} cells in the IAT phase of testing.

Anti-D Duoclone Monoclonal - Standard Grade	740010E	10ml	30 Months
---	---------	------	-----------

All products available in bulk quantities

Lorne Laboratories Limited
Unit 1 Cutbush Park Industrial Estate
Danehill, Lower Earley
Berkshire RG6 4UT United Kingdom

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meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the “final inspection report” of the device. /

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Vacuum Blood Collection Tube

EDTA Tube

Vacuum EDTA Tubes (Lavender or Purple Top) contain EDTA as an anticoagulant and are widely used for determinations in EDTA whole blood clinical hematology tests, especially for a routine blood examination.

Vacuum Blood Collection Tube, EDTA Tube

Material: PET, Glass

Color: Lavender or Purple

Cat. No.	Specification	Volume	Additive	Qty/case (Glass)	Qty/case (PET)
630801	13×75mm	1ml	EDTAK2	1800	1800
630802	13×75mm	2ml	EDTAK2	1800	1800
630803	13×75mm	3ml	EDTAK2	1800	1800
630804	13×75mm	4ml	EDTAK2	1800	1800
630805	13×75mm	5ml	EDTAK2	1800	1800
630806	13×100mm	5ml	EDTAK2	1200	1200
630807	13×100mm	6ml	EDTAK2	1200	1200
630808	13×100mm	7ml	EDTAK2	1200	1200
630809	16×100mm	10ml	EDTAK2	900	900
630901	13×75mm	1ml	EDTAK3	1800	1800
630902	13×75mm	2ml	EDTAK3	1800	1800
630903	13×75mm	3ml	EDTAK3	1800	1800
630904	13×75mm	4ml	EDTAK3	1800	1800
630905	13×75mm	5ml	EDTAK3	1800	1800
630906	13×100mm	5ml	EDTAK3	1200	1200
630907	13×100mm	6ml	EDTAK3	1200	1200
630908	13×100mm	7ml	EDTAK3	1200	1200
630909	16×100mm	10ml	EDTAK3	900	900
631001	13×75mm	2ml	EDTANA2	1800	1800
631002	13×100mm	10ml	EDTANA2	1200	1800