INgezim® FVR Compac

Prod Ref: 13.FVR.K3

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KIT COMPOSITION

<i>REAGENT</i>	<i>2 Plates</i>		<i>5 Plates</i>		<i>10 Plates</i>	
	Uni.	Vol.	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5		10	
Vials with Positive Control	1	4mL	2	4mL	4	4mL
Vials with Negative Control	1	4mL	2	4mL	4	4mL
Vials ready to use (MAb anti N protein- HRPO)	1	30mL	2	30mL	4	30mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL	1	250mL
Bottles with Diluent (DE01-01), ready to use	1	125mL	1	125mL	1	250mL
Bottles with Substrate (ABTS), ready to use	1	30mL	1	60mL	2	60mL
Bottles with Stop Solution, ready to use	1	60mL	1	60mL	2	60mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 ml.
ELISA Reader (405 nm filter)

TECHNICAL BASIS

This kit is based on a blocking enzymatic immunoassay (**Blocking Elisa**). The technique is briefly described below:

The antigen is fixed in a solid support (polystyrene plate). After incubating the sera samples, a RVFV specific monoclonal antibody (MAb peroxidase conjugated) is added. If the sample contains antibodies specific of the virus, they will not allow the binding of the labelled MAb to the antigen whereas if it does not contain specific antibodies, the MAb will bind to the antigen which is coating the plate.

After washing the plate to eliminate all non-fixed material, the presence or absence of labelled MAbs can be detected by adding the substrate which, in presence of the peroxidase, will develop a colorimetric reaction.

If there has been colour development, this means that the conjugate has bound to the antigen, being therefore the sample negative. On the other hand, if the sample contains specific antibodies to RVFV, they will block the binding of the conjugate and there will no colour development will be detected.

The antigen used in this kit is the N recombinant protein from the RVFV. This confers many advantages to the assay. Being produced recombinantly, it therefore totally lacks of infectivity. Plus the N protein is one of the mayor proteins from RVFV, as well as the most antigenic.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used and never return the excess substrate to the bottle.
8. Stop solution can precipitate at refrigeration, in this case, incubate at 37°C before use, to eliminate the precipitation. The stop solution may form foams. Avoid dispensing bubbles on the plate as this may affect the reading of the results.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.

Control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them in aliquots at -20°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry.
- After the last step of washing shake the plate turned over an absorbent filter paper.

PREPARATION OF SAMPLES

Sera samples must be assayed at 1/5 dilution in the diluent provided. This dilution can be done directly in the plate by adding 80 µL of the diluent and then adding 20 µL of sample to each well.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrated washing solution provided in the kit with 24 parts of distilled or deionized water (i.e., 40 mL of concentrate and 960 mL of water). Once prepared, this solution remains stable when stored at +4°C.

Preparation of Controls (+) y (-):

Both controls are ready to use and do not require any further preparation. Dispense 100 µL of each one.

Preparation of conjugate:

This reagent is ready to use. **Do not dilute.**

TEST PROCEDURE

1. Prior to starting the test, bring all reagents to room temperature (20°C to 25°C).
2. **Addition of sera control**
 - Add 80 µL of the diluent and then 20 µL of each sample in every well to be assayed. Add 100 µL of positive and negative controls to their respective wells. We recommend running controls in duplicate.
 - Seal the plate and incubate **45 minutes at room temperature.**
3. Wash 4 times following the procedure previously described.
4. Add 100 µL of the ready to use conjugate to each well. **Seal the plate and incubate for 30 minutes at room temperature.**
5. Wash 4 times following the procedure previously described.
6. Add 100 µL of the substrate solution to each well. Keep the plate **at room temperature for 10 min.**
7. Add 100 µL of the stop solution to each well.
8. Read the OD of each well at 405 nm.

READING AND RESULT INTERPRETATION

When running samples in duplicate, the OD value for each sample should be calculated as the arithmetic mean of both values. In the same way, the arithmetic mean of the values obtained in the two wells of the positive control and the two wells of the negative control must be calculated.

VALIDATION OF THE RESULTS:

The test is considered valid when:

- The OD of the Negative Control (NC) > 0.80
- The OD of the Positive Control (PC) < 0.25

INTERPRETATION OF THE RESULTS

Calculate the inhibition percentage (IP) of each sample as follows:

$$\text{IP of Sample} = 100 - \left(\frac{\text{Sample OD}}{\text{Negative control OD}} \right) \times 100$$

- Samples will be considered **POSITIVE** (i.e. the sample contains specific antibodies to RVFV), when the IP is $\geq 45\%$.
- Samples will be considered **NEGATIVE** (i.e. the sample contains specific antibodies to RVFV) when the IP is $\leq 40\%$
- Samples with an IP between 40 and 45% should be considered as **DOUBTFUL** and be retested by a seroneutralization (SN) assay.

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Designed & manufactured in Spain by:

GOLD STANDARD DIAGNOSTICS MADRID SA
C/Hermanos García Noblejas, 41 2ª planta

28037- MADRID (SPAIN)

Tlf: +34 91368.05.01/04

Fax: +34 91 408.75.98

E-mail: info.spain@eu.goldstandarddiagnostics.com

www.goldstandarddiagnostics.com

9191.INGE
9175.INGE

Distribuido en
Distributed in

por
by