

### SUPRAPHOB REGEN

**Supraphob Regen** is a new generation diffractive-refractive multifocal intra-ocular lens designed to offer EXCELLENT VISION TO THE PATIENTS undergoing CATARACT SURGERY providing spectacle independence.






- Supraphob Regen is made of premium quality biocompatibile hydrophobic acrylic material
  - Fabricated by nano resolution high technology machine with minimum polishing
    - The non-touch preloaded system makes IOL implantation simple
    - Regen is a diffractive-refractive combined multifocal technology
      - The central 4.2 mm diameter has 20 diffractive zones and remaining 1.8 mm (upto 6 mm) has refractive zones
        - The step height gradually declines towards the periphery
          - This distributes the light between the distant, intermediate and near focal points, resulting in a crystal clear visual quality at all distances

# **CLEAR VISION AT ALL DISTANCES**

- Good functional vision at all distances
- Excellent patient satisfaction
- Enjoy the crystal clear vision at all lighting conditions




## **ADD POWER**



- Regen has +3.5D add power. It separates the near and distant images and makes them distinctive.
- With the years of experience in multifocal IOL design, the addition of +3.5D yields the patient, benefit of both near and distant vision at its best level.

# **QUALITY & ACCURACY**

- All lenses are tested for resolution with air force target immediately after lathing
- Each Regen is INDIVIDUALLY TESTED with Optical Bench System called IOLA-MFD for power, quality and MTF for near and far focal points
- Thanks to the advanced technology and material each Regen are assured of higher MTF values
- High MTF values ensure higher contrast at all distances in all light conditions



### NO TO PCO

- Each "Regen" is having 360° square edge posteriorly
- Thanks to the superior quality of the hydrophobic material the posterior capsule securely adheres to the rear surface of the lens
- The anterior placement of diffractive optics protects optical properties from YAG laser
- Together with the material and square edge optics Regen offers excellent barrier to PCO



### SPECIAL FEATURES

- The overall Haptic and Optic design offers excellent centration of IOL irrespective of the capsular bag size
- The negative Asphericity provided in Regen offers better visual quality than any other spherical IOL



SPHERICAL IOL



ASPHERIC OPTICS

- The Regen IOL has ultraviolet filtering chromosphores for the better protection of retina from UV light
- The state of art technology of Regen makes sure that its optical performance is equal to that of young natural crystalline lens
- The unique anterior diffractive design reduces the halos and yields superior vision




OTHER IOL




SUPRA PHOB REGEN IOL

# PRELOADED SYSTEM

• The advanced preloaded system makes the IOL implanation simpler & safer.



Step-1: Please take the preloaded injector as shown



Step-2: Apply the viscoelastic as shown



Step-3: Close the cartridge, the IOL protecting cover will come out automatically



Step-4: Press the cartridge well untill a "CLICK" sound is heard and ensure the cartridge is locked perfectly



Step-5: Apply the viscoelastic as shown in the picture



Step-6: Ready to inject



Overall Length 13.00 mm

Optic Diameter 6.00 mm

Haptic Angulation 0°

Optic Design Equi-biconvex, with diffractive anterior surface.

Diopter Range ! 7.00D TO 30.00D (0.5D increment)

Material Yellow Hydrophobic Acrylic

Refractive Index 1.497

Delivery System i Ready to inject preloaded

Special Features Aspheric Optics with 360° square edge

Manufactured by:

**∕®D** APPASAMY OCULAR DEVICES (P) LTD

R.S. NO. 9/1, 2&3, N.H. 45-A, Villupuram Main Road, Vadamanagalam, Pondicherry - 605 102, INDIA