Cobalt Chromium Coronary Stent System

Yukon[®] CC

*Translumina -
Crossing all lesions*

heart for the heart

translumina[®]

Yukon[®] CC

Micro-porous surface that matters

The Translumina **Yukon[®] CC** represents a new generation stent system technology based on a Cobalt Chromium alloy.

It is the first CoCr Stent with an unique micro-porous surface.

The Yukon CC is in continuation of the Translumina Surface philosophy, improving healing and biocompatibility of the implant.

For detailed information see also literature [1,2,3].


Unique stent surface

The Translumina micro-porous stent surface, called **PEARL Surface**, favours better endothelialisation, which is essential in avoiding thrombosis and restenosis.

Strut thickness

With a low strut thickness of approx. 68µm the stent has excellent healing properties.

SEM (scanning electron microscope)

Surface features

- approx. 1 million pores per cm²
- average micro pore depth of approx. 2 µm
- 100% pore coverage of the surface

References

- ^[1] Dibra et al, Catheterization and Cardiovascular Interventions 2005, 65, 374-380; Results to 6 months, n=100 BMS with rough surface vs. n=100 with polished surface
- ^[2] Palmaz et al, Journal of Vasc. Surg., 2002, 36, 2001-09.
- ^[3] Fitzgerald et al, Expert Opin. Drug Deliv. 2007, 4, 287-295.
- ^[4] Mortier et al, EuroInterv.2008;4,391-396


The new high tech Cobalt Chromium Stent System

A CoCr Stent -
designed for all lesions


Bifurcation stenting

With a cell circumference of up to 18,5mm the Yukon CC allows perfect side branch access which is essential for Bifurcation stenting ^[4].

2-Connector-Design

The new design with only 2 connector struts guarantees maximum flexibility and side branch access.


The suitable catheter...

The luer is designed with an integrated kinking protection.

The flexible tip ensures perfect crossability and trackability.

Technical data & Stent features

Cobalt Chromium Alloy (L605)

Surface containing micro-pores

Crossing profile (Ø 2,5 mm)	0,035" / 0,89 mm
Strut thickness (Ø 2,5 mm)	0,0027" / 68µm (SV), 0,0031" / 79µm (MV)

Metallic surface area	9,1 - 14,9%
Balloon marker material	Platinum / Iridium
Entry profile	0,016" / 0,41 mm
Proximal shaft diameter	1,9 F
Distal shaft diameter	2,7 F
Recommended guide wire	0,014"
Guiding Catheter	min. 5 F

Yukon[®] CC

Ordering information*

Small vessel design (SV)

product matrix

Balloon Ø [mm]	Stent length [mm] & Article number							
	8	12	16	18	21	24	28	32
Ø 2,00	T-CMG2008	T-CMG2012	T-CMG2016	T-CMG2018	T-CMG2021	T-CMG2024	T-CMG2028	T-CMG2032
Ø 2,25	T-CMG22508	T-CMG22512	T-CMG22516	T-CMG22518	T-CMG22521	T-CMG22524	T-CMG22528	T-CMG22532
Ø 2,50	T-CMG2508	T-CMG2512	T-CMG2516	T-CMG2518	T-CMG2521	T-CMG2524	T-CMG2528	T-CMG2532

Medium vessel design (MV)

Balloon Ø [mm]	Stent length [mm] & Article number								
	8	12	16	18	21	24	28	32	40
Ø 2,75	T-CMG27508	T-CMG27512	T-CMG27516	T-CMG27518	T-CMG27521	T-CMG27524	T-CMG27528	T-CMG27532	T-CMG27540
Ø 3,00	T-CMG3008	T-CMG3012	T-CMG3016	T-CMG3018	T-CMG3021	T-CMG3024	T-CMG3028	T-CMG3032	T-CMG3040
Ø 3,50	T-CMG3508	T-CMG3512	T-CMG3516	T-CMG3518	T-CMG3521	T-CMG3524	T-CMG3528	T-CMG3532	T-CMG3540
Ø 4,00	T-CMG4008	T-CMG4012	T-CMG4016	T-CMG4018	T-CMG4021	T-CMG4024	T-CMG4028	T-CMG4032	T-CMG4040

* Please contact our customer service for available sizes

Compliance chart

Ballon Ø [mm]	Inflation pressure [bar or 10 ⁻⁵ Pascal]																
						NP						RBP					
	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20		
Ø 2,00	1,83	1,87	1,90	1,93	1,96	2,00	2,03	2,06	2,10	2,13	2,16	2,20	2,23	2,26	2,29		
Ø 2,25	2,08	2,11	2,14	2,18	2,21	2,25	2,28	2,31	2,35	2,38	2,42	2,45	2,48	2,52	2,55		
Ø 2,50	2,33	2,36	2,40	2,43	2,47	2,50	2,53	2,57	2,60	2,64	2,67	2,70	2,74	2,77	2,81		
Ø 2,75	2,58	2,61	2,65	2,68	2,71	2,75	2,78	2,81	2,85	2,88	2,91	2,94	2,98	3,01	3,04		
Ø 3,00	2,81	2,85	2,89	2,92	2,96	3,00	3,04	3,07	3,11	3,15	3,18	3,22	3,26	3,29	3,33		
Ø 3,50	3,29	3,34	3,38	3,42	3,46	3,50	3,55	3,59	3,63	3,67	3,71	3,76	3,80	3,84	3,88		
Ø 4,00	3,75	3,80	3,85	3,90	3,95	4,00	4,06	4,11	4,16	4,21	4,26	4,31	4,36	4,41	4,46		

CE 0124

© 2010 Translumina GmbH

All Rights Reserved

® YUKON, PEARL Surface and translumina are registered trademarks of translumina GmbH

Please refer to the Instruction for Use supplied with these devices for indications, contraindications, adverse effects, suggested procedures, warnings and precautions.

Translumina GmbH | Neue Rottenburger Strasse 50
72379 Hechingen | Germany | info@translumina.de |
t + 49 7471 98 94 - 0 | f + 49 7471 98 94 - 360

EC Design-Examination Certificate

according to directive 93/42/EEC,
annex II (4)

As a notified body of the European Union, DEKRA Certification GmbH certifies for the manufacturer
translumina GmbH

Neue Rottenburger Straße 50, 72379 Hechingen, Germany

that the design dossier for the product(s) described in the annex complies with the requirements of the directive 93/42/EEC. This certificate is based on the result of the examination of the design dossier according to the directive 93/42/EEC annex II.4 as documented in the report mentioned in the Annex.

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Annex to the EC design examination certificate 50668-23-P2 dated 2016-02-23

Revision status: 0

Date: 2016-03-21

Page 1 of 1

Report number: 50668-P16-02

Product: YUKON CC Stent System

Intended use:

The YUKON CC Stent System is indicated for intra-luminal chronic placement in stenosed coronary artery or aorto-coronary bypass grafts, or in the anterior and posterior tibial artery with a maximum lesion length of 45mm, in order to obtain vessel patency following acute or subacute artery obstruction. It is also indicated for restenosis or arterial dissection after PTCA.

Technical data:

Stent	YUKON CC		
Stent Material	L605 Cobalt-Chromium Alloy with PEARL surface		
Length (mm)	8, 12, 16, 18, 21, 24, 28, 32, 40		
PTCA Catheter	Cathy No. 4		
Balloon Diameter (mm)	2,0 / 2,25 / 2,5 / 2,75 / 3,0 / 3,5 / 4,0		
Balloon Length (mm)	9,5 / 14,0 / 17, 0 / 19,5 / 22,5 / 25,5 / 30,0 / 34,0 / 41,5		
Article Number	T-CMG2008, T-CMG22508, T-CMG2508, T-CMG27508, T-CMG3008, T-CMG3508, T-CMG4008	T-CMG2012, T-CMG22512, T-CMG2512, T-CMG27512, T-CMG3012, T-CMG3512, T-CMG4012	T-CMG2016, T-CMG22516, T-CMG2516, T-CMG27516, T-CMG3016, T-CMG3516, T-CMG4016
	T-CMG2018, T-CMG22518, T-CMG2518, T-CMG27518, T-CMG3018, T-CMG3518, T-CMG4018	T-CMG2021, T-CMG22521, T-CMG2521, T-CMG27521, T-CMG3021, T-CMG3521, T-CMG4021	T-CMG2024, T-CMG22524, T-CMG2524, T-CMG27524, T-CMG3024, T-CMG3524, T-CMG4024
	T-CMG2028, T-CMG22528, T-CMG2528, T-CMG27528, T-CMG3028, T-CMG3528, T-CMG4028	T-CMG2032, T-CMG22532, T-CMG2532, T-CMG27532, T-CMG3032, T-CMG3532, T-CMG4032	T-CMG27540, T-CMG3040, T-CMG3540, T-CMG4040

EC Declaration of Conformity

By the responsible manufacturer:

translumina GmbH
Neue Rottenburger Str. 50
72379 Hechingen
Germany

Zertifikat - Certificate - Certificat

Translumina GmbH declares under its sole responsibility that the below entitled Medical Device complies to the requirements of the Medical Device Directive 93/42/EEC of the European Union.

The Quality Management System is in conformity with the requirements and the Conformity Assessment Procedure was followed in accordance to MDD 93/42/EEC, Annex II.

The device **YUKON CC Stent System** in all variations is classified as Class III as per Annexure IX, Rule 8 of MDD 93/42/EEC.

Issue of this certificate:

YUKON CC Stent System

Ref. No.:

T-CMG2008	T-CMG22508	T-CMG2508	T-CMG27508	T-CMG3008	T-CMG3508	T-CMG4008
T-CMG2012	T-CMG22512	T-CMG2512	T-CMG27512	T-CMG3012	T-CMG3512	T-CMG4012
T-CMG2016	T-CMG22516	T-CMG2516	T-CMG27516	T-CMG3016	T-CMG3516	T-CMG4016
T-CMG2018	T-CMG22518	T-CMG2518	T-CMG27518	T-CMG3018	T-CMG3518	T-CMG4018
T-CMG2021	T-CMG22521	T-CMG2521	T-CMG27521	T-CMG3021	T-CMG3521	T-CMG4021
T-CMG2024	T-CMG22524	T-CMG2524	T-CMG27524	T-CMG3024	T-CMG3524	T-CMG4024
T-CMG2028	T-CMG22528	T-CMG2528	T-CMG27528	T-CMG3028	T-CMG3528	T-CMG4028
T-CMG2032	T-CMG22532	T-CMG2532	T-CMG27532	T-CMG3032	T-CMG3532	T-CMG4032
			T-CMG27540	T-CMG3040	T-CMG3540	T-CMG4040

Certificate Registration Number 50668-23-P2

Notified Body's No. 0124
Notified Body DEKRA Certification
Handwerkstrasse 15
70565 Stuttgart
Germany

This certificate is valid from **2016-03-21** to **2021-03-20**

Hechingen, 2016-03-21

Juergen Ettenhuber
Manager QA/RA

Joerg Weber
Manager R&D

[Vers.2014/01]

This document has to be maintained at least five years after the last production of the concerned product was achieved.

Yukon[®] CC

Cobalt Chromium Stent System

Pre-mounted stent system

Vormontiertes Stent System

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Stent premontato

Sistema de stent pre-montado

Stent system pré-montado

INSTRUCTIONS FOR USE

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MODE D'EMPLOI

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
Table of Contents:

1. Description	3
2. Indications	3
3. Contraindications	4
4. Warnings and Precautions	4
5. Packaging and Storage	6
6. Specifications	6
7. Application	6
8. Potential Complications	10
9. Warranty and Limitation of Liability	11

These instructions are intended for all personnel who will handle the „YUKON® CC“, a premounted stent system, and should be read carefully before the product is used!

Product Label – Symbol Description

	Manufacturer		Batch code
	Date of manufacture		Catalogue number
	Expiry Date		Sterilized using ethylene oxide
	Do not use if packaging is damaged		Consult instructions for use
	Do not re-use		Caution
	Keep away from sunlight		Keep dry
	Nominal diameter of stent		Rated burst pressure
	Nominal length of stent		Nominal pressure
	Do not resterilize		MR Conditional

Note:

These instructions for use pertain to all available YUKON® CC Stent System lengths and diameters.

Information on Sterilization:

The product has been sterilized with ethylene oxide. It is non-pyrogenic. Do not use if the packaging has been opened, damaged or the expiration date is exceeded!

1. Description

The YUKON® CC Stent System consists of two parts: the pre-mounted stent and the PTCA Catheter Cathy No 4 as a delivery system.

1.1 Stent

The balloon-expandable stent used here has been designed for chronic placement in order to keep the vessel open and to stabilize it. The stent is made of an implantable cobalt-chrome alloy (L-605). It has been precision laser cut from one piece of tube.

1.2 PEARL surface

The term PEARL surface stands for a microporous surface developed by Translumina. This surface results in a number of advantages: on the one hand, endothelial cells are to be able to adhere better, which results in quicker growth of the stent into the vascular wall. On the other hand, this specific surface grants the possibility of purposeful administration of medication against restenosis. Various clinical studies have shown that application of medication to a stent surface can clearly reduce the complication rate after implantation, in particular the restenosis rate.

1.3 Delivery System

The Cathy No 4 PTCA catheter has a shaft system with an integrated balloon near the distal tip. The shaft comprises a single-lumen hose, which divides into two lumina in the area of the distal tip; via the short distal lumen, the catheter is inserted over the guide wire, which enables the precise position of the catheter in the area of the stenosis. The balloon is inflated with contrast medium via the continuous lumen. The balloon has been designed to provide an exact defined diameter at a specific pressure (semi-compliant balloon material).

In order to facilitate the positioning of the balloon in the stenosis, a number of markers have been attached. Two radiopaque platinum/iridium markers state the balloon length. Two further markers on the proximal catheter shaft aid positioning of the dilatation catheter in relation to the tip of the guiding catheter. The marker nearer to the Luer-Lock connection of the dilatation catheter is used for femoral application, the other one for brachial application. The delivery system has not been designed for distal injections and measurements, i.e. it does not possess a further lumen for this.

2. Indications

- The YUKON® CC Stent System is indicated for intraluminal chronic placement in stenosed coronary artery or aortocoronary bypass grafts or the tibial anterior and posterior artery of the lower limb with a maximum lesion length of 45mm, in order to obtain vessel patency following acute or subacute artery obstruction. It is also indicated in restenosis or arterial dissection after PTCA procedures.
- Patients being considered for stent implantation should be acceptable candidates for balloon angioplasty.

3. Contraindications

3.1 The coronary application of the YUKON® CC is contra-indicated for:

- Patients with hemorrhagic diathesis or other disorders in which the application of anticoagulant or antiplatelet therapy is contraindicated (e.g. septic ulceration, cerebrovascular complications)
- Patients in whom aortocoronary bypass procedures are contraindicated
- Patients with imminent thrombus formation and alteration of flow after myocardial infarction
- Patients who have experienced a myocardial infarction less than one (1) week prior to the proposed stent implantation
- Patients with proximal atherosclerosis and/or significant vessel tortuosity in whom guide wire access and adequate guiding catheter support is prohibited
- Coronary artery spasm in the absence of a significant stenosis
- Patients with unprotected left main coronary vessel disease
- Lesions that cannot be successfully pre-dilated
- Patients with known allergies to the used materials
- Patients in whom the necessary pharmaceutical treatments are contraindicated
- Severe allergic reaction to the contrast medium, stent and/or delivery system material and/or the necessary pharmaceutical treatment (e.g. aspirin)
- Women during pregnancy

3.2 The peripheral application of the YUKON® CC is contra-indicated for:

- Limited Access to the target vessel
- Diffuse vessel disease
- extremely calcified lesions
- angiographic indication of a massive thrombus amount next to or inside the target lesion
- lesions located in or next to an aneurysm
- extremely tortuous vessels
- extremely angled lesion
- contra-indication against thrombocyte aggregation inhibitors and/or coagulation inhibitors
- Allergy against contrast medium and/or stainless steel, nickel, chrome alloys and/or delivery system material and/or the necessary pharmaceutical treatment

4. Warnings and Precautions

4.1 In General

The medical techniques and procedures described here do not entail all the procedures of PTCA acknowledged in medicine.

- Use the YUKON® CC Stent System before the expiry date noted on the product label.
- Before use, make sure that the packaging has not been damaged or opened.
- The functionality of the stent and the delivery system should be checked prior to use. Ensure that the stent size and length are suitable for the procedure.
- The delivery system is only intended for single use. Do not resterilize and/or reuse it, as this can potentially result in damage and an increased risk of a lack of sterility and contamination.
- The YUKON® CC Stent System should only be used by physicians trained in the techniques of angioplasty and stent implantation.
- Administer appropriate anticoagulation and vasodilation therapy before insertion of the YUKON® CC Stent System (see 7.2). Continue the anticoagulation treatment after the implantation in order to avoid complications and undesired effects.

4.2 Delivery System

- Only use an inflation device with a manometer to inflate the balloon.
- Do not use air or any gaseous medium to inflate the balloon.
- Only use diluted contrast medium.
- The balloon pressure must not exceed the Rated Burst Pressure (RBP) specified on the product label.
- The inflation diameter of the balloon must not exceed the diameter of the artery proximal and distal to the stenosis.
- Do not overtighten the hemostatic valve, as lumen constriction may occur, possibly impairing inflation/deflation of the balloon.
- Do not attempt to straighten a catheter if the shaft is bent or kinked, as it may break. Instead, prepare a new catheter.
- Do not expose the delivery system to organic solvents (e.g. alcohol).

4.3 Stent

- Do not attempt to remove or readjust the stent from/on the delivery system. The stent may not be removed and used on another balloon catheter.
- Do not touch the stent as this could damage the fixation of the stent.
- A stent is an implant and should be handled with maximum care. Avoid contamination by using common sterile working techniques.
- Stenting of patients who are not acceptable candidates for coronary vessel bypass graft surgery can be particularly risky. Therefore, it is necessary to include possibilities for risk minimization, e.g. hemodynamic support during stenting.
- The long-term effectiveness of stent implantation has not been fully established.

4.4 During Implantation

- PTCA should be performed only in hospitals where an emergency surgery of the coronary artery bypass graft can be performed immediately in the event of injuries or life threatening complication.
- Insert the stent or the delivery system into the vessel only under fluoroscopic observation.
- Manipulate the YUKON® CC Stent System forward and backward only when the balloon has been fully deflated under vacuum.
- If the balloon of the delivery system gets caught in the stent during removal of the balloon, it can only be removed with cardiac surgery procedures.
- Never attempt to push the delivery system forwards by force. If you feel a resistance during movement, stop the procedure and only continue when you have found the reason. If this is not possible, then carefully pull the guiding catheter together with the delivery system and stent back as a unit. Do not attempt to retract the delivery system with stent through the guiding catheter as this may cause the stent to come off.
- After the stent has left the guiding catheter, it should not be pulled back into it, since it can become loose and an embolization could occur. If necessary, the stent is pulled back to the distal tip of the guiding catheter and the guiding catheter and the delivery PTCA catheter are removed through the introducer sheath as one unit. The guide wire must be left in its correct position.
- A bursting of the balloon during implantation may cause the migration or deformation of the stent. In such a case, the delivery system with balloon should be retracted carefully. A new balloon catheter with a small profile and semi-compliant balloon material should be inserted over the guide wire to replace it. The expansion can then be completed.
- A stent already implanted proximal to the lesion can be deformed by contact, e.g. with a catheter or a guide wire. In such a case, the correct expansion of the stent must be ensured.
- The re-expansion of a stent can cause dissection. Implantation of another stent can possibly be of help. If a number of stents are implanted, the distal stent should be implanted first.

- Distortions and thus artifacts can be caused in an MRI by the cobalt-chrome alloy.
- Mechanical atherectomy, hyperthermia and the use of laser devices are not recommended in the stented region.

4.5 Compatibility with MRI

Although the YUKON® CC Stent is made from a cobalt-chrome alloy and therefore does not manifest ferromagnetic reactions, it is conditionally MRI-safe.

4.6 Advice due to the single-use

The product „YUKON® CC Stent System“ is designed for the single-use only. The reconditioning for a re-use is prohibited. – The product is exposed high mechanical loads during its use and possibly affected thereafter. Because of that only for the single-use a save application can be warranted.

4.7 Peripheral implantation

The YUKON® CC Stent System is suitable for the treatment of peripheral vessel stenosis of the anterior and posterior tibial artery with lesion lengths of up to 45mm.

For the preparation and procedure please refer to the corresponding literature.

5. Packaging and Storage

Store in a dry, dark and cool place. Do not use it any more after the expiry date (indicated on the product label).

6. Specifications

6.1 Balloon Compliance

Tab.: Ratio: pressure/balloon diameter

The grey areas mark the diameter at Rated Burst Pressure.

Nominal pressure* = 11 [bar, atm or 10⁵ Pa]

Pressure [bar, atm or 10⁵ Pascal]

Balloon ø [mm]	6	7	8	9	10	11*	12	13	14	15	16	17	18	19	20
ø 2.00	1.83	1.87	1.90	1.93	1.96	2.00	2.03	2.06	2.10	2.13	2.16	2.20	2.23	2.26	2.29
ø 2.25	2.08	2.11	2.14	2.18	2.21	2.25	2.28	2.31	2.35	2.38	2.42	2.45	2.48	2.52	2.55
ø 2.50	2.33	2.36	2.40	2.43	2.47	2.50	2.53	2.57	2.60	2.64	2.67	2.70	2.74	2.77	2.81
ø 2.75	2.58	2.61	2.65	2.68	2.71	2.75	2.78	2.81	2.85	2.88	2.91	2.94	2.98	3.01	3.04
ø 3.00	2.81	2.85	2.89	2.92	2.96	3.00	3.04	3.07	3.11	3.15	3.18	3.22	3.26	3.29	3.33
ø 3.50	3.29	3.34	3.38	3.42	3.46	3.50	3.55	3.59	3.63	3.67	3.71	3.76	3.80	3.84	3.88
ø 4.00	3.75	3.80	3.85	3.90	3.95	4.00	4.06	4.11	4.16	4.21	4.26	4.31	4.36	4.41	4.46

* Pressure at which the balloon diameter reaches its nominal value

6.2 Material

- Catheter shaft proximal: PTFE coated hypo-tube
- Catheter shaft distal: Polyamide
- Balloon: Modified Polyamide
- Stent: Implantable cobalt-chrome alloy (L-605), micro porous surface

7. Application

Note:

The physician should consult recent literature on current medical practice on stent implantation, such as published by AHA (American Heart Association) and ACC (American College of Cardiology).

7.1 Accessories

- Introducer sheath, compatible with the guiding catheter
- Femoral or brachial guiding catheter in the appropriate size and configuration to match the coronary artery (5–9 F)
- Hemostatic valve
- 60 % contrast medium diluted 1:1 with heparinized physiological saline
- 20 cc Luer-Lock syringe (optional)
- Inflation device with manometer
- 3-way stopcock for a minimum of 25 bar
- 0.014" Guide Wire (Attention when using a hydrophilic coated guide wire, since the coating can swell and the GW can get stuck in the balloon catheter)
- Guide wire introducer
- Guide wire turning device

7.2 Concomitant Medication

The information provided should be regarded as a recommendation only. It is not intended to dictate medical practice. The conditions under which the stent has been used to date include:

- Administration of aspirin prior to stent implantation.
- Administration of heparin and nitro-glycerin during the procedure.
- Administration of antiplatelets alone or in combination with other anticoagulants up to some months following stent implantation.

Most experience with PTCA procedures has been gained with anticoagulant/antiplatelet therapy up to now. Administration of antiplatelets and anticoagulants depend on the patients' conditions and medical histories. Prior anticoagulant/antiplatelet administration is an essential adjunct to stent implantation. Proper heparin treatment before and after sheath removal is also critical to the overall success of the procedure. Before, during and after stent implantation, specific blood values (e.g. prothrombin time PT, partial thromboplastin time PTT, platelet count etc.) of the patient must be determined in order to create an ideal basis for medication. The selection, duration and dosage of the medication are in the responsibility of the physician treating the patient, as is selection of long-term medication after successful stenting.

7.3 Prior to use

- Please check the sterile packaging of the stent system and all the equipment necessary for this procedure carefully. Do not use, if the packaging is open or damaged.
- Do not use damaged or kinked instruments.
- Please remember that the balloon diameter of the delivery system at nominal pressure should not be larger than the artery diameter above and below the stenosis.
- Check the guide wire size according to the product label.

7.4 Preparation

1. Prepare the inflation device according to the manufacturer's instructions.
2. Prepare diluted contrast medium and sterile saline solution.

3. Fill a 10 or 20 cc syringe with 4 cc of contrast medium.
4. Carefully remove the stent system from the sterile packaging. Use usual sterile technique to do this. Pull the stent system out of the transport tubing. Remove the transport protection.

Caution:

Check the stent system once more, but do not touch the stent. The stent must have close contact to the balloon and it must be positioned precisely in between the balloon markings.

5. Connect the pressure-resistant 3-way stopcock (min. 25 bar) on the Luer connector to the balloon catheter.
6. Connect the syringe to the stopcock. Position the syringe vertically with the piston pointing downwards.
7. Create a strong vacuum by pulling the piston back. Maintain the vacuum for 15 - 20 seconds and make sure that no more bubbles can be seen in the diluted contrast medium.
8. Release the vacuum carefully.
9. Repeat this procedure, if necessary, until no more bubbles can be seen.

Warning

If it is not possible to create or to maintain a vacuum, check the connections. If you can not find the cause, do not use this stent system under any circumstances and return the unit to Translumina with a complaint report.

10. Create another vacuum in the catheter and close the 3-way stopcock.
11. Remove the air from the tubing and connector of the inflation device.
12. Connect the inflation device to the 3-way stopcock of the balloon catheter.
13. Open the stopcock carefully. The inflation device should be in a vertical position in order to avoid air entering the system.
14. Attach a needle to the syringe; insert it into the distal tip of the catheter and flush the guide wire lumen by applying gentle pressure to the syringe.
15. The stent system is now prepared and ready for use.

7.5 Pre-dilation

- We recommend dilating the stenosis with a suitable PTCA dilatation catheter prior to stenting. In general, dilatation of the vessel is required with a balloon to vessel diameter ratio of 1:1.
- For the pre-dilatation of the vessel you must expect atherosclerotic plaque, which may inhibit advancement of the catheter.

7.6 Implantation of the stent

1. Insert a compatible introducer sheath in the femoral or brachial artery using conventional techniques.
2. Insert the guiding catheter into the vessel through the introducer sheath and fix the hemostatic valve to it.
3. Insert the guide wire into the guiding catheter and carefully push it forwards under fluoroscopic control until it is distal to the stenosis. Please pay attention to the admissible guide wire size on the product label.
4. After the guide wire has passed the stenosis, close the hemostatic valve tightly around the guide wire to prevent displacement of the wire.
5. Push the proximal tip of the guide wire back into the distal tip of the stent system, ensuring that the guide wire exits through the side-port located approximately 25 cm from the balloon catheter tip.
6. Open the hemostatic valve and carefully push the YUKON® CC Stent System forwards over the guiding wire to the distal tip of the guiding catheter.
7. Continue the procedure under fluoroscopic control using accepted angioplasty techniques. The radiopaque markers facilitate the positioning of the balloon in the area of the stenosis.

8. When placed in the right position, the stent system needs to be inflated, first with nominal pressure. Keep this pressure for 10 to 90 seconds. Due to the restricted blood flow, the patient can feel heartpain. Adapt the inflation time to the patient's condition.
9. The inflation can be repeated two to three times, but do not exceed the rated burst pressure (RBP). Check the expansion of the stent as long as the balloon catheter has not been retracted after deflation to minimize the risk of the balloon catching in the stent struts (see 7.8)

Caution:

- Do not retract the balloon catheter before you have made sure that the stent has been implanted correctly. A re-expansion with another balloon catheter is possible, but harbors additional risks.
- Do not overtighten the hemostatic valve, since this may restrict the flow of contrast medium into the balloon and thereby increase the inflation/deflation time.
- Do not overstretch the stent, and do not exceed the Rated Burst Pressure (RBP) of the balloon displayed on the product label, as this may cause dissection. Permanently check the correct expansion of the stent with a high resolution fluoroscopic device.
- Expansion of the stent should not be undertaken if the stent has not been properly positioned in the area of the stenosis. If the position of the stent is not optimal, it should be repositioned before inflation.

Please Note:

- It is necessary for the guide wire and/or the balloon catheter to remain in the stenosis until the procedure has been completed.
- Contrast media may effect the inflation/deflation time due to different viscosities.
- Multiple angiographic views and measurements should be taken in order to be sure of a correct positioning of the stent.

7.7 Further dilation in the stented segment

- If the primary angiographic result is not optimal, you can further expand the stent with the help of a semi-compliant balloon catheter.
- If a guide wire has to be pushed through a stent which has already been expanded, its distal tip ought to have an I-shape or should be prepared in this way. In this way, the guide wire can be guided through the cross-section of the stent more easily without getting caught in the struts of the stent.

7.8 Balloon Catheter Withdrawal Procedure

1. Open the hemostatic valve.
2. Hold the guide wire and the hemostatic valve firmly in one hand, while grasping the balloon catheter-shaft with the other hand. Maintain the guide wire's position in the vessel by holding the wire stationary and begin pulling the dilatation catheter carefully out of the guiding catheter.
3. Retract the balloon catheter until the opening in the guide wire lumen is reached. Then proceed in the same way as with an over-the-wire system until the delivery system has been totally removed.
4. Close the hemostatic valve.

Caution:

- Before retracting the balloon catheter, make sure that the balloon has been completely deflated.
- Observe the patient regularly during the first 30 minutes after stent placement, and perform an angiographic examination. If stent placement is associated with the onset of thrombus in the stented region, intracoronary administration of a thrombolytic agent is recommended.

8. Potential Complications

8.1 Potential complications of the coronary application:

The following complications due to PTCA and stent implantations have been reported and may occur. Potential undesired effects are possible, but not limited to the following:

- Dissection, injury, rupture, sudden and total occlusion or other damages to the coronary arteries
- Thrombosis of the coronary arteries
- Unstable angina
- Arterial spasm
- Myocardial infarction
- Cardiac dysrhythmia, including ventricular fibrillation
- Embolism
- Restenosis of the dilated artery
- Hemorrhage or hematoma
- Bleeding from anticoagulation/antiplatelet medication
- Side effects of the concomitant medication
- Allergic reaction to the contrast medium
- Allergic reaction to the stent and/or delivery system material
- Infections
- Arteriovenous fistula
- Emergency surgery of coronary vessel bypass graft
- Deformation of stent symmetry on deployment or on re-crossing with another device (IVUS, balloon catheter, etc.)
- Stent migration
- renal failure
- Death

8.2 Potential complications of the peripheral application:

- Sepsis / infection
- stent migration
- embolism of atherosclerotic or thrombotic material
- acute or sub-acute stent thrombosis
- amputation
- arterio-venous fistula
- Emergency surgery to treat vascular complications
- Hemorrhage or hematoma
- Side effects of the concomitant medication
- Allergic reaction to the contrast medium
- Allergic reaction to the stent and/or delivery system material
- injury of the artery including perforation and dissection
- restenosis of the stented artery
- Thrombus
- tissue necrosis
- total occlusion
- renal failure

9. Warranty and Limitation of Liability

translumina GmbH guarantees that this product has been manufactured, sterilized and packaged with the due care prescribed in EU Directive for Medical Devices 93/42/EEC. Each product is individually tested before it is released for packaging. **translumina GmbH** will replace any device if defects which can be proven to have originated during manufacture or packaging are immediately brought to the attention of **translumina GmbH** or its distributors. This warranty applies in lieu of and ruling out any other verbal or written, express, implicit or other warranty.

As a result of biological differences in individuals, no product is 100% effective under all circumstances. Due to this fact and since we have no control over the conditions, the instruments used, the operation and the diagnosis, **translumina GmbH** and its distributors do not guarantee good effectivity in application; unsatisfactory results also cannot be ruled out. **translumina GmbH** and its distributors shall not be liable for any, damage or losses of any nature arising directly or indirectly from the use of the device. The product liability expires if any damage caused by transport, storage or improper handling by the operator, or other factors outside the control of the manufacturer occurs. **translumina GmbH** and its distributors shall not be liable for damages arising from resterilisation or reuse of the product.

Important Information:

If more than 2 years have elapsed between the revision date stated below and today's date, please contact the manufacturer in order to find out if any changes have occurred.

For further information please contact:

translumina GmbH

Neue Rottenburger Strasse 50
D-72379 Hechingen
Phone: ++ 49 74 71 98 94 - 0
Fax: ++ 49 74 71 98 94 - 3 80
E-mail: info@translumina.de
Internet: www.translumina.de

© 2015 translumina GmbH
All rights reserved.

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® YUKON® and translumina® are registered trademarks of translumina GmbH