Instruction for use

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**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUANTITATIVE DETERMINATION
OF CA15.3 (M12) IN HUMAN SERUM OR PLASMA**

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
Authorized Representative in EU:
Polmed.de
Steinacker 20, D-73773
Aichwald, Germany
e-mail: info@polmed.de

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CONTENT

1. INTENDED USE	10
2. SUMMARY AND EXPLANATION	10
3. PRINCIPLE OF THE TEST	11
4. WARNINGS AND PRECAUTIONS	11
5. KIT COMPONENTS	12
6. SPECIMEN COLLECTION AND STORAGE	13
7. TEST PROCEDURE	13
8. QUALITY CONTROL	15
9. CALCULATION OF RESULTS	15
10. EXPECTED VALUES	15
11. PERFORMANCE CHARACTERISTICS	16
12. LITERATURE	16

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1. McGuire WL, Tandon AK, Allred D, Chamnes GC, Clark GM. How to use prognostic factors in axillary node negative breast cancer patients. J Natl Cancer Inst 1990; 82:1006-7.
2. Nicholson S, Richard J, Sainsbury C, et al. Epidermal growth factor receptor (EGFr): results of a 6-year follow up study in operable breast cancer with emphasis on the node-negative subgroup. Br J Cancer 1991; 63:146-50.
3. Somerville JE, Clarke LA, Biggart JD. C-erb B-2 overexpression and histological type of in-situ and invasive breast carcinoma. J Clin Pathol 1992;45-16-20.
4. Ueronese S, Gambacorta M. Detection of Ki-67 rate in breast cancer. AM J Clin Pathol 1991; 95:30-4.
5. Lotnickier M, Pavesi F, Scarabelli M. Tumor associated antigens CA15-3 and CA125 in ovarian cancer. Int. J. Biolog Markers 1991; 6:115

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF CA15.3 (M12) IN HUMAN SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of CA15.3 (M12) in serum or plasma.

This kit is designed for measurement of CA15.3 (M12) in serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

CA15-3 or MUC1 is a heterogenous glycoprotein with a molecular mass ca. 300-450 kD. Elevation of serum CA15-3 is associated with mammary carcinomas. Quantitative determination of CA15-3 in serum and plasma is helpful in monitoring of patients with such tumours to estimate the course of the disease, effectiveness of its treatment and to reveal recurrence or metastases. However, CA15-3 values obtained should always be interpreted in context of results obtained by other diagnostic and clinical procedures. Besides mammary carcinomas, CA15-3 levels in blood may rise in lung tumours, prostate cancer, ovarian carcinomas, gastro-intestinal tumours. Elevation of CA15-3 level in blood can be also found in benign tumours of the mammary gland and the ovary, endometriosis, hepatitis, liver cirrhosis and lung fibrosis. Pregnancy and lactation may also cause elevation of CA15-3 level in serum.

All CA15-3 test systems are usually not very sensitive (not more than 75% even at stage III mammary carcinoma). Therefore, in monitoring of tumours, we recommend to use this test in conjunction with two other test systems designed by XEMA for diagnostics and monitoring of mammary carcinomas – M20 (analogous to BR27.29 developed by Biomira) and M22 (analogous to MCA, Roche). All three systems should be used to evaluate MUC1 concentration before and after surgery; the test system showing the most pronounced postsurgery decline should be then used for further monitoring.

NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy may develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies gives false results in the present assay. Sera from HAMA positive patients should be treated by depleting adsorbents before assaying.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human CA15.3-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human CA15.3, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5.0% H_2SO_4 . It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	CA15.3 (M12) EIA strips, 8x12 wells	1	pcs		until exp.date
2 CAL 1-5	Calibrator set, 0.8 ml each. The set contains 5 calibrators: 0: 12.5; 50: 125; 250 U/ml	5	pcs	red (C1 - colourless)	2 months
3 CONTROL	Control serum (0.8 ml)	1	pcs	colourless	2 months
4 CONJ HRP	Conjugate, 14 ml	1	pcs	red	until exp.date
5 DIL	EIA buffer 14 ml	1	pcs	blue	until exp.date
6 SUBS TMB	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7 BUF WASH 26X	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate - until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8 STOP	Stop solution, 14 ml	1	pcs	colourless	until exp.date
9 N003	Plate sealing tape	2	pcs		N/A
10 K2261	Instruction CA15.3 (M12) EIA	1	pcs		N/A
11 K226Q	QC data sheet CA15.3 (M12) EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 50–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Microtiter plate shaker. Shaking should be medium to vigorous. Longitudinal shaking approximately 200 strokes/min, oscillations 600–800/min;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL (EIA buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
3	Pipet 50 µl of EIA buffer into each well.
4	Pipet 50 µl of calibrators CAL 1-5, control samples CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at +18...+25 °C and continuous shaking at 600-800 rpm.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
8	Incubate 30 minutes at +18...+25 °C and continuous shaking at 600-800 rpm.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10-20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction.

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus CA15.3 (M12) concentration.

9.3. Determine the corresponding concentration of CA15.3 (M12) in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 U/ml	0.06
CAL 2	12.5 U/ml	0.31
CAL 3	50 U/ml	0.84
CAL 4	125 U/ml	1.68
CAL 5	250 U/ml	2.26

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for CA15.3 (M12). Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, U/ml	
	Lower limit	Upper limit
Males	-	30
Females	-	30
Pregnancy:		
1st trimester	-	55
2nd trimester	5.0	65
3rd trimester	5.0	185
Lactation	-	120

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
CEA	<0.1
CA125	<0.1
CA19.9	<0.1

11.2. Analytical sensitivity

Sensitivity of the assay was assessed as being 0.75 U/ml.

11.3. Linearity

Linearity was checked by assaying dilution series of 5 samples with different CA15.3 (M12) concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery

Recovery was estimated by assaying 5 mixed samples with known CA15.3 (M12) concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

- McGuire WL, Tandon AK, Allred D, Chamnes GC, Clark GM. How to use prognostic factors in axillary node negative breast cancer patients. J Natl Cancer Inst 1990; 82:1006-7.
- Nicholson S, Richard J, Sainsbury C, et al. Epidermal growth factor receptor (EGFr): results of a 6-year follow up study in operable breast cancer with emphasis on the node-negative subgroup. Br J Cancer 1991; 63:146-50.
- Somerville JE, Clarke LA, Biggart JD. C-erb B-2 overexpression and histological type of in-situ and invasive breast carcinoma. J Clin Pathol 1992;45:16-20.
- Ueronese S, Gambacorta M. Detection of Ki-67 rate in breast cancer. AM J Clin Pathol 1991; 95:30-4.
- Lotnickner M, Pavesi F, Scarabelli M. Tumor associated antigens CA15-3 and CA125 in ovarian cancer. Int. J. Biolog Markers 1991; 6:115

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Russian Diagnostic
Manufacturers Association

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RUSSIAN ASSOCIATION
OF MEDICAL LABORATORY
DIAGNOSTICS

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