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ITEM: **15ml conical graduated test tubes - Sterile**

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Sterile test tubes, conical bottom for a complete sample collection. Manufactured in medical PP, particularly transparent. Screw cap which guarantees the perfect leak-proof, with milled rim allowing a safe and easy opening/closing. Accurate graduation, clearly visible. Excellent high and low temperatures, resistance to chemical agents and organic solvents. Resistant to phenol/chloroform. Autoclavable. R.C.F. max: 5,000 x g. Pyrogen and latex free.

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CE Marked product - manufactured in compliance with 98/79/CE Directive and D.lgs 332 dtd 08/09/2000

CARATTERISTICHE PRINCIPALI		TECHNICAL FEATURES
Stato microbiologico	STERILE / STERILE RAGGI BETA / BETA RAYS	Microbiological status
Materiale impiegato provetta	POLIPROPILENE / POLYPROPYLENE	Raw material – Test tube
Materiale impiegato tappo	POLIPROPILENE / POLYPROPYLENE	Raw material - cap
Temperature tollerate provetta	MIN -10°C MAX +121°C	Temperature range - test tube
Temperature tollerate tappo	MIN -10°C MAX +121°C	Temperature range - cap
Scala graduata	MIN. 1ML – MAX 14ML	Graduated scale
Incrementi scala graduata (ml)	0,5	Intervals graduated scale (ml)
Dimensioni provetta (mm)	Ø 17 x 120	Dimensions - test tube (mm)
Dimensione tappo (mm)	Ø 23 x 11	Dimensions – cap (mm)
Peso provetta (gr.)	5,1	Weight - test tube (gr.)
Peso tappo (gr.)	1,2	Weight – cap (gr.)
Validità del prodotto	5 ANNI / YEARS	Shelf life

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Classificazione EDMA > 51091001 - Other containers for samples of human origin

Basic UDI-DI (Global Model Number – GMN): 805577609FT02STETUBE01Y8

*In Vitro Diagnostic Medical Device "specimen receptacles" adapted to contain a human biological sample (blood, urine, faeces, pus, sputum, saliva, semen, cells, nails, hair) to perform in vitro diagnostic analysis laboratory. **For professional use only.***

National classification of medical devices (CND) > W050301020102 (Samples analyses, plastic tubes without additives).

EDMA code > 51091001 - Other containers for samples of human origin

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- For professional use only -

Use only for laboratory diagnostic analysis of human biological samples in healthcare settings.

Device without measuring function.

Any presence of a graduated scale is intended for indicative purposes and without a metrological function.

Do not vary the intended purpose of the device. Please contact the Manufacturer in the event of doubts concerning the intended purpose.

Before use check the perfect closure of the device if applicable.

Handle using personal protective equipment: risk of contamination.

After collection, store the sample as prescribed by the analytical method and deliver to the Laboratory within the maximum time specified by the same.

Do not use the Device in case of evident signs of breakage, injuries, cracks that could compromise its correct use.

Keep out of flame or heat sources which might damage the device.

Do not use after expiry date or if packing is opened.

Do not re-use: Disposable Device. Any reuse could cause sample contamination.

If sterile (ionizing radiations) or aseptic Device: Microbiological status in undamaged pack.

Do not re-sterilize.

Use the Device only with compatible and / or supplied accessories (caps, supports, etc.)

Do not put under direct sun rays; store in a dry, cool place (U.R. max 75% at 26°C)

Store in dry place, temperature range: min. -10°C (14°F) max +50°C (122°F)

Disposal according to applicable regulations, risk of infection.

Before use with particular substances check the resistance/compatibility chart in our catalogue.

Test tubes: centrifuge rotation speed (max 20 minutes, room temperature):

PMMA r.p.m 3.000 (R.C.F. min 1.500 x g) – PS r.p.m 3.000 (R.C.F. min 1.500 x g) – PP r.p.m 5.000 (R.C.F. min 4.000 x g) – TPX r.p.m 4.500 (R.C.F. min 3.200 x g).

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Aptaca S.p.A. Regione Monforte, 30 - 14053 Canelli (Asti) Italy