

SIEMENS



Data Sheet

# Artis one

Floor-mounted system for uncompromised imaging

[siemens.com/healthcare](https://www.siemens.com/healthcare)

# Artis one

## Floor-mounted system for uncompromised imaging

Artis one is an angiography system that leaves the beaten track. It marks a new approach to interventional imaging, because it's different in many ways.

This system delivers proven state-of-the-art technology – like the MEGALIX Cat Plus X-ray tube featuring a unique flat emitter. In addition it offers next-generation tools for uncompromised imaging, such as CLEARstent Live\* to verify stent positioning in real time.

Intelligent operation is enhanced by a configurable Head-up display, allowing you to interact with the system in a completely new, intuitive way. So you can keep your attention where you need it.

And because the solution is so easy to understand and deploy, it will have a positive impact on your whole organization.


Broaden your procedure mix and hit the sweet spot of your business.

Artis one is an angiography system developed for diagnostic imaging and interventional procedures including, but not limited to, pediatric and obese patients.

Procedures that can be performed with the Artis one include cardiac angiography, neuro-angiography, general angiography, rotational angiography, multipurpose angiography and whole body radiographic/fluoroscopic procedures as well as support of procedures next to the table for i.e. patient extremities or hybrid procedures.

Meet a system that is designed differently.

**Artis one. Designed around you.**



\* Option



# Artis one

Floor-mounted system for uncompromised imaging



Head-end position



Right-side table rotated



Left-side table rotated

# Artis one

## Floor-mounted system for uncompromised imaging

### Artis one – designed for uncompromised imaging, intuitive interaction and positive impact

From unique flat-emitter technology to real-time stent enhancement, Artis one offers proven technology with next-generation imaging tools. The CARE+CLEAR packages come standard, the perfect choice for high image quality at the lowest possible dose.

### Optimum system positions?

At the press of a button – or automatically – Artis one smoothly adapts to your procedure-specific needs. Keep the patient's head free for anesthesia and echography or have room to move during an emergency.

Artis one requires less space and provides large returns. Short installation time, easy deployment, and energy savings up to 20 %. Address a broader procedure mix – the smartest way to hit the sweet spot of your business.

### X-ray tube

The MEGALIX Cat Plus tube is the first angiography X-ray tube to offer the unique flat emitter technology. This enables a higher current during fluoroscopy resulting in better image quality especially for obese patients.

### Integrated 3D imaging\*

Parallel processing of the currently examined patient and one additional patient, as well as highspeed C-arm rotation of 60 degrees per second – with integrated 3D.

Two high contrast acquisition modes, quality-optimized with syngo DynaCT Cardiac untriggered, dose optimized with syngo Dyna3D.

**QuickZoom** allows to focus and zoom in the control room or at tableside with just one click and thereby helps to save time and speeds up the workflow.

### CARE+CLEAR

Combined applications to reduce exposure (CARE) help to reduce radiation dose for the operator and the patient.

The CLEAR applications provide outstanding image quality to increase certainty during interventions.

### CLEARstent\* and CLEARstent Live\*

CLEARstent features static stent enhancement allowing you to get a clear picture of implanted stents, whether it was 5 minutes or 5 months ago. CLEARstent Live takes this to real time and provides stent enhancement during an ongoing acquisition – while the stent is being positioned.

### HeartSweep\*

During one single sweep, all necessary angulations required for coronary diagnostics are acquired and help to find the optimal projection for treating a lesion. HeartSweep can follow up to 10 different trajectories with configurable movement speeds and simple one-button operation.

### Unparalleled coverage

Benefit from patient coverage up to 2.10 m. Peripheral run-offs without moving the patient. Imaging the whole body in one go, with ceiling-like imaging workflows on a floor-mounted system.

### Intelligent controls

Artis one features a new, unique heads-up display combined with tactile system operation which allows feeling and operating buttons even when hidden under a cover. This all helps to keep your attention where it's needed, and not on operating the system.

### Artis one display

The new Artis one 30" display offers a 90 % larger image display area than conventional 19" screens. It connects to up to 4 external image sources and users can switch between different layouts which can be freely configured for a specific procedure.

### Artis one Panoramic Display\*

The new Panoramic Display offers an additional 30" screen that is mounted right next to the standard Artis one display. It displays up to 9 additional external image sources in its own configurable layouts, thus keeps even more imaging information visible during the procedure.




\* Option


# Artis one

## Floor-mounted system for uncompromised imaging


### Ideal system positions for every procedure




Transfer position




Head-end position (e.g. for PCI)




Left-side position




Right-side, table-rotated position  
(e.g. for pacemaker implantations)




Peri-position




Left-side, OR position (e.g. for AAA  
treatment)



Extended parking position\*



OR position, right-side\*



Additional transfer position\*

\* Additional system positions, examples

# Artis one

## Floor-mounted system for uncompromised imaging

### Table of contents

<b>Imaging system</b>	<b>7</b>	<b>Connectivity</b>	<b>24</b>
X-ray tube	7	DICOM Functions	24
X-ray generator	8	Networking	24
Detector as30	9	Data export	25
Rotatable collimator	9	Sensis Interface	25
Operating modes	10	Injectors	26
Additional functions	12	Remote Service	26
CARE	13	Security Package	26
Imaging system	15		
CLEAR	15	<b>Room preparation</b>	<b>27</b>
General functions	15	Emergency power supply	27
		Internal line resistance for generator	27
<b>Advanced applications</b>	<b>16</b>	Installation data	28
Quantification	16	Weight	28
CLEARstent	17	Ambient conditions	29
CLEARstent Live	17	Mobile Artis one installations	29
3D Imaging	17	System view	30
LA Segmentation	17	Room layout	31
EP Option	17		
<b>System specifications</b>	<b>18</b>		
Stand	18		
Patient tables	19		
Tabletops	20		
DCS Artis one	21		
Displays	22		
StreamLink	23		
Operation in the examination room	23		
Operation in the control room	23		

# Artis one

## Floor-mounted system for uncompromised imaging

### Imaging System

#### X-ray tube

##### MEGALIX Cat Plus 125/40/90-125GW

- High-performance X-ray tube
- Up to 40% higher fluoro power with unique flat emitter technology
- Increased contrast during fluoroscopy, especially for examinations on obese patients
- Oil/water cooled

Max. exposure voltage (IEC 60613)	125 kV	
Focal spot (IEC 60336)	0.4 <sup>1)</sup>	0.8
Nominal power (IEC 60613) (thermal anode reference power = 300 W)	35 kW	90 kW (IEC 60613:1989)
Nominal power (thermal anode reference power = 0 W)	42 kW	112 kW (IEC 60613:1989)
Nominal radiographic anode input power	38 kW	99 kW (IEC 60613:2010)
Anode angle	8°	
Maximum anode heat content	2,500,000 J (3,375,000 HU)	
Maximum heat content of the X-ray tube assembly	3,600,000 J (4,900,000 HU)	
Maximum cooling capacity of the anode	400,000 J/min. (540,000 HU/min.) / 6667 W	
Continuous heat dissipation of the tube assembly	max. 2900 W	
Anode rotation	150/180 Hz (3-phase current)	
Max. anode current in fluoroscopy	250 mA small focus	
Max. anode current in acquisition	800 mA large focus	
Anode input power	10 min	4000 W
	20 min	3000 W
	> 30 min	2500 W
Total filtration (IEC 60601-1-3)	≥ 2.5 mm Al	
Leakage radiation (IEC 60601-1-3)	< 0.44 mGy/h (at 125 kV in 1 m distance: 2500 W)	
Weight	approx. 36 kg (79.4 lbs.)	

#### Cooling unit

Cooling medium	water (not distilled) with coolant additive
Cooling medium temperature	max. 55 °C
Max. pressure	3.1 bar
Flow rate	3.5 l/min
Weight (cooling system)	< 28 kg (61.73 lbs.) + 6.5 kg (14.33 lbs.) cooling liquid

<sup>1)</sup> With flat emitter technology



# Artis one

## Floor-mounted system for uncompromised imaging

### X-ray generator POLYDOROS A100 Plus

Microprocessor-controlled high-frequency X-ray generator with automatic dose rate control for fluoroscopy and acquisition

Multi-pulse converter frequency	100 kHz
Max. generator power (IEC 60601-2-7 and IEC 60601-2-54)	1000 mA at 100 kV $\triangleq$ 100 kW 800 mA at 125 kV $\triangleq$ 100 kW 800 mA at 100 kV $\triangleq$ 80 kW
Tube current	0.5 mA to 1000 mA in 0.01 mA steps 0.5 mA to < 800 mA <sup>1)</sup> in 0.01 mA steps
Tube current (pulsed fluoroscopy)	15 mA to 250 mA in 0.01 mA steps (small focus)
Pulse frequency	0.5 p/s to 66 p/s or continuous mode
Pulse time	0.5 ms to 800 ms
Max. continuous power in fluoro mode	3000 W
Tube voltage	40 kV to 125 kV in 0.1 kV steps



# Artis one

## Floor-mounted system for uncompromised imaging

### Detector as30 (midsize)

Amorphous silicon flat detector with 39 cm diagonal entrance plane

High-resolution a-Si matrix with 184  $\mu\text{m}$  pixel size and 16-bit digitization depth

Integrated collision sensor	yes
Removable grid	yes
Detector rotation	yes
Detector cooling	air cooled
Imaging size	29 cm x 26 cm
Image display matrix	1560 x 1420 pixels
Size incl. housing and collision protection	430 mm x 378 mm with collision protection
Input fields	39, 32, 25, 20, 16, 11 cm (15.35", 12.6", 9.84", 7.87", 6.3", 4.33")
X-ray conversion technology	a-Si with CsI scintillator
Digitization depth	16-bit (65536 gray scale levels)
Pixel pitch	184 $\mu\text{m}$
Nyquist frequency	2.7 lp/mm
DQE (detective quantum efficiency)	0 lp/mm at 2 $\mu\text{Gy}$ : min. 65 %, typical 70 % (RQA5) 1 lp/mm at 2 $\mu\text{Gy}$ : typical 52 % (RQA5)
MTF (modulation transfer function)	1 lp/mm: min. 53 %, typical 59 % (according to IEC 62220-1-3)
Signal to electronic noise ratio (SENR)	11 dB typical at 5 nGy (RQA5, 1x1 binning, high gain)

### Rotatable collimator

Compact collimator for angiography with rectangular blade and wedge-shaped filter for cardiological applications and DSA

Automatic synchronous rotation of the detector and collimator unit (internal rotation) to compensate for image rotation at different examination positions of the support stand; rotation also possible via remote control enabling upright images of objects or body parts not aligned with the table, e.g. arms (StraightView).

# Artis one

## Floor-mounted system for uncompromised imaging

### Operating modes

#### Fluoroscopy

Digital pulsed fluoroscopy, with 7.5, 10, 15, 30 p/s

Additional fluoroscopy pulse rates from 0.5, 1, 2, 3, 4 p/s (CAREVISION)

Roadmapping (requires DSA option) with automatic pixel shift

Store Monitor: Any image can be stored on the disk

Store Reference: Any image can be stored as a reference image, even during online fluoroscopy

Store Fluoro: 1024 images

Last Image Hold (LIH)

#### Overlay fade

Online superimposing of active fluoro and reference image (overlay reference)

#### Fluoro Loop\*

Storage and display of dynamic fluoro sequences

The maximum fluoro time that can be saved depends on the pulse frequency selected, e.g., 34 s at 30 p/s, 68 s at 15 p/s

#### Roadmap

Individual windowing of vessel map and tool image

#### Cardiac acquisition

Acquisition at 7.5, 10, 15 and 30 f/s, acquisition, display and storage in original matrix, 12-bit

#### DR – 0.5 - 7.5 f/s

Digital radiography with digital real-time filtering, applicable for single images and series with frame rates from 0.5 f/s to 7.5 f/s (to 30 f/s<sup>1)</sup>)

Acquisition, display and storage are performed in original matrix size at a resolution of up to 2.22 megapixels

Time-controlled and manually variable frame rates are included

#### DSA – 0.5 - 7.5 f/s\*

Digital subtraction angiography with digital real-time filtering, applicable for single images and series with frame rates from 0.5 f/s to 7.5 f/s (to 30 f/s<sup>1)</sup>)

Acquisition, display and storage are performed in original matrix size at a resolution of up to 2.22 megapixels

Remask, peak opacification for iodine contrast (MaxOpac) and CO<sub>2</sub> contrast (MinOpac), display of anatomical background (Landmark) from 0 to 100%

Time-controlled and manually variable frame rates are included

#### High-speed acquisition for DR and DSA\*

Acquisition at 10/15/30 f/s

Subtracted display possible only with DSA

# Artis one

## Floor-mounted system for uncompromised imaging

### Operating modes

#### Anatomical background <sup>1)</sup>

Anatomical surroundings visible by fading in the native image

#### Setting new mask <sup>1)</sup>

A new mask can be set with "Move Mask" or "Replace Mask"

#### Pixel shift <sup>1)</sup>

Manual pixel shift, automatic pixel shift, flexible pixel shift (rubber masking)

#### CLEARstent\*

Software for enhanced stent visualization, can be activated from tableside

#### CLEARstent Live\*

Real-time stent enhancement for facilitation of complex procedures

#### 3D Imaging\* for *syngo* DynaCT Cardiac or *syngo* Dyna3D

Allows native and subtracted 3D reconstruction based on digital rotational angiography.

Automatic 3D reconstruction and review in the integrated *syngo* 3D application.

Rotation speed up to 60°/s

Acquisition rate up to 66 f/s

#### PERISTEPPING\*

Peripheral digital angiography stepping of the stand without moving the patient, with a single contrast-medium injection performed while observing the contrast medium bolus

Position-dependent variable frame rates

Fully automatic exposure control

The collimator settings are automatically saved for each stepping increment

# Artis one

## Floor-mounted system for uncompromised imaging

### Operating modes

#### PERIVISION\*

Peripheral digital angiography with stepping of the stand without moving the patient and online subtraction display in one examination procedure with a single contrast-medium injection while observing the contrast medium bolus

One automatically acquired mask image for each individual position

Position-dependent variable frame rates

Fully automatic exposure control

The collimator setting is automatically saved for each stepping increment

#### ECG-triggered fluoroscopy and acquisition\*

ECG-triggered fluoroscopy/acquisition provides a still image of the catheter while compensating for cardiac movement. This enables the use of low pulse frequencies, resulting in a significantly lower dose compared to normal fluoroscopy/acquisition

#### HeartSweep\*

HeartSweep is a dual-axis, rotational angiography following a predefined, configurable trajectory in a single sweep. It supports efficient diagnosis of coronary disease and improves planning of the interventional therapy.

HeartSweep can also cover single-axis rotational angiographies, comparable to Dynavision DR.

Acquisition at 7.5, 10, 15 and 30 f/s, acquisition, display and storage in original matrix, 12-bit

Up to 128 acquisition programs per each mode for flexible adjustment of the X-ray and image processing parameters to the different procedures (selectable in the examination room and in the control room)

### Additional functions

#### ECG recording and storage\*

Recording, storage and display of an ECG waveform

ECG waveform displayed on the display with synchronous image information

# Artis one

## Floor-mounted system for uncompromised imaging

### CARE

Combined applications to reduce exposure (CARE) help to reduce radiation dose for the operator and the patient

#### CAREfilter

Five-level adaptive Cu prefiltration (CAREfilter) for reduction of skin dose; automatic selection control based on the absorption of the object

Filter levels 0.1, 0.2, 0.3, 0.6, 0.9 mm Cu

#### CAREvision

Pulsed fluoroscopy with additional reduced pulse frequencies of 0.5, 1.0, 2.0, 3.0, 4.0 p/s

Pulse frequency can be adjusted to the requirements of each application to significantly reduce radiation exposure, particularly during interventions

#### CAREprofile

Radiation-free positioning of primary and semi-transparent collimators via graphic display in the LIH image on the image display

#### CAREposition

With CAREposition it is possible to perform visually controlled object positioning without radiation

Radiation-free object positioning via graphic display of the central beam and image edges in the LIH image on the image display

When the table is moved, the current positions of the central beam and image edges are superimposed on the LIH image by a graphic overlay

#### CAREwatch

A measurement chamber is integrated into the collimator housing for acquisition of dose area product and reference air kerma / reference air kerma rate

Displayed on the image system display

Different displays can be configured for fluoroscopy and for fluoro pause:

During fluoro: reference air kerma rate

During fluoro pause: accumulated reference air kerma or dose area product or percentage of a configurable dose limit value (total of fluoroscopy and acquisition)

# Artis one

## Floor-mounted system for uncompromised imaging

### CARE

#### **CAREmonitor**

CAREmonitor shows the accumulated peak skin entrance dose according to the current projection in the form of a fill indicator on the live monitor. Any change to the C-arm, table, SID, zoom, or collimator prompts the system to automatically update the calculation.

#### **CAREguard**

CAREguard provides an effective way to control skin dose. Three reference air kerma threshold values can be defined. If the accumulated reference air kerma exceeds a configured threshold, a warning sound is given and a pop-up displays on the system.

#### **CAREreport**

CAREreport is a DICOM structured dose report; it contains all patient demographics, procedure, and dose information. Using commercially available programs or in-house software, this information can be filtered for further processing, such as dose analysis.

#### **Low-dose acquisition**

Low-dose acquisition provides excellent image quality with a dose reduction of up to 67 % in comparison to normal acquisition protocols. One acquisition pedal of the footswitch can be configured as a low-dose acquisition pedal.

#### **Automatic exposure control**

Automatic X-ray control operating five fully independent, self-adjusting, and angulation-driven parameters for optimal dose calculation based on fluoroscopic values

# Artis one

## Floor-mounted system for uncompromised imaging

### Imaging system

High-resolution digital imaging system with outstanding image quality due to real-time image processing

### CLEAR

CLEAR optimizes image quality through real-time processing of the image data.

CLEARcontrol: The new histogram analysis provides a more homogeneous image impression by harmonizing over- and underexposed areas of the image. This is done fully automatically, thus eliminating any further manual user corrections through windowing.

CLEARview: Dose-dependent filtering of the image data efficiently suppresses image noise, enabling clear, sharp images, even for low-dose acquisitions.

CLEARvessel: Every pixel is analyzed in real time, and vessel edges are shown in high contrast without adding noise to the image.

CLEARmotion: Detection of fine structures and effective compensation of motion artifacts.

Fine moving structures, such as small vessels and guidewires, are detected in the image and motion artifacts are suppressed efficiently. The visibility of small moving vessels and guidewires is improved significantly during fluoroscopy.

CLEARchoice: Allows customization of image quality settings according to user preference.

### Image storage capacity

25,000 images in 1k/12-bit matrix

50,000 images in 1k/12-bit matrix\*

100,000 images in 1k/12-bit matrix\*

### General functions

Fast, direct access to all series, single images and reference images, store monitor images, in both the examination room and the control room

#### Changing window values

#### Zooming/Panning

Modification on the fly during postprocessing and pre-configurable for each individual acquisition program

#### Annotation

For inserting predefined or free text and drawing lines, arrows and circles

#### Distance and angle measurement

#### Text functions

Preconfigured image labeling using text modules or free annotation, comment line for image, patient positioning annotation



# Artis one

## Floor-mounted system for uncompromised imaging

### Advanced applications

#### Quantification

##### **QVA – Vascular analysis for vessel diameters of 0.5 mm – 50 mm** (not for coronary analysis)\*<sup>1)</sup>

Measurement program integrated into the imaging system for exact and reproducible vascular analysis

Automatic contour recognition

Stenosis quantification

Automatic and manual determination of reference diameter

Automatic and manual calibration methods

Diameter measurement

##### **LVA – Left ventricular analysis**\*<sup>1)/2)</sup>

Scientific measurement program integrated in the imaging system for evaluating the functional efficiency of the left ventricle

Automatic and manual contour recognition

Calculation of the ejection fraction, volumes and indices (area-length and Simpson methods)

Wall motion (centerline, radial and regional methods)

Automatic and manual calibration

Diameter measurement

##### **QCA – Scientific coronary analysis for vessel diameters of 0.5 mm – 7 mm**\*<sup>1)</sup>

Scientific cardiological vessel analysis with stenosis quantification:

Scientific measurement program integrated into the imaging system for clinically validated, objective, exact and reproducible evaluation of coronary arteries

Automatic contour recognition

Stenosis measurement with geometrical and densitometric calculations

Automatic and manual determination of reference diameter

Automatic and manual calibration methods

Diameter measurement

##### **IZ3D**\*<sup>1)</sup>

IZ3D offers automated detection and 3D analysis of single and bifurcated coronary arteries from 2D angiographic images.

Out-of-plane magnification and foreshortening errors are minimized by calculating true geometric shape in 3D space from two 2D X-ray projections.

In stent planning mode, a virtual stent can be specified. This virtual stent is then displayed in the 3D image and corresponding markers are overlaid onto live fluoro and acquisition.

##### **Angle/length measurement with automatic calibration**

##### **DICOM network connection and syngo user interface**

Remark: Quantitative Coronary Analysis (QCA) is based on the gold standard in coronary analysis:

CAAS II (Cardiovascular Angiography Analysis System Mark II) by Pie Medical, Netherlands. The CAAS II algorithms were developed at Erasmus University in Rotterdam. They have been clinically validated and are internationally recognized for scientific purposes (multi-center studies.)

# Artis one

## Floor-mounted system for uncompromised imaging

### CLEARstent\*

Uses an image-quality optimized algorithm to improve the visibility of the deployed stent during cardiac interventions

Optionally, contrast dye can be given. CLEARstent then calculates a scene alternating between the contrast-filled lumen and the stent-enhanced image (CLEARstent Dynamic).

Resulting images and scenes can be archived in PACS and reviewed on any DICOM viewer

### CLEARstent Live\*

CLEARstent Live algorithm stabilizes moving stent images without noticeable delay

CLEARstent Live supports frame rates up to 15 fps

Works even when the device is moved within the coronary vessel or contrast agent is injected, allowing precise stent positioning relative to previously implanted stents and/or vessel anatomy

Processed images are displayed side by side with original scene on assist segment

The stabilized CLEAR stent Live scenes are automatically saved to scene directory allowing for review of resulting DICOM images on any DICOM viewer

### 3D Imaging\*

Integrated 3D with parallel patient processing

New interface for easier and faster tableside 3D manipulation and viewing

**QuickZoom:** Focus and zoom 3D volumes in the control room or at tableside with just one click

### LA Segmentation\*

One-click segmentation of anatomical structures in 3D image data, especially the left atrium. Structures can be exported to EP mapping systems (e.g. CARTO).

### EP Option

Dedicated measure to improve signal noise in the EP lab. The kit is mounted to the tube and will minimize electromagnetic interference to the other EP recording and EP mapping systems in direct proximity to the system.

# Artis one

## Floor-mounted system for uncompromised imaging

### System specifications

#### Stand

The Artis one angio system is specifically designed to meet the increasing demands of highly flexible imaging for interventional radiology and interventional cardiology.

#### C-arm system

Highly flexible and quick positioning

Single joystick for patient-angle oriented C-arm and detector movements

Integrated computerized collision protection

Programmable positioning	up to 8 system positions, additional 70 user-definable positions and 3 direct positions
Isocenter-to-floor distance	107 cm (42.13")
Focus-to-isocenter distance	75 cm (29.53")
Patient coverage (free floating tabletop, minimum without repositioning)	210 cm (82.68") with table movement
C-arm depth	92.5 cm (36.4")
Stand rotation	motorized programmable positioning
C-arm oblique projections <sup>1)</sup>	± 130° LAO/RAO and + 55°/– 45° CRAN/CAUD at 0° head-end C-arm position; + 81° LAO to – 59° RAO and + 48° CRAN to – 53° CAUD at 35° left-side C-arm position
Angulation speed	variable rotation up to 25°/s with LAO/RAO and 25°/s with CRAN/CAUD; variable rotation, 3D up to 60°/s
Variable focal spot-to-detector distance	approx. 90 cm – 120 cm (35.4" - 47.24"), speed up to 9 cm/s (3.54")
Longitudinal C-arm movement	motorized up to 250 mm/s (9.84"/s)
Transversal C-arm movement	motorized up to 150 mm/s (5.9"/s)

#### Maximum positioning flexibility

Stand rotation for free positioning of system and table relative to one another, for the following positions, in addition to others:

Patient access from the left side

Right-side C-arm positioning 30° relative to the longitudinal axis of the patient and double oblique projections of 55°/69° LAO/RAO and + 45°/– 52° CRAN/CAUD

Stand rotation motorized from ± 160°

Orthogonal system control oriented to the longitudinal axis of the patient

#### Automap stand\*

Automatic stand positioning depending on the reference image selected

#### Automap image\*

Automatic reference image selection depending on the current stand positioning

# Artis one

## Floor-mounted system for uncompromised imaging

### Patient tables

Depending on the diagnostic and therapeutic focus, the patient table enables user-specific application

#### Artis one table

Floor-mounted patient table for angiographic examinations and interventions

Large unobstructed cantilevered tabletop and wide range of rotation enables access to patient from all sides and easy transfer and positioning

Telescoping column with motorized height adjustment

Table control module for operation of all table functions

Table height 75 cm to 110 cm (29.53" to 43.3")

Table width 65 cm (25.6") incl. rails

Table length 284.5 cm (112")

Lift speed 5 cm/s (1.97"/s)

Table rotation  $\pm 120^\circ$  in  $3^\circ$  increments

Manual longitudinal travel  $\pm 62.5$  cm (24.6"); [125 cm (49.2")]

Manual transverse travel  $\pm 17.5$  cm (6.9")

Maximum unobstructed overhang 210 cm (82.7")

Maximum table load 410 kg (903.9 lbs.)

– 250 kg (551.16 lbs.) patient weight

– 100 kg (220.46 lbs.) table accessories

– 60 kg (132.28 lbs.) cardiopulmonary resuscitation (CPR)

# Artis one

## Floor-mounted system for uncompromised imaging

### Free-floating tabletop

#### Tabletop/mattress

Narrow form with recess at head end, e.g., for cardiological applications. The tabletop is tapered in the thorax region for the greatest possible freedom of C-arm angulation.

Tabletop	Length: 279.5 cm (110.04"); width: 48.0 cm (18.9")
Width at foot end	55.5 cm (21.85")
Width at thorax	48.0 cm (18.9")
Max. patient weight	250 kg (551.2 lbs.)
Al equivalent tabletop	≤ 1.4 mm (0.06") at 100 kV, HVL 3.6 mm (0.15") Al
Al equivalent mattress thin	< 0.6 mm (0.02") at 100 kV, HVL 3.6 mm (0.15") Al
Mattress thickness thin	4 cm (1.57")

# Artis one

## Floor-mounted system for uncompromised imaging

### DCS Artis one

Fix mounted display ceiling suspension system DCS Artis one for one 30" display or one 30" and one additional display (up to 30"\*) enables height adjustment and swivel capabilities. Enhanced positioning range and flexibility by double pivot cantilever.

Travel radius	200 cm (78.74")
Vertical lift (height adjustment)	52.2 cm (20.6")
Length of cantilever	100 cm and 100 cm (39.37")
Rotation range between extension and stand	330°, settings every 30°
Rotation range between extension and cantilever	330°, settings every 30°
Rotation range of displays	330°, settings every 30°

### Integrated Data Display

All examination-relevant data of the system and table geometric data, system messages, and dose data with the CAREWATCH option are displayed on the data area on the examination or control room display of the imaging system

# Artis one

## Floor-mounted system for uncompromised imaging

### Display control room

#### 21" Color Display

21" TFT high-contrast color display for flicker-free, distortion-free image display for X-ray diagnostics as well as interventional therapeutic procedures

Light weight, high luminance and contrast values

Ambient light sensor for optimum adaption to the room brightness

Diagonal screen measurement	21" (54 cm)
Image display	1600 x 1200
Pixel size	0.270 mm x 0.270 mm
Calibrated luminance	270 cd/m <sup>2</sup>
Contrast ratio	1500 : 1
Horizontal viewing area	178°
Power consumption	< 48 VA (W)
Power save mode	< 0.5 VA (W)

### Display examination room

#### Display 30"

Selectable display layouts	flexible free
Resolution	2560 x 1600
Pixel size	0.256 mm x 0.256 mm
Contrast ratio	typical 1500 : 1
Maximum luminance	1050 cd/m <sup>2</sup>
Calibrated luminance	400 cd/m <sup>2</sup>
Display area (W x H)	655.36 x 409.6 mm
Dimensions without stand (W x H x D)	730.9 x 484.9 x 80.6 mm
Weight (net)	approx 14.3 kg (31.5 lbs.)
Power consumption	57 W



# Artis one

## Floor-mounted system for uncompromised imaging

### StreamLink\*

Streaming of the examination room display content via IP network. Supports up to two streaming destinations for remote display on a Windows PC, e.g. in a conference or lecture room.

StreamLink also supports recording of the examination room display for later download.

Requires a separate Windows PC <sup>1)</sup> (Windows 7 or 8), with the VLC player plugin for Internet Explorer 10 or higher / Firefox 22 or higher.

### Second 21" color display\*

Second 21" display, examination room, mounted next to the standard Artis one 30" display\*

Third party display, connectable up to 24"

### Artis one Panoramic Display\*

A second 30" display, identical to the standard 30" Artis one display, mounted in the same display ceiling suspension (DCS).

Displays up to 9 additional external image sources such as mapping systems, recording systems, IVUS or PACS in its own configurable layouts.

Individual image sources can have resolutions of up to 1920 x 1200 pixels, and can be shown unscaled (1:1), scaled up or scaled down, depending on the selected layout. Special scaling algorithms ensure optimal image quality also for physiological curves.

User interactions like switching of layouts and image sources and operation of the external user interfaces, are done from the control room via a cloned Panoramic Display. Operation of the Artis one system is unchanged, as it uses its standard displays.

### Operation in the examination room

System operation via modular control elements at the patient table for controlling C-arm movement, patient table and collimators.

Artis one Head-up display and tactile system operation through tableside control for operating the imaging system including post-processing as well as selecting organ programs

Ergonomically designed footswitch for releasing fluoroscopy, acquisition, and table brake, as well as an additional configurable function.

### Operation in the control room

Siemens Healthcare universal *syngo* interface using keyboard and mouse for activating system functions such as post-processing, archiving, and configuring fluoro and acquisition programs

Multi-functional hand switch\* for acquisition control, switching acquisition frame rates and/or step movements (option for PERISTEPPING and/or PERIVISION)

\* Option; <sup>1)</sup> Streaming to *syngo* X Workplace or Artis systems is not supported

# Artis one

## Floor-mounted system for uncompromised imaging

### Connectivity

#### DICOM Functions

##### DICOM Send

Sends images and series to DICOM networks or workstations

##### DICOM StC (Storage Commitment)

Receives archiving confirmation from the image archive

##### DICOM Print\*

Prints image material using virtual film sheets via DICOM print laser camera or network laser printer

##### DICOM Query/Retrieve

Searches for images and series in DICOM networks (Query)

Imports images and series from DICOM networks (Retrieve)

##### DICOM Get Worklist\*

Imports patient and procedure data from a DICOM patient management system

##### DICOM MPPS\* (Modality Performed Procedure Step)

Sends dose data as well as patient examination status to a patient data management system

##### Exam protocol can be sent as DICOM image

##### DICOM SR

Stores quantification results and relevant dose data as DICOM Structured Report and sends it to DICOM network

##### Ready Processed Images

Configurable transfer mode to store and archive overlays and post-processing results in the image pixels

#### Networking

Ethernet interface, full-duplex, gigabit transfer rate

# Artis one

## Floor-mounted system for uncompromised imaging

### Data export

DVD drive for automatic digital image storage (incl. DICOM viewer) on a DVD or CD-R for offline data exchange in DICOM format, JPEG, Bitmap or AVI

---

Scene recorder\* for archiving fluoroscopies and acquisitions on a DVD

---

USB interface to copy images on a memory stick or on an external hard disk

---

### Integration of the Siemens Recording System

#### Sensis Interface\*

Interface to Sensis hemodynamic and electrophysiological recording system for automatic acquisition or transfer of patient demographic data and system parameters (dose report)

---

Sensis can be operated tableside via Artis one Head-up display

---

# Artis one

## Floor-mounted system for uncompromised imaging

### Injectors

For more information and additional injectors, please refer to the accessories catalog

### Standard and optional accessories

Please refer to separate catalog

### Remote Service\*

Preparation for Siemens Remote Service (SRS):

Allows hardware and software remote diagnosis

Allows remote system configuration, e.g., adding a DICOM node

Early warning system to help ensure system operation (Guardian)

### Security Package

#### **syngo Security Package\***

SW option for Artis with expanded security features such as user management and audit trail function

# Artis one

## Floor-mounted system for uncompromised imaging

### Room preparation

#### Emergency power supply\*

##### Emergency power supply\* for the imaging system

Bridging of the imaging system power supply (50/60 Hz) until line voltage is back. In case of power failures of more than 90 seconds the imaging system will be shut down automatically.

Nominal power 2 kVA

##### Emergency power supply\* for all system, table movements and imaging system

Emergency power supply for uninterrupted power supply for all system and table movements, as well as imaging system and monitors for a period of at least 10 min. during a primary power failure.

On-site emergency power supply system is a legal requirement in accordance with IEC 60601-2-43

Nominal power 15 kVA

Line voltage 400 V / 440 V or 480 V; an adaptation to 440/480 V is required.

##### Emergency power supply\* for the entire system incl. emergency fluoro

Emergency power supply for the entire system incl. emergency fluoro for a period of at least 10 minutes during a primary power failure. Uninterrupted power supply for all system and table movements, as well as imaging system and monitors.

Approx. 25 seconds after switching on and restarting the generator, you will be able to work with continuous fluoroscopy in emergency operation mode.

Nominal power 40 kVA

Line voltage 400 V / 440 V or 480 V; an adaptation to 440/480 V is required.

#### Internal line resistance for generator A100 Plus <sup>1)</sup>

Support of hospital emergency power generator (diesel generator) switch to UPS mode (continuous Fluoro, same as with system EPS when emergency power active)

$U_N/P$	80 kW	100 kW
380 V	$\leq 135 \text{ mOhm}$	$\leq 110 \text{ mOhm}$
400 V – 460 V	$\leq 135 \text{ mOhm}$	$\leq 135 \text{ mOhm}$
480 V	$\leq 135 \text{ mOhm}$	$\leq 125 \text{ mOhm}$

<sup>1)</sup> To achieve the full generator power, the measured internal line resistance should not exceed the following values.  
Resistance values in Ohm at  $U_N \pm 10\%$

# Artis one

## Floor-mounted system for uncompromised imaging

### Installation data

Line voltage connection, 3-phase current, TN-S

<b>Generator</b>	<b>POLYDOROS A100 Plus</b>
Nominal voltage <sup>1)</sup> (3 phase)	380 V, 400 V, 420 V, 440 V, 460 V ± 10%, 50/60 Hz ± 1 Hz; 480 V, 60 Hz
Fuse	internal 50 A, external 63 A
Power consumption	0.2 kVA system off 1.1 kVA system off (when tube cooling is running) 1.6 kVA system in stand-by 2.5 kVA system in stand-by (when tube cooling is running) 8.0 kVA for fluoroscopy 160 kVA for acquisition

### System control cabinet

Nominal voltage <sup>1)</sup> (3 phase)	380 V, 400 V, 420 V, 440 V, 460 V ± 10%, 50/60 Hz ± 1 Hz; 480 V, 60 Hz
Fuse	internal 8 A, external 25 A slow-blow fuse
Power consumption	max. 7.2 kVA

<sup>1)</sup> Max. allowable nominal voltage between phases (L1, L2, L3) and PE 300 V

### Weight

Examination room	Stand	approx. 570 kg	(1257 lbs.)
	Display ceiling suspension (DCS) 30"	approx. 100 kg	(221 lbs.)
	Display ceiling suspension (DCS) 30" + 24"	approx. 117 kg	(258 lbs.)
	Display ceiling suspension (DCS) 30" + 30"	approx. 125 kg	(276 lbs.)
	Patient table	255 kg	(563 lbs.)
Control room	Injector wall connection box	approx. 5 kg	(11 lbs.)
	Imaging system	approx. 75 kg	(166 lbs.)
	30" display	approx. 17.5 kg	(38.6 lbs.)
	Display controller	approx. 15 kg	(33 lbs.)
	UPS for image system (option)	13 kg	(29 lbs.)
	Miscellaneous	10 kg	(22 lbs.)
Equipment room	Generator	300 kg	(662 lbs.)
	Cooling system (X-ray tube)	< 39 kg	(< 86 lbs.)
	System control cabinet	270 kg	(595 lbs.)
	Cable cabinet (option)	120 kg	(265 lbs.)

# Artis one

## Floor-mounted system for uncompromised imaging

### Ambient conditions (operation)

General room specifications for:

Examination room, control room and equipment room	Temperature range:	+ 15 °C to + 30 °C (recommended temp. 22 °C)
	Relative humidity:	20 – 75% below dew point
	Max. temperature gradient:	5 °C/h
	Barometric pressure:	70 kPa – 106 kPa

Imaging system	For climatic conditions, see the general room specifications	
	Air flow:	156 m <sup>3</sup> /h
	Max. noise level:	53 dB (A)

Generator	For climatic conditions, see the general room specifications	
	Air flow:	160 m <sup>3</sup> /h
	Max. noise level:	55 dB (A)

System control cabinet	For climatic conditions, see the general room specifications	
	Air flow:	295 m <sup>3</sup> /h
	Max. noise level:	55 dB (A)

Cooling unit (for MEGALIX X-ray tube)	Cooling air:	+ 15 °C to + 30 °C (frost-free room)
	Air flow:	950 m <sup>3</sup> /h
	Max. noise level:	55 dB (A) at 50 Hz; 59 dB (A) at 60 Hz

Stand	Mechanical impact:	max. 10 g/16 ms
	Vibrations:	max. 0.1 g/10-200 Hz

Operation altitude	Less than or equal to 3000 meters (10,000 ft)
--------------------	---

Overvoltage category	II
----------------------	----

Pollution degree	2
------------------	---

Oxygen enriched environment	n/a
-----------------------------	-----

### Mobile Artis one installations\*

The Artis one mobile kit option enables the installation of an Artis one system in a mobile cath-lab (trailer).

Mobile Artis one systems are typically operated by carriers who set up mobile cath-labs temporarily, e.g. while a hospital's own cath-lab is under construction.


\* The Artis one mobile kit option is not commercially available in the United States or some other countries. Due to regulatory reasons their future availability cannot be guaranteed. Please contact your local Siemens organization for further details.



# Artis one


Floor-mounted system for uncompromised imaging

System view (mm)



# Artis one

## Floor-mounted system for uncompromised imaging



On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.

For floor-mounted systems 10848600

#### **Siemens Healthcare Headquarters**

Siemens Healthcare GmbH

Henkestr. 127

91052 Erlangen

Germany

Phone +49 9131 84-0

[siemens.com/healthcare](https://www.siemens.com/healthcare)

#### **Legal Manufacturer**

Siemens Shenzhen Magnetic Resonance Ltd.

Siemens MRI Center,

Gaoxin C. Ave., 2nd,

Hi-Tech Industrial Park

518057 Shenzhen,

China