# Bioshot

Disposable automatic guillotine needle  
for biopsy on soft tissues

Code: OBS

Safety catch

Lower releaser  
for automatic shot

Upper releaser  
for automatic shot


Cursor to  
charge the  
stylet

Cursor to  
charge the  
cannula

Choice between  
two different  
sampling depths

Sliding stopper  
Internal echo-marker  
Depth marks

CE  
0426


biomedical

# Bioshot

Code: OBS

Disposable automatic guillotine needle for biopsy on soft tissues

## Instruction of use


1  
Withdraw the cannula cursor until it stops (hooking click)  
Withdraw the stylet cursor until it stops (hooking click)

2  
When it is charged, both indicators are colored.

3  
Choose the sample's length to take by moving the proper selector on the longer or the shorter position.

4  
Put the safety catch on the LOAD position.

5  
Push the upper shooting button\*

6  
Or alternatively push the lower shooting button\*


7  
To extract the needle: withdraw the cannula cursor until it stops (hooking click) in order to uncover the stylet for the frustule taking

(\*) The device can be used also as semi-automatic needle pushing the shooting button softly:

A  
The first press starts the stylet,


B  
The second one starts the cannula.

If the indicators are not colored the device is uncharged.

To order:

CODE	SIZE	PACKAGE	C. CODE
OBS 1410	14Gx10 cm	20	DARK GREEN
OBS 1415	14Gx15 cm	20	DARK GREEN
OBS 1420	14Gx20 cm	20	DARK GREEN
OBS 1610	16Gx10 cm	20	WHITE
OBS 1615	16Gx15 cm	20	WHITE
OBS 1620	16Gx20 cm	20	WHITE
OBS 1625	16Gx25 cm	20	WHITE
OBS 1810	18Gx10 cm	20	PINK
OBS 1815	18Gx15 cm	20	PINK
OBS 1820	18Gx20 cm	20	PINK
OBS 1825	18Gx25 cm	20	PINK

On request Biomedical is available to provide needles with different lengths and 20G diameter.

## biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy  
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629  
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

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**BIOMEDICAL S.r.l.**

Via G.B. Lulli, 43 - 50144 FIRENZE (FI) - ITALIA

nella sede operativa di  
*in the operational headquarter located in*

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## Conformity declaration

The undersigned company BIOMEDICAL S.r.l., declares, on its own responsibility, that the following disposable sterile medical devices produced and supplied by the company itself:

- are in accordance with the essential requirements of Annex I of the European Directive 93/42/EEC and subsequent amendments and additions.
- are classified as surgical invasive medical device for transient use (class IIa) in accordance to the rule 6 , annex IX of Directive 93/42/EEC and subsequent amendments and additions.
- are manufactured in strict accordance with Italian sanitary laws and with European Directive 93/42/EEC and subsequent amendments and additions - Annex V and Annex VII, UNI EN ISO 9001:2015, UNI ISO 13485:2016, using medical/pharmaceutical grade materials.

This conformity declaration is related to the certificate number 178-02-00-DM, issued by the N.B. 0426 ITALCERT, and to the codes listed below.

Firenze, 28/04/2022

BIOMEDICAL S.r.l.  
Legal representative

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**biomedical** srl

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	αα	DESCRIPTION
0Blxxyy/zzz-αα	da 14G a 18G	da 05 a 99 mm	da 005 a 300 mm		0BI = disposable needle for sternal or iliac aspiration
0BIPxxyy/zzz-αα	da 14G a 18G	da 05 a 99 mm	da 005 a 300 mm		0BIP = pre-assembled disposable needle for iliac or sternal aspiration
0BIExxyy/zzz-αα	da 14G a 18G	da 05 a 99 mm	da 005 a 300 mm		0BIE = medullar disposable needle for removal and perfusion
0BILxxyy/zzz-αα	da 14G a 25G	da 05 a 99 mm	da 005 a 300 mm		0BIL = disposable needle for iliac aspiration with ergonomic handle
0BILExxyy/zzz-αα	da 14G a 25G	da 05 a 99 mm	da 005 a 300 mm		0BILE = disposable needle for iliac aspiration with ergonomic handle for medullar removal
0BLxxyy-αα	da 07G a 18G	da 05 a 25 cm	---		0BL = disposable needle for bone-medullar biopsy
0BLExxyy-αα	da 07G a 18G	da 05 a 25 cm	---		0BLE = disposable needle for medullar removal with lateral holes
0BWxxyy-αα	da 07G a 18G	da 05 a 25 cm	---		0BW = disposable needle for bone-medullar biopsy (patent SOMATEX)
0BLTxxyy-αα	da 07G a 18G	da 05 a 25 cm	---		0BLT = disposable needle for bone-medullar biopsy (Trap System-type)
0BLTxxyy-02	da 07G a 18G	da 05 a 25 cm	---	02 = punta a diamante	0BLT-02 = disposable needle for bone-medullar biopsy (Trap System-type, diamond shaped tip)
0BMxxyy-αα	da 11G a 22G	da 05 a 50 cm	---		0BM = disposable needle for soft tissue biopsy with manual action

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	αα	DESCRIPTION
0BDxxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0BD = disposable needle for soft tissue biopsy with semiautomatic action
0BDSxxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0BDS = disposable guillotine needle for soft tissue biopsy with semiautomatic action (shuttle handle), without removable cannula
0BDSExxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0BDSE = disposable guillotine needle for soft tissue biopsy with semiautomatic action (shuttle handle), without removable cannula, echogenic version (PTFE covered)
0BDPxxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0BDP = disposable guillotine needle for coaxial biopsy
0BDPSxxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0BDPS = disposable guillotine needle for coaxial soft tissue biopsy with semiautomatic action (shuttle handle)
0BExxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0BE = disposable guillotine needle for automatic gun
0BEExxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0BEE = disposable guillotine needle for automatic gun with PTFE covering
0BECxxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0BEC = disposable guillotine needle for automatic gun with coaxial introducer
0BExxyy-01	da 14G a 22G	da 05 a 50 cm	---	01 = versione Promag 2.5	0BE-01 = disposable guillotine needle for automatic gun- Promag 2.5 version
0BECxxyy-01	da 14G a 22G	da 05 a 50 cm	---	01 = versione Promag 2.5	0BEC-01 = disposable guillotine needle for automatic gun- Promag 2.5 version with coaxial introducer
0PGxxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0PG = disposable guillotine needle compatible with Bard Magnum automatic gun
0PGxxyy-01	da 14G a 22G	da 05 a 50 cm	---	01 = punta a diamante	0PG = disposable guillotine needle compatible with Bard Magnum automatic gun, diamond shaped tip
0PGExxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0PGE = disposable guillotine needle compatible with Bard Magnum automatic gun, echogenic version (PTFE covered)

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	αα	DESCRIPTION
0PGCxxyy-αα	da 14G a 22G	da 05 a 50 cm	---		0PGC = disposable guillotine needle with coaxial introducer compatible with Bard Magnum automatic gun
0BSxxyy-αα	da 14G a 20G	da 10 a 25 cm	---		0BS = automatic disposable gun for soft tissue biopsy (BIOSHOT)
0BSExxyy-αα	da 14G a 20G	da 10 a 25 cm	---		0BSE = automatic disposable gun for soft tissue biopsy (BIOSHOT), echogenic version (PTFE covered)
0ITxxyy-αα	da 11G a 23G	da 02 a 50 cm	---		0IT = disposable introducer coaxial needle for introduction of cannula, biopsy devices, hook wires, catheters or other devices, without mandrel and smoothed tip cannula
0ITExxyy-αα	da 11G a 23G	da 02 a 50 cm	---		0ITE = disposable introducer coaxial needle for introduction of cannula, biopsy devices, hook wires, catheters or other devices, without mandrel and smoothed tip cannula, echogenic version (PTFE covered)
0BTxxyy-αα	da 11G a 23G	da 02 a 50 cm	---		0BT = disposable introducer coaxial needle for introduction of cannula, biopsy devices, hook wires, catheters or other devices, with Trocar mandrel and biopsy tip cannula
0RFxxyyyβzzδ	Da 13 a 19G	Da 50 a 250 mm	---		0RF = disposable introducer needle for radiofrequency probes, PTFE covered
0BAxxyy-αα	da 14G a 22G	da 05 a 35 cm	---		0BA = disposable set of hepatic biopsy by aspiration as for Menghini
0BCxxyy-αα	da 14G a 25G	da 05 a 30 cm	---		0BC = disposable needle of histological and cytological biopsy by aspiration as for Menghini modified
0BCAxxyy-αα	da 14G a 25G	da 05 a 30 cm	---		0BCA = automatic disposable needle of histological and cytological biopsy by aspiration as for Menghini modified
0BCACxxyy-αα	da 14G a 25G	da 05 a 30 cm	---		0BCAC = automatic disposable needle of histological and cytological biopsy by aspiration as for Menghini modified with coaxial introducer

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	αα	DESCRIPTION
OMEMxxyy-αα	da 16G a 26G	da 05 a 30 cm	---		OMEM = Chiba disposable needle histological and cytological aspiration or insertion, Menchetti version
OMEFxxyy-αα	da 16G a 26G	da 05 a 30 cm	---		OMEF = Disposable needle for cytological and histological biopsy, Franseen style
OMEWxxyy-αα	da 16G a 26G	da 05 a 30 cm	---		OMEW = Wescott disposable needle for cytological and histological biopsy
OICxxyy-αα	da 18G a 23G	da 05 a 35 cm	---		OIC = disposable needle for alcoholization with closed trocar point
OILxxyy-αα	da 18G a 23G	da 05 a 35 cm	---		OIL = open pointed disposable needle for alcoholization
OMExxxyy-αα	da 14G a 29G	da 05 a 30 cm	---		OME = Chiba disposable needle for cytological aspiration
OMExxxyy-01	da 14G a 29G	da 05 a 30 cm	---	01 = punta triplice affilatura	OME-01 = Chiba disposable needle for cytological aspiration Triple sharpening tip
OMExxxyy-02	da 14G a 29G	da 05 a 30 cm	---	02 = punta Quincke	OME-02 = Chiba disposable needle for cytological aspiration Quincke tip
OMExxxyy-03	da 14G a 29G	da 05 a 30 cm	---	03 = punta Ultra	OME-03 = Chiba disposable needle for cytological aspiration Ultra tip
OMExxxyy-04	da 14G a 29G	da 05 a 30 cm	---	04 = punta triplice affilatura e tappo Quick Release	OME-04 = Chiba disposable needle for cytological aspiration Triple sharpening tip and Quick Release plug
OMExxxyy-05	da 14G a 29G	da 05 a 30 cm	---	05 = punta Quincke e tappo Quick Release	OME-05 = Chiba disposable needle for cytological aspiration Quincke tip and Quick Release plug
OMExxxyy-06	da 14G a 29G	da 05 a 30 cm	---	06 = punta Ultra e tappo Quick Release	OME-06 = Chiba disposable needle for cytological aspiration with Ultra tip and Quick Release plug

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	αα	DESCRIPTION
OMExxyy-07	da 14G a 29G	da 05 a 30 cm	---	07 = punta triplice affilatura, tappo Quick Release e prolunga	OME-07 = Chiba disposable needle for cytological aspiration Triple sharpening tip Quick Release plug and extension
OMExxyy-08	da 14G a 29G	da 05 a 30 cm	---	08 = punta Quincke, tappo Quick Release e prolunga	OME-08 = Chiba disposable needle for cytological aspiration Quincke tip Quick Release plug and extension
OMExxyy-09	da 14G a 29G	da 05 a 30 cm	---	09 = punta Ultra, tappo Quick Release e prolunga	OME-09 = Chiba disposable needle for cytological aspiration Ultra tip, Quick Release plug and extension
OMEpxxyy-αα	da 16G a 26G	da 050 a 300 mm	---		OMEp = disposable needle with peripheral nerve block under echography – without mandrel – with connecting tube
OMEpxxyy-01	da 16G a 26G	da 050 a 300 mm	---	01 = affilatura a 30°	OMEp-01 = disposable needle with peripheral nerve block under echography – without mandrel – with connecting tube (30° angle tip)
OMECxxyy-αα	da 16G a 26G	da 05 a 30 cm	---		OMEC = Chiba disposable needle for cytological aspiration, echogenic cannula and tube
OMEBxxyy-αα	da 14G a 26G	da 05 a 30 cm	---		OMEB = disposable needle for brachytherapy
OIDxx/Dyy-αα	da 18G a 30G	da 20 a 40 mm	---		OIDxx/D = set for galactography with straight needle
OIDxx/Cyy-αα	da 18G a 30G	da 20 a 40 mm	---		OIDxx/C = set for galactography with curved needle
OIDKxx/Cyy-αα	da 18G a 30G	da 20 a 40 mm	---		OIDKxx/C = set for galactography with curved needle, 6 dilators and 1ml syringe
OIDKxx/Dyy-αα	da 18G a 30G	da 20 a 40 mm	---		OIDKxx/D = set for galactography with straight needle, 6 dilators and 1ml syringe
OIMxxyy-αα	da 14G a 25G	da 05 a 35 cm	---		OIM = disposable needle for mammary localization

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	α	DESCRIPTION
0IMXxyyy-αα	da 14G a 25G	da 05 a 35 cm	---		0IMX = disposable needle for mammary localization with double hook
0PNAxyyy-αα	da 14G a 25G	da 05 a 35 cm	---		0PNA = disposable needle for mammary localization, not repositionable; BCA-type handle
0PRxyyy-αα	da 14G a 25G	da 05 a 35 cm	---		0PR = disposable needle for mammary localization; normal handle
0AETxyyyy-αα	da 18G a 30G	da 001 a 200 mm	---		0AET = disposable needle with triple tip for ENS electro stimulation
0AETxyyyy-01	da 18G a 30G	da 001 a 200 mm	---	01 = affilatura a 30°	0AET-01 = disposable needle with triple tip for ENS electro stimulation (30° angle tip)
0AETxyyyy-02	da 18G a 30G	da 001 a 200 mm	---	02 = affilatura a 45°	0AET-02 = disposable needle with triple tip for ENS electro stimulation (45° angle tip)
0AEQxyyyy-αα	da 18G a 30G	da 001 a 200 mm	---		0AEQ = disposable needle with Quincke tip for ENS electro stimulation
0AEQxyyyy-01	da 18G a 30G	da 001 a 200 mm	---	01 = affilatura a 30°	0AEQ-01 = disposable needle with Quincke tip for ENS electro stimulation (30° angle tip)
0AEQxyyyy-02	da 18G a 30G	da 001 a 200 mm	---	02 = affilatura a 45°	0AEQ-02 = disposable needle with Quincke tip for ENS electro stimulation (45° angle tip)
0AECxyyyy-αα	da 18G a 30G	da 001 a 200 mm	---		0AEC = disposable needle with Chiba tip for ENS electro stimulation
0AECxyyyy-01	da 18G a 30G	da 001 a 200 mm	---	01 = affilatura a 30°	0AEC-01 = disposable needle with Chiba tip for ENS electro stimulation (30° angle tip)
0AECxyyyy-02	da 18G a 30G	da 001 a 200 mm	---	02 = affilatura a 45°	0AEC-02 = disposable needle with Chiba tip for ENS electro stimulation (45° angle tip)

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	α	DESCRIPTION
0AEYxxyyy-αα	da 18G a 30G	da 001 a 200 mm	---		0AEY = disposable needle with Thuoy tip for ENS electro stimulation
0AEYxxyyy-01	da 18G a 30G	da 001 a 200 mm	---	01 = affilatura a 30°	0AEY-01 = disposable needle with Thuoy tip for ENS electro stimulation (30° angle tip)
0AEYxxyyy-02	da 18G a 30G	da 001 a 200 mm	---	02 = affilatura a 45°	0AEY-02 = disposable needle with Thuoy tip for ENS electro stimulation (45° angle tip)
0AEPxxyyy-αα	da 18G a 30G	da 001 a 200 mm	---		0AEP = disposable needle with Pencil point tip and lateral hole for ENS electro stimulation
0AEPxxyyy-01	da 18G a 30G	da 001 a 200 mm	---	01 = affilatura a 30°	0AEP-01 = disposable needle with Pencil point tip for ENS electro stimulation (30° angle tip)
0AEPxxyyy-02	da 18G a 30G	da 001 a 200 mm	---	02 = affilatura a 45°	0AEP-02 = disposable needle with Pencil point tip for ENS electro stimulation (45° angle tip)
0AV01xxyyA	da 12G a 20G	da 05 a 25 cm	A = punta aperta		0AV01xxyyA = set for thoracentesis/paracentesis with tap and Verres needle, with open tip mandrel
0AV01xxyyB	da 12G a 20G	da 05 a 25 cm	B = 2 fori		0AV01xxyyB = set for thoracentesis/paracentesis with tap and Verres needle, with two holes mandrel
0AV01xxyyC	da 12G a 20G	da 05 a 25 cm	C = 3 cave		0AV01xxyyC = set for thoracentesis/paracentesis with tap and Verres needle, with 3 pits mandrel
0AV0103yy	03	da 05 a 25 cm	---		0AV0103yy = set for thoracentesis/paracentesis with tap and 3 needles
0AV02xxyyA	da 12G a 20G	da 05 a 25 cm	A = punta aperta		0AV02xxyyA = set for thoracentesis/paracentesis with valve and Verres needle with open tip mandrel
0AV02xxyyB	da 12G a 20G	da 05 a 25 cm	B = 2 fori		0AV02xxyyB = set for thoracentesis/paracentesis with valve and Verres needle, with 2 holes mandrel

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	αα	DESCRIPTION
0AV02xxyyC	da 12G a 20G	da 05 a 25 cm	C = 3 cave		0AV02xxyyC = set for thoracentesis/paracentesis with valve and Verres needle, with 3 pits mandrel
0AV0203yy	03	da 05 a 25 cm	---		0AV0203yy = set for thoracentesis/paracentesis with valve and 3 needles
0AV03xxyyy	da 08 a 21G	da 030 a 300 mm	---		0AV03xxyyy = set for thoracentesis/paracentesis for removal of ascites and drainage needle
0AV0102	---	---	---		0AV0102 = set for drainage with 2 funnel needles, tap and tube 4.1 x 3 mm
0AVAxxyy	da 12G a 20G	da 05 a 25 cm	A = punta aperta		0AVAxxyy = Verres needle with normal slim joint, open point, 16° sharpen, without tap
0AVBxxyy	da 12G a 20G	da 05 a 25 cm	B = 2 fori		0AVBxxyy = Verres needle with normal slim joint, 2 holes, 16° sharpen, without tap
0AVCxxyy	da 12G a 20G	da 05 a 25 cm	C = 3 cave		0AVCxxyy = Verres needle with normal slim joint, 3 pits, 16° sharpen, without tap
0AVDxxyy	da 12G a 20G	da 05 a 25 cm	D = 1 foro		0AVDxxyy = Verres needle with normal slim joint, 1 hole, 16° sharpen, without tap
0AVAQxxyy	da 12G a 20G	da 05 a 25 cm	AQ = punta aperta, affilatura a 30° e rubinetto		0AVAQxxyy = Verres needle with normal slim joint, open tip, 30° sharpen, and tap
0AVBQxxyy	da 12G a 20G	da 05 a 25 cm	BQ = 2 fori, affilatura a 30° e rubinetto		0AVBQxxyy = Verres needle with normal slim joint, 2 holes, 30° sharpen, and tap
0AVCQxxyy	da 12G a 20G	da 05 a 25 cm	CQ = 3 cave, affilatura a 30° e rubinetto		0AVCQxxyy = Verres needle with normal slim joint, 3 pits, 30° sharpen, and tap

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	αα	DESCRIPTION
0AVDQxxyy	da 12G a 20G	da 05 a 25 cm	DQ = 1 foro, affilatura a 30° e rubinetto		0AVDQxxyy = Verres needle with normal slim joint, 1 hole, 30° sharpen and tap
0AVAxxyy-02	da 12G a 20G	da 05 a 25 cm	A = punta aperta	02 = copertura PTFE	0AVAxxyy-02 = Verres needle with normal slim joint, open tip, 16° sharpen without tap and PTFE cover
0AVBxxyy-02	da 12G a 20G	da 05 a 25 cm	B = 2 fori	02 = copertura PTFE	0AVBxxyy-02 = Verres needle with normal slim joint, 2 holes, 16° sharpen without tap and PTFE cover
0AVCxxyy-02	da 12G a 20G	da 05 a 25 cm	C = 3 cave	02 = copertura PTFE	0AVCxxyy-02 = Verres needle with normal slim joint, 3 pits, 16° sharpen without tap and PTFE cover
0AVDxxyy-02	da 12G a 20G	da 05 a 25 cm	D = 1 foro	02 = copertura PTFE	0AVDxxyy-02 = Verres needle with normal slim joint, 1 hole, 16° sharpen without tap and PTFE cover
0AVAQxxyy-02	da 12G a 20G	da 05 a 25 cm	AQ = punta aperta, affilatura a 30° e rubinetto	02 = copertura PTFE	0AVAQxxyy-02 = Verres needle with normal slim joint, open tip, 30° sharpen, tap and PTFE cover
0AVBQxxyy-02	da 12G a 20G	da 05 a 25 cm	BQ = 2 fori, affilatura a 30° e rubinetto	02 = copertura PTFE	0AVBQxxyy-02 = Verres needle with normal slim joint, 2 holes, 30° sharpen, tap and PTFE cover
0AVCQxxyy-02	da 12G a 20G	da 05 a 25 cm	CQ = 3 cave, affilatura a 30° e rubinetto	02 = copertura PTFE	0AVCQxxyy-02 = Verres needle with normal slim joint, 3 pits, 30° sharpen, tap and PTFE cover
0AVDQxxyy-02	da 12G a 20G	da 05 a 25 cm	DQ = 1 foro, affilatura a 30° e rubinetto	02 = copertura PTFE	0AVDQxxyy-02 = Verres needle with normal slim joint, 1 hole, 30° sharpen, tap and PTFE cover
0AVAxxyy-01	da 12G a 20G	da 05 a 25 cm	A = punta aperta	01 = raccordo slim corto	0AVAxxyy-01 = Verres needle with normal slim joint, open tip, 16° sharpen, without tap

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	α	DESCRIPTION
0AVBxxyy-01	da 12G a 20G	da 05 a 25 cm	B = 2 fori	01 = raccordo slim corto	0AVBxxyy-01 = Verres needle with normal slim joint, 2 holes, 16° sharpen without tap
0AVCxxyy-01	da 12G a 20G	da 05 a 25 cm	C = 3 cave	01 = raccordo slim corto	0AVCxxyy-01 = Verres needle with normal slim joint, 3 pits, 16° sharpen without tap
0AVDxxyy-01	da 12G a 20G	da 05 a 25 cm	D = 1 foro	01 = raccordo slim corto	0AVDxxyy-01 = Verres needle with normal slim joint, 1 hole, 16° sharpen without tap
0AVAQxxyy-01	da 12G a 20G	da 05 a 25 cm	AQ = punta aperta, affilatur a a 30° e rubinett o	01 = raccordo slim corto	0AVAQxxyy-01 = Verres needle with normal slim joint, open tip, 30° sharpen and tap
0AVBQxxyy-01	da 12G a 20G	da 05 a 25 cm	BQ = 2 fori, affilatur a a 30° e rubinett o	01 = raccordo slim corto	0AVBQxxyy-01 = Verres needle with normal slim joint, open tip, 30° sharpen, tap and PTFE cover
0AVCQxxyy-01	da 12G a 20G	da 05 a 25 cm	CQ = 3 cave, affilatur a a 30° e rubinett o	01 = raccordo slim corto	0AVCQxxyy-01 = Verres needle with normal slim joint, 3 pits, 30° sharpen and tap
0AVDQxxyy-01	da 12G a 20G	da 05 a 25 cm	DQ = 1 foro, affilatur a a 30° e rubinett o	01 = raccordo slim corto	0AVDQxxyy-01 = Verres needle with normal slim joint, 1 hole, 30° sharpen and tap
0AVAxxyy-03	da 12G a 20G	da 05 a 25 cm	A = punta aperta	03 = raccordo slim corto e copertura PTFE	0AVAxxyy-03 = Verres needle with normal slim joint, open tip, 16° sharpen without tap and PTFE cover
0AVBxxyy-03	da 12G a 20G	da 05 a 25 cm	B = 2 fori	03 = raccordo slim corto e copertura PTFE	0AVBxxyy-03 = Verres needle with normal slim joint, 2 holes, 16° sharpen without tap and PTFE cover

CODE	RANGE XX	RANGE YY(Y)	RANGE Z(ZZ)	αα	DESCRIPTION
0AVCxxyy-03	da 12G a 20G	da 05 a 25 cm	C = 3 cave	03 = raccordo slim corto e copertura PTFE	0AVCxxyy-03 = Verres needle with normal slim joint, 3 pits, 16° sharpen without tap and PTFE cover
0AVDxxyy-03	da 12G a 20G	da 05 a 25 cm	D = 1 foro	03 = raccordo slim corto e copertura PTFE	0AVDxxyy-03 = Verres needle with normal slim joint, 1 hole, 16° sharpen without tap and PTFE cover
0AVAQxxyy-03	da 12G a 20G	da 05 a 25 cm	AQ = punta aperta, affilatura a 30° e rubinetto	03 = raccordo slim corto e copertura PTFE	0AVAQxxyy-03 = Verres needle with normal slim joint, open tip, 30° sharpen without tap and PTFE cover
0AVBQxxyy-03	da 12G a 20G	da 05 a 25 cm	BQ = 2 fori, affilatura a 30° e rubinetto	03 = raccordo slim corto e copertura PTFE	0AVBQxxyy-03 = Verres needle with normal slim joint, 2 holes, 30° sharpen without tap and PTFE cover
0AVCQxxyy-03	da 12G a 20G	da 05 a 25 cm	CQ = 3 cave, affilatura a 30° e rubinetto	03 = raccordo slim corto e copertura PTFE	0AVCQxxyy-03 = Verres needle with normal slim joint, 3 pits, 30° sharpen without tap and PTFE cover
0AVDQxxyy-03	da 12G a 20G	da 05 a 25 cm	DQ = 1 foro, affilatura a 30° e rubinetto	03 = raccordo slim corto e copertura PTFE	0AVDQxxyy-03 = Verres needle with short slim joint, 1 hole, 30° sharpen without tap and PTFE cover