
Intrapur® Inline





B | BRAUN

| | |
|-----------|------------------------|
| en | Instructions for use |
| de | Gebrauchsanweisung |
| bg | Инструкция за употреба |
| cs | Návod k použití |
| da | Brugsanvisning |
| el | Οδηγίες χρήσης |
| es | Instrucciones de uso |
| et | Kasutusjuhised |
| fi | Käyttöohje |
| fr | Mode d'emploi |
| hr | Upute za uporabu |
| hu | Használati utasítás |
| id | Petunjuk penggunaan |
| it | Istruzioni per l'uso |
| lt | Naudojimo instrukcija |

| | |
|-----------|---------------------------|
| lv | Lietošanas instrukcija |
| nl | Gebruikersinformatie |
| no | Bruksanvisning |
| pl | Instrukcja użytkowania |
| pt | Instruções de utilização |
| ro | Instrucțiuni de utilizare |
| ru | Инструкция по применению |
| sl | Navodila za uporabo |
| sk | Návod na použitie |
| sr | Uputstvo za upotrebu |
| sv | Bruksanvisning |
| tr | Kullanım Kılavuzu |
| vi | Hướng dẫn sử dụng |
| zh | 使用说明 |

Intrapur® Inline

en Instructions for use

Device description

IV administration set can be used with devices which are in compliance with EN ISO 80369-7 / ISO 80369-7; ISO 80369-20/ EN ISO 80369-20.

IV administration set is in accordance with EN ISO 8536-4/ ISO 8536-4.

Materials used

Acrylstyrene polybutadiene, Supor®, polyester, polytetrafluoroethylene, polyethylene, polypropylene, polyamide, polystyrene, impact-resistant polystyrene, acrylnitrile-butadiene-styrenecopolymer, polyurethane.

Infusion set with 0.2 µm inline filter

Sterility

Sterilized EO. Please refer to the information on primary packaging.

Do not use if package is damaged.

Visually inspect the product packaging to ensure that the sterile barrier system is intact.

Intended purpose

I.V. administration sets for gravity infusion.

The inline filter of Intrapur® Inline is low protein binding and the membrane of 0.2 µm offers an effective filtration area 10cm² made out of Polyethersulfone which helps to prevent particulate contamination by retaining particles > 0.2 µm and reduces microbiological contamination by retaining bacteria and fungi.

Do not re-use.

Indication

Infusion sets are used for infusion therapy e.g.: Fluid volume maintenance, fluid volume replacement, infusion of medications (e.g. antibiotics, chemotherapeutic agents) with bottles and bags according to the SPC (Summary of product characteristics) of the drugs/solutions.

Patient population

The device can be used for all patients for which infusion therapy is prescribed.

No gender or age-related limitations. Intrapur® Inline can be used for adults, pediatric and neonates.

Intended User

Intrapur® Inline should only be used by healthcare worker who have been adequately trained within the educational background in this technique.

Upon medical assessment and adequate instruction, patients and/or caregivers may also be allowed to take over defined handling steps according to national guidelines. Healthcare workers are responsible for ensuring that patients and/or caregivers are instructed in the correct use of Intrapur® Inline and that the instruction of caregiver and/or patient is documented by the healthcare worker, according to the local documentation guidelines.

Contraindications

IV administration set shall not be used in case of proven incompatibility of the drug/IV solution to one of the used materials. SPCs of the drugs/solutions for contraindications should be considered.

No administration of blood products or blood components.

Intrapur® Inline with 0.2 µm pore size must not be used for the filtration of lipids, suspensions or emulsions and colloidal solutions. The filter should not be used for sterile filtration of solutions with visible contamination by bacteria or pyrogens.

Residual Risks/Side Effects

General risk occurring during infusion therapy including air embolism, local- to systemic infection, particulate contamination, medication error and kinking that may result in decrease or stop of infusion flow. Leakage of cytotoxic or biological fluids can be an additional hazard to patients, visitors or staff.

Precautions

- Visual check of the IV-set for damages before use.
- Do not re-attach the protective cap of the patient connector.
- Do not use if protective caps are missing or loose.
- Mind the general guidelines for parallel infusions
- Please also consider the corresponding SPCs of the drugs/ solutions regarding possible incompatibilities of drugs/ solutions.
- Monitor all connections for tightness over the whole infusion procedure.
- Allow the filter to hang straight down and initiate fluid flow. The filter will automatically prime if left in a vertical position. Do not invert to prime the filter.
- IV filters help to achieve required standards of hygiene, but it is not a replacement for them. Use aseptic techniques in accordance with national guidelines and / or hospital protocols.
- Avoid that the air vent opening of the IV filter gets in contact with e.g. disinfectants or soap from the outside.

Warnings

- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Administration of light sensitive drugs requires UV-protected IV administration sets.
- Always check, if incompatibilities with other drugs or with materials is stated in the SPC (Summary of product characteristics) of the drug.
- Always check, if the use of a 0.2 µm filter is stated in the SPC (Summary of product characteristics) of the drug.
- Do not re-sterilize. Do not use after expiration date.
- Do not use in connection with an infusion pump.

Duration of use

Duration of use is dependent on the intended therapy according to the SPC of the drug or solution. In general, change of administration sets should occur according national guidelines (e.g. CDC guideline) and/or hospital protocols. Intrapur® Inline should be changed after 24h at the latest.

Disposal

Disposal acc. to local guidelines and/or clinical protocols.

Use integrated spike protection at roller clamp for disposal.

Storage and handling conditions

Storage conditions:

- Keep away from sunlight.
- Keep dry.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>

Operating Instructions

☞ Illustrations on page 3.

For Intrapur® Inline re-spike according to hospital protocols and / or guidelines is allowed during a single procedure.

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B. Braun Melsungen AG,