Cerus
Endovascular

NEQSTENT

Coil Assisted Flow Diverter

CE 0297

Instructions For Use

English.....	1
Deutsch (German).....	3
Français (French).....	6
Español (Spanish).....	8
Magyar (Hungarian).....	11
Nederlands (Dutch).....	13
Dansk (Danish).....	16
Italiano (Italian).....	18
Suomi (Finnish).....	21
Norsk (Norwegian).....	23
Svenska (Swedish).....	25
Türk (Turkish).....	27
Português (Portuguese).....	29
Ελληνικά (Greek).....	31
Polski (Polish).....	33


Cerus Endovascular, Ltd.
John Eccles House
Oxford Science Park
Oxford OX4 4GP
United Kingdom
www.cerusendo.com

Customer Service

Healthlink Europe BV
De Tweeling 20-22
5215 MC 's-Hertogenbosch
The Netherlands
+31.13.5479300

Technical Support

Cerus Endovascular, Inc.
47757 Fremont Blvd.
Fremont, CA USA 94538
1 510 651-4000

Emergo Europe B.V.
Prinsessegracht 20
2514AP The Hague
The Netherlands

INSTRUCTIONS FOR USE

Neqstent

Device Description

The Neqstent comprises an implant with pre-attached Detachable Pusher Wire (DPW), and an Introducer. The DPW functions as a guidewire for implant delivery. All devices are provided sterile and nonpyrogenic and are for single patient use only.

The implant is packaged unconstrained in a protective dispenser hoop with an Introducer preloaded onto the DPW shaft just proximal to the implant.

The implant is a self-expanding concave shaped device comprised of a double layer nickel-titanium and platinum wire mesh and a platinum marker. The pre-attached DPW has a composite stainless steel and polymer construction, and a platinum marker. The DPW has a white fluoro-saver mark on its proximal shaft that provides a visual indication of when the implant is approaching but not exiting the distal end of the microcatheter (MC).

The Introducer is a single lumen polymer tube that is used to constrain the implant during introduction into the MC hub. The implant body and proximal marker can be visualized under fluoroscopy.

The implant is delivered to the aneurysm under fluoroscopic guidance using standard endovascular techniques and a commercially available 0.027" or 0.021" inner diameter MC. The implant is electrolytically detached from the DPW using a commercially available detachable coil power supply.

The implant is available in 7 mm, 9 mm, 11 mm and 14 mm diameters to accommodate the range of aneurysm sizes described below.

Figure 1. Packaging hoop.



Figure 2. Device distal end (unconstrained, constrained in Introducer).

Figure 3. Neqstent System.

Indications for Use

The Neqstent is indicated for use in conjunction with conventional embolic coils for endovascular embolization of saccular intracranial aneurysms.

Contraindications

- Allergy to platinum, nickel or titanium

Device Compatibility

The following devices (not supplied) are required for use with the Neqstent:

- Neurovascular Microcatheter (MC):

Neqstent	Any compatible neurovascular microcatheter with an 0.027" ID and up to 160 cm length.
Neqstent 021	Any compatible neurovascular microcatheter with an 0.021" ID and up to 160 cm length.

- Detachable Coil Power Supply: Any compatible electrolytic detachment system
- Neurovascular Microcatheter appropriate to embolization coil delivery through the Neqstent with distal O.D $\leq 0.025"$

Recommended Accessory Products

In addition to the devices listed in Device Compatibility section, the following accessory products are required for the procedure:

- Continuous flush setup, including two rotating hemostasis valves (RHVs), heparinized saline, one 3-way stopcock, and one 1-way stopcock
- Contrast media
- Syringes for contrast media and heparinized saline
- Femoral sheath
- Guide Catheter (minimum 0.071" (1.80 mm) inner diameter is recommended)
- Microguidewire compatible with MC
- Sterile 20 or 22-gauge uncoated stainless steel hypodermic needle to provide electrical ground during implant detachment
- Alcohol-dampened gauze (to clean DPW proximal end before connecting power supply detachment cable)

Recommended Anticoagulation/Antiplatelet Regimen

- Procedural Regimen: Placement of the Neqstent implant in an endovascular procedure requires the use of continuous heparinization. The use of antiplatelet agents during the procedure should be considered on an individual basis.

INSTRUCTIONS FOR USE

Neqstent

- Post-Procedural Regimen: The anti-platelet and anti-coagulation therapy should be administered in accordance with standard medical practice.

Adverse Events

Potential adverse events associated with use of Neqstent, some of which could be fatal or cause severe neurologic deficits, include:

- Aneurysm rupture causing intracranial hemorrhage
- Injury to parent artery causing thrombosis or hemorrhage
- Distal embolization of particles or blood clot causing stroke
- Parent artery vasospasm
- Parent artery dissection
- Aneurysm recanalization
- Infection
- Device migration causing incomplete occlusion, hemorrhage or ischemic stroke

Cerebral angiography carries its own risks such as allergic reaction to contrast media, radiation exposure, and groin femoral artery bleeding/infection

Warnings

- The Neqstent should only be used by physicians who have received appropriate training in neurointerventional techniques and the use of the Neqstent.
- Treatment should not be attempted if the patient's anatomy or physiology is considered unsuitable for endovascular treatment.
- Do not advance or retract the device against unusual resistance without careful assessment of its cause using direct fluoroscopic visualization. Moving device against resistance may result in patient complications, or device malfunctions such as premature implant detachment.
- If unusual resistance is encountered when retracting the implant into the MC, withdraw the MC until the resistance subsides, or remove the MC, implant and DPW from the patient as a single unit.
- Do not rotate the DPW during use. Rotating the DPW may result in premature detachment, which could lead to undesirable implant deployment position or migration.
- Verify that MC is not under stress before implant detachment by gently removing any tension in the MC. Stored forces in the MC could cause the tip to move upon detachment, which could lead to suboptimal implant position and/or aneurysm rupture.

Precautions

- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may compromise the structural integrity of the device or increase the risk of contamination or infection leading to device failure and/or cross-infection and potential patient injury, illness or death.
- Do not use kinked or damaged devices.
- Do not use opened or damaged packages.
- Do not use devices after "Use By" date.
- Use device only in conjunction with fluoroscopic guidance.
- To maintain optimal performance and reduce the risk of thromboembolic events, maintain continuous infusion of heparinized saline solution through the femoral artery sheath, guide catheter, and MC.

Magnetic Resonance Imaging (MRI) Safety Information

Non-clinical testing demonstrated that the Neqstent implant is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, only
- Maximum spatial gradient magnetic field of 2,000-gauss/cm (20 Tesla/m) (extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the First Level Controlled Operating Mode of operation for the MR system

Under the scan conditions defined above, the implant is expected to produce a maximum temperature rise of 2.3°C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant extends approximately 10 mm from the implant when imaged using a gradient echo pulse sequence and a 3 Tesla MR system.


Recommended Procedure Steps

Note: Perform steps below under fluoroscopic guidance as applicable.

1. Continuous Flush Setup (Figure 4)

- Attach first RHV (C1) to the guide catheter (B) hub. Attach second RHV (C2) to first RHV (C1) sidearm. Connect a three-way stopcock (H) to the second RHV (C2) sidearm. Attach a continuous flush line (J) and syringe (I) to the stopcock (H).
- Attach a third RHV (C3) to the Neqstent Delivery MC (D) hub. Attach a continuous flush line (J) to the third RHV (C3) sidearm.
- Insert the Neqstent Delivery MC (D) into the first RHV (C1) connected to the guide catheter (B) hub.
- Attach a fourth RHV (C4) to the coiling MC (F) hub. Attach a continuous flush line (J) to the fourth RHV (C4) sidearm.
- Insert the coiling MC (F) into the second RHV (C2).
- Adjust continuous flush pressure to maintain appropriate infusion of solution during procedure.
- Check all fittings to ensure no air is introduced into system

Figure 4. Flush setup.

A	Femoral Sheath
B	Guide Catheter
C1, C2, C3, C4	RHV
D	Neqstent Delivery Microcatheter
E	Neqstent Detachable Pusher Wire
F	Coiling Microcatheter
G	Coil Detachable Pusher Wire
H	3-Way Stopcock
I	Syringe
J	Line to Flush Solution
K	Double RHV (optional)


2. Implant Size Reference (Figure 5)

- Select appropriate implant size based on Table 1.

Table 1. Implant size reference vs. aneurysm neck dimension

REF (Catalog Number) – Diameter	Aneurysm Neck (mm)
NQS407-XX – 7 mm	3.0 – 5.0
NQS21407-XX – 7 mm	
NQS409-XX – 9 mm	4.0 – 6.0
NQS21409-xx – 9 mm	
NQS411-XX – 11 mm	5.0 – 8.0
NQS21411-XX, 11 mm	
NQS414-XX – 14 mm	7.0 – 10.0

Figure 5: Implant Size Reference

3. Microcatheter Preparation and Placement

- Insert Microcatheter (MC) into sheath and navigate MC over microguidewire through guide catheter to desired vascular location.
- After the MC is in position, remove the microguidewire. Maintain heparinized flush through MC using RHV sidearm.

4. Device Preparation and Placement into Microcatheter

- Detach DPW proximal end from package card.
- Advance DPW into and through package hoop until implant emerges from other end of hoop.
- Grasp DPW distal end and gently pull it out of hoop until entire Introducer is visible.
- Immerse implant and Introducer distal end in saline. Do not tap or apply pressure to implant against the bottom of the saline bath.
- Hold the Introducer at the distal end (close to the implant) and slowly withdraw implant into Introducer.
- Insert Introducer into MC RHV until Introducer distal end is fully seated in MC hub.
- Gently tighten RHV and wait for saline drops to appear at Introducer proximal end. To purge air from space between implant and Introducer, allow continuous flush solution to flow retrograde through the Introducer and out of the Introducer proximal end.
- Advance device through hub until the DPW is completely inside the Introducer.
- Remove Introducer from the wire.
- Maintain continuous flush through RHV.

5. Implant Delivery and Deployment

- Perform all remaining steps under fluoroscopy.
- Confirm MC tip is located in the desired position.
- Carefully advance implant into aneurysm while continually verifying sizing and position.
- If implant placement is not satisfactory, slowly retract into MC by pulling the DPW proximally.
 - The device may be retracted into the microcatheter and redeployed up to 3 times as needed to reposition or remove the device.
- Prior to detachment, verify implant stability:
 - Holding the DPW, keeping the device in its deployed position, pull the MC back approximately 2 cm proximal to the implant;
 - Gently push the DPW to observe a slight bowing of the DPW confirming device stability;
 - The implant should not move into the aneurysm.
- If the implant size is inappropriate, remove and replace with an appropriately sized implant.
 - To facilitate withdrawal for redeployment, ensure the Introducer is securely seated in the MC hub prior to withdrawal (to reload implant into the Introducer for redeployment or to remove it from patient).
- Once satisfactory Neqstent positioning is achieved, advance a coiling microcatheter through the Neqstent mesh, and deliver embolic coils into the aneurysm using standard endovascular techniques.

6. Implant Detachment

- When the final embolic coil has been delivered through the Neqstent into the aneurysm, position MC tip marker proximal to distal end of the DPW to ensure detachment zone has exited MC tip.
- Verify that no tension is present in the MC distal shaft. Tension forces stored in the MC could cause the MC tip or the implant to move upon detachment.

INSTRUCTIONS FOR USE**Neqstent**

- Verify that the RHV is firmly locked around the DPW before attaching the Power Supply to ensure that the implant does not move during the connection process.
- Ensure that the DPW is straight between the RHV and the Power Supply. Straightening this section of the DPW minimizes potential forces on the implant that could affect its position.
- While maintaining DPW and MC position, attach Power Supply and Cable Set as follows:
 - Connect the 2-wire cable plug into the receptacle on the power supply or follow applicable instructions for your power supply and cables.
 - Insert a sterile needle into the skin, and clip the black cable end to the needle.
 - Clip the red cable end to the uncoated proximal section of the DPW.
- Attach Detachment Cable as per Power Supply IFU.
- Actuate Power Supply to initiate implant electrolytic detachment process.
- When Power Supply signals detachment, verify implant detachment by slowly pulling back on the DPW while monitoring the fluoroscopic image. Detachment is indicated by movement of the DPW marker without corresponding movement of the implant marker.
- If the implant is not detached, check and adjust the flush system, flush the system to clear any contrast that may exist around the detachment zone, check and adjust the grounding setup, realign the DPW with the MC, and tighten the RHV. Repeat the detachment process after confirming the implant is in the required position.
- If the implant will not detach, perform the following:
 - Prior to withdrawing the device, advance MC over the implant marker band.
 - Gently pull DPW to retrieve implant into the MC.
- Withdraw MC into guide catheter.
- If the implant was successfully detached, verify implant position and determine need for further actions.

Symbols

	Caution		Lot number
	Catalog number		Magnetic Resonance Conditional, see Instructions For Use
	Contents of package		Manufacturer
	Do not re-sterilize		Nonpyrogenic
	Do not re-use		Sterilized using ethylene oxide
	Do not use if package is damaged		Use by
	Prescription Use Only		Consult Instructions For Use
	Keep dry		Keep away from sunlight

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

Cerus Endovascular warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of Cerus directly affect the device and the results obtained from its use.

The obligation of Cerus under this warranty is limited to the repair or replacement of this device and Cerus shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Cerus neither assumes, nor authorizes any other person to assume for it, any other or additional liability, or responsibility in connection with this device. Cerus assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such device.

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