VIBRANT SOUNDBRIDGE

Vibroplasty Couplers

Incus-Symphonix-Coupler

Incus-LP-Coupler

Incus-SP-Coupler

RW-Soft-Coupler

English | Lietuvių kalba | Latviešu valoda | Eesti keel | Română

English	Instructions for use.....	1
Lietuvių kalba	Naudojimo instrukcija.....	15
Latviešu valoda	Lietošanas instrukcija.....	29
Eesti keel	Kasutusjuhised.....	43
Română	Instrucțiuni de utilizare.....	57

Instructions for use

Vibroplasty Couplers

Introduction

These instructions for use provide information for the following Vibroplasty Couplers:

- Incus-Symphonix-Coupler
- Incus-LP-Coupler
- Incus-SP-Coupler
- RW-Soft-Coupler

Please read these instructions for use carefully and completely so that you are familiar with the use of the Vibroplasty Couplers. Please do not hesitate to contact your MED-EL representative with any additional questions you may have.

Refer to the instructions for use of the Vibrating Ossicular Prosthesis (VORP) of the Vibrant Soundbridge system.


NOTE:

All illustrations are representative only.


Vibroplasty Couplers – Overview

All specified dimensions are shown in millimeters (mm).


Incus-Symphonix-Coupler

Incus-LP-Coupler

Incus-SP-Coupler

RW-Soft-Coupler - 2x


Figure 1: Vibroplasty Couplers – Overview

Description of the device

The Vibroplasty Couplers are middle ear prostheses intended to be used exclusively with the VORP. They are considered part of the Vibrant Soundbridge system.

The types of Coupler for sensorineural hearing loss are:

- Incus-Symphonix-Coupler – left and right versions
- Incus-LP-Coupler – left and right versions
- Incus-SP-Coupler – one version

The type of Coupler for conductive and mixed hearing loss is:

- RW-Soft-Coupler – one version

NOTE:

Peel-off labels are included in the packaging of the Vibroplasty Couplers. These labels are intended to be placed on the Vibrant Soundbridge Implant Registration Card and in the patient's record.

Intended use

The Vibroplasty Couplers are intended to be used in combination with the Vibrant Soundbridge to facilitate the coupling between the FMT and a vibratory structure of the middle ear. The prosthesis type is chosen on the basis of the ossicular remnants once all primary disease has been removed from the middle ear.

- The Incus-Symphonix-Coupler and Incus-LP-Coupler are intended to be placed on the long process of the incus.
- The Incus-SP-Coupler is intended to be placed on the short process of the incus.
- The RW-Soft-Coupler is intended to be placed on the round window (RW) membrane.

Indication

Refer to the instructions for use for the Vibrating Ossicular Prosthesis (VORP).

Contraindication

Refer to the instructions for use for the Vibrating Ossicular Prosthesis (VORP).

The devices shall not be used:

- in the case of a fixed stapes footplate.

The Incus-LP-Coupler and Incus-SP-Coupler shall not be used:

- in the case of a fragile incus.

General information about use

The Vibroplasty Couplers are middle ear prostheses to be used exclusively with the VORP to treat conductive hearing loss, mixed hearing loss and sensorineural hearing loss.

Attaching the FMT to the Coupler

For use, the FMT is first attached to the Coupler. The Coupler is delivered with a holding frame and a retainer which holds the Coupler in place. Attach the FMT to the Coupler while the Coupler is still secured by the holding frame and the retainer. This ensures a safe and correct connection.

For the Incus-Symphonix-Coupler, the Incus-LP-Coupler and the Incus-SP-Coupler, place the FMT on the cage and push it down using surgical tweezers, a needle or a similar tool. Ensure that the conductor link is properly placed into the groove on the holding frame. This can be achieved by pushing the conductor link down using a finger (see Figure 2). Placing the conductor link into the groove ensures proper positioning of the FMT onto the Coupler and the correct orientation of the conductor link. This minimizes unnecessary bending of the conductor link during implantation.


Figure 2: Placing the FMT into the cage

After attaching the FMT to the Coupler, the retainer must be removed from the holding frame and the Coupler by squeezing the two handles of the retainer together (see Figure 3a) and by tilting the retainer into the direction of the holding frame (see Figure 3b).

NOTE:

Remove the retainer before you remove the Coupler with the FMT.

Use surgical tweezers or similar to remove the Coupler with the FMT from the holding frame. Lift the Coupler with the FMT upwards (see Figure 3c). Do not pull at the conductor link.


Figure 3a: Squeezing the two handles together


Figure 3b: Tilting the retainer

Figure 3c: Removing Coupler and FMT from the holding frame

When using the RW-Soft-Coupler, first remove the retainer above the RW-Soft-Coupler: To remove the retainer, squeeze the two handles of the retainer together (see Figure 4a) and tilt the retainer into the direction of the holding frame (see Figure 4b).


Figure 4a: Squeezing the two handles together


Figure 4b: Tilting the retainer

Figure 4c: RW-Soft-Coupler secured by the holding frame

While the RW-Soft-Coupler is still secured by the holding frame, place the FMT onto the Coupler's adhesive pad and press down slightly (see Figure 5). By using this method, the FMT is easily connected to the RW-Soft-Coupler in the right position. Use a microscope to verify the correct position of the FMT on the Coupler. After attaching the Coupler, press the Coupler and the FMT together firmly to ensure a safe connection. If the Coupler is not attached centrally, it is possible to move it again as the adhesive pad does not harden.

NOTE:

Ensure that the RW-Soft-Coupler does not come into contact with bodily fluids before the FMT is attached.

Should the FMT come into contact with bodily fluids, it shall be cleaned with a lint-free tissue before attaching the RW-Soft-Coupler.


Figure 5: Placing the FMT onto the RW-Soft-Coupler

Once a Vibroplasty Coupler is attached to the FMT, the FMT can be connected to various middle ear structures as explained in the following sections.

Incus-Symphonix-Coupler

The Incus-Symphonix-Coupler (see Figure 6) is used to place the FMT onto the long process of the incus via a posterior tympanotomy. Crimp the clip of the Coupler onto the long process of the incus.

Suitable forceps (e.g. MED-EL Forming Forceps) should be used, to mold the Coupler clip in position. Extreme care should be exercised when using the forceps.

Ensure that the FMT is as close as possible or actually touches the stapes. If necessary, bend the connection rod to bring the FMT closer to the stapes. It is possible to bend the connection rod before positioning the FMT with the Coupler in the middle ear.

A left and a right version of the Incus-Symphonix-Coupler are available.


Figure 6: Incus-Symphonix-Coupler placed in a left ear

Incus-LP-Coupler

The Incus-LP-Coupler (see Figure 7) is also used to place the FMT onto the long process of the incus via a posterior tympanotomy. To place the Incus-LP-Coupler, clip the flexible incus clipping structure of the Coupler onto the long process of the incus. Crimping is not necessary.

Ensure that the FMT is as close as possible or actually touches the stapes. If necessary, bend the connection rod to bring the FMT closer to the stapes. It is possible to bend the connection rod before positioning the FMT with the Coupler in the middle ear.

A left and a right version of the Incus-LP-Coupler are available.


Figure 7: Incus-LP-Coupler placed in a left ear

Incus-SP-Coupler

The Incus-SP-Coupler (see Figure 8) is used to place the FMT onto the short process of the incus via an attico-antrotomy. To place the Incus-SP-Coupler, clamp the flexible structure of the Coupler onto the short process of the incus. The two shorter legs of the Coupler shall be placed on the short process of the incus, the two longer legs shall hold the incus body (see Figure 8). The Coupler has to be firmly fixed to the incus.


Figure 8: Incus-SP-Coupler

RW-Soft-Coupler

The RW-Soft-Coupler (see Figure 9) is placed on the FMT by using an adhesive pad. The RW-Soft-Coupler together with the FMT is then positioned on the round window membrane. Cartilage is placed on the contralateral side of the FMT to maintain the position. Ensure that the Coupler does not come off during surgery. When positioning the Coupler and the FMT in the middle ear, be careful not to stick a needle (or similar) between the Coupler and the FMT.


Figure 9: RW-Soft-Coupler

NOTE:

Please refer to the instructions for use for the VORP for information regarding surgical techniques.

Warnings

The warnings of the VORP apply. Carefully read the instructions for use provided for the VORP.

- The Vibroplasty Couplers are shipped sterile and are for single patient use only. Do not re-use or re-sterilize.
- The Vibroplasty Couplers are very small medical devices and should be handled with appropriate care.
- Do not use if the sterile packaging is damaged or the Use-by date (expiration date) has been exceeded.
- The Vibroplasty Couplers must not be altered and may only be used as intended.

Precautions

The precautions of the VORP apply. Carefully read the instructions for use provided for the VORP.

General precautions

- Persons known to be intolerant of the materials used in the Vibroplasty Couplers shall not receive any Vibroplasty Couplers. Please refer to the section **Technical data** at the end of these instructions for use for materials used in the Vibroplasty Couplers.
- The implantation of the VORP together with any Coupler is a sophisticated surgical intervention and should only be performed by experienced and trained surgeons.

Intraoperative precautions

- Accidental bending of the Coupler during removal from its package must be avoided in order to prevent functional damage.
- Always attach the FMT to the Coupler in the holding frame as this guarantees the correct positioning of the Coupler onto the FMT and of the conductor link.
- The middle ear structures must be treated with utmost care.

Postoperative precautions

- When the patient's audio processor is fitted and he/she perceives distortion or noise it may be an indication of slippage of the prosthesis and/or FMT.
- A loss of benefit over time may be the result of unwanted tissue growth that can limit or prevent the conduction of vibrations to the cochlea.

Possible adverse effects

Complications of middle ear surgery also apply for the Vibroplasty Couplers.

Surgery of the middle ear involves manipulation of the fragile ossicular bones and exposes the inner ear to the risk of surgical trauma. Serious complications may arise either during or after surgery that may result in irreparable damage to otologic structures causing irreversible partial or total loss of hearing. Subsequent surgical procedures may be required to correct these conditions, if possible.

Complications that may occur include: sensorineural deafness due to trauma during surgery; granuloma; post-surgery displacement of the implant due to the development of scar tissue; vertigo; post-surgical translocation of the FMT due to trauma or inferior surgical positioning; extrusion of the implant; and post-surgical infection, including otitis media. It is also possible that the device could cause a significant loss in residual hearing.

The responsible physician should ensure that the patient is fully apprised of the possible adverse effects of the implant procedure.

The above-noted additional events are known to be possible adverse events associated with middle ear surgery. Additionally, implant patients are exposed to the normal risk of surgery and general anaesthesia. Major ear surgery may result in numbness, swelling, or discomfort around the ear, the possibility of facial paresis, disturbance of balance or taste, or neck pain. If these occur, they are usually transient and resolve within a few weeks of surgery.

Interference with other equipment

- **Magnetic Resonance Imaging (MRI):** MRI examination with the VORP 503 is only permissible in 1.5 Tesla (T) closed-bore MRI scanners. The following conditions have to be considered:
 - The MRI scanner has to be limited to "Normal Operating Mode"; "First Level Controlled Operating Mode" has to be avoided.
 - Local transmitting RF coils must not be used in the head and neck region. Local receiver coils are not restricted in use.
 - Before patients enter any MRI room, the audio processor must be removed from the head. Audible interference can occur during the scan. Patients shall be advised to indicate any possible discomfort that may arise and to request that the MRI be discontinued if needed. After the MRI examination, the patient shall put on the audio processor only after leaving the MRI room.
 - During MRI examination straight head orientation is required.
 - An artefact of approximately 14cm around the implant will be present on the images.
 - MRI examination with $>1.5T$ will damage the implant and must be avoided.
 - If an MRI examination is needed prior to the first activation of the Vibrant Soundbridge, safety measures shall be taken into account to prevent wound healing complications due to the possible movement of the implant in the strong MR field.
 - When lower extremities are to be examined, it is recommended, but not required, that the patient's legs are positioned in the scanner first.
- The VORP 502 is MR Unsafe by all means.

Miscellaneous

Storage and shipping

The Vibroplasty Couplers shall be stored in their sterile packaging in a cool and dry place. During shipment the temperature shall remain between -20°C and $+60^{\circ}\text{C}$.

Disposal

After use or once the sterile package is opened or damaged, the Vibroplasty Couplers shall be disposed of in a medical waste bin in accordance with federal, state and local regulations. Packaging should be disposed of in accordance with local legislation.

Technical data

Materials

- Incus-Symphonix-Coupler: Titanium Grade 1 (ASTM F67) and Titanium Grade 5 ELI (ASTM F136)
- Incus-LP-Coupler: Titanium Grade 5 ELI (ASTM F136), Titanium Grade 2 (ASTM F67)
- Incus-SP-Coupler: Titanium Grade 5 ELI (ASTM F136)
- RW-Soft-Coupler: Implant Grade Silicone Rubber and Gel

Dimensions

- Refer to Figure 1

Warranty statement

Please refer to the accompanying Warranty Statement for information on our warranty provisions.

Symbols

CE marking, first applied in 2014

Caution

Consult instructions for use

Do not re-use

Manufacturer

Date of manufacture

Use-by date

Catalogue number

Batch code

Sterilised using ethylene oxide

Temperature limit

Please help us to improve the quality of these instructions for use by making any suggestions. For further information regarding the use of this MED-EL product, or to report any problems, please contact:

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Please refer to the accompanying Contact Sheet for your local office.

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MED-EL Elektromedizinische Geräte GmbH