Estradiol

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Estradiol

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of estradiol (estradiol-17 β , E2) in serum, as an aid in the differential diagnosis of amenorrhea, and monitoring of ovulation induction with and without stimulation in assisted reproductive technology (ART).

Catalog Number: **L2KE22** (200 tests),

L2KE26 (600 tests)

Test Code: **E2** Color: **Dark Pink**

Summary and Explanation

Estradiol is a steroid hormone with a molecular mass of 272.3 daltons which circulates predominantly protein-bound. In addition to estradiol, other natural steroidal estrogens include estrone, estriol and their conjugates. Estrogens are hormones secreted principally by the ovarian follicles and also by the adrenals, corpus luteum, placenta and, in males, by the testes.

Estrogenic hormones are secreted at varying rates during the menstrual cycle throughout the period of ovarian activity. The gonadotropins of the anterior pituitary regulate secretion of the ovarian hormones, estradiol and progesterone; hypothalamic control of pituitary gonadotropin production is in turn regulated by plasma concentrations of the estrogens and progesterone. This complex feedback system results in the cyclic phenomenon of ovulation and menstruation.

During pregnancy, the placenta becomes the main source of estrogens. At the menopause, ovarian secretion of estrogens declines at varying rates.

Principle of the Procedure

IMMULITE 2000 Estradiol is a solid-phase, enzyme-labeled chemiluminescent competitive immunoassay. The solid phase (bead) is coated with polyclonal rabbit anti-estradiol

antibody. The liquid phase consists of alkaline phosphatase (bovine calf intestine) conjugated to estradiol.

The patient sample and the reagent are incubated together with the coated bead for 60 minutes. During this time, estradiol in the sample competes with enzyme conjugated estradiol in the reagent for a limited number of antibody binding sites on the bead. Unbound patient sample and enzyme conjugate are then removed by centrifugal washes. Finally, chemiluminescent substrate is added to the reaction tube containing the bead and the signal is generated in proportion to the bound enzyme.

Incubation Cycles: 1 $\times$ 60 minutes

Time to First Result: 65 minutes

Specimen Collection

Lipemia may interfere with the assay. An ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Estradiol has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 25 μ L serum

Storage: 2 days at 2–8°C, or
2 months at –20°C.¹⁸

Dilution: It is recommended that patient samples yielding results greater than 1200 pg/mL be diluted and reassayed. A clinically relevant result is obtained without dilution but the most accurate value is obtained when tested on the sensitive portion of the curve.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.²⁰⁻²²

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Estradiol Bead Pack (L2E212)

With barcode. 200 beads, coated with polyclonal rabbit anti-estradiol. Stable at 2–8°C until expiration date.

L2KE22: 1 pack **L2KE26:** 3 packs

Estradiol Reagent Wedge (L2E2A2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to estradiol, in buffer with preservative. Stable at 2–8°C until expiration date.

L2KE22: 1 wedge **L2KE26:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Estradiol Adjustors (LE2L, LE2H)

Two vials (Low and High), 2 mL each, of estradiol in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KE22: 1 set **L2KE26:** 2 sets

Before running adjustors, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Estradiol Sample Diluent (L2E2Z)

For the on-board dilution of patient samples. One vial containing 25 mL of concentrated (ready-to-use) processed human serum (with preservative) containing undetectable to low levels of estradiol. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2E2Z: 25 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2E2Z: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

E2TC: Estradiol Control (single level)

Also Required

Distilled or deionized water; test tubes; controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of estradiol.

Siemens Healthcare Diagnostics recommends the use of commercially

available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Expected Values

Based on its relationship to IMMULITE Estradiol (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

Reference ranges were generated using IMMULITE Estradiol in a multi-national study involving women in apparent good health (age: 16–44 years), who volunteered to have blood samples drawn, on a daily basis, throughout one complete ovulatory cycle. (See "Menstrual Cycle Graph.")

		Estradiol, pg/mL		
Ovulatory Cycles	n*	Median	Central	95%
Follicular Phase	54 (708)	42	ND	160
Follicular Phase, Days 2 to 3	54 (108)	31	ND	84
Periovulatory, ± 3 days	54 (378)	133	34	400
Luteal Phase	54 (604)	93	27	246

*Number of subjects (total number of results)

		Estradiol, pmol/L		
Ovulatory Cycles	n*	Median	Central	95%
Follicular Phase	54 (708)	154	ND	587
Follicular Phase, Days 2 to 3	54 (108)	114	ND	308
Periovulatory, ± 3 days	54 (378)	489	124	1468
Luteal Phase	54 (604)	343	101	905

*Number of subjects (total number of results)

Another study performed with IMMULITE Estradiol yielded the following results.

Group	n	Estradiol, pg/mL		
		Mean	Median	90% Range
Adult males	50	30.5	29.7	ND–56
Adult females:				
Untreated Postmenopausal	27	ND	ND	ND–30
Treated Postmenopausal	27	ND	ND	ND–93
Oral Contraceptives	61	35.2	24.5	ND–102

ND: not detectable

Group	n	Estradiol, pmol/L		
		Mean	Median	90% Range
Adult males	50	112	109	ND–206
Adult females:				
Untreated Postmenopausal	27	ND	ND	ND–110
Treated Postmenopausal	27	ND	ND	ND–341
Oral Contraceptives	61	129	90	ND–374

ND: not detectable

Pediatric: Reference intervals for the pediatric population (children and adolescents) were established for the IMMULITE Estradiol assay in accordance with CLSI guideline EP28-A3C.²³ Samples were collected prospectively from apparently healthy pediatric subjects, using predefined inclusion criteria. Reference values were generated for subpopulations based on age and Tanner stage subgroups based on physiological development. The study was designed to establish reference values across genders, and to include approximately equal numbers of males and females within each age or Tanner stage subgroup. The subject's Tanner stage was assessed based on pubic hair and genitalia/breast development. The scale proposed by Neinstein and Kaufman was used for the determination of the Tanner stages.²⁴

The reference intervals and Tanner values are based on the central 90% (5th and 95th percentiles). Where sample sizes were insufficient to calculate the 5th or 95th percentile, the minimum or maximum observed values are presented in the Reference Intervals and Reference Values tables.

IMMULITE 2000/2000 XPi Estradiol Pediatric Reference Intervals

Male

Age (Years)	n	Median	Range
		pg/mL	
2–3	10	< 20.0	< 20.0*→ 28.0†
4–9	57	< 20.0	< 20.0–41.9
10–13	119	< 20.0	< 20.0–34.7
14–21	117	26.9	< 20.0–59.2

Age (Years)	n	Median	Range
		pmol/L	
2–3	10	< 73.4	< 73.4*→ 102.8†
4–9	57	< 73.4	< 73.4–53.9
10–13	119	< 73.4	< 73.4–127.3
14–21	117	98.7	< 73.4–217.4

Female

Age (Years)	n	Median	Range
		pg/mL	
2–3	17	< 20.0	< 20.0*→ 34.5†
4–9	44	20.9	< 20.0–40.7
10–11	55	27.5	< 20.0–176.0
12–21	162	60.0	< 20.0–235.2

Age (Years)	n	Median	Range
		pmol/L	
2–3	17	< 73.4	< 73.4*→ 126.6†
4–9	44	76.7	< 73.4–149.4
10–11	55	100.9	< 73.4–645.9
12–21	162	220.2	< 73.4–863.0

*Value presented is the minimum reportable value observed; insufficient sample size to calculate a 5th percentile limit.

† Value presented is the maximum value observed; insufficient sample size to calculate a 95th percentile limit.

**IMMULITE 2000/2000 XPI Estradiol
Pediatric Reference Values by Tanner
Stage**

Male

Tanner Stage	n	Median	Range
		µg/dL	
1	73	< 20.0	< 20.0–38.0
2	63	< 20.0	< 20.0–30.9
3	63	< 20.0	< 20.0–54.2
4	58	26.4	< 20.0–51.8
5	46	30.0	< 20.0–62.3

Tanner Stage	n	Median	Range
		µmol/L	
1	73	< 73.4	< 73.4–139.5
2	63	< 73.4	< 73.4–113.4
3	63	< 73.4	< 73.4–199.0
4	58	96.7	< 73.4–190.2
5	46	110.1	< 73.4–228.7

Female

Tanner Stage	n	Median	Range
		µg/dL	
1	72	< 20.0	< 20.0–54.2
2	46	34.3	< 20.0–177.2
3	64	47.5	< 20.0–203.5
4	46	56.7	< 20.0–227.5
5	50	89.7	< 20.0–363.5

Tanner Stage	n	Median	Range
		µmol/L	
1	72	< 73.4	< 73.4–199.1
2	46	125.7	< 73.4–650.3
3	64	174.1	< 73.4–746.8
4	46	207.9	< 73.4–834.7
5	50	329.0	< 73.4–1334.0

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Pregnancy Samples: Exercise caution when assaying pregnancy samples, since estriol levels may be high enough to interfere.

Neonatal Samples: The assay has not been validated for use on neonatal samples. Crossreacting steroids, including estriol, circulating at high concentrations during this period may cause spuriously elevated results.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this

assay should always be used in combination with the clinical examination, patient medical history, and other findings.

The drug Fulvestrant (FASLODEX®)* may cause falsely elevated estradiol results in immunoassays. For patients being treated with Fulvestrant, an alternate method that is not expected to show cross reactivity to Fulvestrant, such as Liquid Chromatography-Mass Spectrometry (LC-MS), should be used.

With the advent of new steroid based medications (analogues) with similar chemical structures to estradiol, there is the possibility of cross-reactivity and results inconsistent with the patients clinical history. For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings. If the estradiol results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.

* FASLODEX® is a registered trademark of AstraZeneca

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

pg/mL $\times$ 3.671 $\rightarrow$ pmol/L

Reportable Range: 20 to 2000 pg/mL (73–7342 pmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 15 pg/mL (55 pmol/L)

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Precision Profile: In a series of studies involving 10 kit lots, samples spanning the working range were assayed in 10 or 20

replicates per run. (See "Precision Profile" graph.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three estradiol solutions (800, 1440 and 4020 pg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antiserum was generated with estradiol derivative at 6- position, and is highly specific for estradiol. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 μ L/mL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL may interfere with the assay. (See "Lipemia" table.)

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 33 volunteers into plain, heparinized, and Becton Dickinson SST® Vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of estradiol, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Estradiol procedure.

(Heparin) = 0.98 (Serum) + 11 pg/mL
 $r = 0.996$

(SST) = 0.98 (Plain Tubes) + 15 pg/mL
 $r = 0.994$

Means:
476 pg/mL (Serum)
477 pg/mL (Heparin)
482 pg/mL (SST)

Method Comparison: The assay was compared to IMMULITE Estradiol on 97 samples. (Concentration range: approximately 20 to 1800 pg/mL. See graph.) By linear regression:

(IML 2000) = 0.90 (IML) + 18 pg/mL
 $r = 0.990$

Means:

337 pg/mL (IMMULITE 2000)

354 pg/mL (IMMULITE)

References

- 1) Batzer F. Hormonal evaluation of early pregnancy. *Fertil Steril* 1980;34:1–13.
- 2) Bergquist C, Nilius SJ, Wide L. Human gonadotropin therapy: 1. Serum estradiol and progesterone patterns during conceptual cycles. *Fertil Steril* 1983;39:761–5. 3) Erickson GF. Normal ovarian function. *Clin Obstet Gynecol* 1978;21:31–52. 4) Garcia JE, Jones GS, Wright GL. Prediction of the time of ovulation. *Fertil Steril* 1981;36:308–15. 5) Gautray JP, et al. Clinical investigation of the menstrual cycle: 3. Clinical, endometrial, and endocrine aspects of luteal defect. *Fertil Steril* 1981;35:296–303.
- 6) Goebelsmann U, Mishell DR. The menstrual cycle. In: Mishell DR, Davajan V, editors. *Reproductive endocrinology, infertility and contraception*. Philadelphia: Davis FA, 1979: 67–89. 7) Goldstein D, et al. Correlation between estradiol and progesterone in cycles with luteal phase deficiency. *Fertil Steril* 1982;37:348–54.
- 8) Haning RV, et al. Plasma estradiol is superior to ultrasound and urinary estriol glucuronide as a predictor of ovarian hyperstimulation during induction of ovulation with menotropins. *Fertil Steril* 1983;40:31–6. 9) Judd HL, Korenman SG. Effects of aging on reproductive function in women. In: Korenman SG, editor. *Endocrine aspects of aging*. New York: Elsevier Biomedical, 1982:163–97. 10) Landgren BM, Aedo AR, Diczfalussy E. Hormonal changes associated with ovulation and luteal function. In: Flamigni C, Givens JR, editors. *The gonadotropins: basic science and clinical aspects in females*. London: Academic Press, 1982:187–201. 11) March CM, Goebelsmann U, Nakamura RM, Mishell DR. Roles of estradiol and progesterone in eliciting the midcycle luteinizing hormone and follicle-stimulating hormone surges. *J Clin Endocrinol Metab* 1979;49:507–12. 12) Polan ML, et al. Abnormal ovarian cycles as diagnosed by ultrasound and serum estradiol levels. *Fertil Steril* 1982; 37:342–7. 13) Radwanska E, et al. Plasma progesterone and oestradiol estimation in the diagnosis and treatment of luteal insufficiency in menstruating infertile women. *Acta Eur Fertil* 1976;39–47. 14) Radwanska E, et al. Correlation between preovulatory serum estradiol and

midluteal progesterone levels during induction of ovulation with clomid and HCG. *J Reprod Med* 1980; 24:79–82. 15) Rebar RV, Yen SSC. Endocrine rhythms in gonadotropins and ovarian steroids with reference to reproductive processes. In: Krieger, editor. *Endocrine rhythms*. New York: Raven Press, 1979: 259–98. 16) Robertson RD, et al. Assessment of ovulation by ultrasound and plasma estradiol determinations. *Obstet Gynecol* 1979; 54:686–90. 17) Siiteri PK, Febres F. Ovarian hormone synthesis, circulation, and mechanisms of action. In: DeGroot L, et al, editors. *Endocrinology*, vol 3. New York: Grune & Stratton, 1979:1401–17. 18) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: W.B. Saunders, 1994.

- 19) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 20) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377–82, 387–8.
- 21) Clinical and Laboratory Standards Institute (formerly NCCLS). *Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline - Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
- 22) Federal Occupational Safety and Health Administration, *Bloodborne Pathogens Standard*, 29 CFR 1910.1030. 23) Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory: Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Guideline EP28-A3C. 24) Neinstein LS and Kaufman FR, Chapter 1: Normal Physical Growth and Development in Neinstein L.S. *Adolescent Health Care: A Practical Guide*, 4th ed.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (pg/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	89	8.8	9.9%	14	16%
2	180	14	7.8%	19	11%
3	461	20	4.3%	31	6.7%
4	1309	84	6.4%	102	7.8%
5	1800	89	4.9%	127	7.1%

Linearity (pg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	370	—	—
	4 in 8	204	185	110%
	2 in 8	92	93	99%
	1 in 8	43	46	93%
2	8 in 8	756	—	—
	4 in 8	394	378	104%
	2 in 8	193	189	102%
	1 in 8	97	95	102%
3	8 in 8	1849	—	—
	4 in 8	975	925	105%
	2 in 8	466	462	101%
	1 in 8	220	231	95%

Recovery (pg/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	208	—	—
	A	221	238	93%
	B	254	270	94%
	C	407	399	102%
2	—	273	—	—
	A	272	299	91%
	B	329	331	99%
	C	483	460	105%
3	—	614	—	—
	A	612	623	98%
	B	645	655	98%
	C	819	784	104%


Specificity

Compound ¹	ng/mL Added ²	% Cross-reactivity ³
Androstenedione	100,000	0.0004%
Cortisol	100,000	0.000026%
DHEA	4000	0.0006%
17 α -Estradiol	100	0.026%
17 β -Estradiol-3-glucuronide-17-sulfate	50	0.038%
17 β -Estradiol-3-sulfate	50	0.15%
5-Androstene-3 β -17 β -diol	100,000	0.00035%
	10,000	0.00082%
β -Estradiol-17-propionate	3.6	1.25%
β -Estradiol-3-sulfate-17-glucuronide	50	0.04%
β -Estradiol-17-valerate	3.6	1.14%
Estriol	20	0.81%
Estriol-3-sulfate	10,000	0.00053%
Estriol-3-glucuronide	1000	0.0018%
Ethinyl-Estradiol	20	1.382%
	10	1.829%
Estrone	3.6	0.69%
Estrone- β -D-glucuronide	3.6	0.19%
Estrone-3-sulfate	3.6	ND
<i>d</i> -Equilenin	3.6	3.3%
Equilin	3.6	0.38%
Fulvestrant	20	0.31%*
Norgestrel	100	0.029%
	50	0.126%
Progesterone	100,000	ND
Raloxifene Hydrochloride	125	ND
Tamoxifen Citrate	125	ND
Testosterone	100,000	ND

ND: not detectable⁴



*Initial concentration of Estradiol was 21.6 pg/mL. A 286% change in concentration was observed.⁵

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2018 Siemens Healthcare Diagnostics. All rights reserved.