

Passive Middle Ear Implant

Tympanoplasty Prosthesis

AW58347_2.0

hearLIFE

English	2
Deutsch	5
Français	9
Italiano.....	12
Español	15
Português	18
Nederlands	21
Dansk.....	24
Svenska	27
Suomi.....	30
Norsk	33
Polski	36
Türkçe	40
Magyar	43
Čeština.....	46
Slovenčina.....	49
Hrvatski	52
Slovenščina.....	55
Srpski jezik	58
Български език.....	61
Română.....	65
Русский.....	68
Ελληνικά.....	72
Lietuvių kalba	75
Latviešu valoda.....	78
Eesti keel	81
Українська	84
Қазақша.....	87
한국어	91
Tiếng Việt	94
中文.....	97
中文.....	100
العربية.....	105

Instructions for use

Passive Middle Ear Implant Tympanoplasty Prosthesis

Device description

The MED-EL passive middle ear implant – tympanoplasty prosthesis is an ossicular replacement prosthesis made of titanium. The ossicular chain can be reconstructed to restore the conductive hearing mechanism within the middle ear.

Intended use

The tympanoplasty prosthesis is intended to be used as a partial or total replacement of the ossicles of the human middle ear. This passive implant restores the mechanical sound transmission from the tympanic membrane to the oval window.

The partial prosthesis is intended to be used for replacement of the ossicular chain if at least the stapes and its footplate are still present and functional. The total prosthesis is intended to be used for replacement of the ossicular chain if at least the stapes footplate is functional. The tympanoplasty prosthesis is intended to treat patients of all ages.

Indications

The tympanoplasty prosthesis is indicated in case of congenital or acquired defects of the ossicular chain due to e.g.

- chronic otitis media
- traumatic injury
- malformation
- cholesteatoma

The tympanoplasty prosthesis is indicated to treat patients with inadequate conductive hearing from previous middle ear surgery.

Contraindications

The tympanoplasty prosthesis should not be used if:

- the patient is known to be intolerant to the implant materials
- an acute middle ear infection is present

Caution: The physician must fully assess the potential risks and benefits for the patient and his/her realistic expectations with the device prior to the decision to implant. The physician must exercise medical judgement and consider the patient's complete medical history.

Possible adverse effects

Surgery of the middle ear involves manipulation of the fragile ossicular bones and exposes the inner ear to the risk of surgical trauma. Serious complications may arise either during or after middle ear surgery that may result in irreparable damage to otologic structures causing irreversible partial or total loss of hearing. Subsequent surgical procedures may be required to correct these conditions, if possible.

Complications that may occur include: sensorineural hearing loss due to trauma during surgery; vertigo; tinnitus; granuloma; perforation of the tympanic membrane; post-surgery displacement of the prosthesis due to the development of scar tissue; post-surgical translocation of the prosthesis due to trauma or incorrect position; extrusion of the prosthesis; and post-surgical infection, including otitis media. It is also possible that the prosthesis could cause a significant loss in residual hearing. The responsible physician should ensure that the patient is fully apprised of the possible adverse effects of the implant procedure.

Additionally, patients are exposed to the normal risk of surgery and anaesthesia. Major ear surgery may result in numbness, swelling, or discomfort around the ear, the possibility of facial palsy, disturbance of balance or taste, or neck pain. If these occur, they are usually transient and resolve within a few weeks after surgery.

Warnings and precautions

Pre-operative

- The device is intended for single use only. Do not reuse or resterilize.
- Do not use if the sterile packaging is damaged, or if the use-by date (expiration date) has passed.
- Open the package only in a sterile environment to avoid contamination.
- Do not alter the device and only use it as intended.
- Do not use a damaged device.
- To eliminate unwanted electrostatic effects, water the prosthesis with a few drops of sodium chloride solution (0.9%) before you remove the prosthesis from its package.
- Do not implant if the patient is known to be intolerant to the implant material (titanium). Consider patient intolerances for surgical tool selection.
- Do not implant the device in direct contact with passive middle ear implants other than from MED-EL.

Intra-operative

Ossiculoplasty is a sophisticated surgical intervention and should be performed by a trained surgeon. Inadequate surgical technique can lead to irreparable damage.

MED-EL offers a variety of ossicular prostheses. The correct type and length of the prosthesis to be implanted is determined intra-operatively. Care must be taken when choosing the proper type and length of the prosthesis to be implanted to avoid post-operative complications (e.g. damage to the inner ear if the selected prosthesis is too long; insufficient contact between the prosthesis and the foreseen middle ear structures if the selected prosthesis is not in direct contact with the stapes head or stapes footplate).

For prostheses with shaft it is possible to carefully bend the shaft max. three times, max. 45° in any direction. For prostheses with a joint ball connection between head plate and coupling structure it is possible to carefully tilt the head plate, approximately 12° from initial position in any direction.

Avoid accidental bending of the prosthesis in order to prevent functional damage.

Avoid excessive pressure on auditory structures of the middle ear, especially on the stapes head or stapes footplate.

The middle ear structures must be treated with care, particularly the area of the oval window and the stapes footplate. Avoid traumatization of the inner ear due to luxation of the stapes footplate.

A slice of cartilage (recommended: 0.5 mm) must be placed onto the implant, just between the head of the implant and the tympanic membrane.

Post-operative

The patient shall be advised regarding proper hygiene of the operated ear. After surgery, both surgeon and patient must be vigilant for any signs of complications (such as development of a granuloma or loss of benefit over time).

The physician must inform the patient that the patient should avoid situations such as sudden, loud noises or severe ambient pressure variations and that the patient should consult his/her physician with questions regarding post-operative activities.

The Implant Card shall be filled out and handed over to the patient.

MRI safety information

Non-clinical testing has demonstrated that the device is MR Conditional.

A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of up to 7T
- Maximum spatial field gradient of 79T/m (7 900 G/cm)
- Maximum force product of 356T²/m (356 000 000 G²/cm)

In non-clinical testing, the worst-case image artefact caused by the device extends approximately 4.59 mm if imaged with a gradient echo pulse sequence in a 7T MR system.

Storage and disposal

The product must be stored in its sterile packaging. Devices must not be used after the use-by date specified on the package. Packaging should be disposed of in accordance with local requirements.

Warranty statement

Please contact your local MED-EL representative for information on our warranty provisions.

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MED-EL Elektromedizinische Geräte GmbH