Bonaduz, 29th January 2014

Reprocessing guide according to EN ISO 17664

Reprocessing the reusable Airway Adapter

Note This is only a copy of the reprocessing guide of Respironics. For the latest version or further questions contact Respironics directly www.respironics.com

General guidelines

1. Treat all reusable airway adapters in accordance with institutional protocol for single-patient use items. General guidelines:
 - a. To clean the airway adapter rinsing it in a warm soapy solution, then let it soak in a liquid disinfectant consisting of the following substances:
 - Isopropyl alcohol 70%.
 - 10% aqueous solution of chlorine bleach.
 - Gluteraldehyde 2.4% solution such as Cidex® or Steris System 1®.
 - b. Rinse thoroughly with sterile water and dry.
2. The Neonatal reusable Airway Adapters are not intended for use with steam sterilization.
3. The Reusable Airway Adapters may be sterilized using the methods listed below:

Method	Suggested temperature/time
ETO	38°C, 3 hours
Steam autoclave (adult only)	121°C, 20 min
Steam autoclave (adult only)	134°C, 20 min
2% Glutaraldehyde	20°C ± 5°C, 10 hours
PeraSafe™	20°C ± 5°C, 10 hours

Before reusing the adapter, ensure the windows are dry and residue free, and that the adapter has not been damaged.

Cleaning Test Criteria

The adapter testing criteria for cleaning included testing of the physical and dimensional integrity, optical performance, and gas leaks.

Methods tested

The test adapters were cycled 100 times for each method tested.

- Warm water rinse, cold disinfecting with Cidex or Steris Systems, pasteurization and autoclave
- Autoclave at 121°C (250°F), 20 minutes, unwrapped
- Autoclave at 134°C (250°F), 20 minutes, unwrapped

Test Results

There was no significant difference between the baseline data and the data recorded after 100 cycles.

Therefore, the methods described in this design guide and applicable user manuals are the only recommended cleaning and disinfecting methods for the Respironics Novamatrix Reusable CO₂ Airway Adapters.

Sterilization Evaluation

The following methods were evaluated to be effective:

Method	Suggested temperature/time	Tested for:
ETO	38°C, 3 hours	3 cycles at ½ cycle of 1.5 hours with 12 hour aeration; ETO residual with 1 hour extraction required to be ≤ 0.001 mg
Steam autoclave	121°C, 20 min	3 cycles at ½ cycle with gravity displacement autoclave for 10 min with 15 min dry time
Steam autoclave	134°C, 20 min	100 cycles at higher temperature to demonstrate that adult airway adapter could withstand this number of cycles. Since steam sterilization was shown to be effective at 121°C, this testing consisted only of autoclaving and visual inspection for damage after every cycle.
2% Glutaraldehyde	20°C ± 5°C, 10 hours	3 cycles; Cidex [®] Plus – Glutaraldehyde residuals tested to below 5.0 ppm standard.
PeraSafe [™]	20°C ± 5°C, 10 hours	3 cycles

With kind regards

Hamilton Medical AG

Melanie Obert

Product Manager
mobert@hamilton-medical.ch

☎ +41 58 610 18 86

📠 +41 58 610 00 20

Frederike Brühshwein

Senior Manager Regulatory Affairs

Declaration of Conformity

We,

Wir,

Nous,

Hamilton Medical AG, Via Crusch 8, CH-7402 Bonaduz, Switzerland,
(SRN: CH-MF-000013790)

confirm under our sole
responsibility that the
following products

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CEDCL-HAM-C6, Attachment on page 2

comply with:

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EC Medical Device Directive
93/42/EEC, Annex II, Art. 3

All listed products are
classified as class IIb.

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medin Medical Innovations GmbH
Adam-Geisler-Strasse 1
82140 Olching
Germany

Validity:

This declaration is valid for
products manufactured in
2023. Lot numbers are
traceable via
manufacturing protocols.
This declaration is valid in
connection with the final
inspection report.

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General terms & conditions

Patents

The products shown here are not available for purchase by the general public. The information provided here is intended for healthcare professionals only. Always read the labels and follow the product's instructions for use. Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. All images are used for illustrative purposes only and may not accurately represent the product or its use.

© Hamilton Medical All rights reserved.

America

CERTIFICATE

No. U8 058067 0085 Rev. 01

Holder of Certificate: Hamilton Medical AG
Via Crusch 8
7402 Bonaduz
SWITZERLAND

Certification Mark:

Product: Breathing therapy
Lung Ventilator

Model(s): HAMILTON-C6
HAMILTON-C6S

Parameters:

Rated Input Voltage:	100-240 Vac
Rated Frequency:	50/60 Hz
Rated Power:	510 VA
Ingress Protection:	IP22
Protection Class:	I
Applied Part Type:	B
Applied Part Type:	BF for SpO2, CO2 and Aerogen

Remarks:
The certificate is valid for USA and Canada.
Other countries might have other requirements,
which were not part of this certification.

Tested according to: ANSI/AAMI ES60601-1:2005/(R)2012
CAN/CSA-C22.2 No. 60601-1:2014

This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited Certification body.

Test report no.: 713196648

Date, 2021-10-15

(Stefan Reis)