## Declaration of Conformity

According to annex III of the Council Directive 98/79/EC on in vitro diagnostic medical device  
We,

### **EKVITESTLAB LLC**

**Velyka Vasylkivska St. 114, Kyiv, Ukraine, 03150, tel. 0(800)31-89-87; +38 (044)334-89-87**  
e-mail: [info@equitest.com.ua](mailto:info@equitest.com.ua), web-site: [www.equitest.com.ua](http://www.equitest.com.ua)

Declare under our sole responsibility that the following in vitro diagnostic medical devices  
other than those covered by annex II and devices for performance evaluation

#### **EQUI HAV IgM - ELISA kit for the qualitative detection of IgM antibodies to hepatitis A virus, REF EI-031**

Meet the provisions of the Council Directive 98/79/EC concerning medical devices which  
apply to them.

Undersigned declares to fulfill the obligations imposed by Annex III section 2 to 5:

- availability of the technical documentation set in Annex III (section 3), allowing the assessment of conformity of the product with the requirements of the Directive.
- the manufacturer shall take necessary measures to ensure that the manufacturing process follows the principles of quality assurance as appropriate for the products manufactured (Annex III section 4).
- the manufacturer shall institute and keep up to date a systematic procedure to review experience gained from devices in the post-production phase and to implement appropriate means to apply any necessary corrective actions (Annex III section 5).

Conformity assessment was performed according to Article 9 (7) and Annex III, section 3.

Our current Quality System is formatted to international standards:

- **ISO 13485:2016 «Medical devices — Quality management systems — Requirements for regulatory purposes»**

Corporate Contact Information

**EKVITESTLAB LLC**

Velyka Vasylkivska St. 114, Kyiv, Ukraine, 03150

tel. 0(800)31-89-87; +38 (044)334-89-87

e-mail: info@equitest.com.ua

RESPONSIBLE PERSON'S name: Anna Yurchuk

Position: Director

SIGNATURE :

Date : October 25, 2021

Stamp

European Authorized Representative:

Registered Address:

Obelis s.a.

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B-1030 Brussels, Belgium

Phone: 32.2.732.59.54

Fax: 32.2.732.60.03

E-mail: mail@obelis.net

Representative: Mr. Gideon ELKAYAM (CEO)

# CERTIFICATE

MANAGEMENT SYSTEM CERTIFICATION BODY  
«CONFORMITY ASSESSMENT BODY «PROMSTANDART», LLC  
certifies that the enterprise

**EKVITESTLAB**  
Limited Liability Company

registration code 38745936

legal address:

Ukraine, 03150, Kyiv, 114 Velyka Vasylkivska street,

manufacturer's address:

Ukraine, 04212, Kyiv, 60/2 Peremohy Avenue

has established and applies quality management system for  
**development, production, storage and sale**  
**of ELISA kits for in vitro diagnostic**

Audit, № report 2020/015-20.2.1  
confirmed that the requirements

**ISO 13485:2016**

**«Medical devices — Quality management systems —  
Requirements for regulatory purposes»**

are performed.

The control of conformity of the certified quality management system to the requirements of the specified standard is carried out by means of supervisory audit, the periodicity and procedures of which are regulated by the program.

Certificate registration number № UA.QMS.00014-21  
Registered 06 April 2021  
Valid until 05 April 2024

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1. Brillanti S., Masci C., Miglioli M. and Barbara L. Serum IgM antibodies to hepatitis C virus in acute and chronic hepatitis C. // Archives of virology. Supplementum. - 1993. – N8. – P.213-218.
2. Casino C., Lilli D., Rivanera D., et. al. Diagnostic value of anti-hepatitis C virus (HCV) core immunoglobulin M in recurrence of HCV infection after orthotropic liver transplantation. // Journal of Clinical Microbiology. - 1999. – Vol.37, N.8. – P.2726-2728.
3. CDC Testing Recommendations for Hepatitis C Virus Infection, 2015 // <https://www.cdc.gov/hepatitis/hcv/guidelinesc.htm>.
4. Dal Molin G., Tiribelli C. Campello C. A rational use of laboratory tests in the diagnosis and management of hepatitis C virus infection. // Annals of hepatology. - 2003. – Vol.2. – P.76-83.
5. Gupta E., Bajpai M. and Choudhary A. Hepatitis C virus: Screening, diagnosis, and interpretation of laboratory assays // Asian Journal of Transfusion Science. - 2014. - Vol.8. - P.19-25.
6. Lefrere J.J., Guiramand S., et al. Full or partial seroreversion in patients infected by hepatitis C virus // Journal of Infectious Diseases. - 1997. - Vol.175. - P.316–322.
7. Richter S. Laboratory assays for diagnosis and management of hepatitis C virus infection. // Journal of Clinical Microbiology. - 2002. – Vol.40. – N.12. – P.4407-4412.
8. Urdea M.S., Wuestehube L.J., Laurenson P.M., and Wilber J.C. Hepatitis C – diagnosis and monitoring. // Clinical Chemistry. - 1997. – V.43. – N.8. – P.1507-1511.
9. WHO Global Hepatitis Report, 2017 // <http://www.who.int/hepatitis/publications/globalhepatitis-report2017/en/>
10. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU.
11. Закон України «Про відходи» // Відомості Верховної Ради України. - 1998. - №36-37.
12. Наказ МОЗ України №325 від 08.06.2015 «Про затвердження Державних санітарно-протиепідемічних правил і норм щодо поводження з медичними відходами».
13. Постанова КМУ від 02 жовтня 2013р. №754 «Про затвердження технічного регламенту щодо медичних виробів для діагностики in vitro».
14. Hanna Tolonen, Kari Kuulasmaa, Tiina Laatikainen, Hermann Wolf and the European Health Risk Monitoring Project. Recommendation for indicators, international collaboration, protocol and manual of operations for chronic disease risk factor surveys Part 4.Storage and transfer of serum/plasma samples// Finnish National Public Health Institute 2002// [https://thl.fi/publications/ehrm/product2/part\\_iii4.htm](https://thl.fi/publications/ehrm/product2/part_iii4.htm)
15. Surveillance Guidelines for Measles, Rubella and Congenital Rubella Syndrome in the WHO European Region. Annex 3.Collection, storage and shipment of specimens for laboratory diagnosis and interpretation of results//Geneva: World Health Organization; 2012 Dec.

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
---

# HAV IgM

**ELISA kit for the qualitative detection of IgM  
antibodies to hepatitis A virus**

---

Instruction for use

IVD

REF  
EI-031

$\Sigma$  96  
tests

CE

# EQUI HAV IgM

ELISA kit for the qualitative detection  
of IgM antibodies to hepatitis A virus

## 1. INTENDED USE

The «EQUI HAV IgM» is ELISA kit intended for the qualitative detection of IgM antibodies to hepatitis A virus in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) to diagnose acute hepatitis A. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

**Target group:** blood or organ donors; pregnant women and children born to infected mothers; patients with symptoms of liver disease.

**Usage:** ELISA kit is used in clinical diagnostic laboratories, blood transfusion stations, as well as in other institutions working in the field of *in vitro* diagnostics.

## 2. CLINICAL SIGNIFICANCE

One of the most common foodborne infections is hepatitis A. The hepatitis A virus (HAV) causes acute liver disease, which can be mild or severe. Unlike hepatitis B and C, this hepatitis does not become chronic, but can lead to acute liver failure, which is characterized by high mortality.

Hepatitis A virus is a small shellless RNA virus from the *Picornaviridae* family. It is characterized by being highly stable in different environments and can be stored at + 4 ° C for several months, but becomes inactive after 5 minutes at 100 ° C. The virus replicates in liver cells, and then is released through the bile into the environment with the fecal masses of the infected person. The cellular immune response to HAV infection leads to the destruction of hepatocytes, liver dysfunction and the development of symptoms, typical for other types of hepatitis.

Acute hepatitis A, even with the clinical manifestations, does not differ from other viral hepatitis. Therefore, serological markers of infection are used for diagnosis, namely the detection of specific antibodies to HAV antigens. IgM antibodies are detected in the serum 1-2 weeks after infection, with the onset of symptoms or a few days before. In the maximum titer of anti-CAA IgM are detected in the jaundice period, after that their level gradually decreases. In most patients, specific IgM ceases to be detected after 6 months, and may occasionally circulate in the blood for more than a year. IgG antibodies to HAV antigens begin to be released shortly after IgM antibodies and stay in the blood throughout life. Also, specific IgG is produced after a vaccination. Detection of IgG antibodies to CAA indicates the formation of a stable immunity due to infection or immunization.

## 3. ANALYSIS PRINCIPLE

Detection of specific IgM antibodies to hepatitis A virus in the «EQUI HAV IgM» ELISA kit is based on the principle of «IgM capture» of solid-phase ELISA in a two-stage incubation. Monoclonal antibodies specific for human IgM immunoglobulins are adsorbed into the wells of the plate. During the first step of incubation of the

test samples in the wells of the ELISA plate, IgM immunoglobulins, if present in the samples, bind to monoclonal antibodies in the solid phase. The wells are washed to remove unbound components, leaving only specific antibody-antibody complexes. A mixture of HAV antigen and peroxidase conjugate of HAV-specific antibodies that bind to solid-phase immune complexes is then added. Unbound components are removed during washing. Immune complexes are detected by adding a solution of chromogen 3,3', 5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After a 30-minute incubation, the reaction is halted by adding a stop solution. Optical density (OG) in the wells is determined on a spectrophotometer at a wavelength of 450 / 620-695 nm. The intensity of the yellow color is proportional to the number of antibodies in the sample.

## 4. MATERIALS AND EQUIPMENT

### 4.1. Contents of the ELISA kit

<b>STRIPS</b>	1 x 96 wells	<p><b>Microplate</b></p> <p>Each plate well is coated with monoclonal antibodies specific for human IgM immunoglobulins. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months</p>
<b>CONTROL +</b>	1 x 0,25 ml	<p><b>Positive control</b></p> <p>The solution of human IgM immunoglobulins crosslinked with monoclonal antibodies specific for horseradish peroxidase, with a preservative (pink). Store at 2-8 °C</p>
<b>CONTROL -</b>	1 x 0,6 ml	<p><b>Negative control</b></p> <p>Negative human serum with a preservative (yellow). Store at 2-8 °C</p>
<b>DIL SAMPLE</b>	1 x 13 ml	<p><b>Serum dilution solution</b></p> <p>Buffer solution with monoclonal antibodies to human IgG, milk extract, detergent and preservative (brown). Store at 2-8 °C</p>
<b>CONJ 11x</b>	1 x 1,3 ml	<p><b>Conjugate (11x concentrated)</b></p> <p>11-fold concentrate of conjugate of antibodies to hepatitis A virus with horseradish peroxidase in buffer solution with stabilizers (purple).</p> <p>Dilute the conjugate (11x) 1:11 with the conjugate dilution solution before use (eg 100 µl concentrate + 1 ml conjugate dilution solution, enough for 8 wells). Diluted solution should be stored at 2-8 °C for no more than 1 day</p>
<b>DIL CONJ</b>	1 x 13 ml	<p><b>Conjugate dilution solution</b></p> <p>Buffer solution of inactivated hepatitis A virus antigen with detergent and preservative (yellow). Store at 2-8 °C</p>
<b>SOLN TMB</b>	1 x 13 ml	<p><b>TMB solution (ready to use)</b></p> <p>TMB solution, H<sub>2</sub>O<sub>2</sub>, a stabilizer, a preservative (colourless). Store at 2-8 °C</p>

<b>TWEEN WASH 20x</b>	1 x 50 ml	<p><b>Washing solution TWEEN (20x concentrated)</b> 20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days</p>
<b>SOLN STOP</b>	1 x 13 ml	<p><b>Stop Solution (ready to use)</b> 0.5 mol H<sub>2</sub>SO<sub>4</sub> solution (colourless). Store at 2-8 °C</p>

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

## 4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620-695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

## 5. PRECAUTIONS AND SAFETY

### 5.1. Precautions

*Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.*

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI HAV IgM» ELISA kit;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;
- **SOLN|TMB** solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of **SOLN|TMB** with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for conjugate solution and **SOLN|TMB**;
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

## 5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI HAV IgM» ELISA kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of [SOLN|TMB], [SOLN|STOP] and conjugate solution with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

## 5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

## 6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

## 7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by

centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

## 8. REAGENT PREPARATION

*NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!*

### 8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

### 8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

### 8.3. Conjugate solution preparation

Working dilution of the conjugate is prepared as follows: dilute **CONJ|11x** (purple) in a clean vial of solution **DIL|CONJ** (yellow) in the ratio 1:11 (ie, 1 + 10), the solution turns green. For example, for 8 well analysis add to 1 ml **DIL|CONJ** 100 µl **CONJ|11x**. The solution of the conjugate in the working dilution is stable during the day when stored at 2-8 °C.

## 9. ASSAY PROCEDURE

9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame.

Be sure to add control wells in every test run.

9.2. Fill in the sample application plan.

9.3. Prepare the detergent as per para. 8.2.

9.4. Add 90 µL of **DIL|SAMPLE** into each plate well.

9.5. Add 10 µL of controls and test samples into the wells:

**CONTROL|+** – into well A1,

**CONTROL|-** – into wells B1, C1 and D1,

and test samples into the remaining wells.

At the time of adding, the solution changes its colour from brown to blue.

Pipette the mix in the wells carefully to avoid foaming.

9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.

- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
- aspirate the content of all wells into a liquid waste container;
  - add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
  - aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
  - repeat the washing step 4 more times;
  - after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Prepare conjugate solution as per para. 8.3.
- 9.9. Add 100 µL of conjugate solution into each well. Cover the strips with a new piece of adhesive film and incubate for **60 minutes at 37 °C**.
- 9.10. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.11. Add 100 µL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.12. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.13. Add 100 µL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.14. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

*Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added in blank well).*

## 10. CALCULATION AND INTERPRETATION OF RESULTS

### 10.1. Calculation of results

Calculate the average OD of the negative control ( $\bar{Nc}$ ), Cut off (CO) and a sample positivity index ( $IP_{\text{sample}}$ ).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,3$$

$$IP_{\text{sample}} = OD_{\text{sample}} / CO, \text{ where } OD_{\text{sample}} \text{ is the OD sample.}$$

### 10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

CONTROL +

OD ≥ 1,2

CONTROL -

OD ≤ 0,150

CONTROL -

$\bar{N}_c \times 0,5 \leq N_{cn} \leq \bar{N}_c \times 2,0$  where  $N_{cn}$  is the OD for each  $N_c$  run

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and  $\bar{N}_c$  is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

### 10.3. Interpretation of results

$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE*
$IP_{\text{sample}} < 0,9$	NEGATIVE

\* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Analytical performance characteristics

#### Precision of measurement

##### *Intra assay repeatability*

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated in 32 replicates on one series of ELISA kits.

Sample No.	OD <sub>av</sub>	IP <sub>av</sub>	CV, %
37s/2	1,576	4,85	5,3
24s	2,462	7,57	4,3

##### *Inter assay reproducibility*

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD <sub>av</sub>	IP <sub>av</sub>	CV, %
37s/2	1,600	4,78	6,3
24s	2,463	7,36	7,1

#### Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 μmol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

## 11.2. Diagnostic characteristics

The diagnostic characteristics of the ELISA kit were evaluated by examining a set of 30 samples containing IgM antibodies to hepatitis A virus, a set of donor serum samples (188 samples) and samples of PHT202 Anti-Hepatitis A Virus (HAVed) panel. Performance Panel (contains 8 positive and 13 negative samples) - a total of 239 samples - was compared to similar commercial kits. For this set (239 samples) the relative sensitivity of the «EQUI HAV IgM» ELISA kit is 100%, the relative specificity - 100%, the percentage of coincidence - 100%.

## 12. LIMITATIONS OF ASSAY

A positive result in the «EQUI HAV IgM» ELISA kit is the evidence that the patient has IgM antibodies specific for hepatitis A virus. Anti-CAA specific IgM antibodies are usually markers of active replication of hepatitis A virus.

In order to counteract the false-positive results caused by the presence of autoantibodies specific for class G immunoglobulins (rheumatoid factor) in human serum samples, the kit uses a special block component that prevents the formation of immune complexes with anti-human antibodies in the solid phase.

The final diagnosis cannot be established solely on the basis of serological test results. When making a diagnosis the results of a set of laboratory and instrumental studies, as well as clinical manifestations of the disease should all be taken into account.

## 13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<b><i>High background in all wells</i></b>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<b><i>High background in a row of wells</i></b>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer

<b><i>Received OD of the positive control is below the border value</i></b>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<b><i>The colour density of the wells fails to meet the obtained optical density value</i></b>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

#### **14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE**

In case of technical problems, you can obtain assistance by contacting the manufacturer.

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Manufacturer

Authorized Representative in the European Community

In vitro diagnostic medical device

Catalogue number

Date of manufacture

Use by date

Batch code

Temperature limit

Contains sufficient for <n> tests

Caution

Non-Sterile

Consult instructions for use

Keep away from sunlight

Keep dry

Compliance with EU safety requirements

Edition 6, 01.11.2021

For questions and suggestions regarding the ELISA kit contact:

Obelis s.a.  
Bd Général Wahis 53  
1030 Brussels  
Belgium  
Tel: +(32)2 732-59-54  
Fax: +(32)2 732-60-03  
mail@obelis.net

Ekvitestlab LLC  
Velyka Vasylykivska St. 114, Kyiv, Ukraine, 03150  
Tel: 0(800)31-89-87, +38 (044)334-89-87,  
e-mail: info@equitest.com.ua, www.equitest.com.ua

## ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL|SAMPLE] into the wells  
(brown)

Add to 10 µl of controls and samples into the wells:  
A1 – [CONTROL|+], B1, C1, D1 – [CONTROL|-],  
other wells – examined samples  
(change of colour from brown to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN  
(300 µl per well)

Add 100 µl of prepared 1:11 (1+10) conjugate solution into all wells  
(green)

Cover strips with an adhesive film, incubate for **60 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN  
(300 µl per well)

Add 100 µl of [SOLN|TMB] into all wells

Incubate for **30 min in the dark at 18-25°C**

Add 100 µl of [SOLN|STOP] into all wells  
(change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at  
450/620-695 nm

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