Overview

Components

| Item | Description | P/N | Item | Description | P/N |
|------|---|--|---|--|------------------------------------|
| 1 | Central Control Monitor System Configuration Card Service Data Card | 802100 803739 803740 | 17 | Occluder Module | 803480 |
| | | | | Occluder Head | 806455 |
| | Solvice Bala Sala | 0001 10 | | Temperature Module | 802114 |
| 2 | Monitor Mounting Arm | 801441 | | Temperature Sensor (YSI (Yellow Springs Instruments) Series 700 temperature probes, or equivalent) | 24710008 |
| 3 | Roller Pump, 6 in | 801041 | | Pressure Module | 802112 |
| 4 | Roller Pump, 4 in | 801040 | | Pressure Transducer Kit | 16433301 |
| 5 | Centrifugal Control Unit | 801046 | | Pressure Transducer Cable | 806747 |
| 6 | Centrifugal Drive Motor Manual Drive Unit Disposable Pump (8/ctn) | 164267 164268 164275 | | Pressure Adapter Cable | 806748 |
| | | | 18 | Pressure Monitoring Kit | 16066100 |
| 7 | Telescoping Pole with solution rack Mounting Pole 4 ft Short Pole | 16431701 16553401 801407 | 19 | Flowmeter Module Flow Sensor Flowmeter Mounting Bracket Ultrasonic Sensor Gel | 802018 6382 801550 164278 |
| 8 | Crossbar | 16426 | 20 | Electronic Gas System | 801188 |
| 9 | Pole Collar | 801403 | | Oxygen Sensor | 801074 |
| 10 | Crossbar Fitting | 145980 | | Gas Supply Hose Kit (US) Oxygen Hose-green; Air Hose-yellow | 814475 |
| 11 | Bracket, Hanging Bags | 146819 | | Gas Supply Hose Adapter Sets (US) | |
| 12 | Lamp, Long Lamp, Short | 801238 801558 | NCG Hose Adapter Set DISS Hand Tight Hose Adapter Set Ohio Diamond Hose Adapter Set | 144207 144215 144223 | |
| 13 | Hand Crank Bracket (includes two hand cranks) | 802089 801016 | | | |
| 14 | Roller Pump Pole Mount | 801093 | | | |
| 15 | Centrifugal Unit Pole Mount Roller Pump Mntg. Ext. Large Roller Pump Mntg. Ext. Small | 804372 802523 802524 | | Gas Supply Hose Kit (Outside the US) | 814474 |
| | | | | Oxygen Hose-white; Air Hose-white/black | |
| | | | | Hoses do not include adapters for hospital gas outlets. | |
| 16 | Modules in Retaining Rack | | | 95-5 Supply Hose (%O2/%CO2) | 164595 |
| | Air Bubble Detect Module Air Sensor 3/8 x 3/32 | 802110 5773 5785 5791 | | Interface Module for CDI™ 101/CDI™ 100 Serial Cable | 802558 804981 |
| | Air Sensor 1/4 x 1/16 Air Sensor 1/4 x 3/32 | | | Interface Module for CDI™ 500 Serial Cable | 803479 804981 |
| | Cable Assembly Sensor Holder | 149892 149876 | | Interface Module RS-232 Serial Cable | 802113 804982 |
| | | | | Interface Module RS-485 Serial Cable | 803518 810166 |
| | Level Sensor II Module Level Sensor II Alarm Transducer, Red Level Sensor II Alert Transducer, Yellow Level Sensor II Mounting Pads (60 pack) Level Sensor II Gel Pad | 802111 195274 195215 195240 217390 | | Module Cover Right Side Left Side | 804405 804404 |
| | | | | | |



Level Detection System





Level Detect Overview

Purpose of Device

The Level Detect Module serves as the interface between alert and alarm level sensors and the rest of the system. This module is responsible for determining if the blood level has fallen below a specified point in the venous reservoir and for communicating that status to the system.

The level sensor attaches to a rigid reservoir with a disposable adhesive pad. Ultrasonic signals are used to detect the presence or absence of fluid at the sensor position for a range of venous reservoir thicknesses.

Affected Product

This document is applicable to the following products:

- 802111 Level Detect Module
- 195215 Ultrasonic Level II Sensor Yellow
- 195274 Ultrasonic Level II Sensor Red
- 195240 Level Sensor II Pads 60 PK

Occluder System

Occluder System

P/N 803480 Occluder Module

P/N 806455 Occluder Head



Overview

Purpose of the Device

The Occluder System provides a precision, computer-controlled, tube clamping mechanism used to regulate venous return flow in the perfusion circuit. The occluder system consists of an occluder module and occluder head.

The occluder head serves as the tube clamping mechanism. A section of tubing is positioned between the plunger and the occluder cap. The tubing occlusion can be adjusted by controlling the occluder head, using the central control monitor (CCM) controls to move the plunger toward or away from the occluder cap.

The occluder module serves as the interface between the occluder head and the system. The occluder head receives its power and control signals from the occluder module. The occluder head can be unplugged from the occluder module.

The user controls the occluder system from the CCM and can configure it to respond to the arterial pump, either roller or centrifugal types. The responses can be set to open, close, or go to a particular %Flow.

Operating Parameters

Occluder Requirements

Tube Size

The occluder supports tubing sizes in the range 1/4 in (6,4 mm) I.D. x 1/16 in (1,6 mm) wall thickness to 1/2 in (12,7 mm) I.D. x 3/32 in (2,4 mm) wall thickness.

Fluid Temperatures

20°C to 40°C

Occluder Performance

Voltage

24 VDC nominal


12 VDC \pm 5%




Roller Pump Repair

Service Procedure No. 875315 Rev F

| Roller Pump: External Components | | | | | | | |
|----------------------------------|--------------------------------|--|---|--|--|--|--|
| | Part or Sub-Assembly | Disassembly | Reassembly | | | | |
| 1. | Lid Assembly Molded Lid Hinge | a. Remove molded cap from top of occlusion knob.b. Remove lid assembly from race.c. Pull molded lid hinge from assembly. | a. Press buna inserts into lid (if new lid). NOTE: Use P-80 rubber lubricant if needed. b. Install molded lid hinge to lid. c. Install lid assembly onto race. | | | | |
| | Buna Drag Insert | NOTE: Applicable for lids P/Ns 803719 and 803720 only. Hinge is not replaceable for lids P/Ns 871928 and 871929. | d. For lid (P/Ns 871928 and 871929), apply linear lube to hinge o-rings if hinge squeaks. e. Insert cap onto top of occlusion knob. | | | | |
| | Occlusion Cap | NOTE: There is no disassembly required for drag inserts. | | | | | |









EC Certificate

Directive 93/42/EEC Annex II, excluding Section 4 Full Quality Assurance System Medical Devices

Registration No.: HD 60121893 0001

Report No.: 12031336 001

Manufacturer: Terumo Corporation

44-1, 2-chome, Hatagaya Shibuya-Ku, Tokyo 151-0072

Japan

Products: see attachement for products included

Replaces Approval, Registration No.: HD 60077473 0001

Expiry Date: 2022-08-29

The Notified Body hereby declares that the requirements of Annex II, excluding section 4 of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex II, section 5 of the aforementioned directive. For placing on the market of class III devices covered by this certificate an EC design-examination certificate according to Annex II, section 4 is required.

Effective Date: 2017-08-30

Date: 2017-08-25

M.Sc. M. Aihara

Notified Body

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg

TÜVRheinland

TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.



Doc. 1/2, Rev. 0

TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to

Certificate
Registration No.:

HD 60121893 0001

Report No.:

12031336 001

Manufacturer:

Terumo Corporation 44-1, 2-chome, Hatagaya

Shibuya-Ku, Tokyo 151-0072

Japan

Products included:

- Blood Bags
- Blood Donor Set with/without Blood Transfusion Filter
- Blood Transfusion Filter
- Intravenous Catheter
- Intravenous Administration Set
- Hypodermic Syringe
- Winged Needle
- Dental Needle
- Other Medical Needle
- Blood Administration Set
- Lancet

Date: 2017-08-25

Notified Body
TüvRheinland
M.Sc. M. Aihara



Doc. 2/2, Rev. 0

TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate

Registration No.:

HD 60121893 0001

12031336 001

Manufacturer:

Report No.:

Terumo Corporation 44-1, 2-chome, Hatagaya

Shibuya-Ku, Tokyo 151-0072

Japan

Products included:

- Extra-corporeal Membrane Oxygenator
- Cardiopulmonary Bypass Arterial Line Blood Filter
- Heart-Lung Bypass Defoamer
- Cardiotomy Reservoir
- Cardiopulmonary Bypass Blood Reservoir
- Haemoconcentration Filter
- Centrifugal Pump
- Angiographic Catheter
- Balloon Dilatation Catheter
- Catheter Guide Wire
- Guiding Catheter
- Catheter Introducer
- Stents
- Extension Tube
- Temperature Control Unit for Heart-Lung Bypass System Module
- Infusion Pump
- Syringe Infusion Pump
- Clinical Electronic Thermometer

Notified Body

Date: 2017-08-25

M.Sc. M. Aihara







Certificate

No. Q5 037584 0026 Rev. 01

Holder of Certificate: Terumo Medical Products

(Hangzhou) Co., Ltd. M4-9-5 Hangzhou Economic & Technological Development Zone

310018 Hangzhou

PEOPLE'S REPUBLIC OF CHINA

Certification Mark:



Scope of Certificate:

Design and Development, Production and Distribution of Connector, Balloon Catheter, Silicone Balloon Catheter, Silicone Balloon Catheter with Temperature Sensor, Extension Tube, Three-Way Stopcock with Extension Tube, Suction Catheter, P-Type Catheter with Aspirator Collector, Nelaton Catheter, Stomach Tube, Feeding Tube, Connecting Tube, Oxygen Catheter, Three-way Stopcock, Surplug Needle-Free Connector, Intravenous Hyperalimentation, Solution Administration Set for Infusion Pump, Solution Administration Set for Pediatric, Solution Administration Set, Volumetric Solution Administration Set for Infusion Pump, Disposable Volumetric Solution Administration Set, Volumetric Blood Administration Set, Winged Infusion Set, Winged Blood Sampling Set, Surplug Needle-Free Connector with Extension Tube, Surplug Solution Administration Set, Surplug Three-Way Stopcock, Final Filter PS, Automated Peritoneal Dialysis, Blood Administration Set for Infusion Pump, Extra-Corporeal Membrane Oxygenator, Digital Electronic Sphygmomanometer, Electronic Thermometer, Surplug AD, Surplug AD Manifold, Surplug AD Three Way Stopcock, Surplug AD Extension Tube, Surplug AD Administration Set, Intravenous Catheter for Single Use, Nasogastric Tube for Single

Use

Production and Distribution of Cardiotomy Reservoir

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

SH1907321 Report No.:

Valid from: 2019-12-01 Valid until: 2022-11-30

Date, 2019-11-19

Christoph Dicks

Head of Certification/Notified Body





Certificate

No. Q5 037584 0026 Rev. 01

EN ISO 13485:2016 Applied Standard(s):

Medical devices - Quality management systems -

Requirements for regulatory purposes

(ISO 13485:2016) DIN EN ISO 13485:2016

Terumo Medical Products (Hangzhou) Co., Ltd. Facility(ies):

M4-9-5 Hangzhou Economic &, Technological Development Zone,

310018 Hangzhou, PEOPLE'S REPUBLIC OF CHINA