**EN ENGLISH**

## Hydrophilic One-Piece Acrylic Foldable Intraocular Lens

### DESCRIPTION / INTENDED USE

SIDA-LENS hydrophilic one-piece acrylic foldable intraocular lenses (IOLs) are optical devices of high precision, manufactured by high-quality clinical grade UV-absorbing hydroxyl ethyl methyl acrylate (HEMA). They are available either in clear or yellow (with blue-light filtering chromophore bonded to the material) optic and in aspheric or spherical types (consult label on the carton box for the technical characteristics of each model).

These intraocular lenses are indicated for implantation in the posterior chamber of the eye, replacing the natural human crystalline lens that has developed Cataract, thus functioning as a refractive medium for the visual correction of aphakia in adult patients. Their 360° square-edge profile minimizes the chances of PCO, while their biconvex surfaces minimize effects of tilt, decentration and aberration.

### INDICATIONS

SIDA-LENS IOLs are indicated in the following (non-limited) conditions: in conjunction with or following extracapsular cataract extraction or phacoemulsification in the elderly; patients who are not good candidates for contact lenses; who cannot adapt well to cataract spectacles; who are of advanced age.

### CONTRA-INDICATIONS

Surgeons should weigh potential risk vs. potential benefit before implanting a lens. IOL implantation is not advisable in cases that it might worsen an existing condition, interfere with a diagnosis or treatment of pathology, or present a risk to patient's sight.

Some such conditions that have to be respected are: severe optic atrophy; chronic uveitis; microphthalmia; uncontrolled glaucoma; retinal detachment; corneal dystrophy.

### WARNINGS

Surgeon has to evaluate the risk which is involved in IOL implantation, as with all surgical procedures. Complications accompanying cataract surgery may include the following: retinal detachment; corneal endothelial damage; infection (endophthalmitis); hypopyon; corneal edema; persistent glaucoma; pupillary block. The list is indicative only.

The safety and effectiveness of IOL implants has not been verified in patients with preexisting ocular conditions. Alternative methods should be examined before considering the possibility of lens implantation.

Surgeons should continue to monitor patients postoperatively on a regular basis, as long-term effects on IOL implantation have not been determined.

An against-the-guidelines reuse of IOL may cause poor mechanical properties and cross contamination leading to microbial infection.

Some adverse reactions which have been associated with IOL implantation are intraocular infection, hypopyon and secondary surgical interventions, like lens repositioning, wound leak and retinal detachment repair.

### PACKAGING

SIDA-LENS hydrophilic foldable intraocular lenses (IOLs) are supplied sterile by steam, in a hydrated state, in a blister pack which is enclosed in an also sterile autoclaved pouch. It is strictly advised to open the peel pouch under sterile conditions only.

Contents remain sterile till the written expiration date as long as the package is not opened or damaged.

Apart from the lens case, each box contains: the appropriate cartridge and injector; enough product identification stickers for maintaining a record of the IOL implantation; a Patient ID-card, which patients are advised to save after their lens implantation; this informational leaflet.

### PRECAUTIONS

- Store the lens out of direct sunlight, at room temperature <40°C. Do not freeze.
- Do not use if sterile pouch is opened or damaged.
- Pouch should be opened only under sterile conditions and within short time before use.
- Only skilled surgeons with experience should attempt implantation.
- It is strictly advised to leave the Lens in saline or equivalent substance, to avoid dehydration until ready to fold for implant. For obtaining best results, the Lens should be folded and implanted within 3 minutes after removing from the saline.


- Do not soak or rinse lens in solution other than sterile balanced salt solution or sterile normal saline solution or equivalent of such.
- Handle the lens carefully without the use of locking forceps.
- Do not reshape the supporting structures (haptics).
- Do not re-sterilize this lens.

### DIRECTIONS FOR USE

1. In a sterile environment, peel apart pouches of injector and cartridge, and take them out.
2. Hold cartridge in 180° open position for loading the IOL and apply a thin layer of viscoelastic in the lens loading chamber and in the barrel, to serve as a lubricant for IOL's free movement during the injection.
3. Open lens' tray and take it out using forceps. Place it in the cartridge and gently press the lens using flat forceps. At the same time, close the cartridge without pinching the haptic or the edge of the optic.
4. Make sure neither the optic nor the haptic is caught between the shutters (If this happens, then open the cartridge, remove the lens and reload the lens according steps 2-3).
5. Load cartridge into the injector. Ensure that there is no gap between the shutters, before loading the cartridge into the injector.
6. Gently push the plunger until lens moves towards the cartridge tip end, till lens will completely come out of the cartridge. To avoid rotating the lens when it comes out of the cartridge, lens has to be properly loaded in the center of the cartridge and near the barrel.

### NOTES

- The A constant provided on the outer label is an estimated value and given as a guideline for implant power calculation. It is recommended that surgeon derives his/her own constant value based on clinical experience, surgical techniques, measuring equipment and postoperative results.
- SIDA-LENS intraocular lenses are available from -10.0 to +40.0 diopters (D) in steps of 0.5D or 1.0D, depending on the model and the diopter range. Out-of-range diopters are available upon request.
- Refractive Index of lenses: 1.46

### DISCLAIMER

The manufacturer will not be liable for any injury suffered to a patient as a result of: any implantation method or technique used by a surgeon to implant the lens; any prescription selection and use of the lens for any individual patient or patient's condition.

The manufacturer makes no expressed or implied warranties in connection with the sale of this intraocular lens. The manufacturer is liable only for the design and production of the intraocular lens and not for any incidental, consequential, indirect or exemplary damages of any kind, directly or indirectly arising from the purchase or use of this product.

  
**SIDAPHARM**<sup>®</sup>  
 MOCHINTRA DIANA "SIDAPHARM GREECE"  
 21, Stageiriti & 24 Em.Fili str.  
 Thessaloniki GR-54352, Greece, European Union  
 T: +30 2310906660  
 F: +30 2310989846  
 info@sidapharm.gr

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