SPHIGMOMANOMETRO ANEROIDE A-100 E A-50
ANEROID SPHYGMOMANOMETER A-100 AND A-50
TENSIOMÈTRE ANÉROÏDE A-100 ET A-50
ANEROID-BLUTDRUCKMESSGERÄT A-100 UND A-50

MANUALE D'ISTRUZIONE
INSTRUCTION MANUAL
MANUEL D'INSTRUCTION
BENUTZERHANDBUCH

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1. INTRODUCTION

This instruction manual is intended to assist the user for safe and efficient operation of the aneroid sphygmomanometers A-100 (with single head stethoscope included) and A-50 (stethoscope not included) (hereinafter: device).

The device must be used in accordance with the procedures described in this manual. It is important to read and understand the entire manual, especially the section “Information about blood pressure measurement”.

These devices are intended for the non-invasive measurement of systolic and diastolic arterial blood pressure. Always contact the physician if measurement is taken on children or persons with arrhythmia, as errors may occur.

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE
THE DEVICE MUST NOT BE DISASSEMBLED
FOR ANY PROBLEM PLEASE CONTACT THE TECHNICAL SERVICE OF DISTRIBUTOR
FOR SPECIFIC INFORMATION ON YOUR OWN BLOOD PRESSURE PLEASE CONTACT YOUR PHYSICIAN.**

SAFETY RULES

1. Before first use, open the packaging, check all components are enclosed and good working.
2. The medical device must be set-up and used according to information supplied with the present manual;
3. Respect the safety rules and particularly:
 - Children and/or disable persons may only use the device under the supervision of an adult in possession of his/her full faculties.
 - If the device is stored at temperature out by indicated environment temperature limits, do not use it immediately. Leave the device at room temperature for at least one hour.

- Avoid extreme environmental conditions and the exposure to direct sunlight.
 - **NEVER DIP ANY PART OF THE DEVICE AND ACCESSORIES IN WATER, IT ISN'T WATER PROOF.**
 - The device must be used only by adult people.
 - Cleaning the device with soft and dry cloth. Do not use petrol, thinners or similar solvents. Spot on the cuff can be removed carefully with a damp cloth. Do not wash any part of the device.
 - Do not wash or iron the cuff
 - Do not squeeze the air hoses nor the cuff
 - Keep the cuff away from cutting edges and do not twist it
 - Avoid to drop the device or threat it roughly in any way.
4. The device must not be disassembled. For any problem please contact the Technical Service of the Manufacturer. Failure to comply with the above invalidates the safety of the device and the guarantee claims.
 5. This medical device must be used exclusively for the use for which it has been designed and as described in this manual. Any different use must be considered incorrect and therefore dangerous; the manufacturer and the distributor cannot be considered liable for damage caused to people, things or animals by improper, incorrect and/or unreasonable use.
 6. Measurement on people with serious peripheral circulation or renal problems, people who have had an ictus or if they are unconscious and the case of serious arrhythmia may be difficult and bring extra error; it is recommended to consult the physician.
 7. In case of malfunctions contact the technical service of the Manufacturer. If by technical check the device will result opened by user, the guarantee of Manufacturer and Distributor will be expired.


A-100 AND A-50 CHARACTERISTICS

A-100 and A-50 devices are aneroid sphygmomanometers (only A-100 with stethoscope included), that with manual procedure using Riva Rocci/Korotkoff measuring method, can be used to measure with great reliability the systolic and diastolic pressures. These devices have high accuracy and are designed to be very easy to use. Thanks to manometer's dimension and readability, is easy to ready clearly the cuff pressure all measuring time.

WARNING: These devices can not provide reasonable accuracy if used or stored in the temperature or humidity beyond the range stated in the section "Technical Specification" of this manual.


2. INFORMATION ABOUT BLOOD PRESSURE MEASUREMENT

1. It is necessary to know that arterial blood pressure is subjected to sharp fluctuations. The level of the arterial blood pressure depends on many factors. Generally arterial blood pressure is lower in summer and higher in winter. Arterial blood pressure changes with atmosphere pressure and it is affected considerably by many factors, e.g. physical loads, emotional excitability, stress, diet, etc. **Medicines, alcohol, smoking, affect greatly the level of blood pressure.** When blood pressure is measured at the hospital, the value is always higher than that at home. The reason is the stress caused by the ambient, which is known as “White coat effect”. Blood pressure will raise in low temperature, so it is recommended to take the blood pressure at room temperature. If the device was stored at temperatures out of range recommended, it is necessary to leave it at room temperature for at least 1 hour, otherwise malfunction may occur. Blood pressure varies with age and individually, and it is recommended to write down daily the measures, to evaluate with physician the one’s “normal blood pressure”.

The illustration is from British Hypertension Society

2. In case of patients with cardio-vascular diseases, the measurement should be taken under physician’s instruction. Not in any case it is allowed to change the dosage of the medications suggested by the physician without his consultation.
3. An accurate measurement of blood pressure may be difficult in serious arteriosclerosis, weak pulses, or in patients with manifest fluctuation of heart contraction rhythm. It is recommended to consult a physician for a correct reading of the blood pressure measurements.
4. It is necessary to keep quiet during measurement to get accurate readings. Measurement should be conducted in quiet environment at room temperature. Don’t eat or smoke before a measurement. A-100 and A-50 supplied with the


standard cuff which is fit for the arm size 22-32 cm. Care should be taken to ensure that the cuff size is appropriate for the person whose blood pressure is being taken. Children and adults with cuff size fall outside the range 22-32 cm should select special size cuffs. Please contact the dealer to get these special size cuffs.
ATTENTION: Use only original cuff!

3. INSTRUCTIONS

Sit down comfortably and place the device on a table with everything ready for measurement. Make sure you have a stethoscope which is not included in models A-50.


3.1 FITTING THE CUFF

1. Insert the edge of the cuff approximately 5 centimeters into the D-ring as shown. Connect to one air tubes of cuff the aneroid manometer and to other one tube the inflation bulbe with valve.


2. Put the cuff on the left upper arm with the tubes toward the direction of the palm. If measurement on your left arm is difficult, you can use right arm. In this case, it is better to know that the readings may differ from the left one so is important to remember to measure always the pressure at same arm, so the reading can be compared and well evaluated.

3. Wrap cuff around your upper arm with the lower edge of the cuff approximately 2-3 centimeters above the elbow. The mark ARTERY must be over the artery of the arm so that the good arrow printed on cuff for each arm is exactly over the artery during measurement (Right Arm arrow for right arm and Left Arm arrow for left arm).

4. Press the cuff to make sure that it is safety closed. The cuff should not be too tight or too loose. One finger should be easily put in between cuff and upper arm.

5. Sometimes it is difficult to make a regular fit due to the shape of the user's upper arm, the cone-shape assembly of cuff is also acceptable.

6. If the clothes fits closely to the upper arm, or the sleeve has been rolled up take them off to get accurate measurement.

7. When cuff is wrapped around your arm, place under the cuff the head of stethoscope, **in contact with skin**. The head of stethoscope will must be under the cuff and in line with space into the zone between the two "ARTERIAL" arrows printed on the cuff.

Once the cuff is wrapped, make sure that its position is approx at the same level of the heart, and that the forearm is extended naturally on the table.

3.2 CARRY OUT A MEASUREMENT

It is necessary to keep quiet during measurement to get accurate readings. Measurement should be conducted in quiet environment at room temperature. **Do not eat, smoke or exercise before measurement. Measurements are advisable before meals or at least two hours after. In case measurements has to be repeated, at least 10-15 minutes should be waited.**

ATTENTION: Do not use cuffs other than the original cuff contained in this kit. DON'T INFLATE THE CUFF MORE THAN 300mmHg. The measuring time usually shall not exceed 2 minutes.

Before the measurement, take 3~5 minutes of deep breathing and relax. Don't talk or move the arm.

1. Hold the measuring aneroid manometer unit with the hand with cuff or use the hook in the back to attach it at the buttonhole needle on the cuff. The inflation bulb with the hand without cuff. Check that manometer air valve is fully closed (take it between thumb and forefinger and turn on it clockwise).
2. Pump the inflation bulb to increase pressure in the arm cuff. The manometer will show the pressure value within the cuff. Continue to pump until the pressure reaches approx 30~40 mmHg above the regular systolic pressure. If the regular systolic pressure is unknown, pump until 180 mmHg or more until when heartbeat is heard.
3. Stop inflation and wait 3-4 seconds (or 3-4 heartbeats) before you open slowly the air valve such to let the indicator go down by 2-3 mmHg/sec.
4. During pressure decrease, look continuously the manometer and read the values indicated at first heartbeat heard into the stethoscope. This is the Systolic Pressure value (Max pressure). Lower and lower heartbeats will be then heard in the stethoscope. The pressure indicated on the manometer in the moment of last heartbeat is the value of Diastolic Pressure (Min Pressure).
5. The measuring is finished and you can open the air valve of manometer to deflate fully the cuff.

3.3 FAST INTERRUPTION DURING MEASUREMENT

If you do not feel well during measurement or want to stop the measurement for any reason, open the air valve (turn on the valve anticlockwise) to deflate the cuff.

4. MAINTENANCE AND REPAIR

1. It's necessary to protect these devices and their accessories against high moisture and heat, direct sunlight, shock, solvents, alcohol.
2. This device don't needs ordinary maintenance, if necessary contact the technical service of distributor.
3. Keep the cuff and inflation bulb away from sharp subjects and don't extend or twist the cuff.
4. Use only soft and dry cloth to clean the device.
5. The cuffs are delicate and must be handled with care. You can clean the cuff cover with damp cloth.

WARNING: Under no circumstances any part of device and its accessories can be washed. Keep the device and all of its accessories out of the reach of children.

5. THOUBLESHOOTING

SYMPTOM	CHECK POINT	REMEDY
Manometer does not show any pressure during the cuff inflate.	<p>The manometer air valve is not fully closed The air tube is folded, bad connected or damaged. The cuff is damaged and leaks air.</p> <p>The manometer is damaged.</p>	<p>Close fully the air valve on the manometer. Connect correctly the air tube. Change the cuff.</p> <p>Call the technical service of distributor.</p>
The reading is extremely low or high	<p>Is the cuff at the same level as the heart? Is the cuff wrapped right? Did you move your arm during measurement? Did you talk or move your arm (or hand) during measurement?</p>	<p>Make sure that your posture is right Wrap the cuff correctly Relax during measurement Keep quiet and silent during the measurement</p>
Pulse rate is too low or too high	<p>Did you talk or move your arm (or hand) during measurement?</p> <p>Did you make measurement right after exercise?</p>	<p>Keep quiet and silent during the measurement</p> <p>Take measurement again after resting for almost 15 minutes</p>

6. TECHNICAL CHARACTERISTICS

Model	A-100 (REF PR 315100)
Size	Ø 50(D) x 90(H) mm
Weight	Approximately 110 g
Measuring method	Riva Rocci/Korotkoff
Measuring range	40 to 260 mmHg (pressure)
Pressure Range	0 to 300 mmHg
Measuring accuracy	± 3 mmHg for systolic and diastolic pressure
Inflation	Manual
Rapid deflation	Manual Valve / Fast deflation button
Operation temperature and humidity	+10° C to + 40° C , 20 to 85% RH
Storage temperature and humidity	-20° C to + 50° C, 20 to 85% RH
Cuff size	Applicable for arm size 22-32 cm
STETHOSCOPE ENCLOSED IN TO THE PACK	

Model	A-50 (REF PR 315050)
Size	Ø 50(D) x 90(H) mm
Weight	Approximately 110 g
Measuring method	Riva Rocci/Korotkoff
Measuring range	40 to 260 mmHg (pressure)
Pressure Range	0 to 300 mmHg
Measuring accuracy	± 3 mmHg for systolic and diastolic pressure
Inflation	Manual
Rapid deflation	Manual Valve
Operation temperature and humidity	+10° C to + 40° C , 20 to 85% RH
Storage temperature and humidity	-20° C to + 50° C, 20 to 85% RH
Cuff size	Applicable for arm size 22-32 cm
STETHOSCOPE NOT ENCLOSED	

7. SYMBOLS

CE 0123	CE Marking in conformity with EC directive 93/42/EEC and subsequent changes
	Warning!
	Consult the instruction manual
	To keep in cool and dry place
	Storage temperature
LOT	LOT number
REF	REF number

8. WARRANTY

CA-MI warrants its products for **24 months** after purchasing date. In front of this warranty, **CA-MI** will be obliged only to repair or substitute free of charge the products or parts of them that, after verification effected on our factory, or our authorized Service Center, by the Technical Service, results defective. **CA-MI** will substitute the device, only if a **SALE RECEIPT**. We also request to specify the kind of fault, in order to speed up repairing procedures. The product must be accompanied by a description of the defect. The warranty, with exclusion of responsibility for direct and indirect damages, it is thought limited to the solos defects of material or workmanship and it stops having effect when the device results however gotten off, tampered or sheltered out of the Factory or from the Authorized Service center. The commodity always travels to risk and danger of the buyer, without any responsibility of **CA-MI** for damages caused by the transport or dismay from the vector.

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ENG T-Classic / T-Flat, is a Mercury Free thermometer, useful for clinical measurement of body temperature when positioned in axillary, oral or rectal position.

Use and Storage conditions

Measurement principle:	Thermal expansion of a calibrated liquid alloy
Data Display:	Displayed mechanically on a calibrated graduated scale.
Measurement range:	35 - 42°C
Minimum scale interval:	0.10°C
Measurement accuracy (at 37°C and 41°C):	+0.10 to -0.15°C
Environmental operational limits:	+20°C - +30°C
Shipping/storage environmental limits:	0°C - 42°C

- The ideal environment temperature for correct measurement should be between 20 and 30°C. If the thermometer has been kept in an environment with a temperature below 20°C, before using it, it is necessary to keep it for 1-2 hours in an environment with a temperature that is between 20 and 30°C or kept for approximately ten minutes in a room with a temperature between 30 and 35°C.

- If environment temperature exceeds 35°C, immerse the tip of the thermometer in cold water for approximately 10-20 seconds and then shake it in order to lower the reference column and then position it immediately in the site used for measurement.

- Temperature must be below 42 °C when stored and transported in order to avoid irreparable damage.

- After use, always immediately lower the liquid metal column, then correctly sanitise/clean it and place it back in the protective plastic cover, in order to avoid any damage or breakage that may place it in a dangerous situation.

IMPORTANT: To lower the liquid metal column in the T-Classic model, carefully tighten the final part of the thermometer between the thumb and index finger and, paying attention not to crush it and not hitting it against anything, shake the thermometer downwards multiple times, until the liquid column does not go below 36°C. The same procedure can also be used with the T-Flat model that is however equipped with a plastic package that was specifically designed to safely and quickly lower the liquid metal column. If you have purchased T-Flat, it is sufficient to insert the device into its plastic container, tighten the cover well and then tightly hold the plastic bag located at its end between thumb and index finger, allowing the thermometer to swing quickly as if it was a pendulum, until the temperature shown falls below 36°C.

Measurement method

It is important to remember that the body temperature reading depends on where it is measured.

Axillary (positioning it under the arm)

If the thermometer is in axillary position, a surface temperature is obtained that can oscillate approximately between 0.5°C and 1.5°C compared to rectal temperatures in adults. Once the thermometer has been positioned, the patient must keep his arm tight against his body, without applying pressure but ensuring contact between the bulb and the axilla. The average measurement time for this method is approximately 5 minutes. It is worthwhile to note however that it is not possible to obtain an exact reading if the axilla are cooled down. In this case, it is recommended to extend measurement time by approximately 2 minutes in order to obtain the most precise as possible reading that corresponds as much as possible to body temperature.

Oral (positioned in the mouth)

In the mouth there are various hot areas. In general, oral temperature is 0,3° C - 0,8° C lower than rectal temperature. In order to guarantee a reading as precise as possible, position the tip of the thermometer on the left or right of the root of the tongue. The tip of the thermometer must be in constant contact with tissue while reading and must be positioned under the tongue. Keep the mouth closed while reading and breath regularly using nose. Do not eat or drink anything before measurement. Average measurement time is approximately 3 minutes

Rectal (positioning it in the rectum)

This method is the most precise from a medical point of view, since it is the closest to body temperature. The tip of the thermometer is carefully inserted into the rectum a maximum of 0.5 cm. The average measurement time is approximately 3 minutes.

Before proceeding with a measurement, always remember that:

1. The liquid metal column must always be brought under 36°C before starting a measurement.
2. Before each use it is necessary to correctly clean / disinfect/sterilise the thermometer (depending on the type of application, the environment where it is used and single or promiscuous use). Sterilise in medical use alcohol or in plastic sterilisation tubes.
3. In case of rectal measurement, it is recommended to apply a specific lubricant before insertion.
4. It is best to avoid physical activity, as well as taking a bath, eating or smoking at least 30 minutes before carrying out the measurement.
5. Avoid axillary measurement after applying refreshing deodorants or deodorants with particles of silver or other. Their application may interfere in part with a correct measurement.

IMPORTANT PRECAUTIONS

1. The use of clinical thermometers for oral temperature measurement is extremely inadvisable in children that are under 10 years of age;
2. Patients to whom the thermometer is applied orally, must pay careful attention to avoid chewing, talking, playing, running or moving until the thermometer is removed from the mouth;
3. Rectal use of the thermometer on children must be applied in a silent environment, under constant control and sustaining the thermometer with a hand in order to keep it from penetrating too deep or sliding outdoors.
4. Always verify, at least once a year, structural integrity of the device.

Safety precautions

- This thermometer must be used exclusively to read bodily temperature (fever).
- This thermometer was not designed for special applications (for example thermometers for premature babies, ovulation thermometers, etc.);
- The device must be used following the information provided in the attached documents.
- This device must only be used for its destined use and according to instructions supplied in this manual. A different use must be considered incorrect and dangerous.
- Manufacturer and Distributor cannot be considered liable for damages due to an improper use or storage.
- In case of malfunctioning, only contact the distributor's after-sales service or an after-sales service centre authorised by the Manufacturer.
- The device must not come into contact with hot water.
- Do not expose it to high temperatures or direct sunlight.
- Do not let the thermometer fall. It is not impact resistant.
- Do not clean it using solvents, petrol or benzene. Only clean using water or disinfectant.
- The thermometer body is made of glass, a very fragile material that is potentially dangerous in case of breakage. In the same way it is also dangerous if swallowed, the metal liquid contained inside may come out if the glass body were to break. For this reason, never leave the thermometer in children's hands.
- Persistent fever, especially in children, must be treated by a physician. In presence of fever, always consult your physician.

 USING THIS THERMOMETER DOES NOT REPLACE YOUR PHYSICIAN'S ADVICE

DISPOSAL

The body of this device is glass and inside it contains an alloy of liquid metals made of gallium (68% - 0.544g), indium (20% - 0.16g) and tin (12% - 0.096g) and considering its fragility, in order to avoid dispersing the materials, this clinical thermometer must always be handled with extreme care. In particular, one must pay attention to never let the instrument fall and never subject it to treatments or force that may cause damage or breakage. In case of breakage, always proceed by picking up glass and liquid metal that may have come out, using a soft cloth soaked with water and soap. Carefully pick up any residue and throw away even the cloth used for this operation.

Technical Specification

Model References :	TR.100200 – T-CLASSIC TR.100300 – T-FLAP
Size:	Bulb length= 14 ± 3mm Product length = 125±5mm H/W = 12.1/8.8 mm
Weight:	approx 7.2 g without plastic cover
Degree of protection against liquids: :	Waterproof

Legal requirements and guidelines


This product conforms to the European Directive on Medical Devices 93/42/EEC. The device conforms to the requirements foreseen by standard EN12470-1: 2000 + A1: 2009, "Clinical Thermometers - Part 1: Metallic liquid in glass thermometer with maximum degree".

Calibration control

This thermometer does not require any additional periodical or long term calibration.

Warranty

This product is covered by a 2 year warranty starting from the date the purchase document was issued. Damages caused by an incorrect use are not covered by warranty. The warranty is void in case of tampering, if the defect is caused by an improper use or if the damage is not a result of manufacturing defects (accidental fall, careless transport, etc.).

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
digiT-40 (rigid tip) ed digiT-40F (flexible tip) DIGITAL THERMOMETER

INTENDED USE:

Device **digiT-40 / digiT-40F** measures body temperature, which is measured at the metal measuring tip. This digital thermometer can be used under the tongue (oral), in the armpit (axillary) and in the rectum (rectal) to take temperature measuring results. An acoustic signal in "Bi, Bi, Bi..." will beep out once the temperature measuring process is complete. To ensure the best results and long-term satisfaction with this Digital Thermometer, we recommend that you read the instruction carefully before use and every time you will need in the future

SCHEMATIC DIAGRAM:

NOTICE:

- During the measurement, do not walk, run or talk.
- The thermometer should only be used under the supervision of an adult.
- Do not store the thermometer where it will be exposed to direct sunlight, dust or humidity. Avoid extreme temperature.
- Dropping or subjecting your thermometer to strong shocks should be avoided.
- Do not disassemble the thermometer, except to replace the battery.

CAUTION:

- The thermometer contains small parts (e.g. battery), which could be swallowed by children. Therefore, keep the thermometer inaccessible for children.
- If the device will not be used for an extended period of time, make sure you remove the battery when you store it.


OPERATING INSTRUCTION

RECTAL (WAY) MEASUREMENT:

This is the most accurate method from a medical point of view, because it comes closest to the core body temperature. The thermometer tip is inserted carefully into the rectum for a maximum of 2 cm, after having previously applied a water-soluble lubricating solution to the orifice and the tip of the thermometer. Never force the insertion of the thermometer in case of resistance. Keep the thermometer in position until the end measurement beep alarm is emitted. **The usual measuring time is approximately 40 seconds.**

ORAL (WAY) MEASUREMENT:




There are different heat zones in the mouth. To ensure that reading is as accurate as possible, place the thermometer tip to the left or right of the root of the tongue. The thermometer tip must have constant contact with the tissue during the reading and be placed under the tongue in one of the two heat pockets at the back. Keep the mouth closed during the reading and breathe evenly through the nose.

Do not eat or drink anything before the measurement (30 minutes before do the measurement at least). **The usual measuring time is approximately 60 seconds.**

AXILLARY MEASUREMENT:

Placing the thermometer in the armpit provides a measurement of surface temperature. The usual measuring time for this method is **approximately 90 seconds**. It should be noted, however, that an exact reading cannot be obtained if, for example, the armpits have been allowed to cool.

Under this case, we recommend extending the measuring time by around 5 minutes in order to obtain the most precise possible reading that corresponds as closely as possible to the core body temperature (this independently by end measurement beep sounds).

Note: We strongly recommend the rectal method as the most accurate method for identifying the basal temperature, and advise you to extend the measuring time by 3 minutes after the beep (independently by

measurement mode).

OPERATING:

- Press the ON /OFF button to turn on the thermometer. After a short beep indicates that the thermometer is now switched on. Simultaneously, a test display is performed in which all the symbols

to illuminate the display.

- The thermometer is ready to use when the "Lo" and a flashing "° C" appears in the display.

Notes: In an ambient temperature below 32 ° C will show "Lo ° C," appears in an ambient temperature above 42.9 ° C "Hi ° C".

During the measurement, the current temperature will be displayed continuously. During this period, the symbol "° C" blinks. The measurement is ended when a constant temperature is reached. The end of the measurement signal 10 beeps. The highest measured temperature is displayed.

To extend battery life, turn off the thermometer after use by pressing the ON /OFF button again. Otherwise the thermometer will automatically turn off after approximately 8 minutes.

MEMORY FUNCTION:

- Press the ON / OFF button to turn on the thermometer. The last measured temperature is automatically saved is displayed. The temperature appears with a small "M" in the display.
- You can now perform another measurement or turn off the power with the on / off button.
- The stored value is automatically overwritten when a new measurement is started.

CLEANING AND DISINFECTION:

To prevent infections the thermometer should be cleaned thoroughly after each use. We therefore recommend the use of commercially available disinfectants (e.g. a mild solution containing alcohol).

BODY TEMPERATURE:

Normal body temperature varies within a certain range. The following table shows the normal range depends on the measurement method. Therefore, measurement results can be directly compared different measuring points either.

Normal temperature ranges in ° C:

Axilla	34.7 - 37.3
Oral	35.5 - 37.5
Rectal	36.6 - 37.7

However, the normal range of a persons subject toother slight variations. Therefore, it is important to determine the respective normal range for your self. Run the temperature measurement by your self andother family members from - ideally in a healthy state - to determine the normal range for each individual.

BATTERY REPLACEMENT:

Replace the battery when the battery level indicator shows low battery status.

USE ONLY BUTTON BATTERIES 1.5 V DC, type LR41 or AG3


1. Place the thermometer with face back up and pull-out the battery cover by pressing it slightly downwards and at the same time pulling it outwards. CAUTION, during this operation the battery could fall to the ground easily if you apply excessive energy.

2. Take out the battery holder for approx. 18mm, take out the battery by a tool and replace it with a new one

3. Put back the battery cover in the right direction and make sure that it is properly locked in place.

INSTRUCTIONS FOR MAINTENANCE:

This device don't need of particular procedures of maintenance during its expected life except that battery replacement. Life expectation for this device is 3 years from the manufacturing data or 2 years from first use data (without batteries contained and in respect of storage conditions limits). After this period, is recommended to use device after an accuracy check done from qualified third part institute.

SPECIFICATIONS:

Model No.:	digit-40 (rigid tip model)// digit-40F (flexible tip model)
Display:	LCD,3.5digits, 16mmx7.5mm
Temperature measurement:	thermistor
Apply for:	Rectal / Oral / Axillary
Temperature measuring time:	40/60/90 seconds
Measuring range:	32.0 °C ~ 42.9 °C
Accuracy:	+/- 0.1 °C from 35.5 °C to 42 °C +/- 0.2 °C in the rest measuring range
Resolution:	0.1 °C
Memory records:	1 set
Low voltage Indication	YES
Automatically turn-off function	YES (after 8 minutes)
Power source and voltage:	DC 1.5V button batteries, type LR41 or AG3
Battery life:	Approx. 1200 times or 2 years
Dimension:	126mm×20mm×11mm
Unit weight:	Approx. 9.4g (digit-40) – 9.2g (digit-40F) including battery
Work environment:	Temperature 5°C to 40°C; Humidity: 15%RH~85%RH; Atmospheric pressure: 70KPa~106KPa
Transport/Storage condition:	Temperature -20°C~55°C; Humidity: 10%RH~85%RH; Atmospheric pressure: 70KPa~106KPa

*If there is any specification change for the product, these will applied without prior notice.

SYMBOL EXPLANATION:

1		Caution
2		Keep dry
3		Read the instruction before use
4		Type BF Applied Part Device
5	IP54	Degree of protection of an electrical device against accidental or intentional penetration of foreign solid bodies, and protection against the penetration of liquids. 1 st Number - Protected against dust (no harmful deposit) 2° Number - Protected against water sprayed from all directions
6		Low battery status
7		Date of Manufacturer
8		Manufacturer

CA-MI S.r.l. - Via Ugo La Malfa, 13 – Frazione Pilastro, 43013 Langhirano (PR) – Italy

This product is MADE IN CHINA

THE THERMOTER WAS DESIGNED AND MANUFACTURED IN COMPLIANCE WITH THE FOLLOWING EUROPEAN STANDARDS:

IEC 60601-1	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests
ISO 80601-2-56	medical electrical equipment part 2-56: particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement
ISO15223-1	Medical devices--symbols to be used with medical device labels, labelling and information to be supplied lied--part1:general requirements

This product complies with the European Directive for Medical Devices 93/42/EEC and carries the CE mark.

Electromagnetic compatibility:

The device complies with the IEC 60601-1-2 standard for electromagnetic compatibility. Please do not let thermometer become source of the disturbance, also try your best avoid source of the disturbance. Inquire at the Manufacturer for details on this measurement data. We reserve the right to make technical and design changes in the course of continuous product improvement .

WARNINGS!

- This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
- The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

Warranty

This product is covered by warranty for 2 years from the date of purchasing. Damage resulting from incorrect use or abuse is not covered by the warranty. Battery and packaging are excluded from the warranty. If you find that the thermometer is defective and not operating correctly, please firstly check the battery before sending in for repair. Warranty is voided in case of tampering, if the defect is caused by improper use or if the damage is not due to manufacturing faults (accidental drop, inaccurate transportation, etc.). Warranty is valid from date of purchasing (certified from Invoice or receipt).

IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2012/19/UE-WEEE: The symbol on the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, don't dispose it as mixed solid municipal waste, but dispose it referring to a specific collection center located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same. Caution: The wrong disposal of electric and electronic equipment may involve sanctions.

DISPOSAL OF SPENT BATTERIES (Directive 2006/66/EC): This symbol on the product indicates that the batteries should not be treated as household waste. Ensuring that the batteries are correctly disposed of contributes to preventing potential negative consequences for health and the environment, which could otherwise be caused by their inappropriate disposal. Recycling materials helps to preserve our natural resources. Take the spent batteries to collection points for recycling. For more information about collection and recycling of spent batteries or the product, please contact your Municipality, your local waste disposal service or the shop where you purchased the device.

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
The pack includes:

- Main Unit
- Cuff for 22-32cm diameter with air tube
- Travel Bag
- 4 alkaline batteries LR06 type AA - 1.5V
- Instruction Manual

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1. INTRODUCTION

1.1 FEATURES OF THE MY PRESSURE 2.0

This completely automated digital sphygmomanometer uses an oscillometric measurement procedure that permits fast and reliable measurements of the systolic and diastolic pressure and the heart rate.

MY PRESSURE 2.0 uses Fuzzy Logic technology to measure arterial pressure. As the cuff deflates, the inner system takes a pressure reading based on oscillations.

Once the reading has been taken, the rapid deflating process commences and the display shows the systolic and diastolic pressure together with the heart rates.

The high precision has been clinically tested and the system has been specially designed for the outmost ease of use.

Before using, read the contents of this manual carefully and then store it in a safe place. Contact a physician for further queries concerning blood pressure and how to measure it.

ATTENTION!

1.2 IMPORTANT INFORMATION ABOUT SELF-MEASUREMENT

- Do not forget: **self-measurement means Control**, not diagnosis or treatment. Unusual values must always be discussed with the doctor.
Any change in dosage of medications prescribed by the doctor should be discussed with him/her personally.
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE.
THE DEVICE MUST NOT BE DISASSEMBLED.
FOR ANY PROBLEM PLEASE CONTACT THE TECHNICAL SERVICE
OF MANUFACTURER AND/OR DISTRIBUTOR.
FOR SPECIFIC INFORMATION ON YOUR OWN BLOOD PRESSURE
PLEASE CONTACT A DOCTOR.
USE ONLY THE ORIGINAL DC ADAPTER.**

2. IMPORTANT SAFETY INSTRUCTIONS

It is necessary to know that arterial blood pressure is subjected to sharp fluctuations. The level of the arterial blood pressure depends on many factors. Generally arterial blood pressure is lower in summer and higher in winter. Arterial blood pressure changes with atmospheric pressure and is affected considerably by many factors, e.g. physical loads, emotional excitability, stress, meals, etc. Medicines, drinking, smoking affect greatly the level of an individual's blood pressure. Blood pressure does vary with age and individuals, and it is recommended to write down the readings from blood pressure records daily, then you can check with your doctor to find out what is a "normal blood pressure measurement" for you. Please read the instruction manual carefully before using this device, especially **IMPORTANT SAFETY INSTRUCTIONS**, it can help you use the device correctly and safely!

Please keep the instruction manual for future use. For specific information about your own blood pressure, consult your physician.

2.1 WARNINGS

- Consult your physician if you suffer from illnesses prior to using the device.
- The device is not suitable for persons who have electrical implants.
- If you had a mastectomy (breast amputation) do not use this blood pressure monitor on the arm on the side of the mastectomy.
- Pregnant women should only measure their own blood pressure in consultation with their doctor, since the readings may be changed with pregnancy.
- Do not service or maintain the cuff while in use with patient.
- Do not use this blood pressure monitor on any arm where intravascular access or therapy (such as an intravenous drip or a blood transfusion), or an arteriovenous shunt (A-V shunt) is present. The temporary interference to blood flow by the blood pressure measurement could result in injury.
- Do not use the device with other medical electrical (ME) equipment simultaneously.
- Do not use the device in the area the HF surgical equipment, MRI, or CT scanner exists, or in the oxygen rich environment.
- Do not use a mobile phone or other devices that emit electromagnetic fields, near the device. This may result in incorrect operation of the device.
- Never use any accessories or parts from other manufacturers. Using such accessories or parts could cause a hazardous situation for the user or damage to the device.
- Do not modify this equipment without authorization of the manufacturer.
- The batteries used in this device may present a fire or chemical burn hazard if mistreated. Do not disassemble, heat or incinerate.
- Keep equipment away from fire and heat sources to prevent fire or explosion
- Please keep the unit out of reach of infants, children or pets, since inhalation or swallowing of small parts can be dangerous or even fatal.
- Please pay attention to the continuous CUFF pressure due to connection tubing kinking, because it can cause a harmful injury.
- Do not use an extension cord with this device.
- The air tube or the AC adapter cable may cause accidental strangulation in infants.
- Do not put the air tube around your neck - danger of suffocation!
- A device should never be left unattended when plugged in.
- Do not reach for a corded device that has fallen into water. Unplug immediately.
- It is quite normal that two measurements taken in quick succession may produce significantly different results, because too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.

2.2 CAUTIONS

- Use this device under the right environmental conditions as indicated in this user manual. If not, this could affect the performance, life-time of the device and measurement results.
- Only use this device for its intended purpose as described in this user manual.
- Do not confuse self-monitoring with self-diagnosis. This device allows you to monitor your blood pressure. Do not begin, modify or end medical treatment based on the measurement results. Always consult your physician for treatment advice.
- Do not take any therapeutic measures on the basis of a self-measurement. Never change prescribed medication without consulting your physician. Consult your physician if you have any questions about your blood pressure.

- If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure.
- Consult the physician if measurement errors occur in children or persons with arrhythmia.
- The pulse display is not suitable for monitoring the frequency of cardiac pacemakers.
- Common arrhythmias (such as atrial or ventricular premature beats or atrial fibrillation) and peripheral artery disease / arteriosclerosis can affect the accuracy of this blood pressure monitor. Please consult your physician how to best use this blood pressure monitor if you suffer from any of these conditions. Blood pressure measurement is not suitable in cases of serious arteriosclerosis (hardening of the arteries).
- The effectiveness of this blood pressure monitor has not been established in pregnant women.
- Always check the device and cuff before you use it. Do not use the device or cuff if one of them is damaged, because this may cause injury.
- This device is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- Do not attach the cuff on the same arm on which other monitoring medical electrical equipment is attached simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring medical electrical equipment.
- Never attach the cuff on injured skin, an injured arm or an arm under medical treatment as this can cause further injury.
- Do not forcibly crease the arm cuff or the air tube excessively.
- Do not press the air tube while taking a measurement.
- Do not use the device in case of existing polyester or nylon material allergies.
- This device is not suitable for continuous monitoring during medical emergencies or operations.
- This device cannot be used with HF (High Frequency) surgical equipment at the same time.
- This device is not washable. Never immerse the device in water and do not rinse it under the tap.
- This device should keep dry to prevent from moisture.
- The equipment is not AP/APG equipment and is not suitable for use in the presence of a flammable anesthetic mixture with air, with oxygen or nitrous.
- To avoid measurement errors, do not use the device near strong electromagnetic fields, radiated interference signal or electrical fast transient/burst signal. For example magnets, radio transmitters, microwave ovens.
- If this device was stored in low temperature, leave it in right room temperature for at least 1 hour before use it.
- Repeated measurements with an interval of 3 minutes are recommended, so you can calculate the average to get a more accurate measurement. The interlude of 3 minutes can also ensure that the operation of the device does not result in prolonged impairment of the circulation of the blood.
- Atherosclerosis patients may require longer interval (10-15minutes) as elasticity of patient's vessels decreases significantly with the disease. 10-15minutes interval is also applicable for patients suffering from diabetes for a long period of time.
- Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.
- Connecting electrical equipment to MSO effectively leads to creating a ME system, and can result in a reduced level of safety.

2.3 CLASSIFICATION

- ME EQUIPMENT not intended for use in an oxygen rich environment or in the presence of flammable mixers.

- Internally powered equipment.
- Type BF applied part, recognize the cuff as applied part.

3. PUTTING THE BLOOD-PRESSURE MONITOR INTO OPERATION

3.1 INSERTING THE BATTERIES

After you have unpacked your device, first insert the batteries. The battery compartment is located on the bottom of device.

- Remove the battery cover;
- Insert the batteries (4 x size AA 1,5V), thereby observing the indicated polarity and close the cover.

Caution: the life span of the batteries included can be shorter than expected in consequence of storage period.

3.2 BATTERY REPLACEMENT

- During battery replacement the records in your device will not be lost. Regarding clock and date information, will be necessary to set again these because lost during this procedure.
- When weak battery symbol blinks on the display (see below) means that charge of batteries begin to be low and is good to replace batteries, but in case you can use again device. When battery is full discharged the weak battery symbol appear fix on the LCD and device sounds four times before switch off. In this last case replace immediately the battery.

Note: if the device cannot be actuated, please remove the batteries and wait some minutes, then replace the batteries again.

- When weak battery symbol appears on the display please replace with a set of four new batteries. Do not leave exhaust batteries into device compartment.
- Do not use rechargeable batteries. Use only batteries in compliance with specification declared in this manual.
- Remove the batteries from the compartment if the device will be used with DC adapter (optional, not included in the pack) for long periods.

3.3 USE WITH DC ADAPTOR (OPTIONAL ACCESSORY)

This device can operate with DC adaptor accessory an optional accessory to buy separately. To operate with adaptor, you must only connect the DC connector to the plug on the right of the device and connect the adapter to a main power socket adapt to your adapter (Fig. 2).

Please use only the original adapter accessory. In fact only this adapter model is certified for use with you MY PRESSURE 2.0. Remove the batteries from the compartment if the device will be used with DC adapter for long periods.



Fig. 2

3.4 DISPLAY SYMBOLS EXPLANATION

Memory 1
Memory 2

Low battery

Heart beating

A

Average of last three measurements

WHO
(World Health
Organization)
Blood
Pressure
Classification
(to see below
table)

188

Heart beating data

888

Diastolic
pressure data


88:88

Systolic
pressure data

88:88^D

Month/Day display
under reading data

The Pulse mark  also has function to advise, at the end of measurement, in case of Irregular Heartbeat (arrhythmia). When Heartbeat results irregular the Pulse mark  will flash automatically at the end of measurement.

Grade	Systolic	Diastolic
(**) Optimal	>120 mmHg	<80 mmHg
Normal	120 ~ 129 mmHg	80 ~ 84 mmHg
High - Normal	130 ~ 139 mmHg	85 ~ 89 mmHg
Mild Hypertension	140 ~ 159 mmHg	90 ~ 99 mmHg
Moderate Hypertension	160 ~ 179 mmHg	100 ~ 109 mmHg
Severe Hypertension	>180 mmHg	>100 mmHg

(**) **Note:** there is no universally accepted definition of hypotension. However, those having the systolic pressure below 90-100 mmHg are assumed as hypotensive.

3.5 SETTING THE DATE AND TIME

The function provides accurate measuring time for each measurement.

To get accurate date and time, the user should preset the date and time correctly before first use of device.

The operation procedure for presetting Date/Time is follows:

1. Press and don't release button M1, in the meantime press one time only the ON/OFF button, the year number flashes as in Fig. 4 now you can release the M1 button.
2. Press buttons M1 or M2 respectively to decrease or increase the number that blink on LCD.
3. At the end of setting, date and time will be displayed on LCD for 3 minutes; after this time these switch off and LCD will appears empty as in Fig. 5. At any time, pressing once the ON/OFF button, date and time will be displayed on LCD for three minutes.

Fig. 3


Fig. 4

Fig. 5

4. HOW TO MEASURE YOUR BLOOD PRESSURE

4.1 BEFORE THE MEASUREMENT

- Avoid eating, smoking as well as all forms of exercise before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quite atmosphere for about ten minutes before the measurement.
- Remove any garment that fits closely to your upper arm.
- Measure always on the same arm (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

4.2 COMMON SOURCES OF ERROR

Note: comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower (or higher) than the heart, an erroneously higher (or lower) blood-pressure will be measured (Fig. 6)!
- Cuff which is too narrow or too short result in false measurement values. Selecting the correct cuff is very important. The cuff size depends upon the circumference of the arm (measured in the centre). The cuff range is printed on the cuff. If this is not suitable for your use, will be necessary to use other good cuff.


Fig. 6

- **Warning: do not use cuffs other than the original cuff contained in this pack.**
- A loose cuff or a kinked tube may cause false measurement values.
- Please avoid strong electromagnetism interference, (such as mobile telephone, microwave oven, etc). This can cause low accuracy in the measurement.
- Use only alkaline batteries.
- Do not judge the measuring results and adjust your treatment yourself, it may be dangerous. Please give measuring results to physician to make a diagnose. Please take your medications as by physician's instructions.
- This equipment may stop working if there is a strong interference of electrostatic discharge.
- If the device is stored at low temperature, do not use it immediately. Leave the device at room temperature for at least one hour before using it.
- The Arrhythmia may cause mistakes in the measurements, so is good to use with great caution the device on peoples that have arrhythmia symptoms.

4.3 FITTING THE CUFF

- Close the cuff with the Velcro closer in such a way that it lies comfortably and is not too tight. Connect the armband tube to the appropriate opening on the left of the apparatus as shown in the Fig. 7.

Fig. 7

- Pass the end of the cuff (Fig. 8) (with sewn-in rubber stopper) through the metal stirrup so that a loop is formed. The Velcro closer must be facing outwards. (Ignore this step if the cuff has already been prepared.)


Fig. 8

- c. Push the cuff over the left upper arm so that the tube points in the direction of the lower arm.
- d. Lay the cuff on the arm as the Fig. 9. Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and the air tube lays on the inner side of the arm.


Fig. 9 Arteria principale (Brachiale)

Important: ensure that the reference sign (**ARTERY**) is over the Brachial Artery that is placed on the inside of the arm.

- e. Tighten the free end of the cuff and close the cuff with the closer.
- f. There must be no free space between the arm and the cuff as this would influence the result. Clothing must not restrict the arm. Any piece of clothing which does (e.g. a pullover) must be taken off.
- g. Lay the arm on the table (palm upwards) so that the cuff is at the same height as the heart. Make sure that the tube is not kinked (Fig. 10).


Fig. 10

- h. Remain seated quietly for two minutes before starting the measurement.
- i. If the “INDEX” line falls within the **NORMAL (22-32 cm)** range, indicated at the edge of the cuff, this cuff is suitable for the use. If the “INDEX” line falls out of the **NORMAL (22- 32 cm)** range, it could be necessary a special cuff

Note: if it is not possible to fit the cuff to the left arm, it can also be placed on the right one. However all measurements should be made using the same arm.

4.4 MEASURING PROCEDURE

1. Press ON/OFF button once to switch on device (date and time are so displayed on LCD), press it again to start with measurement.
2. The display will switch on with all symbols of LCD for one second.
3. About 2-3 seconds later, device sound two times quickly, appears on LCD the indication 0 and pressure into the cuff begin to increase while pressure value is displayed.
4. This monitor will inflate automatically if the system detects that your body requires more pressure for measurement. The cuff pressure will fall slowly during measuring. Do not shake the cuff or the monitor. Keep quiet. When the cuff pressure rises to 190 mmHg, the inflation movements stop and the indicated number drops down to start measuring.
5. Along with the measurement, besides the values of blood pressure which can detect if your blood pressure has problem or not, the heartbeat rate is also recorded (along with

the measurements, the heart sign will appear on the display, blinking to any heartbeat and detected if you have an irregular heartbeat (arrhythmia) or not. The heartbeat reading values will appear on the display at the end of measurement.


- After finishing measurement device issue a sound beep, the cuff deflates rapidly and automatically and the values of systolic and diastolic pressures and pulse rate will be displayed on the display together an automatic interpretation symbol (as indicated at paragraph 4.4 of this instruction manual) and the heart symbol (blinking or fix as below pictures). When device detects an irregular heartbeat or any excessive body movement during measurement, at the end of measuring will appear on the LCD the symbol 

 BRegular Heart Beat

 Irregular heart beat

Nota: It is recommended to consult a physician if at the end of measuring appear symbol of Irregular Heartbeat (Arytmia) or if values are in red range of colour scale.

- The measuring results will be shown automatically on the LCD until when the ON/OFF button is pushed again or till automatic switch off 3 minutes later. During the same time, the symbols   will flash to remind the user to record the reading.
- When the reading results are displayed, press M1 or M2 buttons to record data in the Memory Set.
- This device has average function. Press M1 or M2 to display the average of last three measurements so to help you to get your blood pressure values more correctly (under symbol of memory appears A indication that means that values on LCD is the average of last three measurements).

4.5 FAST INTERRUPTION DURING A MEASUREMENT

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the ON/OFF button can simply be pressed at any time. The device immediately reduces the cuff pressure.

4.6 MEMORY - DISPLAYING THE LAST MEASUREMENT

The device can store twin 90 sets (total of 180 sets) of recorded data, and each set data includes values of diastolic and systolic pressure, pulse rate and date/time, symbols of readings. After turning on the device by pressing the ON/OFF button, when M1 button is pressed will be displayed symbol  for memory channel M1 and the A letter (average) together average values of last three measurements recorder on Memory 1. To go on to read values of any measurement, press again M1 button, one time for any Regular Heart Beat Irregular heart beat measurement (on the display will be displayed the record number of any measurement

recorded, by more recently 01 to more old 02, 03 etc...). To see the record of memory channel M2 proceed with same procedure used for M1 but pressing M2 button. In this case will be displayed symbol .

Note: when memories are full, the more old values will be erased one to one, starting by more old, when you record a new measurement. To erase all records on one of two memory channels you must press for 5 second at least the button of memory you want erase.

When data are erased on LCD will be displayed Clr indication together the symbol of memory erased ( for M1 and  for M2).

FURTHER INFORMATION

Measurements should not occur soon after each other, since otherwise, the results will be falsified. Wait therefore for several minutes (at least 15 min.) in a relaxed position, sitting or lying, before repeating a measurement.

5. CARE, STORAGE AND REPAIR

1. It's necessary to protect this device against high moisture and heat, direct sunlight, shock, solvents, alcohol.
2. Remove the batteries if the device will be unused for a long time, and keep the batteries out of the reach of children.
3. Keep the cuff away from sharp subjects and don't extend or twist the cuff.
4. Use only soft and dry cloth to clean the device.
5. The cuff is delicate and must be handled with care. You can clean the cuff cover with damp cloth.

Warning: under no circumstances any part of device and its accessories could be washed. Keep the device and all of its accessories out of the reach of children.

6. ERROR MESSAGE

INDICATION	PROBLEM	CAUSE AND SOLUTION
Err	Device doesn't measure.	<p>Low Voltage or interior circuit problem. Remove and change the batteries.</p> <p>If the device is used with DC adapter, please check that the main connection plug is in compliance with rating indicated on the adapter.</p> <p>If after these checks the monitor still does not work, contact the customer service of Manufacturer and/or Distributor.</p>
	Inflation not enough or deflation fault.	<ol style="list-style-type: none"> 1. Cuff is not fastened well, please refasten it again according to the instruction. 2. Cuff is broken, change it with new one. 3. Pump, venting or deflation valve fail, please contact the Manufacturer and/or Distributor.
	Interferences and noises.	<ol style="list-style-type: none"> 1. Move or talk during measurement. 2. Posture, etc. was different for each measurement. 3. User is too nervous. 4. Repeat the measurement following the manual and consult physician if the problem persists.
	Irregular measurements.	<p>When measuring, the value will increase and be shown on the display</p> <ol style="list-style-type: none"> 1. If value shown is lower than 300 mmHg, check the correct position of cuff and repeat the measure. If it still does not work please consult a physician. 2. If value shown is higher than 300 mmHg, please contact the customer service of Manufacturer and/or Distributor.
	Over safety pressure.	<ol style="list-style-type: none"> 1. Air pipe blocked. 2. Inflation controller failure. 3. Pressure sensor fault. <p>Please check and repeat the measurement and/or ask to your Manufacturer and/or Distributor.</p>
	Low Battery.	Replace with four new batteries at the same time.

7. TROUBLESHOOTING

TROUBLE	CHECK POINT	SOLUTION
Though the batteries are installed, there is no indication on the display when the ON/OFF button is pressed.	<ol style="list-style-type: none"> 1. The batteries are not good positioned. 2. The batteries are exhausted. 3. The polarity of batteries is wrong. 4. The contacts of batteries compartment are dirty. 	<ol style="list-style-type: none"> 1. Remove batteries and wait some minutes before install yet these. 2. Replace the batteries. 3. Install the batteries correctly. 4. Clean the battery terminals with dry cloth.
With the DC adapter installed, there is no indication on the display when the ON/OFF button is pressed.	<ol style="list-style-type: none"> 1. The adapter is connected to a main socket with different rating. 2. The adapter is not well connected to the main socket or to the device. 3. The adapter is damaged. 	<ol style="list-style-type: none"> 1. Check if data on the adapter are suitable for the main socket. 2. Check the connections of adapter to the device and main socket. 3. Insert batteries in the device and check if with batteries the device works. If yes, the adapter must be changed because damaged.
The inflation action cannot be done or the air pressure cannot rise.	<ol style="list-style-type: none"> 1. The cuff position is not correct. 2. The cuff connection is not correct. 	<ol style="list-style-type: none"> 1. Check the cuff position and repeat the measurement. 2. Check the cuff connection and repeat the measurement.
The weak battery sign appears on the display.	Battery low.	Replace the batteries.
The blood pressure cannot be taken and on the display appears error message.	<ol style="list-style-type: none"> 1. The cuff is not well fastened. 2. User is nervous. 3. The posture is not correct. 	<ol style="list-style-type: none"> 1. Re-fasten the cuff. 2. Sit down and relax. Keep quiet during the measurement. 3. Keep the cuff and heart at the same height during the measurement. <p>If the patient has severe heart problem, then the blood pressure may not be read correctly.</p>
Under normal measurement circumstances, the reading at home is different from that of the clinics, or each measurement has various reading.	<ol style="list-style-type: none"> 1. The variation is due to the different environments. 2. The blood pressure changes with physiological or psychological status. 	Record everyday value trend and consult it with a physician.

8. TECHNICAL CHARACTERISTICS

Model	MY PRESSURE 2.0 (RE 503000)
Size	105 x 120 x 55 mm
Weight	Approximately 310 gr not including batteries
Measuring method	Oscillometry
Measuring range	40 to 260 mmHg (pressure) 0 to 160 beats/minute (pulse rate)
Measuring accuracy	± 3 mmHg for static pressure value ± 5% of the reading for the pulse rate
Inflation	Automatic
Rapid deflation	Automatic valve
Operating out	Automatic (after 3 minutes)
Batteries	4x1.5V LR06 alkaline batteries type AA
Operating condizions	Room temperature: 10 ÷ 40° C Relative humidity: 0 ÷ 85 % RH Atmospheric pressure: 80 ÷ 106 kPa
Storage condizions	Room temperature: - 20 ÷ 50° C Relative humidity: 0 ÷ 85 % RH Atmospheric pressure: 50 ÷ 106 kPa
Transport condizions	Room temperature: - 20 ÷ 50° C Relative humidity: 0 ÷ 85 % RH Atmospheric pressure: 50 ÷ 106 kPa
Cuff size	Applicable for arm size 22-32 cm

Note: Specification are subject to change without further notice.

9. MANUFACTURER'S DECLARATION

Compliance information for each EMC test	
Electromagnetic Emission (Home Healthcare Environment)	
Emission test (IEC60601-1-2:2014)	Compliance
Conducted and radiated RF emissions	CLSPR 11 Group 1 Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies

Compliance information for each EMC test		
Declaration-Electromagnetic Immunity (Home Healthcare Environment)		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3V 150 kHz to 80 MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz	3V 150 kHz to 80 MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz also meet the requirement of table 9 of 60601-1-2:2014	10 V/m 80 MHz to 2.7 GHz also meet the requirement of table 9 of 60601-1-2:2014
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV,±4 kV,±8 kV,±15 kV air	±8 kV contact ±2 kV,±4 kV,±8 kV,±15 kV air
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines
Surge IEC 61000-4-5	±0.5 kV,± 1 kV line(s) to lines	±0.5 kV,± 1 kV line(s) to lines
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000- 4-11	0% U_T , 0.5 Cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T , 1 Cycle and 70% U_T , 25/30 cycles single phase: at 0° 0% U_T , 250/300 cycles	0% U_T , 0.5 Cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T , 1 Cycle and 70% U_T , 25 cycles single phase: at 0° 0% U_T , 250 cycles
Power frequency (50/60Hz) ma- gnetic field IEC 61000-4-8	30 A/m	30 A/m

NOTE: The EUT is the a.c. mains voltage prior to application of the test level. The following phenomenon is still fulfill the requirement of basic safety and essential performance.

*UT:230V ~/50Hz.The pressure of the EUT is deviation the normal value but the value is still more than 10psi when flow is 4.5l/min.

**UT:230V ~/50Hz.The EUT stop working when adding 0%UT,but the EUT can restore its normal mode automatically.

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment(including peripherals such as antenna cables and external antennas) should be used no closer than 30cm(12 inches to any part of this device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Under the test condition specified in immunity, the product can provide the basic safety and essential performance.
- If the essential performance is lost or degraded,additional measures are necessary,such as reorienting or relocating the device.

10. SYMBOLS

	Device with BF Type Applied Part (cuff)
	CE marking in conformity with EC directive 93/42/EEC and subsequent changes
	Warning!
	Consult the instruction manual
	To keep in dry and cool place Device is not protected from the entry of water/solids
	Storage temperature
	Insulation Class II equipment
	Alternate Current
	Direct Current
	ON/OFF
M1 ou M2	Memory / Recall
	Manufacturer

11. WARRANTY

CA-MI warrants its products for **24 months** after purchasing date. In front of this warranty, CA-MI will be obliged only to repair or substitute free of charge the products or parts of them that, after verification effected on our factory, or our authorized Service Center, by the Technical Service, results defective.

CA-MI will substitute the device, only if a SALE RECEIPT. We also request to specify the kind of fault, in order to speed up repairing procedures.

The product must be accompanied by a description of the defect. The warranty, with exclusion of responsibility for direct and indirect damages, it is thought limited to the solos defects of material or workmanship and it stops having effect when the device results however gotten off, tampered or sheltered out of the Factory or from the Authorized Service center.

The commodity always travels to risk and danger of the buyer, without any responsibility of CA-MI for damages caused by the transport or dismay from the vector.

Software clinically validated in compliance with the international protocol of E.S.H. (European Society of Hypertension).

DISPOSAL

IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2012/19/UE-WEEE:

The symbol on the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, don't dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same.

Caution: The wrong disposal of electric and electronic equipment may involve sanctions.

DISPOSAL OF WASTE BATTERIES (Directive 2006/66/EC)

This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. The recycling of the materials will help to conserve natural resources. At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

my-pressure 2.0

Dispositivo medico di Classe IIa (93/42/EEC e s.m.i.)
Class IIa medical device (MDD 93/42/EEC and subsequent changes)
Dispositif médical Classe IIa (Directive 93/42/EEC et modifications ultérieures)
Geräteklasse IIa (93/42/EWG und nachfolgende Änderungen)

CE 0123

CA-MI

CA-MI Srl - Via Ugo La Malfa 13 - Frazione Pilastro - 43013 Langhirano (PR) Italia
Tel. +39 0521 / 637133 - 631138 - Fax. +39 0521 / 639041
E-mail: vendite@ca-mi.it / export@ca-mi.it www.ca-mi.it

CA·MI

DIAGNOSTICS

T·EASY

EN USER MANUAL

IT MANUALE D'USO

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DE HANDBUCH

THANK YOU FOR PURCHASING T-EASY

Be sure to read this Instruction Manual before using the unit for you to use it safely and correctly. This instruction manual should be well kept for your reference at any time.

- No contraindications.
- You should stop using the device and should consult with your physicians if you experience adverse reactions from the device.
- This device can be applied to all people. No special training is needed, but the user should read the manual carefully before use.

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1. FOREWORD

1.1 WARNINGS

- This unit is only used for human body temperature measurement without use of disease diagnosis; it cannot be used for emergency and continuous measurement in surgery.
- The patients cannot diagnose the disease and get treatment by themselves on the base of measurement result, they must follow the instruction of doctors.
- The main material of the product's body is plastic. Be careful to the potential allergic reactions to these materials.
- Please do not use with infectious wound.
- Please do not touch with your hands or blow infrared sensor with your mouth.
- Please ask professional doctors to explain the measured value of body temperature.
- Do not use mobile phones near the product. Please do not use equipment that generates electromagnetic fields near the product.
- The use of heat and cold producing devices, such as electric heating blankets, heating pads or ice packs, may impair the performance of device and increase the risk of injury to the patient.
- Do not store the unit under sunlight, at a high temperature, in high humidity or dust. Performance may be degraded.
- Please do not disassemble or repair this device by yourself including changing device.
- Please do not clean or maintain the device while using.
- Please wait for 30 minutes before use if the environment has a sudden change(e.g. from under sunlight to an air-conditioned room), or the accuracy may be influenced.
- The effects of degraded sensors and electrodes, or loosened electrodes, can degrade performance or cause other problems.
- DO NOT open the battery cover around any source of ignition, which has the potential to ignite the batteries and cause a fire.
- Please remove the batteries if it is not in use for a long time.
- Please keep the device out of reach of infants, children or pets. Inhalation or swallowing of small parts contained into device can be dangerous or even fatal.


1.2 PRECAUTIONS

- The patient is intended as operator. Under normal circumstances the patient can measure and change battery and maintains the device and its accessories in according to what described the user manual

- The kids should be kept away from this product. The kids until twelve years old, they can use this device only under supervision of their parents or adult people.
- Children under the age of 12 and / or people with disabilities that limit their ability to understand and use products like this, must not have access to the device and can only use it under the supervision of an adult with full mental faculties.
- Please do not bend or stretch the device hard.
- Do not impact or drop the thermometer or expose it to heavy shock. Do not subject the product to sudden and high temperature changes.
- Performance of the device may be degraded if: Operated or stored outside stated temperature and humidity ranges or if the patient's temperature is below the room temperature.
- If the thermometer has been stored at below freezing temperatures, please warm it naturally to room temperature before using.
- Clean the thermometer probe before storing.

2. DEVICE DESCRIPTION

This Infra-red Forehead Thermometer intended to measure human body temperature by measuring forehead. It is developed for displaying body or an object temperature under measuring, by detection the heat radiation emitted from his/its surface.


	Probe: Type BF Applied part
	ON - OFF button symbol
	Measuring / Memory button symbol
	Symbol for "SERIAL NUMBER". This symbol shall be accompanied by the manufacturer's serial number.
	Symbol for "Device reference code"
	Symbol for "LOT production reference"
	DATE OF MANUFACTURE. This symbol shall be accompanied by a date to indicate the date of manufacture.

	<p>Symbol for "MANUFACTURER". This symbol shall be accompanied by the name and the address of the manufacturer.</p>
<p>IP22</p>	<p>Degree of protection of an electrical device against accidental or intentional penetration of foreign solid bodies, and protection against the penetration of liquids. 1st DIGIT: Penetration of solids Protected against the entry of objects with a diameter greater than Ø12mm 2nd DIGIT: Penetration of liquids Protected against drops of water diverted up to 15° inclination</p>
	<p>Refer to user manual</p>
<p>CE 0123</p>	<p>Conformity indication with the essential health and safety requirements set out in European Directives.</p>

3. USING INSTRUCTIONS

3.1 BATTERIES PLACEMENT


- Press battery cover, the battery cover will unhook automatically. Prepare 2x 1.5V AAA batteries, install it into the battery chamber according to correct positive and negative poles.

- Press and release the button once. Backlight will come on and the start-up sequence begins. When the thermometer is ready and correctly working, a horizontal line of dashes ("---") will appear on the screen.

CAUTION: The environmental temperature variation can significantly affect the device measurement accuracy. This device can detect the environmental temperature and if its fluctuation goes beyond the normal range (eg. the device is moved from outside to room in Winter or the other way round), or the environmental temperature goes beyond 15°C-40°C (59°F-104°F), the device would display (“Err”). To avoid such (“Err”), leave the device for 30mins before operation if it is moved from different environment. And make sure the device always works in the correct range of environmental temperature.

3.2 THE SWITCH OF °F OR °C

In the battery slot, press °F/°C change button to switch between Fahrenheit degree (°F) and Celsius degree (°C).


3.3 MEMORY STORAGE FUNCTION

In OFF status, press “” button for 3 second, the product can read and save 30 sets of measurement values in order (as below pictures show).

Device will turn off automatically without operating for 30 seconds or by pressing “Power” button.

Note: Thermometer stores body temperature only, no object temperature.

3.4 MEASURING INSTRUCTIONS

- Place the thermometer on (the device also works in contact) or up to 3 cm away from the center of the forehead, just between and a little over the eyebrows. If the eyebrow area is covered with hair, sweat or dirt, please clean the area beforehand to improve the reading accuracy. It is important to hold


the thermometer and the forehead steady during measurement. Movement and wrong positioning of the device (to a different body area or with a distance greater than 3cm from the forehead) can affect the correct temperature reading.

- Press the ON/OFF button. When the thermometer is placed correctly, the display shows the temperature reading. The appropriate fever light color is displayed on screen and the confirmation beep is heard. Backlight will remain on 5 seconds and auto turn off.

- Visual/audible indications of the temperature readings. If the detected temperature is <math><32^{\circ}\text{C}</math> or $>43^{\circ}\text{C}$, a quick double beep will sound and the “Lo” or “Hi” symbols will appear on the display respecti-

vely and the measurement must then be repeated. In the case of normal body temperature (green backlight) and low fever (yellow backlight), the device will issue a single long beep 1 second long. In the case of body temperature corresponding to high fever (red backlight), 10 short and quick beeps will be issued. Backlight will remain on 5 seconds before than auto turn off.

Backlight Color	Temperature	Meaning
Green	<37.5 °C	Normal
Yellow	≥37.5°C and ≤38.0°C	Low fever
Red	≥38.0 °C	High fever

For their homogeneity of accuracy, between one measurement and the next one, the thermometer requires a technical interval of 5 seconds, during which it is not possible to do measurements. User can see the countdown from 5 to 0 seconds on the display.

To repeat measurement go to step 3.4

3.5 OBJECT MODE: Surrounding environment or object temperature measurement

In ON status, press "Power" button for 3 seconds, the product enters into object mode and LCD will display the "object mode" symbol in the upper right corner. Object mode measures the surrounding environment or object temperature. To ensure the accuracy of measurement, please do not conduct forehead measurement in object mode.

Note: No matter which mode is selected before turning the device off, the default value is always "body mode" after turning it on.

3.6 TURN OFF THE DEVICE

Device can be turned off manually by pressing "Power" button for 5 seconds, or automatically without using it for 30 seconds.

4. CLEANING AND DISINFECTION

4.1 CLEANING AND DISINFECTION

- About thermosensor (probe): if there are dust or other dirt in the mirror or tunnel of the sensor, please clean the probe with a cotton swab dipped denatured alcohol. (Note: please do not place the product directly under the faucet to wash.)
- Device housing: please wipe the product with a soft and dry cloth to avoid scratching of the product. (Note: please do not clean the product directly with water.)

Note: Infrared thermometers adopts highly sensitive technique to detect the temperature of the target object: any dust paste not only can affect the measurement accuracy, but also can cause bacterial infection. We suggest cleaning the device well after every use (as over indicated).

4.2 MAINTENANCE

- Before every use, check the device. Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Storage and use the device at the cool, dry and ventilated environment. Avoid to approach the fire and the heat source, or it will cause the battery explode.
- If you have any problems with the device, such as setting up, maintaining or using, please contact with service personnel of manufacturer or distributor. Don't open or repair the device by yourself.
- Please report to Manufacturer/Distributor if any unexpected operation or events occur.

5. TROUBLE SHOOTING

PROBLEM	POSSIBLE CAUSE	SOLUTION
	Batteries are not in proper direction or out of power	Make sure the batteries are installed in the proper direction or change batteries.
	Low power	Change batteries.
Abnormal results	Probe is not tightened properly or its position is incorrect.	Check if the probe is positioned properly, repeat the measurement procedure again.
	The probe is dirty and the measure time is not adequate	Clean the probe and take the measurement again with the proper measurement time.
	The probe is damaged	Contact service personnel
	Use outside of operating parameters/limits of use (e.g. use from outside to the room in winter or from sunlight to an air-conditioned room).	Wait 30 minutes and repeat the measurement following carefully the instructions and the specified environmental limits of use.
	$T < 32^{\circ}\text{C}$ (89.6°F)	Repeat the temperature measurement, carefully following the instructions for Use.
	$T > 43^{\circ}\text{C}$ (109.4°F)	

Note: Contact your local Distributor if you:

a) Need assistance in device setting up, using or maintaining.

b) Need to report unexpected operation or events.

5.2 SPECIFICATIONS

Product Name	T-Easy (Infra-red Forehead Thermometer)
Model	REF TR 100500/01
Power Consumption	Max. 50 mW in measurement mode
Rating	2×1.5V AAA alkaline batteries
Battery Life	More than 1000 measures (continuous operation)
Auto power-off	30 seconds
Dimension	149(L) x 34(W) x 52(H) mm
Weight	Approx. 67g (not include batteries)
Display Screen	Green <37.5°C, Normal Yellow ≥37.5°C and ≤38.0°C Low fever Red ≥38.0°C High fever
Measurement Range	Body Mode 32°C - 43°C (89.6°F - 109.4°F) Object Mode 0°C - 100°C (32°F -212°F)
Measurement Mode	Body Mode/Object Mode
Minimum Scale	0.1°C/0.1°F
Measurement Accuracy	±0.2°C , for range 35.0°C - 42.0°C ±0.3°C, outside this temperature range At Standard room temperature of 25°C (77.0°F)
Memory	30 memories
Button	Two buttons: - Power /Mode/ Measurement/Memory Button -°C/°F
Alarm	Approx. 5 seconds sound when peak temperature reached
Calibration	No need for calibration before use. But to ensure the measurement accuracy, please use and store the device in appropriate environment described below.
Working Environment	Temperature: 15°C - 40°C (59°F - 104°F) Relative humidity: 15%RH-93%RH Pressure: 70KPa to 106KPa

Storage and Transportation Environment	Temperature: -25°C - 60°C (-13°F - 140°F) Relative humidity: 0%RH-90%RH, non-condensing Pressure: 70KPa-106KPa
Expected Service Life	5 years

MADE IN CHINA

6. DISPOSAL

DISPOSAL OF SPENT BATTERIES **(DIRECTIVE 2006/66/EC)**

This symbol on the product indicates that the batteries should not be treated as household waste. Ensuring that the batteries are correctly disposed of contributes to preventing potential negative consequences for health and the environment, which could otherwise be caused by their inappropriate disposal. Recycling materials helps to preserve our natural resources. Take the spent batteries to collection points for recycling. For more information about collection and recycling of spent batteries or the product, please contact your Municipality, your local waste disposal service or the shop where you purchased the device.

WARNING REGARDING THE CORRECT DISPOSAL OF THE PRODUCT **■ ACCORDING TO THE EUROPEAN DIRECTIVE 2012/19/EU-WEEE:**

After use, the product must not be disposed of with household waste. It can be delivered to specific separate collection centres provided by the municipal administration, or returned to the distributor when purchasing a new device of the same type and having the same functions. Disposing of the product separately avoids possible negative effects on the environment and human health and allows the recovery of materials of which it is composed in order to obtain significant savings in energy and resources. The symbol on the data label indicates the separate collection of electrical and electronic equipment.

Attention: An incorrect disposal of electrical and electronic equipment could lead to fines.

7. ELECTROMAGNETIC COMPATIBILITY

Infra-red Forehead Thermometer has been tested and found to comply with the electromagnetic compatibility (EMC) limits for medical devices. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

CAUTION: do not use this device simultaneously with devices having high EMI levels.

MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

Infra-red Forehead Thermometer is intended for use in the electromagnetic environments specified below. The customer or the user of Infra-red Forehead Thermometer should assume that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Infra-red Forehead Thermometer use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Infra-red Forehead Thermometer is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage function / flicker emissions IEC 61000-3-2	Not applicable	

MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

Infra-red Forehead Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of Infra-red Forehead Thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2kV, ±4kV, ±8kV ±15kV air	±8kV contact ±2kV, ±4kV, ±8kV ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Electrostatic transient / burst IEC 61000-4-4	$\pm 2\text{kV}$ for power supply lines $\pm 1\text{kV}$ for input/output lines	Not applicable	Not applicable
Surge IEC 61000-4-5	$\pm 1\text{kV}$ differential mode $\pm 2\text{kV}$ common mode	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ ($>95\%$ dip in U^T) for 5 sec	Not applicable	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

<p>Infra-red Forehead Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of Infra-red Forehead Thermometer should assure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
	3 Vrms 150kHz to 80MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of Infra-red Forehead Thermometer, including cables, than the recommended separation distance calculated from the equation application to the frequency of the transmitter.

Conducted RF IEC 61000-4-6	6 Vrms in ISM bands	Not applicable	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m 80MHz to 2.5GHz	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p>
	385MHz-5786MHz Test specification for ENCLOSURE POR IMMUNITY to RF wireless communication equipment (refer to table 9 IEC 60601-1-2:2014)	385MHz - 5786MHz Test specification for ENCLOSURE POR IMMUNITY to RF wireless communication equipment (refer to table 9 IEC 60601-1-2:2014)	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$ <p>800 MHz to 2.5 GHz Where p is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m).^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

^a Field strengths from transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Infra-red Forehead Thermometer is used exceeds the applicable RF compliance level above. Infra-red Forehead Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Infra-red Forehead Thermometer.

^b Over the frequency range 150kHz to 80 MHz, field strengths should be less than 3V/m.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE EQUIPMENT OR SYSTEM - FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE - SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and Infra-red Forehead Thermometer

Infra-red Forehead Thermometer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Infra-red Forehead Thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Infra-red Forehead Thermometer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{F_1}\right]\sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1}\right]\sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{F_1}\right]\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations.

Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

For any question, please contact us, CA-MI Srl, first, we will sort the problem for you, please refer to the last page for contact information.

8. MANUFACTURER WARRANTY

Manufacturer guarantees its products for a period of at least 24 months from the purchase date. Based on this warranty the Manufacturer only undertakes to repair or replace the product or its parts that are found to be defective after verification at our premises by the Technical Assistance Service.

The product must be returned accompanied by a description of the defect found and it must be sent together with the Receipt. We always ask for the defect to be specified, in order to perform repair in the shortest possible time.

The warranty, with the exclusion of liability for direct or indirect damage, is limited only to material or workmanship defects and is not valid when the pieces returned turn out to have been disassembled, altered or repaired outside the factory or by authorised service centres. The returned device, even if under warranty, must be shipped CARRIAGE PAID. Goods always travel at the risk of the buyer, without any liability on the part of Manufacturer for damage in transit or loss by the carrier, even if sent carriage paid.

The batteries included in this device at the time of purchase are of a commercial type and their charging features may have deteriorated depending on their date of manufacture, their storage and time in the warehouse, therefore, they are not covered by the standard Manufacturer's Warranty.

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