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1. Christian De Geyter, Maria De Geyter, Peter R. Huber, Eberhard Nieschlag, and Wolfgang Holzgreve – Progesterone serum levels during the follicular phase of the menstrual cycle originate from the crosstalk between the ovaries and the adrenal cortex. Hum. Reprod., Apr 2002; 17: 933 – 939.
2. J. Jaroslav Stern, F. Voss, and C. B. Coulam – Early diagnosis of ectopic pregnancy using receiver – operator characteristic curves of serum progesterone concentrations. Hum. Reprod., May 1993; 8: 775 – 779.
3. B. Gellersen, M. S. Fernandes, and J. J. Brosens – Non-genomic progesterone actions in female reproduction. Hum. Reprod. Update, Jan 2009; 15: 119 – 138.
4. J. Dinny Graham and Christine L. Clarke – Physiological Action of Progesterone in Target Tissues. Endocr. Rev., Aug 1997; 18: 502 – 519.

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF PROGESTERONE IN HUMAN BLOOD SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of progesterone in blood serum or plasma.

This kit is designed for measurement of progesterone in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Progesterone is a gestagen with a MW of 314.5 Dalton. Progesterone is secreted by corpus luteum, adrenals and testis; it plays a role of a precursor for corticosteroids and androgens. Being an estrogen antagonist, Progesterone induces characteristic changes in endometrium necessary for implantation of an impregnated ovum.

During normal menstrual cycle, Progesterone level remains low until LH peak level begins to drop: only slight but statistically significant elevation of Progesterone level occurs concomitantly with LH peak followed by a decrease of Progesterone concentration. During second stage of the cycle, Progesterone and Estradiol levels increase again to complete luteinization. By the end of the cycle, Progesterone level drops again up to levels seen during follicular phase. This quick drop causes menstrual bleeding.

During pregnancy, Progesterone concentration continuously increases, and it induces proliferation and development of mammary glands and inhibits ovulation. During the first trimester, Progesterone is secreted by corpus luteum while from month 3–4 – by mitochondria of the trophoblast. Progesterone level in maternal blood increases rapidly – by week 7–8 it increases 2-fold and continues to increase by week 37–38. Decreased Progesterone levels indicate pathology of pregnancy while elevated levels suggest renal insufficiency.

Elevated Progesterone levels are found in pregnancy, tumours of adrenals or testicles, chorionepithelioma, in lipid tumours of ovaries as well as after intake of preparations of Progesterone or its analogues.

Decreased Progesterone levels are seen in galactorrhea-amenorrhea syndrome, in pregnant women at risk of premature delivery, and in persons taking some drugs such as oral contraceptives, ampicilline, ethynilestradiol.

### 3. PRINCIPLE OF THE TEST

This test is based on competition enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to progesterone-antibodies simultaneously with conjugated Progesterone-peroxidase. Progesterone from the specimen competes with the conjugated Progesterone for coating antibodies. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is inversely related to the quantity of the measured analyte in the specimen.

### 4. WARNINGS AND PRECAUTIONS

**4.1.** For professional use only.

**4.2.** This kit is intended for in vitro diagnostic use only.

**4.3. INFECTION HAZARD:** There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

**4.4.** Avoid contact with stop solution containing 5,0%  $H_2SO_4$ . It may cause skin irritation and burns.

**4.5.** Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

**4.6.** Do not use the kit beyond the expiration date.

**4.7.** All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

**4.8.** Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

**4.9.** Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

**4.10.** Do not mix reagents from different lots.

**4.11.** Replace caps on reagents immediately. Do not swap caps.

**4.12.** Do not pipette reagents by mouth.

**4.13.** Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

**4.14.** Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

**4.15.** The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

| Symbol         | Description                                                                                                             | Qty | Units | Colour code              | Stability of opened/diluted components                                                         |
|----------------|-------------------------------------------------------------------------------------------------------------------------|-----|-------|--------------------------|------------------------------------------------------------------------------------------------|
| 1 SORB MTP     | Progesterone EIA strips, 8x12 wells                                                                                     | 1   | pcs   |                          | until exp.date                                                                                 |
| 2 CAL 1-7      | Calibrator set, 0.8 ml each, zero calibrator C1 – 2 ml The set contains 6 calibrators: 0; 1; 3; 10; 30, 100, 300 nmol/l | 7   | pcs   | purple (C1 – colourless) | 2 months                                                                                       |
| 3 CONTROL      | Control serum (0.8 ml)                                                                                                  | 1   | pcs   | colourless               | 2 months                                                                                       |
| 4 CONJ HRP     | Conjugate, 22 ml                                                                                                        | 1   | pcs   | purple                   | until exp.date                                                                                 |
| 5 SUBS TMB     | Substrate solution, 11 ml                                                                                               | 1   | pcs   | colourless               | until exp.date                                                                                 |
| 6 BUF WASH 21X | Washing solution concentrate 21x, 22 ml                                                                                 | 1   | pcs   | colourless               | Concentrate – until exp.date<br>Diluted washing solution – 1 month at 2...+8°C or 5 days at RT |
| 7 STOP         | Stop solution, 11 ml                                                                                                    | 1   | pcs   | colourless               | until exp.date                                                                                 |
| 8 N003         | Plate sealing tape                                                                                                      | 2   | pcs   |                          | N/A                                                                                            |
| 9 K207I        | Instruction Progesterone EIA                                                                                            | 1   | pcs   |                          | N/A                                                                                            |
| 10 K207Q       | QC data sheet Progesterone EIA                                                                                          | 1   | pcs   |                          | N/A                                                                                            |

**5.2. Equipment and material required but not provided**

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for +37 °C ±0,1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

**5.3. Storage and stability of the Kit**

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

**6. SPECIMEN COLLECTION AND STORAGE**

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

**7. TEST PROCEDURE****7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

**7.2. Procedural Note:**

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

**7.3. Assay flowchart**

See the example of calibration graphic in Quality Control data sheet.

**7.4. Alternative units:**

1 nmol/l = 0.318 ng/ml

**7.5. Assay procedure**

|    |                                                                                                                                                                                                                                |
|----|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | Put the desired number of microstrips into the frame; allocate 16 wells for the calibrators CAL 1-7 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS. |
| 2  | Pipet 25 µl of calibrators CAL 1-7, control samples CONTROL and unknown samples into the wells.                                                                                                                                |
| 3  | Dispense 200 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape (included into the kit).                                                                                                                    |
| 4  | Incubate 120 minutes at 37 °C.                                                                                                                                                                                                 |
| 5  | Prepare washing solution by 21x dilution of washing solution concentrate (BUF WASH 21X) by distilled water. Wash the strips 5 times.                                                                                           |
| 6  | Dispense 100 µl of SUBS TMB into the wells                                                                                                                                                                                     |
| 7  | Incubate 10-20 minutes at +18...+25 °C                                                                                                                                                                                         |
| 8  | Dispense 100 µl of STOP into the wells.                                                                                                                                                                                        |
| 9  | Measure OD (optical density) at 450 nm.                                                                                                                                                                                        |
| 10 | Set photometer blank on air                                                                                                                                                                                                    |
| 11 | Apply lin-log method for data reduction.                                                                                                                                                                                       |

## 8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.


The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

## 9. CALCULATION OF RESULTS

1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.
2. Plot a calibration curve on graph paper: OD versus progesterone concentration.
3. Determine the corresponding concentration of progesterone in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.
4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

| Calibrators | Value      | Absorbance Units<br>(450 nm) x 1000 |
|-------------|------------|-------------------------------------|
| CAL 1       | 0 nmol/l   | 2473                                |
| CAL 2       | 1 nmol/l   | 2315                                |
| CAL 3       | 3 nmol/l   | 2007                                |
| CAL 4       | 10 nmol/l  | 1769                                |
| CAL 5       | 30 nmol/l  | 1427                                |
| CAL 6       | 100 nmol/l | 942                                 |
| CAL 7       | 300 nmol/l | 464                                 |

## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for Progesterone. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

| Sex, age         | Units, nmol/l |             | Units alternative, ng/ml |             |
|------------------|---------------|-------------|--------------------------|-------------|
|                  | Lower limit   | Upper limit | Lower limit              | Upper limit |
| Males            | -             | 4.0         | -                        | 1.27        |
| 12-17 yrs        | 0.3           | 4.3         | 0.1                      | 1.37        |
| Females          |               |             |                          |             |
| 12-17 yrs        | 0.3           | 41          | 0.1                      | 13          |
| post menopausal  | -             | 2.3         | -                        | 0.73        |
| Pregnancy:       |               |             |                          |             |
| 1st trimester    | 36            | 240         | 11.4                     | 76.3        |
| 2nd trimester    | 60            | 240         | 19.1                     | 76.3        |
| 3rd trimester    | 156           | 722         | 49.6                     | 229.6       |
| Menstrual cycle: |               |             |                          |             |
| follicular phase | 0.6           | 4.6         | 0.19                     | 1.46        |
| luteinic phase   | 7.5           | 80          | 2.39                     | 25.4        |
| ovulation        | 11            | 80          | 3.5                      | 25.4        |

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Analytical specificity / Cross reactivity

| Analyte                | Cross-reactivity, % wt/wt |
|------------------------|---------------------------|
| Progesterone           | 100                       |
| 17-Hydroxyprogesterone | 1                         |
| 11-Hydroxyprogesterone | 25                        |
| Corticosterone         | 0.01                      |
| Pregnenolone           | 0.9                       |
| Deoxycorticosterone    | 0.3                       |
| Deoxycortisol          | 0.03                      |
| Cortisol               | 0.002                     |

**11.2.** Analytical sensitivity. Sensitivity of the assay was assessed as being 0.25 nmol/l.

**11.3.** Linearity. Linearity was checked by assaying dilution series of 5 samples with different progesterone concentrations. Linearity percentages obtained ranged within 90 to 110%.

**11.4.** Recovery. Recovery was estimated by assaying 5 mixed samples with known progesterone concentrations. The recovery percentages ranged from 90 to 110%.

## 12. LITERATURE

1. Christian De Geyter, Maria De Geyter, Peter R. Huber, Eberhard Nieschlag, and Wolfgang Holzgreve – Progesterone serum levels during the follicular phase of the menstrual cycle originate from the crosstalk between the ovaries and the adrenal cortex. Hum. Reprod., Apr 2002; 17: 933–939.
2. J. Jaroslav Stern, F. Voss, and C.B. Coulam – Early diagnosis of ectopic pregnancy using receiver – operator characteristic curves of serum progesterone concentrations. Hum. Reprod., May 1993; 8: 775 – 779.
3. B. Gellersen, M.S. Fernandes, and J.J. Brosens – Non-genomic progesterone actions in female reproduction. Hum. Reprod. Update, Jan 2009; 15: 119 – 138.
4. J. Dinny Graham and Christine L. Clarke – Physiological Action of Progesterone in Target Tissues. Endocr. Rev., Aug 1997; 18: 502 – 519.