English Summary:  
Examination of the chemothermal disinfection procedure with  
**Citronensäure-Lösung 50 %®** for its virucidal activity


**Citronensäure-Lösung 50 %®** contains citric acid. **Citronensäure-Lösung 50 %®** was examined for its virucidal activity at 73° C and 83° C and concentrations of 3 and 5 % at 15 min. contact times using parvovirus, strain Haden, and the quantitative virucidal suspension test.

**Citronensäure-Lösung 50 %®** inactivated the highly resistant parvovirus, that is included as testvirus in the EN 14476, in a concentration of three and five percent solution and at 73° C and 83° C and 15 minutes contact time for more than 99.999 % or a five log<sub>10</sub> reduction of infectivity.

Parvovirus is a highly resistant virus and serves in the EN 14476 as surrogate virus for testing chemothermal disinfection procedures. Parvovirus was selected, because its inactivation is proof of inactivation of other clinical relevant viruses, such as hepatitis B virus (HBV), hepatitis C Virus (HCV), or human immunodeficiency virus (HIV)<sup>1)</sup>.

**Citronensäure-Lösung 50 %®** can therefore be recommended for the chemothermal disinfection procedure in a solution of 3 or 5 % at more than 73° C and a contact time of 15 minutes to inactivate clinical relevant viruses, such as parvovirus, enteroviruses, HBV, HCV und HIV.

Luckenwalde, 26. January 2003

PD Dr. Thraenhart  
Director of Eurovir

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D-34209 Melsungen

Date: May 20, 1998  
 Ohg./Br

Test of the sporocidal efficiency  
 of Citric Acid 50 % (120 ml)  
 at 83 °C in the dialysis machine Dialog 710/700/5  
 (B. Braun Melsungen AG)

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Seiten: 6

## 1. Question

It should be tested, whether the Citric Acid 50 % preparation (120 ml) has a sufficient disinfection efficiency at 83 °C in the dialysis machine Dialog 710/700/5.

## 2. Procedure

The disinfection performance of the machine was tested after artificial contamination of the fluid side circuit. A  $\log_{10}$  reduction factor of minimum 5.0 should be reached at the various check points of the water conducting system.

### 2.1 Test Germs

Spores of bacillus subtilis, variatio globigii were applied, which were suspended in ethanol at a concentration of  $3 \times 10^{10}$ /ml.

These test germs were applied, because the standard vegetative germs, which are commonly used for disinfectant agent tests (e.g. e. coli, pseudomonas aeruginosa, candida albican etc.) are not sufficient for the evaluation of highly efficacious procedures. In general streptococcus faecium is the most resistant vegetative (non sporodial) test germ. A reduction factor of 8 for this test germ is equivalent to a value of approx. 5 with b. subtilis spores. A germ number of  $10^9$  (final dilution in the machine) should be applied for the practical tests. Therefore there would be an unrealistic high contamination with biological material.

### 2.2 Artificial Contamination of the Dialysis Machine

Prior to a test run 1 ml of the test germ suspension was injected into the water circuit system, via the mixing chamber. Thus a final concentration of the test germs was realized, which was between  $10^7$  and  $10^8$  /ml of the dialysing fluid.

The artificially contaminated dialysing fluid was recirculated for 10 minutes, to guarantee an uniform distribution of the test germs.

The samples were then taken prior to disinfection at the measurement point 3 (see table 1).

### 2.3 Test Cycles

The disinfection program was started with 120 ml Citric Acid 50 % after the artificial contamination.

Samples were taken after disinfection after completion of the disinfection program at the measurement points 1 through 6 (see table T).

A total of 3 independent test cycles were performed. The tests were performed on April 28, 1998, in the above mentioned institute.

### 2.4 Evaluation

The number of germs were determined, if necessary, via diluted test rows in sterile CPS broth.

The log reduction factor was calculated as a comparative figure for the obtained disinfection efficiency:

$$\log_{10} \frac{\text{Number of Germs prior to Disinfection}}{\text{Number of Germs after Disinfection}}$$

### 3. Results

The germ contamination prior to disinfection for all test cycles was higher than  $10^7$  colony forming units per ml.

The obtained reduction factors are listed in table 2. The figures are arithmetical mean values from double determinations.

### 4 Evaluation

The necessary reduction factors of 5.0 were exceeded at all measurement points, with 120 ml Citric Acid 50 %, at 83 °C.

The tested method is as sufficient as other methods, due to the detected high sporocidal efficiency, for decontamination, e.g. for hepatitis.

(Prof. Dr. med. W. Solbach)

(Prof. Dr. med. H. Ohgke)

Prof. Dr. med. W. Salbach  
Institut für Med. Mikrobiologie und Hygiene  
Med. Universität zu Lübeck  
(Institute for Med. Microbiology and Hygiene  
Med. University of Lübeck, Germany)  
Ratzeburger Allee 160  
23538 Lübeck  
Tel.: 0451/500-2800; Fax: 0451/500-2749

Table 1 for the expertise from May 20, 1998

Measurement points at the Dialysis Machine Dialog 710/700/5

| No. | Location                    |
|-----|-----------------------------|
| 1   | Upline tank                 |
| 2   | In front of balance chamber |
| 3   | Dialyser, arterial inlet    |
| 4   | After balance chamber       |
| 5   | Outlet tubing               |
| 6   | Suction rod                 |

Prof. Dr. med. W. Solbach  
 Institut für Med. Mikrobiologie und Hygiene  
 Med. Universität zu Lübeck  
 (Institute for Med. Microbiology and Hygiene  
 Med. University of Lübeck, Germany)  
 Ratzeburger Allee, 160  
 23538 Lübeck  
 Tel.: 0451/500-2800; Fax: 0451/500-2749

Table 2 for the expertise from May 20, 1998

Test of the Sporocidal Efficiency of Citric Acid 50 % (120 ml) at 83 °C and 15 Minutes  
 Reaction Time in the Dialysis Machine Dialog 710/700/5 from April 28, 1998

| Measurement Point | Test Cycle |       |       |
|-------------------|------------|-------|-------|
|                   | 1          | 2     | 3     |
| 1                 | 6.6        | > 6.6 | > 6.6 |
| 2                 | > 6.6      | 6.6   | 6.3   |
| 3                 | > 6.6      | 6.6   | > 6.6 |
| 4                 | 6.0        | 6.6   | 6.3   |
| 5                 | 6.6        | 6.6   | 6.6   |
| 6                 | > 6.6      | > 6.6 | 5.9   |

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