

Documentation No.	Appendix F Instructions for Use
Creation Date	2022/2/18
Revision No.	1.0
Revision Date	/

Technical File of EC Type Examination Certification

Instructions for Use

This product reference to PPE Regulation (EU) 2016/425 and has been tested and approved in accordance with EN 149: 2001+A1: 2009.

Name: Respirator protective face mask



Model: LMRPFMWE01

Class: FFP2 NR

Applicable standard: EN 149: 2001+A1: 2009

Description: This product provide effective respiratory protection for use in industries where workers will be exposed to dust particles and/or non-volatile liquid particles.

Fitting instructions:



Before fitting respirator, ensure hands are clean.

1. Ensure you are using a clean respirator. Inspect the respirator to ensure there are no visible signs of damage before each use.
2. With the nose clip facing up, open the mask, make the inner side of the mask face to yourself, and pull a mask strap with each hand.
3. Put the mask on the chin.
4. Pull the mask strap behind your ear and adjust it to make it as comfortable as possible.
5. Place the fingers of both hands in the middle of the nose clip and press inwards while moving the fingertips along the clip to the sides until the clip is fully attached to the bridge of the nose. Pinching the nose clip using only one hand may result in improper fit and less effective respirator performance.
6. Before entering the contaminated area, the seal of the respirator on the face should be fit-checked as follows:
 - a) Cover the front of the respirator with both hands being careful not to disturb its fit, exhale sharply.
 - b) If air leaks around the nose, re-adjust the nose clip to eliminate leakage.
 - c) If air leaks at the respirator edges, re-adjust the straps to eliminate leakage.

If you CANNOT achieve a proper fit, DO NOT enter the contaminated area. See your Supervisor.

Documentation No.	Appendix F Instructions for Use
Creation Date	2022/2/18
Revision No.	1.0
Revision Date	/

Technical File of EC Type Examination Certification

Warnings and Use Limitations

1. See information supplied by the manufacturer before use these respirators. Misuse may result in sickness or death. For more information, see your occupational safety supervisor or contact with manufacturer.
2. Failure to follow all instructions and limitations concerning the use of the respirator and/or failure to wear the respirator for the whole duration of exposure to contaminants may seriously reduce the performance characteristics of the respirator and lead to illness, injury or death.
3. This respirator marked “NR”, shall not be used for more than one shift.
4. This respirator helps protect against certain particulate contaminants but does not completely eliminate exposure to the risk of contracting disease or infection.
5. Before use respirator, the wearer must first be trained in the correct use of the respirator in accordance with national regulations and checks prior to use.
 - Prior to use of respirator, always check the shelf life and ensure there are no visible signs of damage.
 - A filtering half mask should never be used if it has become damaged or has exceeded its shelf life.
6. This respirator need to be fit-checked and inspected for damage prior to each use.
7. These respirators do not contain components made from natural rubber latex.
8. Do NOT use the respirator or enter or stay in a contaminated area under the following circumstances:
 - Atmosphere contains less than 19.5% oxygen.
 - Contaminants or their concentrations are unknown or immediately dangerous to life or health.
 - For protection against gases or vapors.
 - In fire or explosive atmospheres.
9. Do NOT use the respirator with beards or other facial hair that interferes with direct contact between the face and the edge of the respirator, or any other conditions that may prevent a good seal.
10. Discard and replace the respirator if it becomes damaged, breathing resistance becomes excessive or at the end of the shift.
11. Leave the contaminated area immediately if dizziness, irritation or other distress occurs.
12. Never alter, modify or repair this respirator.

Documentation No.	Appendix F Instructions for Use
Creation Date	2022/2/18
Revision No.	1.0
Revision Date	/

Technical File of EC Type Examination Certification

Disposal: Contaminated products should be disposed as hazardous waste in accordance with national regulations.

Storage and Transportation: These products have shelf life of 3 years. End of shelf life is marked on the product packaging. Before initial use, always check that the product is within the stated end of shelf life. Products should be stored in clean, dry conditions within relative humidity of < 80% and temperature range of -30°C to +70°C. When storing or transporting this product use original packaging provided.

Company Name: Lyncmed Medical Technology (Beijing) Co., Ltd.

Address: Room 1601, Building No. 2, Zhubang 2000 Business Building, Balizhuang Xili 99, Chaoyang District, 100022 Beijing, PEOPLE'S REPUBLIC OF CHINA

Tel: +86 01086468058

Notified Body: SGS Fimko Oy

Address: Takomotie 8, FI-00380 Helsinki, Finland Notified Body No. 0598

EC REP: Lyncmed technology SRL

Address: Italy Milano(MI) Via Procaccini Giulio Cesare 32 Cap 20154

Tel: +39 3778578323 **E-mail:** Italy@lync2m.com

EU DECLARATION OF CONFORMITY:

The internet address where the EU declaration of conformity can be accessed is: www.lyncmed.com

Explanation of symbols:

