
## INGEZIM BLV Compac 2.0

R.12.BLV.K3

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

CADUCIDAD: **18 meses**  
Conservado a 4°C

## INGEZIM BLV Compac 2.0 R.12.BLV.K3

**INGEZIM BLV Compac 2.0** is an immunoenzymatic assay based in the blocking ELISA technique, which uses a monoclonal antibody (Mab) specific for Bovine Leukaemia (BLV) gp51 protein.

### TECHNICAL BASIS OF THE KIT

1. Plates are supplied coated with inactivated BLV gp51 protein. Sera samples are added and incubated.
2. If samples contain antibodies specific for gp51 protein of BLV, they will bind the antigen.
3. When Mab-PO specific of gp51 is added, it will bind to the protein only if there are no antibodies of the sample blocking the antigen (negative animals). In case that there are antibodies blocking the antigen (infected animals), conjugate will not be able to bind it. Binding is detected by a colorimetric reaction after the addition of the substrate.

### APPLICATION

1. Detection of antibodies specific for BLV gp51 protein.
2. Sera samples from cattle individually assayed or pooled up to 10.
3. Milk samples individually assayed.

### INTERPRETATION OF THE RESULTS

Two cut off are used for results interpretation: positive and negative.

The samples will be considered as **Positive** (they contain antibodies to BLV), when its OD value at 450nm will be equal or lower than the positive cut off.

The samples will be considered as **Negative** (they do not contain antibodies to BLV), when its OD value at 450 nm will be equal or higher than the negative cut off.

The samples which OD value is between both values will be considered as **Doubtful**.

### SENSITIVITY

Using O.I.E. reference sera

O.I.E. reference sera E-4 and E05 were assayed without prior dilution (as strongly positive serum), at 1/10 dilution in negative serum (as weakly positive serum) and at dilution 1/250 in negative milk (as individual milk sample). Positive results were obtained in all cases, concluding that INGEZIM BLV Compac 2.0 maintains the level of sensitivity required by O.I.E. (Manual Standards for Diagnosis Test and Vaccines 6<sup>th</sup> Edition 2008) for individual and pools of 10 sera and for individual milk.

### SPECIFICITY

In order to determine the specificity of the assay, a set of 354 samples of bovine sera from four different herds, were analyzed. Results obtained indicated a specificity of 100% for this study.

### FIELD TRIALS.

Correspondence between INGEZIM BLV and reference technique, AGID, is 100% when field sera are analyzed.

### COMPOSITION OF THE KIT

- Microtitration plates of 96 well
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

REGISTRATION NUMBER 808 RD  
PRODUCT MANUFACTURED BY INGENASA

SELF LIFE: **18 months**  
Stored at 4°C