



ARCHITECT Free T₃

REF 7K63-27
REF 7K63-37
REF 7K63-32



Free T₃
7K63
G56973R03
B7K6W0

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

ARCHITECT Free T₃

INTENDED USE

The ARCHITECT Free T₃ (FT₃) assay is a Chemiluminescent Microparticle Immunoassay (CMIA) for the quantitative determination of free triiodothyronine (Free T₃) in human serum and plasma.

SUMMARY AND EXPLANATION OF THE TEST

3,5,3' Triiodothyronine (T₃) is a thyroid hormone with a molecular weight of 651 daltons¹ and a half-life in serum of 1.5 days.²

T₃ circulates in the blood as an equilibrium mixture of free and protein bound hormone.³

T₃ is bound to thyroxine binding globulin (TBG), prealbumin, and albumin. The actual distribution of T₃ among these binding proteins is controversial as estimates range from 38-80% for TBG, 9-27% for prealbumin, and 11-35% for albumin.⁴

The binding of these proteins is such that only 0.2-0.4% of the total T₃ is present in solution as unbound or free T₃.⁵

This free fraction represents the physiologically active thyroid hormone.³

Free T₃ is typically elevated to a greater degree than free thyroxine (T₄) in Graves' disease.^{6,7}

Occasionally, free T₃ alone is elevated (T₃ thyrotoxicosis) in about 5% of the hyperthyroid population.⁸

In contrast, levels of free T₄ are elevated to a greater degree than free T₃ in toxic multinodular goiter and excessive T₄ therapy.⁹

Serum free T₃ is useful in distinguishing these forms of hyperthyroidism. Free T₃ may also be important in monitoring patients on anti-thyroid therapy where treatment is focused on reducing the T₃ production and the T₄ conversion to T₃. Serum free T₃ may also be useful in assessing the severity of the thyrotoxic state.

The ARCHITECT Free T₃ assay is to be used as an aid in the assessment of thyroid status.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT Free T₃ assay is a two-step immunoassay to determine the presence of free (unbound) T₃ in human serum and plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

1. Sample and anti-T₃ coated paramagnetic microparticles are combined. Free T₃ (unbound) present in the sample binds to the anti-T₃ coated microparticles.
2. After washing, T₃ acridinium-labeled conjugate is added.
3. Pre-Trigger and Trigger Solutions are then added to the reaction mixture.
4. The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is an inverse relationship between the amount of Free T₃ in the sample and the RLUs detected by the ARCHITECT iSystem optics.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Kit Contents

ARCHITECT Free T₃ 7K63

NOTE: Some kit sizes are not available in all countries or for use on all ARCHITECT iSystems. Please contact your local distributor.

REF	7K63-27	7K63-37	7K63-32
	100	500	2000
MICROPARTICLES	1 x 6.6 mL	1 x 27.0 mL	4 x 27.0 mL
CONJUGATE	1 x 5.9 mL	1 x 26.3 mL	4 x 26.3 mL
MICROPARTICLES	anti-T ₃ (sheep) coated Microparticles in MES buffer with sheep IgG stabilizers. Minimum Concentration: 0.085% solids. Preservative: antimicrobial agent.		
CONJUGATE	T ₃ acridinium-labeled Conjugate in citrate buffer with NaCl and Triton X-100 stabilizers. Minimum concentration: 0.33 ng/mL. Preservative: antimicrobial agent.		

Other Reagents

PRE-TRIGGER SOLUTION ARCHITECT Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

TRIGGER SOLUTION ARCHITECT Trigger Solution containing 0.35 N sodium hydroxide.

WASH BUFFER ARCHITECT Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

NOTE: Bottle and volume vary based on order.


Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁰⁻¹³



The following warnings and precautions apply to: MICROPARTICLES / CONJUGATE	
	
WARNING:	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not use reagent kits beyond the expiration date.
- Do not pool reagents within a kit or between kits.
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- Septums **MUST** be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.
 - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
 - Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
 - Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

When stored and handled as directed, reagents are stable until the expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened/Opened*	2-8°C	Until expiration date	May be used immediately after removal from 2-8°C storage.
On board	System temperature	28 days	Discard after 28 days. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.

* Reagents may be stored on or off the ARCHITECT iSystem. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. **If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.** For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The ARCHITECT Free T₃ assay file must be installed on the ARCHITECT iSystem from an ARCHITECT iSystem Assay CD-ROM prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

Alternate Result Units

Edit assay parameter "Result concentration units" to select an alternate unit.

Conversion formula:

$$(\text{Concentration in Default result unit}) \times (\text{Conversion factor}) = (\text{Concentration in Alternate result unit})$$

Default result unit	Conversion factor	Alternate result unit
pg/mL	1,536	pmol/L

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

Validated specimen types to be used with this assay.

Specimen Types	Collection Tubes
Human serum	Serum
	Serum separator tubes
Human plasma	Sodium heparin
	Lithium heparin
	Potassium EDTA

- Other anticoagulants have not been validated for use with this assay.
- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.
- Performance has not been established for the use of neonatal specimens.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated
- For optimal results, serum and plasma specimens should be free of fibrin, red blood cells or other particulate matter.



- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for specimen collection tubes.
- Specimens must be mixed THOROUGHLY after thawing, by LOW speed vortexing or by gently inverting, and centrifuged prior to use to remove red blood cells or particulate matter to ensure consistency in the results.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

Specimen Storage

Specimen Type	Storage Temperature	Maximum Storage Time
Serum/Plasma	2-8°C	≤ 6 days

If testing will be delayed more than 24 hours, remove serum or plasma from the clot, serum separator or red blood cells.

If testing will be delayed more than 6 days, specimens should be frozen at -10°C or colder.

Specimens stored frozen at -10°C or colder for 6 days showed no performance difference.

Avoid multiple freeze/thaw cycles.

Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.

PROCEDURE

Materials Provided

7K63 ARCHITECT Free T₃ Reagent Kit

Materials Required but not Provided

- ARCHITECT Free T₃ Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on www.abbottdiagnostics.com.
- 7K63-02 ARCHITECT Free T₃ Callibrators
- 7K63-12 ARCHITECT Free T₃ Controls
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
 - **Invert the microparticle bottle 30 times.**
 - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.

- **If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
- Once the microparticles have been resuspended, discard the cap and place a septum on the bottle. For instructions about placing septums on bottles, refer to the **Reagent Handling** section of this package insert.

- Load the reagent kit on the ARCHITECT iSystem.
 - Verify that all necessary reagents are present.
 - Ensure that septums are present on all reagent bottles.
- Order tests.
 - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

Maximum number of replicates sampled from the same sample cup: 10

- Priority:
 - Sample volume for first test: 72 µL
 - Sample volume for each additional test from same sample cup: 22 µL
- ≤ 3 hours on board:
 - Sample volume for first test: 150 µL
 - Sample volume for each additional test from same sample cup: 22 µL
- > 3 hours on board: additional sample volume is required. Refer to the ARCHITECT System Operations Manual, Section 5 for information on sample evaporation and volumes.
- If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT Free T₃ Calibrators and Controls.
 - Mix calibrator(s) and controls by gentle inversion before use.
 - Hold bottles vertically and dispense recommended volumes into each respective sample cup.
 - Recommended volumes:
 - for each calibrator: 4 drops
 - for each control: 4 drops
- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Specimen Dilution Procedures

Specimens cannot be diluted for Free T₃ determinations. Specimens which read > 20.00 pg/mL should be reported as such.

Calibration

- Test Calibrators A to F in duplicate. The calibrators should be priority loaded.
 - A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Calibration Range: 0.0 - 30.0 pg/mL.
- Once an ARCHITECT Free T₃ calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used or
 - Controls are out of range.
- For detailed information on how to perform an assay calibration refer to the ARCHITECT System Operations Manual, Section 6.



Quality Control Procedures

The recommended control requirement for the ARCHITECT Free T₃ assay is that a single sample of each control level be tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.

Ensure that assay control values are within the concentration ranges specified in the control package insert.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.

The ARCHITECT Free T₃ assay belongs to method group 2.

RESULTS

Calculation

The ARCHITECT Free T₃ assay utilizes a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration curve.

For information on alternate result units, refer to the INSTRUMENT PROCEDURE, Alternate Result Units section of this package insert.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

Measuring Interval

Measuring interval is defined as the range of values in pg/mL which meets the limits of acceptable performance for both imprecision and linearity. The measuring interval for the ARCHITECT Free T₃ assay is 1.5 (Limit of Quantitation - LoQ) to 20 pg/mL.

LIMITATIONS OF THE PROCEDURE

- Results should be used in conjunction with other data; e.g., symptoms, results of other thyroid tests, clinical impressions, etc.
- If the Free T₃ results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.

EXPECTED VALUES

A normal range of 1.88-3.18 pg/mL (Central 95% interval) was obtained by testing serum specimens from 260 individuals determined as normal by ARCHITECT Anti-Tg, Anti TPO and TSH assays. It is recommended that each laboratory establish its own normal range, which may be unique to the population it serves depending upon geographical, patient, dietary, or environmental factors.

Free T₃ is a secondary indicator of thyroid status. Although the majority of patients with hyperthyroidism will have free T₃ values greater than the upper limit of the euthyroid range, some may have free T₃ values which fall within the normal range.^{14, 15} Specimens from patients described as "sick euthyroids" generally yield values in the low to normal range.^{16, 17}

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The ARCHITECT Free T₃ assay is designed to have a precision of ≤ 10% (total CV). A study based on guidance from Clinical and Laboratory Standards Institute (CLSI, formerly NCCLS) document EP5-A¹⁸ was performed for the ARCHITECT Free T₃ assay. A three member processed human serum based panel was assayed, using two lots of reagents, in replicates of two at two separate times per day for 20 testing days. Data from this study are summarized in the following table.*

Panel Member	Reagent Lot	Instrument	n	Mean Conc.			Within Run		Total	
				Value (pg/mL)	SD	%CV	SD	%CV	SD	%CV
1	1	1	80	3.22	0.096	3.0	0.115	3.6		
1	1	2	80	3.16	0.133	4.2	0.143	4.5		
1	2	1	80	3.60	0.108	3.0	0.141	3.9		
1	2	2	80	3.35	0.113	3.4	0.131	3.9		

Panel Member	Reagent Lot	Instrument	n	Mean Conc. Value (pg/mL)	Within Run			Total	
					SD	%CV	SD	%CV	
2	1	1	80	6.00	0.099	1.7	0.168	2.8	
2	1	2	80	5.88	0.166	2.8	0.184	3.1	
2	2	1	80	6.28	0.154	2.5	0.176	2.8	
2	2	2	80	6.06	0.194	3.2	0.225	3.7	
3	1	1	80	10.50	0.252	2.4	0.481	4.6	
3	1	2	80	10.01	0.289	2.9	0.496	5.0	
3	2	1	80	10.50	0.145	1.4	0.237	2.3	
3	2	2	80	10.12	0.217	2.1	0.265	2.6	

* Representative data; results in individual laboratories may vary from these data.

Sensitivity

The ARCHITECT Free T₃ assay is designed to have a Limit of Quantitation (LoQ) of ≤ 1.5 pg/mL. The LoQ is defined as the lowest concentration at which analyte in a sample can be accurately quantified with precision of ≤ 20% CV.

A study was performed based on guidance from the CLSI document EP17-A2²⁰ with four zero-level samples and 8 samples with T₃ target concentrations ranging from 1.0 to 3.2 pg/mL. The samples were tested over a minimum of 3 days using 2 reagent lots and 6 instruments. In this study, the Limit of Blank (LoB) was 0.94 pg/mL, Limit of Detection (LoD) was 1.07 pg/mL and LoQ was 1.25 pg/mL.*

* Representative data; results in individual laboratories may vary from these data.

Analytical Specificity

The ARCHITECT Free T₃ assay is designed to have a mean analytical specificity of ≤ 0.001% cross reactivity with thyroxine (T₄) at a concentration of 1,000,000 pg/mL.

Interference

The ARCHITECT Free T₃ assay is designed to have a mean potential interference from hemoglobin, bilirubin, triglycerides, and protein of < 10% at the levels indicated below.

Hemoglobin	≤ 500 mg/dL
Bilirubin	≤ 20 mg/dL
Triglycerides	≤ 2000 mg/dL
Protein	≤ 12 g/dL

Accuracy by Correlation

The ARCHITECT Free T₃ 6 point assay is designed to have a slope of 1.00 +/- 0.15 and a correlation coefficient (r) of ≥ 0.90 when compared to the ARCHITECT Free T₃ 2 point assay. A study was performed where specimens were tested using the ARCHITECT Free T₃ 6 point assay and ARCHITECT Free T₃ 2 point assay. Data from this study were analyzed using least squares and Passing Bablok¹⁹ regression methods and are summarized in the following table.*

Abbott ARCHITECT Free T₃ 6 point assay vs. Abbott ARCHITECT Free T₃ 2 point assay

Method	Number of Specimens	Intercept	Slope	Correlation Coefficient
Least Squares Linear Regression	144	-0.32	1.04	0.99
Passing-Bablok Linear Regression**	144	-0.00	0.97	0.99

In this evaluation, serum specimens tested ranged from 1.50 pg/mL to 17.81 pg/mL with the ARCHITECT Free T₃ 6 point assay and from 1.55 pg/mL to 17.66 pg/mL with the ARCHITECT Free T₃ 2 point assay.

* Representative data; variables such as differences in sampling size and sample population may impact the correlation of the assay; therefore, results in individual laboratories may vary from these data.

** A linear regression method with no special assumptions regarding the distribution of the samples and the measurement errors.¹⁹



BIBLIOGRAPHY

- Budavari S, editor. *Merck Index* (11th Ed.). Rahway, NJ: Merck and Co., Inc., 1989:866.
- Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972;21:1073-1092.
- Ekins RP, editor. *Methods for the Measurement of Free Thyroid Hormones*. Amsterdam: Excerpta Medica Foundation. 1979;72-92.
- Robbins J, Rall JE. The Iodine-Containing Hormones. In: *Hormones in Blood* (3rd Ed.). London: Academic Press, 1979;1:632-667.
- DeGroot LJ, Larsen PR, Refetoff S, Stanbury JB. Transport of Thyroid Hormone and Cell Uptake. In: *The Thyroid and Its Diseases*. New York: Wiley and Sons, 1984;62-66.
- Hamburger JL. Evolution of Toxicity in Solitary Nontoxic Autonomously Functioning Thyroid Nodules. *J Clin Endocrinol Metab* 1980;50: 1089-1093.
- Ladenson PW. Diagnosis of Thyrotoxicosis. In: Braverman LE, Utiger RD, editors. *The Thyroid* (6th Ed.). Philadelphia: JB Lippincott Co., 1991:880-886.
- Wahner HW. T3 Hyperthyroidism. *Mayo Clin Proc* 1972;47:938-943.
- Lum SM, Nicoloff JT, Spencer CA, Kaptein EM. Peripheral Tissue Mechanism for Maintenance of Serum Triiodothyronine Values in a Thyroxine-Deficient State in Man. *J Clin Invest* 1984;73:570-575.
- US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Lee LA, Mooney RA, Woolf PD. Clinical Utility of Measuring Free Thyroxin and Free Triiodothyronine in Serum of Critically Ill Patients by Ultrafiltration. *Clin Chem* 1986; 32:797-800.
- Wilke TJ, Eastment HT. Discriminative Ability of Tests for Free and Total Thyroid Hormones in Diagnosing Thyroid Disease. *Clin Chem* 1986; 32:1746:50.
- Kaptein EM, Grieb DA, Spencer CA, et al. Thyroxine Metabolism in the Low Thyroxine State of Critical Nonthyroidal Illnesses. *J Clin Endocrinol Metab* 1981; 53:764-771.
- Merimee TJ, Fineberg ES. Starvation-Induced Alterations of Circulating Thyroid Hormone Concentrations in Man. *Metabolism* 1976; 25:79-83.
- National Committee for Clinical Laboratory Standards. *Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline*. NCCLS Document EPS-A. Wayne, PA: NCCLS, 2001.
- Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. Application of linear regression procedures for method comparison studies in clinical chemistry. Part I. *J Clin Chem Clin Biochem* 1983;21 (11):709-720.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
CONJUGATE	Conjugate
CONTROL NO.	Control Number
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
MICROPARTICLES	Microparticles
PRE-TRIGGER SOLUTION	Pre-Trigger Solution
PRODUCT OF IRELAND	Product of Ireland
REACTION VESSELS	Reaction Vessels
REAGENT LOT	Reagent Lot
REF	List Number
REPLACEMENT CAPS	Replacement Caps
SAMPLE CUPS	Sample Cups
SEPTUM	Septum
SN	Serial number
TRIGGER SOLUTION	Trigger Solution
WARNING: SENSITIZER	Warning: May cause an allergic reaction.
WASH BUFFER	Wash Buffer

ARCHITECT and Chemiflex are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000



DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA.

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised March 2017.

©2017 Abbott Laboratories





ARCHITECT Free T₃ Controls

REF 7K63-10



en

Free T₃
7K63
G5-9120/R05
C7K630

Read Highlighted Changes: Revised June 2015.

INTENDED USE

The ARCHITECT Free T₃ Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT iSystem (reagents, calibrators and instrument), when used for the quantitative measurement of free triiodothyronine (Free T₃) in human serum and plasma. Refer to the ARCHITECT assay-specific reagent package insert for additional information.

CONTENTS

3 Bottles (8 mL each) of ARCHITECT Free T₃ Controls contain T₃ in human serum. Preservative: Sodium Azide.

The following concentration ranges may be used for individual replicate control specifications on the ARCHITECT iSystem:

Control	Target	Range	Target	Range
	Concentration		Concentration	
	(pg/mL)		(pmol/L)	
CONTROL L	3.1	2.02 - 4.09	4.8	3.1 - 6.29
CONTROL M	6.0	4.20 - 7.80	9.2	6.45 - 11.98
CONTROL H	10.5	7.88 - 14.18	16.1	12.10 - 21.77


Each laboratory should establish its own concentration ranges for new control lots at each control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days. Sources of variation that can be expected should be included in this study in order to be representative of future system performance. These may include:

- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules
- Data points collected at different times of the day

These results should be applied to your laboratory's quality control practices.

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.
- The human serum used in the ARCHITECT Free T₃ Controls is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HCV, and anti-HIV-1/HIV-2.


The following warnings and precautions apply to: CONTROL L /	
CONTROL L	CONTROL H
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

STORAGE

- Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.



BIBLIOGRAPHY


1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization, *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI), *Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.



Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONC	Concentration
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL L	Control Low, Medium, High (L,M,H)
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LDT	Lot Number
PRODUCT OF IRELAND	Product of Ireland
RANGE	Range
REF	List Number

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions.

 Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000



DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised June 2015.
©2005, 2015 Abbott Laboratories



ARCHITECT Free T₃ Calibrators

REF 7K63-02



en

Free T₃
7K63
G6-7702 / R01
S7K6W0

INTENDED USE

The ARCHITECT Free T₃ Calibrators are for the calibration of the ARCHITECT iSystem when used for the quantitative determination of free triiodothyronine (Free T₃) in human serum and plasma. Refer to the ARCHITECT assay-specific reagent package insert for additional information.

CONTENTS

6 Bottles (4 mL each) of ARCHITECT Free T₃ Calibrators prepared in human serum. Preservative: Sodium Azide.

The calibrators yield the following concentrations:

Callibrator	Free T ₃ Concentration	
	(pg/mL)	(pmol/L)
CAL A	0.0	0.00
CAL B	1.4	2.15
CAL C	3.5	5.38
CAL D	7.0	10.75
CAL E	17.2	26.42
CAL F	30.0	46.08

STANDARDIZATION

The calibrators are matched to an Abbott internal reference standard. This internal reference standard is manufactured by gravimetric methods based on the Free Triiodothyronine calculation (FT_{3c}) using L-Triiodothyronine (sodium salt, not less than 95% pure by HPLC) and L-Thyroxine (sodium pentahydrate, not less than 95% pure by HPLC) at each concentration level. The FT_{3c} is a calculation of the Free Triiodothyronine hormone concentration, which depends on the amount of Total T₃ and Total T₄ found in the serum as well as the serum's thyroid hormone binding capacity.

PRECAUTIONS

- IVD
- For *In Vitro* Diagnostic Use
- Rx ONLY



- **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human serum used in the ARCHITECT Free T₃ Calibrators is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HCV, and anti-HIV-1/HIV-2.


The following warnings and precautions apply to: CAL A - CAL F	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

STORAGE

- Calibrators are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.



PREPARATION FOR USE

Calibrators may be used immediately after removal from 2-8°C storage.

Prior to each use, mix by gentle inversion.

After each use, tightly close the caps and return the calibrators to 2-8°C storage.

BIBLIOGRAPHY


1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.



Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CAL A	Calibrator (A,B,C,D,E or F)
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF IRELAND	Product of Ireland
REF	List Number
Re ONLY	For use by or on the order of a physician only (applicable to USA classification only).

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions.

 Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000



DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Created September 2016.
©2016 Abbott Laboratories





ARCHITECT Free T₄ Controls

REF 7K65-10



Free T₄
7K65
G5-9239 / R04
C7K650

Read Highlighted Changes: Revised June 2015.

INTENDED USE

The ARCHITECT Free T₄ Controls are for the verification of the accuracy and precision of the ARCHITECT iSystem when used for the quantitative determination of free thyroxine (Free T₄) in human serum and plasma. Refer to the ARCHITECT assay-specific reagent package insert for additional information.

CONTENTS

3 Bottles (8 mL each) of ARCHITECT Free T₄ Controls contain T₄ prepared in human serum. Preservative: sodium azide. The following concentration ranges may be used for individual replicate control specifications on the ARCHITECT iSystem:

Control	Target Concentration		Target Concentration	
	ng/dL	Range	pmol/L	Range
CONTROL L	0.65	0.42 - 0.85	8.4	5.41 - 10.94
CONTROL M	1.2	0.86 - 1.62	15.4	11.07 - 20.85
CONTROL H	2.8	1.82 - 3.78	36.0	23.42 - 48.65

Each laboratory should establish its own concentration ranges for new control lots at each control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days. Sources of variation that can be expected should be included in this study in order to be representative of future system performance. These may include:

- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules
- Data points collected at different times of the day


These results should be applied to your laboratory's quality control practices.

STANDARDIZATION

The controls are matched to an Abbott internal reference standard. This internal reference standard is manufactured by gravimetric methods based on the Free Thyroxine calculation (FT_{4c}) using L-Thyroxine, sodium salt pentahydrate, (HPLC grade), at each concentration level. The FT_{4c} is a calculation of the free thyroid hormone concentration, which depends on the amount of Total T₄ found in the serum and the serum's T₄ binding capacity.

PRECAUTIONS

- IVD
- For *In Vitro* Diagnostic Use

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.
- The human serum used in the ARCHITECT Free T₄ Controls is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HCV and anti-HIV-1/HIV-2.

The following warnings and precautions apply to: CONTROL L / CONTROL M / CONTROL H	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

STORAGE

- Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.



BIBLIOGRAPHY


1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.



Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONC.	Concentration
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL L	Control Low, Medium, High (L,M,H)
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF IRELAND	Product of Ireland
RANGE	Range
REF	List Number

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions.

 Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000



DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised June 2015.
©2006, 2015 Abbott Laboratories





ARCHITECT Free T₄ Calibrators

REF 7K65-02



en

Free T₄
7K65
G5-9241 / R02
S7K6U0

Read Highlighted Changes: Revised June 2015.

INTENDED USE

The ARCHITECT Free T₄ Calibrators are for the calibration of the ARCHITECT iSystem when used for the quantitative determination of free thyroxine (Free T₄) in human serum and plasma. Refer to the ARCHITECT assay-specific reagent package insert for additional information.

CONTENTS

6 Bottles (4 mL each) of ARCHITECT Free T₄ Calibrators prepared in human serum. Preservative: sodium azide.

The calibrators yield the following concentrations:

Callibrators	Free T ₄ Concentration	
	(ng/dL)	(pmol/L)
CAL A	0.0	0.0
CAL B	0.5	6.4
CAL C	1.0	12.9
CAL D	2.0	25.7
CAL E	3.5	45.0
CAL F	6.0	77.2

STANDARDIZATION

The calibrators are matched to an Abbott internal reference standard. This internal reference standard is manufactured by gravimetric methods based on the Free Thyroxine calculation (FT_{4c}) using L-Thyroxine, sodium salt pentahydrate (HPLC grade), at each concentration level. The FT_{4c} is a calculation of the free thyroid hormone concentration, which depends on the amount of Total T₄ found in the serum and the serum's T₄ binding capacity.

PRECAUTIONS

• IVD

• For *In Vitro* Diagnostic Use



- **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.
- The human serum used in the ARCHITECT Free T₄ Calibrators is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HCV and anti-HIV-1/HIV-2.

The following warnings and precautions apply to: CAL A CAL F	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

STORAGE

- Calibrators are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.



BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Calibrator (A,B,C,D,E or F)
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Distributed in the USA by
	Information needed for United States of America only
	<i>In Vitro</i> Diagnostic Medical Device
	Lot Number
	Product of Ireland
	List Number

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions.





Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000



DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

**Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com**

Revised June 2015.
©2013, 2015 Abbott Laboratories





ARCHITECT
Free T₄

REF 7K65-29
REF 7K65-24
REF 7K65-39
REF 7K65-34



en

Free T₄
7K65
G59250R03
B7K6F0

Read Highlighted Changes: Revised April 2017.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

ARCHITECT Free T₄

INTENDED USE

The ARCHITECT Free T₄ (FT4) assay is a Chemiluminescent Microparticle Immunoassay (CMIA) for the quantitative determination of free thyroxine (Free T₄) in human serum and plasma.

SUMMARY AND EXPLANATION OF THE TEST

Thyroxine (T₄) circulates in the blood as an equilibrium mixture of free and serum protein bound hormones. Thyroxine binding globulin (TBG), albumin and pre-albumin bind approximately 75%, 10% and 15% of the total circulating T₄ respectively.¹⁻³ The binding of T₄ by these proteins is such that less than 0.03% is present in the circulation as unbound, free T₄.⁴ This small percentage of the total T₄ represents the physiologically available hormone which is biologically active. Once the free T₄ is absorbed by the target cells, the equilibrium reestablishes circulating free T₄ levels. The equilibrium results in the maintenance of a constant level of free T₄ when alterations occur in either the concentration or affinity of the serum binding proteins. Therefore, in a variety of normal (pregnancy)⁴ and abnormal (Familial Dysalbuminemic Hyperthyroxinemia, FDH)⁵⁻⁷ states, or as a result of the administration of certain drugs (e.g. furosemide^{8, 9} and fenclofenac¹⁰⁻¹²), the target tissues are assured of receiving the required amount of hormone. Free T₄ values may, therefore, provide the best indication of thyroid dysfunction, since free T₄ is less sensitive to changes in the serum binding proteins.

Historically, the diagnosis of thyroid function has involved performing a total T₄ assay^{13, 14} in addition¹⁵ to a Thyroxine Uptake (TU) assay of the same sample. The mathematical combination of these two assays produces a Free Thyroxine Index (FTI) which provides an indirect proportional estimate for free T₄.¹⁶

Alternatively, direct assays have been developed using equilibrium dialysis,^{17, 18} ultrafiltration,^{19, 20} RIA,²¹ and solid-phase EIA technology²² to measure free T₄. In these methods, separation of free and bound tracer is achieved either with a membrane, or by binding free T₄ to a solid phase antibody. This extraction step removes an amount of T₄ which is proportional to the original amount of free T₄ present in the patient sample. Provided that the extracted T₄ is less than approximately 5% of the T₄ in the sample, a true estimation of the free T₄ content can be obtained.

The ARCHITECT Free T₄ assay is to be used as an aid in the assessment of thyroid status.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT Free T₄ assay is a two-step immunoassay to determine the presence of free thyroxine (Free T₄) in human serum and plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

1. Sample and anti-T₄ coated paramagnetic microparticles are combined. Free T₄ (unbound) present in the sample binds to the anti-T₄ coated microparticles.
2. After washing, T₃ acridinium-labeled conjugate is added to create a reaction mixture.
3. Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
4. The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is an inverse relationship between the amount of Free T₄ in the sample and the RLUs detected by the ARCHITECT iSystem optics.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Kit Contents

ARCHITECT Free T₄ 7K65

NOTE: Some kit sizes are not available in all countries or for use on all ARCHITECT iSystems. Please contact your local distributor.

REF	7K65-29	7K65-24	7K65-39	7K65-34
	100	400	500	2000
MICROPARTICLES	1 x 6.6 mL	4 x 6.6 mL	1 x 27.0 mL	4 x 27.0 mL
CONJUGATE	1 x 5.9 mL	4 x 5.9 mL	1 x 26.3 mL	4 x 26.3 mL
MICROPARTICLES	anti-T ₄ (sheep) coated Microparticles in TRIS buffer with sheep IgG stabilizers. Minimum concentration: 0.05% solids. Preservative: sodium azide.			
CONJUGATE	T ₃ acridinium-labeled Conjugate in MES buffer with NaCl and Triton X-100 stabilizers. Minimum concentration: 0.2 ng/mL. Preservative: ProClin.			

Other Reagents

PRE-TRIGGER SOLUTION ARCHITECT Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

TRIGGER SOLUTION ARCHITECT Trigger Solution containing 0.35 N sodium hydroxide.

WASH BUFFER ARCHITECT Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

NOTE: Bottle and volume varies based on order.


Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use



Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.²³⁻²⁶

The following warnings and precautions apply to: CONJUGATE	
	
WARNING:	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.
The following warnings and precautions apply to: MICROPARTICLES	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not use reagent kits beyond the expiration date.
- **Do not pool reagents within a kit or between kits.**
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- **Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the Instructions in this package insert.**
 - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
 - Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
 - Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

When stored and handled as directed, reagents are stable until the expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened/ Opened*	2-8°C	Until expiration date	May be used immediately after removal from 2-8°C storage. Store in upright position.
On board	System temperature	30 days	Discard after 30 days. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.

* Reagents may be stored on or off the ARCHITECT iSystem. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. If the **microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.** For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The ARCHITECT Free T₄ assay file must be installed on the ARCHITECT iSystem from an ARCHITECT iSystem Assay CD-ROM prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

Alternate Result Units

Edit assay parameter "Result concentration units" to select an alternate unit.

Conversion formula:

$$(\text{Concentration in Default result unit}) \times (\text{Conversion factor}) = (\text{Concentration in Alternate result unit})$$

Default result unit	Conversion factor	Alternate result unit
ng/dL	12.87	pmol/L



SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

Validated specimen types to be used with this assay:

Specimen Types	Collection Tubes
Human serum	Serum Serum separator tubes
Human plasma	Sodium heparin Lithium heparin Lithium heparin plasma separator tubes Potassium EDTA

- Other anticoagulants have not been validated for use with this assay.
- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.
- Performance of this test has not been established with neonatal specimens.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.
- For optimal results, serum and plasma specimens should be free of fibrin, red blood cells or other particulate matter.
- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for specimen collection tubes.
- Specimens must be mixed THOROUGHLY after thawing, by LOW speed vortexing or by gently inverting, and centrifuged prior to use to remove red blood cells or particulate matter to ensure consistency in the results.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

Specimen Storage

Specimen Type	Storage Temperature	Maximum Storage Time
Serum/Plasma	2-8°C	≤ 6 days

If testing will be delayed more than 24 hours, remove serum or plasma from the clot, serum separator, plasma separator or red blood cells.

Follow the manufacturer's processing instructions for serum or plasma collection tubes if a removal time of less than 24 hours is specified.

If testing will be delayed more than 6 days specimens should be frozen at -10°C or colder.

Specimens stored frozen at -10°C or colder for 6 days showed no performance difference.

Avoid multiple freeze/thaw cycles.

Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.

PROCEDURE

Materials Provided

7K65 ARCHITECT Free T₄ Reagent Kit

Materials Required but not Provided

- ARCHITECT Free T₄ Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on www.abbottdiagnostics.com.
- 7K65-02 ARCHITECT Free T₄ Calibrators
- 7K65-10 ARCHITECT Free T₄ Controls
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
 - Invert the microparticle bottle 30 times.**
 - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
 - If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
 - Once the microparticles have been resuspended, place a septum on the bottle. For instructions about placing septums on bottles, refer to the **Reagent Handling** section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
 - Verify that all necessary reagents are present.
 - Ensure that septums are present on all reagent bottles.
- Order tests.
 - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
 - Maximum number of replicates sampled from the same sample cup: 10
 - Priority:
 - Sample volume for first test: 95 µL
 - Sample volume for each additional test from same sample cup: 45 µL
 - ≤ 3 hours on board:
 - Sample volume for first test: 150 µL
 - Sample volume for each additional test from same sample cup: 45 µL
 - > 3 hours on board: Additional sample volume is required. Refer to the ARCHITECT System Operations Manual, Section 5 for information on sample evaporation and volumes.
 - If using primary or aliquot tubes, use the **sample volume** to ensure sufficient patient specimen is present.



- Prepare ARCHITECT Free T₄ Calibrators and Controls.
 - Mix calibrator(s) and controls by gentle inversion before use.
 - Hold bottles **vertically** and dispense recommended volumes into each respective sample cup.
 - Recommended volumes:
 - for each calibrator: 4 drops
 - for each control: 4 drops
- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Specimen Dilution Procedures

Samples cannot be diluted for Free T₄ determinations. Samples which read > 5.00 ng/dL should be reported as such.

Calibration

- Test Calibrators A-F in duplicate. The calibrators should be priority loaded.

A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Calibration Range: 0.0 - 6.0 ng/dL.
- Once an ARCHITECT Free T₄ calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used or
 - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Quality Control Procedures

The recommended control requirement for the ARCHITECT Free T₄ assay is that a single sample of each control level be tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.

Ensure that assay control values are within the concentration ranges specified in the control package insert.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.

The ARCHITECT Free T₄ assay belongs to method group 6.

RESULTS

Calculation

The ARCHITECT Free T₄ assay utilizes a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration curve.

For information on alternate result units, refer to the **INSTRUMENT PROCEDURE, Alternate Result Units** section of this package insert.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

Measuring Interval

Measuring interval is defined as the range of values in ng/dL which meets the limits of acceptable performance for both imprecision and linearity.

The measuring interval for the ARCHITECT Free T₄ assay is 0.40 (Limit of Quantitation - LoQ) to 5.00 ng/dL.

LIMITATIONS OF THE PROCEDURE

- Results should be used in conjunction with other data; e.g., symptoms, results of other thyroid tests, clinical impressions, etc.
- If the Free T₄ results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.

EXPECTED VALUES

A normal range of 0.70 ng/dL to 1.48 ng/dL (central 99% interval) was obtained by testing serum specimens from 411 individuals determined as normal by AxSYM Ultrasensitive hTSH II and AxSYM Free T₄ assays. It is recommended that each laboratory establish its own normal range, which may be unique to the population it serves depending upon geographical, patient, dietary, or environmental factors.

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The ARCHITECT Free T₄ assay is designed to have a precision of ≤ 10% (total CV) for concentrations in the range of the low control (0.65 ng/dL), medium control (1.2 ng/dL) and high control (2.8 ng/dL). A study based on guidance from Clinical and Laboratory Standards Institute (CLSI, formerly NCCLS) document EP5-A²⁷ was performed for the ARCHITECT Free T₄ assay. A three member processed human serum based panel was assayed, using two lots of reagents, in replicates of two at two separate times per day for 20 testing days. Data from this study are summarized in the following table.*

Panel Member	Reagent Lot	Instrument	n	Mean Conc. Value (ng/dL)	Within Run		Total	
					SD	%CV	SD	%CV
1	1	1	80	0.69	0.021	3.0	0.032	4.7
	1	2	80	0.67	0.036	5.3	0.041	6.1
1	2	1	80	0.70	0.021	3.0	0.055	7.8
	2	2	80	0.72	0.027	3.7	0.043	6.0
2	1	1	80	1.26	0.048	3.8	0.061	4.8
2	1	2	80	1.22	0.029	2.3	0.044	3.6
	2	1	80	1.25	0.029	2.3	0.066	5.2
2	2	2	80	1.27	0.033	2.6	0.048	3.8
3	1	1	80	2.94	0.084	2.8	0.148	5.1
3	1	2	80	2.87	0.097	3.4	0.151	5.3
	2	1	80	3.03	0.098	3.3	0.191	6.3
3	2	2	80	3.00	0.088	2.9	0.134	4.5

* Representative data; results in individual laboratories may vary from these data.

Sensitivity

The ARCHITECT Free T₄ assay is designed to have a Limit of Quantitation (LoQ) of ≤ 0.4 ng/dL. The LoQ is defined as the lowest amount of analyte in a sample that can be accurately quantitated with precision of ≤ 10% CV.

A study was performed based on guidance from the NCCLS document EP17-A²⁹ with four zero-level samples and 8 samples with Free T₄ concentrations ranging from 0.25 to 1.0 ng/dL.

The samples were tested in at least 5 separate runs over a minimum of 3 days using 2 reagent lots and 6 instruments. In this study, the Limit of Blank (LoB) was 0.22 ng/dL, Limit of Detection (LoD) was 0.28 ng/dL and LoQ was 0.40 ng/dL.*

* Representative data; results in individual laboratories may vary from these data.

Analytical Specificity

The ARCHITECT Free T₄ assay is designed to have a mean analytical specificity of ≤ 0.0035% cross reactivity with triiodothyronine (T₃) at a concentration of 12,000 ng/dL in a sample containing 0.5 ng/dL of Free T₄.



Interference

The ARCHITECT Free T₄ assay is designed to have a mean potential interference from hemoglobin, bilirubin, triglycerides, and protein of < 10% at the levels indicated below.

Hemoglobin	≤ 500 mg/dL
Bilirubin	≤ 20 mg/dL
Triglycerides	≤ 3000 mg/dL
Protein	≤ 12 g/dL

Accuracy by Correlation

The ARCHITECT Free T₄ assay is designed to have a slope of 1.00 ± 0.20 and a correlation coefficient (r) of ≥ 0.90 when compared to the AxSYM Free T₄ assay.

A study was performed where specimens were tested using the ARCHITECT Free T₄ assay and AxSYM Free T₄ assay. Data from this study were analyzed using least squares and Passing Bablok²⁸ regression methods and are summarized in the following table.*

Abbott ARCHITECT Free T₄ vs. Abbott AxSYM Free T₄

Method	Number of Specimens	Intercept	Slope	Correlation Coefficient
Least Squares Linear Regression	675	0.03	0.96	0.953
Passing-Bablok Linear Regression**	675	-0.02	1.00	0.953

In this evaluation, serum specimens tested ranged from 0.52 ng/dL to 3.88 ng/dL with the ARCHITECT Free T₄ assay and from 0.46 ng/dL to 4.14 ng/dL with the AxSYM Free T₄ assay.

* Representative data; variables such as differences in sampling size and sample population may impact the correlation of the assay; therefore, results in individual laboratories may vary from these data.

** A linear regression method with no special assumptions regarding the distribution of the samples and the measurement errors.²⁸

BIBLIOGRAPHY

- Oppenheimer JH. Role of Plasma Proteins in the Binding, Distribution and Metabolism of the Thyroid Hormones. *N Engl J Med* 1968; 278:1153-1162.
- Woeber KA, Ingbar SH. The Contribution of Thyroxine-Binding Prealbumin to the Binding of Thyroxine in Human Serum, as Assessed by Immunoabsorption. *J Clin Invest* 1968;47:1710-1721.
- Tabachnick M, Giorgio NA Jr. *Thyroxine - Protein Interactions Arch Biochem Biophys* 1964;105:563-569.
- DæGroot LJ, Larsen PR, Refetoff S, Stanbury JB. Transport of Thyroid Hormone and Cell Uptake. In: *The Thyroid and Its Diseases*. New York: Wiley and Sons, 1984;62-65.
- Stockigt JR, Topliss DJ, Barlow JW, White EL, Hurley DM, Taft P. Familial Euthyroid Thyroxine Excess: An Appropriate Response to Abnormal Thyroxine Binding Associated with Albumin. *J Clin Endocrinol Metab* 1981;53:353-359.
- Ruiz M, Rajatanavin R, Young RA, Taylor C, Brown R, Braverman LE, et al. Familial Dysalbuminemic Hyperthyroxinemia. *N Engl J Med* 1982;306:635-639.
- Hennemann G, Docter R, Krenning EP, Bos G, Otten M, Vissers TJ. Raised Total Thyroxine and Free Thyroxine Index but Normal Free Thyroxine. *Lancet* 1979;1:639-642.
- Stockigt JR, Lim CF, Barlow JW, Stevens V, Topliss DJ, Wynne KN. High Concentrations of Furosemide Inhibit Serum Binding of Thyroxine. *J Clin Endocrinol Metab* 1984;59:62-66.
- Stockigt JR, Lim CF, Barlow JW, Stevens V, Wynne KN, Mohr VS, Topliss DJ, et al. Interaction of Furosemide with Serum Thyroxine-Binding Sites: In Vivo and In Vitro Studies and Comparison with Other Inhibitors. *J Clin Endocrinol Metab* 1985;60:1025-1031.
- Ratcliffe WA, Hazelton RA, Thomson JA, Ratcliffe JG. The Effect of Fenclufenac On Thyroid Function Tests In Vivo and In Vitro. *Clin Endocrinol* 1980;13:569-575.
- Kurtz AB, Capper SJ, Clifford J, Humphrey MJ, Lukinac L. The Effect of Fenclufenac On Thyroid Function. *Clin Endocrinol* 1981;15:117-124.
- Pearson DWM, Ratcliffe WA, Thomson JA, Ratcliffe JG. Biochemical and Clinical Effects of Fenclufenac In Thyrotoxicosis. *Clin Endocrinol* 1982;16:369-373.
- Larsen PR, Dockalova J, Sipula D, Wu FM. Immunoassay of Thyroxine in Unextracted Human Serum. *J Clin Endocrinol Metab* 1973; 37:177-182.
- Chopra IJ. A Radioimmunoassay for Measurement of Thyroxine in Unextracted Human Serum. *J Clin Endocrinol Metab* 1972; 34:938-947.
- Refetoff S, Hayen SR, Selenkow HA. Estimation of the T₄ Binding Capacity of Serum TBG and TBPA by a Single T₄ Load Ion Exchange Resin Method. *J Nucl Med* 1972;13:2-12.
- Clark F, Horn DB. Assessment of Thyroid function by the Combined Use of the Serum Protein-Bound Iodine and Resin Uptake of 131 I-Triiodothyronine. *J Clin Endocrinol Metab* 1965;25:39-45.
- Ellis SM, Ekins RP. The Radioimmunoassay of Serum Free Triiodothyronine and Thyroxine. In: Pasternak CA, editor. *Radioimmunoassay in Clinical Biochemistry*. London: Heyden, 1975;25:187-194.
- Nelson JC, Tomei RT. Direct Determination of Free Thyroxine in Undiluted Serum by Equilibrium Dialysis/Radioimmunoassay in Serum. *Clin Chem* 1988;34:1737-1744.
- Sophianopoulos J, Jerkunica I, Lee CN, Sgoutas D. An Improved Ultrafiltration Method for Free Thyroxine and Triiodothyronine in Serum. *Clin Chem* 1980;26:159-162.
- Shannon N, Woolf PD. Determination of Free Thyroxine in serum by Ultrafiltration: Validation of a Method and Preliminary Results. *Clin Chem* 1984;30:1770-1773.
- Bayer MF, McDougall IR. Radioimmunoassay of Free Thyroxine in Serum: Comparison with Clinical Findings and Results of Conventional Thyroid-Function Tests. *Clin Chem* 1980;26:1186-1192.
- Sturgess ML, Weeks I, Evans PJ, Mpoko CN, Lalng I, Woodhams JS. An Immunochemiluminometric Assay for Serum Free Thyroxine. *Clin Endocrinol* 1987;27:383-393.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline*. NCCLS Document EP5-A. Wayne, PA: NCCLS; 1999.
- Passing HA, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. *J Clin Chem Clin Biochem* 1983;21:709-720.
- National Committee for Clinical Laboratory Standards (NCCLS). *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. NCCLS Document EP17-A. Wayne, PA: NCCLS; 2004.



Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
CONJUGATE	Conjugate
CONTAINS AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL NO.	Control Number
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
MICROPARTICLES	Microparticles
PRE-TRIGGER SOLUTION	Pre-Trigger Solution
PRODUCT OF IRELAND	Product of Ireland
REACTION VESSELS	Reaction Vessels
REAGENT LOT	Reagent Lot
REF	List Number
REPLACEMENT CAPS	Replacement Caps
SAMPLE CUPS	Sample Cups
SEPTUM	Septum
SN	Serial number
TRIGGER SOLUTION	Trigger Solution
WARNING: SENSITIZER	Warning: May cause an allergic reaction.
WASH BUFFER	Wash Buffer

ARCHITECT, AxSYM, and Chemiflex are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.



Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000



DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised April 2017.

©2015, 2017 Abbott Laboratories





REF 7K64-25
REF 7K64-20
REF 7K64-35
REF 7K64-30



en
Total T3
7K64
G65478R09
B7K640

Revised April 2017.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

ARCHITECT Total T₃

INTENDED USE

The ARCHITECT Total T₃ (TT₃) assay is a chemiluminescent microparticle immunoassay (CMIA) for the quantitative determination of total triiodothyronine (Total T₃) in human serum and plasma.

SUMMARY AND EXPLANATION OF THE TEST

3,5,3' Triiodothyronine (T₃) is a thyroid hormone with a molecular weight of 651 daltons¹ and a half-life in serum of 1.5 days.² T₃ circulates in the blood as an equilibrium mixture of free and protein bound hormone.³ T₃ is bound to thyroxine binding globulin (TBG), prealbumin, and albumin. The actual distribution of T₃ among these binding proteins is controversial as estimates range from 38-80% for TBG, 9-27% for prealbumin, and 11-35% for albumin.⁴ The binding of these proteins is such that only 0.2-0.4% of the total T₃ is present in solution as unbound or free T₃.⁵ This free fraction represents the physiologically active thyroid hormone.³

It has become apparent in recent years that T₃ plays an important role in the maintenance of the euthyroid state. Serum T₃ measurements can be a valuable component of a thyroid screening panel in diagnosing certain disorders of thyroid function as well as conditions caused by iodine deficiency. Clinically, measurements of serum T₃ concentration are especially valuable in diagnosing hyperthyroidism and in following the course of therapy for this disorder.^{2, 6, 7} Under conditions of strong thyroid stimulation, the T₃ measurement provides a good estimation of thyroid reserve.² Recognition of a thyroid dysfunction called T₃-thyrotoxicosis, associated with an increased serum T₃ level but normal thyroxine (T₄), free T₄, and *in vitro* uptake results have further highlighted the importance of serum T₃ measurements.^{2, 8-11} Dietary iodine deficiency results in inadequate production of thyroid hormones despite the presence of normal thyroid tissue. In these cases, the serum T₄ concentration is often low while the Thyroid Stimulating Hormone (TSH) concentration is elevated. Elevated TSH associated with low T₄ is normally indicative of hypothyroidism. However, in iodine deficiency, these results together with normal or slightly elevated serum T₃ are indicative of euthyroid status in most individuals.¹²

T₃ levels are also affected by conditions which affect TBG concentration.¹³⁻¹⁵ Slightly elevated T₃ levels may occur in pregnancy or during estrogen therapy, while depressed levels may occur during severe illness, malnutrition, in renal failure and during therapy with anti-thyroid drugs, propranolol and propylthiouracil and salicylates.^{2, 16, 17} In patients with severe or chronic illnesses, many abnormalities of thyroid hormone balance occur. T₄ production and the extent of serum thyroid hormone binding may be independently abnormal, resulting in a low, normal or high free T₄ estimate. Serum T₃ concentrations are often low; TSH levels may be normal or slightly elevated. Total T₃ measurements may be valuable when hyperthyroidism is suspected and the free T₄ estimate is normal.¹³ The ARCHITECT Total T₃ assay is to be used as an aid in the assessment of thyroid status.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT Total T₃ assay is a two-step immunoassay to determine the presence of Total T₃ in human serum and plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

1. Sample and anti-T₃ coated paramagnetic microparticles are combined. The T₃ present in the sample binds to the anti-T₃ coated microparticles.
2. After washing, T₃ acridinium-labeled conjugate is added to create a reaction mixture.
3. Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
4. The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is an inverse relationship between the amount of Total T₃ in the sample and the RLUs detected by the ARCHITECT iSystem optics.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Kit Contents

ARCHITECT Total T₃ 7K64

NOTE: Some kit sizes are not available in all countries or for use on all ARCHITECT iSystems. Please contact your local distributor.

REF	7K64-25	7K64-20	7K64-35	7K64-30
	100	400	500	2000
MICROPARTICLES	1 x 6.6 mL	4 x 6.6 mL	1 x 27.0 mL	4 x 27.0 mL
CONJUGATE	1 x 5.9 mL	4 x 5.9 mL	1 x 26.3 mL	4 x 26.3 mL
MICROPARTICLES	Anti-T ₃ (sheep) coated microparticles in MES buffer with sheep IgG stabilizers. Minimum Concentration: 0.08% solids. Preservative: ProClin 300.			
CONJUGATE	T ₃ acridinium-labeled conjugate in citrate buffer with NaCl and Triton X-100 stabilizers. Minimum concentration: 0.33 ng/mL. Preservative: ProClin 300.			

Other Reagents

PRE-TRIGGER SOLUTION ARCHITECT Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

TRIGGER SOLUTION ARCHITECT Trigger Solution containing 0.35 N sodium hydroxide.

WASH BUFFER ARCHITECT Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

NOTE: Bottle and volume vary based on order.




Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁸⁻²¹

The following warnings and precautions apply to: MICROPARTICLES / CONJUGATE	
	
WARNING:	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not use reagent kits beyond the expiration date.
- **Do not pool reagents within a kit or between kits.**
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- **Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**
 - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
 - Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
 - Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

When stored and handled as directed, reagents are stable until the expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened/ Opened*	2-8°C	Until expiration date	May be used immediately after removal from 2-8°C storage.
On board	System temperature	30 days	Discard after 30 days. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.

* Reagents may be stored on or off the ARCHITECT iSystem. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. **If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.** For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The ARCHITECT Total T₃ assay file must be installed on the ARCHITECT iSystem prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

Alternate Result Units

Edit assay parameter "Result concentration units" to select an alternate unit.

Conversion formula:

$$(\text{Concentration in Default result unit}) \times (\text{Conversion factor}) = (\text{Concentration in Alternate result unit})$$

Default result unit	Conversion factor	Alternate result unit
ng/mL	1.536	nmol/L
	100.0	ng/dL*

* iSystem Assay CD-ROM version 6.0 and higher will be required to install this alternate result unit (ng/dL).

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

Validated specimen types to be used with this assay:

Specimen Types	Collection Tubes
Human serum	Serum Serum separator tubes
Human plasma	Sodium heparin Lithium heparin Potassium EDTA



- Other anticoagulants have not been validated for use with this assay.
- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.
- Performance has not been established for the use of neonatal specimens.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated
- For optimal results, serum and plasma specimens should be free of fibrin, red blood cells or other particulate matter.
- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for specimen collection tubes.
- Specimens must be mixed THOROUGHLY after thawing, by LOW speed vortexing or by gently inverting, and centrifuged prior to use to remove red blood cells or particulate matter to ensure consistency in the results.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

Specimen Storage

Specimen Type	Storage Temperature	Maximum Storage Time
Serum/Plasma	2-8°C	≤ 6 days

If testing will be delayed more than 24 hours, remove serum or plasma from the clot, serum separator or red blood cells.

If testing will be delayed more than 6 days, specimens should be frozen at -10°C or colder.

Specimens stored frozen at -10°C or colder for 6 days showed no performance difference.

Avoid multiple freeze/thaw cycles.

Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.

PROCEDURE

Materials Provided

7K64 ARCHITECT Total T₃ Reagent Kit

Materials Required but not Provided

- ARCHITECT Total T₃ Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on www.abbottdiagnostics.com.
- 7K64-01 ARCHITECT Total T₃ Calibrators
- 7K64-50 ARCHITECT Total T₃ Manual Diluent
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum

- ARCHITECT Replacement Caps
- Any commercially available controls
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
 - **Invert the microparticle bottle 30 times.**
 - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
 - **If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
 - Once the microparticles have been resuspended, place a septum on the bottle. For instructions about placing septums on bottles, refer to the **Reagent Handling** section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
 - Verify that all necessary reagents are present.
 - Ensure that septums are present on all reagent bottles.
- Order calibration, if necessary.
 - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
 - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

Maximum number of replicates sampled from the same sample cup: 10

- Priority:
 - Sample volume for first test: 75 µL
 - Sample volume for each additional test from same sample cup: 25 µL
- ≤ 3 hours on board:
 - Sample volume for first test: 150 µL
 - Sample volume for each additional test from same sample cup: 25 µL
- > 3 hours on board: additional sample volume is required. Refer to the ARCHITECT System Operations Manual, Section 5 for information on sample evaporation and volumes.
- If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT Total T₃ Calibrators.
 - Mix calibrator(s) by gentle inversion before use.
 - Hold bottles **vertically** and dispense recommended volumes into each respective sample cup.
 - Recommended volumes:
 - for each calibrator: 4 drops
 - Follow the manufacturer's instructions for preparation of commercially available control material.
- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.



- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Specimen Dilution Procedures

Specimens with a Total T₃ value exceeding 8.00 ng/mL are flagged with the code "> 8.00" and may be diluted using the Manual Dilution Procedure.

Manual Dilution Procedure

Suggested dilution: 1:2

It is recommended that dilutions not exceed 1:2.

1. Add a minimum of 75 µL of the patient specimen to 75 µL of ARCHITECT Total T₃ Manual Diluent.
To avoid contamination of Manual Diluent, dispense several drops of Manual Diluent into a clean test tube prior to pipetting.
2. The operator must enter the dilution factor (2) in the Patient or Control order screen. The system will use this dilution factor to automatically calculate the concentration of the sample before dilution. This will be the reported result. The dilution should be performed so that the reported result reads greater than 1.0 ng/mL.

If the operator does not enter the dilution factor, the reported result will be that of the diluted sample. This result (before dilution factor is applied) should be greater than 0.5 ng/mL.

For detailed information on ordering dilutions, refer to the ARCHITECT System Operations Manual, Section 5.

Calibration

- Test Calibrators 1 and 2 in duplicate. The calibrators should be priority loaded.
A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Calibration Range: 0.0 - 8.0 ng/mL.
- Once an ARCHITECT Total T₃ calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used or
 - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Quality Control Procedures

The recommended control requirement for the ARCHITECT Total T₃ assay is a single sample of all control levels tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures. Ensure that assay control values are within the concentration ranges specified in the package insert.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.

The ARCHITECT Total T₃ assay belongs to method group 2.

RESULTS

Calculation

The ARCHITECT Total T₃ assay utilizes a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration curve.

For information on alternate result units, refer to the INSTRUMENT PROCEDURE, Alternate Result Units section of this package insert.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

LIMITATIONS OF THE PROCEDURE

- For diagnostic purposes, results should be used in conjunction with other data; e.g., symptoms, results of other thyroid tests, clinical impressions, etc.
- If the Total T₃ results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.

EXPECTED VALUES

A normal range of 0.58 ng/mL to 1.59 ng/mL (central 95% interval) was obtained by testing serum specimens from 438 individuals determined as normal by AxSYM Ultrasensitive hTSH II and AxSYM Free T₄ assays. It is recommended that each laboratory establish its own normal range, which may be unique to the population it serves depending upon geographical, patient, dietary, or environmental factors.

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The ARCHITECT Total T₃ assay is designed to have a precision of ≤ 10% (total CV). A study based on guidance from Clinical and Laboratory Standards Institute (CLSI, formerly NCCLS) document EP5-A²² was performed for the ARCHITECT Total T₃ assay. A three member processed human serum based panel was assayed, using two lots of reagents, in replicates of two at two separate times per day for 20 testing days. Data from this study are summarized in the following table.*

Panel Member	Reagent Lot	Instrument	n	Mean Conc. Value (ng/mL)	Within Run		Total	
					SD	%CV	SD	%CV
1	1	1	80	0.75	0.021	2.7	0.027	3.6
1	1	2	80	0.73	0.023	3.1	0.030	4.1
1	2	1	80	0.79	0.036	4.6	0.043	5.4
1	2	2	80	0.81	0.047	5.8	0.057	7.0
2	1	1	80	1.50	0.029	1.9	0.110	7.3
2	1	2	80	1.49	0.040	2.7	0.053	3.6
2	2	1	80	1.53	0.031	2.0	0.035	2.3
2	2	2	80	1.54	0.040	2.6	0.049	3.2
3	1	1	80	3.27	0.062	1.9	0.128	3.9
3	1	2	80	3.29	0.107	3.3	0.140	4.2
3	2	1	80	3.55	0.054	1.5	0.071	2.0
3	2	2	80	3.54	0.066	1.9	0.077	2.2

* Representative data; results in individual laboratories may vary from these data.



Recovery

The ARCHITECT Total T₃ assay is designed to have a mean recovery of 100 ± 10% when analyzing samples spiked with known amounts of T₃. T₃ was added to nine normal human serum samples. The concentration of T₃ was determined using the ARCHITECT Total T₃ assay and the resulting percent recovery was calculated.*

Sample	Endogenous T ₃ Concentration (ng/mL)	T ₃ Added (ng/mL)	Observed Total T ₃ Concentration (ng/mL)	% Recovery**
1	2.01	0.77	2.74	94.8
2	0.97	0.78	1.64	85.9
3	1.13	0.79	1.95	103.8
4	0.99	1.54	2.43	93.5
5	0.88	1.53	2.41	100.0
6	0.90	1.54	2.54	106.5
7	1.07	3.03	4.28	105.9
8	1.23	3.04	4.21	98.0
9	0.90	3.03	3.89	98.7

Average Recovery: 98.6%

* Representative data; results in individual laboratories may vary from these data.

$$\text{** \% Recovery} = \frac{\text{Observed Total T}_3 \text{ Conc. (ng/mL)} - \text{Endogenous Total T}_3 \text{ Conc. (ng/mL)}}{\text{T}_3 \text{ Added (ng/mL)}} \times 100$$

Analytical Sensitivity

The ARCHITECT Total T₃ assay is designed to have an analytical sensitivity of ≤ 0.25 ng/mL.

Analytical sensitivity is defined as the concentration calculated as the mean plus two standard deviations of replicates of the ARCHITECT Total T₃ MasterCheck Level 0 (0.0 ng/mL). The analytical sensitivity (low-linearity) is defined in the ARCHITECT Total T₃ assay parameters as 0.25 ng/mL.

Analytical Specificity

The ARCHITECT Total T₃ assay is designed to have a mean analytical specificity of ≤ 0.1% cross reactivity with thyroxine (T₄) at a concentration of 1,100 ng/mL.

Interference

The ARCHITECT Total T₃ assay is designed to have a mean potential interference from hemoglobin, bilirubin, triglycerides, and protein of ≤ 10% at the levels indicated below.

Hemoglobin	≤ 500 mg/dL
Bilirubin	≤ 20 mg/dL
Triglycerides	≤ 2000 mg/dL
Protein	≤ 12 g/dL

Accuracy by Correlation

The ARCHITECT Total T₃ assay is designed to have a slope of 1.00 ± 0.20 and a correlation coefficient (r) of ≥ 0.90 when compared to the AxSYM Total T₃ assay.

A study was performed where specimens were tested using the ARCHITECT Total T₃ assay and AxSYM Total T₃ assay. Data from this study were analyzed using Least Squares and Passing-Bablok²³ regression methods and are summarized in the following table.*

Method	Abbott ARCHITECT Total T ₃ vs. Abbott AxSYM Total T ₃			Correlation Coefficient
	Number of Specimens	Intercept	Slope	
Least Squares				
Linear Regression	1440	0.02	1.04	0.954
Passing-Bablok				
Linear Regression**	1440	-0.08	1.13	0.964

In this evaluation, serum specimens tested ranged from 0.25 ng/mL to 5.83 ng/mL with the ARCHITECT Total T₃ assay and from 0.34 ng/mL to 5.19 ng/mL with the AxSYM Total T₃ assay.

* Representative data; variables such as differences in sampling size and sample population may impact correlation of the assay; therefore, results in individual laboratories may vary from these data.

** A linear regression method with no special assumptions regarding the distribution of the samples and the measurement errors.²³

BIBLIOGRAPHY


- Budavari S, editor. *Merck Index* (11th Ed.). Rahway, NJ: Merck and Co., Inc., 1989:868.
- Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972;21:1073-1092.
- Ekins RP, editor. *Methods for the Measurement of Free Thyroid Hormones*. Amsterdam: Excerpta Medica Foundation, 1979:72-92.
- Robbins J, Rall JE. The Iodine-Containing Hormones. In: *Hormones in Blood* (3rd Ed.). London: Academic Press, 1979:1:632-667.
- DeGroot LJ, Larsen PR, Refetoff S, Stanbury JB. Transport of Thyroid Hormone and Cell Uptake. In: *The Thyroid and Its Diseases*. New York: Wiley and Sons, 1984:62-66.
- Wahner HW, Gorman CA. Interpretation of Serum Tri-iodothyronine Levels Measured by the Sterling Technic. *N Engl J Med* 1971;284:225-230.
- Marsden P, McKerron CG. Serum Triiodothyronine Concentration in the Diagnosis of Hyperthyroidism. *Clin Endocrinol* 1975;4:183-189.
- Ivy HK, Washner HW, Gorman CA. Triiodothyronine (T₃) Toxicosis: Its Role in Graves' Disease. *Arch Intern Med* 1971;128:529-534.
- Hollander CS, Mitsuma T, Nihei N, Shenkman L, Burday SZ, Blum M. Clinical and Laboratory Observations in Cases of Triiodothyronine Toxicosis Confirmed by Radioimmunoassay. *Lancet* 1972;1:609-611.
- Sterling K, Refetoff S, Selenkow HA. T₃ Thyrotoxicosis: Thyrotoxicosis Due to Elevated Serum Triiodothyronine Levels. *JAMA* 1970;213: 571-575.
- Hollander CS, Mitsuma T, Shenkman L, Stevenson C, Pineda G, Silva E. T₃ Toxicosis in an Iodine-Deficient Area. *Lancet* 1972;2:1276-1278.
- Ermans AM. Disorders of Iodine Deficiency. In: Ingbar SH, Braverman LE, editors. *The Thyroid* (5th Ed.). Philadelphia: JB Lippincott Co., 1986:705-721.
- Kaplan MM, Larsen PR, Crantz FR, Dzau VJ, Rossing TH, Haddow JE. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. *Am J Med* 1982;72:9-16.
- Bermudez F, Surks MI, Oppenheimer JH. High Incidence of Decreased Serum Triiodothyronine Concentration in Patients with Nonthyroid Disease. *J Clin Endocrinol Metab.* 1975;41:27-40.
- Oppenheimer JH. Thyroid Function Tests in Nonthyroidal Disease. *J Chronic Dis* 1982;35:697-701.
- Abuid J, Larsen PR. Triiodothyronine and Thyroxine in Hyperthyroidism: Comparison of the Acute Changes During Therapy with Antithyroid Agents. *J Clin Invest* 1974;54:201-208.
- Felig P, Baxter JD, Broadus AE, Frohman LA, editors. *Endocrinology and Metabolism* (2nd Ed.). New York: McGraw-Hill Book Co., 1987:408-416.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline*. NCCLS Document EP5-A. Wayne, PA: NCCLS; 1999.
- Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. Application of linear regression procedures for method comparison studies in clinical chemistry, Part I. *J Clin Chem Clin Biochem* 1983;21(11):709-720.



Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
CONJUGATE	Conjugate
CONTROL NO.	Control Number
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
MICROPARTICLES	Microparticles
PRE-TRIGGER SOLUTION	Pre-Trigger Solution
PRODUCT OF IRELAND	Product of Ireland
REACTION VESSELS	Reaction Vessels
REAGENT LOT	Reagent Lot
REF	List Number
REPLACEMENT CAPS	Replacement Caps
SAMPLE CUPS	Sample Cups
SEPTUM	Septum
SN	Serial number
TRIGGER SOLUTION	Trigger Solution
WARNING: SENSITIZER	Warning: May cause an allergic reaction.
WASH BUFFER	Wash Buffer

ARCHITECT, AxSYM, MasterCheck, and Chemflex are trademarks of Abbott Laboratories in various jurisdictions. Triton is property of its respective owner.

 Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000



DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised April 2017.
©2005, 2017 Abbott Laboratories

