**SV600** 

**Ventilator** 

**Operator's Manual** 



© 2017-2019 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.

For this Operator's Manual, the issue date is December, 2019.

## 4.8 Screen Capture

By pressing this key on the main screen \_\_\_\_\_, the system will capture and save the screen automatically. The screen capture is saved in "bmp" format. The system can store up to 20 screen captures.

#### 4.9 Lock Screen

Press the [Lock] key on the main screen to enter locked status, and the prompt message [Screen locked. Press the Lock button to unlock screen.] will be displayed. During the

period of screen locked, only  $[O_2 \uparrow Suction]$ , and [Lock] keys are enabled. Touch screen, control knob, and other keys are disabled. Press this key a second time to unlock the screen.

# **5** System Settings

## 5.1 Date & Time Settings

- 1. Select the system time field on the main screen to pop up time setup menu.
- 2. Set [Date] and [Time].
- 3. Set [Date Format] to [YYYY-MM-DD], [MM-DD-YYYY] or [DD-MM-YYYY].
- 4. Set [Time Format]: [24 h] or [12 h].

### 5.2 Export to USB

The ventilator's exportation function provides the ability to export some data or settings to USB device.

#### 5.2.1 Export Screen

Screen exportation involves exporting a saved screen capture for the ventilator. The exported file is saved in "bmp" format. This ventilator could save up to 20 screen captures. To export screen capture,

1. Insert the USB device into the USB connector of the ventilator. The key is highlighted on the main screen.

- 2. By selecting the key, the system will open the USB settings interface.
- 3. On the opened interface, select the [Export Screenshot] tab first and then click the [Export Screenshot] key. The system will run a check to verify that there is enough storage space available on the USB device. If there is sufficient space, the system will start to export the screen.
- 4. After exporting is completed, select [Remove USB Device] to remove the USB device.

### 5.2.2 Export Data

Exporting data means to export data from the ventilator, such as patient demographics, current setting parameters, current alarm limits, trend data and so on.

To export data,

- 1. Insert the USB device into the USB connector of the ventilator. The key is highlighted on the main screen.
- 2. By selecting the key, the system will open the USB settings interface.

- 3. On the opened interface, select the [Export Data] tab and then select the [User Export] key. The system will run a check to verify that there is enough storage space available on the USB device. If there is sufficient space, the system will export data including patient information, current parameter settings, current alarm limits, tabular trend, PEEPi measured value, P0.1 measured value, Vtrap measured value, and NIF measured value, etc. The format of the exported data is "html".
- 4. If you need to export calibration data, event logbook and self-check logbook in addition to the above data, select the [Factory Export] tab and enter password. The system will run a check to verify that there is sufficient storage space available on the USD device. If there is sufficient space, the system will start to export data. The exported data is encrypted in the format of "blg".
- 5. After exporting is completed, select [Remove USB Device] to remove the USB device.

#### NOTE

• If you need to check the exported data in format of "blg", please contact the Customer Service Department.

#### 5.2.3 Transfer Settings

You can export or import settings, while unit is in standby.

To export settings,

- 1. Make sure that the machine is in Standby status.
- 2. Insert the USB device into the USB connector of the ventilator. The key is highlighted on the main screen
- 3. By selecting the key, the system will open the USB settings interface.
- 4. Select [Transfer Settings] → Enter system password → [Export Settings] in the opened interface. The system will run a check to verify that there is sufficient storage space available on the USB device. If there is sufficient space, the system will save the current settings and machine defaults to the USB device.
- 5. After exporting is completed, select [**Remove USB Device**] to remove the USB device.

# 6 Start Ventilation

### 6.1 Turn on the System

- 1. Insert the power cord into the power receptacle. Ensure the external power indicator light is lit.
- 2. Press the OO hard key.
- 3. The alarm indicator light flashes yellow and red once in turn, and then the system conducts a self check of the speaker and buzzer once respectively.
- 4. A start-up screen and start-up check progress bar appear. Then the System Check screen is displayed.

#### NOTE

 When the ventilator is started, the system detects whether audible alarm tones and alarm lamp function normally. If yes, the alarm lamp flashes yellow and red successively, and the speaker and the buzzer give check tones. If not, do not use the equipment and contact us immediately.

## 6.2 System Check

# **ACAUTION**

- If the ventilator fails any tests, remove it from clinical use. Do not use the ventilator until necessary repairs are completed and all tests have passed.
- Before running System Check, disconnect the patient from the equipment and ensure that a backup ventilation mode is available for patient ventilation.

To enter the System Check screen,

- The System Check screen is accessed automatically after powering on the system.
- On the non-standby screen, select the [Standby] key and enter the Standby status after your confirmation. Select the [System Check] key in the Standby status to enter the System Check screen.

The system check screen displays the last system check time and total system check result. Select the key to query the last system check information of the ventilator system, including system check items and System Check results.

Connect the gas supply and block the Y piece as illustrated. Then select [Continue] to start System Check item by item.

#### System Check items include:

- Backup Air Supply Test: test the speed of backup air supply.
- $\blacksquare$  O<sub>2</sub> Flow Sensor Test: test the O<sub>2</sub> Insp. Valve and O<sub>2</sub> Flow Sensor.
- Air Flow Sensor Test: test the Air Insp. Valve and Air Flow Sensor.
- Exp. Flow Sensor Test: test the expiratory flow sensor.
- Pressure Sensor Test: test the pressure sensors at the inspiratory and expiratory ports.
- Exp. Valve Test
- Safety Valve Test
- Leakage (mL/min)
- $\blacksquare$  Compliance (mL/cmH<sub>2</sub>O)
- $\blacksquare$  Circuit Resistance (cmH<sub>2</sub>O/L/S)
- O<sub>2</sub> Sensor Test
- Neonatal Flow Sensor Test

#### System Check result can be:

- Pass: indicates that check of this item is completed and is passed;
- Fail: indicates that check of this item is completed but is failed;
- Cancel: indicates that check of this item is cancelled;
- No Gas Supply: indicates that air or O<sub>2</sub> sources are not connected.
- Monitoring Off: indicates that sensor monitoring function may not be switched on when O₂ sensor test or neonatal flow sensor test is being carried out.
- No Sensor: indicates that the neonatal flow sensor is not connected.
- Sensor Reversed: indicates that the neonatal flow sensor is connected reversed.
- Sensor Failure: indicates that the oxygen sensor may not be working.
- High leakage: indicates that there is high leakage from the test tubing, probably because the tubing is disconnected, not properly installed, the safety valve is not closed, or the expiratory valve membrane is not installed.

Total selftest results are listed as follows after all selftest items have been completed:

- Pass: all selftest items successfully pass the seftest.
- Partially Pass: some selftest items fail, but the mechanical ventilation is allowed.
- Fail. Ventilation Disabled: some important selftest items fail, but the mechanical ventilation is not allowed.
- High Leakage, Ventilation Disabled: Exp. Flow Sensor Test, Pressure Sensor Test, Exp. Valve Test, or Safety Valve Test fails, the mechanical ventilation is not allowed.
- Cancel: some selftest items cancelled and other selftest items have been successfully passed.

# **WARNING**

• During the use of the electronic nebulizer, please pay attention to the connection of the nebulizer to prevent the nebulization interruption.

## 10.5 O<sub>2</sub>↑(Oxygen Enrichment)

 $O_2\uparrow$  is also called as  $O_2$  enrichment. It means to deliver oxygen with concentration higher than normal level within the specified time period. The oxygenation magnitude can be set by selecting [**Menu**]  $\rightarrow$  [**Setup**]  $\rightarrow$  [**Ventilation**]. The default oxygen enrichment magnitude is 60% for adult and pediatric patients, and 10% for neonate patients.

Press the  $[O_2 \uparrow Suction]$  key and the ventilator starts oxygen enrichment. At that time, the indicator light for  $[O_2 \uparrow Suction]$  key will be illuminated, and the remaining oxygen enrichment time will be displayed. Oxygen enrichment is active for maximum two minutes. During oxygen enrichment, the currently set oxygen concentration is displayed in the  $[O_2 \%]$  parameter setup quick key field.

When the 2-minute period of oxygen enrichment is up or the  $[O_2 \uparrow Suction]$  key is pressed again, the ventilator terminates oxygen enrichment.

#### NOTE

- The system cannot start O<sub>2</sub>↑ (oxygen enrichment) in the standby, oxygen therapy, or CPRV modes.
- The system cannot start O<sub>2</sub>↑ (oxygen enrichment) in the P-V tool test process.
- When [O₂ Supply Failure] alarm or [No Gas Supply Pressure] alarm is triggered, click [O₂↑ Suction] key, O₂↑ is disabled and prompts [O₂ Supply Failure, O₂↑ disabled].
- If  $O_2 \uparrow$  process triggers  $[O_2 \text{ Supply Failure}]$  alarm or [No Gas Supply Pressure] alarm,  $O_2 \uparrow$  stops.
- Removing the patient tubing during oxygen enrichment will start suction function.
   Refer to 10.6Suction section.

# **B** Product Specifications

The ventilator is already integrated with an expiratory volume monitor, pressure measurement device, pressure release device, built-in gas mixer, alarm system, SpO<sub>2</sub> monitor, O<sub>2</sub> monitor, and CO<sub>2</sub> monitor. Among them:

- The expiratory volume monitor, pressure measurement device, and pressure release device comply with ISO 80601-2-12.
- The alarm system complies with IEC 60601-1-8.
- The gas mixer complies with ISO 11195.
- The SpO<sub>2</sub> monitor complies with ISO 80601-2-61;
- The  $O_2$  monitor complies with ISO 80601-2-55.
- The  $CO_2$  monitor complies with ISO 80601-2-55.
- The gas supply hose assembly complies with ISO 5359.

# **B.1 Safety Specifications**

Type of protection against electric shock	Ventilator: Class I device with internal electrical power supply.  Air compressor: Class I device.
Degree of protection against electric shock	Ventilator: Mixed BF and CF applied part type, with respiratory circuit, and CO <sub>2</sub> being BF type, and SpO <sub>2</sub> being CF type.  Air compressor: No applied part.
Operating mode	Continuous
Degree of protection against	Ordinary equipment, without protection against explosion; not
hazards of explosion	for use with flammable anaesthetics.
Degree of protection against harmful ingress of water	Ventilator: Degrees of protection provided by enclosures(IP Code)—IP21 Air compressor: Degrees of protection provided by enclosures(IP Code)—IPX0 Protection Index according the EN 60529 standard: 2: Protected against solid foreign objects of 12.5 mm diameter and greater 1: Protected against vertically falling water drops 0: no protection
Disinfection and sterilization methods.	Ventilator: The device disinfection and sterilization methods are recommended by manufacturer.  Air compressor: The device does not require disinfection and sterilization.
Equipment type	Mobile

# **D** Alarm Messages

This chapter lists physiological and technical alarm messages.

Note that in this chapter:

- Column P stands for the default alarm level: H for high, M for medium and L for low.
- For each alarm message, corresponding actions are given instructing you to troubleshoot problems. If the problem persists, contact your service personnel.

# **D.1 Physiological Alarm Messages**

#### **D.1.1 Ventilator Parameters**

Alarm Messages	P	Cause and action
		The airway pressure exceeds the set pressure high alarm limit.
		1. Check the patient.
Paw Too High	Н	2. Check the ventilation parameter setup.
		3. Check the alarm limits.
		4. Check the patient tubing for occlusion.
		Airway pressure setting is lower than the low limit of pressure alarm.
		1. Check the patient.
Paw Too Low	Н	2. Check the ventilation parameter setup.
		3. Check the alarm limits.
		4. Check if the patient tubing are leaked or disconnected.
		The inspired O <sub>2</sub> concentration is greater than the FiO <sub>2</sub> high alarm limit
	Н	for at least 30s.
FiO <sub>2</sub> Too High		1. Check air supply.
110 <sub>2</sub> 100 High		2. Check the HEPA filter for occlusion.
		3. If the ventilator uses the $O_2$ cell, calibrate the $O_2$ sensor. If the
		ventilator uses the paramagnetic O <sub>2</sub> sensor, perform the System Check.
		The inspired O <sub>2</sub> concentration has been lower than the FiO <sub>2</sub> low alarm
	Н	limit for at least 30 s or is less than 18%.
FiO <sub>2</sub> Too Low		1. Check air supply.
		2. If the ventilator uses the $O_2$ cell, calibrate the $O_2$ sensor. If the
		ventilator uses the paramagnetic O <sub>2</sub> sensor, perform the System Check.
TVe Too High		The TVe monitored value is greater than TVe high alarm limit for
	M	continuous 3 mechanical ventilation cycles.
		1. Check the ventilation parameter setup.
		2. Check the alarm limits.
TVe Too Low	M	The TVe monitored value is less than TVe low alarm limit for

		continuous 3 mechanical ventilation cycles.
		1. Check the patient.
		2. Check the ventilation parameter setup.
		3. Check the alarm limits.
		4. Check the patient tubing for leakage or occlusion.
		5. Perform System Check to test the leakage
		MVe is greater than MVe high alarm limit.
MVe Too High	Н	1. Check the ventilation parameter setup.
		2. Check the alarm limits.
		MVe is less than MVe low alarm limit.
		1. Check the ventilation parameter setup.
MVe Too Low	Н	2. Check the alarm limits.
		3. Check the patient tubing for leakage or occlusion.
		4. Perform System Check to test the leakage
		The time of failure to detect respiration exceeds Tapnea.
		1. Check the patient.
Apnea	Н	2. Manual breath.
		3. Check apnea time setup.
		4. Check if the patient tubing are disconnected.
		The time of failure to detect respiration exceeds Tapnea. Start apnea
Apnea Vent	Н	ventilation mode.
		Check apnea ventilation parameter setup.
		ftotal is greater than ftotal high alarm limit.
		1. Check the patient.
ftotal Too High	M	2. Check the ventilation parameter setup.
		3. Check the alarm limits.
		ftotal is lower than the ftot low alarm limit.
0.15		1. Check the patient.
ftotal Too Low	M	2. Check the ventilation parameter setup.
		3. Check the alarm limits.
Apnea		
Ventilation	L	This alarm is given when apnea ventilation ends. There is no need to
Ended		process this alarm.

# D.1.2 CO<sub>2</sub> Module

Alarm Messages	P	Cause and action
EtCO <sub>2</sub> Too High	M	The monitored parameter value exceeds the alarm limit.
		1. Check the patient type.
		2. Check the alarm limits.
		The monitored parameter value exceeds the alarm limit.
EtCO <sub>2</sub> Too Low	M	1. Check the patient type.
		2. Check the alarm limits.

		The time of failure to detect respiration by the CO <sub>2</sub> module exceeds
		Apnea Tinsp. Whenever the CO <sub>2</sub> apnea alarm is on, block the [EtCO <sub>2</sub>
		Too High] alarm and [EtCO <sub>2</sub> Too Low] alarm until the alarm is
Apnea CO <sub>2</sub>	M	cleared.
		1. Check the patient.
		2. Check apnea time setup.
		3. Check the connections of CO <sub>2</sub> module sampling device.

# D.1.3 SpO<sub>2</sub> Module

Alarm Messages	P	Cause and action
	M	SpO <sub>2</sub> value is greater than the high alarm limit.
SpO <sub>2</sub> Too High		Check the patient's condition and ventilator settings.
3pO <sub>2</sub> 100 High		2. Check the patient's inspiratory O <sub>2</sub> %.
		3. Check the alarm limits.
		SpO <sub>2</sub> value is lower than the low alarm limit.
SpO <sub>2</sub> Too LOW	M	Check the patient's condition and ventilator settings.
3pO <sub>2</sub> 100 LOW	171	2. Check the patient's inspiratory O <sub>2</sub> %.
		3. Check the alarm limits.
		SpO <sub>2</sub> value is lower than the desaturation alarm limit.
SpO <sub>2</sub> Desat	Н	Check the patient's condition and ventilator settings
Spo <sub>2</sub> Desar	11	2. Check the patient's inspiratory O <sub>2</sub> %.
		3. Check the alarm limits.
	M	PR value exceeds the high alarm limit.
PR Too High		1. Check the patient's condition.
TK 100 mgn		2. Check ventilator settings.
		3. Check the alarm limits.
	M	PR value is lower than the low alarm limit.
PR Too LOW		1. Check the patient's condition.
T K 100 LOW		2. Check ventilator settings.
		3. Check the alarm limits.
	Н	The patient's pulse signal is too weak, and the system cannot perform
No pulse		analysis.
140 puise		1. Check the patient's condition.
		2. Check SpO <sub>2</sub> sensor and measurement site connection

# **D.2 Technical Alarm Messages**

# **D.2.1 Power Board**

Alarm Messages	P	Cause and action
Battery 1 Failure	Н	Battery 1 Charge Failure
02		Contact your service personnel.
Battery 1 Failure	,,	Battery 1 Aging
03	Н	Contact your service personnel.
Battery 1 Failure	11	Battery 1 Comm Error
04	Н	Contact your service personnel.
Battery 1 Failure	,,	Battery 1 Failure
05	Н	Contact your service personnel.
Battery 2 Failure	11	Battery 2 Charge Failure
02	Н	Contact your service personnel.
Battery 2 Failure	11	Battery 2 Aging
03	Н	Contact your service personnel.
Battery 2 Failure	,,	Battery 2 Comm Error
04	Н	Contact your service personnel.
Battery 2 Failure		Battery 2 Failure
05	Н	Contact your service personnel.
Blower Battery		Backup air supply battery failed.
Failure 02	Н	Contact your service personnel.
Blower Battery	7.7	Backup air supply battery failed.
Failure 03	Н	Contact your service personnel.
Blower Battery	11	Backup air supply battery failed.
Failure 04	Н	Contact your service personnel.
Blower Battery	,,	Backup air supply battery failed.
Failure 05	Н	Contact your service personnel.
Battery Temp.		Battery temperature is a bit high during discharge.
High. Connect	M	
Ext.Pwr.		Connect to the external power supply.
Battery Temp		Battery temperature is too high during discharge. The system may be
High. Syst maybe	Н	down.
Down		Connect to the external power supply.
Dattomy in Haa	T	The current system is powered by battery.
Battery in Use	L	Connect to the external power supply.
Low Battery.		The remaining battery power is lower than a threshold.
Connect Ext. Power.	M	Connect to the external power supply.
System DOWN.	Н	Battery power is depleted. The system will shut down in a few
Connect Ext.	L	minutes.

Power.		Connect to the external power supply immediately.
Battery	11	No battery in main unit or backup air supply at present
Undetected	H	Contact your service personnel.
		Power board fan speed abnormal. If it can't be solved, please restart
Fan Failure	M	the machine.
		Contact your service personnel.
Davis a Failure 02	Н	Power Board Selftest Error.
Device Failure 03		Contact your service personnel.

# **D.2.2 Main Control Board**

Alarm Messages	P	Cause and action
Please Reset Date	L	Button cell is available in the system. But the clock is powered down and reset.
and Time		Re-set the date and time.
Key Error	L	Hardkey or rotary encoder is depressed continuously for more than 35s.
		Contact your service personnel.
Device Failure 04	Н	Ctrl Module Init Error.
Device railule 04		Contact your service personnel.
Device Failure 05	Н	Ctrl Module Comm Stop.
Device railule 03		Contact your service personnel.
Device Failure 19	Н	Power Board Comm Stop.
Device Failule 19	11	Contact your service personnel.
Device Failure 20	Н	SpO <sub>2</sub> Module Comm Stop.
Device Famule 20		Restart the ventilator or contact your service personnel.
Device Failure 22	Н	Protecting Module Comm Stop.
Device Failule 22	П	Contact your service personnel.

# **D.2.3 Monitor Board**

Alarm Messages	P	Cause and action
Technical Error	L	Buzzer Failure.
04	L	Contact your service personnel.
Technical Error	М	Atmospheric Pressure Sensor Failure.
05	IVI	Contact your service personnel.
Technical Error	M	3-way Valve Failure.
07	IVI	Contact your service personnel.
Technical Error	M	Nebulizer Valve Failure.
08		Contact your service personnel.
Technical Error	М	Insp. Temp Sensor Failure.
09	IVI	Contact your service personnel.

Technical Error	_	Heating function of the expiration valve is faulty.
10	L	Contact your service personnel.
- · - · · · · · · · · · · · · · · · · ·	Н	Power Supply Voltage Error.
Device Failure 01		Contact your service personnel.
		Memory Error.
Device Failure 02	Н	Contact your service personnel.
		Ctrl Module Comm Stop.
Device Failure 05	Н	Contact your service personnel.
		Ctrl Module Selftest Error.
Device Failure 06	Н	Contact your service personnel.
		Pressure Sensor Failure.
Device Failure 09	Н	Contact your service personnel.
		Safety Valve Failure.
Device Failure 10	Н	Contact your service personnel.
Device Failure 12	Н	Air Insp. Limb Failure.
		Contact your service personnel.
Device Failure 13	Н	O <sub>2</sub> Limb Failure.
		Contact your service personnel.
Device Failure 21	Н	Pressure Sensor Zero Point Error.
		Contact your service personnel.
Device Failure 22	Н	Protecting Module Comm Stop.
		Contact your service personnel.
Device Failure 23	Н	Protection Module Self Check Error.
		Contact your service personnel.
		Monitored PEEP exceeds PEEP + 5 cmH <sub>2</sub> O (PEEP + 10 cmH <sub>2</sub> O for
PEEP Too High	Н	APRV mode) within any fully mechanical ventilation cycle.
		1. Check the ventilation parameter setup.
		2. Check the patient tubing for occlusion.
		Patient's PEEP is less than the setting value to a certain extent.
PEEP Too Low	M	1. Check the patient tubing for leakage.
		2. Perform System Check to test the leakage
Airway		Tube is occluded.
Obstructed?	Н	1. Check and clean the patient tubing.
Obstructed:		2. Check and clean the expiration valve.
Insp. Limb		The patient tubing is bent or occluded in case of $O_2$ therapy.
Airway	M	
Obstructed?		Check if the patient tubing is occluded or bent. If yes, clear it.
		The airway pressure measured by any pressure sensor is greater than
Sustained Airway	Н	the setting PEEP + 15 cm $H_2O$ for 15 s consecutively.
Pressure		1. Check the patient.
1 1088016		2. Check the ventilation parameter setup.
		3. Check the patient tubing for occlusion.
Airway Leak?	L	Tube is leaky.

		1. Check the patient tubing for leakage.
		Perform System Check to test the leakage
Tube		Tube is disconnected.
Disconnected?	Н	
Disconnected?		Re-connect the patient tubing.
		In volume mode or pressure mode when ATRC function is enabled,
		the pressure reaches Paw high alarm limit-5.
Pressure Limited	L	1. Check the patient.
		2. Check the ventilation parameter setup.
		3. Check pressure high alarm limit.
		In pressure mode, delivered gas volume exceeds the set TV high limit.
Volume Limited	L	1. Check the patient.
		2. Check the ventilation parameter setup.
		3. Check the alarm limits.
		Pinsp is lower than the pressure setting value by 3 cmH <sub>2</sub> O or 2/3 of
	v	the pressure setting value, whichever is less.
Pinsp Not		1. Check the patient.
Achieved	L	2. Check TV alarm limits.
7 temeved		3. Check the O <sub>2</sub> supply.
		4. Check the patient tubing for leakage.
		5. Check the HEPA filter for occlusion.
		TVi is less than the TV setting value by more than 10 mL + 10 % of
		the setting value.
		1. Check the patient.
TV Not Achieved	L	2. Check pressure high alarm limit.
I v Not Achieved		3. Check the high-pressure gas supply or the HEPA filter for
		occlusion.
		4. Check the O <sub>2</sub> supply.
		5. Check the patient tubing for leakage or occlusion.
		The pressure reaches Paw high alarm limit-5 in sigh cycle.
Danasana I innita i		1. Check the patient.
Pressure Limited	L	2. Check pressure high alarm limit.
in Sigh cycle		3. Check the patient tubing for occlusion.
		4. Consider to turn off sigh.
		Oxygen supply is not sufficient to support normal ventilator operation.
O <sub>2</sub> Supply Failure	Н	1. Check connection with O <sub>2</sub> supply.
		2. Check O <sub>2</sub> supply pressure.
		Air supply is not sufficient to support normal ventilator operation.
Air Supply Failure	Н	1. Check connection with Air supply.
	**	2. Check air supply pressure
		Both oxygen and air supply are not sufficient to support normal
		ventilator operation.
No Gas Supply	Н	1. Check connection with air and O <sub>2</sub> supply.
Pressure	**	2. Check air and O <sub>2</sub> supply pressure.
		3. For machines with backup air supply configuration, check whether
		3. For machines with backup an suppry configuration, theck whether

1	1	
		the Blower Disabled switch for user maintenance is on.
		4. Check backup air supply for failure.
		In PSV mode, Tinsp exceeds 4s for adult, 1.5s for pediatric, and the
		maximum inspiration time set by the user for neonates for continuous
Tinsp Too Long	L	3 cycles.
Thisp Too Long	L	1. Check the patient.
		2. Check the ventilation parameter setup.
		3. Check the patient tubing for leakage.
Please Check Exp.	Н	Installing the expiratory flow sensor fails.
Flow Sensor		Contact your service personnel.
		The gas temperature exceeds 55°C.
Insp. Gas Temp	Н	1. Disconnect the patient.
Too High		2.Restart the machine. Contact the specified service personnel if the
		issue persists.
Flow Sensor Type		Installation error with air flow sensor or $O_2$ flow sensor.
Error	Н	Contact your service personnel.
21101		Backup air supply fan speed error. If it can't be solved, restart the
Blower Fan		machine.
Failure	M	Please contact your service personnel (turning off backup air supply
1 andic		could also resolve the alarm).
		Backup air supply temperature exceeds the threshold.
		1. Check if the operating ambient temperature of the machine exceeds
Blower		the maximum operating temperature specified by the vendor.
Temperature High	H	2. Check if the fan inlet and outlet are occluded. If yes, clear the
		foreign substance and dust.
		3. Check the rotation of the fan. If it runs abnormally (such as
		abnormal sound or rotation speed), replace the fan.
AMV: Cannot		Cannot meet established MV%
Meet Target	L	1. Check the ventilation parameter setup.
-		2. Check the alarm limits setting.
Technical Error.		Three-way valve failure, only blower gas supply available.
Only Blower Gas	Н	
Supply Available.		Contact your service personnel.
Blower Failure		Three-way valve failure, blower module disabled.
3-way Valve	Н	
Failure		Contact the specified service personnel.
Replace HEPA	L	HEPA filter occluded, resistance increased.
1		1
Filter	L	Contact the specified service personnel.
Filter Blower Technical		Contact the specified service personnel.  Backup air supply Temp Sensor Failure.
	M	
Blower Technical	M	Backup air supply Temp Sensor Failure.
Blower Technical Error 01		Backup air supply Temp Sensor Failure.  Contact your service personnel.

Error 03		Contact your service personnel.
DI E I OI		Insp. Limb valve or flow sensor fails.
	Н	1. Use another device for ventilation.
Blower Failure 01	H	2.Restart the machine.
		3. Contact the specified service personnel if the issue persists.
Blower Failure 02	Н	Insp. Valve Disconnected.
Blower Failure 02		Contact your service personnel.
Blower Failure 03	Н	Backup air supply Temp Too High.
Blower Failure 03	11	Contact your service personnel.
Blower Failure 04	Н	Backup air supply Failure.
Blower Failure 04	11	Contact your service personnel.
O <sub>2</sub> Sensor	L	The $O_2$ sensor is not connected.
Unconnected	L	Connect the O <sub>2</sub> sensor.
Please Replace O <sub>2</sub>	M	The chemical O <sub>2</sub> sensor is expired.
Sensor.	IVI	Please replace the O <sub>2</sub> sensor.
Please calibrate	L	Please calibrate the O <sub>2</sub> sensor.
O <sub>2</sub> sensor	L	Please calibrate O <sub>2</sub> concentration.
Please reset O <sub>2</sub>	M	The oxygen concentration measured by the paramagnetic oxygen
sensor		sensor has a large error.
		Contact your service personnel.
Please perform		Calibrate the pressure sensor.
pressure	Н	
calibration.		Contact your service personnel.
Please perform	Н	Calibrate the flow sensor.
flow calibration.	Н	Please perform flow calibration.

# D.2.4 CO<sub>2</sub> Module

Alarm Messages	P	Cause and action
CO <sub>2</sub> Module Failure 01	М	Sidestream CO <sub>2</sub> module zeroing fails. The gain input signal offset is too large, exceeding the adjustable range.
CO <sub>2</sub> Module Failure 02	M	Contact your service personnel.  CO <sub>2</sub> Init Error. An error occurs to the CO <sub>2</sub> module during initialization.
CO <sub>2</sub> Module Failure 03	M M	Contact your service personnel.  CO <sub>2</sub> self check error. An error occured in the CO <sub>2</sub> module during self check.
		Contact your service personnel.
CO <sub>2</sub> Module Failure 04		CO <sub>2</sub> Hardware Error. Contact your service personnel.
CO <sub>2</sub> Module Failure 05	М	CO <sub>2</sub> Comm Stop, CO <sub>2</sub> Module Failure, CO <sub>2</sub> Comm Error or communication failure reaches 10s.  Contact your service personnel.

CO <sub>2</sub> Module	M	Mainstream CO <sub>2</sub> module zeroing fails.
Failure 06		Contact your service personnel.
CO <sub>2</sub> Sensor High	L	The sensor temperature is too high (above 63°C).
Temp		Contact your service personnel.
	_	Sampling line is faulty or occluded.
CO <sub>2</sub> Sampleline		1. Check the sampling line for occlusion.
Occluded	L	2. Replace the sampling line.
		3. Replace the water trap.
	L	The water trap is disconnected or not connected properly. Check the
CO <sub>2</sub> No Watertrap		water trap.
		Re-install the water trap.
	L	Parameter measured values exceed the measurement range (error
Et CO <sub>2</sub> Overrange		range is included).
		1. Perform CO <sub>2</sub> module zeroing.
		2. Contact your service personnel.
Please Replace	M	The mainstream CO <sub>2</sub> module sensor is faulty.
CO <sub>2</sub> Sensor		Contact your service personnel.
CO. No Sangar	L	The mainstream CO <sub>2</sub> module sensor is not connected.
CO <sub>2</sub> No Sensor		Connect the CO <sub>2</sub> sensor.

# D.2.5 SpO<sub>2</sub> Module

Alarm Messages	P	Cause and action
SpO <sub>2</sub> Sensor Off		Connected SpO <sub>2</sub> sensor became disconnected from patient tubing (e.g.
	L	wire disconnection or short circuit).
		Check SpO <sub>2</sub> sensor and measurement site connection.
Dlaga Daplaga	M	SpO <sub>2</sub> sensor failed (e.g. wire disconnection or short circuit).
Please Replace		1. Replace SpO <sub>2</sub> sensor.
SpO <sub>2</sub> Sensor		2. Contact your service personnel.
	L	Main cable has disconnected from module. Connection between
SpO <sub>2</sub> No Sensor		sensor and main cable has disconnected.
		Check that SpO <sub>2</sub> cable is connected to the module.
	L	The light to which the sensor is exposed is so bright that the sensor's
SpO <sub>2</sub> Too Much		photodetector is absorbing the surrounding light.
Light		Put SpO <sub>2</sub> sensor in a place with lower ambient light levels.
	L	SpO <sub>2</sub> sensor cannot obtain pulse signal (or incomplete signal).
Caro Na Dalas		1. Check the patient's condition.
SpO <sub>2</sub> No Pulse		2. Check SpO <sub>2</sub> sensor and measurement site connection
		3. Replace SpO <sub>2</sub> sensor.
G.O. M. 1.1.	М	SpO <sub>2</sub> module error\SpO <sub>2</sub> initialization error
SpO <sub>2</sub> Module		1. Replace SpO <sub>2</sub> sensor.
Error		2. Contact your service personnel.
SpO <sub>2</sub> Overrange	L	Measured values of parameter SpO <sub>2</sub> exceed the measurement range.

		<ol> <li>Replace SpO<sub>2</sub> sensor.</li> <li>Contact your service personnel.</li> </ol>
PR Overrange	L	Measured values of parameter PR exceed the measurement range.
		1. Replace SpO <sub>2</sub> sensor.
		2. Contact your service personnel.

## D.2.6 Neo. Module

Alarm Messages	P	Cause and action
Reverse the		Neonatal flow sensor connected reversed.
neonatal flow	M	
sensor.		Please reverse the neonatal flow sensor.
Noo Flow Congor		Range of neonatal flow sensor exceeds 32 L/min.
Neo. Flow Sensor	M	Check the patient's condition and ventilator settings
Overrange		2. Change patient type if necessary.
Neo. Flow Sensor	L	Neonatal flow sensor failure.
Failure		1. Replace neonatal flow sensor
ranure		2. Contact your service personnel.
		Serial cable of neonatal flow sensor is not connected. Sampling line of
No Neo. Flow	M	neonatal flow sensor is not connected.
Sensor	IVI	Check the connection of the neonatal flow sensor cable and sampling
		line.
Wrong Neo. Flow	L	Adult proximal flow sensor is used.
Sensor Type		Use neonatal flow sensor.
Neo. Flow Sensor	M	Neonatal flow sensor monitor off in the volume mode.
Monitoring Off		Neonatal flow sensor monitor on.