**SIEMENS**

CE 0197

## HEMASTIX® Reagent Strips

### Test for blood in urine.

**INTENDED USE:**

HEMASTIX Reagent Strips are intended for self-testing by persons and for healthcare professionals in Near Patient/Point of Care settings to detect the presence of blood in urine.<sup>1-4</sup> The reagent strips are for *in vitro* diagnostic use (IVD). Non-Automated. The strips may be read visually. Siemens Strips for urinalysis provide results in semi-quantitative formats. This test is carried out for human urine samples to aid in the assessment of conditions such as kidney disease.

**[!]** It is important to read the entire HEMAStix insert carefully before you begin testing. The directions must be followed exactly. Accurate timing is essential to provide optimal results. The HEMAStix Reagent Strip is ready to use upon removal from the bottle and the Reagent Strip is disposable. No additional laboratory equipment is necessary for testing.

**SPECIMEN COLLECTION AND PREPARATION:**

Wash and rinse your hands with an appropriate cleanser. Collect urine in a clean dry container and test it as soon as possible. If testing cannot be done within an hour after voiding, refrigerate the specimen immediately. Let it return room temperature before testing. Do not let the urine stand at room temperature for an extended period of time before testing or microbial growth could occur and cause false HEMAStix results. Your work area and specimen container should be free of detergents and other substances that may interfere with your HEMAStix test result.

**CAUTION: Ensure that work areas and specimen containers are always free of detergents and other contaminating substances. Some substances can interfere with patient results.**

**DEVICE FOR NEAR-PATIENT TESTING****DEVICE FOR SELF-TESTING**

For Professionals: Do not centrifuge. The use of urine preservatives is not recommended.

**HEMASTIX STORAGE AND HANDLING:**

**[!] DO NOT REUSE:** Each test strip is for single use only.  
**[!] Do not use HEMAStix reagent strips after expiration date. Record the opening date on the bottle label.**

**[!] Store at room temperature between 15–30°C (59–86°F).**

**[!] Do not store the bottle in direct sunlight.**

**[!] IMPORTANT: Protection against ambient moisture, light and heat is essential to guard against altered reagent reaction.**

The Reagent Strips must be kept in the bottle with the cap tightly closed. Do not remove desiccant from bottle. Replace cap immediately and tightly after removing Reagent Strip. Do not touch test pad on the Reagent Strip. The strips should be disposed off after use in the household waste. If visible blood is present in the urine, carefully dispose off the test materials as per your local jurisdictional guidelines.


**[!] WARNING: HEMAStix Reagent Strips are for *in vitro* IVD diagnostic use only. DO NOT SWALLOW.**

**DIRECTIONS FOR TESTING:**

**Material Provided:** Urine Strips  
**Special Materials Required (Not Provided) - Professional use:** Clean dry container, Gloves, Stop Watch/Timer, QC: e.g., Chek-Stix® (REF: 10325654 (1364E)), Paper Towel or Tissue.

**Special Materials Required (Not Provided) - Self-Test use:** Clean dry container, Gloves, Stop Watch/Timer, Paper Towel or Tissue.

Step 1. Check the Expiration Date. If the date has passed, discard and test with strips from a new bottle. Record the opening date on the label. **Use of Reagent Strips beyond the expiration date may yield inaccurate results.** Remove the Reagent Strip from the bottle. Replace the cap immediately and tightly. **Do not touch the pad on the Reagent Strip.**

Step 2. Dip the test end of the Reagent Strip into a fresh, well-mixed, uncentrifuged urine sample and remove it immediately, drawing the edge of the Reagent Strip against the rim of the urine container to remove excess urine (or pass the end of the Reagent Strip through a sample of urine.) **Immediately begin timing.**

Step 3. Wait 60 seconds after removing from the urine, and then match the test pad to the color chart on the bottle. Step 4. Record the results.

**SPECIFIC PERFORMANCE CHARACTERISTICS:**

Specific performance characteristics are based on clinical and analytical studies. In clinical specimens, the sensitivity depends upon several factors: the variability of color perception; the presence or absence of inhibitory factors typically found in urine, the specific gravity, and the pH; and the lighting conditions under which the product is read.

RESULTS: Results with HEMAStix Reagent Strips are obtained in clinically meaningful units directly from the Color Chart comparison. The color blocks represent nominal values; actual values will vary around the nominal values. It is important to dip test pad in urine completely, but briefly, to avoid dissolving out reagents. The test pad may read between 1 and 2 minutes for identifying negative specimens. Read the results carefully within the specified time, in a good light and with the test pad held near the Color Chart on the bottle label.

**IMPORTANT: HEMAStix reagent strip results should never be used as the sole basis for adjusting medication dosage. Always consult your physician or healthcare professional before adjusting your medication.**

**EXPECTED VALUES FOR BLOOD [BLO]:** The test generally detects 0.015–0.062 mg/dL free hemoglobin; however, because of the inherent variability of clinical urines, lesser concentrations may be detected under certain conditions. The sensitivity of the test is reduced in urines with high specific gravity. The test is equally sensitive to myoglobin as to hemoglobin. The appearance of green spots on the reacted reagent pads indicates the presence of intact erythrocytes in urine. The color chart includes examples of trace and moderate non-hemolyzed color blocks.

Reactions ranging from trace to large, with proportionately more numerous spots, may be observed. (Hemoglobin concentration of 0.015–0.062 mg/dL is approximately equivalent to 5–20 intact red blood cells per microliter.)

**CHEMICAL PRINCIPLES OF THE PROCEDURE:** The test is based on the peroxidase-like activity of hemoglobin that catalyzes the reaction of diisopropylbenzoyldihydroperoxide and 3,3',5,5'-tetramethylbenzidine. The resulting color ranges from orange through green; very high levels of blood may cause the color change to continue to blue.

**REAGENTS:** (Based on dry weight at time of impregnation) 6.8% w/w diisopropylbenzoyldihydroperoxide; 4.0% w/w 3,3',5,5'-tetramethylbenzidine; 48.0% w/w buffer; 41.2% w/w nonreactive ingredients.

**BIBLIOGRAPHY:** 1. Elder, T. R. *et al.*: Simplified Test for Hematuria. *Ann. Intern. Med.* 56: 957; 1962.  
 2. Adams, E. *et al.*: Hemolysis in Hematuria. *J. Urol.* 88:427; 1962.  
 3. Leonards, J. R.: Simple test for hematuria compared with established test. *JAMA* 179:807; 1962.  
 4. Henry, J.B. (ed.): Clinical Diagnosis and Management by Laboratory Methods, 19th ed. Philadelphia: Saunders; 1996; pp. 164, 411–456.  
 5. Newall, R.G. and Howell R.: Clinical Urinalysis, The Principles and Practice of Urine Testing in the Hospital and Community. Buckinghamshire, UK: Bayer Corporation; 1990; pp. 25–30.

**TEST INFORMATION:** BLOOD [BLO]: Normally, no haemoglobin is detectable in urine (<100 µg/L or 0.010 mg/dL; 3 RBC/µL). The significance of the trace reaction may vary among patients, and clinical judgment is required for assessment in an individual case. Blood is often, but not always, found in the urine of menstruating females. The test is equally sensitive to myoglobin as to haemoglobin. A haemoglobin concentration of 150–620 µg/L (0.015–0.062 mg/dL) is approximately equivalent to 5–20 intact red blood cells per microliter.

**LIMITATIONS:** All other trademarks are the property of their respective owners.

**PRODUCT NO.:** 10317353 (2816A).

**TECHNICAL ASSISTANCE:** Safety Data Sheets (MSDS/SDS) are available on [siemens-healthineers.com/poc](http://siemens-healthineers.com/poc). For customer support, contact your local technical support provider. [siemens-healthineers.com/poc](http://siemens-healthineers.com/poc).

The summary of safety and performance for this *in vitro* diagnostic medical device is available to the public in the European database on medical devices (EU-DAMED) when this database is available and the information has been uploaded by the Notified Body. The web address of the EU-DAMED public website is: <https://ec.europa.eu/tools/eudamed>.

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

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**SYMBOLS GLOSSARY:**

Symbol	Definition
CE 0197	CE Mark (+notified body ID#)
ΣXX	Contains sufficient for (n) tests
REF	Catalog number
UDI	Unique Device Identification
LOT	Batch code
EC REP	Manufacturer
	Authorized Representative in the European Community
	Indicates changes to the previous revision

© 2011–2023 Siemens Healthcare Diagnostics. All rights reserved.

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