ADVANTAGE PLUS™

Endoscope Reprocessing System





THE MOST ADVANCED AUTOMATED ENDOSCOPE

REPROCESSOR FROM THE MARKET LEADER

REPROCESS A VARIETY OF ENDOSCOPES

- Compatible with all flexible endoscopes
- Validated for reprocessing TEE probes, dilators, cystoscopes and semi-critical rigid scopes

COMPREHENSIVE ENDOSCOPE CHANNEL MONITORING

- Guarantees proper hookup connection
- Monitors endoscope channels for blockage
- Ensures all channels are fully reprocessed

HIGH PERFORMANCE CLEANING CYCLE

- Automated cycle eliminates variability of manual cleaning
- Optimizes endoscope HLD process
- Documented wash cycle

AUTOMATED LEAK TESTING AND ALCOHOL FLUSH

- Enables disinfection of an endoscope prior to repair
- Minimizes leak-related repair costs
- Facilitates drying of endoscope channels

NETWORKING AND REMOTE DIAGNOSTICS

- Network multiple systems together for easy administration and reporting
- Connect to local or off-site server for centralized data collection and convenient back-ups
- Fast technical support via on-line, remote diagnostics
- Integrates with CANEXIS[™] Integrated Workflow Solution

RAPICIDE™ PA HIGH-LEVEL DISINFECTANT

- Proven endoscope compatibility
- Faster endoscope availability
- Reduced worker exposure; environmentally friendly
- Low-temp disinfection and Aldehyde-free





THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, **only Cantel delivers the Complete Circle of Protection**, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

REPROCESS High-level disinfection is the cornerstone of infection prevention. Reprocessing technologies from Cantel are designed to ensure patient safety by adhering to the strictest of standards for high-level disinfection, optimise workflow efficiency, and deliver versatility by supporting a wide range of scopes.

ADVANTAGE PLUS™

Endoscope Reprocessing System | Reprocess





ADVANTAGE PLUS™ ENDOSCOPE REPROCESSING SYSTEM ORDERING INFORMATION - NORTH AMERICA

MODEL NUMBER	DESCRIPTION
ADV-1007	ADVANTAGE PLUS Endoscope Reprocessing System for RAPICIDE™ PA Disinfectant with air compressor
ADV-1008	ADVANTAGE PLUS Endoscope Reprocessing System for RAPICIDE PA Disinfectant without air compressor
CAS-1000	Endoscope Transport Cassette

ADVANTAGE PLUS ENDOSCOPE REPROCESSING SYSTEM ORDERING INFORMATION - INTERNATIONAL

MODEL NUMBER	DESCRIPTION
ADV-2001	ADVANTAGE PLUS Endoscope Reprocessing System for RAPICIDE PA Disinfectant with air compressor
ADV-2002	ADVANTAGE PLUS Endoscope Reprocessing System for RAPICIDE PA Disinfectant without air compressor

Dimensions

 $59 \text{ H} \times 45 \text{ W} \times 36 \text{ D}$ inches (157 H x 114 W x 92 D cm) Lid Height: Open - 54″ (137 cm), Closed - 38″ (97 cm) Depth with air compressor: 37.5″ (95.cm)

Weight

400 lbs (182 kg)

Electrical Requirements

100-240 VAC ± 10%: single phase, 50/60HZ; 1000 Watts

Electrical Safety Certifications

ETL Intertek Listed #75811

UL Standard 61010-1

EN Standard 61010-1

CSA/CAN Standard C22.2, No. 61010-1

Water Supply

Minimum flow rate: 1.32 GPM (5 LPM) at 35 to 87 PSI (2.4 to 6 bar)

 $Water\ Temperature:$

(35°C ± 2°C) - North America; (25°C ± 2°C) - International

Water Hardness: less than 12 gpg (200 ppm)

Pre-filter wall space: 26 H \times 38 W \times 6 D inch (66 H \times 96 W \times 15 D cm)

Basin Size

19 W x 17 L x 7 D inch (48 H x 43 L x 18 D cm)



ADVANTAGE PLUSTM and RAPICIDETM are trademarks of Medivators Inc. CANEXISTM is a trademark of Cantel Medical Corp.

MEDIVATORS IS NOW CANTEL

www.cantelmedical.eu

Medivators BV | p: +31 45 5 471 471 (customer service)

