## DECLARATION DE CONFORMITE CE

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---

## DECLARATION OF EC CONFORMITY

*We, ELITech Clinical Systems SAS, Zone Industrielle 61500 SEES France, hereby certify, under our own responsibility, that the reagents such as listed attached (2 pages), conform to the essential requirements of appendices I and III of European Directive 98/79/EC, relating to in vitro diagnostic medical devices and to the public health code.*

*These devices are classified in the "other device" category since they do not belong neither to list A or list B of annex II nor to self-testing class.*

*This declaration is based on the contents of each technical file and is supported by the certification of our quality system according to the standard NF EN ISO 13485 : 2016 (Certification valid until July 27<sup>th</sup>, 2023).*

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**ELITECH CLINICAL SYSTEMS SAS**  
**Zone Industrielle**  
**61500 SEES FRANCE**

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*Design, production, control and sales of clinical chemistry products intended to be used for in vitro diagnostics. Validation of the combination reagents and analyzers. Distribution of clinical chemistry analyzers and products for in vitro diagnostics.*

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**ELITech Clinical Systems SAS**  
**Zone industrielle - 61500 SEES - FRA**

**est conforme aux exigences des normes internationales**  
*complies with the requirements of the international standards*

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GMED N° 10462-7

Ce certificat est délivré selon les règles de certification GMED / This certificate is issued according to the rules of GMED certification

Renouvelle le certificat 10462-6

**On behalf of the President**  
**Lionel DREUX**  
**Certification Director**

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☞ **BIBLIOGRAFIE**

1. Wu, A.H.B., Clinical guide to laboratory tests, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 66.
2. Dufour D. R., The Liver: Function and Chemical Pathology. Clinical Chemistry: Theory, Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 586 and appendix.
3. Johnson, A. M., Amino Acids and Proteins. Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 286.
4. Doumas, B. T., et al., Albumin standards and the measurement of serum albumin with bromocresol green, Clin Chim Acta (1971), **31**, 87.
5. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
6. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
7. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
8. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
9. Interference Testing in Clinical Chemistry ; Approved Guideline—Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25**(27).
10. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
11. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> edition, AACC Press (1997).
12. Young D.S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> edition, AACC Press (1995).

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8. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).

9. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).

10. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

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

12. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.

13. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).

14. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995)

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**BIBLIOGRAFIE**

- Panteghini, M., Bais, R., *Enzyme, Tietz Fundamentals of Clinical Chemistry*, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317
- Dufour, D.R., *The Pancreas: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 651, and appendix.
- Winn-Deen, E.S., David H., Sigler E. et Chavez R., *Development of a direct assay for  $\alpha$ -amylase*, *Clin. Chem.*, (1988), **34** (10), 2005
- Guder W.G., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma, and serum samples*, *World Health Organization*, WHO/DIL/LAB/99.1 Rev.2, (2002).
- Schumann, G., et al. *IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37°C*, *Clin Chem Lab Med*, (2006), **44** (9), 1146.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI (NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).

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### BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzyme, Tietz Fundamentals of Clinical Chemistry*, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4<sup>th</sup> Ed., (W.B. Saunders Company), (2006), 154
3. Dufour, R., *The liver: Function and chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 586, appendix
4. Ward, M.K., Cockayne, S., *Enzymology. Clinical Chemistry : Concepts and Application*, Anderson, S.C., Cockayne, S. (W.B. Saunders eds. Philadelphia USA), (1993), 238.
5. Schumann, G., et al *Clin Chem Lab Med.*, (2002), **40**, 725.
6. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations*, WHO/DIL/LAB/99.1 Rev.2, (2002), 21.
7. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
8. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
9. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
10. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

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12. Vassault A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686

13. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry ana-lysis of blood samples containing monoclonal immuno-globulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.

14. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).

15. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).

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### BIBLIOGRAFIE

- Higgins, T., et al., *Hemoglobin, Iron, and Bilirubin*, Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds, Philadelphia USA), (2008), 509.
- Dufour, D.R. *The liver: Function and chemical Pathology*, Clinical Chemistry: Theory, Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc.), (2010), 586 and appendix.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
- Wu, H.B., *General Clinical Tests*. Tietz Clinical guide to laboratory tests, 4<sup>th</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 172.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).8.
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
- Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> edition, AACC Press (1997).
- Young D.S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> edition, AACC Press (1995).

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### ☞ BIBLIOGRAFIE

1. Wu, A.H.B., Tietz Clinical guide to laboratory test, 4<sup>th</sup> Ed., (W.B. Saunders Company), (2006), 684.
2. Itani, O., Tsang, R.C., Bone disease, Clinical Chemistry: Theory Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 614 and appendix.
3. Foley, K. F., Boccuzzi, L., Urine Calcium: Laboratory Measurement and Clinical Utility, Labmedicine, (2010), **41**, 683.
4. Bauer, P. J., Anal. Biochem., (1981), 110, 61.
5. Guder, W.G., et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Endres, D.B., Rude, R. K., Disorders of Bone, Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 711.
7. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
8. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
9. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
10. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
11. Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25**(27).
12. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
13. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> edition, AACC Press (1997).
14. Young D.S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> edition, AACC Press (1995).

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



Stabilitatea la bord: 7 zile

**BIBLIOGRAFIE**

- Henderson, A.R., Donald W.M., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 5<sup>th</sup> Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 352.
- Sanhai, W.R., Christenson, R.H., *Cardiac and muscle disease. Clinical Chemistry: Theory, Analysis, Correlation*, 4<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 566 and appendix.
- Neumeier, D., *et al.*, *Clin.Chim.Acta.*, (1976), **73**, 445
- Schumann, G., *et al.*, *Clin. Chem. Lab. Med.*, (2002), **40**, 635.
- Klein, G., *et al.*, *Clin. Chem.*, (2001), **47**, Suppl. A30.
- Tietz N.W., *Clinical guide to laboratory tests*. 3<sup>rd</sup> Ed., W.B. Saunders eds. Philadelphia USA, (1995), 180.
- Vassault, A., *et al.*, *Ann. Biol. Clin.*, (1986), **44**, 686.
- Vassault, A., *et al.*, *Ann. Biol. Clin.*, (1999), **57**, 685.
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995)

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FTRO-CRCO-v23(12/2018)\_PIT-CRCO-4-v23

### BIBLIOGRAFIE

1. Allston, C.A., *Non protein nitrogenous compounds and renal function*. Clinical Chemistry: Concepts and Application, Anderson, S.C., Cockayne, S. (W.B. Saunders eds. Philadelphia USA), (1993), 369.
2. Newman, D.J., Price C.P., *Non protein nitrogen metabolite*. Tietz Fundamentals of Clinical Chemistry, 5<sup>th</sup> Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 414.
3. Butler, A.R., The Jaffe reaction. *Identification of the coloured species*. Clin. Chim. Acta., (1975), 59, 227.
4. Vasiliades, J., *Reaction of alkaline picrate with creatinine. 1. Kinetics and mechanism of formation of the mono-creatinine picric acid complex*. Clin. Chem., (1976), **22**, 1664.
5. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
6. Tietz, N.W. Clinical guide to laboratory tests, 3<sup>rd</sup> Ed, (W.B. Saunders eds. Philadelphia USA), (1995), 186.
7. Vassault A., *et al.*, Ann. Biol. Clin., (1986), 44, 686.
8. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
9. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> edition, AACC Press (1997).
10. Young, D.S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> edition, AACC Press (1995).

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## BIBLIOGRAFIE

- Schreiber, W.E., *Iron and Porphyrin Metabolism. Clinical Chemistry: Theory Analysis. Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 755 and appendix.
- Higgins, T., Beutler, M.D. and Doumas, B.T., *Hemoglobin, Iron and Bilirubin, Tietz Fundamentals of Clinical Chemistry*, 6<sup>th</sup> Ed, Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds. Philadelphia USA), (2008), 509.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
- Wu, H.B., *General Clinical Tests. Tietz Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 634-639.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI (NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25**(27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> edition, AACC Press (1997).
- Young D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> edition, AACC Press (1995).

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

**BIBLIOGRAFIE**

1. Wu, A.H.B., Tietz Clinical guide to laboratory test, 4<sup>th</sup> Ed., (W.B. Saunders Company), (2006), 648.
2. Sanhai, W.R., Eloff, B.C., Christenson, R.H., Cardiac and Muscle Disease, Clinical Chemistry: Theory Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 677, appendix.
3. Panteghini, M., Bais, R., Enzyme, Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
4. Schumann, G. *et al.* Clin Chem Lab Med., (2002), **40**, 643.
5. Guder, W.G., *et al.*, Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. (2002). WHO/DIL/LAB/99.1 Rev.2
6. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
7. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
8. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
9. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
10. Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
11. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
12. Young, D. S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> Ed., AACC Press, (1997).
13. Young, D. S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> Ed., AACC Press, (1995).

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#### BIBLIOGRAFIE

1. Dufour, D.R., *The Pancreas: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 651, and appendix.
2. Panteghini, M., Bais, R., *Enzyme, Tietz Fundamentals of Clinical Chemistry*, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
3. Tietz NW et al. *Lipase in serum - the elusive enzyme: An overview. Clin Chem* (1993); **39**:746-756.
4. Neumann U et al. *New substrates for the optical determination of lipase. EP 207252* (1987).
5. Panteghini, M., *The never-ending search of an acceptable compromise for pancreatic lipase standardisation, Clin Chem Lab Med* 2012; **50(3)**: 419-421
6. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
7. Wu, H.B., *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 634.
8. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI (NCCLS) document EP6-A (2003), **23** (16).
9. *Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
10. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
11. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
12. *Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
13. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263
14. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).
15. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).
16. Glick MR, Ryder KW, Jackson SA. *Graphical comparisons of Interferences in Clinical Chemistry instrumentation. Clin Chem* 1986;**32**:471-475.

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
FTRO-PASL-v21(12/2018)\_PIT-PASL-4-v21

**BIBLIOGRAFIE**

- Henderson, A.R., Donald W.M., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 5<sup>th</sup> Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 352.
- Scherwin, J.E., *Liver fonction. Clinical Chemistry: Theory, Analysis, Correlation*, 4<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 492 and appendix.
- German Society for Clinical Chemistry, *Z. Klin. Chem. Klin. Biochem.*, (1972), **10**, 281.
- Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology, *Scand. J. Clin. Lab. Invest.*, (1974), **33**, 291.
- Guder, World Health Organization, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*, WHO/ DIL/LAB/99.1 Rev.2, (2002), 21.
- Rosalki, S.B., *Clin. Chem.*, (1993), **37**, 648.
- Vassault A., *et al.*, *Ann. Biol. Clin.*, (1986), **44**, 686.
- Vassault A., *et al.*, *Ann. Biol. Clin.*, (1999), **57**, 685.
- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immuno-globulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).

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**BIBLIOGRAFIE**

- Naito, H.K., *Coronary Artery Disease and Disorders of Lipid Metabolism. Clinical Chemistry: Theory, Analysis, Correlation*, 4<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C. (Mosby, Inc. eds. St Louis USA), (2003), 603.
- Expert Panel on Detection, *Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III)*, Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP). *JAMA*, (2001), **285**, 2486.
- Fossati, P., Prencipe, L., *Serum triglycerides determined colorimetrically with an enzyme that produces hydrogen peroxide. Clin. Chem.*, (1982), **28**, 2077.
- Tietz, N.W., *Clinical guide to laboratory tests*, 3<sup>rd</sup> Ed. (W.B. Saunders eds. Philadelphia USA), (1995), 610.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), **23** (16).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline-Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry ; Approved Guideline-Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> edition, AACC Press (1997).
- Young D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> edition, AACC Press (1995).

Referințe:	Compoziția trusei:
TGML-0425	R 6 x 50 mL
TGML-0515	R 6 x 100 mL
TGML-0700	R 4 x 250 mL
TGML-0427	R 6 x 50 mL + Std 1 x 5 mL
TGML-0497	R 1 x 100 mL + Std 1 x 5 mL
TGML-0517	R 6 x 100 mL + Std 1 x 5 mL
TGML-0707	R 4 x 250 mL + Std 1 x 5 mL

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## English - EN

### ❖INTENDED USE

ELITech Clinical Systems URIC ACID MONO SL is an *in vitro* diagnostic reagent intended for the quantitative determination of Uric Acid in human serum, plasma and urine samples on analyzers or semi-automatic analyzers.

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### ❖CLINICAL SIGNIFICANCE <sup>(1-3)</sup>

Uric acid is the major product of the catabolism of endogenous and exogenous purines (adenosine and guanosine). Uric acid is poorly soluble in water. Therefore, in case of high serum concentration, urate crystals can form and deposit in joints, triggering painful inflammations (gout), or can damage kidneys. Increased serum uric acid level can be caused either by increased production (increased intake of purines, increased nucleic acid turn-over particularly in case of certain cancers or following anti-cancer treatments, genetic metabolic disorders such as Lesch-Nyhan syndrome, psoriasis) or by decreased excretion (renal failure, drugs such as diuretics). In the case of pre-eclampsia, serum uric acid can be increased due to both mechanisms. Decreased serum uric acid is more uncommon. It can occur for example in impaired renal elimination such as in Fanconi syndrome or in Hodgkin's disease.



When uric acid concentration is abnormally high in urines, there is a risk of stone formation. In practice, the determination of uric acid in serum is indicated to help diagnose gout, for the follow-up of patients under radiation or chemotherapy, and sometimes to help diagnose pre-eclampsia. The determination of uric acid in urines is indicated to help diagnose gout and in the prevention of recurring kidney stones.

### ❖LIMITATION OF USE

The quantitative assay of Uric Acid alone can not be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

### ❖METHOD & PRINCIPLE <sup>(4)</sup>

Enzymatic / PAP - End Point.

EHSPT = N-Ethyl-N-(2-Hydroxy-3-Sulfopropyl) *m*-Toluidine
4-AAP = Amino-4-antipirine

### ❖COMPOSITION

<b>Reagent: R</b>			
Buffer pH 7.00 (20-25°C)			
EHSPT	0.72	mmol/L	
Amino-4-antipirine	0.37	mmol/L	
Uricase	≥ 150	U/L	
Peroxidase	≥ 12 000	U/L	
Sodium azide	< 0.1	% (w/w)	
<b>Standard: Std (Ref<span> </span>: AUML-0427/0497/0507/0707)</b>			
Uric acid	6	mg/dL	
	357	µmol/L	
Sodium azide	<	0.5	% (w/w)

### MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II

- Normal saline solution (NaCl 9 g/L).
- Analyzers or semi-automatic analyzers.
- General Laboratory equipment (e.g. pipette).
- Do not use materials that are not required as indicated above.

### ❖PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.

- Reagent R and Standard Std contain sodium azide which may react with lead or copper plating to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.

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**On board stability :** The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

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### PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

### ❖SAMPLES

### Specimen <sup>(2)</sup>

- Serum
- Plasma (lithium heparin)
- Urine
- Using any other specimen type should be validated by the laboratory.

### Warnings and precautions

- To prevent urate precipitation, urine samples may be adjusted to pH>8.0 with NaOH.<sup>(2)</sup>

- Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

### Storage and stability <sup>(2)</sup>

*Serum/ plasma*
- 3-5 days at 2-8°C
- 6 months at -20°C
*Urine (Alkalinized)*
- 3 days at room temperature
Do not refrigerate urine samples

### ❖REFERENCE VALUES <sup>(1)</sup>

<i>Serum/plasma</i>	mg/dL	µmol/L	
Men	3.5 - 7.2	208 - 428	
Women	2.6 - 6.0	155 - 357	

<i>Urine (24 h collection)</i>	mg/24h	mmol/24	
	250 - 750	1.48 - 4.43	
	mg/dL*	mmol/L*	
	16.7 - 50	0.99 - 2.97	

\*for a urinary volume of 1.5 L per 24 hours

With a purine-free diet, excretion may decrease by 20 to 25%

*Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

### ❖PROCEDURE

#### Manual Procedure

*Wavelength :* 546 nm
*Optical path :* 1 cm
*Sample/ Reagent ratio :* 1:40
*Temperature :* 37 °C

*Urine samples must be diluted 1:10 with NaCl 9 g/L solution before measurement.*

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#### Automatic Procedure

These reagents may be used on several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

Urine samples must be diluted 1:10 with NaCl 9 g/L solution before measurement. For users of Selectra TouchPro software, urine dilution is performed automatically.

### ❖CALCULATION

<sup>A</sup> Sample	x n	n = Calibrator/ standard
A Calibrator/ Standard		concentration

For the calculation of uric acid concentration in urine, multiply the result by the dilution factor (10). For users of Selectra TouchPro software, the results take the dilution factor into account.

**Conversion factor :** mg/dL x 59.48 = µmol/L
mg/dL x 0.059 = mmol/L

### CALIBRATION

*For the reference AUML-0427/0497/0507/0707 :* ELICAL 2 and Uric Acid Standard 6 mg/dL are traceable to ID-MS (Isotope Dilution - Mass Spectrometry) reference method.

*For the reference AUML-0250/0420/0500 :* ELICAL 2 is traceable to ID-MS (Isotope Dilution - Mass Spectrometry) reference method.

**Calibration frequency :** The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

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Controls have to be performed :
- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements.
Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

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### - Measuring range

*a) Serum/Plasma*
1.50 - 25.00 mg/dL (89 - 1487 µmol/L).
Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 125.00 mg/dL (7436 µmol/L).
Do not report results outside this extended range.

*b) Urine*
5.0 - 250.0 mg/dL (0.30 - 14.87 mmol/L).
Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 800.0 mg/dL (47.59 mmol/L).
Do not report results outside this extended range.

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### Limit of Detection (LoD) and Limit of Quantification (LoQ)

*a) Serum/Plasma*
LoD = 0.09 mg/dL (5 µmol/L)
LoQ = 1.00 mg/dL (59 µmol/L)

*b) Urine*
LoD = 0.6 mg/dL (0.04 mmol/L)
LoQ = 5.0 mg/dL (0.30 mmol/L)

### - Precision

*Serum/Plasma*

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented in the following table.

		Mean	Within-run	Total	
	n	mg/dL	µmol/L	CV (%)	
<b>Low level</b>	80	2.41	143	0.5	2.8
<b>Medium level</b>	80	4.95	294	0.7	2.3
<b>High level</b>	80	6.86	408	0.7	2.2

### b) Urine

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented in the following table.

		Mean	Within-run	Total	
	n	mg/dL	mmol/L	CV (%)	
<b>Low level</b>	80	10.3	0.61	1.8	6.6
<b>Medium level</b>	80	23.9	1.42	1.1	3.8
<b>High level</b>	80	77.9	4.63	1.2	3.3

### - Correlation

*a) Serum/Plasma*

A comparative study has been performed between URIC ACID MONO SL reagent on a Selectra ProM analyzer and a similar commercially available system on 100 human serum samples. The sample concentrations ranged from 1.55 and 23.94 mg/dL (92 and 1424 µmol/L).

The results are as follows :
Correlation coefficient : (r) = 0.999
Linear regression: y = 1.044 x - 0.04 mg/dL (2 µmol/L)

### b) Urine

A comparative study has been performed between URIC ACID MONO SL reagent on a Selectra ProM analyzer and a similar commercially available system on 49 human urine samples.

The sample concentrations ranged from 5.6 and 220.2 mg/dL (0.33 and 13.10 mmol/L).
The results are as follows :
Correlation coefficient : (r) = 0.996
Linear regression: y = 1.061 x + 0.1 mg/dL (0.01 mmol/L)

### - Limitations/Interferences

- Do not use visibly turbid or hemolyzed samples.

#### a) Serum/Plasma

- Studies have been performed to determine the level of interference from different compounds. The following uric acid levels were tested: 2.52 mg/dL and 7.56 mg/dL.

No significant interference is defined by a recovery ≤±10% of the initial value.

**Unconjugated bilirubin:** No significant interference up to 30.0 mg/dL (513 µmol/L).

**Conjugated bilirubin:** No significant interference up to 14.8 mg/dL (253 µmol/L).

**Hemoglobin:** No significant interference up to 50 mg/dL.

**Glucose:** No significant interference up to 500 mg/dL (27.8 mmol/L).

**Triglycerides:** No significant interference up to 2095 mg/dL (23.7 mmol/L).

**Turbidity:** Interference occurs at all levels of Intralipid®
**Ascorbic acid:** Significant interference on samples containing ascorbic acid.

**Methyl-dopa:** No significant interference up to 1 mg/dL.

**CalciumDOBESILATE:** Induces falsely low results on individuals taking calciumDOBESILATE.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.<sup>(6)</sup>

- Results can be falsely lowered by significant levels in the sample of NAC (*N*-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metanzolite.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(6-7)</sup>

#### b) Urine

- Studies have been performed to determine the level of interference from different compounds.

The following uric acid levels were tested: 10.0 mg/dL and 75.0 mg/dL.

No significant interference is defined by a recovery ≤±10% of the initial value.

**Conjugated bilirubin:** No significant interference up to 29.5 mg/dL (505 µmol/L).

**Hemoglobin:** No significant interference up to 300 mg/dL.

**Urea:** No significant interference up to 5000 mg/dL (833 mmol/L).

**Ascorbic acid:** No significant interference up to 20 mg/dL.

**Methyl-dopa:** Induces falsely high results at therapeutic concentrations.

- Results can be falsely lowered by significant levels in the sample of NAC (*N*-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metanzolite.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(6-7)</sup>

### - On board stability/Calibration frequency

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### ❖DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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#### ☛BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFÍA/BIBLIOGRAFIA

- Rifai, N., Warnick, G.R., Remaley, A.T., *Lipids, lipoproteins, apolipoproteins and other cardiovascular risk factors. Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 402.
- Burnett, J.R., *Coronary Artery Disease: Lipid metabolism. Clinical Chemistry: Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 691 and appendix.

- Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), *JAMA*, (2001), **285**, 2486.

- Ahmadraji, T. & Killard, A., *Analytical Methods*, (2013), **5**, 3612.

- Wu, A.H.B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 564.

- Langlois, M.R., *et al.* for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, *Clin. Chem. Lab. Med.*, (2020), **58**, 496.

- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.

- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCC Press, (1997).

- Young, D.S. *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCC Press, (1995).

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- Symbols used are defined on ISO 15223-1 standard, except those presented below.

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❖**CLINICAL SIGNIFICANCE** <sup>(1-3)</sup>

Cholesterol in the serum is derived from dietary sources or is synthesized endogenously, mainly in hepatic and intestinal cells. Cholesterol is an important structural component of cell and organelle membranes. It is a precursor of bile acids, vitamin D and steroid hormones. Cholesterol, being insoluble in water, circulates in association with lipoproteins (HDL, LDL, VLDL and chylomicrons).

HDL transports cholesterol from the cells into the liver where it can be catabolized and eliminated. Low HDL levels have thus been shown to be a risk factor for coronary heart diseases in otherwise healthy people. In practice, HDL cholesterol measurement is requested to assess predisposition of patients to cardiovascular risk as part of a lipid profile. HDL Cholesterol measurment is also important to help diagnose hyperlipoproteinaemias.

❖**LIMITATION OF USE**

In patients where abnormal lipoprotein is observed, the HDL cholesterol results may be discrepant between homogeneous measurement methods. <sup>(4)</sup> The quantitative assay of HDL cholesterol alone can not be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient’s medical history.

❖**METHOD & PRINCIPLE** <sup>(4)</sup>

Phosphate complex inhibition / PAP – End Point.

1<sup>st</sup> step :

While the sample is mixed with the reagent, surfactant and organic and inorganic phosphorus compounds block non-HDL cholesterol.

2<sup>nd</sup> step:

Cholesterol from HDL is released by the action of cholesterol esterase, then undergoes a cascade of enzymatic reactions ending with a Trinder-type reaction peroxidase/ aminoantipyrine:

    
*CHE*
→ free cholesterol + fatty acids
    
*CHO*
→ Cholest-4-en-3-one+H<sub>2</sub>O<sub>2</sub>

    
*Peroxidase*
→Quinone (blue-violet compound)

CHE : cholesterol esterase

CHO : cholesterol oxidase

4-AA : 4-aminoantipyrine

HDAOS : N-(2-Hydroxy-3-sulfopropil)-3,5-dimethoxyaniline, sodium salt

**COMPOSITION**

**Reagent 1: R1**

Buffer, pH 6.85

N-(2-Hydroxy-3-sulfopropil)-3,5-dimethoxyaniline sodium salt (HDAOS) ≥ 0,5 mmol/L

Peroxidase ≥ 8 000 U/L

Inorganic and organic phosphorous compounds

**Reagent 2 : R2**

Buffer, pH 8.15

Cholesterol Oxidase ≥ 2 000 U/L

Cholesterol Esterase ≥ 180 U/L

Peroxidase ≥ 15 000 U/L

4-aminoantipyrine ≥ 2 mmol/L

Surfactant 0,6 %

Sodium azide < 0,1 % (p/p)

Also contains ascorbate oxidase for optimal performance.

**MATERIALS REQUIRED BUT NOT PROVIDED**

- HLCA-0041 HDL LDL CALIBRATOR

- CONT-0060 ELITROL I

- CONT-0160 ELITROL II

- Analyzers or semi-automatic analyzers.

- General Laboratory equipment (e.g. pipette).

- Do not use materials that are not required as indicated above.

**PRECAUTIONS FOR USE AND WARNINGS**

- Reagent R2 contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

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- Use clean or single use laboratory equipment only to avoid contamination.

- Do not interchange reagent vials from different kits. - Consult Safety Data Sheet (SDS) for a proper handling.

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- Serum

- Plasma (lithium heparin)

- Using any other specimen type should be validated by the laboratory.

❖**Warnings and precautions**

- A non-fasting or fasting sample is suitable for the determination of a lipid profile. A repeat of lipid profile on a fasting sample could be performed in cases where a nonfasting Triglyceride result > 400 mg/dL (4.5 mmol/L) or hypertriglyceridemia is known <sup>(6)</sup>

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**Storage and stability** <sup>(1,5)</sup>

- 7 days at 2-8°C

- For longer storage, freeze at -70°C or lower (freeze only once).

❖**REFERENCE VALUES** <sup>(6)</sup>

Most recent publications recommend reporting non-HDL cholesterol in addition to HDL cholesterol and adapting decision-making thresholds for treatments as part of an overall risk assessment. At the laboratory level, the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) recommends that the following concentrations be reported as abnormal :

<b>HDL Cholesterol<span> </span>:</b>			
<i>Serum/plasma</i>	mg/dL	mmol/L	
Men	≤ 40	≤ 1.0	
Women	≤ 45	≤ 1.2	

<b>Non-HDL Cholesterol<span> </span>:</b>			
<i>Serum/plasma</i>	mg/dL	mmol/L	
Non-fasting sample	≥ 150	≥ 3.9	
Fasting sample	≥ 145	≥ 3.8	

***Note :*** *Laboratories should follow locally applicable recommendations for lipid targets if they differ from those reported above.*

**PROCEDURE**

*Manual Procedure*

*Wavelength :* 578 nm
*Optical path :* 1 cm
*Sample/ Reagent ratio :* 1:100
*Temperature:* 37 °C
*Read against distilled water.*

	<b>CALIBRATION</b>	<b>TEST</b>
<b>Reagent R1</b>	900 µL	900 µL
<b>Calibrator</b>	12 µL	-
<b>Sample</b>	-	12 µL

*Mix and read the absorbances (A1) after an incubation of 5 minutes.*

*Then add :*

<b>Reagent R2</b>	300 µL	300 µL
-------------------	--------	--------

*Mix and read the absorbances (A2) after an incubation of 2.5 minutes*

*Automatic Procedure*

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request.

For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

- **Several reagents may be contaminated by HDL CHOLESTEROL reagent.**

- **On Selectra ProM and ProXL : In order to avoid contamination on these instruments, program the following incompatibilities :**

<b>Software</b>	<b>Menu</b>	<b>Parameter</b>
<b>TouchPro</b>	Probe incompatibilities	Link / HDL CHOLESTEROL –Acid Solution
<b>Other</b>	Needle incompatibility	HDL CHOLESTEROL <<HCl

- **For other Selectra instruments: repeat any absurd results after programming a needle wash.**

**CALCULATION**

(A2 - A1 x Fdil) Sample x n

(A2 - A1 x Fdil) Calibrator

n = calibrator concentration

Fdil = Dilution factor = (R1 volume + sample volume) / (R1 volume + R2 volume + sample volume)

*Conversion factor :* mg/dL x 0.0259 = mmol/L

**CALIBRATION**

HDL LDL CALIBRATOR is traceable to the reference method approved by the CDC (Centers for Disease Control and Prevention).

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- Using any other specimen type should be validated by the laboratory.

Controls have to be performed :

- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements.
Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

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**PERFORMANCES**

Performances were obtained on Selectra Pro, following CLSI technical recommendations, under controlled environmental conditions.

**- Measuring range**

5-200 mg/dL (0.13 - 5.17 mmol/L)

Do not report results outside the measuring range.

**Limit of Detection (LoD) and Limit of Quantification (LoQ)**
LoD = 0 mg/dL (0.00 mmol/L)
LoQ = 5 mg/dL (0.13 mmol/L)

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**- Correlation**

A comparative study has been performed between HDL CHOLESTEROL reagent on a Selectra Pro analyzer and a similar commercially available system on 106 human serum samples.

The sample concentrations ranged from 11 to 185 mg/dL (0.28 - 4.78 mmol/L). The results are as follows :
Correlation coefficient : (r) = 0.998
Linear regression: y = 1.031 x – 4 mg/dL (0.10 mmol/L).

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The following HDL cholesterol levels were tested:
31 mg/dL and 54 mg/dL.
No significant interference is defined by a recovery ≤ ±4.2 mg/dL of the initial value at HDL Cholesterol concentration of 31 mg/dL (or 13.6%) and ≤ ±10 % of initial value at HDL Cholesterol concentration of 54 mg/dL
*Unconjugated bilirubin*: No significant interference up to 30.0 mg/dL (513 µmol/L).
*Conjugated bilirubin*: No significant interference up to 29.5 mg/dL (505 µmol/L).
*Hemoglobin*: No significant interference up to 500 mg/dL
*Turbidity*: No significant interference up to 600 mg/dL (6.8 mmol/L) triglyceride equivalent.
*Methyl dopa*: No significant interference up to 2.0 mg/dL.
*Acetaminophen*: No significant interference up to 30 mg/dL.
*Acetylsalicylic acid*: No significant interference up to 20 mg/dL (1190 µmol/L).
*Uric acid*: No significant interference up to 20 mg/dL (1190 µmol/L).
*Ascorbic acid*: No significant interference up to 19.8 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström’s macroglobulinemia) can cause unreliable results.<sup>(7)</sup>

- Results can be falsely lowered by significant levels in the sample of NAC (*N*-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metanzole.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(8,9)</sup>

**- On board stability/Calibration frequency**

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*These performances have been obtained using ELITech Selectra Pro analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.*

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6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature.* Acta Clin Belg., (2004), 59, 263.

7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests.* 2<sup>nd</sup> Ed., AACC Press, (1997).

8. Young D.S., *Effects of drugs on clinical laboratory tests.* 4<sup>th</sup> Ed., AACC Press, (1995).

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- Only for ref. CHSL-0250/0455, used with Selectra TouchPro software.  
**- see § PROCEDURE:** Contamination risk

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2. Burnett, J.R., *Coronary Artery Disease: Lipid metabolism. Clinical Chemistry: Theory, Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., (Mostly Inc. eds), (2010), 691 and appendix.

3. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), *JAMA*, (2001), **285**, 2486.

4. Wu, A.H.B., *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 684.

5. Ahmadraji, T. & Killard, A., *Analytical Methods*, (2013), **5**, 3612.

6. Langlois, M.R., *et al.* for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, *Clin. Chem. Lab. Med.*, (2020), **58**, 496.

7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.

8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACCPress, (1997).

9. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACCPress, (1995).

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## BIBLIOGRAPHIE/BIBLIOGRAPHY

### BIBLIOGRAFIA/BIBLIOGRAFIA

1. Sacks, D.B., *Carbohydrates*. Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A. Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 373.
2. Dods, R.F., *Diabetes Mellitus*. Clinical Chemistry - Theory, Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 729 and appendix
3. Wu, A. H. B., *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 444.
4. World Health Organization (WHO), *Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycemia*, (2006).
5. Neeley, W.E., *Clin. Chem.*, (1972), **18**, 509.
6. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev 2.
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*. *Acta Clin Belg.*, (2004), **59**, 263.
8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACCPress, (1997).
9. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACCPress, (1995).

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- Only for ref. GHSL-0250, used with Selectra TouchPro software.

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

### MÉTHODE & PRINCIPE (5)

Hexokinase - Point final.

G-6-PDH = Glucose-6-phosphate-déshydrogenase

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4. World Health Organization (WHO), *Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycemia*, (2006).

5. Trinder, P., *Ann. Clin. Biochem.*, (1969), **6**, 24.

6. Gunder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.

7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.

8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*. 2<sup>nd</sup> Ed., AAC Press, (1997).

9. Young, D.S., *Effects of drugs on clinical laboratory tests*. 4<sup>th</sup> Ed., AAC Press, (1995).

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## English - EN

### ☛INTENDED USE

ELITech Clinical Systems GLUCOSE PAP SL is an *in vitro* diagnostic reagent intended for the quantitative determination of glucose in human serum and plasma samples on analyzers or semi-automatic analyzers. The standard is intended for the calibration of the reagent.

These *in vitro* diagnostic devices are for professional use only.

### ☛CLINICAL SIGNIFICANCE<sup>(1-3)</sup>

Glucose is the main source of energy for the human body. Glucose is converted either into glycogen or into triglycerides to be stored. Glucose blood level is mainly regulated by two antagonist hormones: insulin and glucagon.

Glycemia disorders appear mostly in type I or type II diabetes as well as in gestational diabetes. They can be also associated to various endocrinal, pancreatic or hepatic disorders, or linked to drugs.

In normal health condition, glucose is filtered then reabsorbed by kidneys and is therefore not present in urine. Elevated concentrations in urine are observed when the blood concentration is high or in case of impaired tubular reabsorption. Glucose measurement in blood is indicated for diabetes, in screening, diagnosis or follow-up of patients. It is also indicated to monitor patients with symptoms of hyperglycemia or hypoglycemia.

### ☛LIMITATION OF USE

For diabetes assessment, the collection conditions and interpretation of serum glucose concentrations should follow local recommendations such as those published by WHO.<sup>(4)</sup>

The quantitative assay of glucose alone can not be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

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## BIBLIOGRAPHIE/BIBLIOGRAPHY

### BIBLIOGRAFIA/BIBLIOGRAFIA

- Endres, D. B. & Rude, R.K., *Disorders of Bone*. Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 711.
- Itani O. & Tsang, C., *Bone Disease. Clinical Chemistry: Theory, Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 614 and appendix.
- Swaminathan, R., *Clin. Biochem. Rev.*, (2003), 24, 47.
- Liamis, G. et al., Am. J. Nephrol., (2013), 38, 50.
- Mann, C.K, Yoe, J.H., *Spectrophotometric determination of magnesium with sodium 1-azo-2-hydroxy-3-(2,4-dimethylcarboxanilido)-naphthalene-1-(2-hydroxybenzene-5-sulfonate)*, *Anal. Chem.*, (1956), 28, 202.
- Guder W.G., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
- Wu, A.H.B. *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 706.
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *(Acta Clin Belg.*, (2004), 59, 263.

9. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACCC Press (1997).

10. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACCC Press (1995)

## SYMBOLES/SYMBOLS/ SIMBOLOS/SIMBOLS

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## English - EN

### INTENDED USE

ELITech Clinical Systems MAGNESIUM XB is an *in vitro* diagnostic reagent intended for the quantitative determination of magnesium in human serum, plasma and urine samples on analyzers or semi-automatic analyzers.

The standard is intended for the calibration of reagent. These *in vitro* diagnostic devices are for professional use only.

#### CLINICAL SIGNIFICANCE <sup>(1-3)</sup>

Magnesium is the fourth most abundant cation in the human body. It is the cofactor for many enzymatic systems including ATP-dependent enzymes and it also plays an active role in bone mineral homeostasis. Less than 1% of the total magnesium of the human body is carried by blood where it is mainly found as free ion, but also protein-bound (mainly albumin) or complexed with various anions. Magnesemia measures total magnesium but only free magnesium is biologically active. Hypomagnesemia can be due to important renal losses (excessive intake of diuretics), gastrointestinal disorders (malabsorption, diarrheas), endocrinal disorders (less often due to inadequate intake. Hypomagnesemia is often associated with hypokalemia and/or hypocalcemia. Hypermagnesemia is rare and is caused almost always by excessive intake (parenteral therapy, excessive intake of magnesium-containing drugs) or renal insufficiency.

The measure of urinary magnesium allows, in presence of magnesium deficiency, to identify a kidney etiology. In clinical practice, measurement of magnesium is indicated to help diagnose the causes of abnormal calcium and/or potassium levels or symptoms suggesting hypo or hypermagnesemia as well as for monitoring magnesium- and/or calcium-based treatments.

#### LIMITATION OF USE

Because hypoalbuminemia may trigger pseudo-hypomagnesemia, serum magnesium results must be interpreted in regards of serum total protein and/or albumin levels. <sup>(4)</sup>

The quantitative assay of MAGNESIUM XB alone can not be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

#### METHOD & PRINCIPLE <sup>(5)</sup>

Xylydyl Blue - End Point

Xylydyl blue in the reagent combines with the magnesium from the sample to form a red-purple chelate. EGTA is used to complex calcium and thus prevents it from interfering with the test.

The simultaneous increase in absorbance at 505-510 nm and decrease of the 620-630 nm absorbance are proportional to the magnesium concentration in the sample.

#### COMPOSITION

**Reagent: R**  
AMP buffer, pH 11.2  
Xylydyl Blue 120 µmol/L  
Sodium azide < 0.1 %(w/w)  
Also contains EGTA and surfactants for optimal performance

AMP: 2-Amino-2-methyl-1-propanol

**Standard: Std**  
Magnesium 2.0 mg/dL  
823 µmol/L

#### MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2  
- CONT-0060 ELITROL I  
- CONT-0160 ELITROL II  
- Normal saline solution (NaCl 9 g/L).  
- Analyzers or semi-automatic analyzers.  
- General Laboratory equipment (e.g. pipette).  
- Do not use materials that are not required as indicated above.

#### PRECAUTIONS FOR USE AND WARNINGS

- Reagent R contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.  
- Consult Safety Data Sheet (SDS) for a proper handling.

#### STABILITY

**Store at 2-8 °C and protect from light. Do not freeze.** Do not use after expiration dates indicated on the vial labels.

The standard should be immediately and tightly capped to prevent contamination and evaporation.

*On board stability:*

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

#### PREPARATION

The reagent and standard are ready to use.

#### PRODUCT DETERIORATION

- These products should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

#### SAMPLES

**Specimen** <sup>(1,6)</sup>

- Serum.  
- Plasma (lithium heparin).  
- Urine.  
- Using any other specimen type should be validated by the laboratory.

**Warnings and precautions** <sup>(1,2)</sup>

- Serum must be separated from cells as rapidly as possible. <sup>(1,3)</sup>

- Samples must be free from hemolysis <sup>(1,2)</sup>

- After collection, urine specimens should be acidified with hydrochloric acid 6N to a pH < 2 to prevent magnesium salt precipitation. <sup>(1)</sup>

Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

**Storage and stability** <sup>(6)</sup>

Serum/plasma  
Serum/plasma  
- 7 days at room temperature.  
- 7 days at 2-8°C.  
- 1 year at -20°C.  
Urine (acidified)  
- 3 days at room temperature  
- 3 days at 2-8°C  
- 1 year at -20°C

#### REFERENCE VALUES <sup>(7)</sup>

<i>Serum/plasma</i>	mg/dL	mmol/L
	1.5 - 2.6	0.63 - 1.05

<i>Urine (24h - collection)</i>	mg/24h	mmol/24 h
	73 - 122	3.0 - 5.0

	mg/dL*	mmol/L*
	4.9 - 8.1	2.0 - 3.3

\* for a urinary volume of 1.5 L per 24 hours

*Note* : *The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

#### PROCEDURE

*Manual Procedure*

*Wavelength:* 505 - 625 nm  
*Optical path:* 1 cm  
*Sample/ Reagent ratio:* 1:100  
*Temperature:* 37 °C  
*Urine samples must be diluted 1:5 with NaCl 9 g/L urine before measurement.*  
*Read against reagent blank/distilled water.*

	CALIBRATION	TEST
<b>Reagent R</b>	2 ml	2 ml
<b>Standard/ Calibrator</b>	20 µl	-
<b>Sample</b>	-	20 µl

*Mix and read the absorbances (A) after an incubation of 5 minutes.*

*Automatic Procedure*

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

Urine samples must be diluted 1:5 with NaCl 9 g/L urine before measurement. For users of Selectra TouchPro software, urine dilution is performed automatically.

**Important set-up information:**  
**MAGNESIUM XB reagent can be weakly contaminated by CHOLESTEROL SL on Selectra ProM and ProXL.**  
**In order to avoid contamination on these instruments, program the following incompatibilities:**

<i>Software</i>	<i>Menu</i>	<i>Parameter</i>
<b>TouchPro</b>	Probe incompatibilities	incompatibility/ CHOLESTEROL - MAGNESIUM
<b>Other</b>	Needle incompatibility	CHOLESTEROL: MAGNESIUM

#### CALCULATION

<sup>ΔA</sup> Sample x n = Calibrator/ standard  
<sup>ΔA</sup> Calibrator/ Standard concentration

For the calculation of magnesium concentration in urine, multiply the result by the dilution factor (5). For users of Selectra TouchPro software, the results take the dilution factor into account.

**Conversion factor :** mg/dL x 0.41 = mmol/L

#### CALIBRATION

For the reference MAGX-0600 : ELICAL 2 and Magnesium XL standard are traceable to the atomic absorption reference method  
For the reference MAGX-0250 : ELICAL 2 is traceable to the atomic absorption reference method.

*Calibration frequency* : The calibration is specific for each analyzer. § PERFORMANCE DATA).

#### QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :  
- prior to assaying patient samples,  
- at least once per day,  
- after every calibration,  
- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

#### WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

#### PERFORMANCE DATA

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

#### Measuring range

*a) Serum/Plasma*

0.20 - 5.00 mg/dL (0.08 - 2.06 mmol/L)

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 25 mg/dL (10.28 mmol/L).  
Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

*b) Urine*

1.0 - 20.0 mg/dL (0.41-8.23 mmol/L)

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 100 mg/dL (41.14 mmol/L).  
Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

**Limit of Detection (LOD) and Limit of Quantification (LoQ)**

*a) Serum/Plasma*

LoD = 0.03 mg/dL (0.01 mmol/L)  
LoQ = 0.20 mg/dL (0.08 mmol/L)

*b) Urine*

LoD = 0.1 mg/dL (0.04 mmol/L)  
LoQ = 1.0 mg/dL (0.41 mmol/L)

#### Precision

Imprecision data has been obtained on 2 Selectra Pro analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below.

*a) Serum/Plasma*

	Mean	Within-run	Total		
	n	mg/dL	mmol/L	CV (%)	
<b>Low level</b>	80	1.55	0.64	1.0	3.5
<b>Medium level</b>	80	2.51	1.03	1.3	3.7
<b>High level</b>	80	3.84	1.58	0.7	3.0

*b) Urine*

	Mean	Within-run	Total		
	n	mg/dL	mmol/L	CV (%)	
<b>Low level</b>	80	1.2	0.49	2.9	9.2
<b>Medium level</b>	80	5.0	2.06	0.7	4.0
<b>High level</b>	80	16.0	6.58	0.8	3.3

#### Correlation

*a) Serum/Plasma*

A comparative study has been performed between MAGNESIUM XB reagent on a Selectra ProM analyzer and a similar commercially available system on 102 human serum samples.

The sample concentrations ranged from 0.27 to 4.99 mg/dL (0.11 - 2.05 mmol/L).  
The results are as follows :  
Correlation coefficient : (r) = 0.996  
Linear regression: y = 1.041x - 0.01 mg/dL (0.00 mmol/L).

*b) Urine*

A comparative study has been performed between MAGNESIUM XB reagent on a Selectra ProM analyzer and a similar commercially available system on 80 human urine samples.  
The sample concentrations ranged from 1.7 to 19.6 mg/dL (0.70 - 8.06 mmol/L).  
The results are as follows :  
Correlation coefficient : (r) = 0.999  
Linear regression: y = 1.011x + 0.1 (0.00 mmol/L).

#### Limitations/Interferences

*a) Serum/Plasma*

- Studies have been performed to determine the level of interference from different compounds.  
The following analyte levels were tested: 1.50 mg/dL and 3.90 mg/dL

No significant interference is defined by a recovery  $\pm$ 10% of the initial value.

**Triglycerides** : No significant interference up to 3000 mg/dL(33.90 mmol/L).  
**Unconjugated bilirubin** : No significant interference up to 30.0 mg/dL (513 µmol/L).  
**Conjugated bilirubin** : No significant interference up to 29.5 mg/dL (504 µmol/L)  
**Hemoglobin** : No significant interference up to 500 mg/dL  
**Calcium** : No significant interference up to 20.0 mg/dL (4.99 mmol/L)  
**Ascorbic acid**: No significant interference up to 30 mg/dL.  
**Acetaminophen** : No significant interference up to 30 mg/dL.

**Acetylsalicylic acid** : No significant interference up to 200 mg/dL

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results. <sup>(8)</sup>

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. <sup>(9-10)</sup>

*b) Urine*

- Studies have been performed to determine the level of interference from different compounds.  
The following magnesium levels were tested: 1.2 mg /dL and 15 mg/dL.

No significant interference is defined by a recovery  $\pm$ 10% of the initial value.  
**Conjugated bilirubin** : Aucune interférence significative jusqu'à 29.5 mg/dL (504 µmol/L)  
**Hemoglobin** : No significant interference up to 500 mg/dL.  
**Calcium** : No significant interference up to 60.0 mg/dL (14.97 mmol/L)  
**Uric Acid** : No significant interference up to100 mg/dL (5.95 mmol/L).  
**Urea** : No significant interference up to 5000 mg/dL (832.50 mmol/L).  
**Ascorbic acid** : No significant interference up to 19.8 mg/dL.  
**pH**: No significant interference for pH values ranging between 2.5 à 6.0

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. <sup>(9-10)</sup>

**- On board stability/Calibration frequency**  
**On Board Stability:** 14 days  
**Calibration frequency:** 7 days

Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

*These performances have been obtained using ELITech Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.*

#### DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

#### TECHNICAL ASSISTANCE

Contact your local distributor or ELITech Clinical Systems SAS (CCsupport@elitechgroup.com).

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**BIBLIOGRAFIE**

1. Wu, A.H.B., Tietz Clinical guide to laboratory test, 4<sup>th</sup> Ed., (W.B. Saunders Company), (2006), 470
2. Panteghini, M., Bais, R., Enzyme, Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
3. Dufour, R., The liver function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 586, appendix.
4. Schumann, G., et al., IFCC Primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 6. Reference procedure for the measurement of catalytic concentration of  $\gamma$ -Glutamyltransferase. Clin. Chem. Lab. Med., (2002), **40**, 734.
5. Guder, W.G., et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples, WHO/DIL/LAB/99.1 Rev.2, (2002).
6. Kytzia H-J, Reference intervals for GGT according to the new IFCC 37°C reference procedure, Clin. Chem. Lab. Med., (2005), **43**, A69.
7. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
8. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
9. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
10. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

11. Interference Testing in Clinical Chemistry ; Approved Guideline—Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).

12. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.

13. Young, D. S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> Ed., AACC Press, (1997).

14. Young, D. S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> Ed., AACC Press, (1995).

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4. Bretauière, J.P., *et al., Clin. Chem.*, (1976), **22**, 1614.

5. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.

6. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACCPress (1997).

7. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACCPress (1995).

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- Only for ref. **URSL-0250/0455**, used with Selectra TouchPro software.

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**BIBLIOGRAFIE**

1. Scherwin, J.E., *Liver function. Clinical Chemistry: Theory, Analysis, Correlation*, 4<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 492 and appendix.
2. Tietz, N.W., *Clinical guide to laboratory tests*, 3<sup>rd</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 518.
3. Doumas, B.T., et al., *Clin. Chem.*, (1981), **27**, 1642.
4. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
5. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI (NCCLS) document EP6-A (2003), **23** (16).
6. *Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
7. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
8. *Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
9. *Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
10. Vassault, A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686.
11. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
12. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).
13. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).

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**BIBLIOGRAFIE**

- Panteghini, M. & Bais R., *Enzymes. Tietz Fundamentals of Clinical Chemistry*, 6th Ed, Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2008), 317
- Sanhai, W.R. & Christenson R.F., *Protein Isoforms: Isoenzymes and Isoforms. Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 677 and appendix.
- Apple, F. S., *Cardiac function, Tietz Fundamentals of Clinical Chemistry*, 6th Ed, Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2008), 614.
- Boot, E., de Haan, L., *Massive increase in serum creatine kinase during olanzapine and quetiapine treatment, not during treatment with clozapine. Psychopharmacology (Berl)* (2000);150 :347-348.
- Schumann, G., et al. *Clin Chem Lab Med.*, (2002), **40**, 635-42.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline.* CLSI (NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline.* CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition.* CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition.* CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).

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Les Hautes Rives  
02160, Maizy, France

# CALCIUM CPC method

Reagent for quantitative determination of calcium  
in human plasma or serum, or urines.

REF 80004 R1 1 x 200 mL R2 1 x 200 mL R3 1 x 10 mL

## TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

Fax: (33) 03 23 256 256

IVD IN VITRO DIAGNOSTIC USE

## CLINICAL SIGNIFICANCE (1) (2)

Total Calcium exists in 3 physiochemical states in plasma, of which approximately 50 % is free or ionised calcium, 40 % is bound to plasma proteins, and 10 % are bound with small anions.

The level of serum calcium may be affected by intestinal malabsorption, by alterations in plasma proteins level, especially albumin, which should be measured concurrently with calcium.

Hypercalcemia is found in hyperparathyroidism, multiple myeloma, bone and parathyroid neoplasms and in states with bones demineralisation.

Hypocalcemia is encountered in hypoparathyroidism and in several cases of necrosis and acute pancreatitis.

## PRINCIPLE (4)

Moorehead and Briggs derived CPC (O-Cresol Phtalein Complexone) method allows to determinate total Calcium concentration in serum, plasma or urines.

In alkaline solution CPC reacts with calcium to form a dark-red coloured complex which absorbance measured at 570 nm is proportional to the amount of calcium in the specimen.

## REAGENTS

**Vial R1** **BUFFER** (Concentration in the Test)

Amino-2-methyl-2-propanol-1 1.70 mol/L  
pH 11.0 at 20°C  
Hydrochloric acid 210 mmol/L

Xi, R36/38: Irritating to eyes and skins

R52/53: Harmful to aquatic organisms. May cause long term effects to aquatic environment

S36/37/39: Wear suitable protective clothing, gloves and eye/face protection

**Vial R2** **CHROMOGEN** (Concentration in the Test)

O-cresolphthalein complexone (CPC) 78 µmol/L  
Hydroxy-8-Quinoline 3.36 mmol/L  
Hydrochloric acid 25 mmol/L

**Vial R3** **STANDARD**

Calcium 10 mg/dL (2.5 mmol/L)

## REAGENTS PREPARATION

Working reagent preparation:

In a carefully cleaned container with HCl 0.1 N and well rinsed with demineralised water, mix 1 volume of vial R1 and 1 volume of vial R2.

Reagents may also be added separately.

## SAFETY CAUTIONS

BIOLABO reagents are designated for professional, in vitro diagnostic use.

- Verify the integrity of the contents before use.
- Use adequate protections (overall, gloves, glasses).
- Do not pipette by mouth.
- In case of contact with skin and eyes, thoroughly wash affected areas with plenty of water and seek medical advice.
- Material Safety Data Sheet is available upon request.
- Waste disposal: Respect legislation in force in the country.

All specimens should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions. Respect legislation in force in the country.

## STABILITY AND STORAGE

**Store at 18-25°C, well cap in the original vial and away from light.**

- **Standard (vial R3):** Transfer requested quantity, recap and store at 18-25°C.
- Reagents are stable until expiry date stated on the label of the kit when stored and used as described in the insert and free from contamination.
- Working reagent is stable for 1 day when free from contamination.
- Discard reagent if cloudy or if reagent blank measured at 570 nm is > 0.400.
- Don't use working reagent after expiry date stated on the label of the Kit.

## SPECIMEN COLLECTION AND HANDLING (1)

Serum or heparinised plasma. Do not use citrate, oxalate or EDTA. Blood obtained on fasting patient with minimal venous occlusion and without exercise or after restoring circulation at least for 1 minute.

24 h Urines: Acidify before assay with 20 to 30 mL HCl 6 N to dissolve calcium salts.

Mix well then dilute (1 + 2) with demineralised water before performing the test.

Calcium is stable in serum for:

- at least 7 days at 2-8°C.
- 6 months at -20°C.

Plastic and glass may adsorbed calcium during storage, especially with diluted solution.

Long-term freezing may lead to associated evaporation, lyophilisation or co precipitation with fibrin (i.e. heparinised plasma) or lipids.

## INTERFERENCES (3)

Haemolysis, icterus, lipemia, paraproteins: Perform bichromatic or multi-wavelengths analysis or specimen blank to reduce positive or negative interferences. Magnesium interference is reduced by Hydroxy-8-quinoline.

Handle with care specimens, calibrators and controls to avoid contamination by environmental calcium. Use disposable tubes and cuvettes and glassware cleaned with HCl 0.1 N, well rinsed with demineralised water.

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

## MATERIAL REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Normal and pathological control sera.
3. EDTA Solution 10 mM to realize serum blank

## CALIBRATION (6)

- Standard enclosed in the kit (vial R3) or BIOLABO Multicalibrator [REF] 95015 traceable to SRM 909b.
- Or any calibrator traceable to a reference method or material.

The calibration frequency depends on proper instrument functions and on the preservation of the reagent.

It is recommended to calibrate in the following cases:

1. When changing batch of reagent.
2. After maintenance operations on the instrument.
3. If control values are out of range, even after using a new vial of fresh serum.

## QUALITY CONTROL

- BIOLABO EXATROL-N Level I [REF] 95010.
- BIOLABO EXATROL-P Level II [REF] 95011.
- Assayed control sera referring to the same method.
- External quality control program.

It is recommended to control in the following cases:

- At least once a run.
  - At least once within 24 hours.
  - When changing vial of reagent.
  - After maintenance operations on the instrument.
- If control is out of range, apply following actions:
1. Repeat the test with the same control.
  2. If control is still out of range, prepare a fresh control serum and repeat the test.
  3. If control is still out of range, use a new vial of calibrator or a fresh calibrator and repeat the test.
  4. If control is still out of range, calibrate with a new vial of reagent.
  5. If control is still out of range, please contact BIOLABO technical support or your local Agent.

## EXPECTED VALUES (2)

### CALCIUM, total

#### in serum:

Population	mg/dL	mmol/L
Premature	6.2-11.0	[1.55-2.75]
0-10 days	7.6-10.4	[1.90-2.60]
10 days –24 months	9.0-11.0	[2.25-2.75]
24 months –12 years	8.8-10.8	[2.20-2.70]
12 years -18 years	8.4-10.2	[2.10-2.55]
18-60 years	8.6-10.0	[2.15-2.50]
60-90 years	8.8-10.2	[2.20-2.55]
> 90 years	8.2-9.6	[2.05-2.40]

**in 24 h Urines:** 100-300 mg/24 h (2.5-7.5 mmol/24 h).

Each laboratory should establish its own normal ranges for the population that it serves.

## PERFORMANCES CHARACTERISTICS

	Within run n = 20	Medium level	High level	Between run n = 20	Medium level	High level
Mean mg/dL	8.69		13.94	Mean mg/dL	10.35	13.06
S.D. mg/dL	0.06		0.15	S.D. mg/dL	0.116	0.212
C.V. %	0.7		1.1	C.V. %	1.12	1.62

Detection limit: approximately 0.16 mg/dL

Sensitivity for 1 mg/dL: 0.105 Abs. at 570 nm.

Comparison study with commercially available reagent:

$$y = 1.0249x - 0.237 \quad r = 0.9985$$

## LINEARITY (5)

The assay is linear between 7.5 and 17.5 mg/dL (1.88 to 4.37 mmol/L). Under 7.5 mg/dL, increase the specimen volume and re-calculate the result (example: if  $V = 50 \mu\text{L}$ , multiply the result by 0.51).

Above 17.5 mg/dL, dilute the specimen with demineralised water and re-assay taking into account the dilution factor to calculate the result. Linearity limit depends on specimen/reagent ratio.

## MANUAL PROCEDURE

Temperature should be held constant as the absorbance of the dye is temperature sensitive.

Pipette into well identified test tubes:	Blank	Standard	Assay
Reagent	1 mL	1 mL	1 mL
Demineralised water	25 $\mu\text{L}$		
Standard		25 $\mu\text{L}$	
Specimen			25 $\mu\text{L}$

Mix well. Incubate for 5 minutes at room temperature.  
Read absorbance at 570 nm (550-590) against reagent blank.  
The coloration is stable for 1 hour away from light.

### Notes:

- ✓ Cloudy or icteric sera: after reading absorbance according to this procedure, perform a specimen blank by adding a drop (25  $\mu\text{L}$ ) of EDTA solution 10 mM in assay and blank reagent tubes, mix and re-read. Deduct this value from the previously reading obtained for the specimen.
- ✓ Bichromatic analysis: the 2<sup>nd</sup> wavelength may be 500, 650 or 700 nm.
- ✓ Specific procedures are available upon request for automated instruments. Please contact BIOLABO technical support.

## CALCULATION

Calculate the result as follows:

**Serum or plasma:** 
$$\text{Result} = \frac{\text{Abs(Assay)}}{\text{Abs(Standard)}} \times \text{Standard concentration}$$

**Urinis diluted (1 + 2):** multiply the result by 3 (dilution factor).

## REFERENCES

- (1) TIETZ N.W. *Text book of clinical chemistry*, 3<sup>rd</sup> Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 1395-1406, p.1435-1439.
- (2) *Clinical Guide to Laboratory Test*, 4<sup>th</sup> Ed., N.W. TIETZ (2006) p. 202-207
- (3) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4<sup>th</sup> Ed. (1995) p. 3-115 to 3-124
- (4) MOOREHEAD W.R., BRIGGS H.G., *Clin. Chem.*, (1974), 20, p.1458-1460
- (5) W. L. CLARK, E L BAGINSKI, S S MARIE, et B. ZAK, *Spectrometric Study of a direct determination of Serum Calcium*, *Microchem. J.*, 20, (1975), p.22-32.
- (6) SRM:Standard Reference Material®

Manufacturer

Use by

In vitro diagnostic

Temperature limitation

Catalogue number

See insert

Batch number

Store away from light

sufficient for

dilute with

**BIOLABO**  
www.biolabo.fr

**MANUFACTURER:**  
**BIOLABO SAS,**

Les Hautes Rives  
02160, Maizy, France

# UREA Colorimetric Method

Reagent for quantitative determination of urea  
in human serum and plasma or urines.

REF 80221	R1 1 x 125 mL	R2 1 x 1.25 mL	R3 1 x 31 mL	R4 1 x 10 mL
REF 80321	R1 1 x 500 mL	R2 1 x 5 mL	R3 1 x 125 mL	R4 1 x 10 mL

## TECHNICAL SUPPORT AND ORDERS

Tel : (33) 03 23 25 15 50

Fax: (33) 03 23 256 256

**IVD** IN VITRO DIAGNOSTIC USE

## CLINICAL SIGNIFICANCE (1) (5)

More than 90% of urea is excreted through the kidneys in urines. Measurement of the plasma or serum urea concentration is widely regarded as a test of renal function. However, a number of nonrenal factors also influence the circulating urea concentration: Urea increased level occurs when proteins catabolism is accelerated, burns, stress, myocardial infarction... Urea is decreased in acute liver destruction and is accompanied with increased ammonium level. Urea level is generally studied in conjunction with creatinine level (urea/creatinine ratio) to refine post-renal or pre-renal diagnosis.

## PRINCIPLE (4)

Enzymatic and colorimetric method based on the specific action of urease which hydrolyses urea in ammonium ions and carbon dioxide. Ammonium ions then form with chloride and salicylate a blue-green complex. This coloration, proportional to urea concentration in the specimen, is measured at 600 nm.

## REAGENTS COMPOSITION

<b>Vial R1</b>	<b>SALICYLATE</b>	
Salicylate		31 mmol/L
Nitroprussiate		1.67 mmol/L
<b>Vial R2</b>	<b>UREASE</b>	
Urease		≥ 15 KUI/L
<b>Vial R3</b>	<b>BASE</b>	
Sodium hypochlorite		7 mmol/L
Sodium hydroxide		62 mmol/L
<b>Before dilution: Xi, R36/38, Irritating to eyes and skin</b>		
<b>S24-25-26-28: Avoid contact with skin and eyes. After contact with skin, rinse immediately with plenty of water. After contact with eyes, rinse immediately with plenty of water and seek medical advice</b>		
Once diluted: None		
<b>Vial R4</b>	<b>STANDARD</b>	Urea 40 mg/dL (6.66 mmol/L)

## SAFETY CAUTIONS

BIOLABO reagents are designated for professional, in vitro diagnostic use.

- Verify the integrity of the contents before use.
- Use adequate protections (overall, gloves, glasses).
- Do not pipette by mouth.
- In case of contact with skin or eyes, thoroughly wash affected areas with plenty of water and seek medical advice.
- Reagents contain sodium azide (concentration < 0.1%) which may react with copper and lead plumbing. Flush with plenty of water when disposing.
- Material Safety Data Sheet is available upon request.
- Waste disposal: Respect legislation in force in the country.

All specimens should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions. Respect legislation in force in the country.

## REAGENTS PREPARATION

Add contents of vial R2 (Urease) into vial R1 (Salicylate). Mix gently by inversion.

### Alkaline (vial R3):

Procedure n°1 and n°2 (manual): Dilute (1 + 3) with demineralised water

Procedure n°3 (manual or automatic): Ready to use

**Standard (vial R4):** transfer the requested quantity, recap and store at 2-8°C.

## STABILITY AND STORAGE

**Store at 2-8°C, well recap in the original vial and away from light.**

- Unopened, reagents are stable until expiry date stated on the label when stored and used as described in the insert.
- Once reconstituted, working reagent (R1+R2) is stable for 1 month when free from contamination.
- Once opened, Base (vial R3) diluted ¼ is stable for 3 months when free from contamination.
- Once opened, the contents of vial R4 is stable for at least 3 months when free from contamination.

Discard any reagent if cloudy or if absorbance of blank against water at 600 nm > 0.100.

Don't use working reagent or diluted contents of vial R3 after expiry date stated on the label of the kit.

## SPECIMEN COLLECTION AND HANDLING (2)

Unhemolysed serum or heparinised plasma. Avoid fluoride or ammonium as anticoagulant which interfere with the assay.

Urea is stable in serum or plasma for:

- 24 h at room temperature.
- several days at 2-8°C.
- at least 2-3 months frozen.

24h Urine: diluted (1+19) with demineralised water before assay.

Urea is stable in urines for:

- 4 days at 2-8°C.

Add antibacterial agent as Thymol to improve the stability.

## INTERFERENCES (3)

No interference of assayed substances (ascorbic acid, bilirubin, haemoglobin, triglycerides) with the test.

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

## MATERIAL REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Normal and pathological control sera

## CALIBRATION

- Standard enclosed in the kit (vial R4) or BIOLABO Multicalibrator [REF] 95015 traceable to SRM 909b.
- Or any calibrator traceable to a reference method or material.

The calibration frequency depends on proper instrument functions and on the preservation of reagent.

It is recommended to calibrate in the following cases :

- When changing batch of reagent.
- After maintenance operations on the instrument.
- If control values are out of range, even after using a new vial of fresh control serum.

## QUALITY CONTROL

- BIOLABO EXATROL-N Level I [REF] 95010.
- BIOLABO EXATROL-P Level II [REF] 95011.
- Other assayed control sera referring to the same method.
- External quality control program.

It is recommended to control in the following cases:

- At least once a run.
- At least once within 24 hours.
- When changing vial of reagent.
- After maintenance operations on the instrument.

If control is out of range, apply following actions:

- Repeat the test with the same control.
- If control is still out of range, prepare a fresh control and repeat the test.
- If control is still out of range, use a new vial of calibrator or a fresh calibrator and repeat the test.
- If control is still out of range, calibrate with a new vial of reagent.
- If control is still out of range, please contact BIOLABO technical support or your local Agent.

## EXPECTED VALUES (2)

In serum and plasma	mg/dL	[mmol/L]
In cord	45-86	[7.5-14.3]
Premature	6-54	[1.1-8.9]
< 1 year	9-41	[1.4-6.8]
Children	11-39	[1.8-6.4]
18-60 years	13-43	[2.1-7.1]
60-90 years	17-49	[2.9-8.2]
> 90 years	21-66	[3.6-11.1]

**In urines** 26-43 g/24 h [0.43-0.71 mol/24 h]

To calculate blood urea nitrogen (BUN): multiply the value of urea (mg/dL) by 0.467.

Each laboratory should establish its own normal ranges for the population that it serves.

## PERFORMANCE CHARACTERISTICS (Procedure n°1)

Within run N = 20	Normal level	High level	Between run N = 20	Normal level	High level
Mean mg/dL	40	141.6	Mean mg/dL	35	111
S.D. mg/dL	0.76	1.66	S.D. mg/dL	1.58	3.44
C.V. %	1.89	1.17	C.V. %	4.5	3.1

Detection limit: approximately 10 mg/dL.

Comparison with commercially available reagent:  $y = 0.9816 + 0.87$

$r = 0.9961$

	Sensitivity for 100 mg/dL at 600 nm
Procedure n°1	Approx. 0.400 abs
Procedure n°2	Approx. 0.800 abs
Procedure n°3	Approx. 0.700 abs

Note: Sensitivity is higher at upper wavelength and lower at inferior wavelength.

## LINEARITY

Procedure n°1 and n°3:

The reaction is linear up to 250 mg/dL (41.7 mmol/L).

Procedure n°2: The reaction is linear up to 125 mg/dL (20.9 mmol/L).

Above, dilute the specimen with saline solution and reassay taking into account dilution factor to calculate the result. Linearity limit depends on specimen/reagent ratio.

## MANUAL PROCEDURE

Let stand reagents and specimens at room temperature.

### Procedure n°1

Pipette into test tubes	Blank	Standard	Assay
Working reagent (R1+R2)	1 mL	1 mL	1 mL
Deminerilised water	5 µL		
Standard		5 µL	
Specimen (Note 1)			5 µL
Mix and wait for 4 minutes at room temperature or 2 minutes at 37°C			
Base (vial R3) diluted ¼	1 mL	1 mL	1 mL
Mix. Let stands for 8 minutes at room temperature or 5 minutes at 37°C. Read absorbance at 600 nm (590-610) against blank (Note 3). Reaction coloration is stable for 2 hours.			

### Procedure n°2

Pipette into test tubes	Blank	Standard	Assay
Working reagent (R1+R2)	1 mL	1 mL	1 mL
Deminerilised water	10 µL		
Standard		10 µL	
Specimen (Note 1)			10 µL
Mix and wait for 4 minutes at room temperature or 2 minutes at 37°C			
Base (vial R3) diluted ¼	1 mL	1 mL	1 mL
Mix. Let stands for 8 minutes at room temperature or 5 minutes at 37°C. Read absorbance at 600 nm (590-610) against blank (Note 3). Reaction coloration is stable for 2 hours.			

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Pipette into test tubes	Blanc	Etalon	Dosage
Working reagent (R1+R2)	1 mL	1 mL	1 mL
Deminerilised water	5 µL		
Standard		5 µL	
Specimen (Note 1)			5 µL
Mix and wait for 4 minutes at room temperature or 2 minutes at 37°C			
Base (vial R3) PUR	250 µL	250 µL	250 µL
Mix. Let stands for 8 minutes at room temperature or 5 minutes at 37°C. Read absorbance at 600 nm (590-610) against blank (Note 3). Reaction coloration is stable for 2 hours.			

- Serum, plasma or urines diluted (1+19) with deminerilised water.
- Specific procedures are available upon request for automated instruments. Please contact BIOLABO technical support.
- The test may be performed at **578 nm**. In this case, Procedure n°2 is linear up to 300 mg/dL.

## CALCULATION

Calculate the result as follows:

**Serum and plasma:**

$$\text{Result} = \frac{\text{Abs (Assay)}}{\text{Abs (Standard)}} \times \text{Standard concentration}$$

**Urines diluted (1+19):**

Multiply the result by 20 (dilution factor).

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## **BIOLABO LES HAUTES RIVES 02160 Maizy - France**

est conforme aux exigences de la norme suivante :  
is in compliance with the requirements of the following standard:

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*DESIGN, MANUFACTURING AND SALE OF IN VITRO DIAGNOSTIC MEDICAL DEVICES.  
TECHNICAL SUPPORT AND SUPPORT SERVICES.*

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
This certificate is valid until its expiry date unless further notice, provided that the compliance and implementation of the Quality Management System are found to be satisfactory at follow-up audits and that AB Certification contract rules are fulfilled.

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
*DESIGN, MANUFACTURING AND SALE OF IN VITRO DIAGNOSTIC MEDICAL DEVICES.  
TECHNICAL SUPPORT AND SUPPORT SERVICES*

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
**dispositifs médicaux de diagnostic in vitro relevant de la directive n°98/79/CE**  
*in vitro diagnostic medical devices covered by Directive 98/79/EC*

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*Section to be completed by the applicant*

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*Device(s) category: Reagents & Instruments for Medical Biology*

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*The name of the device(s) appears on the EC declaration(s) of conformity of the manufacturer or the authorized representative*

**Classification du(des) dispositif(s) :**

*Classification of the device(s) :*

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*device of list A annex II*

**autotest hors annexe II**

*device for self-testing not listed in annex II*

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*other device (all devices except annex II and self-testing devices)*

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*Name and address of the manufacturer or the authorized representative:*

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*Name and address of Production site (optional):*

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*I the undersigned Isabelle, Oget, Director of Regulatory Affairs declare that the information above-mentioned is correct and the in vitro diagnostic medical devices on the EC declaration(s) of conformity are CE marked under my responsibility within the meaning of the European directive n°98/79/EC and fulfil the essential requirements of health and safety.*


Date : 30/08/2018

**PARTIE RESERVEE A LA CCIR PARIS IDF**

*Section reserved for the administration*

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
*The in vitro diagnostic medical devices CE marked in conformity with the directive 98/79/EC can be placed on the French market and in the other Member states of the European Union and part of the European Free Trade Association, and be exported in the non-EC Member States. This free sale certificate is valid until the maintenance, by the manufacturer of the concerned devices, of an CE declaration of conformity (other devices) together with when appropriate, the certificates delivered by a notified body (devices of list A and B, annex II, devices for self-testing not listed in annex II). This free sale certificate can only be used for exportation outside European Union.*

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PARIS  
2016

**EC DECLARATION OF CONFORMITY**  
for *in vitro* diagnostic medical device (IVD) - Directive 98/79/EC

The undersigned Sentinel CH. SpA, having its premises in Milan, Italy - Via Robert Koch 2, manufacturer of the family of devices named as "kits for clinical chemistry, immunochemistry, coagulation, molecular biology and rapid tests for immunology and serology" declares, under its own responsibility that the device

REF: **17642H**Description: **Bilirubin Total Liquid**EDMA: **11.02.01.03**1° LOT placed on the market: **10149**

complies with all the essential requirements listed in Annex I of the 98/79/EC Directive, as prescribed in Annex III of such Directive and its Italian transposition (legislative decree nr. 332/2000).

It therefore declares and assures, under its own responsibility, that the device:

1. complies with the applicable provisions of the Directive
2. is not included in the list A and B of Annex II of the Directive
3. is designed, manufactured and placed on the market according to the company certified quality system, in compliance with the current quality regulations as expressed in Annex III of the Directive.

Furthermore, the manufacturer declares to:

1. keep and make available for the Competent Authority the product technical file, as specified in Annex III of the 98/79/CE Directive, as well as to retain the batch records for a period of at least ten years after the production date of the last lot
2. have instituted and keep up to date an adequate procedure to guarantee the market surveillance requested by the Directive.

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Date / Data

28/03/2011

**Certified Quality Management System ISO 9001:2008 - ISO 13485:2003 and ISO 13485:2003 CMDCAS**

<b>IN</b>	<b>REF 17642H</b>	<b>Colorimetric determination of total bilirubin (BIL-T) concentration in serum and plasma according to Jendrassik-Grof mod. method</b>	<b>IVD</b>
	<b>Bilirubin Total Liquid</b>		<b>REAGENT 1a: 7 x 45 mL - REAGENT 1b: 1 x 20 mL</b>
<b>STANDARD/CALIBRATOR:</b> the term refers to the standard / the calibrator			
<b>REAGENT:</b> the term refers to the single reagent <b>CONTROL:</b> the term refers to the control			

**SUMMARY**

Total and direct bilirubin levels are used for the diagnosis and control of hepatic (hepatitis, cirrhosis), haemolytic and biliary disorders. In particular, high levels of total bilirubin can be indicative of excessive haemolysis or hepatic disorders of the hemoglobin catabolism.

**PRINCIPLE**

Total bilirubin, in an acid solution containing caffeine, reacts with diazotized sulfanilic acid forming a coloured azocompound (azobilirubin) which colour intensity is proportional to the total bilirubin concentration in the sample.

**REAGENTS**

The components of the kit, stored at 2-8 °C in unopened vials, are stable up to the expiry date indicated on the package.

Components of the kit and concentration of reactive ingredients:

**REAGENT 1a**

sulfanilic acid 10 mmol/L, hydrochloric acid 0.8 mol/L, caffeine ≤ 0.6 mol/L

**REAGENT 1b**

sodium nitrite ≤ 50 mmol/L

**PREPARATION OF REAGENT SOLUTIONS**

**REAGENT 1a and REAGENT 1b:** to be used in the preparation of the **R1 Solution**.

**Preparation of R1 Solution:** add to one vial of REAGENT 1a 2.5 mL of Reagent 1b. Mix for a few minutes.

**REAGENT 1a contains urea as additive.** In the sample order setting, do **not input urea test immediately after total bilirubin** in "random access" automatic analyzers.

**STABILITY**

On Board: **30** days, if contamination is avoided.

Calibration: **14** days. Repeat the calibration at any variation in the reagent lot.

**CALIBRATION**

For the calibration, use the following material:

**Clin Chem Cal**      **REF 16550**      **4x3 mL**

Multiparametric lyophilised calibration serum. For use, follow the instructions contained in the kit.

**QUALITY CONTROL**

Use the following control materials to verify test accuracy:

**Clin Chem Control 1**      **REF 16150**      **6x5 mL**

Multiparametric lyophilised control serum. For use, follow the instructions contained in the kit.

**Clin Chem Control 2**      **REF 16250**      **6x5 mL**

Multiparametric lyophilised control serum. For use, follow the instructions contained in the kit.

**STANDARDIZATION**

Jendrassik-Grof mod.

**SAMPLE**

Serum, plasma (heparin). **Do not use emolyzed samples.** Collect samples in accordance with the NCCLS procedure reported in the bibliography.<sup>1</sup>

**STABILITY OF THE SAMPLE**

7 days at 2-8 °C or 6 months at -20 °C. It is recommended to store samples **protected** from direct light.

**WARNING AND PRECAUTIONS**

- For in vitro diagnostic use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.
- Safety Data Sheets are available at [www.sentinel diagnostics.com](http://www.sentinel diagnostics.com) or contact your local representative.

-  **CAUTION:** This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and be handled in accordance with the OSHA Standard on Bloodborne Pathogens<sup>2</sup>, Biosafety Level 2<sup>3</sup> or other appropriate biosafety practices<sup>4,5</sup> should be used for materials that contain or are suspected of containing infectious agents.

**REAGENT 1a**

**DANGER:** contains **CITRIC ACID, CAFFEINE, HYDROGEN CHLORIDE, SULFANILIC ACID**

Hazard statements:

- H314 Causes severe skin burns and eye damage.

Precautionary statements:

- P280 Wear protective gloves/protective clothing/eye protection/face protection.
- P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
- P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower.
- P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
- P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P310 Immediately call a POISON CENTER/doctor/...
- P501 Dispose of contents/container in accordance with applicable regulations.

**REAGENT 1b**

**DANGER:** contains **SODIUM NITRITE**

Hazard statements:

- H314 Causes severe skin burns and eye damage.

Precautionary statements:

- P280 Wear protective gloves/protective clothing/eye protection/face protection.
- P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
- P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower.
- P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
- P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P310 Immediately call a POISON CENTER/doctor/...
- P501 Dispose of contents/container in accordance with applicable regulations.

**INSTRUMENTATION AND MATERIALS REQUIRED BUT NOT PROVIDED**

General laboratory equipment.

**ANALYTICAL PROCEDURE**

Wavelength:	570 nm (main) 700 nm (reference)
Pathlength:	1 cm
Temperature:	37 °C
Sample/R1 Solution:	1/15
Reaction:	end point (increase)

Allow reagents to reach working temperature before using. A proportional variation of the indicated reaction volumes does not change the results.

**EXAMPLE OF ANALYTICAL PROCEDURE ON AUTOMATED INSTRUMENTS****Time 0**

Calibrator/Controls/Sample = 17 µL

R1 Solution = 250 µL

**After 600 sec**

Reading

**CALCULATION**

- Plot a calibration curve on a graph paper by tracing absorbance (y axis) according to corresponding mg/dL concentration (x axis) for each Calibrator.
- Indicate on the calibration curve the absorbance value obtained for Samples and Controls.
- Extrapolate the mg/dL value for Samples and Controls from the calibration curve.

**CONVERSION FACTOR**

BIL-T [ mg/dL ] x 17.1 = BIL-T [ µmol/L ]

**REFERENCE VALUES**

<b>Adults:</b>	0.3 - 1.2 mg/dL
<b>Newborns (0-1 day):</b>	1.4 - 8.7 mg/dL
<b>Newborns (1-2 days):</b>	3.4 - 11.5 mg/dL
<b>Newborns (3-5 days):</b>	1.5 - 12.0 mg/dL

It is recommended that each laboratory establish its own expected range. For diagnostic purposes, results obtained should always be evaluated taking into consideration the patient's history and all other clinical findings.

**PERFORMANCES AUTOMATIC ANALYZER****SENSITIVITY:**

**0.05 mg/dL.** Sensitivity was calculated on 20x3 replicates x run of normal saline and reported as the "mean zero value + 3 SD".

**INTRA-ASSAY PRECISION:**

was determined on 20 replicates on two controls (L1/L2) and on a human serum sample (L3). The results were as follows:

mg/dL	L1	L2	L3
mean	1.29	3.96	2.60
SD	0.01	0.01	0.02
CV%	1.1	0.4	0.7

**INTER-ASSAY PRECISION:**

was determined from 5x2x4 tests (day x run x rep) on each control (3 levels - L1/L2/L3). The results were as follows:

	mean	Total Imprecision		Between Days		Repeatability	
		mg/dL	SD	CV%	SD	CV%	SD
L1	1.60	0.06	3.6	0.05	3.2	0.02	1.1
L2	3.25	0.13	4.0	0.11	3.4	0.02	0.6
L3	4.97	0.20	4.1	0.20	3.9	0.02	0.5

**ACCURACY:**

this test (y) was compared with a commercially available method (x). The results were as follows:

**N = 60, r = 0.96, y = 1.01x + 0.05**

**INTERFERENCES:**

the test is not affected by the presence of ascorbic acid up to **100 mg/dL**, haemoglobin up to **50 mg/dL** and lipids up to **1000 mg/dL**. In haemolysed samples, the presence of haemoglobin can cause erroneously low bilirubin results; it is therefore recommended **not to use haemolysed samples**.

**ANALYTICAL RANGE: 0.05 - 20 mg/dL.**



Samples with values higher than **20 mg/dL** must be diluted 1:10 with normal saline and the result multiplied by 10.

**WASTE MANAGEMENT**

Reagents must be disposed of in accordance with local regulations.

**BIBLIOGRAPHY**

- NCCLS Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard - Fifth Edition (H3-A5). Wayne, PA: The National Committee for Clinical Laboratory Standards, 2003.
- US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne Pathogens.
- US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories, 5th ed. Washington, DC: US Government Printing Office, January 2007.
- World Health Organization. Laboratory Biosafety Manual, 3rd ed. Geneva: World Health Organization, 2004.
- Sewell DL, Bove KE, Callihan DR, et al. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline — Third Edition (M29-A3). Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
- Pesce, A.J., Kaplan, L.A.: "Methods in Clinical Chemistry", Mosby Ed. (1987).
- Burtis C.A., Ashwood E.R.: "Tietz Textbook of Clinical Chemistry", W.B. Saunders Company Ed. (3<sup>rd</sup> edition, 1999).
- Guder W.G.: "The Quality of Diagnostic Sample". Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and the German Society for Laboratory Medicine. (1<sup>st</sup> Edition - 2001).
- Jakobs, D.S., Kasten, Jr., B.L., DeMott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2<sup>nd</sup> Edition - 1990).

<b>Explanation of symbols</b>	
	In vitro Diagnostic Medical Device
	Catalogue number
	Batch code
	Contents of kit
	Caution, consult accompanying documents
	Consult instructions for use
	Use by (last day of the month)
	XXX Contains sufficient for <n> tests
	Temperature limitation
	Manufacturer

**Note:** changes in comparison to the previous version are indicated by a vertical bar in the text margin.

# Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 9001:2015

This is to certify that:

SENTINEL CH. S.p.A.  
Via Robert Koch, 2  
Milano (MI)  
20152  
Italy

Holds Certificate No:

**FM 504784**

and operates a Quality Management System which complies with the requirements of ISO 9001:2015 for the following scope:

The design, development, manufacture and sale of in-vitro diagnostic test kits and reagent for the determination of the clinical chemistry, immunochemistry, coagulation, molecular biology, rapid tests for immunology and serology; systems and instruments for clinical diagnostic and related servicing activities.

For and on behalf of BSI:

Andrew Launn, EMEA Systems Certification Director

Original Registration Date: 2006-07-04

Latest Revision Date: 2021-05-10

Effective Date: 2021-06-11

Expiry Date: 2024-06-10

Page: 1 of 1

003

...making excellence a habit.™

# Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016 & EN ISO 13485:2016

This is to certify that:

SENTINEL CH. S.p.A.  
Via Robert Koch, 2  
Milano (MI)  
20152  
Italy

Holds Certificate Number:

**MD 504781**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 & EN ISO 13485:2016 for the following scope:

The design, development, manufacture and sale of in-vitro diagnostic test kits and reagent for the determination of the clinical chemistry, immunochemistry, coagulation, molecular biology, rapid tests for immunology and serology; systems and instruments for clinical diagnostic and related servicing activities.

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For and on behalf of BSI:

\_\_\_\_\_  
Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2006-06-19

Latest Revision Date: 2021-06-15

Effective Date: 2021-06-19

Expiry Date: 2024-06-18

Page: 1 of 1

...making excellence a habit.™

<b>EN</b>	<b>REF 17643H</b>	<b>Colorimetric determination of direct bilirubin (BIL-D) concentration in serum and plasma according to Jendrassik-Grof mod. method</b>	<b>IVD</b>
	<b>Bilirubin Direct Liquid</b>	<b>REAGENT 1a: 4 x 45 mL - REAGENT 1b: 1 x 15 mL</b>	<b>CE</b>
<b>STANDARD/CALIBRATOR:</b> the term refers to the standard / the calibrator			
<b>REAGENT:</b> the term refers to the single reagent <b>CONTROL:</b> the term refers to the control			

**SUMMARY**

Total and direct bilirubin levels are used for the diagnosis and control of hepatic (hepatitis, cirrhosis), haemolytic and biliary disorders. In particular, high levels of direct bilirubin (or "conjugated bilirubin"), which are normally absent in blood or present in negligible quantities, are detectable in the following cases: extrahepatic biliary disorders (e.g. choledoch duct's calculus, pancreas tumour); intrahepatic biliary disorders (e.g. cirrhosis, hepatitis, liver tumours).

**PRINCIPLE**

Direct bilirubin, in acid solution, reacts with diazotized sulfanilic acid forming a coloured azocompound (azobilirubin) which colour intensity is proportional to the direct bilirubin concentration in the sample.

**REAGENTS**

The components of the kit, stored at 2-8 °C in unopened vials, are stable up to the expiry date indicated on the package.

Components of the kit and concentration of reactive ingredients:

**REAGENT 1a**

sulfanilic acid 30 mmol/L, hydrochloric acid 0.25 mol/L

**REAGENT 1b**

sodium nitrite ≤ 10 mmol/L

**PREPARATION OF REAGENT SOLUTIONS**

**REAGENT 1a and REAGENT 1b:** to be used in the preparation of the **R1 Solution**.

**Preparation of R1 Solution:** add to one vial of REAGENT 1a **2.5 mL** of REAGENT 1b. Mix for a few minutes.

**STABILITY**

On Board: **30** days, if contamination is avoided.

Calibration: **14** days. Repeat the calibration at any variation in the reagent lot.

**CALIBRATION**

For the calibration, use the following material:

**Clin Chem Cal**      **REF 16550**      **4x3 mL**

Multiparametric lyophilised calibration serum. For use, follow the instructions contained in the kit.

**QUALITY CONTROL**

Use the following control materials to verify test accuracy:

**Clin Chem Control 1**      **REF 16150**      **6x5 mL**

Multiparametric lyophilised control serum. For use, follow the instructions contained in the kit.

**Clin Chem Control 2**      **REF 16250**      **6x5 mL**

Multiparametric lyophilised control serum. For use, follow the instructions contained in the kit.

**STANDARDIZATION**

Jendrassik-Grof mod.

**SAMPLE**

Serum, plasma (heparin). **Do not use emolyzed samples.** Collect samples in accordance with the NCCLS procedure reported in the bibliography.<sup>1</sup>

**STABILITY OF THE SAMPLE**

**7** days at 2-8 °C or **6** months at -20 °C. It is recommended to store samples **protected** from direct light.

**WARNING AND PRECAUTIONS**

- For in vitro diagnostic use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.

- Safety Data Sheets are available at [www.sentinel diagnostics.com](http://www.sentinel diagnostics.com) or contact your local representative.

-  **CAUTION:** This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and be handled in accordance with the OSHA Standard on Bloodborne Pathogens<sup>2</sup>, Biosafety Level 2<sup>3</sup> or other appropriate biosafety practices<sup>4,5</sup> should be used for materials that contain or are suspected of containing infectious agents.

**REAGENT 1a****DANGER:**

contains **HYDROGEN CHLORIDE, SULFANILIC ACID**

**Hazard statements:**

- H314 Causes severe skin burns and eye damage.

**Precautionary statements:**

- P280 Wear protective gloves/protective clothing/eye protection/face protection.
- P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
- P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower.
- P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
- P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P310 Immediately call a POISON CENTER/doctor/...
- P501 Dispose of contents/container in accordance with applicable regulations.

**REAGENT 1b**

**DANGER:** contains **SODIUM HYDROXIDE**

**Hazard statements:**

- H314 Causes severe skin burns and eye damage.

**Precautionary statements:**

- P280 Wear protective gloves/protective clothing/eye protection/face protection.
- P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
- P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower.
- P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
- P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P310 Immediately call a POISON CENTER/doctor/...
- P501 Dispose of contents/container in accordance with applicable regulations.

**INSTRUMENTATION AND MATERIALS REQUIRED BUT NOT PROVIDED**

General laboratory equipment.

**ANALYTICAL PROCEDURE**

Wavelength:	570 nm (main) 700 nm (reference)
Pathlength:	1 cm
Temperature:	37 °C
Sample/R1 Solution:	1/15
Reaction:	end point (increase)

Allow reagents to reach working temperature before using. A proportional variation of the indicated reaction volumes does not change the results.

**EXAMPLE OF ANALYTICAL PROCEDURE ON AUTOMATED INSTRUMENTS****Time 0**

Calibrator/Controls/Sample = 17 µL  
R1 Solution = 250 µL

**After 600 sec**

Reading

**CALCULATION**

- Plot a calibration curve on a graph paper by tracing absorbance (y axis) according to corresponding mg/dL concentration (x axis) for each Calibrator.
- Indicate on the calibration curve the absorbance value obtained for Samples and Controls.
- Extrapolate the mg/dL value for Samples and Controls from the calibration curve.

**CONVERSION FACTOR**

BIL-D [ mg/dL ] x 17.1 = BIL-D [ µmol/L ]

**REFERENCE VALUES**

**Adults:** 0.0 - 0.2 mg/dL

It is recommended that each laboratory establish its own expected range. For diagnostic purposes, results obtained should always be evaluated taking into consideration the patient's history and all other clinical findings.

**PERFORMANCES AUTOMATIC ANALYZER****SENSITIVITY:**

**0.05 mg/dL.** Sensitivity was calculated on 20x3 replicates x run of normal saline and reported as the "mean zero value + 3 SD".

**INTRA-ASSAY PRECISION:**

was determined on 20 replicates on two controls (L1/L2) and on a human serum sample (L3). The results were as follows:

mg/dL	L1	L2	L3
mean	0.81	2.29	1.43
SD	0.01	0.01	0.01
CV%	1.0	0.4	0.4

**INTER-ASSAY PRECISION:**

was determined from 5x2x4 tests (day x run x rep) on each control (3 levels - L1/L2/L3). The results were as follows:

	mean	Total Imprecision		Between Days		Repeatability	
	mg/dL	SD	CV%	SD	CV%	SD	CV%
L1	0.81	0.02	2.5	0.02	2.1	0.01	1.0
L2	1.43	0.04	2.5	0.03	1.9	0.01	0.4
L3	2.28	0.05	2.0	0.04	1.7	0.01	0.4

**ACCURACY:**

this test (y) was compared with a commercially available method (x). The results were as follows:

**N = 60, r = 0.998, y = 0.80x + 0.13**

**INTERFERENCES:**

the test is not affected by the presence of ascorbic acid up to **10 mg/dL**, haemoglobin up to **50 mg/dL** and lipids up to **200 mg/dL**. In haemolysed samples, the presence of haemoglobin can cause erroneously low bilirubin results; it is therefore recommended **not to use haemolysed samples**.

**ANALYTICAL RANGE:**

**0.05 - 20 mg/dL.**

Samples with values higher than **20 mg/dL** must be diluted 1:10 with normal saline and the result multiplied by 10.

**WASTE MANAGEMENT**

Reagents must be disposed of in accordance with local regulations.

**BIBLIOGRAPHY**

- NCCLS Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard - Fifth Edition (H3-A5). Wayne, PA: The National Committee for Clinical Laboratory Standards, 2003.
- US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. Bloodborne Pathogens.
- US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories, 5th ed. Washington, DC: US Government Printing Office, January 2007.
- World Health Organization. Laboratory Biosafety Manual, 3<sup>rd</sup> ed. Geneva: World Health Organization, 2004.
- Sewell DL, Bove KE, Callihan DR, et al. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline — Third Edition (M29-A3). Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
- Pesce, A.J., Kaplan, L.A.: "Methods in Clinical Chemistry", Mosby Ed. (1987).
- Burtis C.A., Ashwood E.R.: "Tietz Textbook of Clinical Chemistry", W.B. Saunders Company Ed. (3<sup>rd</sup> edition, 1999).
- Guder W.G.: "The Quality of Diagnostic Sample". Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and the German Society for Laboratory Medicine. (1<sup>st</sup> Edition - 2001).
- Jakobs, D.S., Kasten, Jr., B.L., DeMott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2<sup>nd</sup> Edition - 1990).

**Explanation of symbols**

IVD

In vitro Diagnostic Medical Device

REF

Catalogue number

LOT

Batch code

Cont.

Contents of kit

Caution, consult accompanying documents

Consult instructions for use

Use by (last day of the month)

XXX

Contains sufficient for &lt;n&gt; tests

Temperature limitation

Manufacturer

**Note:** changes in comparison to the previous version are indicated by a vertical bar in the text margin.

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