# EXCELL MCDSe

ELETTROBISTURI PER CHIRURGIA MONOPOLARE E BIPOLARE  
ELECTROSURGICAL UNITS FOR MONOPOLAR AND BIPOLAR USE

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The equipments of the series Excell MCDSe System are electrosurgical units for major surgery, suitable for all operations where the monopolar and the bipolar techniques are used, as well as the monopolar mode with the Argon gas flow. The units are available in 5 different models.

The models 400 MCDSe, 350 MCDSe, 250 MCDSe, and 200 MCDSe are indicated for electrosurgery. By the contrary, the model 400/A MCDSe is suitable both for electrosurgery and electrosurgery with Argon gas. It is equipped with all the necessary functions for this type of application.

## AUTO-CHECK AND SELF-REGULATION OF OUTPUT POWERS MCDSE SYSTEM

A new system with double microprocessor *Master Sleeve System* ensures a wide range of functions and performances at a very high level.

**Auto-Check:** complete auto-control of the functioning, through:

- **Main Auto-Check**, main auto-diagnosis procedure at the switching-on and **Standard Auto-Check**, continuous auto-diagnosis during the functioning.
- **General Error Control and Output Error Control:** the auto-check systems interrupt the functioning and inform the operators through specific **Error Codes** in case of damages or errors in the functioning/switching on, or in case of problems that could provoke a sudden delivery of power, which is higher than the selected one.

### ADC System-constant power e Starting Impulse Control.

All the currents are automatically and dynamically self-regulated, according to a continuous feed-back between the unit and the tissues on which the operator intervenes, and according to the electrode used and its rapidity of sliding on the tissues themselves.

All the cutting currents are characterized by a starting pulse, which is controlled and automatically optimized to obtain the maximum efficacy. These two systems guarantee the best results in all the operating fields, by using the minimum possible power (the interventions under liquid included).

## 10 PROGRAMMABLE MEMORIES

The MCDSe models have the possibility to memorize 10 functioning programs that are selectable through a simple button. In this way, the different performances can be personalized and then reproduced.

## HIGH FREQUENCY LEAKAGE CURRENTS CONTROL CIRCUIT

A specific circuit verifies that the HF leakage currents to hearth always remain within the safety limits, in order to avoid all risks of burns, even in dangerous situations like the contact between the patient and a small part of the operating table.

## SAFETY CIRCUIT OF THE NEUTRAL ELECTRODE NPCC SYSTEM

The safety circuit of the neutral electrode guarantees its complete control.

With the *twin section electrodes Split* type, it fully checks the connection and the quality of the contact between the electrode and the patients' tissues.

If the contact is not optimum (approx. 80% of a standard electrode for adults), it gives the first luminous indication. If the contact decreases to approx. 50% of a standard electrode for adults, it gives another luminous indication, and automatically reduces the power at a safety level of 200 W. If the contact is dangerous or the electrode cable is broken / not connected at all, it interrupts the delivery of the power, through a luminous and acoustic alarm.

With the *one single section electrodes Non Split* type, it only controls if the cable is integral and correctly connected.

## PRACTICALNESS AND EASINESS

On the front panel, which is completely smooth and easily washable, there are all the control devices with selection buttons for the functions and adjustable buttons for the powers; there are also luminous activation signals or alarm signals, as well as the powers control displays. The outputs are characterized by luminous indications that light up according to the usage modes that have been selected, in order to favour the connection and the activation of the accessories.

## GREAT RELIABILITY, TECHNICAL ASSISTANCE AND TECHNOLOGICAL EVOLUTION

The Excell MCDSe models have been designed with interchangeable boards, and they are equipped with a new generator *RF Mosfet Powered Damped Oscillator* with alimentation device *switching* that does not cause any problem of heating; the test procedures also include a "stress test" of 100 hours of functioning at max. power. The Auto-Check systems memorize the last 32 **Error Codes** for the

future controls, while the functioning and the auto-check softwares can be easily updated, through a serial port for a connection to a PC.

## WIDE RANGE OF CURRENTS AND PERFORMANCES

All the Excell models are equipped with 13 currents to satisfy all possible operating needs; these currents become 14 into the model 400/A MCDSe.

**Four for the monopolar cut**, with sinusoidal wave, which can be used even with the Argon gas flow into the model 400/A MCDSe in order to reduce the smokes and the smells. They ensure the maximum rapidity and precision in all the operating fields, the procedure under liquid included.

- **Pure:** Pure cut without any coagulating effect.
- **Blend 1:** Coagulating cut with normal haemostatic effect.

The two currents are characterized by special self-regulation softwares, which are particularly useful for the TUR interventions, the vaporisation, etc.

- **Blend 2:** Coagulating cut, with strong haemostatic effect, *spray* type.
- **Endo:** Coagulating cut, with cut phases alternated to coagulation phases and a special self-regulation software for its usage in flexible endoscopy.

**Four for the monopolar coagulation**, which allow to obtain different coagulation effects.

- **Soft:** Very delicate, with soft superficial effect and strong deep action.
- **Fulg Forced:** With strong superficial and deep effect.
- **Pin Point Contact:** Similar to the previous one, but softer.
- **Spray:** Without any contact and with a very strong superficial effect.

### One for the Argon gas coagulation, into the model 400/A MCDSe

This current consists of a special coagulation without any contact, which can be obtained by combining the spray coagulation with an Argon gas flow through specific electrodes that guarantee a very rapid superficial effect, a very small deep action (3mm max.), an absolute absence of the sticking effect of the electrode on the tissues, and a minimum risk of perforation.

**Two for the bipolar cut**, which allow to use this technique with scissors, hooks for laparoscopy, etc.

- **Pure:** Pure cut with minimum coagulating effect.
- **Blend:** Coagulating cut with strong coagulating effect.

### Three for the bipolar coagulation.

- **Micro:** Very delicate, *Micro Precise* type, with a minimum sticking effect of the tissues on the tips of the forceps.
- **Micro Auto:** The same as Micro, but with *Impedance Sensing* automatic *Auto Start / Auto Stop* device.
- **Macro:** *Standard* type, which is very rapid and efficacious, ideal for forceps with bigger section (i.e. laparoscopy).

## TWO MONOPOLAR AND ONE BIPOLAR OUTPUTS

**The two monopolar outputs, which can be used at the same time by two different operators, can be activated by three modes:**

- Both of them by the hand-switch electrodes holder handles, both of them by the double pedal foot-switch or one by the hand-switches and one by the double pedal foot-switch.

Into the model 400/A MCDSe, when also the functioning with Argon gas is selected, the outputs are activated alternatively:

- Both of them by the hand-switch electrodes holder handles or both of them by the double pedal foot-switch, one for the normal electrosurgery and one for the Argon gas enhanced electrosurgery.

**The bipolar output, which can be used independently and at the same time as the monopolar ones, can be activated by two modes:**

- For the cut or the coagulation, by the double pedal foot-switch.
- For the coagulation, by selecting the current *Micro Auto* with *Impedance Sensing* automatic *Auto Start / Auto Stop* device with *Start Delay* adjustable from 0 to 5sec. This represents a very useful system, which allows obtaining coagulation by the hand-switch activation without using any special forceps.

## THE PEDAL FOOT-SWITCHES

The Excell units can be used by one double foot-switch pedal for the activation of the monopolar or bipolar functions and by two double foot-switch pedals, one for the monopolar functions and the other for the bipolar functions.

## ELIMINATION AND FILTERING OF THE ELECTROSURGICAL SMOKES

In order to eliminate the bacterial/viral contamination, the Excell units can be used together with electrosurgical smoke evacuators, which are equipped with suction functions expressly studied for the laparoscopic procedures and for solving the problem of a perfect visual field while operating.

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## TECHNICAL FEATURES

**Electronic generator, compliant with:** IEC 601-1 and IEC 60601-2-2 Standards.

**Class and EC Certificate (93/42+2007/47):** IIB, CE0051-nr. 187/MDD.

**Electromagnetic compatibility:** Unit is compliant with IEC 60601-1-2, Test report IMQ nr.80SF00709/1, 80SF00709/2.

**Classification and type IEC 601-1:** Class I, Type CF.

Typical LF leakage currents: on the patient: 4  $\mu$ A= 0,004mA, into the enclosure: 1  $\mu$ A= 0,001mA, to earth: 30  $\mu$ A= 0,3mA.

**Output circuit according to IEC 60601-2-2:**

"Floating" insulated at low and high frequencies, protected against the use of the defibrillator.

Typical high frequency leakage currents: 100 mA.

**Monopolar and bipolar working frequency:** 440kHz.

**Self-Regulation System with double micro-processor Master Sleeve System:**

**Auto-Check**, complete auto-control of the functioning, through:

- **Main Auto-Check:** main auto-diagnosis procedure at the switching-on;
- **Standard Auto-Check:** continuous auto-diagnosis during the functioning;
- **General Error Control and Output Error Control:** the auto-check systems interrupt the functioning and inform the operators through specific Error Codes in case of errors in the activation, damages in the functioning or in case of problems that could provoke a delivery of power, which is higher than the selected one.

**ADC System-constant power e Starting Impulse Control**

All the currents are automatically self-regulated, according to a continuous feed-back in real time (7,000 checks / sec) between the unit and the tissues on which the operator intervenes, and according to the electrode used and its rapidity of sliding on the tissues themselves.

The initial pulse of the cut currents is optimized automatically in order to obtain the maximum efficacy.

**Functioning memorization:** 10 programs.

**Control of high frequency leakage currents:** by a specific circuit with alarm signal and automatic decrease of output powers.

**Outputs, activations, and pedal foot-switches**

There are two monopolar outputs that can be used at the same time by two operators:

- Both of them by the hand-switches, both of them by the double foot-switch pedal or one by the hand-switches and one by the double foot-switch pedal.

Into the model *MCDSe* 400/A, when also the functioning with Argon gas is selected:

- Both of them by the hand-switches or by the double foot-switch pedal, one for the normal electro-surgery and the other for the electro-surgery with Argon gas.

There is one bipolar output, which can be used independently and at the same time as the monopolar ones:

- For the cut or the coagulation, by the double foot-switch pedal.
- For the coagulation, by selecting the current *Micro Auto* with *Impedance Sensing* automatic *Auto Start / Auto Stop* device with *Start Delay* adjustable from 0 to 5sec.

The Excell units can be equipped with:

- One double foot-switch pedal for the activation of the monopolar or bipolar functions.
- Two double foot-switch pedals (one for the monopolar functions and the other for the bipolar functions).

The pedals are compliant with the Standard IEC 60601-2-2, they are waterproof (IPX7) with low voltage 12VDC supply for medical safety.

**Micro/macro adjusting powers by push-buttons with low/fast changing and steps:**

Monopolar: 0/30W=1W, 30/100W= 2W, 100/200W=5W, da 200W=10W.

Bipolar: 0/10W=0,5W, 10/30W=1W, 30/100=2W, da 100 W=5W.

**Powers control:** by displays.

**Safety circuit of the neutral electrode NPCC System**

Control of the connection/contact of the neutral electrodes *Split* - twin section type, as well as of the electrodes *Non Split* - one single section type.

With the *Split* type, it fully checks the connection and the good contact of the electrode.

If the contact is not optimum (approx. 80% of the surface of a standard electrode for adults), it gives the first luminous indication. If the contact decreases to approx. 50%, it gives another luminous indication, and automatically reduces the output power at max. 200W. If the contact is dangerous or the electrode cable is broken / not connected at all, it interrupts the delivery of the power and gives a luminous (red) and acoustic (loud, intermittent, not adjustable) alarm.

With the *Non Split* type, it only controls if the cable is integral and correctly connected.

**The acoustic and luminous signals are compliant with IEC 60601-2-2**

Activation of the Cut (yellow light - loud sound); Activation of the Coagulation (blue light - acute sound); Neutral Electrode (red light - loud, intermittent and not adjustable sound).

**Mains:** 230/115V $\pm$  10% ~ 50/60Hz.

**Absorption at 230V:** Max power 3,6A=828VA, Stand-by 0,4A= 92VA.

**Enclosure:** protected against the penetration of liquids, according to IEC 60601-2.

**Cooling:** by convection, without fan.

**Equipotential Connection:** standard plug DIN 42801.

**Dimensions and Weight (LxDxH):** 400/A *MCDSe*: 38x38x16cm. – 16kgs.

400 *MCDSe*, 350 *MCDSe*, 250 *MCDSe*, 200 *MCDSe*: 38x35x16cm. – 15kgs.

**Argon gas section (only into the model Excell 400/A *MCDSe*)**

**Supply:** Both with 1 or 2 cylinders of 5lt. and with centralized system.

**Max. Flow and Gas Pressure:** 15lt/min. At the entrance: 2.5atm; during the functioning: 1atm.

**Control of the flow through *Constant Flow System*:** from 1 to 15lt/min, through an electronic sensor with adjusting buttons, visual check on the bar-led, automatic auto-compensation according to the used electrode, and alarm just in case of absence of gas.

**Control of the pressure into the circuit *Safety Gas System*:** Reducer of the pressure on the cylinder; internal reducer of the pressure with safety valve; pressure sensor connected to the electronic control system with Auto-Check when the gas section is switched on.

**Protection of the delivered gas flow:** The gas output is equipped with antibacterial filter.

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Insulated bipolar scissors: **BITec** with standard blades - **BITecD** with delicate blades (all insulated except the tips)

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Insulated bipolar forceps: **PMC - PBC** standard, **PMC/L - PBC/L** with irrigation, **PMC/ns - PBC/ns** with no-stick tips.

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# THE UNITS AND THE ACCESSORIES

## UNITS AND STANDARD SETS OF ACCESSORIES

B900	<b>EXCELL 400 MCDSe</b> , without accessories
B910	<b>EXCELL 350 MCDSe</b> , without accessories
B920	<b>EXCELL 250 MCDSe</b> , without accessories
B930	<b>EXCELL 200 MCDSe</b> , without accessories
B940	<b>EXCELL 400/A MCDSe</b> , without accessories
B610/A	<b>SET OF STANDARD ACCESSORIES</b> as follows:
1 DS/E	Double waterproof foot-switch pedal
1 NP/A	Stainless steel neutral electrode for adult - cable 2.5 mt. long
2 MPE/E	Sterilizable electrodes-holder handle, cable 3.5 mt. long
1 SEL/E	<b>Composed by:</b> <b>2 E1</b> - Straight knife (blade) electrode, <b>2 E5</b> - Thick needle electrode, <b>1 E7</b> - Thin needle electrode, <b>1 E12</b> - Straight ball electrode Ø 2.5 mm, <b>2 E14</b> - Straight ball electrode Ø 4 mm
B610/B	<b>SET OF STANDARD ACCESSORIES</b> same as B610/A but with flexible neutral electrode made of conductive rubber for adult NP/GA
B610/P	Ditto, but with paediatric neutral electrode NP/GP

## GENERAL LIST OF ACCESSORIES

### HANDLES FOR USE BY FOOT-SWITCH PEDALS

MPE/E	Handle, cable 3.5 mt. long, sterilizable
MPE/E5	Ditto, cable 5 mt. long

### HAND-SWITCH HANDLES

MPE/CMS	Double push-button handle with knife electrode, cable 3 mt. (100 times sterilizable)
MPE/CMS5	Ditto, cable 5 mt.
MP/CM	Ditto, disposable type, cable 3 mt. (*)

### ACTIVE ELECTRODES, made of stainless steel, insulated stem Ø 2.3÷2.4 mm. sterilizable

#### SHORT TYPE, STEM 70 mm. LONG

E1	Straight knife (blade) electrode
E1/I	Totally insulated except final 5 mm. straight knife (blade) electrode
E3	Angled knife (blade) electrode
E1/L	Straight lancet electrode
E3/L	Ditto, angled
E5	Straight thick needle electrode
E6	Ditto, angled
E7	Straight thin needle electrode
E7/I	Totally insulated except final 5 mm. straight thin needle electrode
E8	Angled thin needle electrode
E10	Very fine needle electrode, Ø 0.40 mm.
E12	Straight ball electrode Ø 2.5 mm.
E13	Ditto, angled
E14	Straight ball electrode Ø 4 mm.
E15	Ditto, angled
E16	Straight ball electrode Ø 6 mm.
E17	Ditto, angled
E18	Wire diamond loop electrode, 5x10 mm.
E19	Ditto, 10x10 mm.
E21	Round wire loop electrode Ø 5 mm.
E23	Ditto, Ø 10 mm.
E25	Ditto, Ø 15 mm.
E23/N	Round ribbon loop electrode Ø 10 mm.
E25/N	Ditto, Ø 15 mm.
E26	Plate electrode
EXT/15	Extension, long 15cm. for all the electrodes with stems Ø2.3÷2.4 mm.

#### LONG TYPE, FLEXIBLE STEM 130 mm. LONG

E40	Straight knife (blade) electrode
E40/I	Totally insulated except final 5 mm. straight knife (blade) electrode
E41	Straight thick needle electrode
E42	Straight thin needle electrode
E42/I	Totally insulated except final 5 mm. straight thin needle electrode
E43	Round wire loop electrode Ø 5 mm.
E44	Ditto, Ø 10 mm.
E45	Ditto, Ø 15 mm.
E46	Straight ball electrode Ø 2.5 mm.
E47	Straight ball electrode Ø 4 mm.
E47/6	Straight ball electrode Ø 6 mm.

### ELECTRODES FOR MICROSURGERY, STERILIZABLE

MID	Reducer for needles (for all the electrode holders)
SAD	Series of 12 needles Ø 0.10 mm.
SAD/1	Ditto, Ø 0.15 mm.
SAD/2	Ditto, Ø 0.20 mm.
SAD/3	Ditto, Ø 0.40 mm.

### RE-USABLE NEUTRAL ELECTRODES

NP/A	Stainless steel electrode for adult (cm. 25x13), cable 2.5 mt. long
NP/GA	Electrode for adults (25x12cm.), flexible conductive rubber, cable 2.5 mt.
NP/GA5	Ditto, cable 5 mt. long
NP/GP	Ditto, pediatric type (18x8 cm.), cable 2.5 mt. long
NP/GP5	Ditto, cable 5 mt. long
FGE	Fixing belt for electrodes with 2 buttons (*)

### DISPOSABLE, ADHESIVE, NEUTRAL ELECTRODES

CMS/E	Re-usable connection cable 2.5 mt. long
CMS/E5	Ditto, cable 5 mt. long
EIP/DA	Adhesive, no-Split single section, adult type (25pcs, each pack) (*)
EIP/SA	Ditto, Split double section, adult type (25pcs, each pack) (*)
EIP/DP	Adhesive, no-Split single section, paediatric type (25pcs. each pack) (*)
EIP/SP	Ditto, Split double section, paediatric type (25pcs. each pack) (*)

### INSULATED MONOPOLAR FORCEPS FOR COAGULATION, WITHOUT CONNECTING CABLES TO THE UNIT, STERILIZABLE

PIC/1	Straight forceps (Cushing/Potts-Smith) ("grasping" tips 1 mm. - L. 18 cm.)
PIC/1-25	Ditto, L. 25 cm.
PIC/2	Straight forceps (Cushing/Potts-Smith) ("grasping" tips 2 mm. - L. 25 cm.)

### INSULATED MONOPOLAR FORCEPS FOR COAGULATION, WITH CONNECTING CABLES TO THE UNIT, STERILIZABLE

CPI	Connecting cable for PMI, 3.5 mt.
CPI/5	Ditto, 5 mt. long
PMI/1	Straight forceps (Cushing/Potts-Smith) ("grasping" tips 1 mm. - L. 18 cm.)

PMI/1-20	Ditto, 20 cm. long
PMI/1-25	Ditto, 25 cm. long
PMI/2	Straight forceps (Cushing/Potts-Smith) ("grasping" tips 2 mm. - L. 25 cm.)
PMI/B	Bayonet forceps (Jansen/Yasargil) ("grasping" tips 2 mm. - L. 20 cm.)

### HAND-SWITCH INSULATED MONOPOLAR FORCEPS FOR COAGULATION, WITH CONNECTING CABLE TO THE UNIT, STERILIZABLE

PMI/BJ21	Bayonet forceps, cable 3.5 mt. long ("grasping" tips 2 mm. - L. 21 cm.)
PMI/PJ21	Straight forceps, cable 3,5 mt long ("grasping" tips 2 mm. - L. 21 cm.)
PMI/PJ24	Ditto ("grasping" tips 2 mm. - L. 24 cm.)

### Monopolar accessories for laparoscopy, please ask for specific detail.

### CONNECTING CABLES FOR INSTRUMENTS FOR LAPAROSCOPY

CPE	Connecting cable for instruments with male / female connector Ø 4 mm, 3,5 mt.
CPE/5	Ditto, L. 5 mt.

### Bipolar scissors for surgery with their cables, please ask for specific detail.

### CONNECTING CABLES FOR BIPOLAR FORCEPS/ELECTRODES AND FOR BIPOLAR HOOKS, FORCEPS, SCISSORS FOR LAPAROSCOPY, STERILIZABLE

CPB/E	Connecting cables, 3 mt.
CPB/E5	Ditto, L. 5 mt.

### INSULATED BIPOLAR FORCEPS AND RIGID ELECTRODES, STERILIZABLE

#### Standard forceps for bipolar coagulation

PMC/JR	Straight forceps (Jeweler) (straight tips 0.5 mm. - L. 11.5/12 cm.)
PMC/JC	Ditto, angled tips
PMC/RS	Straight forceps (Cushing/Potts-Smith) (straight tips 0.7 mm. - L. 15.5/16 cm.)
PMC/CS	Ditto, angled tips
PMC/R	Straight forceps (Cushing/Potts-Smith) (straight tips 1 mm. - L. 20 cm.)
PMC/C	Ditto, angled tips
PBC/R	Straight forceps (Cushing/Potts-Smith) (straight tips 2 mm. - L. 20 cm.)
PBC/C	Ditto, angled tips
PMC/R25	Straight forceps (Cushing/Potts-Smith) (straight tips 1mm. - L. 25 cm.)
PMC/C25	Ditto, angled tips
PBC/R25	Straight forceps (Cushing/Potts-Smith) (straight tips 2 mm. - L. 25 cm.)
PBC/C25	Ditto, angled tips
PMC/RSB	Bayonet forceps (Jensen/Yasargil) (straight tips 0.7 mm. - L. 16.5/17 cm.)
PMC/B	Bayonet forceps (Jensen/Yasargil) (straight tips 1 mm. - L. 20cm.)
PMC/BCD	Ditto, angled tips down
PMC/BCU	Ditto, angled tips up
PBC/B	Bayonet forceps (Jensen/Yasargil) (straight tips 1 mm. - L. 20cm.)
PBC/BCD	Ditto, angled tips down
PBC/BCU	Ditto, angled tips up
PMC/B25	Bayonet forceps (Jensen/Yasargil) (straight tips 1 mm. - L. 25 cm.)
PBC/B25	Bayonet forceps (Jensen/Yasargil) (straight tips 1 mm. - L. 25 cm.)

### Forceps with irrigation for bipolar coagulation, please ask for specific detail.

### Forceps with no-stick tips for bipolar coagulation, please ask for specific detail.

### Rigid electrodes for bipolar coagulation of turbinal or larynx, please ask for specific detail.

### Bipolar hooks, forceps and scissors for laparoscopy, please ask for specific detail.

### Accessories and cables for flexible endoscopy, please ask for specific detail.

### ADAPTORS FOR USING CABLES WITH NON ALSA STANDARD PLUG

RD5	For monopolar cables with plug Ø from 2 to 8 mm or Martin standard.
RD/BF	For bipolar cables with double plug Ø 4 mm (international standard) or with plug Valleylab /Commed standard.
RD/BF1	For bipolar cables with coaxial plug Ø 12,5 mm (Erbe/Storz standard)
RD/BF2	For bipolar cables with coaxial plug Ø 8 mm (Martin/Bertchold standard)

### TROLLEYS, FOOT-SWITCH PEDALS, BOXES FOR ACCESSORIES

H10/AB	Economic trolley with 2 shelves 45x40x82 cm. antistatic wheels, 2 with brakes
H23/SE	Trolley with 3 shelves. Dimensions: 50x50x80 cm., antistatic wheels, 2 with brakes
H26	Trolley for electrosurgical unit and smoke evacuator, with seats for Argon gas tanks. Dimensions: 52x58x80 cm., Antistatic wheels, 2 with brakes
DS/E	Double foot-switch pedal, water-proof (IPX7)
DS/B	Ditto, for bipolar activation only
BOX-TE	Stainless steel round box for electrodes
BOX-RA	Stainless steel rectangular box for accessories, 21x10x5 cm.

## ACCESSORIES FOR ARGON GAS SURGERY

### For model 400/A MCDSe only

### HAND-SWITCH HANDLE AND ELECTRODES FOR SURGERY / LAPAROSCOPY

AC/HANDLE	Double push-button electrode-holder, sterilizable, cable 3.5 mt.
AC/E25-C	Rigid electrode for coagulation, 25 mm. long, sterilizable
AC/E100-C	Ditto, 100 mm. long
AC/E320-C	Ditto, 320 mm. long (for laparoscopy only)
AC/E 320-H	L Hook rigid electrode, 320 mm. long, sterilizable (for laparoscopy only)
AC/E40-A	Rigid needle electrode, 40 mm. long, sterilizable
AC/E100-A	Idem, 100 mm. long
AC/E40-L	Rigid knife electrode, 40 mm. long, sterilizable
AC/E100-L	Ditto, 100 mm. long

### CABLE AND FLEXIBLE ELECTRODES / PROBES FOR FLEXIBLE ENDOSCOPY

AC/CABLE+	Connecting cable for flexible electrodes/probes, sterilizable, 3,5mt. long
-----------	--

### ELECTRODES / FLEXIBLE PROBES: please ask for specific detail.

### TROLLEY, GAS CYLINDER, PRESSURE REDUCER, BACTERIAL FILTER

H26	See above
B5	Argon gas cylinder, capacity 5lt (RD/P - Pressure reducer for cylinder)
ESU/TG	Tube for gas input with quick-action coupling (for B5)
ESU/F	Bacterial filter for gas output (ESU/FC - Metallic connector for filter)

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**POLIVAC B4/SLT SERIES**  
for high surgery / vacuum-extraction

**GENERALS FEATURES**

Built in accordance with the latest International Standards IEC 60601-1 and EN 10079-1 and the EC Certificates, the series of suction units POLIVAC B4/SLT, POLIVAC B4/SLT-EXP allow the final user to meet all his needs.

All the models, characterized by a big power suction pumps and high suction speed, are very quite. They have many common characteristics, but different suction group and standard accessories.

**POLIVAC B4/SLT 30 2**

It is the basic model of the series. It is suitable for the operating theatre and other medical suction interventions. Its technical characteristics are as follows:

- Class IEC/CEI 60601-1: I BF.
- EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
- Usage for unlimited time / Noise: less than 50 db.
- Power supply / Absorption: 230V - one phase, 50Hz / 250VA.
- Rotary oil bath pump-unit (max vacuum 730 mmHg "0,96 atm, 0,96 bar", flow rate 50lt/min, engine: from 1/3 of HP).
- Painted metal case and antistatic wheels, two with castors.
- Two bottles of 4 lt. made of unbreakable and autoclavable polycarbonate, with internal over-flow devices, and a quick system for the aspiration changing between bottles; silicone autoclavable tubing (Ø mm 7x12).
- Vacuum regulator and vacuum check device (0-76 cm/Hg).
- Main switch and pedal switch (pneumatic, very safe, water-proof, anti-explosion).
- Dimensions (width/depth/height), weight: cm 40x35x75 Kg 24.

**It can be equipped with antibacterial filter!!:**  
See note 1)  
**It can be equipped with many components!!:**  
See note 2)

**POLIVAC B4/SLT 30 1**

It is the same as model POLIVAC B4/SLT 30 2, but with 1 bottle 4 lt. only.

**POLIVAC B4/SLT 30 2-EXP**

It is the basic model of the series. It is suitable for the operating theatre and other medical suction interventions. Its technical characteristics are as follows:

- Class IEC/CEI 60601-1: I BF.
- EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
- Usage for unlimited time / Noise: less than 50 db.
- Power supply / Absorption: 230V - one phase, 50Hz / 400VA.
- Piston suction-group, maintenance free (max vacuum 700 mmHg "0,9 atm, 0,90 bar", flow rate 70 lt./min, engine: from 1/2 of HP).
- Painted metal case and antistatic wheels, two with castors.
- Two bottles of 4 lt. made of unbreakable and autoclavable polycarbonate, with internal over-flow devices, and a quick system for the aspiration changing between bottles; silicone autoclavable tubing (Ø mm 7x12).
- Vacuum regulator and vacuum check device (0-76 cm/Hg).
- Main switch and pedal switch (pneumatic, very safe, water-proof, anti-explosion).
- Dimensions (width/depth/height), weight: cm 40x35x75 Kg 24.

**It can be equipped with ant bacterial filter!!:**  
See note 1)  
**It can be equipped with many components!!:**  
See note 2)

**POLIVAC B4/SLT 30 1-EXP**

It is the same as model POLIVAC B4/SLT 30 2EXP, but with 1 bottle of 4 lt. only.

**POLIVAC B4/SLT 30 P-EXP**

It is the portable model of the series. It is suitable for the operating theatre and other medical suction interventions. Its technical characteristics are as follows:

- Class IEC/CEI 60601-1: I BF.
- EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
- Usage for unlimited time / Noise: less than 50 db.
- Power supply / Absorption: 230V - one phase, 50Hz / 400VA.
- Piston main suction-group, maintenance free (max vacuum 700 mmHg "0,9 atm, 0,90 bar", flow rate 70 lt./min, engine: from 1/2 of HP).
- Painted metal case, handle and seat for bottles, which is already set for the usage of disposable bottles of diam. from 9 to 13 cm.
- One bottle (Two bottles, on request) of 2 lt. made of unbreakable and autoclavable polycarbonate, with internal over-flow device, silicone autoclavable tubing (Ø mm 7x12).
- Vacuum regulator and vacuum check device (0-76 cm/Hg).
- Main switch and, on request, pedal switch STOP/PN (pneumatic, very safe, water-proof, anti-explosion).
- Dimensions (width/depth/height), weight: cm 38x35x28 Kg 17.

**It can be equipped with antibacterial filter!!:**  
See note 1)  
**It can be equipped with many components!!:**  
See note 3)

**POLIVAC B4/SLT 30 1VE-EXP**

It is the most suitable model both for the operating theatre and other medical suction interventions as well as the Vacuum-Extraction procedures. It is equipped also with the vacuum regulator pedal switch.

- Its technical characteristics are as follows:
- Class IEC/CEI 60601-1: I BF.
  - EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
  - Usage for unlimited time / Noise: less than 50 db.
  - Power supply / Absorption: 230V - one phase, 50Hz / 400VA
  - Piston suction-group, maintenance free (max vacuum 700 mmHg "0,9 atm, 0,90 bar", flow rate 70lt/min, engine: from 1/2 of HP).
  - Painted metal case and antistatic wheels, two with castors.
  - Two bottles of 4 lt. made of unbreakable and autoclavable polycarbonate, with internal over-flow devices, and a quick system for the aspiration changing between bottles; silicone autoclavable tubing (Ø mm 7x12).
  - Vacuum regulators on the panel and by the pedal switch (it stops by itself and it maintains the position set by the operator); vacuum check device (0-76 cm/Hg).
  - Main switch and pedal switch (pneumatic, very safe, water-proof, anti-explosion).
  - Dimensions (width/depth/height), weight: cm 40x35x75 Kg 24.

**It can be equipped with antibacterial filter!!:**  
See note 1)  
**It can be equipped with many components!!:**  
See note 2)

**POLIVAC B4/SLT 30 P VE-EXP**

It is the portable model, suitable both for the operating theatre and other medical suction interventions as well as the Vacuum-Extraction procedures. It is equipped also with the vacuum regulator pedal switch. Its technical characteristics are as follows:

- Class IEC/CEI 60601-1: I BF.

- EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
- Usage for unlimited time / Noise: less than 50 db.
- Power supply / Absorption: 230V - one phase, 50Hz / 400VA.
- Piston suction-group, maintenance free (max vacuum 700 mmHg "0,9 atm, 0,90 bar", flow rate 70 lt./min, engine: from 1/2 of HP).
- Painted metal case, handle and seat for bottles, which is already set for the usage of disposable bottles of diam. from 9 to 13 cm.
- One bottle (Two bottles, on request) of 2 lt. made of unbreakable and autoclavable polycarbonate, with internal over-flow device, silicone autoclavable tubing (Ø mm 7x12).
- Vacuum regulators on the panel and by the pedal switch (it stops by itself and it maintains the position set by the operator); vacuum check device (0-76 cm/Hg).
- Main switch and, on request, pedal switch STOP/PN (pneumatic, very safe, water-proof, anti-explosion).
- Dimensions (width/depth/height), weight: cm 38x35x28 Kg 17.

**It can be equipped with antibacterial filter!!:**  
See note 1)  
**It can be equipped with many components!!:**  
See note 3)

**POLIVAC B4/SLT 50 E**

It is the most complete model of the series for the operating theatre; it is equipped with two independent suction-groups: the main one and the anaesthesia one. Its technical characteristics are as follows:

- Class IEC/CEI 60601-1: I BF;
- EC Certificate 93/42: EC 0051
- Usage for unlimited time / Noise: less than 50 db;
- Power supply / Absorption: 230V - one phase, 50Hz / 450VA
- Rotary oil bath main suction group (max vacuum 730 mmHg "0,96 atm, 0,96 bar", flow rate 50lt/min, engine: from 1/3 of HP);
- Piston anaesthesia suction-group, maintenance free (max vacuum 700 mmHg "0,9atm, 0, 96 bar", flow rate 30 lt./min, engine: from 1/5 of HP);
- Painted metal case, with antistatic wheels, two with castors;
- Two bottles of 4 lt. for the main suction-group, with quick system for the aspiration changing between bottles; one bottle of 1 lt. for the anaesthesia group. All the bottles are made of unbreakable and autoclavable polycarbonate, with internal over-flow devices. Silicone autoclavable tubing (Ø mm 7x11);
- Two vacuum regulators (one for each suction-group) and vacuum check devices (0-76 cm/Hg) ;
- Main switch and pedal switch (pneumatic, very safe, water-proof, anti-explosion);
- Dimensions (width/depth/high), weight: cm 38x35x28 - Kg 30.

**It can be equipped with antibacterial filter !!:**  
See note 1)

- Note 1) It can be equipped with a disposable antibacterial filter, which works also as supplementary over-flow device;
- Note 2) Available pieces: external over-flow device, seat for the usage of disposable bottles of diam. from 9 to 16 cm (also for wall), cannulae-handle, stainless steel box for the accessories;
- Note 3) Available pieces: external over-flow device, cannulae-handle, and stainless steel box for the accessories.

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## ACCESSORIES

GP/TPE- External over-flow device.  
 C/H- Autoclavable transparent cannulae-handle.  
 STOP/PN- Pedal switch (for the activation).  
 STOP/VR- Pedal switch (for the vacuum regulation).  
 H10/AB- Mobile trolley with two shelves (cm 40x40), and antistatic wheels, two of which with castors.  
 GP/FM- Disposable anti-bacterial filter group. (\*)  
 R/FM- Spare filter for the group GP/FM. (\*)  
 VMLT4- Autoclavable, unbreakable jar of 4 lt. made of polycarbonate.  
 VMLT2- Autoclavable, unbreakable jar of 2 lt. made of polycarbonate.  
 TP/PG- Autoclavable cap with over-flow device for VMLT4 and VMLT2, made of nylon.  
 PT/7- Connector for tubing TS/7.  
 TS/7- Autoclavable standard tubing made of silicone - Ø mm 7x12.  
 TS/A- Autoclavable standard tubing made of silicone - Ø mm 7x12 - of 2mt. long, and with connector PT/7.

*JARS FOR DISPOSABLE BAGS OF LT. 2, ask for specific details.*

### CANNULAE FOR SURGERY

1-C/SAC- Series of 4 cannulae with handle: Yankhauer curved, Yankhauer soft, Universal curved and Universal soft (terminal part that does not attach to the tissue).  
 2-C/SAC2- Series of 2 cannulae with handle and finger switch: Yankhauer curved, Yankhauer soft.  
 3-C/SC2- Series of 2 cannulae: straight, curved.

### CANNULAE FOR GENERAL USAGE AND ENT

#### OTO-RHINO-LARYNGOLOGIE

6-C/AS- Curved, saliva sucking.  
 6-C/AM- Angled, mucus sucking.  
 7-C/MS6- Soft, made of silicone, Ø mm 6 with finger switch.  
 7-C/MS8- See above, Ø mm 8  
 8-C/UM6- Universal, curved, with manual finger switch.  
 8-C/UM8- See above, straight, soft.  
 9-C/F1- "Frazier / Ferguson", Ø mm 1,5 with finger switch.  
 9-C/F2- See above, Ø mm 2  
 9-C/F3- See above, Ø mm 3  
 9-C/F4- See above, Ø mm 4

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10-C/J1- "Jackson", cm 35 x Ø mm 3  
 10-C/J2- See above, cm 45  
 10-C/J3- See above, cm 25 x Ø mm 4  
 10-C/J4- See above, cm 35  
 10-C/J5- See above, cm 45  
 10-C/J6- See above, cm 55  
 10-C/J6- See above, cm 60 x Ø mm 4

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## VORTEX-S AS-200 SERIES for medium surgery, Vacuum-extraction

### GENERAL FEATURES

Built in accordance with the latest International Standards IEC 60601-1 and EN 10079-1 and the EC Certificates, the series of suction units VORTEX AS200 and VORTECO AS100 are composed of different models, which allow to meet every kind of need in the field of small and medium surgery (i.e. the small operating theatres and the Vacuum-extraction procedures).

### VORTEX-S AS-200/20LT

It is the basic model of the series. It is suitable for: Emergency Room, hospital, endoscopy, ENT, tracheotomized people, gynaecology - endouterine treatments, small and medium surgery, small operating theatres, etc. Its technical characteristics are as follows:

- EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
- Class IEC 60601-1: I BF.
- Usage for unlimited time (8 hours) / Noise: less than 50 dB.
- Power supply / Absorption: 230V - one phase, 50Hz/138VA.
- Piston pump-unit, maintenance free (max. vacuum: 700 mm Hg "0.9 atm, 0.90 bar", flow rate: 20lt/min, engine: 1/5 HP).
- Painted metal case, handle and seat for the jars, already prepared for the usage of disposable jars (Ø from 9 to 12 cm).
- One (two, on request) collection jar of 2 lt. made of unbreakable and autoclavable polycarbonate, with internal overflow device, and silicone autoclavable tubing (Ø mm 7x12).
- Vacuum regulator and vacuum check device (0-76 cm/Hg).
- Main switch and, on request, pedal switch STOP/PN (pneumatic, very safe, water-proof, anti-explosion).

- Dimensions (width/depth/height), weight: cm 33.5x30x25 - Kg 7.5.
- Dimensions (width/depth/height), weight with trolley H5: cm 50x50x85 - Kg 17.5.
- It can be equipped with antibacterial filter!!: See note 1)
- It can be equipped with many components!!: See note 2)
- It can be equipped with trolley!!: See note 3)

### VORTEX-S AS-200/30LT

It is the same as model VORTEX AS200/20LT, but with a more powerful pump-unit, which assures a higher aspiration speed.

The extra technical characteristics are as follows:

- Power supply / Absorption: 230V - one phase, 50Hz / 140VA.
- Piston pump-unit, maintenance free (max. vacuum 700 mm Hg "0.9 atm, 0.90 bar", flow rate: 30lt/min, engine: 1/5 HP).

### VORTEX-S AS-200/20LT VE

It is the model of the series suitable both for Emergency Room, hospital, endoscopy, gynaecology - endouterine treatments -, small and medium surgery, small operating theatres, as well as the Vacuum-Extraction procedures, because it is equipped also with the pedal switch for the vacuum regulation. Its technical characteristics are as follows:

- EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
- Class IEC 60601-1: I BF.
- Usage for unlimited time (8 hours) / Noise: less than 50 dB.
- Power supply / Absorption: 230V - one phase, 50Hz / 138VA.
- Piston pump-unit, maintenance free (max. vacuum: 700 mm Hg "0.9 atm, 0.90 bar", flow rate: 20lt/min, engine: 1/5 HP).
- Painted metal case, handle and seat for the jars, already prepared for the usage of disposable jars (Ø from 9 to 12 cm).
- One (two, on request) collection jar of 2 lt. made of

- unbreakable and autoclavable polycarbonate, with internal overflow device, and silicone autoclavable tubing (Ø mm 7x12).
- Panel and pedal vacuum regulator (self-blocking, it maintains the initial position set by the operator), and vacuum check device (0-76 cm/Hg)
- Main switch and, on request, pedal switch STOP/PN (pneumatic, very safe, water-proof, anti-explosion)
- Dimensions (width/depth/height), weight: cm 33.5x30x25 - Kg 7.5.
- Dimensions (width/depth/height), weight with trolley H5: cm 50x50x85 - Kg 17.5.
- It can be equipped with antibacterial filter!!: See note 1)
- It can be equipped with many components!!: See note 2)
- It can be equipped with trolley!!: See note 3)

### VORTEX-S AS-200/30LT VE

It is the same as model VORTEX AS200/20LT-VE, but with a more powerful pump-unit, which assures a higher aspiration speed.

The extra technical characteristics are as follows:

- Power supply / Absorption: 230V - one phase, 50Hz / 140VA.
- Piston pump-unit, maintenance free (max. vacuum: 700 mm Hg "0.9 atm, 0.90 bar", flow rate: 30lt/min, engine: 1/5 HP).

- Note 1) It can be equipped with a disposable antibacterial filter, which works also as supplementary overflow device;
- Note 2) Available pieces: external overflow device, cannulae-handle, stainless steel box for the accessories;
- Note 3) Available also with trolley model H5. In this case, the unit does not have any seat for the usage of collection jars (they are on the trolley, and they are already prepared for the usage of disposable jars of diam. from 9 to 16 cm (also for wall). The unit is fixed on the trolley.

## VORTECO AS-100 SERIES for small and medium surgery, Vacuum-extraction, emergency

### VORTECO AS-100/EMERGENCY(\*)

It is the model of the series that can be supplied both with AC mains and with chargeable emergency battery. It is suitable for: private usage, Emergency Room, hospital, the emergency trolleys, etc. Its technical characteristics are as follows:

- EC Certificate (93/42+2007/47): CE0434.
- Class IEC of the unit: I B.
- Usage: 45 min. (even with battery) / Noise: less than 65 dB.
- Power supply / Absorption: 12VDC / 690VA.
- Class IEC / Power supply / Absorption: Class II / 230 V / 990 VA.
- Battery: hermetic and chargeable 12VDC 3 Ah.
- Piston pump-unit, maintenance free (max. vacuum: 600 mmHg "0.8 atm, 0.80 bar", flow rate: 20lt/min).
- Painted metal case, handle and seat for the jars.
- One collection jar of 1 lt. made of unbreakable and autoclavable polycarbonate, with internal overflow device, and silicone autoclavable tubing (Ø mm 7x12).
- Main switch, vacuum regulator and vacuum check device (0-76 cm/Hg).
- Dimensions (width/depth/height), weight: cm 37x26x10.5 - Kg 6.
- The basic version with 1 jar of 1 lt. can be equipped also with a bag for the transport, and a cable for the connection to the car lighter socket.

It can be equipped with antibacterial filter!!: See note 1)

### VORTECO AS-100/15LT

It is the smallest model of the series. It is suitable for: Emergency Room, hospital and small surgery. Its technical characteristics are as follows:

- Class IEC/CEI 60601-1: I BF.
- EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
- Usage: 01 hour / Noise: less than 65 dB.
- Power supply/Absorption: 230V one-phase, 50 Hz/180VA.
- Piston pump-unit, maintenance free (max. vacuum: 650 mmHg "0.86atm, 0.86 bar", flow rate: 15lt/min, engine: 1/4 H/P).
- Painted metal case, handle and seat for the jars, already prepared for the usage of the disposable jars (Ø from 9 to 12 cm).
- One (two, on request) collection jar of 1 or 2 lt. made

- of unbreakable and autoclavable polycarbonate, with internal overflow device, and silicone autoclavable tubing (Ø mm 7x12).
- Main switch, vacuum regulator and vacuum check device (0-76 cm/Hg).
- Dimensions (width/depth/height), weight: cm 33.5x30x25 - Kg 6.8.

It can be equipped with antibacterial filter!!: See note 1)  
It can be equipped with many components!!: See note 2)

### VORTECO AS-100/20LT

It is the model of the series suitable both for Emergency Room, hospital, endoscopy, ENT, tracheotomy, gynaecology - endouterine treatments -, small and medium surgery, small operating theatres, etc. Its technical characteristics are as follows:

- Class IEC 60601-1: I BF.
- EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
- Usage for unlimited time (6 hours) / Noise: less than 50 dB.
- Power supply/Absorption: 230V - one phase, 50Hz/138VA.
- Piston pump-unit, maintenance free (max. vacuum 700 mmHg "0.9atm, 0.90 bar", flow rate: 20lt/min, engine: 1/5 HP).
- Painted metal case, handle and seat for the jars, already prepared for the usage of disposable jars (Ø from 9 to 12 cm).
- One (two, on request) collection jar of 1 or 2 lt. (2 lt. for the endouterine function) made of unbreakable and autoclavable polycarbonate, with internal overflow device, and silicone autoclavable tubing (Ø mm 7x12).
- Main switch, vacuum regulator and vacuum check device (0-76 cm/Hg).
- Dimensions (width/depth/high), weight: cm 33.5x30x25 - Kg 7.5.
- Dimensions (width/depth/height), weight with trolley H5: cm 50x50x85 - Kg 17.5.
- It can be equipped with antibacterial filter!!: See note 1)
- It can be equipped with many components!!: See note 2)
- It can be equipped with trolley!!: See note 3)

### VORTECO AS-100/30LT

It is the same as model VORTECO AS100/20LT, but with a more powerful pump-unit, which assures a higher aspiration speed. The extra technical characteristics are as follows:

- Power supply/Absorption: 230V - one phase, 50Hz/140VA.
- Piston pump-unit, maintenance free (max. vacuum: 650/700 mmHg "0.9atm, 0.9 bar, 900mbar", flow rate: 30lt/min, engine: 1/5 HP).

### VORTECO AS-100/30LT VE

It is the model of the series suitable both for Emergency Room, hospital, endoscopy, gynaecology - endouterine treatments -, small and medium surgery, small operating theatres, as well as the Vacuum-Extraction procedures, because it is equipped also with the pedal switch for the vacuum regulation. Its technical characteristics are as follows:

- Class IEC 60601-1: I BF.
- EC Certificate (93/42+2007/47): CE0051 - nr. 187/MDD.
- Usage for unlimited time (6 hours) / Noise: less than 50 dB.
- Power supply/Absorption: 230V - one phase, 50Hz/140VA.
- Piston pump-unit, maintenance free (max. vacuum: 700 mmHg "0.9atm, 0.90 bar", flow rate: 30lt/min, engine: 1/5 HP).
- Painted metal case, handle and seat for the jars, already prepared for the usage of disposable jars (Ø from 9 to 12 cm).
- One (two, on request) collection jar of 1 or 2 lt. (2 lt. for the endouterine function) made of unbreakable and autoclavable polycarbonate, with internal overflow device, and silicone autoclavable tubing (Ø mm 7x12).
- Main switch, vacuum regulator and vacuum check device (0-76 cm/Hg).
- Panel and pedal vacuum regulator (self-blocking, it maintains the initial position set by the operator), and vacuum check device (0-76 cm/Hg).
- Panel switching device.
- Dimensions (width/depth/height), weight: cm 33.5x30x25 - Kg 7.5.
- Dimensions (width/depth/height), weight with trolley H5: cm 50x50x85 - Kg 17.5.
- It can be equipped with antibacterial filter!!: See note 1)
- It can be equipped with many components!!: See note 2)
- It can be equipped with trolley!!: See note 3)

- Note 1) It can be equipped with a disposable antibacterial filter, which works also as supplementary overflow device;
- Note 2) Available pieces: external overflow device, cannulae-handle, stainless steel box for the accessories;
- Note 3) Available also with trolley model H5. In this case, the unit does not have any seat for the usage of collection jars (they are on the trolley, and they are already prepared for the usage of disposable jars of diam. from 9 to 16 cm also for wall). The unit is fixed on the trolley.

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## ACCESSORIES

GP/TPE- External overflow device.  
 C/H- Transparent and autoclavable cannulae-handle.  
 STOP/PN- Pedal switch (for the activation).  
 STOP/VR- Pedal switch (for the vacuum regulation).  
 H10/AB- Trolley with two shelves (cm 40x40) and antistatic wheels, two with castors.  
 H5- Trolley. (\*)  
 GP/FM- Disposable antibacterial filter group. (\*)  
 R/FM- Spare filter for the group GP/FM. (\*)  
 VMLT2- Jar 2lt. made of autoclavable and unbreakable polycarbonate.  
 TP/PG- Cap made of nylon with overflow device for VMLT2, autoclavable.  
 PT/7- Connector for tubing TS/7.  
 TS/7- Standard autoclavable silicone tubing - Ø mm 7x12.  
 TS/A- See above, 2mt long with connector PT/7.  
 VMLT1- Graduated jar of 1lt. made of autoclavable and unbreakable polycarbonate.  
 TP/B1- Cap made of rubber with overflow device for VMLT1, sterilizable.  
 JARS FOR DISPOSABLE BAGS OF LT. 1 and 2, ask for specific details.

### ONLY FOR VORTECO AS 100/EMERGENCY

V/VE- Graduated jar of 1lt made of sterilizable unbreakable polycarbonate. (\*)  
 T/VE- Cap with overflow device for V/AE, sterilizable. (\*)  
 B/VE- Handle for transport. (\*)  
 C/VE- Connecting cable for the lighter socket. (\*)

### CANNULAE FOR SURGERY

1-C/SAC- Series of 4 cannulae with handle: Yankhauer curved, Yankhauer soft, Universal curved and Universal soft (the final part does not attach).  
 2-C/SAC2- Series of 2 cannulae with handle and finger switch: Yankhauer curved, Yankhauer soft.  
 3-C/SC2- Series of 2 cannulae: straight, curved.

### CANNULAE FOR GENERAL USE AND ENT OTO-RHINO-LARYNGOLOGIE

6-C/AS- Curved saliva sucking cannula.  
 6-C/AM- Angled mucus sucking cannula.  
 7-C/MS6- Silicone soft cannula, Ø mm 6 with finger switch.  
 7-C/MS8- See above, Ø mm 8.  
 8-C/UM6- Universal curved cannula with handle and finger switch.  
 8-C/UM8- See above, soft straight cannula.  
 9-C/F1- "Frazier/Fergusson", Ø mm 1.5 with finger switch.  
 9-C/F2- See above, Ø mm 2  
 9-C/F3- See above, Ø mm 3  
 9-C/F4- See above, Ø mm 4

### CANNULAE FOR ENDOSCOPY (with finger switch)

10-C/J1- "Jackson", cm 35xØ mm 3  
 10-C/J2- See above, cm 45  
 10-C/J3- See above, cm 25xØ mm 4  
 10-C/J4- See above, cm 35  
 10-C/J5- See above, cm 45  
 10-C/J6- See above, cm 55  
 10-C/J6- See above, cm 60x Ø mm 4

### DISPOSABLE CATHETERS


11-R/CH- Connection for catheters. (\*)  
 11-CH/14- cm 60 x Ø mm 4.7 (\*)  
 11-CH/16- cm 60 x Ø mm 5.4 (\*)  
 11-CH/18- cm 60 x Ø mm 6 (\*)

### CANNULAE FOR GYNAECOLOGY

12-SET/E- Set for endouterine aspiration, made of: 5 rotating a-traumatic cannulae (Ø mm 14/12/10/8/6), handle with manual aspiration regulation, tubing (in silicone, Ø mm 12x22) and connectors. All the pieces are autoclavable.  
 13-SET/E-MONO- Set for endouterine aspiration with disposable cannulae, made of: 6 adaptors (Ø mm 4/5/6/8/10/12), handle with manual aspiration regulation, tubing (in silicone, Ø mm 12x22) and connectors. All the pieces are autoclavable.  
 TS/12- Silicone autoclavable tube, Ø mm 12x22.  
 TS/E- See above - 2mt. long, with connector PT/12.

### CUPS FOR VACUUM-EXTRACTION

14-VE/MC- Set of metallic autoclavable cups, made of: 3 metallic cups (Ø mm 30/40/50) with membrane, handle, tyrants and tubing.  
 15-VE/SC60- Autoclavable silicone soft suction cup (Ø mm 60) with handle, tyrant and manual vacuum regulation.  
 15-VE/SC70- See above (Ø mm 70).  
 TS/7- Autoclavable silicone standard tubing, Ø mm 7x12.  
 TS/A- See above - 2mt. long, with connector PT/7.

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
## EC CERTIFICATE

Certificate No 187/MDD

### Full Quality Assurance System Approval Certificate

On the basis of our examination carried out according to Annex II, excluding section 4, of the Directive 93/42/EEC and its revised version, we hereby certify that:

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a quality assurance system ensuring the conformity of the following products:

**Surgical aspirators**

**Accessories for medical suction units and vacuum-extractors**

**Accessories for electrosurgical units**

**Electrosurgical smoke evacuators**

**Electrosurgical unit**

**Non electrical trolleys for electrosurgical equipments, electrosurgical smoke evacuators and medical suction units**

series and type refs in the Annex

with the relevant essential requirements of the aforementioned directive (from design to final inspection and testing) and it is subject to surveillance as specified in section 5 of Annex II. For class III devices, this certificate is valid only with the relevant EC Design-Examination Certificate of Annex II.4.

Reference to IMQ files Nos:

10A9800405; 10A9900240; 10AA00040; 10AA00016; 10AA00060; 10AA00156; 10AA00262; 10AB00169; 10AB000170; 10AC00210; 10AE00050; 10AF00011; 10AF00032; 10AF00192; 10AG00183; 10AH00132; 10AH00244; 10AI00034; 10AI00127; 10AI00122; 10AJ00093; 10EL00005; COMEDCONMHDM110021912-01; 10SN00082; DM16-0007485-01; 10AN00148; 10EO00012; DM16-0007485-01; DM17-0017739-01; DM18-0029069-01; DM18-0029585-01; DM18-0031303; DM19-0046045-01.

**This Approval Certificate is issued by IMQ S.p.A. as Notified Body for the Directive 93/42/EEC and its revised version. Notified Body notified to European Commission under number: 0051.**

Date: 1999-05-10  
 Updated: 2019-11-25  
 Substitution Date: 2018-11-07  
 Expiry Date: 2022-03-13

DocuSign

IMQ

## EC CERTIFICATE

Certificate No 187/MDD

### Annex

#### Surgical aspirators

Series: VACUMSOL Type ref. VACUMSOL AS-60; VACUMSOL AS-65; VACUMSOL AS-65/H.

Series: VORTEX-S Type ref. VORTEX-S AS-200/30LT; VORTEX-S AS-200/20LT VE; VORTEX-S AS-200/30LT VE; VORTEX-S AS-200/20LT.

Series: VORTECO Type ref. VORTECO AS-100/15LT; VORTECO AS-100/20LT; VORTECO AS-100/30LT; VORTECO AS-100/30LT VE.

Series: POLIVAC Type ref. POLIVAC B4/SLT 30 1; POLIVAC B4/SLT 30 2; POLIVAC B4/SLT 30 1-EXP; POLIVAC B4/SLT 30 2-EXP; POLIVAC B4/SLT 30 P-EXP; POLIVAC B4/SLT 50 E; POLIVAC B4/SLT 30 1VE-EXP; POLIVAC B4/SLT 30 PVE-EXP.

#### Accessories for medical suction units and vacuum-extractors

Type ref. As to annexed document "Elenco-AccAspir#A" dated 29.03.2005 attached; valid only if provided with IMQ stamp.

#### Accessories for electro-surgical units

Type ref. as to annexed document "AccElett-Elenco#D" dated 30.03.2008; valid only if provided with IMQ stamp.

#### Electrosurgical smoke evacuators

Series: CS Type ref. CS200; CS200LC; CS900LC

#### Electrosurgical unit

Series: ALSATOM Type ref. ALSATOM SU 50-MPC; ALSATOM SU 100-MPC; ALSATOM SU 140-MPC; ALSATOM SU 140/D-MPC; ALSATOM SU 140/BD-MPC; ALSATOM SU 140/D-MPC (100V-60Hz).

Series: MCDSe Type ref. EXCELL 200 MCDSe; EXCELL 250 MCDSe; EXCELL 350 MCDSe; EXCELL 400 MCDSe; EXCELL 250/A MCDSe; EXCELL 400/A MCDSe.

Series: NHP Type ref. EXCELL NHP 250/D; EXCELL NHP 350/D; EXCELL NHP 400/D; EXCELL NHP 250/DA; EXCELL NHP 400/DA; EXCELL NHP ENDOMED.

Type ref. MB 140-PRO  
Trade mark ALCYON

Type ref. BIPOLAR HF 100; BIPOLAR HF 140.  
Trade mark BUTTERFLY ITALIA

Series: NHP/T Type ref. EXCELL NHP/T-400; EXCELL NHP/T-200; EXCELL NHP/TA-400; EXCELL NHP/TA-200.

Type ref. PWT-400  
Trade mark TONTARRA

Date: 1999-05-10  
Updated: 2019-11-25  
Substitution Date: 2018-11-07  
Expiry Date: 2022-03-13

  
\_\_\_\_\_  
IMQ DocuSign

# EC CERTIFICATE

Certificate No 187/MDD

## Annex

### Electrosurgical unit

Series: TopTom Type ref. TopTom SU 50; TopTom SU 100; TopTom SU 140; TopTom SU 140/D.


Series: AMNOTOM Type ref. AMNOTOM 160 BASE; AMNOTOM 400 COMPACT; AMNOTOM 400 PREMIUM; AMNOTOM 400 PREMIUM A.

Trade mark AMNOTEC

### Non electrical trolleys for electrosurgical equipments, electrosurgical smoke evacuators and medical suction units

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VIA C. BONAZZI 16 - 40013 CASTEL MAGGIORE (BO)

E' CONFORME ALLA NORMA / IS IN COMPLIANCE WITH THE STANDARD

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*THE USE AND THE VALIDITY OF THE CERTIFICATE SHALL SATISFY THE  
REQUIREMENTS OF THE RULES FOR CERTIFICATION OF MANAGEMENT SYSTEMS*

<b>DATE:</b>	<b>PRIMA CERTIFICAZIONE</b>	<b>EMISSIONE CORRENTE</b>	<b>SCADENZA</b>
	<i>FIRST CERTIFICATION</i>	<i>CURRENT ISSUE</i>	<i>EXPIRY</i>
	1998-05-05	2019-10-31	2022-12-19

IMQ S.p.A.- VIA QUINTILIANO, 43 - 20138 MILANO ITALY  
Management Systems Division - Flavio Ornago

*IQNet, the association of the world's first class certification bodies, is the largest provider of management System Certification in the world.  
IQNet is composed of more than 30 bodies and counts over 150 subsidiaries all over the globe.*

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THE INTERNATIONAL CERTIFICATION NETWORK

# CERTIFICATE

*CISQ/IMQ has issued an IQNet recognized certificate that the organization:*

## **ALSA APPARECCHI MEDICALI SRL**

VIA C. BONAZZI 16 - 40013 CASTEL MAGGIORE (BO)

*has implemented and maintains a*

*Quality Management System*

*for the following scope:*

***Design, manufacture, placing on the market, technical service and contract manufacturing of electrosurgical units with or without ARGON gas and related accessories, electrosurgical aspirators and related accessories, electrosurgical smoke evacuators and related accessories***

*Further clarifications regarding the applicability of ISO 9001:2015 requirements may be obtained by consulting the organization*

*which fulfills the requirements of the following standard:*

## **ISO 9001:2015**

*Issued on: 2019 - 10 - 31*

*Expires on: 2022 - 12 - 19*

*This attestation is directly linked to the IQNet Partner's original certificate and shall not be used as a stand-alone document*

**Registration Number: IT - 1231**

*Alex Stoichitoiu*  
*President of IQNET*

*Ing. Claudio Provetti*  
*President of CISQ*

**IQNet Partners\*:**

AENOR Spain AFNOR Certification France APCER Portugal CCC Cyprus CISQ Italy  
CQC China CQM China CQS Czech Republic Cro Cert Croatia DQS Holding GmbH Germany FCAV Brazil  
FONDONORMA Venezuela ICONTEC Colombia Inspecta Sertifointi Oy Finland INTECO Costa Rica  
IRAM Argentina JQA Japan KFQ Korea MIRTEC Greece MSZT Hungary Nemko AS Norway NSAI Ireland  
NYCE-SIGE México PCBC Poland Quality Austria Austria RR Russia SII Israel SIQ Slovenia  
SIRIM QAS International Malaysia SQS Switzerland SRAC Romania TEST St Petersburg Russia TSE Turkey YUQS Serbia  
IQNet is represented in the USA by: AFNOR Certification, CISQ, DQS Holding GmbH and NSAI Inc.

www.imq.it

CERTIFICATO N. 9124.ALS2  
CERTIFICATE N.

IQNet, the association of the world's first class certification bodies, is the largest provider of management System Certification in the world.  
IQNet is composed of more than 30 bodies and counts over 150 subsidiaries all over the globe.

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*Design, manufacture, placing on the market, technical service and contract manufacturing of electrosurgical units with or without ARGON gas and related accessories, electrosurgical aspirators and related accessories, electrosurgical smoke evacuators and related accessories*

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THE INTERNATIONAL CERTIFICATION NETWORK

# CERTIFICATE

*CISQ/IMQ has issued an IQNet recognized certificate that the organization:*

## **ALSA APPARECCHI MEDICALI SRL**

VIA C. BONAZZI 16 - 40013 CASTEL MAGGIORE (BO)

*has implemented and maintains a*

*Quality Management System*

*for the following scope:*

***Design, manufacture, placing on the market, technical service and contract manufacturing of electrosurgical units with or without ARGON gas and related accessories, electrosurgical aspirators and related accessories, electrosurgical smoke evacuators and related accessories***

*Further clarifications regarding the applicability of ISO 13485:2016 requirements may be obtained by consulting the organization*

*which fulfills the requirements of the following standard:*

### **ISO 13485:2016**

*Issued on: 2019 - 10 - 31*

*Expires on: 2022 - 12 - 19*

*This attestation is directly linked to the IQNet Partner's original certificate and shall not be used as a stand-alone document*

**Registration Number: IT - 126301**

*Alex Stoichitoiu*  
*President of IQNET*

*Ing. Claudio Provetti*  
*President of CISQ*

**IQNet Partners\*:**

AENOR Spain AFNOR Certification France APCER Portugal CCC Cyprus CISQ Italy  
CQC China CQM China CQS Czech Republic Cro Cert Croatia DQS Holding GmbH Germany FCAV Brazil  
FONDONORMA Venezuela ICONTEC Colombia Inspecta Sertifointi Oy Finland INTECO Costa Rica  
IRAM Argentina JQA Japan KFQ Korea MIRTEC Greece MSZT Hungary Nemko AS Norway NSAI Ireland  
NYCE-SIGE México PCBC Poland Quality Austria Austria RR Russia SII Israel SIQ Slovenia  
SIRIM QAS International Malaysia SQS Switzerland SRAC Romania TEST St Petersburg Russia TSE Turkey YUQS Serbia  
IQNet is represented in the USA by: AFNOR Certification, CISQ, DQS Holding GmbH and NSAI Inc.