

STERISEPT PLUS

Disinfection and cleaning agent for surfaces of medical devices

- free from phenols, aldehydes and phosphates
- universal cleaning disinfectant for all water-tolerant surfaces
- broad spectrum of microbiological activity
- good cleaning properties and material compatibility

STERISEPT PLUS is a liquid concentrate for disinfection and cleaning of non-invasive medical devices and other water-tolerant surfaces. STERISEPT PLUS combines both cleaning and disinfection, which makes it a universal cleaning disinfectant for all surfaces in medical areas.

Due to specially selected active ingredients and surfactants, it has good material compatibility with good cleaning properties. STERISEPT PLUS is suitable for use in hospitals, food and pharmaceutical industries, veterinary areas, laboratories, dental clinics and in other medical institutions.

STERISEPT PLUS has a broad spectrum of bactericidal (incl. MRSA, VRE), (incl. Mycobacterium avium and Mycobacterium terrae), fungicidal and yeasticidal activity. The product is also active against viruses (incl. HBV, HIV, HCV, Vaccinia and Rota).



USAGE: Recommended dosage: 0,25 - 0,5 % depending on efficacy required (see table)

STERISEPT PLUS is concentrated. Prepare working solution by diluting the required concentration in water (at room temperature). As an example, dilute 5 ml of STERISEPT PLUS in 995 ml of tap water, to get 0,5% working solution. Apply a sufficient amount of ready-made working solution to treated surfaces by using wipe or cloth. Ensure complete coverage. The time of effect depends on microorganisms.

STERISEPT PLUS cannot be used together with aldehydes-based preparations, therefore avoid using products which contain aldehydes before and after disinfecting with STERISEPT PLUS.

STERISEPT PLUS is especially suitable to be used together with CP DRY WIPES DISPENSER SYSTEM for a convenient and user-friendly application. Impregnated wipes in the dispenser system are usable for 28 days.

MICROBIOLOGY:

EFFICACY	TEST METHOD	CONTACT TIME / CONCENTRATION		
		5 MIN	15 MIN	30 MIN
Bactericidal (incl. MRSA, VRE)	EN13727, EN14561		0,25%	
Yeasticidal	EN13624, EN14562		0,25%	
Bacteria and yeast	EN13697		0,25%	
Mycobactericidal	EN14348, EN14563		0,5%	
All enveloped viruses	EN14476	0,5%	0,25%	
Rotavirus	EN14476	0,5%		

INGREDIENTS: 100 g of product contains active ingredients: Didecyl-Dimethyl-Ammonium Chloride (DDAC) 12,5g; N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine 12,5g.

STANDARDS:

EN 13624	EN 13697	EN 13727	EN 14348	EN 14476	EN 14561	EN 14562	EN 14563
----------	----------	----------	----------	----------	----------	----------	----------



1000 ml / 5000 ml