Lyophilised, adsorbed anti-rabies nucleocapsid conjugate

REF 3572112

FOR THE HISTOLOGICAL DIAGNOSIS OF RABIES FROM BRAIN SMEARS BY IMMUNOFLUORESCENCE

Lyophilised: 4 x Q.S. 3 ml vials

I. COMPOSITION

IgG Immunoglobulins are:

- obtained by immunising rabbits with purified nucleocapsids from the fixed Pasteur strain
- purified by ion exchange chromatography
- coupled with fluorescein isothiocyanate (FITC)
- adsorbed (after titration) with 10% (final concentration) normal mouse brain suspension
- clarified by centrifugation and lyophilised

II. INTENDED USE

With the Rabies-specific labelled polyclonal Antibodies contained in the product which are targeted against the ribo-nucleoprotein complex of the Rabies virus, all currently known Lyssavirusgenotypes can be detected with the direct immunofluorescencetest (IFT) irrespective of the animal species tested.

Polyclonal antibodies, which are produced by the means of the Pasteur strain (PV) are known for their good specificity.

To ensure optimum performances, it is recommended to determine beforehand the most appropriate working solution (i.e. use reagent as undiluted or diluted).

The sensitivity of the anti-Rabies nucleocapsid conjugate may be reduced when testing Lyssavirus Genotype 3 (Mokola Virus), 5 (European Bat Lyssavirus, EBLV-1) and 6 (EBLV-2) strains. In this case the samples should be examined more carefully.

BIO-RAD

III. DIRECTION FOR USE

Reagent Preparation

Reconstitute each vial with 3 ml distilled water and centrifuge at 1500 rpm for 5 minutes for clarification.

Procedure

1. Fix the smears in acetone at -20°C for 30 minutes.
2. Add, on each smear, a sufficient quantity of clarified conjugate (i.e. 0.1 ml).
3. Incubate at 37°C for 30 minutes in a moist chamber.
4. Wash in two successive phosphate buffer (PBS) baths for 5 minutes each.
5. Apply a few drops of glycerine buffer.

Note: Counter staining in Evans Blue may improve the contrast and is sometimes useful for better fluorescent observation. A final 1/2000 dilution in Evans Blue is usually satisfactory.

- 1% Evans Blue solution - 20 µl (1 drop)
- clarified conjugate - 400 µl

Low positive samples may however be difficult to detect using Evans Blue.

Reading

Examine slides with a fluorescence microscope.

Controls

For easy diagnosis, we recommend to perform each day:

- a negative control with a normal mouse brain smear
- a positive control with a rabid mouse brain smear.

Other Reagents Required for the I.F. Reaction

- pH 7.2 phosphate buffer (10 X concentrated)
- glycerine buffer for I.F.
- Evans Blue: 1% stabilised solution
- Acetone resistant I.F. slides; 2 x 6 wells (7 mm)

Ask your Bio-Rad representative for more information on the availability of these reagents.

Available as

Kit 4 x q.s. 3 ml vial, lyophilised conjugate (code 3572112).

Each vial contains a quantity equivalent to 3 ml of adsorbed labelled immunoglobulins.

Storage

Storage at +2-8°C until the expiration date marked on the packaging.

After reconstitution the conjugate can be stored at +4°C protected from light for 15 days.

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