FTCE-CALI2-5-v5 (07/2020)

Français - FR

USAGE PRÉVU

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English - EN

INTENDED USE

This *in vitro* diagnostic device is intended for the calibration of ELITech Clinical Systems quantitative tests listed in the value sheet.
This *in vitro* diagnostic device is for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.
- Sodium azide < 0.1 % (w/w)
- Concentrations for each analyte to test are lot-specific.

MATERIALS REQUIRED BUT NOT PROVIDED

- ELITech Clinical Systems reagents listed in the value sheet.
- General Laboratory equipment.

TRACEABILITY

The traceability is indicated in the value sheet enclosed in the kit.

PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved or cleared in compliance with European Directive 98/79/EC, Annex II, List A. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- These products contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- Consult Safety Data Sheet (SDS) for a proper handling.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PREPARATION

- Carefully open the vial, avoiding the loss of lyophilizate.
- Add in exactly 3 mL of distilled or deionized water.
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particulates would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY

Prior to reconstitution :

- Stored at 2-8 °C and protected from light.
- Do not use after expiration date indicated on the vial label.

After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.

- Stability of the components:
Between 15 and 25 °C : 8 hours
Between 2 and 8 °C : 2 days
Between -25 and -15 °C : 4 weeks (when frozen once)

Exceptions:

- Stability of total bilirubin (when stored protected from light):
Between 15 and 25 °C : 6 hours
Between 2 and 8 °C : 1 day
Between -25 and -15 °C : 2 weeks (When frozen once)
- Stability of direct bilirubin (when stored protected from light):
Between 15 and 25 °C : 3 hours
Between 2 and 8 °C : 8 hours
Between -25 and -15 °C : 2 weeks (when frozen once)

Note: Store vials tightly capped and protected from light.

PROCEDURE

To use ELICAL 2, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

LIMITATIONS

Using ELICAL 2 has been validated with the ELITech systems (Analyzers and reagents listed in the value sheet).

Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit.
For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

ASSISTANCE TECHNIQUE

Contact your local distributor or ELITech Clinical Systems SAS (CCsupport@elitechgroup.com).

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❖BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Schreiber,W.E., *Iron and Porphyrin Metabolism*, Clinical Chemistry: Theory Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 755 and appendix.
- Higgins, T., Beutler, E. and Dumas, B.T., *Hemoglobin, Iron and Bilirubin*, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 509.
- Wu, A.H.B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 634.
- Higgins, T., *Clin. Chem.*, (1981), 27, 1619.
- Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), 59, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

**SYMBOLS/SYMBOLS/
SÍMBOLOS/SÍMBOLOS**
- Les symboles utilisés sont décrits dans la norme ISO-15223-1 hormis ceux présentés ci-dessous.
- Symbols used are defined on ISO-15223-1 standard, except those presented below.
- Los símbolos utilizados son descritos en la norma ISO-15223-1 a la excepción de los presentados a continuación.
- Os símbolos utilizados são definidos na norma ISO-15223-1, exceto os apresentados abaixo.

[CONT]	Contient Content Contiene Conteúdo
[R1]	Réactif R1 Reagent R1 Reactivo R1 Reagente R1
[R2]	Réactif R2 Reagent R2 Reactivo R2 Reagente R2
[Std]	Standard Standard Estándar Padrão
	Modification par rapport à la version précédente Modification from previous version Modificación con respecto a la versión anterior Modificação relativamente à versão anterior
	Conformité Européenne European Conformity Conformidad Europea Conformidade Europeia

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English - EN

❖INTENDED USE

ELITech Clinical Systems IRON FERENE is an *in vitro* diagnostic reagent intended for the quantitative determination of iron in human serum samples on analyzers or semi-automatic analyzers.
The standard is intended for the calibration of reagent. These *in vitro* diagnostic devices are for professional use only.

❖CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

In the body, 65 to 70 % of total iron is circulating in erythrocyte hemoglobin, 20 to 25% is stored in cells mainly as an iron-ferritin complex, less than 0.1% is transported by transferrin, the remaining iron being located in other proteins, including myoglobin. Serum iron levels are increased in hemochromatosis or liver damage. Lowered serum iron levels can be associated with inflammation, increased needs, a dietary deficiency or gastro-intestinal disorders.
Iron measurement is indicated to evaluate the body's iron content and to help diagnose iron deficiency or iron overload.



❖LIMITATION OF USE

The quantitative assay of iron alone can not be used to diagnose a disease or a specific pathology.
The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

❖METHOD & PRINCIPE ⁽⁴⁾

Ferene - End Point.

Iron is released from transferrin in acidic pH as a ferric ion Fe³⁺. It is then reduced by the ascorbic acid into ferrous ion Fe²⁺ and eventually forms a colored complex with ferene. The 578 nm absorbance of the iron-ferene complex is proportional to the iron concentration of the sample.

❖COMPOSITION

Reagent 1 : R1
Acetate buffer (pH 4.5) 1 mol/L
Reagent 2 : R2
Ferene 3 mmol/L
Ascorbic acid 240 mmol/L
Standard: Std (ref : FEFE-0600)
Iron 100 µg/dL
17.9 µmol/L

Also contains thiourea for optimal performance.

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8. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).

9. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).

10. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

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13. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).

14. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995)

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BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzyme, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders Company), (2006), 154
3. Dufour, R., *The liver: Function and chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 586, appendix
4. Ward, M.K., Cockayne, S., *Enzymology. Clinical Chemistry : Concepts and Application*, Anderson, S.C., Cockayne, S. (W.B. Saunders eds. Philadelphia USA), (1993), 238.
5. Schumann, G., et al *Clin Chem Lab Med.*, (2002), **40**, 725.
6. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations*, WHO/DIL/LAB/99.1 Rev.2, (2002), 21.
7. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
8. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
9. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
10. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

11. Interference Testing in Clinical Chemistry ; Approved Guideline—Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).

12. Vassault A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686

13. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immuno-globulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.

14. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).

15. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

- Higgins, T., et al., *Hemoglobin, Iron, and Bilirubin*, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds, Philadelphia USA), (2008), 509.
- Dufour, D.R. *The liver: Function and chemical Pathology*, Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc.), (2010), 586 and appendix.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
- Wu, H.B., *General Clinical Tests*. Tietz Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 172.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).8.
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd edition, AACC Press (1997).
- Young D.S., *Effects of drugs on clinical laboratory tests*, 4th edition, AACC Press (1995).

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BIBLIOGRAFIE

1. Panteghini, M. & Bais R., *Enzymes. Tietz Fundamentals of Clinical Chemistry*, 6th Ed, Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2008), 317
2. Sanhai, W.R. & Christenson R.F., *Protein Isoforms: Isoenzymes and Isoforms. Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 677 and appendix.
3. Apple, F. S., *Cardiac function, Tietz Fundamentals of Clinical Chemistry*, 6th Ed, Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2008), 614.
4. Boot, E., de Haan, L., *Massive increase in serum creatine kinase during olanzapine and quetiapine treatment, not during treatment with clozapine. Psychopharmacology (Berl)* (2000);150 :347-348.
5. Schumann, G., et al. *Clin Chem Lab Med.*, (2002), **40**, 635-42.
6. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline.* CLSI (NCCLS) document EP6-A (2003), **23** (16).
7. *Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline.* CLSI (NCCLS) document EP17-A (2004), **24** (34).
8. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition.* CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
9. *Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition.* CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
10. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
11. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
12. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
13. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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Stabilitatea la bord: 7 zile

BIBLIOGRAFIE

- Henderson, A.R., Donald W.M., *Enzymes*, Tietz Fundamentals of Clinical Chemistry, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 352.
- Sanhai, W.R., Christenson, R.H., *Cardiac and muscle disease*. Clinical Chemistry: Theory, Analysis, Correlation, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 566 and appendix.
- Neumeier, D., et al., *Clin.Chim.Acta.*, (1976), **73**, 445
- Schumann, G., et al., *Clin. Chem. Lab. Med.*, (2002), **40**, 635.
- Klein, G., et al., *Clin. Chem.*, (2001), **47**, Suppl. A30.
- Tietz N.W., *Clinical guide to laboratory tests*. 3rd Ed., W.B. Saunders eds. Philadelphia USA, (1995), 180.
- Vassault, A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686.
- Vassault, A., et al., *Ann. Biol. Clin.*, (1999), **57**, 685.
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995)

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FTRO-CRCO-v23(12/2018)_PIT-CRCO-4-v23

BIBLIOGRAFIE

1. Allston, C.A., *Non protein nitrogenous compounds and renal function*. Clinical Chemistry: Concepts and Application, Anderson, S.C., Cockayne, S. (W.B. Saunders eds. Philadelphia USA), (1993), 369.
2. Newman, D.J., Price C.P., *Non protein nitrogen metabolite*. Tietz Fundamentals of Clinical Chemistry, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 414.
3. Butler, A.R., The Jaffe reaction. *Identification of the coloured species*. Clin. Chim. Acta., (1975), 59, 227.
4. Vasiliades, J., *Reaction of alkaline picrate with creatinine. 1. Kinetics and mechanism of formation of the mono-creatinine picric acid complex*. Clin. Chem., (1976), **22**, 1664.
5. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
6. Tietz, N.W. Clinical guide to laboratory tests, 3rd Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 186.
7. Vassault A., et al., Ann. Biol. Clin., (1986), 44, 686.
8. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
9. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd edition, AACC Press (1997).
10. Young, D.S., Effects of drugs on clinical laboratory tests, 4th edition, AACC Press (1995).

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BIBLIOGRAFIE

1. Lamb, E.J., & Price C. P., *Creatinine, Urea, and Uric Acid*. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R. & Bruns D. E. (W.B. Saunders eds. Philadelphia USA), (2008), 363.
2. Wu, A.H.B., Tietz Clinical guide to laboratory tests, 4th Ed, (W.B. Saunders eds. Philadelphia USA), (2006), 1096.
3. Kaplan, J.M. & First, M.R., *Renal Function*. Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds St Louis USA), (2010), 567 and appendix.
4. Bretaudière, J.P., et al., *Direct Enzymatic Determination of Urea in Plasma and Urine with a Centrifugal Analyzer*. Clin. Chem., (1976), 22, 1614.
5. Fawcett, J.K., Scott, J.E., *A rapid and precise method for the determination of urea* J Clin Pathol., (1960), 13:156-9.
6. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI (NCCLS) document EP6-A (2003), **23** (16).
7. *Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
8. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).

Referințe:

URSL-0250	8 x 25 mL
URSL-0455	4 x 55 mL
URSL-0400	2 x 62,5 mL
URSL-0420	4 x 62,5 mL
URSL-0500	5 x 125 mL

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R1 8 x 20 mL + R2 8 x 5 mL
R1 4 x 44 mL + R2 4 x 11 mL
R1 2 x 50 mL + R2 1 x 26 mL
R1 4 x 50 mL + R2 2 x 26 mL
R1 5 x 100 mL + R2 1 x 127 mL

FTRO-URSL-v20(12/2018)_PIT-URSL-4-v20

9. *Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

10. *Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP7-A2 (2005), **25** (27).

11. Vassault A., et al., *Ann. Biol. Clin.*, (1986), 44, 686.

12. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.

13. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd edition, AACC Press (1997).

14. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th edition, AACC Press (1995).

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