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
GE Healthcare

<http://www.gehealthcare.com/>

SIGNA Pioneer 3.0T MR 30

CONFIDENTIAL

**SECRET DE
AFACERI**

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I. SIGNA Pioneer 3.0T MR 30 Fighter EMEA

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**SECRET DE
AFACERI**

1.4 MR workflow concept

The MR system's role in the workflow for an exam is to acquire data and create images for analysis of the patient's condition. The operator needs to be able to start with the patient's requested procedure, execute the set of acquisitions and post-processing activities to fulfill that request and then transfer/export the resulting data.

Figure 1-4 MR exam workflow


Table 1-12 MR exam workflow

Description
<p>Prepare the patient Prior to bringing the patient into the scan room, screen the patient.</p>
<p>Patient handling (1) The patient is brought into the MR scan room and prepared for the exam.</p>
<p>Patient Registration (2a and 2b) The patient is entered into the Worklist Manager through either a RIS system or manually entered and a protocol is attached to the patient's exam.</p> <p>The scan data is acquired (2c) A protocol must be attached to the exam before you can start exam.</p>
<p>View and analyze images (3) After the scan data is acquired the images can be viewed and filmed (3a) and analyzed (3b) using one of several post processing applications.</p>
<p>Store images (4) The exam can be networked to be analyzed on an AW workstation (4a), stored and analyzed on a PACs system (4b) or simply stored on a CD/DVD (4c).</p>

CONFIDENTIAL

SECRET DE AFACERI

- Preventing patient warming is one of the most important safety measures you must take into consideration as you prepare a patient for an MR exam. Appropriate RF padding and proper patient positioning are the most effective means of preventing injury related to RF heating. The following are a few golden rules to remember as you position and pad your patients:
 - Only use GE HealthCare-approved RF padding.
 - Use non-conductive padding that is at least 0.25 inches (0.635 cm) thick between the patient's skin and the magnet bore.
 - Appropriate padding must be used EVERY time without exception
 - Sheets and gowns are not a substitute for approved RF padding.
 - Never allow your patient's skin to come in direct contact with the scanner bore or any surface coil or cable.
 - Never allow skin-to-skin contact.
 - If a patient does not fit in the MR scanner bore with the required padding, another modality should be used to scan the patient.
 - While some of these rules may seem a little tough to follow at times, remember that RF injury, which can in extreme cases include burns such as the one you see below, can happen very quickly and your patient may not have time to warn you in time to prevent an injury.
4. Only use approved RF coils that are not damaged.
- Labels such as the one in the figure below, provide warnings about working with RF coils.

Figure 2-42 Warning label

- Before using the coil, check the integrity of the electrical insulation of the components or accessories of the device.
5. Keep electrically conductive material that must remain in the magnet bore from directly contacting the patient by placing insulation between the conductive material and the patient.
- Place a clean cotton sheet over the coil and comfort pad so the patient's skin does not come in contact with the coil or the comfort pad.
6. Position RF cables down the center and directly out of the bore, without looping or crossing the cables.
- Route the cables so there are no loops in any cables in the magnet. Cable holders are located on both side of the cradle near the edges. Use the appropriate gating cable for surface coil imaging.

- Use only MR system recommended monitoring equipment, ECG leads, wires, electrodes, and other components and accessories.
 - For medical devices that are labeled as MR Safe or MR Conditional consult the device manufacturer's documentation.
 - Follow all instructions for the proper operation of physiologic monitoring or other equipment provided by the manufacturer of the device.
7. Test the patient intercom.
 - Make sure that the patient can hear you and you can hear the patient.
 8. Enter the correct patient weight.
 - Correct weight entries maximize performance and help prevent excessive RF exposure.
 9. Turn on the bore light and fan.
 - Lights turned on inside the bore can help alleviate feelings of claustrophobia.
 - A fan inside the magnet bore provides adequate air movement for the patient. Keep the fan on at all times.
 10. Show the patient how to use the Patient Alert System.
 - Patients experiencing uneasiness or concern can squeeze the Patient Alert bulb.
 - When the Patient Alert bulb is squeezed, an alarm emits a signal.
 11. Respond to bore temperature direction messages throughout the procedure.
 - The temperature messages are located in the message window on your console.
 - If the patient reports feeling warm, discontinue the procedure.

2.13.5 Patient padding

Figure 2-43 Elbow RF burn

- The following are a tips that will assist you in properly positioning and applying RF padding to your patients. Should you need more information on prevention of patient warming than what is provided here, refer to your surface coil and refer to Tissue Heating in this manual. If you need help beyond the documentation please do not hesitate to reach out to your local Applications Specialists.
- Although the photos in the following sections are from a Discovery system, the safety padding guidelines apply to all MR systems.

Whole body padding

- An important consideration when padding your patients is that you will need to double check the position of the pads once the patient is in the bore. Table movement may dislodge padding and expose skin to the scanner bore.

Figure 2-44 Padding between patient and bore. 1 = bore pads

- Notice that padding is positioned not only at the patient's sides to prevent their arms from touching the bore, but that padding is also placed between the hands and thighs and between knees and ankles to prevent forming conductive loops.

Figure 2-45 Patient padding**Surface coil padding**


- Padding with a surface coil presents different challenges from a patient RF padding perspective.
 - First rule of thumb is to remember to use all manufacturer provided padding to prevent motion and the patient's skin from coming in contact with the coil, and to also use additional padding if appropriate to secure an opposing extremity to prevent contact with the coil which could also lead to burns or motion artifacts.
 - Just as with the whole body RF padding demonstration, you'll need to make certain that the patient's skin does not come into contact with the scanner bore and that padding is placed between the hands and thighs to prevent conductive loops.

Figure 2-46 Extremity padding


- A final safety consideration for surface coils is to ensure that the patient does not come into contact with the coil cable, therefore you may need to use additional RF padding to protect the patient.
 - Care should also be taken to ensure the cable is not looped in the bore and that it is routed down the center of the scanner bore.

Figure 2-47 Coil cable with no loop**Cardiac coil padding**

- Follow your basic padding recommendations to prevent contact with the scanner bore and prevent conductive anatomical loops, but there are a couple of additional steps you'll need to take to ensure patient safety
 - The cardiac coil does not require additional RF padding to be placed between the patient and the anterior coil component, but you should use the manufacturer's pad on the posterior component of the coil for patient RF protection. You should also cover the patient with their gown before placing the anterior component of the coil and make certain both the anterior and posterior elements are in alignment.

Figure 2-48 1 = coil pad aligned with coil, 2 = sheet to cover pad, 3 = anterior and posterior coil elements aligned


- Secure the coil snugly, but comfortably with the straps.

Figure 2-49 Coil secured with straps

- As is the case of all surface coils ensure that the cables do not come in contact with the patient and that they are not looped and routed down the center of the bore. As you can see there is significantly more cable that we need to isolate from the patient, so be sure to use as much padding as needed.

Figure 2-50 1 = Pad placed between cable and patient skin

- If you are using the cardiac coil, it's likely you are also using the ECG leads and cable. The rules for the ECG cable are the same as the coil cable. Route the ECG cable down the center of the bore, do not loop the ECG cable and do not allow it to come in contact with the coil cable.

Figure 2-51 1 = ECG cable with no loops

2.13.6 Protect the patient's eyes and ears

About this task

Hearing protection is required. During a scan, all persons in the scan room are required to wear proper hearing protection to avoid possible hearing impairment. Acoustic levels may exceed 99dBA. Properly worn hearing protection is required for all people in the magnet room during a scan, including any MR workers.

Patients must close his or her eyes when the alignment light is on during positioning. Follow these guidelines to ensure proper eye and ear protection for your patient.

Step-by-step instructions

1. Provide the patient with hearing protection.
 - Earplugs or a headphone system with stereo music. For details, see [Acoustic Noise on page 189](#).
 - Earplugs reduce the intensity of the sound, while allowing your patient to hear normal conversations.
 - Headphone systems soften acoustic noise, but may impede verbal communication with patients while the system is operating.


Table 2-44 Disposable ear protection

Description	dB
E8801BA EAR Disposable Foam Earplugs	29
E8801BB EAR Taperfit2 Foam Earplugs	32
E8801BC Max-Lite Foam Earplugs	30

⚠ WARNING

Hearing protection is required for all people in the magnet room during a scan to prevent hearing impairment. Acoustic levels may exceed 99 dB(A)

2. Make sure that the hearing protection device is worn properly.
 - Earplugs should be comfortable for the patient and inserted fully. Pliable earplugs compress when they are rolled between the fingers and conform to the ear after they are inserted.

Figure 2-52 Hearing protection

- The headphone system should be audible and comfortable for the patient.
3. Instruct the patient to close his or her eyes when the alignment light is on.
 - The Laser Alignment Lights for patient positioning can cause eye injury.

⚠ CAUTION

Turn off the laser light after positioning the patient.

lassitude, loss of strength, asthma, nausea, vomiting, headache, fever or hypothermia, anorexia, loss of sense of smell, diarrhea, anuria, liver damage, jaundice and convulsions.

Animal experiments have resulted in observable birth defects.

Cancers of the lung and nasal sinuses in Nickel workers have been known for more than 50 years to be associated with nickel refining, nickel plating, and nickel polishing.

First Aid:

- If this chemical gets in the eyes, immediately flush the eyes with large amounts of water, occasionally lifting the lower and upper lids. Get medical attention immediately. Contact lenses should not be worn when working with this chemical.
- If this chemical gets on the skin, immediately wash contaminated skin with soap or mild detergent and water. If this chemical soaks clothing, immediately remove clothing and wash contaminated skin with soap or mild detergent and water. Get medical attention promptly.
- When this chemical has been swallowed and person is conscious, immediately give the person large quantities of water or milk. Remove by gastric lavage unless patient is vomiting. Do not make an unconscious person vomit. Get medical attention immediately.

Spill, leak, and disposal procedures

Observe all federal, state and/or local regulations when storing or disposing of this substance. Contact local and/or state environmental authorities to insure proper compliance.

This substance does not meet the definition of a hazardous waste as defined by the Resource Conservation and Recovery Act (RCRA) (40CFR260).

4.12.2 Scan the phantom

About this task

The DAQA procedure provides a means for you to track overall system or RF coil functionality. The application supports all GE HealthCare coils that have the coil ID feature.

Important

To obtain meaningful, reproducible data for a given RF coil, consistency and repeatability is key. Always use the same phantom, the same positioning of the phantom in the coil and the same landmark at the same location on the phantom/coil.

Table 4-5 Phantoms used in DAQA tests

Phantom name
Head TLT
Head TLT Phantom Loader
Body TLT phantom sphere
Body TLT Loader

Step-by-step instructions

1. Press the **Alignment light** button on the Magnet control panel.

CAUTION

Exposing eyes to laser alignment lights may result in eye injury.

Do not stare directly into the laser beam.

Instruct patients to close their eyes to avoid eye exposure to the alignment light.

Closely monitor all patients and prevent them from accidentally staring into the beam. Do not leave the laser beam on after you position the patient.

2. Press the table movement buttons to advance the cradle until the axial alignment light rests at the desired landmark. Confirm centering with the sagittal alignment light.
3. Press **Landmark** button from the Magnet control panel. The cradle position reads zero.
4. Make sure all health lines are long enough to accommodate movement and then press **Advance to Scan** to move the cradle to magnet isocenter.
 - The alignment lights automatically turn off.
5. Adjust the in-bore light and fan.
 - For details see [5.11 Adjust the magnet fan or light levels on page 309](#).
6. Leave the scan room and enter the console room to begin scan prescription.
 - Close the scan room door during the acquisition to prevent RF leaks.
 - An open door prevents scanning.

6.5 Patient return to landmark

About this task

Once you have started scanning, you can remove the patient from the magnet bore and then return the patient to the scan position without losing the landmark. This may be done to prep the patient for another phase of the exam, give the anxious or claustrophobic patient a break between acquisitions, etc.

Step-by-step instructions

1. When the system is in between scans, press the **Out** button (fast or slow) to bring the patient out of the magnet bore. For table movement details, see:
 - [Magnet controls considerations on page 307](#)
2. When the patient is ready to be placed back into the magnet, press the **Back to Landmark** button.

6.6 Patient alert

Your MR system has a patient alert system that enables the patient to alert you at the console by squeezing a bulb. It is comprised of a squeeze bulb and a control box, typically located at the operator console.

About this task

Figure 6-9 Patient alert system



Table 6-2 Patient alert system image legend

Number	Description
1	Patient alert bulb attachment area
2	Control box

Squeezing the Patient Alert bulb causes the control box to light up and emit an audible signal. A switch on the control box allows you to set the signal for intermittent or constant light and sound.

Your MR system also has an intercom system that enables you to maintain verbal contact with the patient throughout the examination.

CAUTION

Provide all patients with the Patient Alert bulb. This can be especially important for procedures that require the concerted attention of the technologist/operator at the MR or Advantage Workstation (AW) operator console, e.g., fMRI sequences.

The patient alert system is comprised of a squeeze ball that is connected to the PAC and a control box. The control box is typically located at the operator console.

The patient alert port is a Type BF applied part.

Use the following steps to set the patient up with the alert bulb.

Step-by-step instructions

1. Give the patient the alert bulb.

The bulb is a rubber product and not latex.

2. Instruct the patient to loosely hold the alert bulb and to only squeeze the bulb to bring you into the scan room to consult with the patient and attend to his needs.

Figure 6-10 Patient holding alert bulb


3. If the patient squeezes the rubber ball end of the alert system, a loud sound is heard in the control room. Follow these steps.
 - a. Press **RESET** to stop the alarm and reactivate it.
 - b. Go into the scan room and attend to the patient's needs.

Figure 6-11 Patient Alert control box: Steady/Pulse button on left and Reset button on right

Figure 7-17 The invalid state is only confirmed after the series have been saved in RXd state

- If the AutoBind task displays an **INV** (invalid) state and a message displays that the AutoBind requirements are not met, then the series will not bind. The invalid state is only confirmed after the series have been acquired, thus it is very important that all series meet the AutoBind requirements.
10. Click **Save Rx** for each series.
 11. Click **Scan** to acquire the data for each series. Alternatively, click **Auto Scan > Scan**.
 12. When all the series are completed, and the AutoBind post-process task is in the Done state, select the bound series. From the Workflow Manager control panel, click **View**.

7.1.10 Scan with Auto Paste

Before you begin

- A scan session is open.
- Localizers and other series have been acquired.

CAUTION

Images labeled as water may include signal from fatty tissue, and images labeled as fat may include signal from water. This error may occur in regions of high magnetic field variation, in spatially isolated tissue, due to patient or tissue motion, due to phase wrap artifacts, due to similar water and fat intensity in T2-weighted Flex scans, due to TE values beyond recommended limits, and/or in images with low signal-to-noise ratios. The presence of fat tissue in images labeled as water, or vice versa, may occur within single images or throughout an entire stack of slices. By default, both sets of images (labeled fat and labeled water) will be reconstructed and inserted into the database for review. Proper calibration and center frequency selection will reduce the occurrence of this error. Complete elimination of this error may not be possible and thus interpretation of MR images must be completed by trained personnel.

Results**Image annotation**

M3D/VIBRANT/flip angle

Water and Fat images are synthesized from collected in-phase and out-of-phase images; thus TE values for Water and Fat images are the average of IN-phase and Out-of-Phase TEs. In-phase and out-of-phase gets annotated according to the TE and TE2 values on screen.

8.3 Scan Anatomy: Cardiac

8.3.1 Cardiac introduction

Cardiac Gating/Triggering allows the same data to be collected at the same points of the cardiac cycle for each repetition of the pulse sequence allowing motion artifacts to be reduced, or sometimes eliminated. The benefit is that the motion is always the same and is essentially “frozen,” allowing motion artifacts to be reduced or, sometimes, eliminated. Gating and triggering of MR images can be performed with ECG or Peripheral gating techniques.

Table 8-18 Cardiac applications

Scan plane	Cardiac application
Short Axis	Assesses two chambers, either both atria or both ventricles.
2-Chamber Long Axis	Visualizes left atrium mitral valve and left ventricle.
4-Chamber Long Axis	Assesses septal wall and visualizes the four chambers, mitral, and tricuspid valves in a single plane.
3-Chamber View (left ventricular out-flow tract)	Assesses left atrium, mitral valve, left ventricle, aortic valve, and aortic root in one plane.

8.3.1.1 Cardiac Gating/Triggering

Cardiac gating/triggering is used for:

- imaging the heart's structure and function

CONFIDENTIAL**SECRET DE
AFACERI**

- imaging in the thorax
- Cine-PC arteriography to examine flow

There are two types of gating/triggering:

- **ECG or VCG:** Uses the heart's electrical activity as reported by the ECG. Electrode placement is critical to a successful ECG scan. Adjust the electrodes until a strong signal is received. If the signal is lost during the acquisition, it aborts the scan. For details, see [Set up the ECG gated exam on page 852](#).
- **PPG or Peripheral:** Uses a photopulse sensor that picks up the mechanical action of blood pulsing through the body. For details, see [Set up the peripheral gated exam on page 863](#).

Figure 8-313 Cardiac gated heart

Peripheral Gating is used for head, cervical, and thoracic spines to minimize pulsatile motion of CSF on T2*- or T2-weighted scans. It is safe for patients due to the fiber-optic technology.

There are three types of cardiac gated sequences:

1. **Prospectively gated retrospectively reconstructed multi-cardiac phase Cine or retrospective Cine**

The prospective gating means that new k-space lines are acquired every time the system receives a new trigger. Retrospective reconstruction means that the number of TRs per RR interval (heartbeat) can vary and data is sorted in reconstruction (retrospectively or after the scan is done). Data acquisition is done continuously.

Figure 8-314 ECG waveform with Retrospective CINE


Table 8-19 Image legend

Legend	Description
RR	R-peak to R-peak of ECG cycle
ARW	Arrhythmia Rejection Window
1	Yellow bar represents data acquisition

For retrospective Cine, when ECG or VCG gating is used, the cardiac phases are reconstructed from R peak. For peripheral gating or third-party gating, reconstruction is from trigger detection. For details, see [9.7.10.15 2D CINE on page 1519](#).

2. To open/close a toolbar menu, click the menu item button to open the menu, and then click the button to close it. For example, click **Shim** to open the shim menu, and click **Shim** again to close the shim menu.
3. To close the Graphic Rx toolbar, click the  in the upper-right corner of the toolbar.

9.4.2 AIR x™

AIR x™ uses deep learning algorithms that automatically identify anatomical structures to graphically prescribe a slice range for scans based on specified anatomical references. It is a scan feature that encapsulates automated calculation and management of geometry data from localizer images. It stores that computed geometry data and organizes it by anatomical reference markers and scan planes.

AIR x™ uses the 2D 3-Plane localizer acquired as the first series to automatically position slices for subsequent series based on the orientation and location defined in the protocol.

AIR x™ benefits

- Increased productivity by simplifying workflow steps and thus reducing scan prescription time.
- Improved prescription consistency from one MR technologist to the next.
- Automatic adjustment of slice prescription (reference line and scan plane) for various patient positions, anatomy, and pathology.
- Improved consistency with patients that are likely to have repeated scans over a time period.

Considerations

- Anatomical Reference options can be selected for a task and then linked with other series.
- Anatomical Reference options are reset to None when the following occurs:
 - Copy Rx is selected or linked to copy geometry data from some other task.
 - The graphic prescription computed by AIR x™ is manually adjusted. Note that the AIR x™ prescription can be reset by re-selecting the option from the Anatomical Reference menu.
 - An error occurs during AIR x™ processing.
- Oblique plane must be selected from the Scan Plane menu located on the scan parameter screen.

Figure 9-106 Scan Plane menu

9.4.2.1 AIR x™ brain

Use these steps to acquire an AIR x™ brain scan.

About this task

Consider the following before acquiring an AIR x™ brain scan:

- AIR x™ calculations may be compromised for patients with an artificial eye.
- The patient position must be supine.

Step-by-step instructions


1. From the *Protocol* screen, select a head category, and then select a head protocol.
2. Acquire a whole brain, 2D 3-plane FGRE, or SSFSE localizer. SSFSE localizer is preferred to ensure the best results. The 3-plane localizer must meet the following criteria:
 - Frequency and phase resolution ≥ 192
 - Slice thickness ≤ 15 mm
 - Slice coverage for all axial, sagittal, coronal planes ≥ 10 cm. When the HNU comfort tilt is used, or the patient's head is rotated to the left or right, please use slice coverage ≥ 16 cm to ensure best results.
 - Slice gap ≤ 5 mm
 - TE for SSFSE ≤ 80 ms
3. From the Workflow Manager, select **Task 2** and click **Setup** to review the scan parameters.
4. To view the *Details* screen, click the  0:40  of the scan parameter screen.
5. From the Anatomical Reference menu on the *Details* tab, select an option:
 - None
 - ACPC
 - Temporal Lobe
 - IAC
 - OML
 - Orbits
 - Optic Nerves
 - Pituitary
 - Circle of Willis

Figure 9-107 Anatomical Reference menu


6. From the Anatomical Plane menu on the *Details* tab, select an option. Note that the scan plane options are based on the selected Anatomical Reference.
 - Axial
 - Sagittal
 - Sagittal Left (only appears if Optic Nerves is selected from Anatomical Reference menu)
 - Sagittal Right (only appears if Optic Nerves is selected from Anatomical Reference menu)
 - Coronal

9.5.32 Multi Station

Multi Station is a GE HealthCare protocol that is an alternative to QuickSTEP. It provides automatic table movement and switching of coils between stations for peripheral vascular run-offs. It allows you to prescan at multiple stations to optimize image quality and it properly annotates image locations based on the landmark. Therefore, this protocol allows imaging of multiple phases of contrast enhancement at multiple anatomical locations typically with one contrast injection. It can be performed in conjunction with a bolus triggering technique, such as SmartPrep on page 1176 or Fluoro Trigger on page 1222.

Figure 9-257 Multi Station runoff post processed in Add Subtract

The Multi Station sequence is acquired following the optional acquisition of a localizer meta-series. A meta-series consists of multiple series linked together. The localizer meta-series is a 3-Plane localizer that acquires several series, at different locations, across the anatomy of interest. The localizer images are used for defining the scan locations for the sequence. It is easier to save the localizers as your sites protocols when they are separate series, i.e., not linked. The meta-series can be prescribed from any localizer series. You may choose to acquire a TOF acquisition as a localizer rather than a 3-Plane Localizer acquisition.

The mask, arterial or venous sequences are also meta-series, and include one to four 3D (Fast TOF GRE or Fast TOF SPGR) series with unique scan locations. In the meta-series, the Number of Stations is defined at the Multi Station screen. The Number of Stations indicates the number of separate series that will be scanned during the series. When the Multi Station screen is saved in protocol edit session, the Workflow Manager lists additional new series identical to the series just saved.

Once all series within the meta-series are prescribed and saved, all the series are then prescanned. The prescan values are saved by the system until the station associated with those prescan values is ready to be scanned. When prescan has been completed for any of these stations, the Workflow Manager changes the state from RxD to PSCD. If they have not been prescanned, they will be listed as RXD as long as they have been saved.

At the time of contrast injection, the meta-series is initiated with the Scan enabled. For Fast 3D TOF sequences, the first station can be initiated using SmartPrep or Fluoro Trigger. When the bolus is detected, the 3D volume at station 1 is acquired. Upon completion of the first station, the table moves and the volume at station 2 is acquired. The system automatically moves the table and initiates data acquisition until all stations are acquired.

With Multi Station:

- T2 FLAIR results in bright signal suppressed from CSF on T2-weighted image structures adjacent to fluid filled structures, therefore becoming more apparent.
- PROPELLER T2 FLAIR imaging provides robust T2-weighted imaging with CSF suppression. When CSF is suppressed, you get better visualization of gray and white matter.
- PROPELLER Brain T2 FLAIR reduces patient motion artifact compared to traditional T2 FLAIR with comparable scan time.

9.6.56 T2 Prep

T2 Prep applies a sequence of non-slice selective 90° and 180° RF pulses to suppress cardiac muscle tissue.

T2 Prep with MRCP improves suppression of signal from background tissues.

It is not necessary to program a TI time for the inversion time – it is automatically calculated.

9.6.57 Tailored RF

TailoredRF is used to stabilize the echo amplitudes in the FSE sequences (Spin Echo and Inversion Recovery). Tailored RF produces images with the following characteristics: less blurring, slightly less SNR, flatter contrast for T2s, and slightly more slices per TR.

Figure 9-432 Tailored RF


Table 9-129 Image legend

Number	Description
1	Tailored RF On
2	Tailored RF Off

By varying the flip angle of the first two refocusing pulses and keeping the remaining refocusing pulses the same, the echo amplitudes can be made fairly constant. The result is a reduction of blurring in the image.

g. **Synthetic b-value list:** type and enter the desired b-value for each number. The Synthetic b-value list does not display under the following conditions.

- When the # of Synthetic b-values = 0.
- When the # of T2 images = 0 and the # of b-value = 1.
- When there is 1 b-value less than 2500, for example b= (0,1000), b= (100,2000).

3. Click **Save Rx > Scan**.

- The source and synthetic DWI images, and ADC and eADC maps prescribed as post-process tasks are displayed in the Patient List.
- The MAGiC DWI series number is the original series # x 100.

Figure 9-456 Patient List with DWI source prospective series selected

Series	Type	Images	Description
1	PROSP	27	3Plane Loc SSFSE
2	PROSP	12R	Cal RHRRRAIN
3	PROSP	72	Ax DWI
300	DIFFUSION	36	DW_Synthetic: Ax DWI
350	DIFFUSION	12	ADC (10^-6 mm^2/s):Ma...
351	DIFFUSION	12	eADC:May 03 2016 13-2...


4. To post-process the images, see Synthetic DWI workflow.

9.7.4.5 Diffusion Tensor EPI

DW EPI Tensor or DTI is a technique that produces image contrast proportional to the local diffusion coefficient of water. Both the diffusion coefficient and its directional dependence can be measured using DTI. Data can then be used to image the directional dependence of the local diffusion coefficient in the tissue.

DTI pulse sequence consists of a 90° excitation pulse and a 180° RF. It also includes a pair of DW gradients, one placed before the 180° RF, the other immediately after. This pair of gradients is applied not only for the slice gradient axes, but the readout and the phase gradient axes as well.

Figure 9-457 DW EPI Tensor PSD

Use DTI for brain applications. It is most commonly used for white matter tract visualizations. Visualizing the white matter tracts in the brain can be useful in monitoring patients with stroke, epilepsy, various brain trauma, neonatal brain development, as well as disease management. Up to 50,000 images can be acquired within a single DTI series.

Figure 9-539 Raw k-space data and transformed data


Table 9-183 Image legend

Number	Description
1	Raw k-space data
2	Transformed data (right)

Data Processing

PROPELLER's ability to effectively reduce motion and susceptibility artifacts hinges in large part on the amount of data collected during a PROPELLER scan.

The redundant data collected in the center of k-space makes it possible for PROPELLER to perform several correction steps before final image reconstruction. After the initial signal is obtained, PROPELLER executes phase correction. Then the program performs three additional correction steps: rotation correction, translation correction, and correlation weighting.

All this additional data requires considerably more processing. PROPELLER uses intensive multi-channel image reconstruction and processing techniques, since it uses five times more processing steps than a conventional DWI acquisition.

9.7.7.1 PROPELLER scan parameters selections



Consider this information when modifying the scan parameters for all types of PROPELLER scans. For specific scan parameter values, select a protocol from your GE HealthCare or Site library.

Typically use PROPELLER to acquire brain, spine, liver, pelvis, shoulder, knee and wrist images.

CONFIDENTIAL

SECRET DE AFACERI

Table 11-6 Remote/Local value explanation (Table continued)

Remote/Local examples	Description
<p>Figure 11-30 No images transferred</p> 	<p>No images transferredIn this example, a "0" in the series Local column indicates that no images transferred for that series. A "0" in the exam Local column indicates that no images from the exam transferred.</p>
<p>Figure 11-31 Partial image transfer</p> 	<p>Partial image transferIn this example, the value in the Remote column is less than the value in the Local column, which indicates the following:</p> <ul style="list-style-type: none"> • 4 of the total 23 images in the series transferred. • 4 of the total 174 images in the exam transferred.

11.6 Archive

The archive system only allows a networked device such as a PACS for archiving images. Images stored on these devices are off-line, meaning that in order to display the images you must first restore them to the system disk. Only DICOM images can be archived and restored on your system.

There is a difference between archiving images and simply saving them. An archive device is a network device set up as an archive destination by your service personnel. When an image is archived, an

CONFIDENTIAL

SECRET DE AFACERI

Table 11-15 Enterprise Authentication / Authorization (Table continued)

Specification	Input
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	Microsoft Active Directory / Novell eDirectory
Ports (default)	N/A – outgoing communication initiation only
Traffic characterization and Bandwidth Requirements	On demand, local user initiated. The bandwidth is dependent on the local site.
Latency max	N/A

11.17.3 Required characteristics and configuration for network support of MR scanner

The network must meet the specific requirements above for all traffic flows associated with the subset of features, use cases and workflow required by the responsible organization's users.

In addition, the network must be "flat" (limited to a single IP broadcast domain).

- Connection of the MR System to a network that includes other equipment could result in the loss of supported DICOM functionality if network does not meet the defined specifications.
- The hospital should ensure proper network configuration.
- Subsequent changes to the network could introduce new risks and require additional analysis.
- Changes to the network includes:
 - Changes in the network configuration.
 - Connection of additional items to the network.
 - Disconnecting items from the network.
 - Update of equipment connected to the network.
 - upgrade of equipment connected to the network.

11.17.4 Remote host parameters

The Network function supports DICOM networking. When adding or updating a remote list, there are some new parameters needed. All of the following information, except for Comments, needs to be provided in order to set up a remote host:

- The Host name to be entered is the name of the device. If the device is DICOM, the name must match exactly to the name given to the device.
- The Network address of the device is a static IP address that is provided by the institution's network administrator.
- The Network protocol is DICOM. If the MR system will be sending to this device, the device must be DICOM and the DICOM network protocol must be selected.
- The Port number is unique to the device. If the device is an Advantage Windows workstation or another GE HealthCare MR system, the number will be 4006.
- The AE Title is unique to the device. If the device is an Advantage Windows workstation or another GE HealthCare Systems system, the AE Title will be the same as the Host name.

CONFIDENTIAL

**SECRET DE
AFACERI**

- The Comment field allows you to input a comment.
- The Archive Node refers to the archiving responsibility of the device:
 - If Auto is selected, the MR system will automatically check to see if the device is a Storage Commitment Provider.
 - If Yes is selected, the device will be responsible for archiving images. When the device has received and saved the images, a notification message will be displayed on the scanner console and the Archive status for the exam will be “A” for archived.
 - If No is selected, the device will not be responsible for archiving.
 - **NOTE**
The device must be a Storage Commitment Provider in order for remote archive node to function.
- Access to the local host refers to the device's ability to access the MR system. Select Yes if you want the device to be able to send to and/or query the MR system.
- The Custom search feature enables the Custom search dialog box to be automatically displayed when you select receive from the remote browser. If Yes is selected, the feature is enabled. If No is selected, the Custom search dialog box will not automatically be displayed. You can, however, get to the search feature once the remote browser is displayed, by simply selecting Search, on the remote browser.

11.17.5 Network compatibility

The GE HealthCare MR system image formats are DICOM. This image format may only be transferred between systems using a DICOM network protocol.

11.18 CD/DVD/USB

The CD/DVD/USB allows you to store DICOM images and a DICOM Viewer to a CD-R or DVD-R or USB media that can be played back on a PC or laptop. Images stored on a CD-R, DVD-R or USB can be restored to the AW system or MR EXCITE HD system.

The CD/DVD interchange system is not considered part of the archive system but rather an ancillary storage device designed for sending select images or an entire exam to referring physicians. When an exam/series/image is burned to a CD-R, DVD-R or USB, the system does not consider the exam/series/or image as archived. Therefore, the archive column on the Patient List does not change its state based on CD/DVD/USB interchange activity.

11.18.1 Compatibility

The table below lists the compatible CD and DVD manufacturers.

- To export data, the CD or DVD must be R+W.
- For all other archive activities, the DVD must be DVD-R.

Table 11-16 MR-compatible CDs and DVDs

CD-R	DVD-R
Verbatim 1X-52X DataLife Plus CD-R	Verbatim 16X DataLife Plus DVD-R
Verbatim 4X-52X Media-Disc CD-R	Verbatim Media-Disc DVD-R

- To set the anonymization level, see 13.10.4 Anonymize Patient settings on page 2385.

12.3 AutoBind

About this task

Use these steps to bind axial images from multiple stations into a single sagittal or coronal view.

An intensity, normalization correction algorithm is applied station-by-station, slice-by-slice to the axial bound datasets. This algorithm is applied to minimize a band-like appearance that results from different intensity levels between stations. The images from the different stations may have different pixel sizes. Stations with different pixel sizes are bound together by upscaling the slices, that have lower resolutions than the slice with the maximum resolution.

Do not use the bound images as a single source for any diagnosis.

Step-by-step instructions


- From the header area of the screen, click the  to open the Image Management work area.
- From the Patient List, select AutoBind appropriate axial series from an exam that is whole body with multiple stations.
 - Screen saved images are not allowed.
 - Axial MR series with variable L-R and A-P FOV centers are supported.
 - Variable FOV of up to 50% are supported.
 - Variable matrix sizes between stations are supported.
 - A four slice axial overlap is required for good intensity correction between stations. Fewer overlap slices may result in a band-like appearance between stations.

Figure 12-6 Example of an AutoBind appropriate series

Exam	Station name	Name	Date	Description	Modality	Archived	Transferred
772	Bay7B	0119, Whole Body	Oct 21 16	Whole Body	MR		
383	Bay19	0086, Autobind	Oct 20 16	Autobind	MR		
208	dvmrbay	0005, DW-WBL, Vol3	Oct 06 16	WBL	MR		
176	dvmrbay	0001, WBIAutobind	Sep 23 16		MR		

Series	Type	Images	Description	Modality	Manufacturer	Archived	Transferred to
6	PROSP	72	1 DWI	MR	GE MEDICAL SYS.		
7	DIXON	72	WATER: 1 LAVA-Flex	MR	GE MEDICAL SYS.		
8	PROSP	128	CeJ C Spine+Neck 36	MR	GE MEDICAL SYS.		
9	PROSP	72	2 DWI bilat shim	MR	GE MEDICAL SYS.		
12	DIXON	72	WATER: 2 LAVA-Flex	MR	GE MEDICAL SYS.		
13	PROSP	128	Cal Chest 36 AA	MR	GE MEDICAL SYS.		
14	PROSP	72	3 DWI	MR	GE MEDICAL SYS.		
15	DIXON	72	WATER: 3 LAVA-Flex	MR	GE MEDICAL SYS.		
16	PROSP	156	Cal Body 24 AA2	MR	GE MEDICAL SYS.		
17	PROSP	72	4 DWI	MR	GE MEDICAL SYS.		
18	DIXON	72	WATER: 4 LAVA-Flex	MR	GE MEDICAL SYS.		

Today's diagnostic imaging technology provides an amazing breadth of information in a wide variety of ways. However, analyzing it, compiling what's relevant, and reporting your findings can put a serious crimp in your workflow, causing your operating costs to soar, and possibly even affecting your diagnostic confidence. AW VolumeShare 7 offers an elegant solution to your workflow challenges.

Overview

AW VolumeShare 7 is a multi-modality image review, comparison, and processing workstation with simplicity and power at its core. This workstation features 64-bit technology that allows processing of today's increased slice count studies. With significant improvements in workflow processes and clinical features that may help increase your productivity, AW VolumeShare 7 can elevate your practice to a whole new level.

Highlights

- Contemporary, simplified 2D Viewer user interface is harmonized with new Volume Viewer application. new hardware and operating system with expanded memory.
- Supports new applications in CT, MR and Vascular specialties.
- The Postfetch tool retrieves historical exams for comparison.
- Preference Sharing enables exchange of settings and protocols among users.
- The Key Images feature allows you to flag images and create key objects.
- Usage monitoring dashboard tracks and displays exportable analytics per application.
- Supports Chinese, Japanese, Korean and Russian languages for certain applications
- Supports enhanced security including DICOM TLS

Visit us:

www.gehealthcare.com/aw/platforms/aw-workstation/

CONFIDENTIAL

SECRET DE AFACERI

Features

- Simple drag and drop action for networking, media interchange, and filming.
 - Multi-tasking capabilities.
 - Active directory integration allows enterprise level user authentication. Directories supported include Microsoft® Active Directory® (MSAD) and other LDAP authentication services
 - Patient list to manage images from local workstation or remote systems.
 - Postfetch feature, designed to optimize your reading workflow by automatically gathering a patient's prior exams according to your pre-set criteria.
 - Support for Key Images.
 - End Review automates routine filming and networking tasks with just one click.
 - Search Advantage for fast and easy search of a patient's exam history on PACS or any other DICOM®-compliant device.
 - Enhanced Quick Filters of the Patient List filters studies by Modality, Date, End Review status or Exam Description.
 - DICOM CD/DVD/USB Creation Tool.
 - High-performance database management system.
- 2D Viewer for image display, manipulation, annotation, review.
 - Integrated Filmer with enhanced flexibility to perform all filming and data exporting tasks.
 - Access to a wide variety of applications for greater diagnostic flexibility.

Industry Standards

AW VolumeShare 7 complies with a wide variety of industry standards to facilitate adoption of features and performance improvements as the computing and medical imaging industry evolves.

DICOM Conformance Standards:

- DICOM 3.0 Storage Service Class for RT, CT, MR, CR, X-ray (Angio and R&F), Digital X-ray (DX), MG, NM, PET, U/S, Secondary Capture, Secondary Capture Color DICOM Image Objects. (Service Class User (SCU) for image send and Service Class Provider (SCP) for image receive).
- DICOM 3.0 Query/Retrieve Service Class (SCU and SCP).
- DICOM 3.0 Storage Commitment Service Class (SCU).
- DICOM Print (Color and B&W).
- DICOM Media Interchange (CD-R, DVD+R(W)).

Filming Protocols

- DICOM Print (Color and B&W) and Adobe® Postscript (Color and B&W) for supported printers.

Indications for Use

AW VolumeShare 7 is a review workstation, which allows easy selection, review, processing and filming of multi-modality DICOM images from a variety of diagnostic imaging systems. When interpreted by a trained physician, filmed or displayed images on the AW monitor may be used as a basis for diagnosis, except in the case of mammography images.

Regulatory Compliance

Rx Only

This product complies with European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC.

This product or its features may not be available in all regions. Please contact your sales associate.

CONFIDENTIAL

**SECRET DE
AFACERI**

AW VolumeShare 7

The AW VolumeShare 7 workstation gives you powerful capabilities aimed at increasing your diagnostic confidence, reducing your operating costs, and improving your departmental productivity.

Available Configurations:

VolumeShare is a multi-modality advanced visualization workflow solution that helps to enhance diagnostic confidence and productivity. The AW VolumeShare 7 workstation features software that is optimized for 64-bit technology and multicore processor hardware to provide leading edge performance. It includes AW's premier 3D image analysis package, Volume Viewer, with a suite of volumetric visualization and analysis tools for CT, MR, 3D X-ray and PET. In addition, this package includes the multi-modality 2D Viewer with dedicated functions for review of CT, MR, X-Ray, Angio, DX, U/S and PET images. It also includes the Filmer, a multimedia export tool for creating electronic films, filming and exporting to an internal web server, CD, or DVD.

AW VolumeShare 7 is also available for purchase without Volume Viewer for users who require only 2D image applications. Volume Viewer dependent features will not be supported on this configuration.

Modalities

Supported DICOM modalities include:

- Computed tomography (CT),
- Magnetic resonance imaging (MR),
- Radiofluoroscopy (RF),
- Xray angiography (XA),
- Computed radiography (CR),
- Digital radiography (DX),
- Mammography (MG),
- Nuclear medicine (NM),
- Positron emission tomography (PET),
- Ultrasound (US)
- Structured Reporting (SR),
- Key Objects (KO)

Volume Viewer Application

- VolumeShare 7's 3D image analysis package with a broad range of multimodality image analysis tools including multiplanar reformat, curvilinear reformatting, volume rendering and lumen navigation.

- For detailed description of Volume Viewer features, please refer to the Volume Viewer Product Data Sheet.

License Management

- Application Usage Monitor tracks how often a particular license was used, when, and by whom. You can compile and view usage reports for a specific period of time. This helps you understand how your systems contribute to departmental productivity and to plan for future license requirements based on actual usage data.
- For floating licenses, a secondary license server can be assigned for redundancy. If the primary license server becomes unavailable, the secondary can be contacted so that work can continue.

User interface

- Displays on one or two 1280x1024 monitors with scroll or optical mouse.
- Simple drag-and-drop mechanism for networking, media interchange, filming and data export options.
- Automatic help message display briefly describes function currently under the mouse pointer.
- Shortcut keys and programmable window/level function keys for accelerated control.
- Multi-tasking capabilities that allow Patient List, 2D Viewer, Filmer and Volume Viewer to run simultaneously with Fast Switch capabilities between applications.

Security

- Full support of Enterprise Authentication with Microsoft® Active Directory® and Novell® eDirectory™. Users can access their AW Workstation with existing Single Sign On credentials. This may streamline HIPAA compliance.
- If Enterprise Authentication is not used, an administrator-defined UNIX login and password is required per workstation. A lock screen feature enables you to lock the workstation display when leaving it momentarily, preventing unauthorized access to patient data.
- Main actions performed on the AW (save, film, network, export) are flagged and saved to an audit trail which includes username, ID, exam ID, date, and time.
- Multiple users can have different local UNIX logins, use the same AW workstation, and yet keep their user preferences unique.

GE imagination at work

AW VolumeShare 7 -1

- AW VolumeShare 7 is secured by Product Network Filters (PNF), which keeps unauthorized users from accessing open ports over the network. You can configure this feature to provide only specific remote devices access to AW ports.
- While not recommended, AW VolumeShare 7 does support ClamAV[®] anti-virus protection. This program is not activated by default in order to guarantee optimal performance of the AW workstation. If needed, your GE Healthcare field engineer can activate the program. However, it is the customer's responsibility to provide a safe internet connection from the workstation and to keep virus definitions current.
- AW VolumeShare 7 provides de-identification and encryption capabilities to limit privacy risks to sensitive information. The patient data exported during clinical workflow may be encrypted in flight by using either the SFTP or DICOM TLS protocol depending on the required dataflow.

Workflow management

AW VolumeShare 7 provides features designed to optimize your reading workflow:

- The Patient List provides tools that enable, you to sort and filter imaging studies stored on the local workstation or on remote systems.
- Postfetch allows you to retrieve a patient's prior related DICOM exams from a remote DICOM host. Retrieval can be triggered by the arrival of a new patient study on the workstation from the network or supported media.
- Support for IHE Key Image Notes profile, allowing you to mark certain images in the 2D or 3D viewer as key images. Key images are displayed as a separate series in the exam list and accessed or sent to an IHE KIN-supporting PACS. Specific tags and Key Image Notes may be created using the 2D Viewer application.
- Support for external USB disks to serve as DICOM storage media. The USB disk can be accessed for reading and writing via the media button, the same way you would access a CD or DVD. The maximum number of images that can be stored on a USB device is restricted to 300,000 across all exams. USB media must be formatted using the FAT32 file system. See the AW VolumeShare 7 Basic Display/Viewer/Filmer User Manual for more details on USB storage.
- Preference sharing among users. Sharable preferences include: Filmer layouts, Volume Viewer custom protocols, and viewer preferences. A standard USB flash drive may also be used to share preferences between workstations.

- One-Touch protocols enable you to define an application or preset protocol to launch automatically based on DICOM elements.
- DICOM Query/Retrieve Storage Class User (SCU) and Storage Class Provider (SCP) provide seamless network integration.
- DICOM Storage Commitment SCU lets you know when exams have been archived successfully on DICOM devices, such as PACS, that support Storage Commitment SCP.

End Review

End Review automates the routine tasks required at the end of reviewing each exam. The "End Review" flag in the patient list allows you to mark exams as "Done" after post-processing has been completed. End Review allows you to automatically perform one or several of these actions:

- Print pages prepared in the Filmer to the default printer (DICOM or Postscript filmers supported), with the option to automatically clear the Filmer after printing.
- Save electronic films in the workstation's attached DICOM database.
- Automatically push the entire exam or only the series created on the AW to one or more remote hosts.

Remote Network Host Management

- AW VolumeShare 7 Patient List provides a control panel to select a remote host or destination for network transfer. Icons signify different device types (acquisition, post processing, PACS, etc.) that are accessible to the workstation via the network.
- You can initiate a network transfer simply by dragging an exam, patient, series, or image(s) to the remote host icon. You can activate a remote browser to display detailed exam information from the remote host, including filtering information if supported.

Search Advantage

- This advanced search engine lets you search a patient's exam history on PACS or other DICOM device with a few mouse clicks.

Worklists and Filters

- You can filter the patient list by modality, date, end review status, or exam description. Most filters are available on the remote host patient list as well.
- You can filter the patient list further by choosing one or any of these parameters:
 - Modality.
 - Patient name or patient ID.
 - Exam location (hospital name) and exam description.
 - Series description

CONFIDENTIAL

SECRET DE
AFACERI

- Date and time of day: today, or specified date or date ranges with a specified time or time ranges.
- Radiologist's or referring physician's name.
- Accession number.
- Quick Sort (ascending or descending) and Quick Access (entry field) are available for any of the fields displayed at exam level (patient name, patient ID, exam location and description, radiologist's and referring physician's name, date and time).

Queues

- AW VolumeShare 7 manages three queues: network transfer; media interchange, filming. A menu on the patient list gives you easy access to queue status, and an animated icon gives you continuous network queue status information.

DICOM CD/DVD/USB creation tool

- The DICOM media creation tool offers you the flexibility to save and retrieve from supported CD, DVD, or USB storage devices.
- Add selected exam, series, or images by dragging and dropping to the pre-mastering window.
- CD and DVD usage percent is shown before you launch the writing process.
- You can optionally down sample certain 1024²x512² X-ray angiographic images during CD creation.
- CD/DVD composer lets you manage large amounts of data. If data size is larger than medium capacity, you are prompted to supply additional media storage space.
- You can specify the number of copies created during a CD/DVD save session.
- You can embed a DICOM viewer on DICOM CD/DVD so data can be reviewed on a PC running a Windows operating system.
- You get lossless JPEG reading of any DICOM CD/DVD media. Lossless writing is available for X-ray exams.
- CD/DVD drives operate at 16x write speeds and read speeds up to 48x depending on choice of media.
- For limitations on USB storage, refer to the AW VolumeShare 7 Basic Display/Viewer/Filmer User Manual.

Database management

- High-performance database management system classifies data according to patient folder description of the DICOM standard: Patient/Exam; Study/Series; Sequence/Images.
- Automatically deletes images on first-in first-out (FIFO) basis, which can be toggled on/off.
- Lock Exam feature protects specified exams from deletion.
- Permanent display of available disk space facilitates disk management.

- Study anonymization tool modifies DICOM elements to remove protected health information.

2D Viewer

The 2D Viewer is an application used to display, manipulate, annotate, and review 2D images by a trained physician for diagnostic interpretations. This tool has been designed to provide greater user efficiency, providing a set of intuitive, easy-to-use interface features.

Display customization

Display customization allows the user to manage layouts to display data. Key benefits include:

- Allows modification of number of exams or series displayed to facilitate either single or multi-exam reviews.
- Standard layouts give you flexibility to tailor the image display from 1x1 to 8x8. The current layout is retained if the images are from the same modality and orientation.
- Annotation levels allow selection of the image information fields to display.
- Toolbar customization puts you in control of which buttons are displayed to meet your needs.

Study Navigation

- The navigator lets you assign a series to a view on the fly.
- Cine mode also provides temporal, spatial, or manual playback loops.
- With two exams or series loaded, Cine mode supports a side-by-side display format with synchronized playback loops for more efficient comparisons.

Image Review

- The initial image window and level setting is based on the DICOM header. Once displayed, several methods of adjusting image window and level are provided.
- Provides routine image manipulation features:
 - Flip/Rotate, Zoom, Pan, Magnifying Glass
 - Inverse Video: inverts grayscale color map.
 - Display normal: lets you return the image to its default viewing parameters.

You can access the following features using a single mouse click directly on the image. Having these direct manipulation tools easily at hand gives you fast access with less distraction from your review task by eliminating the need to return to the graphical user interface controls: 2D distance, Angle, Report cursor, Box ROI, elliptical ROI and free-hand ROI

SECRET DE
AFACERI

CONFIDENTIAL

- The program continuously updates statistics on the fly.
- The Annotation feature lets you highlight areas of interest by adding text and line/arrow pointers to any image feature.
 - The Copy/Paste/Erase feature lets you copy/paste/remove any text or graphic placed on an image.
 - Cross-reference indicates the position of the current slice over the scout or localizer or any other non-parallel series.
 - Save lets you store a copy of the image as it appears on the screen for future review.
 - The Key Images feature allows you to flag images and create key objects.

Batch Filming

The Print Series feature lets you automatically batch film an entire series with a single keystroke.

Advanced X-ray Analysis

- You can apply shutters to X-ray images to focus on specific areas within the image. Shutters are black, opaque overlays with elliptical or rectangular cutouts that you can size as desired. Once a shutter is applied, you can move the underlying image to shift the viewable area. Shutters can also be applied to CT and MR images with the Image Matte feature.
- You can apply several levels of edge enhancement filters to an image.
- **Image Subtraction:** A graphical user interface lets you select a mask and subtraction of all images in an associated sequence. Subtracted images may be saved as a sequence with one mouse click.
- **Landscape:** With this tool, you can introduce a percentage of the mask image into the subtracted image for anatomical reference. A graphical user interface lets you specify the percentage of mask to be reintroduced.
- **Pixel Shift:** Auto pixel shift optimizes a mask to image registration. You can also move a mask manually to optimize subtraction quality in a specific region of interest.
- **Split Pixel Shift:** You can split the screen horizontally or vertically for pixel shift.
- **Maximum/Minimum Opacifications:** You can integrate selected images to provide a resulting Max. Op./Min. Op. image.

Filmer

The integrated Filmer enhances the efficiency of the review station and gives you greater filming and exporting flexibility.

This feature supports two modes: the Mini Filmer mode and Full-Screen mode, which provides the ability for full customization and film layout template creation. Three key mechanisms give you exceptional flexibility:

1. Free Format Filming
2. DICOM Structured Reporting (SR)
3. Data Export (HTML/PDF and JPEG/PNG/MPEG/AVI/QTVR)

With the Filmer you can easily extract significant images from any AW application (2D Viewer, Volume Viewer, READY View, etc.).

- You can film images individually by dragging and dropping to the on-screen Filmer, or by using the F1 keystroke.
- Multiple Image formatting lets you film multiple images in a single page frame with the F2 keystroke.
- With Film MID you can send multiple images to a single filmer frame with the F3 keystroke.
- Batch Filming is supported from applications that provide that capability (Print Series in the Viewer, Batch Film Protocols in Volume Viewer).

Communication between applications and the Filmer is accomplished in the mini Filmer mode, which provides a minimized footprint. The mini Filmer mode provides the following:

- Compatibility with Batch Filming from Volume Viewer (may be optional in your configuration)
- Store and position images transferred from an application.
- Type of export (film, media, database).
- Rapid switching between Full-Screen Filmer and application.

The resultant electronic films become a quick summary of the patient study and radiological interpretation to be reviewed by clinicians and physicians. They can also be saved as independent files (DICOM SR and Secondary Captures) for teaching purposes.

The Filmer works on the WYSIWYG (What You See Is What You Get) approach, and the electronic film can contain one or several pages with specific layouts for each page or all pages.

A flexible Edit Mode provides the ability to easily add, manipulate, format, or delete images from the film. Images can contain text and graphics from measurements and your annotations, and may be window/leveled, magnified, flipped, rotated, or cine. You can add additional annotation to the image in edit mode.

The Preview Mode displays the film as it will be printed or exported, taking into account the layout applied to each page and the compression level specified for non-DICOM exporting (JPEG/PNG and MPEG/AVI/QTVR).

Printing

AW VolumeShare 7 includes Network DICOM Print (B&W and color). For supported printers, AW VolumeShare 7 includes Network PostScript capability (B&W and color). For PostScript and

GE imagination at work

SECRET DE
AFACERI

CONFIDENTIAL

DICOM printers, any printing format created in the Filmer is supported (e.g. non-square matrix formats for rectangular images such as CT run-offs).

- AW VolumeShare 7 postscript printing has been tested on the following devices: Codonics 1660M, 1660MD or Horizon, Lexmark Optra 1650N, 1855N, SC1275N, C710N, C720N, T612 or T614, Seiko 1720D, Kodak DMI3600, Quantum GL2101HD, declared with film/thick paper or with plain paper, Tally T8106, HP LaserJet and Xerox Phaser
- Digital cameras and analog cameras using a 3M-952 protocol (including DASM interface) are not supported on AW VolumeShare 7.

Exporting

Data export is integrated in the Filmer, thus providing all image processing tools needed for multimedia image export. The Filmer exports any electronic film to DICOM SR, PDF/HTML, JPEG, PNG, MPEG, AVI or QTVR format. AW VolumeShare 7 enhanced security of the export functionality with the use of SFTP connectivity. Non-DICOM data can also be saved on a multi-session CD/DVD. The Data Export capability is intended only for publishing and communication, not for diagnostic purposes. Its simplicity is reflected in the different export mechanisms available:

- CD/DVD removable media
- USB flash drive
- Network HTTP and FTP protocols

Cardiac Review and Export

Processing and review of CT, MR and PET cardiac exams with manual oblique reformatted protocols can be exported as a multi-phase Cine movie that allows referring physicians to review exams in a dynamic mode.

System Components

Workstation Configuration

- HP Z4 G4 Workstation
- HELIOS 7 operating system
- Intel® Xeon® W-2245 Eight physical core/ sixteen logical cores 3.9 GHz CPU
- 64GB (4x16GB) DDR4 2933 MHz or higher ECC Registered DIMM
- NVIDIA Quadro NVS P620, 2GB Graphics card
- 1 x 512GB M.2 NVMe Drive for OS and Applications
- 1 x 1TB M.2 NVMe Drive for image cache.
- Cache, which is subject to overhead can store approximately:
 - 6,400,000 256*256*16 bits uncompressed images OR
 - 1,600,000 512*512*16 bits uncompressed images OR
 - 400,000 1024*1024*16 bits uncompressed images OR
 - 100,000 2048*2048*16 bits uncompressed images.
- Internal DVD Writer drive for read/write of DICOM CD/DVD media, read/write of Data Export CD/DVD data and service use (DVD Install)

- Integrated dual Ethernet 10/100/1000 Mbit/s ports.
- Additional single Ethernet 10/100/1000 Mbit/s port
- 1 USB QWERTY (or regional) Keyboard and Mouse

Footprint

- Height 38.60 cm (15.2 in.)
- Width 16.89 cm (6.65 in.)
- Depth 44.47 cm (17.5 in.)
- Approximate Weight 17 kg (38 lbs)

Operating Environment

- Temperature: +5°C to +35°C (40° to 95° F)
- Humidity: 10% to 85% (relative non-condensing)
- Altitude: 0 to 5000m (16,404 ft.)
- Acoustics: LWAd less than 4.3 Bels
- Shock: 40 G peak, half-sine, 2-3 ms

Non-Operating Environment

- Temperature: -40°C to +60°C (-40° to 140° F)
- Humidity: 10% to 90% (relative, non-condensing)
- Altitude: 0 to 12192m (37,000 ft.)

Monitors

- (2) 19" color Flat Panel LCD monitors
- DICOM Part 14 factory calibrated
- Native resolution 1280 x 1024 (5:4 aspect ratio)
- Approximate Weight: 5.3 Kg (11.68 lbs)
- AC 100 - 240 V : 50 / 60 Hz

Image Networking

- Standard 10/100/1000 Base-T Ethernet for DICOM
- 1000 Base-T dedicated network for optimal Direct Connect performance
- Protocols supported:
 - DICOM 3.0 Storage SCU/SCP and Query/Retrieve SCU/SCP
 - RSVP
 - TCP/IP network layer
- SdCNet is no longer supported
- AW VolumeShare 7 does not support the AdvantageNET network protocol.
- AW VolumeShare 7 does not support DICOM images from GE Healthcare Signa™ version 5.x 1.5T MR systems
- AW VolumeShare 7 software is supported on the previous generation AW VolumeShare 5 HP Z820 and Z440 workstations.

CONFIDENTIAL

SECRET DE
AFACERI

About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems. Our "healthymagination" vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality around the world. Headquartered in the United Kingdom, GE Healthcare is a unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employees are committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com.

GE Medical Systems SCS
283, rue de la Minière
78530 Buc, FRANCE

©2014-2022 General Electric Company.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. Timing and availability remain at GE's discretion and are subject to change and applicable regulatory approvals. Contact your GE representative for the most current information.

GE, the GE Monogram, Centricity, SIGNA and imagination at work are trademarks of General Electric Company.

Adobe is either a registered trademark or trademark of Adobe Systems Incorporated in the United States and/or other countries.

ClamAV is a registered trademark of Cisco Systems.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

Intel and Xeon, are either registered trademarks or trademarks of Intel Corporation in the United States and/or other countries.

Microsoft, Active Directory, and Windows, are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Novell and eDirectory are registered trademarks of Novell, Inc. in the United States and other countries.

Oracle and Java are registered trademarks of Oracle and/or its affiliates.

All other product names and logos are trademarks or registered trademarks of their respective companies.

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SECRET DE
AFACERI

Volume Viewer

GE Healthcare's premium 3D Advanced Visualization and Image Processing platform brings you powerful capabilities to help you visualize and interpret your multimodality imaging data with confidence and ease.

AutoLaunch

Preprocessing enabler

Imaging modalities provide a wealth of diagnostic information, but also present a challenge to your Radiology workflow as image volumes continue to increase while your time for reading and reporting decreases.

Volume Viewer provides you with a rich 3D image processing toolset aimed at creating and displaying the views you need with little user input and streamlining interpretation and reporting by providing the visualization tools you need with minimum clicks.

With an intuitive, modern user interface, Volume Viewer helps a wide diversity of users to learn and master the expanding portfolio of sophisticated tools and applications that it provides. This customizable user interface maximizes the real estate allocated for displaying images and provides access to the users' favorite tools directly from the image viewport. In addition, a customizable toolbar provides one click access to commonly used tools to facilitate a productive review.

Overview

Volume Viewer provides excellent 3D visualization and processing capabilities for reading and comparing CT, MR, 3D X-ray, PET, PET/MR and PET/CT datasets. Volume Viewer also features a broad portfolio of high-performance analysis tools, automating routine tasks and helping to make 3D image processing a stress-free component of your routine workflow.

Volume Viewer is available on VolumeShare 7, a multi-modality advanced visualization workflow solution that helps to enhance diagnostic precision and productivity.

Highlights

- **Customizable protocol-driven workflows** supporting CT, MR, PET and Innova 3D XA modalities, designed to help enhance your clinical review productivity
- **Compare** multiple exams from **different modalities, at loading² and on the fly²**
- **Adaptive layouts** that take advantage of the dual monitor capabilities in landscape and portrait orientations
- **Easy point-and-click segmentation tools** for segmenting structures of interest
- **CT Dual Energy** and GSI⁴ support with dedicated workflows and tools
- **Advanced 4D PET** support with Q.Freeze⁶ tool and 4D PET Dynamic workflow
- **Interactive Summary Table** collecting measurements as they are deposited on the images.
- **Rich set of export tools** including export for 3D printing³.

<https://www.gehealthcare.com/en/products/advanced-visualization/all-applications/volume-viewer>

CONFIDENTIAL

SECRET DE
AFACERI

General Features

- **Unified user interface**, with one click access to your favorite tools and more space for displaying images.
- **Customizable protocol-driven workflows** designed to help enhance your clinical review productivity for a wide range of care areas, with 3D and 4D real time interaction.
- **"Smart layout"** protocol, which adapts the screen layout to the type of loaded exams.
- **High resolution, real-time rendering modes:** MIP/Min IP/Average; standard volume rendering and enhanced volume rendering with Volume Illumination¹ on the fly; Navigator view for endoluminal fly-through; lumen and curved views.
- **Compare mode** for reviewing multiple exams and modalities, simplified with **Dynamic Load to load and register²** new series into the current review session.
- Rich set of **2D/3D ROI** tools, including **AutoContour** for 3D semi-automatic contouring of structures of interest (CT, MR, PET), and **ROI color coded** for tissue classification based on voxel values.
- **Advanced tools**, to take advantage of full 3D capabilities: **AutoSelect**, for easy point-and-click segmentation; One or two clicks **Quick Vessel Trace** to analyze all vessels, in curved reformat, lumen, or MPR view; Advanced 3D processing tools.
- **Summary Table** extends reporting capabilities, while collecting measurements as they are deposited on the images.

- **Multiple Export capabilities**, in different DICOM and non-DICOM formats, including **Key Image Notes**.
- **3D Suite³:** STL, OBJ, VRML and 3MF file formats to export single or multi mesh objects, for 3D printing. 3D Suite is not intended for clinical use.
- **Save state**, which allows to save and restore to pause and resume the review at a later time or to facilitate collaboration workflows.

Modality Features

- **Dual energy and GSI⁴ GE Healthcare CT images** supported with dedicated review protocols.
- **Bone VCAR⁵** combines spine labelling on the fly for all CT review and a dedicated layout to improve spine review.
- **MR review support:** Recognition of PSD name, and specific MR parameters at image loading, which simplifies protocol customization.
- Dedicated review workflows for **PET/CT** and **PET/MR** images, including **PET SUV measurements**.
- **Q.Freeze⁶** combines the quantitative benefits of **4D PET gated** imaging into an image that provides both frozen patient motion and reduced image noise.
- Dedicated protocol-driven workflows for **Innova 3D XA images⁹**.

Platform features

- **"Smart Compression"** technology for fast and responsive client performance on AW Server.
- Fast automated **Pre-processing⁷**

- Multi-sessions switch on AW workstation with **AutoLaunch⁸**.
- Support of **Chinese, Japanese, Korean and Russian** languages.

System Requirements

- AW Server 3.2 and above, with recommended monitor resolution up to dual 2MP or single 3MP.
- AW 4.7 Workstation and above (some performances limitations can occur with Z800)
- CentricityTM Universal Viewer

Intended use / Indications for use

Volume Viewer is a medical diagnostic software that allows the processing, review, analysis and communication of 3D reconstructed images and their relationship to originally acquired images from CT, MR, X-Ray Angiography and PET Scanning devices. The combination of acquired images, reconstructed images, annotations and measurements performed by the clinician are intended to provide to the referring physician clinically relevant information for diagnosis, surgery and treatment planning.

Regulatory compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC. This product or its feature may not be available in some other countries or regions. Please contact your sale associate.

Rx Only

¹Requires Volume Illumination license
Volume Illumination may not be available in all countries or regions.
Please contact your sale associate.

²Requires Integrated Registration license

³Requires 3D Suite license

⁴Requires Volume Viewer GSI license

⁵Requires Bone VCAR license

⁶Requires Q.Freeze license

⁷Preprocessing is compatible with the following purchasable options: AutoBoneTM Xpress, CardIQ Xpress Reveal, CardIQ Xpress Function, CardIQ Xpress Process, Advanced CTC Pro 3D EC, Colon VCAR, and CT Perfusion 4D and FastStroke. ColonVCAR is not for sale in the US.

⁸Requires AutoLaunch license for AW workstation

⁹Requires Innova Volume Viewer license

SECRET DE
AFACERI

CONFIDENTIAL

Volume Viewer


General Features

Unified user interface

- The user interface provides one click access to favorite tools and more space for clinical images.
- The page layout is organized to facilitate navigation through the review steps. Guided protocols walk the unfamiliar user through a clear workflow of creating and manipulating volumetric data, helping to reduce the learning curve by providing a consistent workflow.
- The customizable toolbar can be docked at the top, bottom, left, or right, by a simple drag & drop. This preferred location is saved for future use for each user.
- An extensive set of advanced tools are grouped into categories on the toolbar to facilitate quick retrieval when needed.
- Each category and tool are easily configurable by a simple drag & drop in the Toolbar so that they are only one click away for your future reviews
- To reduce long mouse trips, up to 6 favorite tools can be accessible from each image via a right mouse menu.
- Customization of the toolbar and tool behaviors with a user's preferences is provided in a single intuitive Preferences menu.
- Keyboard shortcuts and quick access tools for on-viewport controls

The following mouse modes are available to help you manipulate the images:

- Interactive VR adjustment lets you adjust VR opacity interactively with the mouse as Window Level/Window Width
- Direct Paging (continuous or contiguous)
- Free hand rotation on 3D and oblique views
- Percentage fusion, to easily change the transparency of objects on fused images. There are several fusion modes available for your use.

- Mag glass displays zoomed-up imagery around the cursor.

Protocol management and loading

- Customizable protocol-driven workflows for a wide range of clinical applications, including support for multiple review steps. You are guided through the review steps by a Navigation bar at the top of the User Interface.
 - Review protocols created by user with Volume Viewer can be saved as favorite protocols and accessed easily in routine.
 - Progressive Load enhances image loading performance. This feature has two modes:
 - Images are loaded in sequential mode (e.g., Reformat)
 - Images are loaded in interlaced mode (e.g., 3D/VR)
- In either case, image review can begin as soon as the first image is displayed.
- Navigation through series and exams can be performed without exiting the patient list.

Smart layout

- Volume Viewer General Review protocol presents an adaptive layout for single and multi-modality PACS-like reading of data. This protocol distributes and displays all loaded series of a single exam in an intelligent layout to maximize use of available monitor space. With multiple exams, the protocol automatically detects and distributes "like" series for optimal comparison.

GE imagination at work


CONFIDENTIAL

SECRET DE AFACERI

Advanced visualization capabilities

Volume Viewer offers various capabilities to display advanced rendering modes:

- **3D Volume Rendering / Volume Illumination¹ images**
- **MPR views with different slice thickness. The following rendering modes are available for all Thick Slab: MIP, Min IP, Volume Rendering, Average.**
- Navigator views which display interactive endoluminal views exportable as movies.
- Curvilinear reformatting allowing display of curved, lumen and cross section views for various structures (vessels, spine, etc).
- Merge multiple volume rendered models into a single view.

Volume Rendering

Volume Illumination

Predefined cut planes

Volume Viewer lets the user define cut planes to isolate specific structures in the VR model, in 16 configurations, including Left / Right / Inferior / Superior / Anterior / Posterior / Front. You can then display the voxels values on each plane.

Multiple Volume Rendering and Volume Illumination¹ objects

- Volume Viewer allows you to merge up to 8 Volume Rendered models from the same series into a single 3D view.
- You can adjust independently the threshold, colors, and transparency of each Volume Rendering or Volume Illumination


Auto-center fly through with smart cursor

- In fly-through studies (airways, colon, angiography navigation along the centerline of the structure of interest can be recorded step by step and can be exported as a movie.
- In addition to the regular 180° or less Navigator view, Volume Viewer enables a Fish-Eye View with any wide camera angle value from 180° to 360°. This provides a view of structures both in front of and behind the users' virtual location on the same image.

Lumen view

- Lumen View provides an unfolded 3D view around a user defined centerline. The lumen view can be interactively adjusted (rotation around the centerline, width, field of view).

3D and 4D capabilities

- Volume Viewer allows 3D cursor synchronization in any orientation on the fly.

GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

- Real time reformat in oblique planes is available for your exams, as well as simple and double oblique interactive modes.
- Quick access to cross reference display gives you the ability to display the current slice location or all slices, or other planes.
- Real time interaction in 4D mode, by using Cine tool on all planes and 3D views for CT, MR, and PET multiphase data.

Compare mode


- With Volume Viewer several volumes can be loaded either from the same exam or from different exams. This is useful for exam comparison (Compare Mode) or for Multiphase examinations (of the liver for example). The Compare Mode is also available for PET/CT examinations and takes advantage of a dual screen configuration.
- Multi-phase CT and MR data are listed as separate phases in the Series Selection screen, allowing selection of all phases or a subset of phase data for review. The Series Selection Panel displays the multi-phase data split into arterial and portal phases, for example.
- Dynamic Load, compatible with both Volume Viewer and Fusion protocols, lets you drag & drop 3D volumes from CT, MR, PET, and 3D XA modalities into a desired viewport. Together with the Integrated Registration² option, Dynamic Load allows new volumes to be registered and loaded on the fly. Save State series can be restored as a separate session. For dual monitor configurations, a Save State series can be displayed on the right monitor together with a current session on the left monitor. For single monitor configurations, a user may simply switch between the two sessions.
- The zoom and pan functions are propagated to all images displayed in the same orientation.

2D/3D ROI for quantitative measurements


- Volume Viewer contains a set of basic 2D/3D tools: distance, angle, report cursor, arrow tool, annotation, free hand ROI tool, that can be placed on fused images as well.
- Easy deposition and labelling of measurements.
- All 3D ROIs can be customized, and color coded to display statistics computed on different intensity range and help assessing tissue classification based on voxel values.
- All measurements are considered bookmarks and can be accessed by navigating through the list of measurements.

AutoContour tool

- AutoContour provides consistent tools for 3D contouring ROI, with the same workflow used for CT, MR, and PET images.
- For CT and MR images, the contouring tool will suggest contour boundaries after defining 1 point or a diameter on a structure of interest, which can be edited prior to validation.

- For PET images, the contouring tool has been adapted for standard SUV calculations.

Advanced 3D tools

Segmentation tools

- AutoSelect segmentation tool helps adding or removing any continuous structures of interest by simple point and click.
- You can also access advanced segmentation tools such as Scalpel, Threshold, Dilate, Erode, Subtraction, Intersection, Addition, Filter Floaters.

Quick Vessel Trace (CT and MR images)

- You can perform vessel analysis with one or two points clicks from any CT or MR Angiography images.
- Quick Vessel Trace automatically extracts the vessel after user clicks and launches vessel tracking for fast review in curved reformat, cross-sectional, best L-section, lumen, and MRP view.

GE imagination at work

CONFIDENTIAL

SECRET DE AFACERI

- Tracking can be extended proximally or distally for a full view of the vessel.

Summary table

- The Summary Table collects measurements and organizes them according to finding and exam date.
- The Summary Table allows interactive navigation through findings in the image viewports and can be docked or expanded.
- The Summary Table can be copied into the client clipboard and then pasted in email, reports, office software etc. (on AWS only)

Export capabilities

Volume Viewer contains multiple standard options for exporting the results of a review session:

- Save images to the database, as new DICOM series.
- Save the significant images as Key Image Note objects in the database. Key Image Notes and End Review allow you to flag images of interest as Key Image Notes (IHE profile) and push them to the archiving system when you exit the application.
- Save State object are used to save the work in a new series containing all the post processing data (3D Model, displays, measurements, annotations, etc), for future review.
- Save STL, OBJ, VRML or 3MF file from 3D Volume Rendering / Volume Illumination¹ viewport for 3D printing purposes through 3D Suite³

Additional export tools are available to create specific sequences of images to be exported:

- The Batch tool allows creating a sequence of reformatted images or a sequence of rotating 3D views.
- The Movie tool creates a comprehensive movie including different rotations, zooms, and pan of the image, which can be exported as DICOM series or mpeg file.
- Quick Export: Exports in a single click a full batch of contiguous images at the displayed thickness for 2D images, or a batch of rotations of a 3D View.
- Cardiac Review and Export: Processing and reviewing cardiac exams for CT, MR and PET with manual oblique


reformatted protocols can be exported into a multi-phase Cine movie that allows the referring physician to review the exam in a dynamic mode.

- Capture and send images from the viewport to your computer clipboard on the fly (Ctrl+C/Ctrl+V on AW Server only).

Modality Features

CT image quality


- CT filters to enhance or smooth CT images to help improve image quality on the fly.
- The lung filter helps enhance contours of images reconstructed in standard mode for excellent visualization of lung structures

Standard recon w/o filter

Standard recon with Lung filter

- A smoothing filter reduces noise while maintaining high image contrast.

No filter

Smooth filter

No filter

Smooth 3D+ filter

CT Dual energy protocols

- Three protocols are available to enable review of dual energy images acquired on compatible GE Healthcare scanners and enhancing pixels with specific ratio: **DE Calcium** ($HU_{80}/HU_{140} \geq 1.25$), **DE Uric Acid** ($HU_{80}/HU_{140} \leq 1.25$) and **DE Custom** (custom threshold on HU_{Low}/HU_{High}).

GE imagination at work

CONFIDENTIAL

SECRET DE AFACERI

- In addition, all Volume Viewer standard protocols are compatible with GSI⁴ GE Healthcare CT images.

Bone VCAR⁵


- Bone VCAR is available in Volume Viewer toolbar and can be launched on the fly
- Vertebrae are automatically recognized and labelled (deep learning algorithm)
- A dedicated protocol provides a specific layout of the spine including curved views

Review of MR studies

Volume Viewer offers support for review of MR studies, through dedicated protocols.

- The MR-specific parameters (PSD name, image weighting, scan plane, etc) are recognized at loading and are used to customize the layout display.
- Anatomy-based protocols for standard review (e.g. spine, brain stroke...), with dedicated and guided review steps, to help review standard MR exams.
- Dedicated MR Cardiac Viewer for single and multiple Cine Review. Enables comparison between cardiac series such as Time Course and Myocardial Delayed Enhancement (MDE) images.
- Support of MR Multi Echo, Multiphase, Diffusion series.
- Support of parametric series (e.g., ADC, MTT, and other parametric maps) created in READY View, allowing direct functional measurements.


- 2D and 3D ROI propagation to other phases / series / exams allowing you to easily correlate information from multiple sources.
- Improved MR annotation consistency.
- In combination with Integrated Registration², enables direct access to MR image registration (inter/intra exam), image fusion and Whole Body MR Review protocols
- Ability to automatically bind several MR axial series corresponding to sequential axial locations into a single series. This is useful for consolidating multi-stage acquisitions for displaying Whole Body MR scanning.

SUV measurements on PET images

- Volume Viewer supports Standardized Uptake Values (SUV) for image display and measurements. Several SUV scales are available like SUVbw, SUV lbm, SUVbsa, as well as SUV Peak.
- Window/level presets may be user-defined for PET images based on SUV values.
- The SUV values are available in all the basic 2D/3D ROI tools of Volume Viewer and in the AutoContour tool.

PET/CT and PET/MR dedicated protocols

- Volume Viewer contains a list of predefined protocols allowing the review of PET/CT and PET/MR data. The screen layouts contain fused views between the morphological and functional images, as well as 3DMIP rendering of PET images.
- Factory protocols are customizable for your own review.
- PET 4D protocol allows users to load and display gated and dynamic PET series, as well as summing or reframing these series.

GE imagination at work

CONFIDENTIAL

SECRET DE AFACERI

PET IQ improvements

- Enhanced 3D visualization algorithm to fully support Q.Clear PET images.
- Q Freeze⁶ option provides a respiratory motion-corrected Gated PET volume, statistically similar to conventional static PET with significantly reduced or eliminated blurring effects due to patient respiration. It provides also a corrected Gated PET series where counts are summed back to each respiratory phase.

Dedicated features for 3D XA images⁹

The optional Innova 3D XR application produces XA images (512x512 images or 256x256 images depending on your settings). This format opens the Volume Viewer world to Innova 3D XA images:

- Single exam: Innova Navigator, Curved Reformat, Multi Oblique, MPVR 3mm Average, Segment Structure
- Compare exam: Compare Axial, Compare Coronal, Compare Sagittal, Reformat
- Layouts and protocols customization
- The Innova options provide access to Workflow booster tools, easy multi-modality review "on-the-fly" dedicated Innova review protocols with dedicated Review Steps.

Smart compression

- The Smart Compression technology automatically displays images with full fidelity when the image is still, then uses the selected compression level for increased interaction speed during user interaction. This allows for diagnostic reads on full fidelity static images with responsive dynamic display even at low bandwidth. On-image visual indicators notify the user when compression is in effect.

Pre-processing⁷ features on AW Server

- Preprocessing is a Volume Viewer feature that automatically performs routine processing tasks and saves the results so they are ready for you. When new exams are transferred to the database, Volume Viewer recognizes user-defined keywords in the Series Description and launches the appropriate pre-processing. The results are stored in a Save State object with the original study so you can load it when you are ready for review.

AutoLaunch⁸ features on AW workstation

- AutoLaunch is a Volume Viewer feature that automatically preloads exams, so they are ready for your review.

- When you are processing a study and a new exam is transferred to the AW Workstation, this feature automatically launches the new exam with an appropriate Volume Viewer protocol in the background. When you are ready, one click in the AutoLaunch window instantly brings up the exam in the Volume Viewer application, eliminating waiting time and extra steps to load the exam into computer memory for processing.
- AutoLaunch is compatible with CT, MR and PET single volume protocols of Volume Viewer.
- When combined with optional applications⁸, this feature gives access to data already preprocessed, that you can review by just clicking on AutoLaunch window.

Summary of Operation

Volumetric models are loaded by selecting the exam or series. The user can select a protocol category from an anatomical selector or go directly to a Review Layout. In either case, images are loaded progressively in the background; this gives control to the user in just a few seconds after selecting the images. Selecting a Review Layout launches a volumetric display protocol with predefined layout preferences. Review Layouts may be customized and combined to suit the user's workflow. Selecting a protocol category unlocks a variety of visual protocols that include the layout, threshold, rendering mode and filming formats. Some of these protocols direct the user through the process providing capabilities to interactively view and manipulate the model, increasing productivity and consistency for all modalities.

References

¹Requires Volume Illumination license

Volume Illumination may not be available in some other countries or regions.

Please contact your sales associate

²Requires Integrated Registration license

³Requires 3D Suite license

⁴Requires Volume Viewer GSI license

⁵Requires Bone VCAR license

⁶Requires Q.Freeze license

⁷Preprocessing is compatible with the following purchasable options: AutoBone™ Xpress, CardIQ Xpress Reveal, CardIQ Xpress Function, CardIQ Xpress Process, Advanced CTC Pro 3D EC, Colon VCAR, and CT Perfusion 4D and FastStroke. ColonVCAR is not for sale in the US.

⁸Requires AutoLaunch license for AW workstation

⁹Requires Innova Volume Viewer license

CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

GE Medical Systems SCS
283 rue de la Miniere
78530 Buc France

DOC1460639 rev9

General Electric reserves the right to make changes in specifications and features, or discontinue the product or service described at any time, without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. Illustrations are provided for informational purposes, and your configuration may differ.

This information does not constitute legal, financial, coding, or regulatory advice in connection with your use of the product or service. Please consult your professional advisors for any such advice. Operation of GE Healthcare products should neither circumvent nor take precedence over required patient care, including human intervention of healthcare providers. GE Healthcare products and services do not code medical procedures. Accurate coding is the responsibility of the provider or billing professional.

GE, the GE Monogram, imagination at work, Centricity and AutoBone are trademarks of General Electric company.

General Electric Company, by and through its GE Healthcare division.

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
General Electric Company, prin intermediul diviziei sale GE Healthcare.

©2019 General Electric Company

[Scrieți aici]

GE imagination at work

CONFIDENTIAL

SECRET DE
AFACERI

Publications

5819086-1EN
Revision 1

User Guide

Part 1: Overview of Volume Viewer Dynamic Shuttle

CE 0459

© 2019 General Electric Company. All rights reserved.

CONFIDENTIAL

**SECRET DE
AFACERI**

Note: The contours are not displayed on slab thick views but only on views with the minimal slice thickness.

Statistics

Right click on the bookmark → Statistics

List of available statistics for PET ROI (deposited via <Insert> key or via AutoContour tool): Volume, Minimum, Maximum, Average, Standard Deviation, Relative Deviation, TLG (SUV only), Peak, Threshold

Note: A maximum of three statistics can be selected and displayed.

Note: The units of display of the VOI will always match the WL units displayed in the lower left of the viewport (VOI statistics in a fused viewport will be displayed in SUV, not HU).

Note: The TLG is the product between the volume of the tumor and its SUV average. TLG is only available for VOIs that are automatically thresholded.

Note: The negative pixel values from PET images are not taken into account in the calculation of VOI statistics.

Contextual options

Right clicking a ROI/VOI when it is active will enable you to:

- Show/Hide statistics in the ROI/VOI (will show/hide statistics for all ROI's/VOI's in all viewports),
- Lock/Unlock ROI/VOI position, size and comments,
- Delete ROI/VOI,
- Clone ROI/VOI.
- Show/Hide max value in the VOI,
- Send cursor to max value.

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
Brain View

Beyond conventional brain anatomical images.

Advanced Magnetic Resonance Imaging (MRI) methods such as Spectroscopy, Diffusion Weighted Imaging (DWI) and Diffusion Tensor Imaging (DTI), Time Course series, and Functional Imaging have significantly improved the ability to image conditions of the brain. Beyond conventional anatomic and structural data, these advanced techniques provide physiological information on metabolism and hemodynamics.

Highlights

- Guided workflows with intelligent display based on smart layout to help analyze MR data.
- Adapt your application to fit personal and institution requirements for more standardized analysis and improved productivity
- Provides additional clinical information through relative ROI measurements curves and color parametric images.
- Enables fusion of color parametric images with anatomical 2D or 3D images with simple "drag and drop" method.
- Provides an adaptive protocol for brain multi-parametric data processing.
- Enables MR to MR image registration to reduce patient motion effects.
- Accessible from PC, laptop, PACS/RIS workstation for streamlined workflow.

Overview

Integrated with the READY View platform, Brain View offers you advanced techniques designed for easy and confident analysis of information from a variety of MR brain-specific imaging data sets.

SECRET DE
AFACERI

CONFIDENTIAL

Visit us:

www.gehealthcare.com

Features

Brain View offers four advanced protocols:

- **Arterial Spin Labeling (ASL)** for contrast-free assessment of cerebral blood flow. ASL uses water in arterial blood as an endogenous contrast media to visualize tissue perfusion and help evaluate vascular-deficient or vascular-rich brain regions.
- **BrainStat** with Gamma Variate Fitting (GVF).
- **BrainStat** with deconvolution of arterial input function (AIF) and time series 1-Click motion correction. It generates hemodynamic information such as
 - relative Cerebral Blood Flow (rCBF)
 - relative Cerebral Blood Volume (rCBV)
 - Mean Transit Time (MTT & MTT-SVD)
 - Time To Peak (TTP)
 - Bolus Arrival Time (BAT)
 - Time to maximum value of the residue function (Tmax)
 - T1 Leakage & T2* Leakage indicator maps
- **FiberTrak** for tracking White Matter fibers in all orientations & restore them with Save State. FiberTrak provides detailed information on

tissue microstructure by reconstructing fiber pathways from diffusion tensor series.

- Brain View offers an adaptive protocol to process multi-parametric brain data. This factory protocol is adjustable and can be customized to fit your personal requirements in terms of
 - display & workflow (layouts, review steps)
 - parameters and settings
- Easy-to-use slide bars let you segment parametric images in real time.
- Display and export ROI statistics from the Summary table
- Export graph values as csv file.
- Save State let you save and restore the state of the processed images at any stage.
- Contextual help pages that give general assistance about the image processing algorithms.
- Save all generated parametric images in one click.

System Requirements

- AW Server 3.2 and above
- AW 4.7 Workstation and above
- Color Monitors
- Centricity™ Universal Viewer

Software Requirements


- Brain View requires READY View as prerequisite.

Product Description

Brain View, is an image analysis software option of READY View that allows the user to view and process Magnetic Resonance (MR) images of the brain.

Regulatory Compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC

© 2016 General Electric Company.
All rights reserved. Data subject to change.
GE, the GE Monogram, imagination at work and Centricity are trademarks of
General Electric Company.
DOC1915506
GE Medical Systems SCS
283 rue de la Miniere
78530 Buc France

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