

# ONESPRAY WIPES ALCOHOL BASED FAST ACTING DISINFECTANT WIPES FOR MEDICAL INSTRUMENTS



EN 13624  
EN 13727  
EN 14476  
EN 14348  
EN 16615

## USAGE

OneSpray Wipes is used for the disinfection of pre-cleaned, non-invasive, smooth and non-porous medical devices (Ultrasound probes, etc.) which can not be disinfected by the immersion method. It is ready to use.

The entire surface of the device is wiped with enough wipes.

It is based on a combination of ethanol, propanol and quaternary ammonium and does not fade or dye the material used. Has a nice fragrance.

## PROPERTIES

Ready to use product

Aldehyde and phenol free.

Does not do color fading or dyeing on the material used

Pleasant-smelling.



**%99.9**  
Effective

## MICROBIOLOGICAL EFFECT

Bactericidal, Fungicidal, and Virucidal effects.

## COMPOSITION

%30 Ethanol, %10 2-Propanol, %0.25 Didecylmethylpoly(oxethyl) Ammonium Propionate, Excipients, Essence, Deionized Water.

## WARNINGS AND STORAGE CONDITIONS

Keep out of reach of children.

Keep away from sunlight and high temperatures.

Store in tightly closed original packaging/container, dry and cool area.

Keep away from heat, sparks and fire.

Under normal conditions, shelf life is 2 years.

Store locked up.



Alcohol Based



Ready to use



Aldehyde and phenol free



Pleasant-smelling



Fast and effective disinfection

KAF G35 - 1	100 Piece
KAF G35 - 2	150 Piece
KAF G35 - 3	200 Piece

	Test Method	Microorganism	Effect Duration
Bactericidal	EN 13727	Staphylococcus aureus	30 sec
		Pseudomonas aeruginosa	
		Enterococcus hirae	
Fungicidal / Yeast-cidal	EN 13624	Candida albicans	5 min
		Aspergillus niger	
Tuberculocidal / Mycobactericidal	EN 14348	Mycobacterium avium	5 min
		Mycobacterium terrae	
Virucidal	EN 14476	Poliovirus Type 1	1 min
		Human Adenovirus Type 5	
		Murine Norovirus	
Bactericidal/ Yeast-cidal	EN 16615	Pseudomonas aeruginosa	1 min
		Staphylococcus aureus	
		Enterococcus hirae	