STERICOLOR TEST

BIOLOGICAL INDICATOR OF THE STEAM STERILIZATION PROCESS

Stericolor-test is a biological indicator consisting of a suspension of *Bacillus stearothermophilus* ATCC 7953 spores in a cultural ground containing a carbohydrate and a pH indicator.

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The thermoresistance of the spores present in the Stericolor Test vials is such that they survive after the treatment in autoclave for 5 minutes at $121^{\circ}\text{C} \pm 0,5^{\circ}\text{C}$ (1 atm), while they die after 12 minutes at $121^{\circ}\text{C} \pm 0,5^{\circ}\text{C}$ (1 atm).

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Stericolor Test is used to check the sterilization capability of steam autoclaves according to the recommendation of UPS XX. *Bacillus stearothermophilus* ATCC 7953 has bacteriological and biological characteristics such as to render it a biological indicator of absolute reliability and safety. The stock is not pathogenic, is not toxic, has no pyrogen activity and, unlike other microorganisms, it is hardly found as laboratory contaminant.

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WAY OF USE

Mod de folosire

To check the sterilization process, Stericolor Test should be placed inside the autoclave in those places where the sterilization conditions are critical (e.g.: at the autoclave bottom where air stagnation may be suspected, inside air stagnation of clothes, inside baskets, etc.)

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To check the sterilization process of material contained in vessels having dimensions greater than that of the indicator vial, the best results are obtained by inserting Stericolor Test inside the containers themselves (bottles etc.)

The number of vials to be used depends on the practicability condition of the autoclave and The

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Valabilitate
18 months

Preservation

Keep Stericor test at ambient temperature
Pastrati Stericor test la temperatura camerie .

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