Liquichek Sedimentation Rate Control Levels 1 and 2

REF	514 Level 1	4 x 9 mL
	515 Level 2	4 x 9 mL
	514X Bilevel MiniPak	2 x 9 mL

LOT	74370	Level 1	74371
		Level 2	74372

ENGLISH

INTENDED USE

Liquichek Sedimentation Rate Control is intended for use as a quality control to monitor the precision of erythrocyte sedimentation rate (ESR) laboratory testing procedures.

SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. Two levels of control are available to allow performance monitoring within the clinical range.

For customers in Germany: Quality control materials are required for assessment of laboratory performance as described in the "Guideline for Quality Assurance of Medical Laboratory Examinations following the German Medical Association" (Rili-BÄK regulation).

REAGENT

This product is composed of stabilized human red cells suspended in a buffered fluid and preservatives.

STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened at 18 to 30°C. Once opened, this product will be stable for 31 days when stored tightly capped at room temperature (18 to 30°C). Avoid prolonged exposure of opened vials to light. Vials should be kept tightly closed after use to avoid evaporation.

DO NOT FREEZE. DO NOT EXPOSE TO EXCESSIVE HEAT.

PROCEDURE

This product should be treated and analyzed the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit, or reagent being used.

1. This product must be prepared with a fresh tube each time.
2. Gently invert vial until packed cells have been resuspended. Continue mixing for 30 seconds. Avoid foaming. DO NOT VORTEX.
3. Follow the manufacturer's directions for filling sample tubes. The classical Westergren procedure (for example, Chase Disposable Westergren tubes) DOES NOT require prefilling of control material.
4. After each use, wipe any residual material from threads of vial and inside cap. Promptly replace the cap and store as described in the Storage and Stability Section.
5. Dispose of sample tubes; DO NOT re-use sample tubes.

Safety Precautions

- Caution should be used when handling this product to prevent splashing. Wear appropriate eye/face protection when using this product to protect from splashes.
- If, during the use of this device or as a result of its use, a serious incident occurs, please report it to Bio-Rad Laboratories and to your national health authority, as required.
- Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.

LIMITATIONS

1. This product should not be used past the expiration date.
2. This product is not intended for use as a standard.
3. Atypical physical appearance may be observed with this product. If the recovered values are not within the expected ranges, retest with a new bottle of this product. If the same result is obtained, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.
4. For Laboratory/Professional Use Only.

ASSIGNMENT OF VALUES

The mean values printed in this insert were derived from replicate analyses and are specific for this lot of product. The tests listed were performed by using manufacturer supported reagents and a representative sampling of this lot of product. Individual laboratory means should fall within the corresponding acceptable range; however, laboratory means may vary from the listed values during the life of this product. Variations over time and between laboratories may be caused by differences in laboratory technique, instrumentation calibrating method and reagents. It is recommended that each laboratory establish its own means and acceptable ranges and use those provided only as guides. Refer to www.qcnet.com for insert update information.

SPECIFIC PERFORMANCE CHARACTERISTICS

This product is a stabilized liquid product manufactured under rigid quality control standards. To obtain consistent vial-to-vial assay values, this product requires proper storage and handling as described.

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ENGLISH

Biological source material. Treat as potentially infectious.

Each human donor unit used to manufacture this product was tested as required by FDA accepted methods. Tests results were non-reactive or negative for evidence of infection due to Human Immunodeficiency Virus (HIV), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV). This product may also contain other human source materials for which there are no approved tests. In accordance with good laboratory practice, all human source material should be considered potentially infectious and handled with the same precautions used with patient specimens. Safety Data Sheet (SDS) available for professional users on www.bio-rad.com.

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Bio-Rad Laboratories comprehensive line of quality controls and QC data management solutions.

Autoimmune Controls
Blood Gas Controls
Cardiac Assessment Controls
Chemistry Controls
Coagulation Controls
Congenital/Pediatric Disease Controls

Diabetes/Hemoglobin Controls
Hematology Controls
Hepatitis & Retrovirus Controls
Immunoassay Controls
Immunology/Protein Controls
Molecular Controls
Sexually Transmitted Disease Controls

Specialty Infectious Disease Controls
Therapeutic Drug Monitoring Controls
Toxicology: Drugs-of-Abuse Controls
Toxicology: Specialty Controls
Urinalysis Controls
QC Data Management Solutions
External Quality Assurance Services (EQAS)

BIO-RAD

USA: 1-800-224-6723 | www.bio-rad.com/qualitycontrol

GLOSSARY	GLOSSAR	GLOSSAIRE	GLOSSARIO	GLOSARIO	GLOSSÁRIO	ORDLISTA	ORDLISTE
PARAMETER Erythrocyte Sedimentation Rate	PARAMETER Blutsenkungsgeschwindigkeit	PARAMETRE Vitesse de sédimentation globale	PARAMETRO Velocità di sedimentazione eritrocitaria	PARÁMETRO Velocidad de sedimentación globular	PARÁMETRO Taxa de sedimentação de eritrócitos	PARAMETER Sänkingsreaktion, erythrocyter	PARAMETER Blodsänkning (SR)
TERMS Glass Mean Plastic Range Type Units	BEGRIFFE Glas Mittelwert Kunststoff Bereich Typ Einheiten	TERMES Verre Moyenne Plastique Plage de valeurs Type Unités	TERMINI Vetro Media Plastica Intervallo Tipo Unità di misura	TÉRMINOS Vidrio Media Plástico Rango Tipo Unidades	TERMOS Vidro Média Plástico Limites Tipo Unidades	TERMER Glas Medelvärd Plast Område Typ Enheter	ORDLISTE Glas Middelvärdi Plastic Område Type Enheder

METHOD

Methode // Méthode // Metodo // Método // Método // Metod // Metode // Yöntem // 測定法

	Type	Units	Level 1 - 74371		Level 2 - 74372	
			Mean	Range	Mean	Range
ERYTHROCYTE SEDIMENTATION RATE						
Diesse Mini-Cube (1)	Plastic	mm/hr	11	± 5	89	± 25
Diesse MiniVes (Europe)	Plastic	mm/hr	17	± 5	88	± 25
Diesse Ves-matic 5 (1)	Plastic	mm/hr	5	± 5	53	± 25
Diesse Ves-matic 10 (1)	Plastic	mm/hr	11	± 5	81	± 25
Diesse Ves-matic 20 / 20 Plus New (1)	Plastic	mm/hr	10	± 5	89	± 25
Diesse Ves-matic 30 / 30 Plus (1)	Plastic	mm/hr	13	± 5	108	± 25
Diesse Ves-matic 60 / 60 Plus (1)	Plastic	mm/hr	14	± 5	96	± 25
Diesse Ves-matic Cube 30	Plastic	mm/hr	10	± 5	94	± 25
Diesse Ves-matic Cube 30 Touch (1)	Plastic	mm/hr	8	± 5	57	± 25
Diesse Ves-matic Cube 80 / 200	Plastic	mm/hr	6	± 5	76	± 25
Diesse Ves-matic Easy	Plastic	mm/hr	10	± 5	75	± 24
Diesse Ves-matic Easy (Europe)	Plastic	mm/hr	8	± 5	72	± 25
Greiner Bio-One SRS100/II SRT10/II (1)	Plastic	mm/hr	8	± 5	72	± 19
Kimble-Chase Disposable Westergren Tubes (2)	Glass	mm/hr	4	± 4	43	± 15
Polymedco Sedimat	Plastic	mm/hr	8	± 4	47	± 10
Polymedco Sediplast	Plastic	mm/hr	6	± 4	50	± 15
Sarstedt S-Sedivette	Plastic	mm/hr	8	± 5	68	± 14
Streck ESR-Auto Plus	Glass	mm/hr	§		§	
Vital Diagnostics Eocyte 20 (1)	Glass	mm/hr	9	± 6	84	± 12
Vital Diagnostics Eocyte 40 (1)	Glass	mm/hr	9	± 6	87	± 12
Vital Diagnostics Eocyte M / 10 (1)	Glass	mm/hr	9	± 6	81	± 12
Vital Diagnostics Eocyte Mini (1)	Glass	mm/hr	8	± 6	77	± 12
Vital Diagnostics MICROsed-System (1)	Glass	mm/hr	7	± 4	81	± 18
Vital Diagnostics Monitor Series (1)	Glass	mm/hr	8	± 4	81	± 18
Westergren (Modified)	Plastic	mm/hr	5	± 4	56	± 18

FOOTNOTES // Fussnoten // Notes de bas de page // Note a piè pagina // Notas al pie de página // Notas de rodapé // Fotnoter // Fodnoter // Dipnotlar // 脚注

ENGLISH

- (1) TEMPERATURE CORRECTION: Means and Expected Ranges are corrected to 18°C in accordance with the Manley table.
 - (2) METHOD VARIATION: Tubes do not contain Sodium Citrate (Diluent).
 - § The data required to establish the means and acceptable ranges for this assay were not obtained due to limited assignment participation. If your facility is interested in participating in the Value Assignment Program for this assay, please contact your local Bio-Rad office.
- ❖ INTERNATIONAL USE ONLY - The following section contains data for methods that are not available for diagnostic use in the United States.

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- INTERNATIONAL USE ONLY -

The following section contains data for methods that are not available for diagnostic use in the United States. ❖

INTERNATIONAL USE ONLY -

The following section contains data for methods that are not available for diagnostic use in the United States. ❖

METHOD

Methode // Méthode // Metodo // Método // Método // Metod // Metode // Yöntem // 測定法

	Type	Units	Level 1 - 74371		Level 2 - 74372	
			Mean	Range	Mean	Range
ERYTHROCYTE SEDIMENTATION RATE						
Linear Chemicals IRIA / THERMA Analyzers	Plastic	mm/hr	12	± 11	87	± 31
Linear Chemicals LENA 16/32 EDTA Analyzer	Plastic	mm/hr	4	± 4	53	± 16
Linear Chemicals THERMA NE / ERLINE AR / LENA Analyzers	Plastic	mm/hr	8	± 7	58	± 19
Wiener lab e32	Plastic	mm/hr	7	± 4	46	± 14
YHLO Vision / Vision Pro ESR Analyzer	Plastic	mm/hr	§		§	

FOR INTERNATIONAL USE

NOTES

MANUFACTURERS LISTED

Herstellerliste // Liste des fabricants // Elenco Produttori // Lista de fabricantes // Fabricantes enuncidos // Lista över tillverkare // Liste over producenter // Üretici Listesi // 製造元一覧

Diesse Diagnostica Senese, Siena, Italy

ELITECHGroup Inc., Lincoln, Rhode Island

Kimble-Chase Life Science and Research Products LLC., Vineland, New Jersey

Linear Chemicals, Spain

Polymedco, Cortlandt Manor, New York

SARSTEDT, Nümbrecht, Germany

Shenzen YHLO Biotech Co., Ltd., Shenzhen, China

Streck Laboratories, La Vista, Nebraska

Thermo Fisher Scientific Inc., Waltham, Massachusetts

Wiener Laboratories SAIC, Rosario, Argentina.

U.S. Patents: 4,968,629; 5,863,799; 5,888,822; 6,017,764; 6,051,433; 6,331,435; 6,342,391

TRADEMARK // Marke // Marque Déposée // Marchio Commerciale // Marca Comercial // Marca Comercial // Varumärke // Varemærke // Ticari Marka // 商標

ENGLISH

Liquichek is a trademark of Bio-Rad Laboratories, Inc., in certain jurisdictions.

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Liquichek é uma marca registrada da Bio-Rad Laboratories, Inc. em algumas jurisdições.
UNITED STATES, Bio-Rad Laboratories Inc.