



235 x 120 x 45

235 x 120 x 45

intouch™ POWDER FREE NITRILE EXAMINATION GLOVES

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ISO 374-5: 2016

VIRUS

EN POWDER FREE NITRILE EXAMINATION GLOVES • Non sterile • Single Use • Recommended for medical use

Medical Device Regulation (EU) 2017/745 & Personal Protective Equipment Regulation (EU) 2016/425 (Category III) • Medical Devices Regulations 2002 (SI 2002 No. 618, as amended) (UK MDR 2002) • EU Notified Body for (Module B) and Module C2: SATRA Technology Europe Limited (2777) • UK Approved Body for (Module B) and Module C2: SATRA Technology Center Limited (AB0321) • EU Declaration of Conformity • UK Declaration of Conformity is accessible at www.intouchcares.com • This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. • The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical used in a mixture. • It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. • When used, protective glove may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. • Before usage, inspect the gloves for any defect or imperfections. • EN ISO 374-4:2019 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. • The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimens. • Donning – Hold glove by the bead with one hand. Align the glove thumb with your other hand thumb and slide your hand into the glove, one finger into each glove finger. Pull by the glove palm to get a good fit. Don the other glove by the same procedure. • Doffing – Hold glove bead and pull toward the finger until the glove come off. • Where relevant, a list of the substances contained in the glove which are known to cause allergies, per listed in Annex G of EN ISO 21420:2020, shall be supplied on request. • Components used in glove manufacturing may cause allergic reactions in some users. If allergic reactions occur, seek for medical advice immediately. • In storage, avoid excessive heat. Open box should be shielded from exposure to direct sun or fluorescent lighting. Gloves are packed in dispenser which is suitable for transport. Keep the gloves in box when not in use.

Minimum Breakthrough Duration Time > 240min	Chemotherapy Drug & Concentration
	Cisplatin 1.0 mg/mL, Cyclophosphamide (Cytoxan) 20.0 mg/mL, Cytarabine 100 mg/mL, Doxorubicin (DTIC) 10.0 mg/mL, Doxorubicin Hydrochloride 2.0 mg/mL, Etoposide 20.0 mg/mL, Fluorouracil 50.0 mg/mL, Ifosfamide 50.0 mg/mL, Methotrexate 25.0 mg/mL, Mitomycin C0.5 mg/mL, Mitoxantrone 2.0 mg/mL, Paclitaxel (Taxol) 6.0 mg/mL, Vincristine Sulfate 1.0 mg/mL
	Fentanyl citrate and concentration
	Fentanyl citrate injection (100 mcg/2mL)

30.2min Thiotepa (10.0 mg/mL)
10.1min Carmustine (BCNU) (3.3 mg/mL)

Warning: Not recommended for use with Carmustine and Thiotepa.

Tested for use with Chemotherapy Drugs using ASTM D6978.

Resistance to Permeation by Chemicals	Level	Mean Degradation (%)
*4% Chlorhexidine Digluconate	6	19.0
40% Sodium Hydroxide (K)	6	-42.9
10-13% Sodium Hypochlorite	6	14.7
50% Sulphuric Acid	6	-20.5
5% Ethidium Bromide	6	-3.4
50% Glutaraldehyde	6	27.4
0.1% Phenol	6	33.8
1.5% Methanol in Water	6	21.9
3% Povidone-iodine	6	33.7
10% Sodium Percarbonate	6	15.4
10% Acetic Acid	4	66.7
37% Formaldehyde (T)	3	5.0
30% Hydrogen Peroxide (P)	2	22.8
65% Ammonium Hydroxide (O)	0	-52.0
65% Nitric Acid (M)	0	97.6
70% Isopropanol	0	62.2
35% Ethanol	0	38.8
99% Acetic Acid (N)	0	93.9

* Permeation rate 7µg/ cm²/min

Permeation Performance Level & Measured Breakthrough Time (minutes)
0 > * 1 > 10 min, 2 > 30 min, 3 > 60 min, 4 > 120 min, 5 > 240 min, 6 > 480 min

*Indicates that the glove falls below the minimum performance level as stated in EN ISO 374-1:2016 +A1:2018 for the given individual hazard.

G-LF-XXXXXX-F-R0
Last update in 2022

FR POWDER FREE NITRILE EXAMINATION GLOVES • Non sterile • Single Use • Recommended for medical use

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RO POWDER FREE NITRILE EXAMINATION GLOVES • Non sterile • Single Use • Recommended for medical use

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