EN

INSTRUCTIONS FOR USE

INDICATIONS FOR USE:

The CERENOVUS NIMBUS™ Device (the Device) is intended to be used to restore blood flow in patients experiencing an acute ischemic stroke due to a large vessel neurovascular occlusion. The Device is designed for use in the anterior and posterior neurovasculature in vessels of diameter 1.5 mm to 5.0 mm, such as the internal carotid artery, the M1 and M2 segments of the middle cerebral artery, the A1 and A2 segments of the anterior cerebral artery, the basilar, the posterior cerebral and the vertebral arteries.

The Device should only be used by physicians trained in neurointerventional catheterization and the treatment of ischemic stroke.

DEVICE DESCRIPTION:

The Device is composed of a three-dimensional nitinol cage assembly at the distal end of a tapered nitinol shaft as shown in Figure 1 below. The Device is supplied pre-loaded in an Insertion Tool.


Figure 1: Device Illustration
1. Shaft/Push Wire 2. Proximal Radiopaque Coil
3. Outer Cage 4. Mid-Markers (X2) 5. Distal Markers (X2)

Table 1: Device Length details (from Figure 1.)

Catalogue Number	A	B	C	D	E
	Distal Outer Cage Diameter	Working Length	Proximal Outer Cage Diameter	Proximal Coil Length	Overall Length
GCE4528	4.5 mm	28 mm	2.25 mm	20 mm	1927 mm

COMPATIBILITY:

Microcatheter: The Device should be introduced and delivered through a size “18” microcatheter with an internal diameter of 0.021” and an overall length of not more than 155 cm. The Device is compatible with a Prowler® Select™ Plus microcatheter. Performance of the Device may be impacted if a different microcatheter is used.

Guide Catheter / Intermediate Catheter: The Device should be used in conjunction with a Neurovascular Guide or Sheath with an internal diameter of 0.075” or greater. A Balloon Guide Catheter may be used. An Intermediate Catheter may also be used.

CONTRAINDICATIONS:

- Allergy or hypersensitivity to Nickel-Titanium.
- Excessive vessel tortuosity that may prevent device delivery.
- Patients with angiographic evidence of carotid dissection.

POTENTIAL COMPLICATIONS:

Possible complications include but are not limited to: air embolism, hematoma or hemorrhage at puncture site, infection, distal embolization, vascular spasm, neurological deterioration including stroke and death, thrombosis, dissection or perforation, device malfunction or fracture, vascular occlusion, ischemia, intracranial hemorrhage, pseudo-aneurysm formation, or other minor or major complications related to cerebral angiography.

WARNINGS & PRECAUTIONS:

- Store in a cool, dry place.
- Do not use the Device if the product sterilization barrier or its packaging is compromised.
- Do not use after the expiration date printed on the product label.
- Do not deploy in vessels of less than 1.5 mm diameter.
- Do not use damaged or kinked Devices.
- Do not advance the Device through the microcatheter against significant resistance.
- Do not torque or rotate the Device.
- Do not advance the Device after it has been deployed or partially deployed from the microcatheter.
- Do not withdraw the Device against significant resistance. Assess the cause of resistance using fluoroscopy and if necessary advance the microcatheter over the Device to resheath or partially resheath to aid withdrawal. Resheathing a Device with captured clot may result in loss of clot and distal embolization.
- This Device is for use in one patient only. Do not resterilize or reuse in more than one patient. Reusing the Device in another patient could result in compromised Device performance, cross infection and other safety related hazards.
- Do not use any CERENOVUS NIMBUS™ Revascularization Device more than 3 times.
- Do not perform more than 3 retrieval attempts in any one vessel.
- Do not withdraw the deployed Device through a previously stented vessel. If it becomes necessary to cross a deployed stent to access a clot, first ensure that the stent can be crossed with the Guide or Intermediate Catheter. The Device should then be fully pulled back into Guide or Intermediate Catheter prior to retrieval through the stent.
- Users shall take all necessary precautions to limit X-radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible. The risks of angiographic and fluoroscopic X-ray radiation include alopecia, burns ranging in severity from skin reddening to ulcers, cataracts and delayed neoplasia that increase in probability as procedure time and number of procedures increase.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

PROCEDURAL INFORMATION:

The Device is delivered endovascularly under fluoroscopic guidance in a similar manner to that of other neurovascular clot retrieval systems. The use of antiplatelet and anticoagulation agents is recommended at the discretion of the treating physician.

Preparation & Delivery:

1. Using standard interventional techniques, access the arterial system and using angiography, determine the location of the occluded vessel.
2. Advance an appropriate Guide Catheter, Sheath or Balloon Guide Catheter as close to the occlusion site as possible. Connect a rotating hemostasis valve (RHV) to the proximal end of this catheter and connect to a continuous flush system.
3. Select an appropriate microcatheter (see “Compatibility” section). Connect an RHV to the proximal end of the microcatheter and connect to a continuous flush system.
4. With the aid of a suitable Guidewire, and using standard catheterization techniques and fluoroscopic guidance, advance the Microcatheter up to and across the occlusion so that the distal end of the Microcatheter is positioned distal to the occlusion.
5. Remove the guidewire from the microcatheter and, if desired, gently infuse contrast media through the microcatheter to visualize the distal end of the occlusion.
6. Remove the Insertion Tool with the preloaded Device from the packaging hoop, taking care not to accidentally retract or advance the Device from the Insertion Tool while doing so.
7. Insert the distal end of the Insertion Tool through the RHV of the microcatheter and wait until fluid is seen exiting the proximal end of the Insertion Tool, confirming that the Device is flushed. Advance the Insertion Tool until it is fully seated in the hub of the microcatheter. Fully tighten the RHV to hold the Insertion Tool securely in position.
8. Confirm that the Insertion Tool is fully seated in the hub of the RHV before proceeding to advance the Device. Advance the Device until at least half of the shaft length has been inserted into the microcatheter, at which point the Insertion Tool may be removed.
9. Do not discard the Insertion Tool as it will be needed if additional device passes are required during the procedure.

Positioning and Deployment:

1. Continue to advance the Device until the radiopaque distal markers of the Device approach the distal region of the Microcatheter.
Note: There are two silver bands positioned on the proximal portion of the Device shaft to assist in minimizing the amount of fluoroscopic exposure required during Device insertion. For example, if using a standard microcatheter (total length of 155 cm and a 7 cm RHV), then when the first silver band on the shaft approaches the RHV, the distal end of the device will be approximately 8 cm from the distal end of the microcatheter. When the second silver band on the Device shaft approaches the RHV, the distal end of the device will be nearing the distal end of the microcatheter.
2. The device should be positioned in the clot ideally such that the end of the proximal radiopaque coil (Figure 2) is aligned with the proximal face of the clot.


Figure 2

3. While stabilizing the Device shaft, retract the microcatheter to deploy the Device within the clot. The Device will be fully deployed once the distal tip of the Microcatheter has been retracted to a position over the proximal radiopaque coil of the Device.
4. After deployment, the device may be left to embed in the clot for up to 5 minutes prior to withdrawal.

Retrieval:


1. A torque device may be advanced over and tightened against the Device shaft to achieve additional control of the Device during retrieval.
2. If using a Balloon Guide Catheter, inflate the balloon to occlude the vessel as per the Balloon Guide Catheter IFU.
3. Prior to clot retrieval re-advance the microcatheter to the clot while holding the CERENOVUS NIMBUS™ push wire static. Do not continue to advance against significant resistance.
4. Withdraw the Device and microcatheter slowly and carefully as a single unit to the Guide or Intermediate Catheter while aspirating through the guide or Intermediate Catheter maintaining microcatheter and CERENOVUS NIMBUS™ position relative to each other during the withdrawal step.
5. Applying vigorous aspiration (by syringe), withdraw the Device and microcatheter into the Guide or Intermediate Catheter and continue to aspirate until the Device reaches the proximal RHV.
Note: If withdrawal in the guide or intermediate catheter is difficult (as may be the case with a large clot burden) then deflate the balloon (if applicable) and withdraw the guide or intermediate catheter, microcatheter and Device together through the introducer sheath.

6. Disconnect the RHV from the Guide Catheter and remove the Device, microcatheter and RHV together from the Guide Catheter.
7. Use a syringe to further aspirate the Guide and/or Intermediate Catheter to ensure it is clear of any thrombus material.

Cleaning and Re-use:

1. The Device may be used for up to three retrieval attempts.
2. If an additional pass is to be made with the Device then carefully remove any captured thrombus from the Device and clean the Device in heparinized saline, rubbing gently from proximal to distal to remove any residual thrombus material. Inspect the Device carefully and if any damage or deformity is observed do not use the Device – instead use a new Device for any subsequent passes. If using the same Device, replace the Insertion Tool onto the proximal shaft and retract the Device into the Insertion Tool until it is fully resheathed. Do not retract the device into the Insertion Tool against significant resistance – ascertain the cause of the resistance and if required use a new device for any subsequent passes.
3. Do not attempt more than three retrieval attempts in one vessel.

INFORMATIVE SYMBOLS FOR CERENOVUS NIMBUS™:

	Caution		Batch Code
	Consult instructions for use		Catalogue Number
	Do not use if package is damaged		Sterilized using ethylene oxide
	Do not reuse		Contents
	Keep dry		Minimum microcatheter inner diameter (ID)
	Keep away from sunlight		Use-by date
	Manufacturer		

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Neuravi LTD