HbA1c Control P

Package size

4×1 mL

Intended use

HbA1c Control P is used in routine chemistry analytes quality control by monitoring accuracy and precision of Mindray BS measurement system and test ability of clinical laboratory.

Summary

HbA1c Control P contains the following analytes:

Hemoglobin A1c (HbA1c)

Mindray BS measurement system is composed of Mindray BS-series chemistry analyzers, Mindray reagent kits, calibrators and control materials.

Components

HbA1c Control P is a lyophilized control based on human blood. The concentrations or activities of the control components are lot-specific and almost at abnormal or pathological levels for the methods used.

Warnings and Precautions

1. For in vitro diagnostic use only.
2. Professional or trained person is required. Please take the necessary precautions for handling all laboratory reagents.
3. The control values are lot-specific with the matched models listed in the value sheet.
4. The control should be run daily in parallel with the patient samples and after calibration or reagent lot change.
5. Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.
6. The control was tested with Elisa methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen¹.
7. Material safety data sheet is available on request.
8. Disposal of all waste material should be in accordance with local guidelines.

Preparation

1. Take the vial out of the freezer and allow reaching room temperature.
2. Tap the vertically positioned vial gently and ensure that the lyophilized

- material is at the bottom of the vial.
3. Remove the screw cap and rubber stopper carefully, avoiding any loss of lyophilizate.
 4. Reconstitute by adding 1.0 mL distilled/deionized water to the side of the vial slowly and exactly.
 5. Carefully replace the rubber stopper and leave to stand for 30 minutes at room temperature.
 6. During the standing time, mix contents by inverting several times and swirling gently to ensure that all components of lyophilizate are dissolved. Avoid the formation of foam.
 7. Dispense the required volume into a sample cup, and then run QC test you have requested or have it analyzed in the same way as patient samples.

Storage and Stability

The control is stable up to expiration date indicated on the label when stored in unopened vial at 2~8°C and protected from light. Once reconstituted, please store control tightly capped when not in use. Once thawed to room temperature, the control should not be re-frozen.

Once opened, the reagents are stable for 15 days when refrigerated on the analyzer or refrigerator.

Materials required but not provided

1. Distilled/deionized water.
2. Mindray BS-series chemistry analyzer, calibrator and reagent kit.
3. General laboratory equipment.

Assay Values

The control assay values (the target value and range) determined by Mindray standard transfer procedure and routine method are listed in the value sheet. The target value was obtained from Mindray measurement system and the range was calculated as the target values ± 3 standard deviations. The control values are specific for different lot numbers and different model of chemistry analyzers, so please check the lot number and model of analyzer before use.


Quality control

The control result should be within the defined range as shown in the control value sheet. If the control result falls outside the range, the measurement system should be checked. For instance, position of reagent or sample placed on the analyzer; expiration date or storage condition of the calibrator, reagent and control; the parameter setting in software; performance of the analyzer; or the calibration process. Each laboratory should establish its own internal quality control scheme and procedure for corrective action.

References

1. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.

Graphical symbols

© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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