# Glu(GOD)

## Glucose Kit (GOD-POD Method)

### Order Information

| Cat. No.      | Package size              |
|---------------|---------------------------|
| 105-000849-00 | R1: 4×35 mL + R2: 2×18 mL |
| 105-000888-00 | R1: 4×40 mL + R2: 2×20 mL |
| 105-001624-00 | R1: 6×57 mL + R2: 3×32 mL |

### Intended Purpose

In vitro test for the quantitative determination of Glucose (Glu) concentration in human serum and plasma on Mindray BS series chemistry analyzers. It is intended to be used for diagnosis and monitoring therapeutic effect of diabetes.


### Summary

Carbohydrates supply the body energy with glucose, which is the most important monosaccharide in blood, and it is an indispensable energy supplier for cellular function<sup>1</sup>.

Measuring blood glucose is used for the diagnosis of carbohydrate metabolism disorders and monitoring of treatment in diabetes mellitus, neonatal hypoglycemia, idiopathic hypoglycemia, pharmlc hypoglycemia and insulinoma<sup>2, 3</sup>.

### Assay Principle

Glucose oxidase-Peroxidase (GOD-POD) method

By the catalysis of GOD (Glucose oxidase), glucose is oxidated to yield H<sub>2</sub>O<sub>2</sub>, and then in the presence of POD (Peroxidase), H<sub>2</sub>O<sub>2</sub> oxidates 4-Aminoantipyrine with p-Hydroxybenzoic acid sodium to form a colored dye of quinoneimine. The absorbency increase is directly proportional to the concentration of glucose.

### Reagents Components

|            |                                     |            |
|------------|-------------------------------------|------------|
| <b>R1:</b> | Phosphate buffer                    | 70 mmol/L  |
|            | Sodium azide                        | 0.08%      |
|            | Peroxidase                          | >4.7 KU/L  |
|            | N-ethyl-N-(2-hydroxy-3-sulfopropyl) |            |
|            | -3-methylaniline, sodium salt       | 1.5 mmol/L |

|            |                   |            |
|------------|-------------------|------------|
| <b>R2:</b> | Phosphate buffer  | 70 mmol/L  |
|            | Sodium azide      | 0.08%      |
|            | 4-Aminoantipyrine | 2.5 mmol/L |
|            | Glucose oxidase   | >100KU/L   |

### **Storage and stability**

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

On board in use, the reagents are stable for 30 days when refrigerated on the analyzer.

Contamination must be avoided.

Do not freeze the reagent.

### **Specimen collection and preparation**

#### **■ Specimen types**

Serum, heparin and NaF plasma are suitable for samples.

#### **■ Preparation for Analysis**

1. Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.
2. Centrifuge samples containing precipitate before performing the assay.
3. Specimens should be tested as soon as possible after sample collection and pre-analytical treatment.

#### **■ Sample Stability**

Serum/heparin plasma:

8 hours at 15-25°C

72 hours at 2-8°C

NaF plasma<sup>4</sup>: 24 hours at 15-25°C

The stability of glucose in specimens is affected by storage temperature, bacterial contamination, and glycolysis. Plasma or serum samples without preservative (NaF) should be separated from the cells or clot within half an hour of being drawn. When blood is drawn and permitted to clot and to stand uncentrifuged at room temperature, the average decrease in serum glucose is ~ 7 % in 1 hour (0.28 to 0.56 mmol/L or 5 to 10 mg/dL). This decrease is the result of glycolysis. Glycolysis can be inhibited by collecting the specimen in fluoride tubes<sup>3</sup>.

Sample stability claims were established by manufacturer and/or based upon references, each laboratory should establish its own sample stability criteria.

**Reagent Preparation**

R1 and R2 are ready to use.

Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

**Materials required but not provided**

1. General laboratory materials: NaCl solution 9 g/L (saline), distilled/deionized water.
2. Calibrator and Control: Please check the section of reagent instruction of Calibration and Quality Control.
3. Mindray BS series chemistry analyzers and General laboratory equipment.

**Assay procedure**

| Parameters Item                                                                                               | BS-2000 chemistry analyzers |
|---------------------------------------------------------------------------------------------------------------|-----------------------------|
| Assay type                                                                                                    | Endpoint                    |
| Wavelength (Primary/Secondary)                                                                                | 505/660 nm                  |
| Reaction direction                                                                                            | Increase                    |
| R1                                                                                                            | 160 $\mu$ L                 |
| Sample or Calibrator                                                                                          | 2 $\mu$ L                   |
| Mix, incubate at 37°C for 5 min, read the absorbance A1, then add:                                            |                             |
| R2                                                                                                            | 40 $\mu$ L                  |
| Mix thoroughly, incubate at 37°C for 10 min, read the absorbance A2,<br>Then calculate $\Delta A = (A2 - A1)$ |                             |

Parameters may vary in different chemistry analyzers, may adjust in proportion if necessary. For Mindray BS series chemistry analyzers, Reagent Parameters is available on request. Please refer to the appropriate operation manual for the analyzers.

**Calibration**

1. It is recommended to use the Mindray Calibrator (Multi Sera Calibrator: 105-001144-00 or other suitable calibrators) and 9 g/L NaCl (saline) for two-point calibration. Traceability of the Mindray Multi Sera Calibrator can refer to the calibrator instructions for use of Mindray Company.
2. Calibration frequency

Calibration is stable for approximately 30 days on BS-2000 chemistry analyzers. The calibration stability may vary in different instruments, each laboratory should set a calibration frequency in the instrument parameters appropriate to their usage pattern.

Recalibration may be necessary when the following occur:

- As changed reagent lot.
- As required following quality control procedures or out of control.
- As executes specific maintenance or troubleshooting procedure of chemistry analyzers.

3.The calibrator values are lot-specific with the matched models listed in the value sheet.

## Quality control

- 1.It is recommended to use the Mindray Control (ClinChem Multi Control: 105-009119-00, 105-009120-00 or other suitable controls) to verify the performance of the measurement procedure; other suitable control material can be used in addition.
- 2.Two levels of control material are recommended to analyze each batch of samples. In addition, the control should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or trouble shooting procedures as detailed in the appropriate system manual.
- 3.Each laboratory should establish its own internal quality control scheme and procedures for corrective action if control doesn't recover within the acceptable tolerances.

## Calculation

$\Delta A = (A_2 - A_1)$  sample or calibrator.

$C \text{ sample} = (\Delta A \text{ sample} / \Delta A \text{ calibrator}) \times C \text{ calibrator}$ .

The BS series chemistry analyzer detects the change of absorbance( $\Delta A$ ) and calculates the Glu concentration of each sample automatically after calibration.

Conversion factor:  $\text{mg/dL} \times 0.0555 = \text{mmol/L}$

## Dilution

If the value of sample exceeds 25 mmol/L, the sample should be diluted with 9 g/L NaCl solution (saline) (e.g. 1+4) and rerun; the result should be multiplied by 5.

## Expected values

| Sample Type               |        | mg/dL  | mmol/L  |
|---------------------------|--------|--------|---------|
| Serum/Plasma <sup>5</sup> | Adults | 70-110 | 3.9-6.1 |

The expected value above is provided from reference, and Mindray has verified it by 168 samples of people from China.

The reference intervals according to Tietz<sup>4</sup> are as follow:

| Sample Type  |             | mg/dL  | mmol/L  |
|--------------|-------------|--------|---------|
| Serum/Plasma | 1 day       | 40-60  | 2.2-3.3 |
|              | >1 day      | 50-80  | 2.8-4.4 |
|              | Children    | 60-100 | 3.3-5.6 |
|              | Adults      | 74-106 | 4.1-5.9 |
|              | 60-90 years | 82-115 | 4.6-6.4 |
|              | >90 years   | 75-121 | 4.2-6.7 |

Each laboratory should establish its own reference intervals based upon its particular locale and population characteristics since expected values may vary with geography, race, sex and age.

## Performance Characteristics

### Analytical Sensitivity

The Glucose Kit has an analytical sensitivity of 0.3 mmol/L on BS-2000. Analytical sensitivity is defined as the lowest concentration of analyte that can be distinguished from a sample that contains no analyte. It is calculated as the value lying 3 standard deviations above that of the mean from 20 replicates of an analyte-free sample.

### Measuring range

The Mindray BS series systems provide the following linearity range:

| Sample Type  | Units         |
|--------------|---------------|
| Serum/Plasma | 0.3-25 mmol/L |

A high Glu concentration sample (approximately 25 mmol/L) is mixed with a low concentration sample (<0.3 mmol/L) at different ratios, generating a series of dilutions. The Glu concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient  $r \geq 0.990$ . The reportable range is 0.3-125 mmol/L.

### Precision

Precision was determined by following CLSI Approved Guideline EP05-A3<sup>6</sup>, each sample was assayed 2 times per run, 2 runs per day, a total of 20 days. The precision data of controls and human samples on BS-2000 are summarized below\*.

| Specimen Type<br>(N=80) | Mean<br>(mmol/L) | Repeatability  |      | Within-Lab     |      |
|-------------------------|------------------|----------------|------|----------------|------|
|                         |                  | SD<br>(mmol/L) | CV % | SD<br>(mmol/L) | CV % |
| Control Level 1         | 5.23             | 0.02           | 0.42 | 0.06           | 1.16 |
| Control Level 2         | 14.46            | 0.06           | 0.38 | 0.18           | 1.24 |
| Serum 1                 | 2.86             | 0.01           | 0.39 | 0.06           | 1.98 |
| Serum 2                 | 7.32             | 0.03           | 0.37 | 0.05           | 0.74 |
| Serum 3                 | 10.54            | 0.04           | 0.34 | 0.09           | 0.89 |

\*Representative data, results in different instruments, laboratories may vary.

### Analytical Specificity

The samples with different concentration interfering substance were prepared by addition of interferent to human serum pools, and recovers are within  $\pm 10\%$  of the corresponding control value to be considered as no significant interference.

## Glu(GOD)

**mindray**

No significant interference was observed when the following substances were tested for interference with this methodology. The data of interference studies on BS-2000 are summarized below.

| Interfering Substance | Interferent Concentration (mg/dL) | Analyte Concentration (mmol/L) | Relative Deviation (%)* |
|-----------------------|-----------------------------------|--------------------------------|-------------------------|
| Ascorbic acid         | 30                                | 4.40                           | -1.92                   |
| Hemoglobin            | 500                               | 4.36                           | -1.00                   |
| Bilirubin             | 40                                | 4.38                           | +3.82                   |
| Intralipid            | 500                               | 4.40                           | -2.82                   |

*\*Representative data, results in different instruments, laboratories may vary.*

In very rare cases gammopathy, in particular type IgM, may cause unreliable results<sup>7</sup>.

### ■ Method Comparison

Correlation studies were performed using CLSI Approved Guideline EP09-A3<sup>8</sup>. The Mindray System (Mindray BS-2000/Mindray Glu Reagent) (y) was compared with comparison system (Hitachi 7180/Roche Glu Reagent) (x) using the same serum specimens. The statistical data obtained by linear regression are shown in the table below\*:

| Regression Fit     | Correlation Coefficient (r) | Sample (N) | Concentration Range (mmol/L) |
|--------------------|-----------------------------|------------|------------------------------|
| $y=0.9953x+0.0339$ | 0.9996                      | 100        | 2.26-24.07                   |

*\*Representative data, results in different instruments, laboratories may vary.*

### Result interpretation

In high altitude areas with low oxygen content, or in samples with extremely high concentration, falsely low results may occur in patient samples.

The results could be affected by drugs, disease, or endogenous substances<sup>7,9</sup>. When the reaction curve is abnormal, it is recommended to retest and check the result.

### Warnings and precaution

1. For in vitro diagnostic use only. For laboratory professional use.
2. Please take the necessary precautions for handling all laboratory reagents.
3. Please confirm the integrity of the package before use. Do not use the kits with damaged packages. The reagents avoid direct exposure to sunlight and freezing. The results can't be assured when stored at inappropriate condition.

4. If unintentionally opened before used, store the reagents tightly capped at 2-8°C and protected from light, and the stability is equally to in-use stability.
5. Do not mix reagents with different lots and bottles.  
Do not use the reagents beyond the expiration date and the in-use date.  
Do not mix fresh reagents with in-use reagents.  
Avoid the formation of foam.
6. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration or controls do not meet the insert and/or the Mindray System criteria.
7. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
8. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
9. When the reagents accidentally enter the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
10. Safety data sheet is available for professional user on request.
11. Disposal of all waste material should be in accordance with local guidelines.
12. All human material should be considered potentially infectious.
13. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable.
14. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
15. This kit contains components classified as follows in accordance with the Regulation (EC) No 1272/2008:

**Warning**

|      |                                |
|------|--------------------------------|
| H315 | Causes skin irritation.        |
| H319 | Causes serious eye irritation. |

**Prevention:**

|      |                                                                                  |
|------|----------------------------------------------------------------------------------|
| P280 | Wear protective gloves, protective clothing, eye protection and face protection. |
| P264 | Wash all exposed external body areas thoroughly after handling                   |

| Response:          |                                                                                                                                  |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------|
| P305+P351<br>+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P337+P313          | If eye irritation persists: Get medical advice/attention.                                                                        |
| P302+P352          | IF ON SKIN: Wash with plenty of water.                                                                                           |
| P332+P313          | If skin irritation occurs: Get medical advice/attention.                                                                         |
| P362+P364          | Take off contaminated clothing and wash it before reuse.                                                                         |

## References

1. Marshall WJ, Lapsley M, Day AP, Ayling RM, eds. Clinical Biochemistry: Metabolic and Clinical Aspects. 3rd ed. Churchill Livingstone 2014;273-275.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 131-137.
3. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;709-722.
4. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 444.
5. Shang H, Wang YS, Shen ZY, eds. National Guide to Clinical Laboratory Procedures, 4th ed. National Health Commission of the People's Republic of China, 2014; 233.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007,45(9):1240-1243.
8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-349,3-371.

## Graphical symbols

**Glu(GOD)**

**mindray**

In Vitro Diagnostic  
medical device

Unique device  
identifier

European  
Conformity

Consult Instructions  
For use

Use-by  
date

Authorized representative in  
the European Community

Batch Code

Temperature  
limit

Manufacturer

Catalogue  
number

Keep away from sunlight

Indicates a medical device that needs protection from light sources

© 2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

**E-mail Address:** service@mindray.com

**Website:** www.mindray.com

**Tel:** +86-755-81888998; **Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH(Europe)

**Address:** Eiffestraße 80, Hamburg 20537, Germany

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1. Marshall WJ, Lapsley M, Day AP, Ayling RM, eds. Clinical Biochemistry: Metabolic and Clinical Aspects. 3rd ed. Churchill Livingstone 2014;273-275.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 131-137.
3. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;709-722.
4. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 444.

5. Shang H, Wang YS, Shen ZY, eds. National Guide to Clinical Laboratory Procedures, 4th ed. National Health Commission of the People's Republic of China, 2014; 233.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007;45(9):1240-1243.
8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-349,3-371.

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## Referencias

1. Marshall WJ, Lapsley M, Day AP, Ayling RM, eds. Clinical Biochemistry: Metabolic and Clinical Aspects. 3rd ed. Churchill Livingstone 2014; 273-275.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 131-137.
3. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012; 709-722.
4. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 444.
5. Shang H, Wang YS, Shen ZY, eds. National Guide to Clinical Laboratory Procedures, 4th ed. National Health Commission of the People's Republic of China, 2014; 233.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007, 45(9):1240-1243.
8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-349, 3-371.

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**Bibliografia**

1. Marshall WJ, Lapsley M, Day AP, Ayling RM, eds. Clinical Biochemistry: Metabolic and Clinical Aspects. 3rd ed. Churchill Livingstone 2014;273-275.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Francoforte: TH-Books Verlagsgesellschaft 1998; 131-137.
3. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;709-722.
4. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 444.
5. Shang H, Wang YS, Shen ZY, eds. National Guide to Clinical Laboratory Procedures, 4th ed. National Health Commission of the People's Republic of China, 2014; 233.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007;45(9):1240-1243.
8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-349,3-371.

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**Indirizzo:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, Repubblica popolare cinese

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## Referanslar

1. Marshall WJ, Lapsley M, Day AP, Ayling RM, eds. Clinical Biochemistry: Metabolic and Clinical Aspects. 3rd ed. Churchill Livingstone 2014;273-275.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 131-137.
3. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders

2012;709-722.

- 4.Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 444.
- 5.Shang H, Wang YS, Shen ZY, eds. National Guide to Clinical Laboratory Procedures, 4th ed. National Health Commission of the People's Republic of China, 2014; 233.
- 6.CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- 7.Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007;45(9):1240-1243.
- 8.CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
- 9.Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-349,3-371.

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**mindray**

**Üretici:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Adres:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. Çin

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## Références

- 1.Marshall WJ, Lapsley M, Day AP, Ayling RM, eds. Biochimie clinique : Aspects métaboliques et cliniques. 3rd ed. Churchill Livingstone 2014;273-275.
- 2.Thomas L. Clinical Laboratory Diagnostics (diagnostics cliniques de laboratoire) : Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 131-137.
- 3.Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;709-722.
- 4.Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 444.
- 5.Shang H, Wang YS, Shen ZY, eds. National Guide to Clinical Laboratory Procedures, 4th ed. National Health Commission of the People's Republic of China, 2014; 233.
- 6.CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- 7.Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007,45(9):1240-1243.
- 8.CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
- 9.Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-349,3-371.

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© 2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés
**Fabricant** : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
**Représentant en Europe** : Shanghai International Holding Corp. GmbH (Europe)