Description
This medical device is a monovisual intraocular lens. Its sterile and is in single use.
This intraocular lens is indicated for the replacement of a cataractous natural crystalline lens for visual correction of aphakia after cataract surgery. It is a one-piece and posterior chamber foldable monovisual intraocular lens. It is designed for implantation in the capsular bag and to provide the function of the natural crystalline lens by its dioptric power and its UV (Ultraviolet) filter.

Intraocular lens properties
Each intraocular lens is obtained by machining a hydrophilic acrylic material with a hydration rate of 25% that incorporates a filter to protect the retina from UV radiation.
The UV cut-off wavelength is 340 nm (transmission < 10%).
CRISTALENS Industrie recommends using intraocular lenses with equivalent transmission spectra in both of a patient's eyes.
Optics: monofocal, square edges.
Modulation Transfer Function (MTF) > 0.45 (measured @ 100 c/mm for a 3-mm aperture).
Reflective index: 1.46.
Available power range: refer to the technical documentation.
Intraocular lens dimensions: refer to labelling.

Content
The protective packaging contains the sterile product, these instructions for use, the patient card and self-adhesive traceability labels.
Each intraocular lens, packed individually in a double packaging, is sterile and sealed.
The external blister contains the internal blister including the intraocular lens held in place on a holder, the whole hydrated by saline solution.


Indications & claimed clinical performances
This posterior chamber monovisual intraocular lens is indicated for the replacement of a cataractous natural crystalline lens for visual correction of aphakia after cataract surgery. It is designed for implantation in the capsular bag and to provide the natural crystalline lens function.

Special considerations in patient selection
This intraocular lens is contraindicated in newborns, premature babies, infants and children.
Intraocular lens implantation may not be suitable in patients with any of the following conditions (non-exhaustive list):
Non-age-related cataract (e.g. traumatic cataract).
History of ocular or refractive surgery.
Use of systemic or ocular drugs that might affect the vision.
Intraocular or capsular inflammation.
Intraocular infection.
Glaucoma.
Retinal disorders (e.g. macular degeneration, diabetic retinopathy, history of retinal detachment, cystoid macular oedema, macular hole).
Corneal disorders comprising visual acuity (e.g. corneal endothelial diseases, corneal dystrophy, history of keratoconus, keratoectasia).
Chorioidal haemorrhage or any other intraocular haemorrhage.
Intraocular hypertension.
Cataract or corneal anomalies (e.g. zonular laxity, zonulodysplasia) which may affect the postoperative centring or the postoperative tilt of the intraocular lens.
Posterior capsule rupture or wide capsulorhexis (compromised intraocular lens stability).
Anterior chamber collapse.
Narrow anterior chamber.
Microphthalmia.
Anisometropia.
Through preoperative evaluation and careful clinical analysis have to be performed by the surgeon to rigorously assess the benefit/risk ratio before an intraocular lens is implanted in such patients.

LUCIS[®] CLARÉ[®] CRISTAL[®]

INTRAOCULAR LENS LENTILLE INTRAOCULAIRE INTRAOKULARLINSE LENTE INTRAOCULAR LENTE INTRAOCULARE ІНТРАОКУЛЯРНА ЛІНЗА

25% hydrophilic acrylic
Acrylique hydrophile 25%
25% hydrophiles Acrylat
Acrílico hidrófilo al 25%
Acrílico idrofílico 25%
25% Гідрофільний акрил

Precautions for use
Do not use after the expiration date.
Do not use if the packaging has been damaged or opened. Sterility is only guaranteed if the packaging shows no sign of damage.
Do not use if the intraocular lens is damaged or shows any defect.
Do not reuse.
Various surgical techniques can be used to implant an intraocular lens. It is up to the surgeon to select the most appropriate.
Avoid the use of alcohol (e.g. 70% alcohol), particularly when used in the surgical treatment of retinal detachment, may come in contact with the intraocular lens of the posterior capsule of the crystalline lens is not.
In patients with atypical eyes (including unusually large eyes as in common in high myopia), the surgeon may need to make appropriate additional preoperative measurements (e.g. of the posterior cornea).
Optification or clarification that may lead to intraocular lens exchange have been observed with various brands of various manufacturers of hydrophilic acrylic intraocular lenses following procedures using intracapsular or vitreous aspiration of air or gas of SF₆, C₂F₆ or C₃F₈ type.
CRISTALENS Industrie recommends that surgeons avoid using hydrophilic acrylic intraocular lenses in patients with any retinal or corneal pathology that may require either at the time of the cataract procedure or in the future, surgery using any type of gas, vitreous procedure such as vitrectomy, retinotomy, especially in the event of retinal detachment, endothelial keratoplasty (DSEK, DSAEK, DMEK, ...) (non-exhaustive list).
1. Contact CRISTALENS Industrie without delay (Email: marketing@crislen.fr).
2. Return the device with all the information necessary for traceability, as specified at the end of the text.

Instructions for use
The device should be kept at a temperature of about 22°C for at least three hours before implantation. 25% hydrophilic acrylic intraocular lens stored at a temperature below 10°C may any truly become slightly opaque as a result of thermal shock at the moment of implantation (although this will disappear within 30 minutes).
1. Check labelling, especially model, power and expiration date.
Check the packaging condition. The packaging is damaged, use another intraocular lens.
2. Open the external blister and take out the internal blister using aseptic technique to maintain sterile conditions.
Open the internal blister on the sterile area and carefully remove the intraocular lens. If the intraocular lens is not present in the holder, use another intraocular lens.
3. Meticulously examine the intraocular lens. If any damage or anomaly is noticed, use another intraocular lens.
The intraocular lens hydrated state requires careful handling as well as sterile instruments and intraocular solutions suitable for use. Immediate use is necessary upon opening the packaging in order to avoid dehydration of the intraocular lens.
The holder is designed so that fingers can be used to take the implant out of it. Once removed from the holder, the intraocular lens may be handled with an injector or forceps for implantation.
After implantation, ensure that the intraocular lens is inside the capsular bag, and in the right position by checking the reflexes.
At the end of the procedure, dispose of the packaging and holder.

Possible complications and adverse events
As with any surgical procedure, there are a risk, possible complications and adverse events following cataract surgery and/or implantation of an intraocular lens inside (non-exhaustive list):
Dislocation or luxation of the intraocular lens;
Ocular or intraocular inflammation;
Intraocular infection;
Corneal edema;
Macular oedema;
Intraocular pressure rise;
Retinal detachment;
Glaucoma;
Pupillary block;
Capsular block;
Capsular rupture;
Wound leak;
Hypotony;
Dysphotopsia;
Precipitates on the intraocular lens surface;
Reflective error.
Damage to the intraocular lens (scratch, crack, break of the optic, scratch, crack, deformation, break of the haptics).
Secondary surgical interventions including, but not limited to: intraocular lens repositioning, intraocular lens replacement, vitreous aspiration or induction for pupil dilation, wound leak repair, retinal detachment procedure.

Interruptions
The intraocular lens passed the ISO 10975 laser exposure test according to the ISO 11079-5 standard.
The intraocular lens contains no metal.
No direct interactions with drugs are known. However, current or previous treatments with alpha-1-adrenergic receptors antagonists may increase the risk of complications associated with cataract surgery (Intraoperative Floppy Iris Syndrome (IFIS)).
No adverse events related to interactions with other investigational or specific treatment devices are known, other than those mentioned in these instructions for use.

Patient information
Before surgery, the patient should be given appropriate information about intraocular lenses and about the contraindications and possible complications and adverse events to this type of surgery and this type of implant.
Following the operation, the surgeon must warn the patient that in the event of decreased performance of the operated eye, the patient must consult an ophthalmologist without delay.

Calculation of intraocular lens power
The intraocular lens power must be calculated preoperatively for each patient based on keratometric and biometric data, using formulas available from the literature, the surgeon's experience, and estimated A-constants indicated on the box. These estimates serve as a starting point for the power calculation. They will have to be optimized by each surgeon based on clinical experience, surgical techniques, measuring equipment and achieved postoperative results.
The dioptric power, expressed in diopters, of each intraocular lens is displayed on the labelling with the symbol "D".

Expiration date
The expiration date indicated on the box determines the limit for use. It is indicated in the format YYYY-MM-DD.
The intraocular lens must not be implanted after the expiration date indicated on the packaging.

Storage and transport conditions
Protect from sunlight, moisture (humidity < 65%) and frost. Avoid shocks and crushing on the box.
Specific storage and transport conditions; see symbols in these instructions for use and on the box.

Disposal of the device
The device (holder, packaging and intraocular lens) must be disposed of in adapted special collector in accordance with the rules and recommendations in force at each health care facility.
The intraocular lens constitutes an infectious clinical waste.
The box and instructions for use can be recycled.

Disclaimer
CRISTALENS Industrie cannot be held responsible for any damage incurred by a patient resulting from:
Surgical technique or implantation method used by the surgeon.
An integrity fault, a deterioration of the intraocular lens, noticed during the surgery but not leading to the removal of the damaged intraocular lens.
Choice or prescription of the intraocular lens.
Reuse of the device, reoperation, reusing use of a damaged or already opened device (compromised sterility), or use after its expiration date. The identified risks are: damage to the medical device, contamination, infection, endophthalmitis, inflammation, injury, illness, organ or eye loss.
Non-compliance of these instructions for use.

Guarantee
CRISTALENS Industrie guarantees its intraocular lenses against any manufacturing defect.

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