Cardioprotective Haemodialysis

FX classix

High-Flux Dialysis for Improved Survival

Cardioprotective Haemodialysis **SPOT**

**FRESENIUS
MEDICAL CARE**

Fresenius Medical Care Deutschland GmbH
61357 Bad Homburg • Else-Kröner-Str. 1 • Germany

i.v. Nikolay Kapustin

i.v. Oleksandr Zaika

Protect your Patient

Cardioprotective Haemodialysis

The reduction of risk factors for cardiovascular diseases (CVD) is core to the development of dialysis systems and products at Fresenius Medical Care. Outstanding cardioprotection must be reflected in all levels of product development and application.

Wide-ranging cardioprotection

There have been tremendous improvements in the quality and efficacy of haemodialysis (HD) therapy in recent years. Despite this, cardiovascular diseases (CVD) remain the leading cause of death for patients with end-stage renal disease (ESRD). Moreover,

SFP

Services

Over 30 years' experience in dialysis at your service.

- Project planning and consulting
- Training and Education
- Technical Services
- Water Quality Service (WQS)
- Medical Information Services


Products

State-of-the-art technologies enable advanced cardioprotective therapies.

- CorDiax Product line:
 - 5008 CorDiax and 5008S CorDiax
 - FX CorDiax dialysers
 - BCM-Body Composition Monitor
- Classix Product line:
 - 4008S classix
 - FX classix dialysers
- Therapy Data Management System (TDMS)
- Online Purification Cascade® (OPC)

overall and cardiovascular mortality is markedly greater in ESRD patients than in the general population. This is why we put Cardioprotective Haemodialysis on the SPOT. A comprehensive approach that includes services, products and therapies is needed to achieve

the best therapeutic performance – meaning improved clinical outcomes and better quality of life, enhanced control of therapy costs, and simpler, safer handling.

Outcomes

Achieving better outcomes with cardioprotective therapies.

- Reduced mortality risk
- Fewer cardiovascular complications
- Optimised use of resources

Therapies

Cardioprotective therapies designed by the world market leader in haemodialysis.

- High-Flux dialysis
- ONLINE HDF
- Advanced Fluid Management

Protect your Patient

Cardioprotection at the heart of long-term haemodialysis

Chronic kidney disease (CKD), as well as the effects of dialysis itself, can lead to cardiovascular diseases (CVD) such as atherosclerosis and left ventricular hypertrophy (LVH), the largest causes of death in haemodialysis patients.¹

Fresenius Medical Care's mission is to enable nephrologists to provide the best possible therapy for their long-term haemodialysis patients in order to minimise the risk of CVD.

In addition to the efficient removal of uraemic toxins, protecting patients through a high level of membrane biocompatibility and endotoxin retention is crucial in Cardioprotective Haemodialysis.

Therefore, Fresenius Medical Care has developed a new class of dialyser, which opens the door to cardioprotective renal replacement therapy – the FX classix:

FX classix – highest level of biocompatibility


- INLINE steam sterilisation enables the production of sterile and pyrogen-free dialysers and ensures high biocompatibility.²

FX classix – maximum endotoxin retention

- The Helixone® membrane has a high endotoxin retention capacity, which minimises the risk of inflammation.³

FX classix – cost saving potential

- FX classix dialysers provide an additional cost saving potential thanks to the lower rinsing volumes enabled by INLINE steam sterilisation as well as the lower weight of the dialysers, which could result in lower waste management costs.

References

1. de Jager D. et al., JAMA (2009); 302: 1782–1789.
2. Müller T. F. et al., Nephron (1998); 78: 139–142.
3. Weber V. et al., Blood Purif (2003); 21: 365.


FX classix – high performance

- Performing High-Flux dialysis has advantages over Low-Flux dialysis: thanks to the larger pores on the inner surface of the innovative Helixone® membrane, High-Flux dialysers also remove middle molecules such as β_2 -microglobulin while preventing the loss of essential blood components such as albumin. In addition, the permeability to water is much higher than in Low-Flux dialysers.

These benefits reduce the risk of CVD and help to improve the long-term outcomes of your patients.

SPOT on:

- Highest biocompatibility due to INLINE steam sterilisation.²
- High endotoxin retention of Helixone® membrane.
- Cost saving potential due to lower rinsing volumes.


Protect your Patient

Clinical benefits of High-Flux dialysers

Improved survival

The Membrane Permeability Outcome Study (MPO) revealed superior survival rates in high-risk patients when treated with High-Flux membranes compared to Low-Flux membranes. For patients with hypoalbuminaemia (≤ 4 g/dL of serum albumin) or diabetes mellitus, a reduction in the relative risk of death of up to 37% was observed.¹

Kaplan-Meier survival curves for the population of patients with serum albumin levels ≤ 4.0 g/dL (log-rank test $p=0.032$).¹

(Graph adapted from original publication)


Up to 86% of dialysis patients worldwide have a serum albumin level ≤ 4 g/dL, underlining the relevance of these risk factors in dialysis.²

During the first 4 years of the MPO study, one in eleven events of death was prevented when hypoalbuminaemic patients were treated with High-Flux dialysers instead of Low-Flux dialysers.

Guidelines recommend High-Flux dialysers

As a consequence of the results of the MPO study, High-Flux membranes are now recommended by the European Renal Best Practice Advisory Board for all haemodialysis patients:

“Guideline 2.1: Synthetic High-Flux membranes should be used to delay long-term complications of haemodialysis therapy ... even in low-risk patients...”³


Improved anaemia management

In patients with ESRD, it is often necessary to administer EPO to treat anaemia. In addition to this, inflammation often contributes to EPO hypo-responsiveness.⁴ It was shown that High-Flux membranes improved control of anaemia while allowing a progressive reduction in the exogenous EPO dose by 25 to 45%.⁵

Hence, High-Flux membranes offer the potential to reduce EPO costs.

SPOT on:

- Improved patient survival with High-Flux dialysis.¹
- Relative risk reduction of 37 % for hypoalbuminaemic patients with High-Flux dialysis.¹
- Improved anaemia control through High-Flux dialysis.⁵

Recovery of haemoglobin (Hb) levels was significantly better after 6 months for patients treated with High-Flux vs Low-Flux membranes. Further, in this patient group the mean EPO dose was significantly lower.⁵
(Graph adapted from original publication)

References

1. Locatelli F. et al., Journal of American Society of Nephrology (2009); 20: 645–654.
2. The DOPPS report 2004; http://www.dopps.org/pdf/dopps_report_2004.pdf.
3. Tattersall J., Nephrol Dial Transplant (2010); 25: 1230–1232.
4. Gunnell J. et al., Am J Kidney Dis (1999); 33(1): 63–72.
5. Ayli D. et al., J Nephrol (2004); 17: 701–706.

Protect your Patient

Proven benefits of the FX-class® design

Fresenius Polysulfone® has long been the »gold standard« in dialysis membranes. For over 30 years, Fresenius Polysulfone® has stood for outstanding safety and performance. Derived from established Fresenius Polysulfone® technology, the Helixone® membrane is at the core of FX-class® dialysers.

The new FX classix dialysers are part of the FX-class® series. More than 177 million treatments have been performed with FX-class® dialysers proving the record of success of the Helixone® membrane.


The unique design of the FX-class® dialyser is based on refined and optimised performance and handling. Several state-of-the-art technologies have been combined to offer distinctive benefits:

Helixone® membrane – optimised performance

- Optimised membrane permeability enables efficient removal of low molecular weight substances and middle molecules
- Minimal loss of essential blood components
- Produced with Nano Controlled Spinning (NCS™) technology

Optimised haemodynamics

- Homogenous blood flow in the dialyser header through lateral blood-inlet port
- Fewer stagnation zones in the header region
- Risk of bloodline kinking is diminished

Optimised dialysate flow for higher clearances


- 3-dimensional microwave structure of the fibres, together with a higher packing density, ensures a homogenous distribution of dialysate over the entire cross-section of the dialyser
- Radial flow of the dialysate around each fibre within the bundle

Kind to the environment

- Usage of ecologically-friendly plastics
- Lower carbon-footprint as a result of fewer materials, less packaging and less fuel for transport

SPOT on:

- Proven and trusted Helixone® membrane.
- Optimised haemodynamics.
- Optimised dialysate flow.
- Environmental-friendliness.

Protect your Patient

Purity ensured – with steam


INLINE steam sterilisation

Product safety means patient safety. In the manufacturing process of our dialysers, we comply strictly with highest quality standards. Thus, all the FX classix dialysers pass through the unique INLINE steam sterilisation process specifically developed by Fresenius Medical Care.

The blood and dialysate compartments of the dialysers are rinsed with hot steam > 121°C for 15 minutes. Following this, all dialysers undergo a fibre leakage test to ensure the integrity of every single fibre.

No chemical residuals


INLINE steam sterilisation reduces potential hazards from residuals. The basic principle of this method is extensive rinsing with hot steam – without the need for chemicals or gamma sterilisation processes. Gamma irradiation may induce the degradation and alteration of the material chemistry and generate cytotoxic substances.¹ INLINE steam sterilisation therefore leads to highly purified dialysers free from chemical, cytotoxic and carcinogenic residuals and with excellent haemocompatibility.


INLINE steam sterilisation process and integrity test

SPOT on:

- Rinsing with hot steam leads to highly purified dialysers.

INLINE steam sterilisation process.


Integrity test: air pressure is applied to the fibre bundle from one side while the other side contains sterile water. If any leakages were present in the membrane, air would pass through the membrane and create bubbles.

References

1. Shintani H. et al., Journal of Analytical Toxicology (1989); 13: 354–357.

Protect your Patient

Greater protection through active prevention

Superior endotoxin retention


Endotoxins are large molecules from the outer membrane wall of gram-negative bacteria. They are able to enter dialysis fluid, and thus the bloodstream, via the microbial contamination of water or fluid conduits. Once in the patient's blood, endotoxins can induce inflammatory responses and – in the longer term – complications such as amyloidosis or accelerated atherosclerosis.

Membranes, such as Helixone[®], which have a high endotoxin retention capacity, protect the patient from inflammation, particularly when ultrapure dialysate is not available.¹ Therefore, it is crucial to prevent endotoxins entering the bloodstream by adopting the following hygiene regime:

- Use of dialysis membranes with high endotoxin retention capacities, such as Helixone[®], to protect the patient from inflammation.
- Use of dialysis fluid filters to create ultrapure dialysis fluid free from residual endotoxins.
- Overall hygiene of water supply system.

SPOT on:

- High endotoxin retention of Helixone[®] membrane.
- Improved patient protection through ultrapure dialysis fluid.

Endotoxin adsorption per cm² membrane surface area after 120min in-vitro dialysis with contaminated dialysate (endotoxin from bacterial culture filtrates; initial concentration 50EU/mL).¹

(Graph adapted from original publication)

References

1. Weber V. et al., Blood Purif (2003); 21: 365.

The new FX classix

Optimised use of resources


Lower weight – reduced costs

The reduced weight of the FX classix dialysers – due to less packaging and decreased use of processed materials – allows cleaner, more cost-effective waste management and thus conserves valuable resources.

Less rinsing – reduced costs

Since FX classix dialysers are INLINE steam sterilised, rinsing volumes of only 500 mL are needed per treatment. Consequently, it is possible to apply the dialysers quickly with decreased preparation time.


Moreover, the reduced rinsing volume represents an average cost saving of 50 % for the priming fluid.

The new FX classix

Performance data

FX classix High-Flux dialysers	Molecular weight (Dalton)	FX 50 _{classix}	FX 60 _{classix}	FX 80 _{classix}	FX 100 _{classix}
Clearance ($Q_b = 300$ mL/min)					
Cytochrome c	12,230	55	74	89	100
Inulin	5,200	72	95	113	122
Vitamin B ₁₂	1,355	137	162	185	201
Phosphate	132	204	225	244	253
Creatinine	113	224	243	259	264
Urea	60	253	266	279	280
Clearance ($Q_b = 400$ mL/min)					
Cytochrome c	12,230	-	76	92	105
Inulin	5,200	-	99	119	129
Vitamin B ₁₂	1,355	-	175	202	222
Phosphate	132	-	252	279	291
Creatinine	113	-	277	300	309
Urea	60	-	312	334	336
Ultrafiltration coeff. (mL/h x mmHg)		27	38	53	68
Sieving coefficients					
Albumin	66,500			< 0.001	
Myoglobin	17,053			0.1	
β_2 -microglobulin	11,731			0.7	
Inulin	5,200			1	
<i>In vitro performance: $Q_D = 500$ mL/min, $Q_F = 0$ mL/min, $T = 37^\circ\text{C}$ (EN 1283, ISO 8637). Ultrafiltration coefficients: human blood, Hct 32%, protein content 6%.</i>					
Membrane material				Helixone®	
Sterilisation method				INLINE steam	
Housing material				Polypropylene	
Potting compound				Polyurethane	
Units per box				24	
Effective surface (m ²)		1.0	1.4	1.8	2.2
K ₀ A Urea		866	1,068	1,394	1,429
Priming volume (mL)		53	74	95	116
Article number		F00002385	F00002386	F00002387	F00002388

Protect your Patient

High-Flux dialysis – improved survival. Better outcomes

Almost one in two patients with ESRD dies as a result of cardiovascular disease. That is why Cardioprotective Haemodialysis is a core principle at Fresenius Medical Care, as we work and learn to solve the challenges of modern dialysis. Each step we take is focused on minimising cardiovascular risks and extending patients' lives. In recent years, several studies have demonstrated that patients show improved long-term survival when treated with High-Flux dialysers. Hence, the new FX classix dialysers are a fundamental component of our SPOT programme and help you to protect your patient – day by day.

State-of-the-art technologies enable advanced cardioprotective therapies.

✓ Efficient removal of uraemic toxins

✓ Higher clearances by design

✓ Optimised haemodynamics

✓ INLINE steam sterilisation

✓ Optimised use of resources

SPOT on:


- High-Flux dialysis and design benefits.

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Oleksandr Zaika
i.v. Oleksandr Zaika

Cardioprotective Haemodialysis

Despite significant improvements in the quality and efficacy of haemodialysis therapy in recent years, cardiovascular disease (CVD) remains the leading cause of death for dialysis patients. Today, almost every other dialysis patient dies from cardiovascular complications.

Fresenius Medical Care is supporting nephrologists world-wide in reducing their patients' risks for cardiovascular morbidity and mortality.

Innovative membranes like Fresenius Polysulfone® and Helixone®, modern monitoring devices like the Blood Volume Monitor, the Blood Temperature Monitor and Online Clearance Monitoring (OCM®), ultra-pure dialysis fluid prepared with DIASAFE® *plus* and modern ONLINE haemodiafiltration systems support the reduction of CVD risk factors.

Moreover, one of our major goals in coming years is the development and implementation of innovative new therapies and products that further improve the cardiovascular prognosis of dialysis patients.

The Dialysis Fluid Filter DIASAFE®*plus*

The quality and purity of the dialysis fluid are of major concern in modern-day renal replacement therapies, as large volumes of dialysis fluid come into contact with the patient's bloodstream during each treatment.

Bacterial endotoxins present in contaminated dialysis fluid may elicit undesirable acute reactions and influence the long-term outcome of patients on chronic haemodialysis.

Although water used for the production of dialysis fluid is treated by a series of purification steps, it still may not meet the stringent requirements on bacterial contamination levels laid down by various regulatory bodies.

By the application of special filters that are highly efficient in retaining bacterial contaminations, the required purity grades of dialysis fluid can be achieved easily.

The DIASAFE®*plus* filter, located at the end of the water treatment chain, ensures the safe production of ultrapure dialysis fluid. This is attributed to the excellent endotoxin-retention capabilities of its Fresenius Polysulfone® fibres.

DIASAFE®*plus* is an integral part of contemporary dialysis machines. Only three handling steps are necessary to install or exchange DIASAFE®*plus* (Fig. 1):

- Open the locks of the filter holder
- Slide DIASAFE®*plus* filter into the guide grooves
- Close the locks – DIASAFE®*plus* is ready to use


Fig. 1: DIASAFE®*plus* can be connected with 3 rapid handling steps only. The DIAFIX™ lock system ensures a safe and hygienic connection.

The use of DIASAFE®*plus* is a key step towards Good Dialysis Practice.

Dialysis Fluid Purity

Dialysis fluids may contain microbial impurities such as endotoxins derived from bacterial fragments. Endotoxins are known to cause acute adverse reactions and promote long-term complications in haemodialysis patients ^(1, 2).

The toxic properties of endotoxins can be ascribed mainly to their lipid A component, which is not exposed by intact bacteria, but released only during growth or lysis of gram-negative bacteria ^(3, 4).

Endotoxin fragments may have molecular weights well below 2000 Da. These fragments are small enough to pass across both, low- and high-flux haemodialysis membranes into the patient's blood-stream (Fig. 2).

With respect to endotoxin permeability, significant differences exist between the various types of dialysis membranes, thereby offering variable degrees of safety during haemodialysis ⁽³⁾.

In order to avoid endotoxin-related complications during routine haemodialysis, the European Best Practice Guidelines for Haemodialysis (EBPG) ⁽⁵⁾ advise the usage of water having a purity level in compliance with the recommendations of the European Pharmacopoeia. However, the usage of ultrapure water for conventional high-flux dialysis is strongly recommended by the EBPG (Table 1).

Ultrapure water or dialysis fluid can easily be achieved through the application of special dialysis fluid filters – such as DIASAFE[®] plus.

	Pure water	Ultrapure water
Microbial contaminations (CFU/mL)	≤ 100	< 0.1
Bacterial endotoxins (IU/mL)	< 0.25	< 0.03

Table 1: The different purity levels of pure and ultrapure water according to the EBPG.


Fig. 2: Fragments of bacterial endotoxins enter the patient's blood-stream and activate leucocytes, thereby leading to acute and chronic complications in haemodialysis patients.

Clinical Advantages of Using Ultrapure Dialysis Fluid

Endotoxins can activate immune-competent cells in a number of ways, thereby contributing to chronic inflammation that is present in all haemodialysis patients ⁽⁶⁾ (Fig. 3). Recent evidence demonstrates that chronic inflammation is a major risk-factor for progressive atherosclerotic cardiovascular disease (CVD) ⁽⁷⁾.

Besides the application of haemodialysis membranes with a high biocompatibility, the usage of ultrapure dialysis fluid, in particular, has been shown to reduce markers of chronic inflammation in haemodialysis patients ⁽⁸⁾. Therefore, it is suggestive that ultrapure dialysis fluid has a beneficial effect on inflammatory diseases such as atherosclerotic CVD ⁽⁹⁾.

Moreover, oxidative stress – a situation, in which the normal balance between production of reactive oxygen species (ROS) and antioxidant activity is tilted in favour of ROS – is increased by several treatment-related stimuli, including bacterial endotoxins derived from the dialysis fluid ^(10, 11).

As oxidative stress is associated with the progression of malnutrition, anaemia and inflammatory diseases such as atherosclerosis, it appears desirable to reduce dialysis-induced oxidative mechanisms, e. g. through the usage of biocompatible membranes and ultrapure dialysis fluid ^(10, 11).

The importance of ultrapure dialysis fluid in routine haemodialysis treatments is emphasized by the finding that endotoxins act in synergy with advanced glycation end-products (AGE), which enhance inflammation and oxidative stress ⁽¹²⁾. Furthermore, the use of ultrapure dialysis fluid has been shown to reduce the plasma levels of the AGE compound pentosidine ⁽¹³⁾.

Finally, ultrapure dialysis fluid has also been shown to improve iron utilization and the response to erythropoietin; thus, ultrapure dialysis fluid could be beneficial for anaemia treatment allowing for a reduction in erythropoietin dosage, while maintaining optimal haemoglobin levels ^(14, 15).


Fig. 3: Endotoxins (LPS) stimulate the release of pro-inflammatory cytokines, reactive oxygen species and lipid-mediators from immune-competent cells.

DIASAFE® plus in the ONLINEplus™ System

Ultrapure dialysis fluid prepared with the DIASAFE® plus dialysis fluid filter, together with haemodialysers containing endotoxin-retaining membranes (Fresenius Polysulfone® or Helixone®) are the main building blocks for a high-quality haemodialysis treatment.

The ONLINEplus™ system takes the quality standards of convective treatment modalities as haemodiafiltration/haemofiltration (HDF/HF) one step further: using two DIASAFE® plus dialysis fluid filters in series, an

extremely high microbiological safety is achieved by double filtration of the substitution fluids used in ONLINE HDF/HF therapies ⁽¹⁶⁾.

Besides improving hygiene and safety of convective therapy modalities, the ONLINEplus™ option also offers additional treatment features and adds to ease of handling.


Fig. 4: Schematic flow chart of ONLINE haemodiafiltration with the ONLINEplus™ system

References

1. Dasgupta MK: Biofilms and infection in dialysis patients. *Seminars in Dial* 15: 338-346, 2002.
2. Brunet P and Berland Y: Water quality and complications of haemodialysis. *Nephrol Dial Transplant* 15: 578-580, 2000.
3. Lonnemann G: Chronic inflammation in hemodialysis: the role of contaminated dialysate. *Blood Purif* 18: 214-223, 2000.
4. Golenbock DT et al.: Lipid A-like molecules that antagonise the effects of endotoxins on human monocytes. *J Biol Chem* 266: 19490-19498, 1991.
5. European Best Practice Guidelines for Haemodialysis (Part 1), Section IV – Dialysis fluid purity. *Nephrol Dial Transplant* 17 (Suppl. 7): 45-62, 2002.
6. Ward DM: Hemodialysis water: an update on safety issues, monitoring and adverse clinical effects. *ASAIO J* 50 (6): XIII-XIX, 2004.
7. Yao Q et al.: Inflammation as a cause of malnutrition, atherosclerotic cardiovascular disease and poor outcome in hemodialysis patients. *Hemodial Int* 8: 118-129, 2004.
8. Schiff H et al.: Effects of ultrapure dialysis fluid on nutritional status and inflammatory parameters. *Nephrol Dial Transplant* 16: 1863-1869, 2001.
9. Lonnemann G: When good water goes bad: how it happens, clinical consequences and possible solutions. *Blood Purif* 22: 124-129, 2004.
10. Locatelli F et al.: Oxidative stress in end-stage renal disease: an emerging threat to patient outcome (consensus paper). *Nephrol Dial Transplant* 18: 1272-1280, 2003.
11. Ward RA and McLeish KR: Oxidant stress in hemodialysis patients: what are the determining factors? *Artif Organs* 27: 230-236, 2003.
12. Reznikov LL et al.: Effect of advanced glycation end products on endotoxin-induced TNF- α , IL-1 and IL-8 in human peripheral blood mononuclear cells. *Clin Nephrol* 61: 324-336, 2004.
13. Izuwara Y et al.: Ultrapure dialysate decreases plasma pentosidine, a marker of carbonyl stress. *Am J Kidney Dis* 43: 1024-1029, 2004.
14. Sitter T et al.: Dialysate related cytokine induction and response to recombinant human erythropoietin in hemodialysis patients. *Nephrol Dial Transplant* 15: 1207-1211, 2000.
15. Hsu PY et al.: Ultrapure dialysate improves iron utilisation and erythropoietin response in chronic haemodialysis patients – a prospective cross-over study. *J Nephrol* 17: 693-700, 2004.
16. Weber C et al.: Novel online infusate-assisted dialysis system performs microbiologically safely. *Artif Organs* 24: 323-328, 2000.

Technical Data

Membrane material	Fresenius Polysulfone®
Effective Surface (m ²)	2.2
Weight (g)	170
Housing material	Polypropylene
Potting material	Polyurethane
Sealings	Silicone
Connection to machine	DIAFIX™ Lock System
Filtration rate	5 mL/min mm HG (3.75 L/min bar; max. 2 bar)
Operating time	Standard HD: max. 12 weeks ONLINE HF/HDF, ONLINE priming/rinsing: max. 12 weeks or 100 treatments
Disinfection	Puristeril® 340 or Puristeril® <i>plus</i> (peracetic acid) Diasteril® (hydroxyacetic acid) or Citrosteril® (citric acid) Sporotal® 100 (sodium hypochlorite) max. 11 times
Article number	500 820 1

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Hygiene

Citrosteril

Heat Disinfectant for Haemodialysis Machines
with Recirculation

Fresenius Medical Care Deutschland GmbH
61352 Bad Homburg • Else-Kröner-Str. 1 • Germany

Woz
i.v. Nikolay Kapustin

Aleksandr Zaika
i.v. Aleksandr Zaika

Citrosteril

Heat Disinfectant for Haemodialysis Machines with Recirculation

Disinfectant for chemo-thermal disinfection of haemodialysis machines¹.

- pH value 1.7 to 2.0
- excellent removal of limescale
- disinfection and decalcification in one process
- active ingredients composed of natural substances
- biodegradable
- odourless
- free from colouring additives

Action

The synergistic effect of its components makes Citrosteril a potent disinfectant solution.

Citrosteril at 84°C has a broad spectrum of micro-biocidal activity and works bactericidal and virucidal² including HBV/HCV/HIV.

Specification

100 g Citrosteril contains:

21 g citric acid 1-hydrate; lactic acid, malic acid


Ordering information


Unit	Language combination	Art. No.
1 x 5L (single canister)	multilingual	F00005157
2 x 5L (carton)	multilingual	F00005158
6 x 2L (carton)	multilingual	508 536 1

Literature:

1. Solbach W, Universität zu Lübeck; Verification of the sporicidal efficacy of the product Citrosteril at 85°C in the Fresenius 5008 dialysis machine, 20.12.2002
2. Labor Dr. Merk & Kollegen, Ochsenhausen: Antiviral efficacy of Citrosteril against bovine Parvovirus, 9.9.2005

Further information is available on request.

Haemodialysis

Diasteril®

For Heat Disinfection of Haemodialysis Machines
with Recirculation

**FRESENIUS
MEDICAL CARE**


Fresenius Medical Care Deutschland GmbH
61352 Bad Homburg • Else-Kröner-Str. 1 • Germany

Fresenius Medical Care

Waz
I.V. Nikolay Kapustin

Oleksandr Zaika
I.V. Oleksandr Zaika

Diasteril®

For Heat Disinfection of Haemodialysis Machines with Recirculation

The alternative heat disinfectant is based on hydroxyacetic acid.


In nature, this physiological acid is present in grapes, sugar cane and in fat of sheep's wool.

- Delimiting and cleaning effect.
- Broad spectrum of action, bactericidal, virucidal, inactivates Hepatitis B/C and HIV
- Odourless and non-volatile
- Biological degradable
- Lowest chemical oxygen requirement of the biological degradable waste water
- No crystallizing
- Economical and environmentally friendly procedure
- Optimised consumption of disinfectant
- Test for disinfectant residues carried out by using pH indicator strips (pH-Fix 3.6–6.1, Art. No. 628 816 1)

Specifications

	Volume	Languages	Art. No.
Diasteril®	6 litres	GB//E//P//F//GR	508 565 1
Diasteril®	6 litres	D//NL//F//SK//CZ	508 564 1
Diasteril®	6 litres	RU//SK//RO//PL//CZ//GB	508 566 1

Expert's opinion

Prof. Dr. med. W. Solbach, Luebeck;

- 1) Sporocidal effect (*B. subtilis*) machine at 84°C, July 15, 1997
- 2) Sporocidal effect (*B. subtilis*) Port 1 at 84°C, February 18, 1998
- 3) Disinfectant residues, May 12, 1998

Further evaluations upon request!

Fresenius Medical Care

Fresenius Medical Care Deutschland GmbH · 61346 Bad Homburg v. d. H. · Germany · Phone: +49 (0) 6172-609-0 · Fax: +49 (0) 6172-609-2191
 Head office: Else-Kröner-Straße 1 · 61352 Bad Homburg v. d. H.
www.fmc-ag.com

Haemodialysis

Puristeril® 340

Cold Disinfectant for Haemodialysis Machines

**FRESENIUS
MEDICAL CARE**

Fresenius Medical Care Deutschland GmbH
61352 Bad Homburg • Else-Kröner-Str. 1 • Germany

Fresenius Medical Care

i.v. Nikolay Kapustin


[Signature]
i.v. Oleksandr Zaitsev

Puristeril® 340

Disinfecting Agent for Haemodialysis Machines Based on Peracetic Acid

Superior efficacy

- Puristeril® 340 shows the superior efficacy of a peracetic acid-based disinfectant.
- Peracetic acid is widely used for disinfection due to its exceptionally broad spectrum of microbiocidal activity at low concentrations and short exposure times.
- Puristeril® 340 decomposes in a non-toxic way. The following degradation reactions take place:

- After use Puristeril® 340 is easily removable by rinsing with water.
- Due to the low pH value, the necessary decalcification of haemodialysis machines is easily achieved.
- Puristeril® 340 is designed for cold disinfection. In principle it can be used for all haemodialysis systems like haemodialysis machines, water treatment devices and circuit pipes.

Specification

Puristeril® 340 contains peracetic acid and hydrogen peroxide.

Disinfection

Puristeril® 340 is bactericidal, fungicidal, sporicidal, virucidal (incl. HBV/HCV/HIV).

Use Puristeril® 340 in accordance with the instructions provided by the manufacturer of the machine. Puristeril® 340 can be used in all Fresenius Medical Care 2008, 4008 and 5008 haemodialysis machines as well as the GENIUS® system.

For disinfection of water treatment systems, proceed according to the manufacturer's instruction.

Testing for residual disinfectant

For safety reasons, a test to show the absence of residual disinfectant residues must be performed after the completion of the disinfection procedure. The absence of Puristeril® 340 is detectable by potassium iodide starch paper (art.no. 508 521 3).

Stability and storage

Properly stored, the disinfectant remains fully effective for 18 months after production. Keep container sealed at all times and store in an upright position. If possible, store in well-aired rooms at 5 to 25°C. Do not expose to direct sunlight.

Order information

Article	Quantity	Art. No.
Puristeril® 340	1 × 5 kg	508 562 1 (multilingual)
Puristeril® 340	1 × 10 kg	508 563 1 (multilingual)
Puristeril® 340 GENIUS®	1 × 3 kg	508 567 1 (multilingual)

Evaluations are available on request.

Fresenius Medical Care

Fresenius Medical Care Deutschland GmbH · 61346 Bad Homburg v. d. H. · Germany · Phone: +49 (0) 6172-609-0 · Fax: +49 (0) 6172-609-2191

Head office: Elsa-Kröner-Straße 1 · 61352 Bad Homburg v. d. H.

www.fmc-ag.com